# evo fill™ duo

Cordless Warm Vertical Compaction & Backfill Obturation Device

# User Manual

# User Manual

### Manufacturer

 DiaDent Group International  
16, Osongsaengmyeong 4-ro, Osong-eup,  
Heungdeok-gu, Cheongju-si,  
Chungcheongbuk-do, 28161, Korea  
Tel : 82-43-266-2315  
Fax : 82-43-262-8658  
www.diadent.co.kr  
diadent@diadent.co.kr

### Distributor

DiaDent Group International Inc.  
11-3871 North Fraser Way, Burnaby  
B.C., Canada, V5J 5G6  
  
Tel : 1-604-451-8851  
Fax: 1-604-451-8865  
www.diadent.com  
diadent@diadent.com

### Distributor

DiaDent Europe B.V.  
Antennestraat 70, 1322AS  
Almere, the Netherlands  
  
Tel : 31-36-549-8607  
Fax: 31-36-536-7317  
www.diadent europe.com  
diadent@diadent europe.com

REF Dec.01, 2015 (Rev.0)


# **eVO fill™ duo**

*Cordless Warm Vertical Compaction &  
Backfill Obturation Device*

---

	Page
<b>ENGLISH</b>	04
<b>FRANÇAIS</b>	36

	Page
Precautions for Use	5
Intended Use	6
Purpose & Benefits	7
Product Information	7
Parts & Components	8
EvoPack Component Names & Functions	10
EvoFill Component Names & Functions	12
Connection Part	15
<hr/>	
EvoPack Operation Method	16
EvoPack Directions for Use	17
EvoFill Operation Method	18
EvoFill Directions for Use	20
EvoFill Duo Duty Cycle	21
EvoFill Duo Charger Function & Charging Method	22
<hr/>	
Checklist & Maintenance	23
Cleaning Sterilization & Maintenance	24
Guidance DiaDent's Declaration - Electromagnetic emissions	30
Troubleshooting & Product Warranty	33
Specifications	34

This manual provides detailed operating and setup instructions and is designed to help users prevent unexpected problems and avoid damages.

 <b>Caution</b>	Noncompliance with the manual can cause body injury or property damage.
--	---

- ▶ This product is a medical device designed for dental professional use only.
- ▶ Completely insert the power plug into the outlet. Unstable connection may cause a fire.
- ▶ To prevent fire or electric shock, do not use loose electrical outlet or damaged power cord.
- ▶ To prevent an electric shock, never touch the power plug with wet hands.
- ▶ Do not plug in multiple devices to one outlet at the same time. An overloaded outlet can cause a fire as it gets hot.
- ▶ Do not place the power cord or electrical plug near heating devices. Molten cable sheaths can cause a fire or an electric shock.
- ▶ Keep the battery away from heat, sparks, and flames to prevent an explosion or a fire.
- ▶ Keep the battery away from a live terminal to prevent an electric shock and possible explosion.
- ▶ This product is for indoor use only. Do not use it outside.
- ▶ Operating conditions  
Temperature: 10~40 °C, Atmospheric pressure: 60~101kPa, Relative humidity: 30~75%
- ▶ Do not use this device in dusty or humid places
- ▶ Check the power voltage before using to prevent any product damage.
- ▶ When unplugging the power cord, hold the plug head and never pull the cable. Broken cables can cause a fire.
- ▶ Gutta Percha Cartridge (EvoFill) and Heating Condenser (EvoPack) get hot during operation. Handle with care to prevent burns.
- ▶ In case of device failure, contact the manufacturer or nearest distributor/dealer for repair. Do not disassemble the device.
- ▶ Keep the device clean and properly stored at all times.
- ▶ Physical impacts can cause damages. Be careful not to drop the product.
- ▶ To prevent a direct or indirect breakdown of the device, leave the charger unplugged in case of thunder and lightning.
- ▶ This product package includes dangerous objects such as Heating Condenser(EvoPack). Keep out of reach of children.
- ▶ Do not modify this product without approval of the manufacturer.
- ▶ This product is not waterproof. Do not soak in water or any liquid. (IPX0)

This product has been tested according to the IEC 60601-1-2 and is designed to provide reasonable protection against harmful interference to general operation and installation of medical devices.

This product generates, uses, and can radiate radiofrequency energy and thus cause interference to other devices nearby. In case of harmful interference, user is recommended to try one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation distance between the interference-causing equipment and receiver.
- Consult a local dealer or the manufacturer.

## Intended Use

Intended for use by professionally qualified licensed dentists, endodontists and clinicians to heat gutta percha and place it into the previously prepared root canals of human teeth, in order to provide a quick and complete obturation of the canal.

- 1) Point of body contact : Direct contact of tooth
- 2) Duration of contact : Temporary use (<60 min)

## Purpose & Benefits

### ► EvoPack

- EvoPack is an endodontic device that uses heat to cut, soften, downpack and compact the root canal filling material. It is an integral part of the warm vertical condensation technique.
- EvoPack effectively and tightly compacts and seals all canals including lateral canals.
- EvoPack features LED light-guided plugger, instant heating, and digital temperature control.

### ► EvoFill

- EvoFill is an endodontic device used to obturate root canals with Gutta Percha.
- Hygienic one-time use Gutta Percha Cartridge (not included in the package) has built-in specially-formulated Gutta Percha that can deliver fast, precise, and direct injection of softened Gutta Percha into the root canal.
- With the innovative electric motor, EvoFill features precisely controlled motorized gutta percha extrusion, fast heating of 360° swivel cartridge and digital control panel.
- Save valuable time with hygienic single-use gutta percha cartridge that provides reliable result every time.

## Product Information

### ■ EvoPack

As an integral part of the vertical condensation technique, Heating Condenser (EvoPack) can transfer heat in a temperature-controlled setting to soften, spread, cut, downpack, and compact Gutta Percha.

### ■ EvoFill

EvoFill uses hygienic one-time use Gutta Percha Cartridge to deliver fast, precise, and direct injection of softened Gutta Percha into the root canal. Even though it is possible to obturate the canal using EvoFill only, it is recommended to use EvoPack first to downpack to the apex and backfill the rest of the space using EvoFill to prevent over-filling or under-filling.

### ■ EvoFill Duo Charger (Charging Station)

Place EvoFill and EvoPack on the charging station.

The LED light on the charging station will display a purple light while the device is charging.

Once charging is complete, it will show a solid blue LED light.

## Parts & Components

▪ EvoPack


▪ Heating Condenser (EvoPack):  
Expired date N/A (Reusable)

▪ EvoFill

▪ Safety Cap

▪ Charger

## Parts & Components

▪ Adapter

▪ Power Cord

▪ Battery

▪ Multi-tool

▪ Scrubbing Brush

Accessories are to be used for DiaDent products only and not to be used with other manufacturer devices.


All accessories have the DiaDent's logo. **DiaDent**<sup>®</sup>

▪ Gutta Percha Cartridge (EvoFill): Expiration: 4 years (Disposable)  
20pcs. in package (ITEM #4005-2010)

(\* Not included in EvoFill Duo package)

## EvoPack Component Names & Functions

**Warning :**  
Hot Surface

Heating Condenser (EvoPack) is intended to heat up as this product is used to cut and melt root canal filling materials. Heating Condenser (EvoPack) temperature can reach up to 250 °C, operate with care to prevent burns.

### 1. Power Button

Press and hold the button for a few seconds to turn on and off the device. The device will turn on or off with an audible alert.

### 2. Menu Button

(1) Short Press (less than 1 second)

- Select the heating temperature.
- Pressing the menu button changes the temperature in the following order **170 °C → 200 °C → 250 °C → 170 °C. (The factory setting is at 170 °C.)**

## EvoPack Component Names & Functions

(2) Press and hold (long press - over 1 second)

- Heating Condenser (EvoPack) set the heating time.
- Each long press changes the heating time in the following order **S (3 Seconds) → L (5 Seconds) → S (3 Seconds).**

※ All setting information is stored even if the device is turned off.

### 3. Operation Button

- Press the operation button to heat up Heating Condenser (EvoPack). A beep will sound to signal. Holding on to the operation button will make the beep continue to sound
- The LED light will come on when the device is being heated for the selected heating time. When the button is released or once the heating time is over, the device will automatically stop heating. The LED light will turn off automatically after 5 seconds.

### 4. Display

The display screen shows the current temperature, remaining battery, heating time, and error status.

### 5. Heating Condenser (EvoPack)

- The very end of the Heating Condenser (EvoPack) heats up to cut a Gutta Percha point (Master Cone) placed inside the root canal. It is also used for melting and compacting Gutta Percha.
- The error indicator will flash on the display screen if the tip is not firmly and properly inserted into the EvoPack hand-piece
- Once the tip is properly inserted, the error message will disappear on the screen.

### 6. LED Light

When the operation button is pressed, the LED light will work as a flashlight within the oral cavity. When the operation button is released or when the operation time is over, the LED light will automatically go off after 5 seconds.


\* IEC62471 test completed: Be careful not to point the LED light directly into eyes.

### 7. Battery Pack


The battery pack can be removed or replaced. Use the Multi-Tool provided in the package to remove and replace the battery.

[Label of battery pack]


- ① Power button
- ② Menu button
- ③ Operation button
- ④ Display
- ⑤ GP Window
- ⑥ Gutta Percha Cartridge (accessory not included, purchase separately.)
- ⑦ Safety Cap
- ⑧ Battery pack

- ① Battery indicator
- ② Error Status
- ③ Heating temperature
- ④ Moving direction of the piston

GP Window (Piston forward)

GP Window (Piston backward)

EvoFill is intended to heat up for root canal obturation purposes.  
The device can heat up to 200°C. Operate with care to prevent burns.

1. Power Button

- Press and hold the button for a few seconds to turn the power on and off. An audible alert will sound to signal.
- When the power is turned on, EvoFill automatically starts to pre-heat to the selected temperature. (The factory setting is at 180°C.)
- Preheat the device at a high temperature for a short time in the beginning of the procedure to quickly soften the Gutta Percha.

2. Menu Button

(1) Short Press

Pressing the menu button changes the temperature in the following order:  
160°C → 180°C → 200°C → 160°C.

- ※ All setting information is stored even if the device is turned off.
- (2) The piston automatically moves backwards on a longer press.

3. Operation Button

- When the operation button is pressed, the heat-softened Gutta Percha is extruded from the Gutta Percha cartridge.
- The operation button works only when preheating process is complete. (The temperature indicator on the display screen will stop flashing once the preheating process is complete).
- ※ Preheating enables quick warming and softening of Gutta Percha when the initial heating temperature is low.

- If preheating starts at 120°C or lower, the temperature indicator on the display screen will flash as it heats up. Once it has reached the selected temperature, the temperature indicator will stop flashing and beep to signal.
- When the device is in use or when the preheating temperature is 120°C or higher, the device is in operation mode, not preheating mode.
- Once the Gutta Percha Cartridge has been replaced, wait at least 1 minute before operating the device again to soften the Gutta Percha in the cartridge.

## EvoFill Component Names & Functions

### 4. Display

The display screen shows the current temperature, remaining battery, moving direction of the piston, and error status.

### 5. GP Window

GP Window displays the amount of Gutta Percha available in the Gutta Percha Cartridge (EvoFill).

### 6. Safety Cap

In order to prevent burns, the safety cap must be used at all time.

### 7. Battery Pack

The battery pack can be removed or replaced by users. Use the Multi-Tool provided in the package to remove or replace the battery.


[Label of battery pack]

## Connection Part

### Connection part

EvoFill: Gutta Percha Cartridge


EvoFill Connecting Part = (A) + (B)  
 (A) Gutta Percha Cartridge  
 (B) Heating Chamber  
 (first 15mm at the front-end of the device)

- The front part of the device is equipped with a heating chamber to heat-soften Gutta Percha. Therefore, the connecting part of A and B may become hot during operation. Always use the EvoFill cover to prevent burns.

- The heating chamber and Gutta Percha Cartridge (EvoFill) may still remain hot after use. Use the Multi-Tool provided in the package to prevent burns. Depending on the temperature setting, the Gutta Percha Cartridge (EvoFill) temperature can reach up to 200 °C.

### EvoPack: Heating Condenser (EvoPack)

※ EvoPack Connecting Part = Heating Condenser (EvoPack)

- Do not touch the Heating Condenser (EvoPack) during operation to prevent burns. (Depending on the temperature setting, the Heating Condenser (EvoPack) temperature can reach up to 250 °C.)

## EvoPack Operation Method

### 1. Insert the heating condenser (EvoPack)

- Insert the Heating Condenser (EvoPack) into the front part of the EvoPack handpiece. Ensure the Pack Tip is completely and properly inserted.

\* To prevent burns, check if the Heating Condenser (EvoPack) has cooled down before replacing it.

### 2. Turn on the power.

- Press and hold the power button to turn on the device. An audible alert will sound as the device turns on.

- Press the menu button to change the temperature setting (170 °C, 200 °C, and 250 °C).  
- Press and hold (long press) on the menu button can change the heating time (S: 3 sec and L: 5 sec).

### 3. Press the operation button.

- Pressing the operation button starts heating the tip and turns on the LED light with an audible alert.  
- The device stops heating when the button is released or when heating is completed.  
The LED light automatically goes off after 5 seconds.

- When the Heating Condenser (EvoPack) is not properly inserted, the error indicator will flash.  
If the user presses the operation button in this case, the device will signal a warning sound.

\* Do not touch the Heating Condenser (EvoPack) while in use to prevent burns.

### 4. Turn off


- Press and hold the power button to turn off the device. An audible alert will sound as the device turns off.

※ When the battery level is low, the battery indicator on the display screen will flash and the device will go off automatically. Make sure the device is fully charged before using it.

## EvoPack Directions for Use

### Directions for Use

#### Downpack Obturation with EvoPack




## EvoFill Operation Method

1. Insert the Gutta Percha Cartridge (EvoFill). Gutta Percha Cartridge not included in the package.

Purchase separately.

- Insert the Gutta Percha Cartridge (EvoFill) into the handpiece and turn the cartridge nut clockwise to securely connect the two parts. Turn counterclockwise to release Gutta Percha Cartridge (EvoFill).
- Put on the EvoFill Safety Cap at the front part of the device.
- Bend the tip part of the Gutta Percha Cartridge (EvoFill) using the Multi-Tool provided.

2. Turn on the device.

- Press and hold the power button to turn on the device. It automatically starts to preheat to the selected temperature.
- If the initial heating temperature is low, preheat Gutta Percha at a higher temperature for a short time to soften it quickly.
- The temperature indicator flashes when the device is preheating and stops flashing with an audible alert when preheating is complete.

3. Select the temperature.

- Choose the desired temperature. Press on the menu button to change the temperature in the following order (160 °C, 180 °C, and 200 °C).
- Press and hold the menu button to move the piston backward.

## EvoFill Operation Method


4. Press the operation button.

- If the user presses the operation button before preheating is complete (when the temperature indicator is still flashing), a warning sound will signal and the device will not operate.
- Once preheating is complete, press the operation button to extrude the softened Gutta Percha. Continue extruding while withdrawing the tip from the canal. Use a hand plugger to compact the obturation material.
- When the Gutta Percha in the Gutta Percha Cartridge is used up, the piston will automatically move back. It takes 25-30 seconds for the piston to fully move back.
- Once the piston is moved back completely and has stopped moving, replace the Gutta Percha Cartridge using the Multi-Tool.
- The Gutta Percha Cartridge and heating chamber may still be hot. Handle with care to prevent burns.
- In order for the Gutta Percha inside the cartridge to soften, wait at least 1 minute before operating the device again.
- When the sensor or wire is disconnected, the error indicator will flash with an audible alert.

\* Caution

Use the Multi-Tool to disassemble the tip.

- Do not use the device on any body parts other than a tooth.
- Do not use the Gutta Percha Cartridge on a tooth for more than 20 seconds. The heated Gutta Percha Cartridge, which can reach up to 200 °C, may cause a burn.
- The Gutta Percha Cartridge (EvoFill) is not reusable. The Gutta Percha Cartridge (EvoFill) must be sterilized with rubbing alcohol before use.
- Before replacing the Gutta Percha Cartridge, make sure the device has cooled down completely.
- Be careful while using the device. The motor may still operate without the Gutta Percha Cartridge connected when the device reaches the desired temperature.

5. Turn off the device.

- Press and hold the power button to turn off the device. An audible alert will sound as the device turns off.

※ When the battery level is low, the battery indicator will flash with an audible alert and the device will automatically turn off. Make sure it is fully charged before using.


## EvoFill Directions for Use

### Obturing (Backfill)


- Fit the master Gutta Percha cone that is dipped in sealer and take a radiograph to ensure that it is fitting well.
- Increase the preset temperature of EvoFill, so the Gutta Percha inside the cartridge is heated.
- Using EvoFill's tip, make contact with the upper part of Gutta Percha Point that is already filled in the root canal, transfer heat for 3 to 5 seconds.
- Extrude the melted Gutta Percha Obturator to the root canal and continue extruding while withdrawing the tip from the canal.
- Use an Endodontic hand plugger to compact obturation material.
- Repeat as needed until canal is filled without voids and take a radiograph.

### ROOT CANAL OBTURATION TECHNICAL SEQUENCE

#### Backfill Obturating with EvoFill

- The device will switch to standby mode after the final use.
- The device in standby mode (LCD display screen on) will automatically switch to sleep mode after 5 minutes.
- If a user presses a button during standby mode, the device will go back to operation mode.
- If sleep mode (LCD display screen off) lasts more than 5 minutes, the power will turn off automatically. If a user presses a button during sleep mode, the device will automatically turn back to operation mode.

### EvoFill Duo Operation Mode Cycle


#### [EvoFill]

- A new Gutta Percha Cartridge can be used for 1 to 70 seconds.
- Once the Gutta Percha inside the cartridge has been all used up, the piston moves backwards. It is recommended that users start using the device again after at least 90 seconds, in the meantime, the cartridge can be replaced, and the device can be preheated.

#### [EvoPack]


- Depending on the mode of operation chosen (S or L), the heating condenser can be used for 3 to 5 seconds.
- After heating time (3~5sec.), users are recommended to use remaining heat for 2 to 3 seconds during the waiting time and press the button to re-start and re-use.

## EvoFill Duo Charger Function & Charging Method

### 1. EvoFill Duo Adaptor


- Connect the power adaptor to the back side of the charger.
- Information on adaptor
  - ▶ Manufacturer: BridgePower Corp. (Model name: BPM030S09N03)
  - ▶ Input: 100/240V $\varnothing$ , 50 ~ 60Hz, 1A
  - ▶ Output: 9V $\equiv$ , 3A

\* Only use the adaptor provided by the manufacturer. Using any other adaptor may cause fire or electric shock.

### 2. EvoFill Duo Charger

- To charge the battery, place EvoFill (on the left side) and EvoPack (on the right side) on the charging station.

### 3. LED Display

- Charging : Purple LED light
- Fully charged : Blue LED light

\* If the LED lights do not come on, check if the devices are properly placed on the charging station.

### 4. Battery can be fully charged in 90 minutes.

#### ■ (Battery Pack)

## EvoFill Duo Charger Function & Charging Method

### ▶ Battery replacement

### \* Warning

- Only use the battery provided by the manufacturer.
  - To take the battery out, press the hook at the bottom of the handpiece with the Multi-Tool and pull the battery out.
  - Insert a new battery pack into the handpiece until it snaps into place.
  - Do not apply any physical force or electrical shock to the battery pack. It may cause a fire.
- ▶ Battery must be inserted inside EvoFill and EvoPack in order to charge.  
Contact a local dealer or the manufacturer to purchase a new battery pack.  
Only use the charger provided by the manufacturer.

## Checklist & Maintenance

### 1. Checklist

- In order to prevent contamination, the handpiece must be cleaned thoroughly before and after use.
- Make sure the handpiece and other parts are not damaged or defective.
  - Make sure the Gutta Percha Cartridge (EvoFill) and reusable EvoFill cover are intact without cracks, foreign substance, or dirt.
  - Make sure the connecting part of the Gutta Percha Cartridge (EvoFill) and the heating chamber is clean and free of residual Gutta Percha.
  - Check the Heating Condenser (EvoPack) and make sure it is intact without cracks, foreign substance, or dirt.
  - Make sure the connecting part of Heating Condenser (EvoPack) and the heating chamber is clean.

### 2. Monthly Inspection

- Tighten all the connecting parts.
- Separate the Gutta Percha Cartridge (EvoFill) or Heating Condenser (EvoPack) from the handpiece. Wipe the inner and outer parts of the handpiece.

### 3. Yearly Inspection

- It is important to check if the device is performing safely and properly. Inspect the sheaths of the power cord and cables and make sure they are intact.

## Cleaning, Sterilization & Maintenance

### 1. How to clean the charger


- Remove the handpiece from the charger.
- Disconnect the charger from the power adaptor and outlet.
- Wipe the outer surface of the charger with a soft cloth dampened with a pH neutral disinfecting solution or rubbing alcohol.
- Do not soak the charger in water or any liquid.
- Clean the charger pin with a cotton swab and rubbing alcohol.

### 2. How to clean the handpiece

- Remove the handpiece from the charger.
- Wipe the outer surface of the handpiece with a soft cloth dampened with a pH neutral disinfecting solution or rubbing alcohol.
- Do not soak the handpiece in water or any liquid. (This product is not waterproof. IPX0)
- Clean the handpiece base with cotton swabs and rubbing alcohol.

### 3. How to clean the inner part of the EvoFill

#### 1) Take the EvoFill Safety Cap off.

#### 2) Turn the EvoFill Cartridge Nut and remove the Gutta Percha Cartridge (EvoFill).

- Users must use the Multi-Tool to separate the Gutta Percha Cartridge from the handpiece after use.

Press and hold the menu button for over a second so that the piston inside the device is moved to the very back.

Confirm the position of the piston through the GP window.


## Cleaning, Sterilization & Maintenance

Press the operation button so that the piston inside the device moves to the very front.


Confirm the position of the piston through the GP window.

#### \*Precaution:

Users must be careful as the piston automatically reverses to the back if the button is pressed for too long. (It takes 25-30 seconds for the piston to completely move to the back.)

### 3) Inside of EvoFill

- Lay the handpiece as shown above and remove all residual gutta percha (GP) using the scrubbing brush. If the residues are too hard to remove, turn on the device and melt them at 160 °C and try again.

### 4. EvoFill Safety Cap

- Cleaning method (Refer to Table 2)
- Wet-heat autoclave before and after use.  
(Wet-heat Autoclave at 121 °C for 20 minutes and dry at room temperature for 30 minutes.)

### 5. Heating Condenser (EvoPack) Sterilization (Refer to Table 2)

- Turn the device off.

- To remove the Heating Condenser (EvoPack), pull out the bent part covered with silicone.
- Cleaning method (Refer to Table 2)
- Wet-heat autoclave before and after use.  
(Wet-heat Autoclave at 121 °C for 30 minutes and dry at room temperature for 30 minutes.)

## Cleaning, Sterilization & Maintenance

### 6. Scrubbing Brush Sterilization

- Clean with alcohol and air dry.

### 7. How to reassemble

- Refer to operation Method.

### 8. Disposal

#### A. Disposal of the handpiece and charger

Remove the battery pack before disposal.

Dispose of electronics as electronic wastes.

Dispose of all other parts according to the relevant electronic product waste regulations.

#### B. Disposal of the battery pack

Dispose of the battery pack as a chemical waste.

## Table 1 (Safety Cap: EvoFill)

WARNING	
Limitation on reprocessing	Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use

INSTRUCTIONS			
Point of use:	Remove foreign substances with disposal cloths or papers.		
Containment and transportation:	No particular requirements		
Preparation before cleaning	Enzymatic Detergent should be diluted with water		
	Material		Enzymatic detergent, Ultrasonic cleaner, Nylon brush, dust absorbent, Alcohol, Tap water, Purified water
	Step	Procedure	Descriptions
	1	Pre wash	Wash the contaminated EvoFill Safety Cap to cold softened tap water for 30 seconds.
Cleaning: Automated	2	Enzyme soak	Soak diluted medizyme (complex enzyme detergent) for 3 minutes.
	3	Rinse	Rinse the EvoFill Safety Cap in tap water for 30 seconds and repeat once more.
	4	Detergent wash	Wash the EvoFill Safety Cap with diluted Medizyme at 40 °C ~ 45 °C for 3 minutes in ultrasonic cleanser. And then rinse the EvoFill Safety Cap with tap water for 30 seconds three times.
	5	Purified Rinse	Wash the EvoFill Safety Cap with purified water at 40 °C ~ 45 °C for 3 minutes in ultrasonic cleanser. This cleaning process is repeated two times more.
	6	Dry	Dry at 90 °C (±10 °C) for 20 minutes with hot air.
		Material	
	Step	Procedure	Descriptions
Cleaning: Manual	1	Remove Soil	Completely submerge the test sample (EvoFill Safety Cap) in diluted Medizyme (complex enzyme detergent) and allow to soak for 20 minutes. Scrub using a soft-bristled or nylon brush until visible soil (microorganisms) has been removed.
	2	1 <sup>st</sup> Rinsing	Rinse the EvoFill Safety Cap with tap water for 3 minutes to remove Medizyme(complex enzyme detergent)
	3	Ultrasonic cleaning	Wash the EvoFill Safety Cap with diluted medizyme(complex enzyme detergent) for 10 minutes in ultrasonic cleanser
	4	2 <sup>nd</sup> Rinsing	Wash the EvoFill Safety Cap with purified water for 3 minutes in ultrasonic cleanser. This cleaning process is repeated once more
	5	Alcohol disinfection	Soak the EvoFill Safety Cap in alcohol (70%) for 30 seconds, and the wipe EvoFill Safety Cap with dust absorbent.

**Table 1 (Safety Cap : EvoFill)**

Inspection and Function Testing	<ol style="list-style-type: none"> <li>1. Check the cleaned EvoFill Safety Cap holder.</li> <li>2. Check for any foreign substances that became attached on to the EvoFill Body.</li> <li>3. Check the connection between the body and the Safety Cap.</li> <li>4. Check the EvoFill Safety Cap for possible cracks or damages as a defective Safety Cap may cause injuries to patient.</li> <li>5. Check the connection between the protector and the Safety Cap.</li> </ol>
Packaging	No particular requirements. Safety Cap can be easily cracked or corroded.
Sterilization	Autoclave, minimum 121 °C for 30 minutes, Dry 20 minutes.
Storage	Keep away from any place where products may be damaged.
Additional Information	No particular requirements

**Table 2 (Heating Condenser : EvoPack)**

<b>WARNING</b>	
Limitation on reprocessing	Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use

<b>INSTRUCTIONS</b>			
Point of use:	Remove foreign substances with disposal cloths or papers.		
Containment and transportation:	No particular requirements		
Preparation before cleaning	Enzymatic Detergent should be diluted with water		
	<b>Material</b>	Enzymatic detergent, Tap water, Purified water, Automatic Cleaner	
	<b>Step</b>	<b>Procedure</b>	<b>Descriptions</b>
	1	Pre wash	Wash the contaminated Heating Condenser (EvoPack) to cold softened tap water for 30 seconds.
Cleaning: Automated	2	Enzyme soak	Soak diluted medizyme (complex enzyme detergent) for 3 minutes.
	3	Rinse	Rinse the Heating Condenser in tap water for 30 seconds and repeat once more.
	4	Detergent wash	Wash the Heating Condenser with diluted Medizyme at 40 °C ~ 45 °C for 3 minutes in ultrasonic cleanser. And then rinse the Heating Condenser with tap water for 30 seconds three times.
	5	Purified Rinse	Wash the Heating Condenser with purified water at 40 °C ~ 45 °C for 3 min. in ultrasonic cleanser. This cleaning process is repeated two times more.
	6	Dry	Dry at 90 °C (±5 °C) for 30 minutes with hot air.

**Table 2 (Heating Condenser : EvoPack)**

	<b>Material</b>	Enzymatic detergent, Alcohol, Tap water, Purified water, Ultrasonic Cleaner	
	<b>Step</b>	<b>Procedure</b>	<b>Descriptions</b>
Cleaning: Manual	1	Remove Soil	Completely submerge the Heating Condenser (EvoPack) in diluted Medizyme (complex enzyme detergent) and allow to soak for 20 min. Scrub using a soft-bristled or nylon brush until visible soil (microorganisms) has been removed.
	2	1 <sup>st</sup> Rinsing	Rinse the Heating Condenser (EvoPack) with tap water for 3 minutes to remove Medizyme(complex enzyme detergent)
	3	Ultrasonic cleaning	Wash the Heating Condenser (EvoPack) with diluted medizyme (complex enzyme detergent) for 10 minutes in ultrasonic cleanser
	4	2 <sup>nd</sup> Rinsing	Wash the Heating Condenser (EvoPack) with purified water for 3minutes in Ultrasonic Cleaner. This cleaning process is repeated once more
	5	Alcohol disinfection	Soak the Heating Condenser (EvoPack) in alcohol (70%) for 30 seconds, and the wipe Heating Condenser (EvoPack) with dust absorbent.

Inspection and Function Testing	<ol style="list-style-type: none"> <li>1. Check the cleaned Heating Condenser (EvoPack) holder.</li> <li>2. Check for any foreign substances that became attached on to the EvoPack Body.</li> <li>3. Check the connection between the body and the Heating Condenser (EvoPack).</li> <li>4. Check the Heating Condenser (EvoPack) for possible cracks or damages as a defective Heating Condenser (EvoPack) may cause injuries to patient.</li> <li>5. Check the connection between the protector and the Heating Condenser.</li> </ol>
Packaging	No particular requirements Heating Condenser (EvoPack) can be easily cracked or corroded.
Sterilization	Autoclave, minimum 121 °C for 30 minutes, Dry 30 minutes
Storage	Keep away from any place where products may be damaged.
Additional Information	No particular requirements

## ■ Guidance and DiaDent's declaration - electromagnetic emissions

### 1.1 Electromagnetic emissions

EvoFill Duo is intended for use in the electromagnetic environment specified below.

The customer or the user of EvoFill Duo should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment
RF emissions CLSPR 11	Class B	The EvoFill Duo is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

### 1.2 Electromagnetic immunity

EvoFill Duo is intended for use in the electromagnetic environment specified below.

The customer or the user of EvoFill Duo should assure that it is used in such as environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment
Electrostatic discharge(ESD) IEC61000-4-2	±6kV : Contact ±8kV : Air	±6kV : Contact ±8kV : Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%

### 1.3 Electromagnetic immunity 2

EvoFill Duo is intended for use in the electromagnetic environment specified below.

The customer or the user of EvoFill Duo should assure that it is used in such as environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment
Radiated RF IEC61000-4-3	3V/m 80MHz to 2.5GHz	3V/m 80MHz to 2.5GHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the EvoFill Duo, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance  <math>d = 1.2\sqrt{P}</math>  <math>d = 1.2\sqrt{P}</math> 80 MHz to 800 MHz  <math>d = 1.2\sqrt{P}</math> 800 MHz to 2.5GHz                      where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1. At 80MHz to 800MHz, the separation distance for the higher frequency range applies.  
2. These guidelines may not apply in all situations. electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio(cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which EvoFill Duo is used exceeds the applicable RF compliance level above, the FA should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the EvoFill Duo.  
b Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3 V/m.

1.3 Recommended separation distances between portable and mobile RF communications equipment and EvoFill Duo. EvoFill Duo is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of EvoFill Duo can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and EvoFill Duo as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 1.2\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1.0	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer.

Note 1. At 80MHz to 800MHz, the separation distance for the higher frequency range applies.  
2. These guidelines may not apply in all situations. electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## ■ Guidance and DiaDent's declaration - electromagnetic emissions

### Warning

Potential electromagnetic or other interference between medical equipment and other devices being operated together in the same environmental may exert an adverse influence on functioning of the medical equipment. Non-medical equipment not in compliance with the requirements of EN60601- 1 and EN60601-1-2 should not be used together in the same environmental as the medical equipment.

This equipment has been tested and found to comply with the limits for medical devices in IEC60601-1-2:2001. This limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

The guide of disposal procedure of old electric appliance

1. When this crossed out wheeled bin symbol is attached to a product it means the products is covered by the European Directive 2002/96/EC.
2. All electrical and electronic products should be disposed of separately from the municipal waste stream via designated collection facilities appointed by the government or the local authorities.
3. The correct disposal of your old appliance will help prevent potential negative consequences for the environment and human health.
4. For more detailed information about disposal of your old appliance, please contact city office, waste disposal service or the shop where you purchased the product.

## ■ Symbols

Symbol	Description	Symbol	Description
	Direct Current		"ON"/"OFF"(push-push) NOTE Each position, "ON"/"OFF". is a stable position
	Type B Applied Part		Refer to instruction manual or booklet
	Class II equipment		Humidity limitation
	Caution		Atmospheric Pressure limitation
	Sitting prohibited		Warning: Hot surface
	Manufacturer		Manufacture date
	Serial No		DO NOT REUSE
	EU Representative		e-Waste

## ■ Troubleshooting

EvoPack or EvoFill does not operate:	<ul style="list-style-type: none"> <li>- Press and hold the power button for 2 seconds.</li> <li>- Check if the battery is properly placed inside.</li> <li>- Charge or replace the battery.</li> </ul>
The operation button does not work with the error indicator flashing:	<ul style="list-style-type: none"> <li>- Check that the power is on.</li> <li>- Check that the Heating Condenser (EvoPack) is inserted properly.</li> <li>- Replace the Heating Condenser (EvoPack) and try again.</li> <li>- If the operation button still does not work, please contact your local dealer or manufacturer.</li> </ul>
Gutta Percha does not come out when the operation button is pressed	<ul style="list-style-type: none"> <li>- Check that the power is on.</li> <li>- Check if the device is preheating. Use after preheating is complete.</li> <li>- Check if the motor works.</li> <li>- Check if the Gutta Percha Cartridge is placed inside.</li> </ul>
The error indicator is on:	<ul style="list-style-type: none"> <li>- Turn off the power and restart.</li> <li>- If the same problem persists, please contact a local dealer or manufacturer.</li> </ul>
The power does not turn on:	<ul style="list-style-type: none"> <li>- Fully charge the battery before using the device.</li> <li>- Detach the battery pack and insert it again.</li> <li>- If the same problem persists, please contact a local dealer or manufacturer.</li> </ul>

These examples are general guidelines to the most common issues.

Please consult with your representative or contact a local service center about problems that are not listed here.

## ■ Product Warranty

1. This product is warranted by DiaDent for 1 year from the date of purchase against manufacturing defects (unless otherwise stated).  
(Note: Consumable parts, such as the EvoFill Safety Cap, Heating Condenser (EvoPack), and battery, are not eligible for warranty service.)
2. DiaDent will not accept any liability or provide any repair service in case of :
  - Non-compliance with the operation manual provided
  - Operation, handling, cleaning, and storage without care
  - Damages caused by natural disasters
  - Damages resulted from improper source of electricity

## ■ Specifications

ITEM		Specification	Remark
EvoPack	SIZE (mm)	182 x 25 x 29mm	L x W x H
	Weight (g)	102 ± 1g	Including battery
	Working Voltage (V)	3.6V $\overline{=}$	
	Operating Temperature (°C)	170 °C ~ 250 °C	
	Standby / Auto power off	5min. / 10min.	
EvoFill	SIZE (mm)	222 x 25 x 29mm	L x W x H
	Weight (g)	190 ± 1g	Including battery
	Working Voltage (V)	3.6V $\overline{=}$	
	Operating Temperature (°C)	160 °C ~ 200 °C	
	Standby / Auto power off	5min. / 10min.	
Charger	SIZE (mm)	178 x 135 x 82mm	L x W x H
	Input	9.0V $\overline{=}$	
Battery Pack	SIZE (mm)	78 x 25 x 26mm	L x W x H
	Working Voltage (V)	3.6V $\overline{=}$	
	Battery Capacity	2500mAh	Li-ion Battery
Adapter	Model Name	BridgePower / BPM030S09N03	
	Input	100~240V, 50/60Hz, 1A	
	Output	9.0V $\overline{=}$ , 3A	
	SIZE (mm)	95 x 54 x 32mm	L x W x H
Operating Conditions	Temperature (°C)	10 °C ~ 40 °C	
	Atmospheric Pressure (kPa)	60 ~ 101Kpa	
	Humidity (%)	30 ~ 85%	
Transportation and Storage Conditions	Temperature (°C)	-10 °C ~ 40 °C	
	Atmospheric Pressure (kPa)	60 ~ 101Kpa	
	Humidity (%)	30 ~ 85%	

**EC REP** Authorized representative in the European Community

### DiaDent Europe B.V.

Antennestraat 70, 1322 AS Almere, The Netherlands

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