

Dia Dent Conduction & Backful

Cordless Warm Vertical Compaction & Backfill Obturation System

Revolutionize your practice with the new EvoFill[™] Duo cordless obturation system. With an innovative electric motor that prevents hand fatigue, EvoFill Duo offers fast heating and controlled gutta percha extrusion with precise 3D fills and reliable results.

EvoFill[™] Backfill Obturation System

- Controlled, precise motorized gutta percha extrusion
- Fast heating 360° swivel cartridge needle
- 3 different temperature settings to control extrusion speed
- Hygienic single-use, fast loading gutta percha cartridge
- No mess, easy clean up

EvoPackTM Vertical Compaction Device

- Innovative LED light-guided condenser
- Save time with instant heating
- Tightly seals all canals including lateral canals
- Adjustable temperature setting and time
- Simple one-button operation

Revolutionary Evolution ...

Evolution: a process of formation or growth; continuous development of change, progression, metamorphosis.

Ask your personal dealer representative for detailed pricing and product information. Demo unit is available upon request. Instructional video is available on **www.diadent.com** DiaDent Group International www.diadent.com | 1.877.342.3368 Follow us on FACEBOOK




DiaDent

ROOT CANAL OBTURATION TECHNICAL SEQUENCE

💿 Vertical Condensation with 🖭 🛛 Pock 🏾






- 1a. Shape the root canal with endo files and determine the working length. Select a master cone that closely matches the apical size.
- 1b. Insert Heating Condenser (EvoPack) into the root canal to the point 4~7mm short of the working length. Move and set the endo stop accordingly.
- 2. Insert the selected master Gutta Percha cone that is coated with sealer.
- 3a. Place Heating Condenser (EvoPack) in the orifice of the canal and activate EvoPack by pressing the Operation Button. Using the activated Heating Condenser, cut the middle of Gutta Percha cone and downpack the soften material to the point 6~9mm from the apex for about 2 seconds only.
- 3b. Release the Operation Button to allow the Heating Condenser to cool down, while at the same time continue pushing and holding the cooled Heating Condenser to the apex for 8 seconds. This ensures all lateral and accessory canals are filled. Activate the Heating Condenser again by pressing and holding the Operation Button for about 1.5 seconds, compact the Gutta Percha cone to the point 4~7mm of working length.
- 4. Condense the remaining material using a hand plugger.

ORDER INFORMATION

4005-1001 4005-1201 4005-2010 4005-2020 4005-2020 4005-2050 4005-2070 4005-2071 4005-2072 4005-2073 4005-2074 4005-2075 4005-2075	EvoFill Duo Regular Kit EvoFill Handpiece EvoPack Handpiece Gutta Percha Cartridge (20/box) Charger (for EvoFill & EvoPack) Power Adapter and Cord Battery Heating Condenser XF Heating Condenser F Heating Condenser FM Heating Condenser M Heating Condenser ML Heating Condenser ML Heating Condenser Scombo F/FM (2/box) Heating Condensers Assorted XF/ML (5/box)
1 3 73	5



- 🕤 🛛 Backfill Obturation with 🖉 🖉 🖓 🖓 🖉 🖓 🔊
 - 5. Using EvoFill's activated Gutta Percha Cartridge Tip, make contact with the upper part of filled master cone in the canal, transfer heat for 3 to 5 seconds to soften the material.
 - 6a. Extrude the melted Gutta Percha into the root canal.
 - 6b. Continue flowing while withdrawing the tip from the canal.
 - 7. Use a hand plugger to compact the obturation material.
 - EvoFill[™] Duo Complete Kit (1YEAR WARRANTY)

Contains: 1 EvoFill Handpiece, 1 Safety Cap, 1 EvoPack Handpiece, 2 Heating Condensers, 1 Dual Charger, 1 Power Adapter & Cord, 1 Scrubbing Brush, 1 Multi-Tool, User Manual.

• Gutta Percha Cartridge (EvoFill) •

Contains: 20 tips/box

SPECIFICATIONS

EvoPack	Size (mm) L x W x H	182 x 25 x 29 mm
	Weight (g)	102 +- 1 g (incl. battery)
	Standby/Auto Power Off	5 min./10 min.
EvoFill	Size (mm) L x W x H	222 X 25 X 29 mm
	Weight (g)	190 +- 1 g (incl. battery)
	Standby/Auto Power Off	5 min./10 min.
Charger	Size (mm) L x W x H	178 x 13 x 82 mm
	Charging time	Approx. 90 minutes
Battery	Size (mm) L x W x H	78 x 25 x 26 mm
Adapter	Input	100~240V, 50/60Hz, 1A
Operating Condition	Temperature	10°C~ 40°C
	Atmospheric Pressure	60~101 kPa
	Humidity	30~85%