

Arthrometers Contribute to Accuracy in Orthopedic Assessment

Orthopedics is about precision. We use a variety of devices to collect objective data on/from patients. These include goniometers, hand-held dynamometers, isokinetic machines, pulse-oximeters, and blood pressure cuffs, to name a few. These data contribute to other clinical techniques to assist the practitioner in the identification of pathology. A lesser utilized device to quantify joint mobility is an arthrometer. Unlike a goniometer which measures angles in degrees of motion, i.e., flexion or extension (osteokinematic motion), arthrometers measure linear translation of a joint in millimeters of motion (arthrokinematics/accessory motion). This accessory (gliding/sliding) motion that occurs between two joint surfaces contributes to the osteokinematic motion. Without this accessory motion, joint surfaces would fail to maintain optimal contact as concave and convex surfaces move on each other.

TYPES OF ARTHROMETERS

In the 1980's, Med Metrics Corporation developed a knee arthrometer called the KT1000. For the first time, injury to the anterior cruciate ligament (ACL) could be quantified. The device paved the way for a plethora of research studies and enhanced understanding of the magnitude of translation that defined pathology. Around 1990, the Rolimeter was developed and brought to market (Aircast Europe). Currently, neither of these devices are commercially available. However, several other arthrometers such as the KLT (Karl Storz, Tuttlingen, Germany), KiRA (I+, Italy), and GNRB (Genourob, France) entered the market. In 2020, the Automatic Knee Arthrometer (Peking University) was developed and reported to be able to successfully determine the degree of knee laxity in ACL injuries.¹ Yet, this device is not commercially available. In 2021, the Mobil-Aider arthrometer (Therapeutic Articulations, LLC) entered the market and has begun to populate the literature with several studies.²⁻⁸ The Mobil-Aider and Telos are the only arthrometers that can be used on multiple joints, including the upper extremity. However, the Telos is used for positioning in conjunction with stress radiographs, so it is not truly a stand alone nor a portable device.

KNEE ARTHROMETER DATA

Magnetic resonance imaging (MRI) has been the gold standard for ACL diagnosis, but MRI's are expensive and are static images. Arthrometers can provide dynamic data. There have been numerous studies to evaluate the data obtained from arthrometers to assess anterior translation of the tibia on the femur (Lachman test). A study by Bach⁹ using the KT1000 reported a 0.2 ± 1.6 mm difference between normal knees but 4.8 ± 3.7 mm for acute ACL injuries and 5.5 ± 4.5 mm in chronic ACL injuries. However, Barcellona et al¹⁰ found the KT1000 to over-estimate sagittal knee joint laxity. Other studies¹¹⁻¹⁴ suggested less than 3 mm difference between knee translations was normal but more than 3 mm implied instability.

Klasan et al¹⁵ compared the KT1000 and the GNRB. Maximal anterior translation of healthy knees measured 5.7 mm with the KT1000 and 4.4 mm with the GNRB with an increase of 3 mm or more indicating laxity. Whereas, Michel et al¹⁶ reported a 4 mm difference was consistent with laxity using the GNRB. Runer et al¹⁷ reported maximal anterior tibial translation as 5.8 mm for the Rolimeter, 6.0 mm for KLT, 12.3 mm for KiRA, and 7.7 mm for the KT1000. ICC (intra-rater) ranged from 0.7 to 0.9. Side-to-side differences in healthy knees were also reported with the Rolimeter, KLT, KiRA, and KT1000 as 0.8 mm, 1.1 mm, 2.4 mm, and 1.1 mm, respectively. The authors set equivalence boundaries to ± 1 mm and thus determined results were equivalent for the Rolimeter, KLT, and KT1000 but not for the KiRA. Murgier et al¹⁸ also compared four arthrometers. However, they reported the threshold for side-to-side differences of sagittal knee laxity were 4.28 ± 2.25 mm for the GNRB, 6.05 ± 3.81 mm for the Telos, 4.15 ± 1.76 mm for the KT1000, and 2.82 ± 1.40 mm for the Rolimeter.

As one can see, most of the studies performed on these different arthrometer are comparing them to one another. If a device has not been validated, why would it be used as a comparison? Why aren't devices being compared (validated) to a gold standard? This validation process was performed on the Mobil-Aider arthrometer via bench research using the Zeis Smartzoom (ICC = 0.997).^{2,3} In addition, a clinical study compared healthy knees to those with partial and complete ACL tears as confirmed by an MRI.⁴ The double-blinded data revealed side-to-side differences in sagittal translation (Lachman test) of ± 0.18 mm in the

healthy knees, 2.05 mm in partial ACL tears, and more than 3.38 mm in complete ACL tears with the Mobil-Aider.

ANKLE ARTHROMETER DATA

For the ankle, there are only three arthrometers available: Instrumented Ankle Arthrometer (Blue Bay), Telos SD900 (Telos Medical) and the Mobil-Aider (Therapeutic Articulations). A study validating the Instrumented Ankle Arthrometer reported approximately 2 mm differences between the control and ankle laxity group.¹⁹ This was consistent with the work of Kerkhoff²⁰ and Kovaleski.²¹ Kerkhoff²⁰ reported sectioning the anterior talofibular ligament in cadavers resulted in an increase of 2 mm of anterior laxity. Further, sectioning the calcaneal fibular ligament added another 2 mm of laxity.

A study performed with the Mobil-Aider arthrometer compared healthy ankles to those with first degree (1°) and second degree (2°) lateral ankle sprains.⁵ When testing the anterior talofibular ligament with the Mobil-Aider using an anterior drawer technique, a 1.11 mm side-to-side difference was consistent with a 1° sprain and 2.16 mm difference indicated a 2° sprain. These ratios correspond with a study by Lee et al¹⁴ using stress radiographs. When a ratio of ATFL stress to resting was 1.1, 1.3, and 1.4, the sprain was classified as a 1°, 2°, and 3°, respectively. The ability to have objective values of anterior laxity and the corresponding ratios can significantly impact the treatment decision.

UPPER EXTREMITY ARTHROMETER DATA

The Mobil-Aider is the only arthrometer available for upper extremity testing. When analyzing shoulder posterior glides, the Mobil-Aider arthrometer was found to have a 0.83 correlation with electro-magnetic motion analysis.⁸ Two studies were also performed on shoulder and wrist accessory motions.^{6,7} Intra-rater reliability was reported to be 0.691, 0.771, and 0.789 for posterior shoulder glides and 0.904, 0.917, and 0.932 for wrist volar glides. Whereas, inter-rater reliability for shoulder posterior glides were 0.546 and wrist volar glides were 0.462.

CLINICAL APPLICATION

Testing with arthrometers provide objective data . By placing a numeric value on the magnitude of motion, one can objectively compare right to left joints. While developing normative values may not be practical, data on differences between injured and uninjured joints can be clinically significant. As in other manual techniques, intra-rater reliability has been found to be good to excellent but inter-rater reliability is poor to fair with or without an arthrometer. There are many factors that contribute to this issue. Maximal manual force application has shown low agreement with dependency on investigator gender.²² Because linear translation can be influenced by the force applied during excursion, side-to-side differences are recommended for clinical assessment.^{17,23} Furthermore, many studies describe clinician testers as novice or experienced. Although years of clinical experience may be important, familiarity with the arthrometer is much more important.¹⁷ Easy to operate arthrometers used in a way that is consistent with standard testing techniques can reduce the learning curve and enhance reliability.

In summary, there is technology available to provide objective data on joint laxity that fills the void between subjective assessment and expensive imaging. Clinicians should consider incorporating the use of arthrometers into practice to provide objective data to drive clinical decisions and to populate the literature about the efficacy of manual techniques.

REFERENCES

1. Niu X, Mai H, Wu T, et al. Reliability of a novel automatic knee arthrometer for measuring laxity after anterior cruciate ligament ruptures. *Orthopaedic Journal of Sports Medicine*. 2022,10(2):1-7
2. Gulick DT. Proof of Concept: Taking the guessing out of assessing knee stability. *International Journal of Sports and Exercise Medicine*. June 2019
3. Gulick DT. Reliability and validity of an innovative device for ACL testing: The Mobil-Aider™. *Journal of Sport Rehab*. 2019 May 16 ePub; 2020 Feb 1;29(2):257-261.
4. Hammoud S, Palombaro K, Gulick DT. Use of a New Arthrometer to Assess Knee Pathology. *Global Journal of Orthopedic Research*. January 2022.

5. Taweel NR, Gulick DT, Palombaro KM. Assessing Lateral Ankle Sprains with a New Arthrometer. *Foot & Ankle Specialist*. August 2022
6. Wise CH, O'Donohue JM. Intra-rater and inter-rater reliability of the Mobil-Aider® device for measurement of linear translation: Implications for clinical practice and teaching. *World Journal of Physical & Rehabilitation Medicine*. January 2022.
7. O'Donohue JM, Wise CH. Measurement of accessory motion of the glenohumeral and radiocarpal joints: Intra-rater reliability of the Mobil-Aider® device for measurement of linear translation. *Ann Physiother Clin*. October 2021; 3(1):1014.
8. Tuzson A, Tarleton G. Validating the Mobil-Aider to measure joint accessory motion in healthy adult shoulders. *Journal of Health, Science & Medicine*. October 2021
9. Bach BR Jr, Warren RF, Flynn WM, Kroll M, Wickiewicz TL. Arthrometric evaluation of knees that have a torn anterior cruciate ligament. *J Bone Joint Surg Am*. 1990;72:1299-1306
10. Barcellona MG, Christopher T, Morrissey MC. Bench testing of a knee joint arthrometer. *Orthopedics*. 2013;36(8):e1000-e1006.
11. Arneja S, Leith J. Review article: Validity of the KT-1000 knee ligament arthrometer. *J Orthop Surg (Hong Kong)*. 2009;17(1):77-79.
12. Daniel DM, Stone ML, Sachs R, Malcom L. Instrumented measurement of anterior knee laxity in patients with acute anterior cruciate ligament disruption. *Am J Sports Med*. 1985;13:401-407.
13. Rangger C, Daniel DM, Stone ML, Kaufman K. Diagnosis of an ACL disruption with KT-1000 arthrometer measurements. *Knee Surgery Sports Traumatology & Arthroscopy*. 1993;1(1):60-66
14. Lee KT, Park YU, Jegel H, Park JW, Choi JP, Kim JS. New method of diagnosis for chronic ankle instability: comparison of manual anterior drawer test, stress radiography and stress ultrasound. *Knee Surgery Sports Traumatology & Arthroscopy*. 2014;22:1701-1707.
15. Klasan A, Putnis SE, Kandhari V, Oshima T, Parker DA. Anterior knee translation measurements after ACL reconstruction are influenced by the type of laximeter used. *Knee Surgery, Sports Traumatology, Arthroscopy : Official Journal of the ESSKA*. 2020 Nov;28(11):3639-3646

16. Michel E, Jordan E, Canovas F, Bouchdoug K, Dagneaux L, Gaillard F. Influence of residual anterior laxity on function outcomes after anterior cruciate ligament reconstruction. *Orthopedic Traumatology Surgery & Research*. 2022 May;108(3):103264.
17. Runer A, di Sarsina TR, Starke V, et al. The evaluation of Rolimeter, KLT, KiRA and KT-1000 arthrometer in healthy individuals shows acceptable intra-rater but poor inter-rater reliability in the measurement of anterior tibial knee translation. *Knee Surg Sports Traumatol Arthrosc*. 2021;29(8):2717-2726.
18. Murgier J, Béranger JS, Boisrenoult P, Steltzlen C, Pujol N. Prospective comparative study of knee laxity with four different methods of anterior cruciate ligament tears. *International orthopaedics*. 2018 Aug; 42 (8):1845-1851
19. Disanto TJ, Swanik B, Swanik KA, Straub SJ, Needle AR. Concurrent validity of the anterior drawer test and an arthrometer in evaluating ankle laxity. *Athletic Training & Sports Health Care*. 2011;3(1):15-20
20. Kerkhoffs GM, Blankevoort L, van Poll D, Marti RK, van Dijk CN. Anterior lateral ankle ligament damage and anterior talocrural joint laxity. *Clinical Biomechanics*. 2001;16(8):635-643.
21. Kovalski JE, Hollis J, Heitman, Gurchiek LR, Pearsall AW 4th. Assessment of ankle subtalar joint complex laxity using an instrumented ankle arthrometer: An experimental cadaveric investigation. *Journal of Athletic Training*. 2002;3:467-474
22. Klasan A, Putnis SE, Kandhari V, Oshima T, Fritsch BA, Parker DA. Healthy knee KT1000 measurements of anterior tibial translation have significant variations. *Knee Surg Sports Traumatol Arthroscopy*. 2020 Jul;28(7):2177-2183
23. Wordeman SC, Paterno MV, Quatman CE, Bates NA, Hewett TE. Arthrometric curve-shaped variables to assess anterior cruciate ligament deficiency. *Clinical Biomechanics*. 2012;27:830-836