# Optimized Constants for AF-1 / Nanex™ Preloaded IOLs

(based on clinical results published on IOLCon and optical biometry mostly from the Zeiss IOLMaster)

| Model name | Specifications     | Nominal | Barrett | Haigis         |        |       | Hoffer Q | Holladay-1 | SRK/T | Castrop |   |   |
|------------|--------------------|---------|---------|----------------|--------|-------|----------|------------|-------|---------|---|---|
|            |                    | Α       | LF      | a <sub>o</sub> | $a_1$  | $a_2$ | pACD     | sf         | Α     | С       | Н | R |
| AF-1       |                    |         |         |                |        |       |          |            |       |         |   |   |
| PY-60AD    | 3-Piece, asph.     | 118.4   | 1.67    | -0.093         | -0.023 | 0.208 | 5.30     | 1.54       | 118.6 | -       | - | - |
| PC-60AD    | 3-Piece, asph.     | 118.4   | 1.67    | -0.093         | -0.023 | 0.208 | 5.30     | 1.54       | 118.6 | -       | - | - |
| PY-60R     | 3-Piece, spherical | 118.4   | 1.62    | 1.060          | 0.400  | 0.100 | 5.24     | 1.48       | 118.5 | -       | - | - |
| PC-60R     | 3-Piece, spherical | 118.4   | 1.62    | 1.060          | 0.400  | 0.100 | 5.24     | 1.48       | 118.5 | -       | - | - |
| 251        | 1-Piece, asph.     | 118.4   | 1.62    | -0.542         | 0.161  | 0.204 | 5.30     | 1.52       | 118.5 | -       | - | - |
| 250        | 1-Piece, asph.     | 118.4   | 1.62    | -0.542         | 0.161  | 0.204 | 5.30     | 1.52       | 118.5 | -       | - | - |
| 255        | 1-Piece, asph.     | 118.4   | 1.62    | -0.542         | 0.161  | 0.204 | 5.30     | 1.52       | 118.5 | -       | - | - |
| 254        | 1-Piece, asph.     | 118.4   | 1.62    | -0.542         | 0.161  | 0.204 | 5.30     | 1.52       | 118.5 | -       | - | - |
| 151        | 1-Piece, spherical | 118.4   | 1.62    | -0.542         | 0.161  | 0.204 | 5.30     | 1.52       | 118.5 | -       | - | - |
| 150        | 1-Piece, spherical | 118.4   | 1.62    | -0.542         | 0.161  | 0.204 | 5.30     | 1.52       | 118.5 | -       | - | - |

As the data set for model 250/251 is larger on Ocusoft/ULIB than that submitted to IOLCon, the recommendation is to use the Manufacturer/ULIB constants rather than the optimized IOLCon values.

### Nanex™

| NY1-SP | 1-Piece, asph. | 119.2 | 1.94 | -0.2676 | 0.2382 | 0.1993 | 5.715 | 1.904 | 119.112 | - | - | - | ** |
|--------|----------------|-------|------|---------|--------|--------|-------|-------|---------|---|---|---|----|
| NC1-SP | 1-Piece, asph. | 119.2 | 1.94 | -0.2676 | 0.2382 | 0.1993 | 5.715 | 1.904 | 119.112 | - | - | - | ** |

<sup>\*</sup> http://ocusoft.de/ulib/c1.htm (as of Oct 31, 2016)



Constants for Holladay-2, Hill-RBF or KANE formulas are not published on Ocusoft or IOLCon. Optimized constants for these formulas can be calculated in some biometer devices based on the above-mentioned SRK/T values. The Olsen formula can be used for IOL power calculation in the Oculus Pentacam AXL and Haag-Streit Lenstar / Eyestar devices, as all required IOL parameters are included in the corresponding software.

These optimized constants for the calculation of intraocular lens power are based on actual surgical data and are provided by Ocusoft or IOLCon as a starting point for individual constant optimizations. The information available on these websites are based on data originating from other users and not by HOYA Surgical Optics ("HSO"). HSO therefore does not warrant the correctness, completeness and correctness of the contents on the said websites.



<sup>\*\*</sup> https://iolcon.org (as of May 4, 2023)

# Optimized Constants for Vivinex™ Preloaded IOLs

(based on clinical results published on IOLCon and optical biometry mostly from the Zeiss IOLMaster)

| Model name    | Specifications | Nominal<br>A | Barrett<br>LF | a <sub>o</sub> | Haigis<br>a <sub>1</sub> | a <sub>2</sub> | Hoffer Q<br>pACD | Holladay-1<br>sf | SRK/T<br>A | С      | Castrop<br>H | R      |           |
|---------------|----------------|--------------|---------------|----------------|--------------------------|----------------|------------------|------------------|------------|--------|--------------|--------|-----------|
| Vivinex™      |                |              |               |                |                          |                |                  |                  |            |        |              |        |           |
| XY1<br>XY1-SP | 1-Piece, asph. | 118.9        | 2.0           | -1.0453        | 0.2582                   | 0.2292         | 5.716            | 1.944            | 119.226    | 0.3249 | 0.1267       | 0.1548 | NEW<br>** |
| XC1<br>XC1-SP | 1-Piece, asph. | 118.9        | 2.0           | -1.0453        | 0.2582                   | 0.2292         | 5.716            | 1.944            | 119.226    | 0.3249 | 0.1267       |        | NEW<br>** |

## **Vivinex Impress™**

As there is currently no data set available for Vivinex Impress™ on IOLCon, the recommendation is to use the optimized constants for the Vivinex™ monofocal IOL also for Vivinex Impress™.

#### Vivinex™ Toric

As IOLCon data set is less than 100 cases for Vivinex™ Toric, the recommendation is to use the optimized constants for the Vivinex™ monofocal IOL also for Vivinex™ Toric.

### Vivinex™ Gemetric™ (Toric) / Vivinex™ Gemetric™ Plus (Toric)

<sup>\*\* &</sup>lt;a href="https://iolcon.org">https://iolcon.org</a> (as of May 4, 2023)



Constants for Holladay-2, Hill-RBF or KANE formulas are not published on Ocusoft or IOLCon. Optimized constants for these formulas can be calculated in some biometer devices based on the above-mentioned SRK/T values. The Olsen formula can be used for IOL power calculation in the Oculus Pentacam AXL and Haaq-Streit Lenstar / Eyestar devices, as all required IOL parameters are included in the corresponding software.

These optimized constants for the calculation of intraocular lens power are based on actual surgical data and are provided by Ocusoft or IOLCon as a starting point for individual constant optimizations. The information available on these websites are based on data originating from other users and not by HOYA Surgical Optics ("HSO"). HSO therefore does not warrant the correctness, completeness and correctness of the contents on the said websites.

