DE

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Akupressurmatte+ Kissen

Nutzung:

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Acupressure mat+ pillow

Usage:

The acupressure mat + pillow are suitable for lying or sitting on it to relax. The acupressure mat can be used with or without a T-shirt. If you have sensitive skin, it is recommended to use it with a T-shirt.

The acupressure mat can be placed on the floor, on a chair or on the sofa. The effect is more intense on harder surfaces (e.g. wooden floor). Therefore, for beginners it is recommended to use it on a soft surface, for advanced users it is recommended to use it on a hard surface.

Duration of use:

- Beginner: 10-15min
- Intermediate users: 20-45min
- Advanced users: 45min+

Washing Instructions:

- Remove lining from cover to clean
- Hand wash the cover
- Do not machine dry cover

Precautions

When using the acupressure mat, users must monitor the effect of the product on their body to avoid the risk of injury. Among other things, if you feel discomfort or pain, you should stop using the acupressure mat and consult a doctor.

Products should be kept safe from children.

Contraindications:

- Oncological diseases
- Coagulation disorders (hemophilia) or use of anticoagulants
- Epilepsy
- Significant skin damage in areas of application (eczema, psoriasis, shingles)
- Open wounds
- Pregnancy
- Fever or acute illness
- Diseases of the heart
- Severe form of various diseases
- Thrombophlebitis

For any simple or complex diagnosis, consult your doctor before using the acupressure mat.

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