

# Testsimplets<sup>®</sup>

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### *Procedure for differential blood count*

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
*Capillary blood  
take up a small drop of blood with the  
cover slip*

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*Venous blood  
apply directly to the cover slip*

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*Place the cover slip with the small drop  
of blood on the prestained area of the  
slide*

schnell  
*fast*

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*Ready for evaluation after only  
15 minutes*

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## Testsimplets®

- Ready-to-use slides prestained with standard dyes for time- and cost-saving differential diagnostic microscopy
- Simple and clean procedure for differential blood counts. No blood film preparation or staining procedures needed
- Tried and tested method for office practice, laboratory or hospital ward, particularly suitable for single specimens or small series

## Applications

- Routine evaluation of the differential blood count
  - differentiation of the normal WBC count and its reactive changes
  - visualization of all normal leukocyte forms and platelets
  - simultaneous detection of reticulocytes
- Cancer cytology
- Urinary cytology
- Cerebrospinal fluid cytology
- Nasal cytology
- Spermatozoal staining

## Test principle

- The stain is composed of the two dyes new methylene blue N and cresyl violet acetate in chemically pure form
- Constant mixing ratio and even distribution ensure reliably good staining quality
- Reaction principle analogous to that of the panchromatic (panoptic) stain

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### Characteristics

The simple and clean handling and the fast staining with standardized dyes make Testsimplets® ideal for use in hospitals and doctors' surgeries. Testsimplets® are especially suitable for examination of single samples and small series such as those that occur particularly in doctors' practices or in the laboratories of small hospitals. In addition, a differential count with Testsimplets® provides rapid results on the ward or during night duty in any hospital. In general practice the doctor can have the result of the differential count available before the end of the consultation.

### Clinical aspects

Many diseases are characterised by typical changes in the number and morphology of the blood cells. The differential count, in addition to other tests, provides useful information about the extent and nature of these changes. The differential count is thus clinically relevant in:

- infectious processes (reactive shifts)
- suspected allergic reactions (eosinophilia)
- suspected parasitic infection (eosinophilia)
- checking for drug side-effects (agranulocytosis, thrombocytopenia)
- lymphadenopathy of unknown origin (e.g. lymphocytosis)
- generalised bleeding of unknown origin (thrombocytopenia)
- suspicion of leukaemic diseases (pathological shift to the left)

The differential blood count stain with Testsimplets® also shows up reticulocytes. The reticulocytes in the peripheral blood reflect the regenerative capacity of erythropoiesis in the bone marrow. The reticulocyte count thus provides information about the red cell turnover.

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Detection of reticulocytosis is important, e.g.

- if a haemolytic disorder is suspected
- for monitoring treatment of iron-, folic acid- or vitamin B12 deficiency anaemia

## Procedure

Apply a small drop of capillary or venous blood (about 3 µl) to the middle of a cover slip and place on the pre-stained area of the slide in such a way that the blood spreads evenly in all directions. To obtain even distribution of the sample it can be useful to press gently on the centre of the cover slip with a pencil. Allow 15 minutes for staining and then differentiate the preparation at x800 to x1000 magnification with oil immersion. The preparation remains stable for at least 4 hours at room temperature.

If microscopic differentiation of the preparation cannot be performed until after 4 to 24 hours the preparation should be refrigerated immediately.

## Evaluation

For evaluation, select a site where the cells are spread thinly enough to allow easy identification. Count 100 leukocytes, scanning back and forth across the thin layer with the “battlement tracking technique”.

While scanning the slide it is advisable to keep adjusting the focus of the leukocytes by “playing” with the fine adjustment knob of the microscope. In doubtful cases – such as can occur with the identification of monocytes, band neutrophils and basophils – it is particularly important to scan the entire cell bringing individual cell regions into sharp focus.

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## Leukocyte characterisation

The leukocyte forms occurring in the peripheral blood of healthy individuals have the following characteristic features:

### Neutrophilic granulocytes (neutrophils)

A fine, colourless or faintly violet granulation can be seen in the cytoplasm. The nucleus stains bright purple and thus stands out clearly from the granulated cytoplasm.

### Band neutrophils

The nucleus of the band neutrophils is shaped like a band or rod, i.e. without segmentation.

### Segmented neutrophils

The shape of the nucleus depends on the degree of maturity of the cell and corresponds to familiar picture.

### Eosinophilic granulocytes

The granules are characteristic. The cytoplasm is completely or almost completely filled with clearly defined, bright yellow-green to green granules. They are larger than the granules of the basophils and neutrophils, transparent and appear swollen. If achromates are used a greenish colour prevails.

In the mature forms the nucleus often consists of two light purple segments and is thus clearly distinguishable from the granules on account of its colour.

### Basophilic granulocytes

The basophils are characterised by the different-sized, usually sharply delineated granules which stain purple, dark violet or orange to brown. "Playing" with the fine adjustment reveals orange-brown rings at the outer perimeter of the violet granules.

### Monocytes

Monocytes appear considerably smaller than on conventional blood films. The cytoplasm stains more or less purple. It appears indistinctly granulated without individual granules being distinguishable. After long staining times it can become just as purple as the nucleus. The contours of the nucleus form a sharp and clearly identifiable line which separates the nucleus distinctly from the cytoplasm.

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### Lymphocytes

The cytoplasm is more or less purple and contains no granules. In the small lymphocytes it forms a narrow, scarcely perceptible rim around the nucleus.

The nucleus is purple. The chromatin appears clumped, increasingly so at the nuclear membrane. Distinct nucleoli can often be seen. The nucleus is usually round.

Functional forms of the lymphocytes (reactive lymphocytes) can have nuclei of different shapes which closely resemble the nuclei of monocytes. This considerably complicates evaluation of the functional lymphocyte forms.

### Early leukocyte forms

Early leukocyte forms occurring in the peripheral blood can be seen but are difficult to differentiate.

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## Characterisation of other blood cells

### Platelets

The platelets appear as small purple discs. They are usually in clumps but can also occur singly. Marked thrombocytopenia can be identified.

### Erythrocytes

The erythrocytes stain a pale yellowish brown and are almost always swollen. This deformation is caused by influences of the dye on the vital cell. **Evaluation of erythrocyte morphology is therefore not recommended.**

It is not possible to estimate the haemoglobin content of the erythrocytes. Forms without haemoglobin formation can be confused with lymphocytes on cursory inspection.

### Reticulocytes

The reticulocytes are also deformed. They stain the same colour as the erythrocytes. The typical granular and filamentous structures (reticulin) stain a more or less intense purple.

## Reticulocyte counting

To determine the reticulocyte count, count all the erythrocytes (including reticulocytes) in a visual field and note down the number. Then count the number of reticulocytes in the same visual field. Repeat this procedure in different areas of the film until a total of 1000 erythrocytes have been counted and express the reticulocyte count in ‰.

## Anticoagulants


EDTA is recommended as anticoagulant for venous blood.

EDTA blood should not have been standing for more than 3 hours at room temperature before it is used (do not refrigerate).


Use of citrated blood can cause overstaining.

Oxalate and fluoride are unsuitable, likewise heparinate and blood from patients on heparin therapy.


Normale Leukozytenformen des peripheren Blutes  
*Leukocyte forms normally occurring in peripheral blood*


Stabkerniger  
*Band neutrophil*


Segmentkerniger  
*Segmented neutrophil*

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


---

## Pathologische Leukozytenformen des peripheren Blutes

### *Pathological leukocyte forms in the peripheral blood*

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Testsimplets® can be used in cancer cytology for acute testing in hospitals and for preliminary testing in the doctor's surgery.

Testsimplets® are ideal for rapid staining of cell material suspected of being cancerous, e.g. in bronchial cytology. Early morning sputum samples, aspirates of sputum and bronchial contents obtained by endoscopy or cell material of other origin (pleural, ascites, pericardial or pulmonary aspirates, gastric washings) can be examined under the microscope after staining for only a few minutes.

Testsimplets® thus permit extremely rapid diagnostic evaluation. This means that a cytological result can be obtained while the endoscopy is still in progress so that further diagnostic measures can be initiated immediately if necessary.

### Procedure

Wash native sputum using the Kitasato-Mulder method. Then apply the sputum pellet to the Testsimplets® slide using a platinum loop with a diameter of 5 mm, spread the material carefully and cover with a cover slip. To ensure that the preparation is evenly distributed it can be useful to apply gentle pressure to the middle of the cover slip with a pencil. Allow 15 minutes for staining and then examine the preparation under the microscope. Particularly viscid sputa require longer staining times.


Material obtained by aspiration of body cavities is concentrated by cytocentrifugation. Bronchial contents obtained by endoscopic bronchial aspiration are collected in physiological saline and then concentrated by cytocentrifugation. A small drop of the sediment is applied to the prestained area of the Testsimplets® slide and covered with a cover slip.

After a staining time of 2–4 minutes the preparation can be examined under the microscope.

The stained preparation remains stable for 4 hours at room temperature. If refrigerated, it may still be evaluated after up to 24 hours, or after 48 hours in the case of very viscid sputa.

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
*Native sputum: alveolar histiocytes and abundant undifferentiated lymphocyte-like tumour cells*  
*Histology: undifferentiated bronchial carcinoma*

## Evaluation


The Testsimplets® preparations are evaluated according to the criteria used for conventionally stained preparations (Papanicolaou classes I to V).

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*Endoscopic bronchial aspirate: several polymorphous tumour cells alongside normal ciliated epithelial cells and a goblet cell*  
*Histology: polymorphous bronchial carcinoma*

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
*Pleural effusion: large tumour cells with foamy plasma, one mitosis*

*Histology: large-cell, partially mucinous bronchial carcinoma*

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
**Histologie:** Ovarialkarzinom

*Pleural effusion: aggregate of adenocarcinoma cells*  
*Histology: ovarian carcinoma*

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

*Pleural effusion: aggregate of adenocarcinoma cells, metastases of an ovarian carcinoma*

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

*Native sputum: squamous cell carcinoma cells with some leukocytes lying alongside*

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
*Endoscopic aspirate: moderately keratinising squamous cell carcinoma*

Pleurapunktat: primärer Tumor, Pankreaskopfkarzinom

*Pleural aspirate: primary tumour, pancreas head carcinoma*


Pleuraerguß: metastasierendes Mammakarzinom

*Pleural effusion: metastatic breast cancer*

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

*Native sputum: cells of a dedifferentiated carcinoma, histologically confirmed*

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
*Sputum: undifferentiated carcinoma*

Bronchialsekret: normale Flimmerzellen

*Bronchial secretions: normal ciliated cells*

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*Endoscopic bronchial aspirate: normal ciliated epithelial cells and goblet cells*

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Testsimplets® allow extraordinarily good differentiation and classification of nuclear structures, nuclear membranes, chromatin and nucleoli as well as cytoplasmic structures.

Urinary cytology is useful in:

- all forms of microscopic and macroscopic haematuria
- refractory cystitis
- cystalgia
- dysuria of unclear origin
- monitoring after urothelial tumour surgery
- early detection and monitoring of bladder cancer

The criteria for tumour grading are the same as those used with conventional staining (Papanicolaou and Pappenheim).

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
### Procedure

The processing of the urine is similar to that for preparation of a normal urinary sediment: centrifuge the urine for about 10 minutes at 3000 rpm. Pour off the supernatant carefully and completely. Resuspend the sediment in 3 drops of physiological saline. Apply one drop of this suspension directly to the centre of one of the cover slips provided and then place the cover slip on the prestained area of the slide. To ensure even distribution of the sample it can be useful to press gently on the centre of the cover slip with a pencil. After about 5 to 10 minutes the staining process is completed and the preparation can be examined; at room temperature the preparation is stable for up to 5 hours after staining.


### Evaluation

Begin by surveying the preparation at a magnification of x100 and then examine any striking cells or cell formations at x200 to x400 magnification or at x1000 with oil immersion.


On account of the vital stain the stained cells appear in their natural size. In contrast to other methods, details especially of the nuclear structures and nuclear membrane are therefore extremely well visualized.


Normale Urothelzelle umgeben von Erythrozyten  
*Normal urothelial cell surrounded by erythrocytes*


Blasenkarzinom G 3, entdifferenzierte Zellen mit großen Kernen, Vakuolen und Kernkörperchen  
*Bladder carcinoma G3, dedifferentiated cells with large nuclei, vacuoles and nuclear bodies*

Blasenkarzinom G 2, hyperchromatische Kerne, zahlreiche Kernkörperchen, destruktive Veränderungen im Zytoplasma, mehrere Leukozyten in der Umgebung  
*Bladder carcinoma G2, hyperchromatic nuclei, numerous nuclear bodies, destructive changes in the cytoplasm, several leukocytes in the surrounding area*

Blasenkarzinom G 3, teilw. hyperchromatische Kerne, kleine Zellen, multiple Nukleolen, Kernwandhyperchromasie  
*Bladder carcinoma G3, nuclei partially hyperchromatic, small cells, multiple nucleoli, nuclear wall hyperchromasia*

Blasenkarzinom G 2, deutliche Kerngrößenvarianz, multiple Nukleolen, Kernplasmarelation zu Gunsten der Kerne verschoben  
*Bladder carcinoma G2, marked variation of nuclear size, multiple nucleoli, nucleus-to-cytoplasm ratio shifted in favour of the nucleus*

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Testsimplets® stain all cells occurring in the cerebrospinal fluid. Unlike conventional techniques, e.g. sedimentation chamber and cyto centrifugation, Testsimplets® is suitable for routine cytological testing of CSF:

- rapid processing and isolation of the CSF cells
- small volume of CSF (20 µl native CSF)
- detects even normal leukocyte counts (3–4 cells per µl) without mechanical concentration
- considerably less danger of cell lysis and cell denaturation
- greater proportion of cells suitable for differentiation

### Procedure

After mixing, apply about 20 µl of the cerebrospinal fluid sample to the prestained area of the Testsimplets® slide. Then immediately place a cover slip on the slide in such a way that the CSF is evenly spread in all directions. To ensure even distribution of the sample it can be useful to press gently on the centre of the cover slip with a pencil.

After a staining time of 10 to 15 minutes the preparation is ready for differentiation. Differentiation must be performed within one hour after preparing the slide.

If the slide is stored in a moist chamber at +4°C (Petri dish) it can be kept for up to 2 hours before differentiation. If stored at –20°C the CSF preparation is stable for at least 2 days.

### Evaluation

Evaluate the preparation at a magnification of x800 – x1000 with oil immersion. In preparations with low leukocyte counts 25–50 leukocytes should be counted, in preparations with high leukocyte counts 100–200. The lower detection limit of the method is 3–4 leukocytes per microlitre, i.e. within the normal range.

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## Characterisation of the leukocytes in cerebrospinal fluid

### Granulocytes

Neutrophilic granulocytes have violet-brown nuclei or nuclear segments with clumped chromatin. The cytoplasm is reddish to pink with yellowish-brownish granules.

Basophilic granulocytes have typical dark violet granules.

### Monocytes

Cells with a large eccentric nucleus which is oval, kidney-shaped or lobulated and stains red-brown to violet with a honeycomb-like structure. The nucleus is surrounded by abundant pink-violet, grainy cytoplasm.

### Lymphocytes

Cells with a round, sometimes slightly indented nucleus, which stains red-violet to brown-black, with distinct coarse chromatin clumping, surrounded by a narrow to moderately wide rim of pink to violet-brown cytoplasm without or with scarcely perceptible granules.

In contrast to the lymphocytes with nuclei of different sizes, transformed lymphocytes have abundant pink to red-brown cytoplasm with a few darker staining granules.

### Plasma cells

Cells with a large, roundish-oval marginal nucleus that stains dark brown with a cartwheel pattern of lighter areas and abundant red-brown cytoplasm with honeycomb-like inclusion structures.

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#### Erythrozyten

### Further CSF cells that can be shown with Testsimplets®

#### Monocyte-like cells

Monocytoid cells, macrophages and giant cells

#### Cells lining CSF spaces

Choroid plexus cells, ependymal cells, arachnoid mesothelial cells and leptomeningal cells

#### Glia cells


#### Tumour cells and cells suspicious for cancer

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



Testsimplets® can be used as a simple and reliable screening method for routine rhinological evaluation. In the differential diagnostic evaluation of rhinitis the cytology permits a distinction between:

- bacterial
- viral
- atrophic
- allergic and pseudoallergic inflammation

### Procedure

Collect material from the inferior and/or middle nasal meatus using a cotton wool swab. With children, material obtained by blowing the nose can be used for evaluation. Apply the specimen material to the Testsimplets® cover slip by spreading or rolling with the cotton swab and then press the cover slip gently onto the prestained area of the slide. The time between collection of the material and pressing the cover slip onto the slide should be as short as possible as otherwise there will be insufficient moisture for staining.

### Evaluation

After a staining time of 15–20 minutes cell structures and crystals can be evaluated.

### Interpretation

#### Normal result

A healthy nasal smear contains epithelial cells with a relatively small amount of cytoplasm and large nucleus and otherwise mainly ciliated epithelial cells and goblet cells.

#### Neutrophilic granulocytes

Microscopic detection of increased numbers of neutrophils is proof of rhinitis which can be of bacterial, viral or allergic superinfected origin.

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### **Eosinophilic granulocytes**

The presence of 'nests' of eosinophils or more than 30% eosinophilic granulocytes indicates an anaphylactic allergy or analgesic intolerance.

### **Basophilic cells**

Cells containing basophilic granules are less abundant than eosinophils. The presence of basophilic cells, either basophilic granulocytes or mast cells, suggests rhinitis of allergic origin. Basophilic cells also occur in obstructed nasal ventilation, e.g. due to tumours or polyps, and in patients who have undergone tracheotomy or laryngectomy.

### **Lymphocytes**

In acute viral infections there are increased numbers of lymphocytes.

### **Giant cells and ciliary loss**

Multinuclear giant cells indicate an acute viral infection. Degenerative changes in the ciliated cells can lead to ciliary loss.

### **Epidermoid cells**

Small nucleus, large cytoplasm, typically found in an atrophic rhinitis.


### **Crystals**

Breakdown of eosinophils produces the spindle-shaped Charcot-Leyden crystals which are found in allergic rhinitis and analgesic intolerance. Fern-like crystals are found in nasal secretions with an increased sodium-potassium content accompanied by reduced colloid levels. This reduction of the protective film is a sign of beginning mucosal atrophy in the late stage of chronic rhinitis.

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In routine semen analysis Testsimplets® can be used in place of the time-consuming and costly Papanicolaou staining. Testsimplets® is a labour-saving method for staining the fresh ejaculate.

### Procedure

The amount of liquefied ejaculate needed is dependent on the spermatozoa concentration:

- 5 µl (spermatozoa concentration > 50 mill/ml)
- 10 µl (spermatozoa concentration > 10–50 mill/ml)
- 15 µl (spermatozoa concentration < 10 mill/ml)

Apply a drop of the liquefied ejaculate to the middle of a cover slip and place on the prestained area of the slide in such a way that the ejaculate spreads evenly in all directions. To ensure even distribution of the sample it is useful to press on the middle of the cover slip with a pencil.

After 30 minutes a weak staining is attainable, however the preparation should be left to stand at room temperature for about 2 hours until the majority of the spermatozoa have become immobile and well stained.

Then differentiation into normal and pathological forms is undertaken. The preparation is stable for at least 4 hours at room temperature. If microscopic differentiation cannot be performed until after 4 to 24 hours the slide should be refrigerated immediately after preparation. After 20–24 hours stored in a refrigerator at 4 °C the spermatozoa are still well stained.

Fixation using Eukitt® is not advisable due to the subsequent destruction of the spermatozoa.

### Evaluation


Evaluation is performed at a magnification of x800 to x1000 with oil immersion. Usually 200–500 spermatozoa per preparation are evaluated and documented in percent according to the criteria for normal and pathological forms.

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
*Human spermatozoa with high proportion of morphologically normal spermatozoa*

## Characterization of the spermatozoa

The spermatozoa stain very well. The head, the mid-piece and the tail can all be clearly identified. The nucleus lying within the head stains dark purple. The plasma appears orange to pale violet depending on the focus. The mid-piece also appears pale violet if it is pathologically enlarged. The tail is only lightly stained but is clearly visible.

## Normal findings

More than 60% morphologically normal spermatozoa.

Normal geformtes Spermatozoon

- a) Akrosom
- b) Kern
- c) Mittelteil
- d) Schwanz

*Morphologically normal spermatozoon*

- a) acrosome
- b) nucleus
- c) mid-piece
- d) tail

Art. Nr.	Bezeichnung	Packungsinhalt
191574	Testsimplets®	50 gebrauchsfertige, farbbeschichtete Objektträger 50 staubfreie Deckgläser (24 x 36 mm)
Cat. No.	Name	Pack content
191574	Testsimplets®	50 ready-to-use, prestained slides 50 dust-free cover slips (24 x 36 mm)

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*Differential diagnostic  
microscopy with  
**Testsimplets®** –  
well tried in  
practice and laboratory*

Waldeck GmbH & Co. KG  
Havixbecker Str. 62 · D-48161 Münster  
Tel.: +49 (0) 25 34/970-0  
Fax: +49 (0) 25 34/970-258  
E-Mail: [info@waldeck-online.com](mailto:info@waldeck-online.com)  
Internet: [www.waldeck-online.com](http://www.waldeck-online.com)

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Diagonal GmbH & Co. KG