



Techline Technologies, Inc.
Device Specifications PTT-B54
Injection in a Box™



PTT-B54

Injection in a Box™

*This product is covered by one or more of the following patents:
US10,217,380, US10,535,282, US10,586,470, US11,195,434,
US11,195,435, and/or other(s) pending.*

TECHLINE® Device Specifications
March 11, 2022

TECHLINETRAUMA.COM

Description-

The Injection in a Box™ task trainer is for training students in intramuscular injection. The simulated SOFTTECH® skin provides realistic feedback to the trainee, while the injected fluid is collected by the removable sponge, located in the bottom of the case.

Construction-

The representative skin and tissue of the PTT-B54 is constructed of high-quality silicone rubber, with layers of coloration that provide a depth and realism to the product. In addition, the SOFTTECH® material has a realistic feel, and is very durable and UV resistant.

Operation-

The task trainer is contained in its own storage case. The Injection in a Box™ can be used for administration of intermuscular injections or inject fluid, to simulate IV administration. The moisture collecting sponge allows for the collection of the injected fluid in the cavity. Can be used inside or outside of the case.

Maintenance-

The device should be rinsed in clean, warm water and allowed to dry before storage. The simulation may be removed from the storage case for cleaning. The device may be washed with mild detergents like dish soap or cleaned with a disinfecting agent like isopropyl alcohol or a diluted bleach solution, as it is durable, and chemically resistant. Allow the device and sponge to dry fully before storage to prevent mold or mildew growth.



PTT-B54 Injection in a Box™
Light, Medium, and Dark Skin Tones Available