

# journey SoftSecure

ROTATING SHOWER CHAIR

## PRODUCT MANUAL



ROTATING SHOWER CHAIR

 journey  
HEALTH & LIFESTYLE®

## Welcome to Journey Health & Lifestyle!

Thank you for purchasing your SoftSecure Rotating Shower Chair. You've invested in a high-quality, premium materials product that will facilitate comfort and safety.

### **Safety instructions**

- Read these instructions for use carefully.
- All components should be checked for damage and a secure fit prior to use.
- Approved for indoor use only.
- Clothing or body parts can be pinched when sitting down or standing up, as well as when raising the armrests and locking them into place.
- Be careful on the seat; be aware of wet and slippery surfaces.
- Prevent unauthorized use, for example, by children.
- Take note of the maximum weight permitted.

## II. PRODUCT DESCRIPTION

### Parts of the product



1. Backrest
2. Seat
3. Flip up armrest
4. Rotating plate with knob
5. Push button

### III. INSTRUCTIONS FOR USE

**Step 1-**

Attach the legs to the frame by depressing push buttons into the tubes located on each side of the legs.

Make sure that push buttons “snap” correctly through holes and are stable and adjusted to the same height.



**Step 2-**

Attach backrest to the frame, by depressing push buttons and sliding the backrest into the tubes located on each side and on the back of the seat.



**Step 3-**

Check that the seat is stable and level prior to use.

- (1) Pull the knob and turn the seat to face the user, then lock the seat into position by locking the knob and keep swiveling the seat until it locks in position.

### III. INSTRUCTIONS FOR USE



- (2) Unlock the seat by pulling the knob.
- (3) When user is seated on the chair, position the seat to ensure that it is locked into position and no longer rotates freely.
- (4) Lock the shower chair when the user faces the open side of the bath. The user can now stand up carefully and slowly.
- (5) To flip up the armrest, push the button shown below.



<b>Product Specifications</b>	
<b>Product Dimension</b>	(18"-19") x (18"-19") x (30"-34")
<b>Weight Capacity</b>	300 lbs
<b>Packaging Dimensions</b>	25" x 18" x 16"
<b>Net Weight</b>	21 lbs
<b>Material on product</b>	Aluminum

## **Warranty**

Journey Health & Lifestyle warrants the SoftSecure Rotating Shower Chair frame to be free from defects in materials, workmanship and assembly for twelve (12) months from the original purchase date. The warranty does not extend to non-durable components such as the rubber tips.



 **journey**  
HEALTH & LIFESTYLE®

 **journey**  
**SoftSecure**

ROTATING SHOWER CHAIR

If you have any questions, call us at our toll free number:

**1-800-958-8324**

or visit [www.shopjourney.com](http://www.shopjourney.com)