

INFLATION INSTRUCTIONS

FitPAWS® Balance Products are designed for canine rehabilitation and conditioning exercises but they are not puncture-proof and like any other plastic inflatable under pressure, when pushed past their limits, can burst. If you choose your exercises wisely and follow a few simple guidelines, the FitPAWS® Balance Products are designed to function properly.

PLEASE READ CAREFULLY

1. Always inspect your FitPAWS® products for damage. Creases and folds are normal, as it will come folded, but any snags or tears will prohibit you from inflating properly.
2. Allow some time for the inflatable products to get to room temperature because if the material is too hot or too cold it will affect how well it inflates. We strongly recommend inflating and letting the product sit for 24 hours to stretch and adapt to the environments temperature. You may find that after 24 hours you can inflate slightly more.
3. **DO NOT OVER INFLATE.** FitPAWS® Balance Products should be filled according to size, but there can be variances.

For best performance, you should be able to depress the surface easily, approximately 2" with your fingertips, soft to the touch. Any inflatable that is stretched too tight has the potential to burst. Please pay attention to the surface

tension of the products. If it is too tight to the touch, reduce the air and let the product sit for 24 hours before using.

Maximum sizes noted are ONLY a guideline. If the products, for any reason, appear to be too tight and does not inflate to maximum height, or is hard to the touch, then STOP - DO NOT inflate further. Release air and allow the product to sit overnight before use.

If you have questions, please contact us at info@fitpawsusa.com.

Use tape or a pencil, mark the maximum height printed on your peanut or egg on a wall so that you can measure the ball during inflation. Measure from the highest point of the ball, not the saddle. Use the enclosed pump or air compressor inflate ball up to 80-85% of the maximum height. Plug the ball and let rest for 24 hours before inflating to desired height.

Product	Inches*	CM*	85% Inflation (Inflate Wait 24 Hrs)	95% Inflation (Recommended Max Inflation)
80 cm Peanut	31.5" H, 62" L	80 cm H, 157 cm L	27" or 68 cm	30" or 76 cm
70 cm Peanut	27.5" H, 50.75" L	70 cm H, 128 cm L	24" or 60 cm	27" or 67 cm
60 cm Peanut	23.5" H, 44" L	60 cm H, 111 cm L	20" or 51 cm	23" or 57 cm
50 cm Peanut	19.5" H, 36.25" L	50 cm H, 92 cm L	17" or 43 cm	19" or 48 cm
40 cm Peanut	15.5" H, 32" L	40 cm H, 81 cm L	13" or 34 cm	15" or 38 cm

*Do Not exceed maximum diameter printed on the peanut.
You can always use a ball, less inflated, depending on the type of balance challenge desired.