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

Mölnlycke®
Mepilex® Border Ag
WITH SAFETAC® TECHNOLOGY


STERILE EO Sterilized using ethylene oxide


CE 2797 Marking CE accordance with MDD 93/42/EEC; 2797 is the number on the body of the Mepilex® Border Ag product

Manufacturer
 Mölnlycke Health Care AB
 Gamlestadsvägen 3C, Box 13080, SE-402 52 Göteborg, Sweden

Mepilex® Border Ag


Mepilex® Border Sacrum Ag

- Single use
- See instructions for use
- Keep dry
- Keep away from sunlight
- 25°C / 77°F
Upper temperature limit
- Do not use if package is damaged

- For low exuding wounds
- For moderately exuding wounds
- For highly exuding wounds

Mepilex® Border Ag with Safetac® Technology
 Self-adherent antimicrobial soft silicone foam dressing.

Product description
 Mepilex® Border Ag consists of a Safetac® soft silicone wound contact layer, an absorbent polyurethane foam pad containing silver sulphate and activated carbon, a layer with super absorbent polyacrylate fibres, a non woven and a vapour permeable and waterproof film.

Mepilex® Border Ag is a soft silicone foam dressing that absorbs exudate, maintains a moist wound healing environment and has antimicrobial properties.

Mepilex® Border Ag contains silver sulphate which when in contact with fluid will release silver ions to create an effective bacterial barrier and inactivate a wide range of wound related pathogens (bacteria and fungi), as shown *in vitro*. By reducing the number of microorganisms, Mepilex® Border Ag may also reduce odour.

Mepilex® Border Ag has been shown to inactivate wound related pathogens up to 7 days *in vitro*.

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Safetac® is a unique adhesive technology that minimises pain to patients and trauma to intact skin or wounds.

Dressing material content :

Polyurethane, silicone, polyacrylate, cotton, viscose, polyester, polyolefin, silver sulphate.

Intended use

Mepilex® Border Ag is designed for the management of medium to high exuding wounds such as leg and foot ulcers, pressure ulcers, malignant wounds, partial thickness burns, traumatic and surgical wounds where a moist environment, exudate handling, gentle fixation and an antimicrobial action is indicated.

Mepilex® Border Ag may be used on infected wounds as part of a treatment regimen under supervision of a qualified health care professional.

Mepilex® Border Ag can be used under compression bandaging.

Instructions for use

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Remove the release films and apply the adherent side to the wound. Do not stretch.
4. For best results, **Mepilex®** Border Ag's wound pad should overlap the wound edges by at least 2 cm.

Mepilex® Border Ag is intended for short-term use up to 4 weeks. For continued use, a re-assessment by a physician is recommended.

Frequency of change

Mepilex® Border Ag may be left in place for up to seven days depending on the patient, condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

Initially it might be necessary to change **Mepilex®** Border Ag more frequently. This is due to a change with the treatment regime, which can often result in an initial increase in exudation.

Precautions

- **Mepilex®** Border Ag should be used under the supervision a qualified health care professional.
- Do not use on patients with a known sensitivity to silver or any other content of the dressing.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- **Mepilex®** Border Ag may cause transient discoloration of the wound bed and surrounding skin.
- In the event of clinical infection **Mepilex®** Border Ag does not replace the need for systemic therapy or other adequate infection treatment.
- Do not use **Mepilex®** Border Ag during radiation treatment or examinations e.g X-ray, ultrasound, diathermy or Magnetic Resonance Imaging.
- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).
- Do not use **Mepilex®** Border Ag together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- Other than saline solution or water, the interaction of cleansing agents in combination with **Mepilex®** Border Ag has not been demonstrated.
- The interaction of **Mepilex®** Border Ag with topical treatments has not been demonstrated.
- For external use only.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.
- If the product is used after the expiry date product properties cannot be ensured.

Storage and disposal

- **Mepilex®** Border Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight.
- Disposal should be handled according to local environmental procedures.

Presentations

Mepilex® Border Ag 7,5 x 7,5 cm	Box, 5 pcs
Mepilex® Border Ag 10 x 10 cm	Box, 5 pcs
Mepilex® Border Ag 10 x 20 cm	Box, 5 pcs
Mepilex® Border Ag 10 x 25 cm	Box, 5 pcs
Mepilex® Border Ag 10 x 30 cm	Box, 5 pcs
Mepilex® Border Ag 12,5 x 12,5 cm	Box, 5 pcs
Mepilex® Border Ag 15 x 15 cm	Box, 5 pcs
Mepilex® Border Ag 15 x 20 cm	Box, 5 pcs
Mepilex® Border Ag 18 x 18 cm (Sacrum)	Box, 5 pcs
Mepilex® Border Ag 20 x 20 cm (Sacrum)	Box, 5 pcs
Mepilex® Border Ag 23 x 23 cm (Sacrum)	Box, 5 pcs

KEMENKES RI AKL

Mepilex® and **Safetac®** are registered trademarks of Mölnlycke Health Care AB

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