

SYS-6010 Infusion Pump

Service Manual



Version 2.0 English

Intellectual Property and Statement

Intellectual Property Rights

The intellectual property right of this product and its Service Manual belongs to SHENZHEN MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD. (hereinafter short as MEDCAPTAIN).

©2014-2015 All rights reserved. SHENZHEN MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.

Without prior approval from MEDCAPTAIN in writing, this Service Manual shall not be photocopied, modified or translated, fully or partially, by any individual or organization.  , MEDCAPTAIN and  are registered trademarks or trademarks of MEDCAPTAIN.

Statements

MEDCAPTAIN reserves the right for final interpretation of this Service Manual.

MEDCAPTAIN reserves the right to modify the contents of this service manual for a more accurate and effective service quality. The modified contents should be reflected in the newly published service manual version.

MEDCAPTAIN is responsible for safety, reliability and performance of this equipment only on the condition that:

- All disassembly, replacement, test, modification and repair are conducted by qualified personnel approved by MEDCAPTAIN.
- All replacement parts, supporting accessories and consumables during the maintenance are provided by MEDCAPTAIN;
- Maintenance records for product are reserved.

Version Information

content	language	version	Publish time	code
SYS-6010 Infusion Pump Service Manual	English	V2.0	Sep,2016	

Contents

1	IMPORTANT INFORMATION	1
1.1	SERVICE PERSONNEL	1
1.2	THE LIMITATIONS OF THE SERVICE MANUAL	1
1.3	DEVICE MAINTENANCE AND INSPECTION	1
1.4	QUALITY CONTROL	2
1.5	SAFETY INSTRUCTION	2
1.6	LIST OF ABBREVIATED	2
1.7	LIST OF SYMBOLS	3
1.8	CONTACTS	3
2	SYSTEM INTRODUCTION	4
2.1	INTRODUCTION	4
2.2	THE COMPOSITION OF THE INFUSION PUMP	4
2.3	APPEARANCE	5
2.4	CIRCUIT BLOCK DIAGRAM	8
2.5	THE SOFTWARE	19
2.6	SPECIFICATIONS	19
3	MAINTENANCE	21
3.1	CLEANING AND STERILIZING	21
3.2	PERIODIC MAINTENANCE	21
3.3	REPLACING THE BATTERY	23
4	SOFTWARE UPGRADE	24
4.1	HOW TO UPGRADE THE DRIVER SOFTWARE	24
4.2	UI SOFTWARE UPGRADE	28
5	ALARM, FAULT AND TROUBLESHOOTING	30
5.1	ALARM INFORMATION	31
5.2	COMMON FAULTS AND TROUBLESHOOTING	33
6	THE DISASSEMBLING AND INSTALLING	34
6.1	ATTENTIONS	34
6.2	ASSEMBLY TOOLS	35
6.3	DISASSEMBLY	35

6.4	INSTALLATION	43
7	THE SAFETY TEST	43
7.1	ELECTRIC LEAKAGE TEST OF ENCLOSURE	43
7.2	ELECTRIC LEAKAGE TEST OF GROUND	44
7.3	ELECTRIC LEAKAGE TEST OF PATIENT.....	44
8	MAINTENANCE SPARE PARTS LIST	45

Important information

1 Important Information

1.1 Service personnel

This service manual is for the information of service personnel only. Service may only be conducted by personnel who

- Has basic knowledge of electronic circuit and mechanical.
- Has basic knowledge of medical devices and clinic.
- Has received proper training of maintenance service of MEDCAPTAIN and has certain knowledge of device principle, structure, performance and operation.
- Has the necessary equipment and instruments.
- The permission of MEDCAPTAIN.

1.2 The Limitations of the Service Manual

This service manual describes all performance and configurations of the device. Some other devices performance and configurations are not included here. There are differences between the published manual and the actual status of the device, especially when device has been modified, a corresponding maintenance information is needed. Therefore, this service manual may be used together with some follow-up complementary information.

MEDCAPTAIN will complement the relevant information timely according to the actual situation of the device modification.

Service personnel should also refer to the Operation Manual.

1.3 Device Maintenance and Inspection

- The power supply of this device can reach to 100-240V. Maintenance against the requirements of the service manual may cause an electric shock hurt or death.
- Maintenance against these requirements of this service manual may seriously damage the device.
- The service personnel must be trained and permitted by MEDCAPTAIN.
- The maintenance must be of electrostatic discharge (ESD) protection conditions. Do not touch PCBA or semiconductors by hand without any protective measure.
- Do not operate on the display by using sharp objects. Otherwise, the display may be damaged.
- Do not disinfect the infusion pump by using the high-pressure steam sterilization method.
- Before internal battery operation, check the battery to ensure that sufficient power is available. Recharge, if required.
- Liquid intrusion into the AC power socket, USB or nurse call socket may cause short-circuit. While connecting the power cable, check if the connecting parts are dry. If liquid spills on the infusion pump, clean the pump with a dry cloth. And then check.

- Use the maintenance parts and accessories provided by MEDCAPTAIN to replace and maintain the device.
- Do not use the infusion pump in a flammable environment.
- High-frequency surgical equipment, mobile phone, wireless device and defibrillator may have interference on the infusion pump. Keep away from them while operating.
- After maintenance, have a safety test and cleaning according to this manual.

1.4 Quality control

MEDCAPTAIN satisfies ISO9001 and ISO13485 Certificate of Quality System. The relative products have also obtained the CE Certification, in line with the requirements of MMD instruction.

1.5 Safety instruction

Warnings:

- The power supply of this device can reach to 100-240V. Maintenance against the requirements of this Service Manual may cause an electric shock hurt or death.
- Maintenance against the requirements of this Service Manual may seriously damage the device.

1.6 List of abbreviated

ESD	Electro-Static discharge
PCBA	Printed Circuit Board Assembly
ISO	International Standardization Organization
CE	Council of Europe
MDD	Medical Device Directive
IPX2	Level of Protection from Liquid instruction
ON/OFF	ON/OFF key
HOME	Main Menu Key
CPU	Central Processing Unit
N/A (Not Applying)	Not Applying

1.7 List of symbols

Symbol	Description	Symbol	Description
	CAUTION! Read the accompanying document.		Type CF equipment
	Alternating current	IPX2	Level of protection from liquid instruction
	Direct current		ON/OFF
	Manufacturer		Date of manufacture
	HOME		

1.8 Contacts

- a) When the infusion pump is used and maintained normally and regularly, we provide 18 months free maintenance (replacement of battery and appearance parts are not included). For failure caused by other factors, part replacement will be charged accordingly.
- b) If you have any questions when using the infusion pump, please contact local distributor or directly contact us at any time.

The after-sales service contact details of Medcaptain Medical Technology Co., Ltd. are as follows:

Address: 12th Floor, Baiwang Research Building, No.5158 Shahe West Road, Xili, Nanshan District, Shenzhen, P.R.China

Telephone: 400-809-3369

Fax: 0755-26001651

Postal: 518055

Website:<http://www.medcaptain.com>

E-mail:MC.Service@medcaptain.com

System Introduction

2 System Introduction

2.1 Introduction

2.1.1 Purpose

This product is intended for hospitals to infuse liquid at constant speed or liquid medicine continuously through the veins of patients.

2.1.2 Product Features

MEDCAPTAIN SYS-6010 is a micro-continuous infusion pump. It ensures constant infusion speed and accurate dosing volume during long time infusion.

- Support all disposable IV sets conform to the standard.
- You can customize other IV sets conform to the standard.
- Three occlusion levels can be selected from three levels. The pressure status of tubing also display.
- The maximum infusion rate can be set to 1200mL /h
- Calibration function for infusion accuracy
- Safety design by monitoring infusion status.
- Multiple modes of infusion
- Multi-channel IV delivery system, realizing relay infusion function
- Touch screen, providing quick and convenient man-machine interface.
- Display night mode, reducing light interference to patients and environment.
- Connection to bar code scanner function.
- Providing three types of power supply: AC power supply, DC power supply, and internal lithium battery. The lithium battery can power the infusion pump for no less than 5 hours (when infusion is run at 25mL/h).
- Double CPU and redundancy design for key units.
- Two-way alarm for monitoring the main control circuit and motor drive circuit at any time.
- Independent motor driving CPU and motor subdivided drive chip design.
- Setting and automatic prompt of maintenance interval.
- Modular installation design enables multi-channel pumps among pumps.
- Handle, pole clamp, bar code scanner, relay infusion function are optional, depending on the user's need.

2.2 The Composition of the Infusion Pump

The SYS-6010 infusion pump mainly consists of the pump shell, display and operating system, monitoring system, alarm system, motor drive system, tubing peristaltic module, power supply system, drop sensor, the handle (optional), and pole clamp (optional).

The infusion pump adopts dual processor structure, controls the motor precisely,

drives the peristaltic sheet to infuse through the mechanical drive device, monitors the sensors and infusion process, and provides sound and light alarms.

2.3 Appearance




Figure 1

1 – Touchscreen

2 – [HOME] key

3 – [ON/OFF] key


4 – [OPEN] key

5 – Pump door

6 – Shell

7 – Alarm indicator

8 – Place shelf of drop sensor



- | | | |
|---------------------|----------------------------|--------------------------------|
| 1 –IV groove | 2 –Peristaltic Pump Tablet | 3 –Lighting lamp |
| 4 –Depressor | 5 –Door catch | 6 –Air bubble sensor |
| 7 – Pressure sensor | 8 –Anti-Free-Flow clamp | 9 –Anti-Free-Flow Clamp Button |





Figure 2

- | | | |
|---------------------|---------------------|--------------------------|
| 1 –Cover of Battery | 2 –Pole clamp Hole | 3 – Loud speaker |
| 4 - Auxiliary Alarm | 5-external inlet 1 | 6 –AC power supply Inlet |
| 7 –External Inlet 2 | 8 –External Inlet 3 | 9 – Auxiliary Alarm |
| 10 – Shell | | |

2.3.1 Drop Sensor



1 – Button

2 – Drop Hole

3 – Cable

4 – Socket

2.3.2 Pole Clamp




Figure 4

1 – Mounting Screw

2 - Mounting Knob of Infusion Stand

2.3.3 Nurse Pager




Figure 5

1 – Button

2 – Cable

3 – Socket

2.4 Circuit Block Diagram

The whole circuit block diagram of SYS-6010 is shown as Figure 6. The parts are as follows:

V2.0




Figure 6

2.4.1 AC/DC Power Board

a. Circuit Introduction:

The AC/DC power board has adopted the medical level power supply with low leakage current and high-insulation voltage. It satisfies the requirements of EMC. The AC/DC power board has an input alternating voltage ranging from 90 to 240V, 50/60Hz, direct output voltage of 13.6V and the maximum output current of 2A. The AC/DC power board mainly consists of AC inlet, fuse, EMC filter circuit, rectifier, high voltage DC filter capacitor, switching tube, pulse width modulator, transformer, rectifier and absorbing circuit, DC filter circuit, stabilized light decoupling feedback control circuit, etc.

b. Main Testing Point

No.	Testing point	Description	Range	Unit	Remark
1	TP9	DC Output Voltage	13.5-13.7	V	

c. Circuit Socket Number and Definition

No.	Socket No.	Description	Pin-out Description	Range	Explanation
1	J1	DC output voltage inlet	1: DC output	13.0~13.7V	Connect J1 in the DC/DC Power board
			2: Voltage sampling	2.78~2.94V	
			3: Ground	Output DC ground	
2	J2	AC input inlet	1: Live line	AC L line input	Connect AC power input inlet
			2: NC	/	
			3: Neutral line	AC N line input	
3	J3	Ground connecting inlet	1: Ground	/	Connect AC power input inlet

d. Caution:

- The AC/DC power board can reach to a high voltage of 400V. The workbench is required to be isolated and avoid touching the high voltage circuit parts. High voltage may cause electric shock hurts.
- While using a multimeter or oscilloscope to test, the pens cannot short out the live parts in the circuit board, or the circuit board would be damaged.
- While using an oscilloscope or multimeter powered by mains to test the circuit of mains in the AC/DC power board (primary non-isolated side circuit), the isolated differential voltage probe must be used to measure, or it will damage the board and test instrument.
- While using an oscilloscope multi-channel probe to test the circuit, all ground loops of the probe in the oscilloscope must be connected to the same potential point, or it will damage the board and test instrument.

2.4.2 The DC/DC Power Board

a. The Circuit Introduction

The DC/DC power board is from the AC/DC power board with main functions as follows:

(1) Management circuit of battery charge: the charging current of lithium battery is limited at around 200mA. As charging voltage boots, charging current gradually decreases. A fully charge for 1500mA battery generally needs 8-10 hours.

(2) 18V stepper motor driving power supply: 18V DC power supply is provided to the stepping motor. The maximum current is 600mA.

(3) 5V DC power supply: provides power supply to the digital circuit, CPU, the sensor circuit. A maximum input current is 3A.

(4) 3.7V maintaining power supply: provides maintaining power supply to CPU and

clock circuit in the shutdown situation. The input is from battery.

(5)The stepper motor driving circuit: adopts the specialized stepper motor controlling circuit. Under the control of CPU, the speed, subdivision, current ratio and damping ratio can be adjusted, depending on requirements.

b. The Main Testing Points

No.	Testing point	Description	Range	Unit	Remark
1	TP5	4.0V DC output voltage	4.0-5.0	V	
2	TP6	5V DC output voltage	4.9-5.2	V	
3	TP7	18V DC output voltage	17.5-18.5	V	

c. Number and Definition of the Circuit Board Sockets

No.	Socket	PIN No.	Description	Explanation
1	J1	3	Socket to inlet the DC output voltage	Connect J1 in the AC/DC board
2	J2	14	Socket to connect the power board	Connect J5 in the power board
3	J8	14	Socket to connect the power board	Connect J8 in the power board
4	CON2	8	Socket to connect the interface board	Connect CON2 in the interface board

2.4.3 The Power Board

a. The Circuit Introduction

The power control board has following main functions:

(1) Control the CPU circuit: controls on/off and stepper motor, receives and processes signals from the compression bar sensors, and communicates with interface CPU of SYS02 circuit.

(2)A function of signal switch board: collects signals from every sensor and all kinds of DC power, distributes the corresponding signal to SYS02 circuit board, DC/DC power board, stepper motor, alarm light and speaker, RS485 serial port and external interface, etc.

(3)Charging circuit and boosted circuit of the super-capacitor.

(4)MCU control circuit of buzzer alarm.

b. The Main Testing Points: none.

c. Number and Definition of the Circuit Board Sockets

No.	Socket	PIN No.	Description	Explanation
1	J2	24	Socket to connect the	Connect J2 in the SYS02 board

			SYS02 board	
2	J3	8	Socket to connect the pressure&speed detection board of the infusion pump	Connect J2 in the pressure &speed detection board of the infusion pump
3	J4	7	Socket to connect the pressure&speed detection board of the infusion pump	Connect J1 in the pressure &speed detection board of the infusion pump
4	J5	14	Socket to connect the DC/DC board	Connect J2 in the DC/DC board
5	J7	4	Socket to connect the pressure pressure&speed detection board of the infusion pump	Connect J7 in the pressure&speed detection board of the infusion pump
6	J8	14	Socket to connect the DC/DC board	Connect J8 in the DC/DC board
7	J9	3	Not used	Apply in the infusion pump
8	J10	3	Socket to connect the bubble sensor	Connect bubble sensor
9	J16	3	Socket to connect the power supply	Connect the battery pack plug
10	CON1	8	Socket to connect the interface board	Connect in the CON1 interface board

2.4.4 Battery Pack

a. Introduction

SYS-6010 adopts the 11.1V/1500mA lithium-ion battery pack. Connect battery pack to the power board with a connector. Disconnecting the connector, you can separate the battery pack and connecting wire.

b. The Main Testing Points

No.	Testing point	Description	Range	Unit	Remark
1	Black-black wire	Short circuit test	The two wires should be shorted		
2	Red-black wire	Battery pack voltage	9-12.5	V	Battery level influence the voltage

Note: The quality of battery pack is not judged by its voltage. Connecting the infusion pump to test the actual operating time is needed.

c. Number and Definition of the Sockets

No.	Socket	PIN	Description	Explanation
-----	--------	-----	-------------	-------------

		No.		
1	/	3	Socket to connect the power board	Connect J16 in the power board

2.4.5 Interface Board

a. Circuit Introduction

The interface circuit has following main functions:

(1) Provides the interface to connect SYS-6010 to peripherals, including the switching of power supply, correspondence, drop sensor signal, bar code scanner and the external power supply.

(2) Provides the speaker and buzzer driving circuit.

(3) Transfers the program recording signal of the power board CPU.

b. The Main Testing Points:none

c. Number and Definition of the Circuit Board Sockets

No.	Socket	PIN No.	Description	Explanation
1	J2	4	Socket to connect speaker and buzzer	Connect speaker and buzzer
2	CON1	8	Socket to connect power board	Connect CON1 in the power board
3	CON2	8	Socket to connect DC/DC power board	Connect CON2 in the DC/DC power board
4	BUS1	5	USB2.0 socket	The interface to connect peripherals including power input/output, correspondence and drop sensor; Connect MP-80 workstation, or individually connect the external power supply and drop sensor; Program update of CPU in the SYS02 board.
5	BUS2	10	USB3.0 socket	The interface to connect peripherals including power input/output, correspondence, drop sensor, and nurse call. Program update of CPU in the SYS02 board.
6	BUS3	10	USB3.0 socket	The interface to connect peripherals including power input/output, correspondence, drop sensor, and nurse call. Program update of CPU in the power board and SYS02 board.


No.	Socket	PIN No.	Description	Explanation
1	J2	4	Socket to connect speaker and buzzer	Connect speaker and buzzer
2	CON1	8	Socket to connect power board	Connect CON1 in the power board
3	CON2	8	Socket to connect DC/DC power board	Connect CON2 in the DC/DC power board

2.4.6 Interface Socket

a. Socket Introduction


(1)USB1 socket (BUS1): USB2.0 socket interface. To connect the external DC power supply, bar code scanner of RS485 serial port, power supply of the intravenous workstation with RS485 correspondence, drop sensor signal, program update of SYS02 board CPU. The external power supply should be used with the matched cables. The DC input power is 9-15V and the current is 1A.

The definition of socket signal



(2)USB2 socket (BUS2): USB3.0 socket interface. To connect the external DC power supply, bar code scanner of RS485 serial port, drop sensor signal, nurse call, and program update of SYS02 board CPU. The external power supply should be used with the matched cables. The DC input power is 9-15V and the current is 1A.

The Definition of Socket Signal




Remark: When the voltage of external power supply connecting to the USB1 and USB2 socket is over 13.6V, the external battery in SYS-6010 can be fully charged. If the voltage


is under 13.6V, the battery may not be fully charged or not charged.

(3)USB3 socket (BUS3):USB3.0 socket interface. To connect the external 5V DC power supply, bar code scanner of RS485 serial port, drop sensor signal, nurse call, program update of the power board and SYS02 board CPU.

The Definition of Socket Signal



b. The Main Testing Points: none



2.4.7 SYS02 Board

a. Circuit Introduction

The SYS02 board has follows main functions:

(1)The man-machine interactive operation and management: the touch screen and keyboard input processing; interface and menu display; connecting to the power board serial ports, to obtain the infusion pump parameters, alarm information, the sensor signal display and alarm output; output the infusion parameters set by the user to power board to make a driving operation.

(3)2.4 inch TFT color display: with a resistance to touch. The screen resolution is 320X240 with FPC welded at the circuit board.

b. The main testing points: none

c. Number and definition of the sockets

No.	Socket	PIN No.	Description	Explanation
1	J1	24	Socket to connect the alarm light board	Connect J1 in the alarm light board
2	J2	24	Socket to connect the power board	Connect J2 in the power board
3	J3	8	Socket to connect the button cover	Connect the button cover

2.4.8 Button cover

a. Button cover introduction

(1) Buttons: [ON/OFF] for switching on/off,[HOME] for returning to main screen,[OPEN] for opening the pump door.

(2) Indicator light: 2 LED lights of green and yellow are built in the cover. While the machine is power off and AC inputting, yellow light works. While the power is on, green light works. The indicator lights works as follows:

State of SYS-6010	State of yellow light	State of green light
Power off and no AC	OFF	OFF
Power off and AC	ON	OFF
Power on	OFF	ON

Remark: "ON" means light is on, "OFF" means light is out.

b. The main testing points:none

c. Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	FPC block	8	Socket to connect SYS02 board	Connect J3 in the SYS02 board

2.4.9 Alarm light board

a. The circuit introduction

3 groups of LED lights are installed in the alarm light board. The colors are green, yellow and red, each group of 2 LED lights, to indicate the states of working, high, middle and low level alarm. The states of indicators are as follows:

State of indicators	Meaning of lights indicate
Green light flashing	The infusion pump is normally working, flash frequency indicates infusion speed
Yellow light illuminating	Lower level alarm
Yellow light flashing	Middle level alarm

Red light flashing	High level alarm
--------------------	------------------

- b. The main testing points: none
- c. Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	J1	24	Socket to connect SYS02 board	Connect J1 in the SYS02 board

2.4.10 The pressure& speed detection board

- a. The circuit introduction

Main functions are as follows:

(1)Process the pressure sensor signals: amplify and filter signals from the pressure sensor, output to the power board.

(2)Collect and transfer signals: the peristaltic pump tablets signal, the stepper motor speed signal, the pump door-open signal, the pump door-off signal and the hall switch signal are all collected to the pressure& speed detection board through two cables to the power board.

(3)Switch the driving signal: switch the driving signals from the door-open motor and infusion motor through controlling the relay.

Note: the pressure sensor and the infusion pump pressure& speed detection board should be adjusted and tested together, therefore, they need to be replaced one to one correspondence, or there be a new debugging.

- b. The main testing points: none
- c. Number and definition of the sockets

No.	Socket	PIN No.	Description	Explanation
1	J1	7	Socket to connect the power board	Connect J4 in the power board
2	J2	8	Socket to connect the power board	Connect J3 in the power board
3	J3	4	Socket to connect the pressure sensor	Connect the pressure sensor
4	J4	4	Socket to connect the stepper motor speed board	Connect the stepper motor speed board
5	J5	3	Socket to connect the pump tablets testing board	Connect the pump tablets testing board
6	J6	6	Socket to connect the hall switch board	Connect the hall switch board
7	J7	4	Socket to connect the power board	Connect J7 in the power board
8	J8	4	Socket to connect the	Connect the open motor

			open motor board	
9	J9	4	Socket to connect the infusion motor board	Connect the infusion motor

2.4.11 The infusion motor

- The motor introduction: adopts a stepper motor with 42mm diameter and 48 steps, 2 phases driving. The working voltage is 18V.
- Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	/	4	Socket to connect the infusion pump pressure & speed detection board	Connect J8 in the infusion pump pressure & speed detection board

2.4.12 Stepper motor speed board

- The circuit introduction
Motor speed function: two light couplings collect the guiding and shading light signals A and B from the optical disc connecting on the stepper motor. You can judge the speed and rotating direction of the stepper motor according to the order of A and B and if there are A and B.
- The testing points: none
- Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	J1	4	Socket to connect the infusion pump pressure & speed detection board	Connect J4 in the infusion pump pressure & speed detection board

2.4.13 The pump tablets position testing board

- The circuit introduction
The peristaltic pump position testing function: the light coupling collects the guiding light signals from the optical disc connecting on the peristaltic pump and judge the position of the peristaltic pump during the rotation.
- The main testing points: none.
- Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	J1	3	Socket to connect the infusion pump pressure & speed detection board	Connect J5 in the infusion pump pressure & speed detection board

2.4.14 The pressure sensor

- The pressure sensor introduction
The pressure sensor adopts the resistive Wheatstone bridge strain gauge pressure sensor. When strain gauge pressure changes, the bridge will be out of balance, the

voltage will output and the pressure signal is received.

Note: The pressure sensor is a very sensitive and precision device. Do not touch the surrounding parts or pull the connecting lines while dismounting.

- b. The main testing points: none.
- c. Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	/	4	Socket to connect the infusion pump pressure & speed detection board	Connect J3 in the infusion pump pressure & speed detection board

2.4.15 Bubble sensor

- a. Bubble sensor introduction

The bubble sensor adopts the ultrasonic testing method. Its testing is not affected by the color of infusion line. It has high-sensitivity (greater than or equal to 30ul bubbles) and a fast response speed.

- b. The main testing points: none.
- c. Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	/	3	Socket to connect the power board	Connect J10 in the power board

2.4.16 Hall switch board

- a. The circuit function introduction

(1) To output the open and close signals of the pump door through the two micro switch opening or closing. Then according to the distance between hall parts and hall switch to output signals and eventually output to the power board so as to judge the open or close of the pump door.

(2) Press the micro switch button to control the open and close the anti-free-flow clamp.

- b. The main testing points: none.
- c. Number and definition the socket

No.	Socket	PIN No.	Description	Explanation
1	/	3	Socket to connect the infusion pump pressure & speed detection board	Connect J6 in the infusion pump pressure & speed detection board

2.4.17 The open motor

- c. The motor introduction: adopts a stepper motor with 25mm diameter and 24 steps, 2

phases driving. The working voltage is 24V.

d. Number and definition of the socket

No.	Socket	PIN No.	Description	Explanation
1	/	4	Socket to connect the infusion pump pressure & speed detection board	Connect J9 in the infusion pump pressure & speed detection board

2.5 The software

The SYS-6010 consists of 3 softwares.

The power board software: to obtain all kinds of sensor signals, control and calculate of the stepper motor, obtain self-testing information and process the alarm information. It corresponds with the UI software and power failure alarm software.

The UI software: to be used in the monitor interface display, button operation, function setting and operating of the human-computer interaction, data storage, nurse call and voice function. It corresponds with the power board.

The power failure alarm software: under the situation of built-in battery failure and external power is cut-off, to control the device to output a light and voice high-level alarm. It corresponds with the power board.

2.6 Specifications

Power supply	AC power supply: AC 100-240V,50/60 Hz, power, consumption 25 VA External DC power supply: DC 12 V Internal battery: lithium battery 11.1 V 1500 mAh Time of continuous use: no less than 5 hours (for infusion at 25 mL /h with a new battery)
Compatible IV sets	All disposable infusion sets of 20d/mL and 60d/mL conform to the standards.
Infusion mode	6 infusion modes: rate, drop, time, standard weight, weight, multi-rate
Infusion setting range	0.1-1200.0mL /h or (1-400d/min) (20d/mL IV set) 0.1-200.0mL /h or (1-200d/min) (60d/mL IV set) Least increment of 0.1 mL /h(or 1d/min)
Total volume display	0.0-99999.9mL
Accuracy	±5%
Purge operation	400.0mL /h (20d/mL IV set) 200.0mL /h (60d/mL IV set)
Bolus operation	35.0-1200.0mL /h(20d/mL IV set)

	35.0-400.0mL /h(60d/mL IV set) Auto-calculating the speed by the bolus volume
KVO rate	0.1-5.0mL /h
Occlusion level	11 levels are available 225mmHg~975mmHg, 11 level 1:225mmHg 2:300mmHg 3:375mmHg 4:450mmHg 5:525mmHg 6:600mmHg 7:675mmHg 8:750mmHg 9:825mmHg 10:900mmHg 11:975mmHg
Alarm	Near End, Infusion End, Occlusion Alarm, Low Battery, Battery Empty, No Battery, No Power Supply, Door open, Air Bubble, No Drop Sensor, No Drop Signal, Drop Error, Reminder Alarm, and alarms during self-test and running.
Special function	Repeat alarming: If there is still alarm after mute alarm sound, it will alarm again in 2 minutes. Event recording: can store and playback 1000 events maximum Alarm volume setting: 5 levels of alarm voice are available. Power supply switching: When AC/DC power supply is cut off, the infusion automatically switches to internal battery supply. Barcode scanning: Input the patient information by barcode scanning.
Operating conditions	Ambient temperature: 5°C-40°C Relative humidity: 15%-95%, nocondensation Atmospheric pressure: 70.0 kPa-106.0 kPa
Storage conditions	Ambient temperature: -20°C-60°C Relative humidity: 10%-98%, no condensation Atmospheric pressure: 22.0.0 kPa-107.4.0 kPa
Classification	Class I, Type CF, IPX2

Maintenance

3 Maintenance

3.1 Cleaning and Sterilizing

- Before cleaning the pump, be sure to turn off the power and disconnect the AC or DC power cables.
- If any solution spills on the pump or the pump gets heavily soiled, wipe it with wet soft cloth dampened with cold or lukewarm water.
- Use a piece of dry soft cloth to clean the AC power supply socket, USB socket or the nurse call socket, ensure that the socket is dry before using it.
- If the pusher and clutch need to be dismantled and cleaned, contact your distributors or manufacturer.
- Do not use organic solvent such as alcohol or thinner.
- If disinfection is necessary, using the common disinfectors such as Chlorhexidine gluconate and Benzalkonium chloride. After using the agent with a soft cloth, wipe off it with a soft cloth dampened with water or warm water. When using the disinfecting agent, follow the caution of each agent.
- When sterilize the system, use EOG(Ethylene Oxide Gas) at the temperature of 50°C or less and relative humidity of 60% or less. After sterilize, ventilate the room for over 24 hours or place system in aerator for 8 hours. The above is only a guideline, use adequate method to check sterilization results.
- The infusion pump cannot be autoclaved.
- Never use a dryer or similar device to dry the infusion pump.
- If liquid spills onto the pump, check whether the pump still functions normally. Test the insulation and leakage current when necessary.
- Do not soak the infusion pump into water.

3.2 Periodic Maintenance

Perform a periodic maintenance inspection to ensure safe operation and the longest possible life of the infusion pump, and check the infusion pump once every six months. You can maintain some items by yourself and contact your local distributor to maintain some other items.

3.2.1 Checking the appearance (before every use)

- Appearance checking: There are no cracks and damages.
- Key operations: Keys response smoothly and effectively.
- No abnormal noise during the peristaltic plate working.

3.2.2 Checking the power cable (before every use)

- Check the appearance of the power cable. If the appearance is damaged and the plug and the socket are in poor contact, contact the distributor for replacement in time.
- Connect the infusion pump to the AC/DC power and there is no indication of powering on, contact the distributor for maintenance in time.

3.2.3 Checking the infusion rate (every 6 months)

Check the infusion flow once every 6 months by the graduate and timer.

Checking method: Use the preset brand IV set with a 60 mL /h rate, fill the infusion bag with distilled water and access the graduate. Start, run 10 minutes at a 60mL /h rate and observe the liquid volume in the graduate where 9.5-10.5mL is qualified.

3.2.4 Alarm (every 6 months)

1) Air bubble in tubing

Checking method: Simulate the infusion process and create a bubble by tilting the drip chamber during the infusion. There should be a visible alarm [Air Bubble] on the screen, red alarm light flashing, and audible alarm cycling every 15 seconds like beep-beep-beep ...beep-beep...beep-beep-beep ...beep-beep....

2) Occlusion

Checking method: Use the preset brand IV set with a 25 mL /h rate, set the occlusion 検 level to be P2, and then close the roller clamp of the IV set. Start the infusion. Within 1 minute, there should be a visible alarm [Occlusion] on the screen, red alarm light flashing, and audible alarm cycling every 15 seconds like beep-beep-beep ...beep-beep...beep-beep-beep ...beep-beep....

3.2.5 Electric and mechanical safety (every 12 months)

To ensure safety, test the insulation voltage, leakage current and earthing resistance according to the IEC 60601-1.

3.2.6 Checking the internal battery (every 6 months)

Perform the following inspections on the battery every 6 months:

- Connect to the AC power supply to recharge the battery for over 10 hours.
- Turn on the power.
- Set the infusion rate to 25 mL/h and start the infusion. Record the start time.
- Operate the system until it stops infusing due to low battery alarm. Record the finish time.
 - (1) If the time from the start of the infusing to end of operation is 4 hours or more, the battery condition is good.
 - (2) If the time from the start of the infusing to end of operation is 1 to 1.5 hours, the battery condition is reaching its service life.
 - (3) If the time from the start of the infusing to end of operation is less than 1 hour, the battery has reached its service life. Replace the battery. You are advised to contact the distributor to replace the battery.
- When the battery lever check is complete, recharge the battery for next use.

3.3 Replacing the Battery

- **Remove an internal battery.**

- (1) Use a screwdriver to loosen the battery cover fixing screws at the bottom of the pump.
- (2) Remove the battery cover.
- (3) Disconnect the battery cable connector.
- (4) Remove the battery.

- **Install the internal battery**

- (1) Insert the connector of the battery cable into the battery.
- (2) Insert the new battery into the battery compartment.
- (3) Attach the battery cover.
- (4) Use a screwdriver to tighten the screws securing the battery cover.

Software Upgrade

4 Software Upgrade

SYS-6010's software upgrade involves two parts, being the Driver and UI upgrade.

4.1 How to upgrade the driver software


➤ Upgrade Tool:





➤ Connection: pay attention to the cable color(Blue color cable corresponds to a white triangle) and USB port




➤ Open the software  "MPLAB IDE v8.86" on the desktop of windows;




- In the menu bar, click on the "Configure" -> "Select Device", as shown in figure choice chip PIC16F1527, click OK to exit;



- In the menu bar, click on the "Programmer" -> "Select Programmer" -> Select "PICKit 3" connection PICKit 3;



- Click ok




- PICkit 3 connection is successful, if not connected successfully check the connection is correct or not;

```

MPLAB IDE v8.86 - [Output]
File Edit View Project Debugger Programmer Tools Configure Window
Build Version Control Find in Files PICkit 3

PICkit 3 detected
Connecting to PICkit 3...
Firmware Suite Version..... 01.27.20
Firmware type.....Enhanced Midrange
PICkit 3 Connected.
PK3Err0045: You must connect to a target device to use PICkit
3.
Device ID Revision = 00000005
  
```

- In the menu bar, click on the "Programmer" => "Settings"=> "power", Click the item "power target circuit from PICkit 3"; then click "OK"





- After the success of the set display device ID;

```


MPLAB IDE v8.86 - [Output]
File Edit View Project Debugger Programmer Tools Configure Window
Build Version Control Find in Files PICkit 3

PICkit 3 detected
Connecting to PICkit 3...
Firmware Suite Version..... 01.27.20
Firmware type.....Enhanced Midrange
PICkit 3 Connected.
PK3Err0045: You must connect to a target device to use PICkit
3.
Device ID Revision = 00000005
  
```


- From the menu bar "File" - > "Import" hex File Import products. Burn MP - 30, for example, find the corresponding version in the folder, and open;




➤ Import success;



➤ Click the icon "program", then begin to download;




➤ The download is complete, the software prompts




4.2 UI software Upgrade


- Upgrade tool: please use long calbe to connect pump



- Connect the upper USB port of pump;




- Open the software "Programer. Exe"  Programer ; Turned off pump;



- Turn on the machine, then it would show “connect success”




- In the Programmer software, click "add", then open the upgrade file;




- Select the Com and upgrade file , click start download;





➤ Then it would show “download complete”;



5 Alarm, fault and troubleshooting

Warning!

- There are 90-240V AC high-voltage current inside the MP-60A and 400V DC high-voltage current produced by the switching power supply. A careless touching may cause serious personal injury to the service engineer when they are conducting the checks and troubleshooting. Please be strict in accordance with the requirements of the rules and regulations to maintain the device by professional maintainers.
- When using equipment to check power board or other parts, short-circuit or incorrect measurement may cause serious personal injury or a severe damage to the device.
- After every disassembling and maintenance, check the safety and function of the device. See details in section 7.

The basic instruments of maintaining and checking

- PC (to download software and fill in forms)
- Software download connector
- To be used infusion set
- Digital multimeter
- 20-50mL graduate (accurate to 0.2mL)
- Stopwatch or timing device

Note!

- **The graduate to inspect the infusion precision, is just used for general maintenance.**
- **The special materials for disassembling should refer to the requirements in the disassembling and installation chapter.**
- **The safety project inspection needs specialized testing equipments. After disassembling and installing the device, the safety of the pump must be tested.**

5.1 Alarm information

During normal operation, the pump generate visible and audible alarm to remind users to deal with the alarms in time.

The alarm can be categorized into clinical alarm and technical alarm. Clinical alarm is generated according to the operation state during the general using of the device. It is not equipment failure but to remind users to deal with the device working in time, such as when an infusion is finished or a reminder alarm. A clinical alarm is sounded along with concrete wording of specific information on the screen.

A technical alarm is caused by abnormal data inside the device during normal use of the device. It can be an equipment failure or not, showing as ERR-XXXX.

5.1.1 The Alarm Levels

Alarm	Sound	Light
Low-level alarm	Give out a beep at the time of every alarm.	The yellow indicator illuminates.
Mid-level alarm	Give out three beeps at intervals of 15 seconds.	The yellow indicator flashing.
High-level alarm	Give out a series of beeps at intervals of 15 seconds.	The red indicator flashing.
Prompt message	Give out two short beeps.	The corresponding indicator flashing.

5.1.2 The Clinical Alarms and Troubleshooting

Alarm Symptoms	Alarm level	Causes	Troubleshooting
No Power Supply	Low	No external AC/DC power supply is connected.	Immediately connect the AC power supply or the external DC power supply.
No Battery	Middle	No internal battery or the internal battery operates	Install a (new) battery.

Alarm Symptoms	Alarm level	Causes	Troubleshooting
		abnormally.	
Low Battery	Middle	The internal battery is running critically low.	Immediately connect an AC power supply or an external DC power supply.
Near End	Low	The infusion will finish within 3 minutes.	Wait until the infusion finishes.
Occlusion Alarm	High	1. The tubing is occluded. 2. Low occlusion level has been set during high viscosity.	Click [STOP] to clear the alarm, troubleshoot, and then continue the infusion.
Door Open	High	Door is open during infusion.	Click [STOP] to clear the alarm. Check whether the pump door has been closed.
Air Alarm	High	1. Air bubble in the infusion line. 2. The flatten tube is fixed inside the air bubble detector.	Click [STOP] to clear the alarm. Check the tubing to release the air bubble.
Infusion End	High	The VTBI or the infusion time set finishes.	Click [STOP] to clear the alarm.
Reminder Alarm	Middle	Operation forgotten (no operation on any key for 2 minutes after the IV set loaded.	Click any key to clear the alarm.

5.1.3 Technology alarms and troubleshooting

No.	Trouble	Possible Cause	Troubleshooting
------------	----------------	-----------------------	------------------------

1	ERR-0003	Door position fault	Replace the hall switch board
2	ERR-0004	Pressure sensor is falling or abnormal pressed.	Reinstall the tube again or Replace the pressure sensor and the infusion pressure speed auto driving board.
3	ERR-0005	Stepper motor speed error	Replace the hall switch board or Replace the stepper motor
4	ERR-0006	The 5V power supply voltage or 18V motor drive voltage is abnormal.	Replace DC/DC power board
5	ERR-0007	Pump tablets fault	Replace the peristaltic pump tablets position test board
6	ERR-0008	SYS02 board and power board correspondence fault	SYS02 board damaged or the cable disconnection
7	ERR-0009	System fault	Update software
8	ERR-0010	EEPROM fault	Replace the SYS02 board

5.2 Common faults and Troubleshooting

The device fault checking:

NO.	Failure phenomenon	Cause	Solution
1	Can not turn on pump	No AC power	Check AC power, panel indicator, replace battery
		Key panel failure	Replace panel
		FPC cable failure	Replace FPC cable
		Power board failure	Replace power board
2	Black screen and alarm indicator is off	No AC power	Check AC power, panel indicator, replace battery
		FPC cable failure	Replace FPC cable
		SYS02 board failure	Replace SYS02 board
3	Black screen and alarm indicator is on	FPC cable failure	Replace FPC cable
		SYS02 board failure	Replace SYS02 board
		Touch screen failure	Replace touch screen
4	White screen	Touch screen failure	Replace touch screen
5	Touch screen is no	Touch screen not calibrate	Calibrate touch screen

	response when clicking	The surface of screen is dirty	Clean screen
		Touch screen install incorrectly	Install screen correctly
		FPC cable failure	Replace FPC cable
		Touch screen failure	Replace screen
6	Air bubble alarm	IV Set with air bubbles	Get rid of air bubbles
		IV Set damaged	Replace IV Set
		Air bubble sensor failure	Replace air bubble sensor
7	Alarm of air bubble of No IV set when starting	IV set with air bubbles	Get rid of air bubbles
		IV set install error	Install IV set correctly
		Air bubble sensor failure	Replace air bubble sensor
8	Occlusion alarm	Blocking in IV set	Clear blocking
		Select wrong brand	Select right brand
		IV set of Custom brand not calibrate	Calibrate pressure value

The disassembly and installation

6 The disassembling and installing

6.1 Attentions

Please do the disassembling in an anti-static environment. To avoid damage to the components, please wear wrist strap or other anti-static equipment.

Before each disassembling, please carefully record the fixed way and layout of the board connectors and their cables to ensure consistent before and after the disassembling.

During reinstallation, use screws removed from the original position, otherwise, it can cause a weak installing or damage to the position, involve other positions to be weakly installed or damaged and eventually damage the device; And during use after the reinstallation, the screws and components might fall off, thus causing unpredictable damage or injury to the product or human body.

Disassemble the device in a correct order. Incorrect order or forcible assembly might cause irreversible damage to the device.

Before separating the components, please be sure that all connections are plugged off. During separation, take care not to break the connection wire and damage the connectors.

Please store the removed screws and other parts by category for easy reach during

reinstalling, and avoid dropping, contaminating or losing them.

During disassembly and reinstallation, pay attention to protecting the waterproof seal and do not pull the seal strip hard. Please do not miss or damage it during installation.

After assembling, check if there are extra screws or other components, If there are, a reassembling is needed.

Warnings!

- **Before and after each assembling, be sure to cut off the external AC power supply to avoid causing a human body injury.**
- **During assembling, avoid touching the circuit board and connection cables by sharp instruments and avoid shorting and breaking circuits, or it may cause a serious damage to the device.**

6.2 Assembly Tools

6.2.1 General Tools

- Philips Screwdriver
- M3 inner hexagon screwdriver M3
- Pliers
- M3 nut spanner M3
- Diagonal cutting nippers
- Anti-static wrist strap
- Several Nylon cable ties

6.2.2 Special instruments and materials

- Several screw cover caps

6.3 Disassembly

6.3.1 Remove the battery




Figure 6.3-1 Note:

Code name

- 1 –Cover of battery
- 2 – PB2.3 tapping screw
- 3 –battery pack

1) Use phillips screwdriver to take out the 4 PB2.3 tapping screws (figure 6.3-1/code name 2) in the cover of battery (figure 6.3-1/code name 1).

2) Remove the cover of battery.

3) Take out the battery pack lightly (figure 6.3-1/code name 3)

Attention: pull out the battery pack lightly for there are connection wires between the battery pack and the device.

4) The connection wires connect the battery and the device by a connector .Press the block button in the connector and separate the plug and socket.

6.3.2 Remove the rear cover




Figure 6.3-2

Figure 6.3-2 Notes:

Code Name

1 –Rear Cover

2 – PB2.3 PB2.3 tapping screw

3 – PT3.0 tapping screw

1) Use phillips screwdriver to take out the 6 PB2.3 tapping screws (figure 6.3-2/code name 2) in the rear cover (figure 6.3-2/code name 1).

2) Use phillips screwdriver to take out the 6 PT3.0 tapping screws (figure 6.3-2/code name 3) in the rear cover.

3) Balance the rear cover with left hand, then use screwdriver with right hand to stick the inner side of the AC power socket. Use a uniform force to inner cover and separate the rear cover carefully.

Attention: be careful to remove the rear cover, avoid pulling out other components inside the device. After removing the rear cover, please take a picture first to know well the circuit connections and wire bindings.

6.3.3 Remove the Baffle




Figure 6.3-3

Figure 6.3-3 Notes

Code Name

1 –PCBA stand component

2 –PCBA stand board

3 –PB2.3 tapping screw

4–FRC flat cable

5 –Middle plate

Attention: please connect the power supply, press [OPEN] to open the pump door, then remove the battery.

1)Use the short edge of 2# inner hexagon spanner to push-out a certain space between baffle-axis (figure 6.3-3/code name 2) and round decorative cover (figure 6.3-3/code name 3), take out the round decorative cover.

2)Use pliers tip to clamp the pushed-out baffle-axis, pull it out completely.

3)Remove the baffle (figure 6.3-3/code name 1).

6.3.4 Remove the Bracket for PCB




Figure 6.3-4 Notes:

Code Name

- 1 –Round Decorative Cover
- 2 – PB2.3 Tapping Screw
- 3 –Square Decorative Cover
- 4 –Cover Board
- 5 –Panel
- 6 – FPC Flat Cable
- 7 – Medium plate
- 8 –Flat Cable Sealing Groove
- 9 – PCB Bracket for PCB

1)Use pliers to take down the 9 round decorative covers (figure 6.3-4/code name 1) and 2 square decorative covers (figure 6.3-4/code name 3) in the cover board (figure 6.3-4/code name 4).

Attention: use pliers tip to insert between the decorative covers and stick them up, avoid breaking the cover board.

2)Use phillips screwdriver to remove the 7 PB2.3 tapping screws.

3)Take out the cover board lightly.

Attention: Do not miss the magnets in the cover board.

4)Remove the FRC flat cable (figure 6.3-4/code name 6) inside the cover board (figure 6.3-4/code name 5).

5)Remove the flat cable sealing groove (figure 6.3-4/code name 8) in the middle board (figure 6.3-4/code name 7).

6)Pull out the bracket for PCB (figure 6.3-4/code name 9) lightly, and pull out the FPC flat cable at the same time.

Attention: take care the FPC flat cable not to be hung.

6.3.5 Remove the pump door components




Figure 6.3-5 Notes:

Code Name

1 –Spring

2 –Door Shaft

1)Use pliers to take out the spring (figure 6.3-5/code name 2) in the door shaft (figure 6.3-5/code name 2).

Attention: move the door shaft for a certain distance in the direction of spring to ensure the pliers tip can stick in and stick up the spring. Take care not to scratch the surface of the door shaft, once the surface of the door shaft is scratched, it is easy to be stuck in the shaft sleeve.

2)Use pliers to take out the door shaft.

Attention: move the door shaft for a certain distance in the direction of spring, stick the pliers tip into the groove of the door shaft and pull it out. Take care not to scratch the surface of the door shaft.

3)Remove the cover board (figure 6.3-5/code name 5).

6.3.6 Remove Components of the Peristaltic Pump




Figure 6.3-6

Figure 6.3-6 Notes:

Code Name

- | | |
|------------------------|---------------------------------------|
| 1 –Pressure sensor | 2 –M4 Machine Screw |
| 3 – M4 Spring Washer | 4 – M4 Flat Washer |
| 5 –Middle board Label | 6 –M3.0 Machine Screw |
| 7 – M3.0 Spring Washer | 8 –Components of the Peristaltic pump |
- 1) Use phillips screwdriver to take out the M4 machine screws (figure 6.3-6/code name 2) in the pressure sensor (figure 6.3-6/code name 3).
 - 2) Remove the pressure sensor.
 Attention: Do not break off the connection wire of the pressure sensor and take down the pressure from the board.
 - 3) Use pliers to tear down the middle board label (figure 6.3-6/code name 5).
 Attention: Do not break the surface.
 - 4) Use phillips screwdriver to take out the 6 M3.0 machine screws (figure 6.3-6/code name 6). Remove the M3.0 spring washer (figure 6.3-6/code name 7).
 - 5) Remove components of the peristaltic pump (figure 6.3-6/code name 8).

6.3.7 Remove the Door Drive and Components of the Middle board




Figure 6.3-7

Figure 6.3-7 Notes:

Code Name

- | | |
|------------------------------------|------------------------|
| 1 –Alarm lamp board | 2 –PB2.3 tapping screw |
| 3 –SYS02 module and display screen | 4 – Lamp shade |
| 5 –Lamp board | 6 –Baffle |
| 7 –FPC flat cable | 8 –Panel |
| 9 –M3 machine screw | 10 –M3 spring washer |
| 11 –Door catch | 12 –Baffle spring |
| 13 –Baffle block | 14 –Membrane |
| 15 –Baffle shaft | |

1) Use phillips screwdriver to remove 2 M3 machine screws (figure 6.3-7/code name 9) in the baffle block (figure 6.3-7/code name 13), take down the 2 M3 spring washer (figure 6.3-7/code name 10).

Attention: hand press on the baffle, avoid the baffle spring pushes out it.

2) Take down the baffle block.

3) Turn up the baffle (figure 6.3-7/code name 6), take down the 2 baffle springs (figure 6.3-7/code name 12).

4) Use pliers to take down the baffle shaft (figure 6.3-7/code name 15).

Attention: use the pliers tip to push out it backward for a length and pull put it.

5) Take down the baffle.

6) Use phillips screwdriver to remove the 2 PB2.3 tapping screws (figure 6.3-7/code name 2) in the lamp shade (figure 6.3-7/code name 4).

7) Take down the lamp shade and lamp board (figure 6.3-7/code name 5).

8) Take down the FPC flat cable (figure 6.3-7/code name 7).

9) Remove the connection wire in the alarm lamp board (figure 6.3-7/code name 1) from

SYS02 module board (figure 6.3-7/code name 3) lightly.

10) Use phillips to remove the 2 PB2.3 tapping screws in the alarm lamp board and take down the alarm lamp board and lamp shade lightly.

11) Take out the SYS02 module board and display screen lightly.

Attention: put the display screen upward, avoid scratching or pollute it.

12) Use phillips to remove the 4 machine screws in the 2 door catches (figure 6.3-7/code name 11) and take down the 4 M3 spring washers.

Attention: do not miss the plastic chips under the door catch, 2 pieces of each door catch.

13) Take down the door catch.

14) Tear down the membrane (figure 6.3-7/code name 14) from the panel (figure 6.3-7/code name 8).

Note: the installation operation method: Reverse operation according to the steps above.

The Safety Test

7 The Safety Test

Warning!

- The electric and safety test is used to check the abnormal electric hazards. If the hazards can not be discovered timely, it might cause injury to the patient and operator.
- The equipment of the electric and safety test adopts the safety analyzer, maintainers should master the usage.
- The electric and safety test completely complies with the standard of IEC/EN60601-1.
- Please conduct the electric and safety test in a normal circumstance (including temperature, humidity and atmosphere).
- The electric and safety test in the chapter takes the example of 601 safety analyzer. Different brands of the safety analyzer might have a different analysis.
- If the electric and safety test fails, please contract our after-sales service engineers.

7.1 Electric Leakage Test of Enclosure

1. Connect the 601 safety analyzer to 264 VAC/60 Hz power source.
2. Connect the applied part of the tested equipment by using the applied part connection tooling. The SUM end of the connection tooling for applied part is connected to the RA end of the safety analyzer.

3. Via power cord, connect the tested equipment to the auxiliary power output port on 601 safety analyzer.
4. Connect one end of the red testing wire to the “Red input terminal” on the safety analyzer, and clip another end tightly onto the metal foil on the surface of tested equipment enclosure.
5. Turn on the power source of 601 safety analyzer and press “5—Enclosure leakage” on the panel of 601 safety analyzer to enter the enclosure leakage current testing interface.
6. The enclosure leakage current shall not exceed 100 μ A under normal conditions and shall not exceed 300 μ A in event of a single fault.

7.2 Electric Leakage Test of Ground


1. Connect the 601 safety analyzer to 264 VAC/60 Hz power source.
2. Connect the tested applied part to the the RA end of the safety analyzer.
3. Connect the 601 safety analyzer to 264 VAC/60 Hz power source.
4. Turn on the power source of 601 safety analyzer and press “4—Ground leakage” on the panel of 601 safety analyzer to enter the ground leakage current testing interface.
5. The enclosure leakage current shall not exceed 500 μ A under normal conditions and shall not exceed 1000 μ A in event of a single fault.

7.3 Electric Leakage Test of Patient

1. Connect the applied part of the tested equipment by using the applied part connection tooling. The SUM end of the connection tooling for applied part is connected to the RA end of the safety analyzer.
2. Via power cord, connect the tested equipment to the auxiliary power output port on 601 safety analyzer.
3. Turn on the power source of 601 safety analyzer and press “6—Patient leakage” on the panel of 601 safety analyzer.
4. Keep pressing “APPLIED PART” key to select AC or DC measuring. If DC is selected, “DC” is displayed behind the limit value.
5. The patient leakage current shall not exceed 10 μ A under normal conditions and shall not exceed 50 μ A in event of a single fault.

Maintenance Spare Parts List

8 Maintenance Spare Parts List



Code Name

No.	Part number	Description	No.	Part number	Description
1	60Z0020261	Cover of Battery	a	7901000004	Interface Board
2	7404000006	Battery	b	9114002321	Pressure&speed detecting board
3	9114002441	Rear Cover	c	9114002252	Power Control Board
4	9108008001	PCB Bracket Component	d	9114002242	DC-DC Power Board
5	60Z00204911	Sliding Block	e	9114002231	AC-DC Power Board
6	60Z0020531	Cover of Anti-Free-Flow Clamp	f	9114002321	Pressure Sensor
7	60Z0020191	Anti-Free-Flow Clamp	g	9114002282	Alarm Lamp Board
8	9114002480	Middle Board	h	9116000011	SYS02 Module Board
9	60Z0020131	Cover Board	i	9114002242	Peristaltic Pump Control Board
10	9108011001	Panel	j	9114006001	Peristaltic Tablets Test Board
11	9114002480	Baffle Block		9114002262	Interface Board
12	9114002480	Baffle		9114004001	Hall Board
13	9114002481	Door Driver Component		7901000003	Air Bubble Sensor

MADCAPTAIN TECHNOLOGY CO., LTD

The after-sales service address: 12th Floor, Baiwang Research Building, No.5158
Shahe West Road, Xili, Nanshan District, Shenzhen, P.R.China

Telephone: 0755-26953369

Fax: 0755-26001651

Postal: 518055

Website:<http://www.medcaptain.com>

E-mail:info@medcaptain.com

2014 All rights reserved. SHENZHEN MEDCAPTAIN MEDICAL
NOLOGY CO., LTD.