


# **Environmental Remediation Equipment inc.**

T: 514.326.8852 F: 514.326.8961 sales@ereinc.com 8605, Champ d'Eau. Montréal, Qc. H1P 3B8



# **Alconox**®



**Cleaning method:** Soak, brush, sponge, cloth, ultrasonic, flow through clean-inplace. Will foam—not for spray or machine use.

**Directions:** Make a fresh 1% solution (2 1/2 Tbsp. per gal., 1 1/4 oz. per gal. or 10 grams per liter) in cold, warm, or hot water. If available use warm water. Use cold water for blood stains. For difficult soils, raise water temperature and use more detergent. Clean by soak, circulate, wipe, or ultrasonic method. Not for spray machines, will foam. For nonabrasive scouring, make paste. Use 2% solution to soak frozen stopcocks. To remove silver tarnish, soak in 1% solution in aluminum container. RINSE THOROUGHLY—preferably with running water. For critical cleaning, do final or all rinsing in distilled, deionized, or purified water. For food contact surfaces, rinse with potable water. Used on a wide range of glass, ceramic, plastic, and metal surfaces. Corrosion testing may be advisable.

Convenient sizes	Alconox cat#
Case 9 x 4 lb. Boxes	AL-1104
25 lb. Carton	AL-1125
50 lb. Carton	AL-1150
100 lb. Drum	AL-1101
300 lb. Drum	AL-1103
Case 12 Box 50 x 1/2 oz.	AL-1112

1 lb. makes 13 gal. cleaning solution

### **Powdered Precision Cleaner**

- Concentrated to save you money
- · Biodegradable and readily disposable
- · Penetrating wetting power to save you time
- · Replaces corrosive acids and hazardous solvents
- Free rinsing to give you reliable results and no interfering residues
- Use to pass your cleaning validation tests for lab accreditation and plant inspection approval

**Used to clean:** Healthcare instruments, laboratory ware, vacuum equipment, tissue culture ware, personal protective equipment, sampling apparatus, catheters, tubing, pipes, radioactive contaminated articles, optical parts, electronic components, pharmaceutical apparatus, cosmetics manufacturing equipment, metal castings, forgings and stampings, industrial parts, tanks and reactors. Authorized by USDA for use in federally inspected meat and poultry plants. Passes inhibitory residue test for water analysis. FDA certified.

**Used to remove:** Soil, grit, grime, buffing compound, slime, grease, oils, blood, tissue, salts, deposits, particulates, solvents, chemicals, radioisotopes, radioactive contaminations, silicon oils, mold release agents.

**Surfaces cleaned:** Corrosion inhibited formulation recommended for glass, metal, stainless steel, porcelain, ceramic, plastic, rubber and fiberglass. Can be used on soft metals such as copper, aluminum, zinc and magnesium if rinsed promptly. Corrosion testing may be advisable.







Physical data	Typical value
pH of 1% solution	9.5
Flash Point (degrees F)	None
Phosphate Content (as Phosphorus)	7.3%
Organic Carbon (% calculated w/w)	11%
Fragrance Content	0%
Surface Tension 1% Sol'n (Dyne/cm)	32
Percent active ingredients	100%

Color: White and cream colored flakes

Form: Powder

Solubility in Water: To 10% (w/w) at ambient temperature

Hard Water Effectiveness: Highly Effective

Biodegradability: Biodegradable Foam Tendency: High Foaming

Shelf Life: Two years from the date of manufacture

## **Health Safety Information**

OSHA Hazardous Ingredients: None RCRA Hazard Class: Non-hazardous

Flammability: Non-flammable

Latex Content: None in detergent, packaging materials or

adhesives.

Oral Toxicity: LD50 > 5000 mg/kg oral rat

No ingredient defined as an oral toxicant by OSHA Eve Irritation: Mild to Moderate eve irritant if not rinsed Inhalation Toxicity: Non-irritating solution, powder a

potential irritant **VOC Content: 0%** 

Carcinogenicity: NTP = No IARC = No OSHA = No

All ingredients in Alconox are listed in TSCA inventory. USDA NSF cat A1

### **Chemical Description**

Alconox consists primarily of a homogeneous blend of sodium linear alkylaryl sulfonate, alcohol sulfate, phosphates, carbonates. Alconox is anionic in nature.

## **Cleaning Validation Methods**

Test a parameter of rinse water before and after rinsing the cleaned surface. No significant change in the parameter indicates no detectable residue. Parameters measured include: pH, conductivity, UV, TOC, HPLC, sodium concentration, phosphorus concentration, anionic surfactant concentration using inexpensive detergent water testing kits, surface tension, and surface analysis. Pharmaceutical Cleaning Validation references are available in the tech info section of www.alconox.com.

#### **Precautions**

No special precautions other than good industrial hygiene and safety practices employed with any industrial chemical (see Directions). A Material Safety Data Sheet is available at www.alconox.com or by calling fax-on-demand at 914-948-4040 and following the prompts. Use fax document # 0311. Contact Alconox, Inc. for purchase specifications. Typical data is not a specification. While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

#### AL-1104-1

4 lb box

## AL-1104

Case 9 x 4 lb boxes

#### AL-1112-1

Dispenser box (50 x 1/2 oz)



