



Owner's Manual

Table of Contents

Introduction

Practice before Operating	2
---------------------------------	---

Safety Instructions

Safety Considerations	4
Driving Outdoors	7
Use caution when driving on inclines	9

Familiarize Yourself with your Scooter

The Batteries	11
Feature Diagram	12
Terminology	13
Disassembly and Assembly Base step of the scooter	14

Operating your Scooter

Controls	17
Seat	19
Tiller Angle Adjustment	21
Free Wheel Mode	21

Technical Specifications	22
---------------------------------------	-----------

ICE Symbols	23
--------------------------	-----------

Warranty

Limited Warranty	24
------------------------	----



Welcome aboard your new scooter, and thank you for choosing our product. Please read this manual carefully, and follow all instructions before attempting to operate your scooter for the first time. If there is anything in this manual that you do not understand, or if you require additional assistance for setting up your scooter, please contact your local dealer.

This latest model is designed for specific practical user needs, combining solid, rugged construction, and modern high-tech electronics, to enhance safety and performance.

With a state-of-the-art, programmable electronic control system, your scooter can be programmed and adjusted within a given range of its performance characteristics, to suit your individual needs. The controller is set up at the factory to give the scooter nominal operating performance characteristics.

After becoming familiar with the basic operation of the scooter, you may wish to customize the settings to fit your own personal preferences. A wide range of customization options can be adjusted such as acceleration, deceleration, maximum speed, turning speed, safety controls, and so on. Contact your local dealer for advice on additional equipment you may need.

Having your scooter checked regularly by your local dealer is the best way to ensure smooth operation and safety.

This manual provides users practical tips and information on safety issues, operation, and maintenance. Please read it very carefully to ensure your maximum enjoyment and to fully benefit from your independence and mobility.

Whenever special advice or attention is needed, please contact your local dealer, who has the tools and know-how to provide expert servicing for your scooter.

Your satisfaction and opinions are highly valued by both your local dealer and Merits. Please be sure to fill out the enclosed warranty registration form, and return it to your local dealer. The information is necessary for providing you with the best service, and to be sure all of your needs are met.


Failure to follow these instructions may result in damage to the scooter or serious injury.

■ Practice Before Operating


Find an open area such as a park and have an assistant to help you practice until you have confidence operating this vehicle.

Make sure the unit is OFF before getting on or off. Set the speed control knob according to your driving ability.

We recommend that you keep the speed at the slowest position (fully counter-clockwise) until you are familiar with the driving characteristics of this vehicle.



Stop, forward, and reverse operation practice



Pull the lever on the right side towards you to move **FORWARD**.

Pull the lever on the left side toward you to move **BACKWARD**.

Getting familiar with this vehicle



**First, practice moving forward.
Be sure to set the speed to the lowest setting.**



**After becoming familiar with moving forward,
practice making "S" turns.**



**Once you are familiar with "S" turns, practice
moving in reverse. Note that at any speed control
setting, the vehicle moves more slowly in reverse
than forward**


■ Safety Considerations

DO NOT do any of the following



NO!

Do not carry any passengers



NO!

Do not drive across a slope



NO!

Do not drink and drive

Consult your physician to determine if your medications impair your ability to control this vehicle



NO!

Do not tow a trailer



NO!

Do not turn on or use hand-held personal communication devices such as citizens band (CB) radios and cellular phones

This vehicle has an immunity level of 20 V/m which should protect it from Electromagnetic Interference(EMI) from Radio Wave Sources. The rapid development of electronics, especially in the area of communications, has saturated our environment with electromagnetic (radio) waves that are emitted by television, radio and communication signals. These EM waves are invisible and their strength increases as one approaches the source. All electrical conductors act as antennas to the EM signals and, to varying degrees, all power wheelchairs and scooters are susceptible to electromagnetic interference (EMI). This interference could result in abnormal, unintentional movement and/or erratic control of the vehicle. The United States requires the preceding statement be incorporated into the user manuals for all electric powerchairs and scooters.

Powered wheelchairs and electric scooters (in this text, both will be referred to as powered wheelchairs) may be susceptible to electromagnetic interference(EMI), which is interfering electromagnetic energy emitted from sources such as radio stations, TV stations, amateur radio (HAM) transmitters, two-way radios and cellular phones. The interference (from radio wave sources) can cause the powered wheelchair to release its brakes, move by itself or move in unintended directions. It can also permanently damage the powered wheelchair's control system. The intensity of the EM energy can be measured in volts per meter (V/m). Each powered wheelchair can resist EMI up to a certain intensity. This is called the "immunity level." The higher the immunity level, the greater the protection. At this time, current technology is capable of providing at least 20 V/m of immunity level which would provide useful protection against common sources of radiated EMI.

Following the warnings listed below should reduce the chance of unintended brake release or powered wheelchair movement that could result in serious injury:

- 1) Do not turn on hand-held personal communication devices such as citizens band (CB) radios and cellular phones while the powered wheelchair is turned on.
- 2) Be aware of nearby transmitters such as radio or TV stations and try to avoid coming close to them.

- 3) If unintended movement or brake release occurs, turn the powered wheelchair off as soon as it is safe.
- 4) Be aware that adding accessories or components, or modifying the powered wheelchair, may make it more susceptible to interference from radio wave sources. (Note: there is no easy way to evaluate their effect on the overall immunity of the powered wheelchair).
- 5) Report all incidents of unintended movement or brake release to the powered wheelchair manufacturer, and note whether there is a radio wave source nearby.

TURN OFF YOUR SCOOTER AS SOON AS POSSIBLE WHEN EXPERIENCING ANY OF THE FOLLOWING:

1. Unintentional motions.
2. Unintended or uncontrollable direction.
3. Unexpected brake release.

The FDA has written to the manufacturers of power wheelchairs, asking them to test their new products to ensure they provide a reasonable degree of immunity against EMI. The letter states that powered wheelchairs should have an immunity level of at least 20 V/m, which provides a reasonable degree of protection against the more common sources of EMI. The higher the level, the greater the protection.

■ Driving Outdoors

When you are on the road, please pay attention to the following:



NO!
Do not drive in traffic.



NO!
Do not drive beside a river, port, or lake without a fence or railing.



NO!
If possible, do not drive during the rain.



NO!
If possible, do not drive during or on snow.



NO!
Do not drive off-road or on any uneven surfaced roads.



NO!
If possible, do not drive at night.



NO!

Make sure there are no obstacles behind you when in reverse.

We recommend to set the speed knob at the lowest speed for reversing.




NO!

Do not make sudden stops, weave erratically, or make sharp turns.



NO!

Keep your arms on or inside the armrests and feet on the footrest at all time.



NO!

Do not attempt to climb curbs greater that 1.5"(4cm).



NO!

Do not attempt to cross over a gap greater than 3"(7.5cm).

■ Use caution when driving on inclines

Driving on inclines is more dangerous than on level surfaces. If you fail to heed these warnings, a fall, tip-over or loss of control may occur and cause severe injury to the vehicle user or others.



NO!

Do not attempt to climb an incline greater than 6°



NO!

Do not reverse while driving up an incline.

Forward only. If you reverse while moving up an incline, it may cause the vehicle to tip over.



NO!

Do not attempt to drive across a sloping surface greater than 3°

Driving across a slope greater than 3° is very dangerous and may cause the vehicle to tip over.



NO!

Use caution when driving over soft, uneven or unprotected surfaces such as grass, gravel and decks.



NO!

Do not get on or off on an incline.

Always stop on a level surface to get on or off of the vehicle.



NO!


Use low speed while driving down an incline.

When going down an incline, the tiller will become harder to reach and handle. When braking while moving down an incline, the scooter will take longer to come to a complete stop.



NO!

Do not load or carry heavy items in the basket while driving down an incline.



YES!

Always climb or descend gradients perpendicular to the slope or ramp.

■ The Batteries

Your scooter is equipped with a removable battery pack, which contains two maintenance free, sealed lead-acid batteries. These batteries require no maintenance other than ensuring they are properly charged. If other batteries are used, check with your battery supplier for proper battery care and maintenance instructions.

Charging the Batteries

Because your batteries may only have a partial charge when you first receive your scooter, you may not experience full riding time until you have fully charged them. Your scooter is equipped with either an on-board or an off-board battery charger. Charging your batteries as specified below will ensure maximum life, power, and range.



WARNING: Use of a non-grounded receptacle could result in an electric shock

Off-Board Battery Charger

1. Turn off the power on the scooter.
2. Insert the battery charger plug in the socket on the battery pack.
3. Do not switch on the battery charger until all the plugs are in position.
4. When starting its charge, the battery charger will show an amber light. The charge is complete when the charging light turns to green.



Off-board battery charger

Feature Diagram

In this section, we will acquaint you with the many features of your micro portable scooter and how they work. Upon receipt of your micro portable scooter, inspect it for any damage. Your micro portable scooter consists of the following components.

S730



S740



- | | |
|-------------------------------|----------------------------|
| 1. Tiller | 7. Shroud |
| 2. Tilt Angle Adjustment Knob | 8. Removable Battery Pack |
| 3. Seat Ass'y | 9. Basket |
| 4. Anti-Tip Wheels | 10. Armrests |
| 5. Pin Key | 11. Pu tires (Front wheel) |
| 6. Wig-Wag Control Lever | 12. Pu tires (Rear wheel) |

■ Terminology


1. **Tiller** - Similar in operation to the handlebars of a bicycle, this device is used to steer the vehicle. It also contains the controls necessary to propel and control the vehicle.
2. **Tiller angle Adjustment Knob** - When loosened, allows the tiller angle to be set for optimum driver comfort. During operation of the vehicle, this lever must be fully tightened.
3. **Operator's Seat** - This adjustable seat can be swiveled to ease mounting / dismounting to / from the vehicle and can be removed for easier transportation of the vehicle.
4. **Anti-tip Wheels** - Not visible in figure but positioned one each side of the vehicle at the rear, these wheels prevent a tip-over backwards.
5. **Pin Key** - "Plug in" is on, "Pull up" is off .
6. **Wig-Wag Control Lever** - this is used to control the speed and direction of drive.
7. **Manual brake Release** - This lever must be in the lower position to drive the vehicle. In this position the electric brake is controlled by the vehicle controller and will be released for drive and engaged when stopped. On occasions when it is desirable to push the vehicle, this lever is moved to the upper position to permanently hold the brake release. When the Manual Brake Release lever is in the upper position, the vehicle can not be driven under power.
8. **Removable Battery Pack** - This pack holds the two 12 volt batteries and the battery charger. This module can be removed from the vehicle to reduce the vehicle weight when you intend to lift the vehicle. i.e.- into the trunk of a vehicle.
9. **Removable Front Basket** - This basket clips to the front of the tiller and is used to carry small personal items.
10. **Arm-Rests** - are width and angle adjustable. (see seating adjustments later in this manual).
11. **Anti-tip Wheels** - Not visible in figure but positioned one each side of the vehicle at the rear, these wheels prevent a tip-over backwards.

■ Disassembly and Assembly of the scooter

This scooter can be separated into four parts, the frame, basket, battery pack, and the seat. Follow the procedures below to disassemble your scooter and transport it with care.

1. Frame
2. Battery Pack
3. Seat
4. Basket

S730



S740



Disassembly of the scooter

This scooter can be separated into four parts, the frame, basket, battery pack, and the seat. Follow the procedures below to disassemble your scooter and transport it with care.

Step 1: Loosen the knob of the frame.




Step 2: Remove the seat.




Step 3: Disconnect the power plug.



Step 4: Remove the battery pack.



Step 5: Remove the basket.



Transportation:


Lift the front wheel off the ground and the anti-tip wheel remains grounded, allowing a single person to drag the scooter. Or, one person can lift the front end of the scooter while another lifts the rear end to carry the scooter.

■ Controls

Only drive within your control limitations. Loss of control of your scooter could result in serious injury to yourself or others. If your speed becomes difficult to control, release the speed engager lever and your scooter will come to a complete stop. Only use the on/off switch to stop your scooter in an emergency.


Pin Key :

- Fully insert the key into the pin key to power up (turn on) your scooter.
- Remove the key from the pin key to power down (turn off) your scooter.



Speed Control:


Speed Dial regulates the speed of the scooter. Start at the slowest speed until you feel confident controlling your scooter safely. Turn the speed dial counter-clockwise to decrease the speed. Turn the speed dial clockwise to increase the speed.



Forward / Reverse:


To move forward, pull the lever on the right side toward you.

To move backwards, pull the lever on the left side toward you.




Battery Gauge:

Indicates the charge level of the batteries. If the needle goes into the red area, the batteries need to be charged as soon as possible. If the needle is all the way to the right side of the green area, the batteries are fully charge. As the needle moves to the red area it indicates the depletion level.




Horn:

The horn is activated by pressing the horn button. The horn is useful to warn people or animals that you are coming towards them. You may also find it helpful to use when rounding blind corners or reversing.



Power eye:

When the scooter is on, and all conditions are normal, the Power Eye will be on. When there is some special situation that needs attention, the light will flash. See the reference chart below for the meanings of the different flashing patterns.



Number of Flashes	Meaning	Number of Flashes	Meaning
1	Battery needs recharging	6	Not in neutral at power up
2	Battery voltage too low	7	Speed pot error
3	Battery voltage too high	8	Motor volts error
4	Current limit time out	9	Other internal error
5	Brake fault		


Brakes and Throttle control lever: Whenever the speed engager lever is moved out of the neutral position, the electromagnetic brake will automatically release and your scooter will move. When the speed engager lever is released, it will return to the neutral position and the scooter will decelerate and come to a complete stop. The parking brake will then engage preventing further movement of your scooter. Your unit is equipped with a programmable controller that has a high peddle disable safety feature. This will prevent unexpected acceleration of the scooter, if the speed engage lever is activated the same time you turn the key "ON". To reset the controller, release the speed engage lever and turn the key "OFF" for a couple of seconds and then turn it back "ON".

! If your scooter ever moves in an unexpected manner, release the speed engager lever and turn off the power.

■ Seat


Backrest angle:

- Backrest folds forward for ease of removal / storage.



Seat Height adjustment:

- (1) Loosen the knob(A).
- (2) Remove the retaining pin(B).
- (3) Adjust the pin to the right height and insert the pin into the hole.
- (4) Tighten the knob.




Armrest angle:

- Flip the armrest up to expose the adjustable bolt.
- Turn the bolt in to lower the angle of the armrest to the desired angle.



Armrest width:

- Locate the width socket knob
(follow the armrest support down to the base of the seat).
- Loosen the socket knob by turning counter-clockwise.
- Pull the armrest in or out to reach your desired position.
- Tighten the width socket knob.



■ Tiller Angle Adjustment:

The tiller angle adjustment knob allows you to position the tiller closer or further away for better access to the controls.

- Locate angle adjusting knob located on the right side of the tiller.
- Hold the tiller with one hand and loosen the hand knob by turning counter-clockwise with the other hand.
- Adjust the angle of the tiller to a comfortable position.
- Turn the hand knob clockwise to tighten, test the tiller is set in position by attempting to push back and for to push backward and forward.




■ Free Wheel Mode

Make sure the free-wheel lever is in the drive position to prevent the scooter from rolling. never sit in your scooter if it is in "Free wheel" mode.

Manual free-wheel mode:

Your scooter features a "FREE WHEEL" mode for manual operation. To activate manually, Pull up the pin key OFF and locate the free-wheel lever on the right side of the rear section. Pull the lever to the "FREE WHEEL" position to disengage the brake. Push the lever to the "DRIVE" position to engage the brake. When the scooter is in manual free-wheel mode, you will have no brakes. You will be unable to operate the scooter. When you wish to push your scooter for a short distance, you may put it into Manual Free-Wheel mode.

free wheel position



drive position

Getting in and out:

Your scooter is designed to make getting in and out of the scooter as easy as possible. Make sure the scooter is on a level surface and the key switch is turned OFF. If necessary, raise the armrest to give you maximum space to transfer in or out of the seat. Once transfer is complete, return the armrest before operating the scooter.






! Never operate your scooter without your feet being placed on the scooter platform. Driving your scooter without your feet on the platform could cause serious bodily injury.

■ Technical Specifications - S730


Model	S730
Overall Length	95.5 cm (37.6")
Overall Width	54.3 cm (21.4")
Maximum Rider Weight	114kg (250 lbs)
Batteries	10AH x 2pcs
Battery Pack Weight	7.6kg (16.7 lbs)
Vehicle Weight	21 kg (46.2 lbs)
Seat Weight	6.6 kg (14.5 lbs)
Ground Clearance	5 cm (2")
Total Weight	35 kg (77 lbs)
Maximum Range Up To	9.6 km (6 miles)
Maximum Speed	6 kph (4 mph)
Front Wheel Size	8"PU
Rear Wheel Size	8"PU

■ Technical Specifications - S740

Model	S740
Overall Length	95.5 cm (37.6")
Overall Width	54.3 cm (21.4")
Maximum Rider Weight	100 kg (220 lbs)
Batteries	10AH x 2pcs
Battery Pack Weight	7.6kg (16.7 lbs)
Vehicle Weight	21 kg (46.2 lbs)
Seat Weight	6.6 kg (14.5 lbs)
Ground Clearance	5 cm (2")
Total Weight	35 kg (77 lbs)
Maximum Range Up To	9.6 km (6 miles)
Maximum Speed	6 kph (4 mph)
Front Wheel Size	7"PU
Rear Wheel Size	8"PU

ICE SYMBOLS	
	Caution, attention or consult accompanying documents.
	Alternating Current
	Type BF Equipment
	Double Insulation
	No Smoking or Open Flames

Degree of protection against ingress of water is rated as IPx0.

Serialization format for products	
	
1	The first digit is the last one digit of the year for manufacture.
2	The second and third digits are the month for manufacture.
3	The fourth to seventh digits are counting of how many units were manufactured during the month.

■ Merits Limited Warranty

Merits Corporation warrants to the original purchaser of this wheelchair product that it is free of defect in material and workmanship and that, when operated within the guidelines and restrictions of this manual, will remain so free of defect in material and workmanship for a period of One (1) year from the original date of purchase.

Excluded from this warranty is failure due to negligence, abuse, accident, operation outside of rated limits, commercial or institutional use, damage / wear to upholstery or tires and improper maintenance or storage. The batteries for this wheelchair product are not supplied by Merits Corporation; contact the battery manufacturer / supplier if warranty replacement is requested.

This wheelchair product must not be modified in any way without the express written consent of Merits Corporation. Any such unauthorized modification could cause unreliable and / or unsafe operation and will void this warranty.

Where a failure occurs within the 1- year warranty period that is not excluded above, the failed components will be replaced with similar new or reconditioned components at Merits sole option. Merits Corporation will not be responsible for labor and / or shipping charges.

The foregoing warranty is exclusive and in lieu of all other warranties expressed or implied including, but not limited to, the implied warranty of merchantability and fitness for a particular purpose. Merits Corporation will not be liable for any consequential or incidental damages whatsoever.

■ **Warranty Registration**

**MERITS HEALTH PRODUCTS INC.
WARRANTY REGISTRATION**

MODEL NO. _____

SERIAL NO. _____

DATE PURCHASED _____

NAME _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

DEALER NAME _____

STAMP

RETURN ADDRESS

We wish you a safe and comfortable riding experience!





MERITS HEALTH PRODUCTS, INC.

www.meritshealth.com 730 N.E. 19th Place, Cape Coral, FL 33909
TEL: 1-800-963-7487, 1-239-772-0579 FAX: 1-239-574-2661