DE

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EN

Specifications:

13 AM pelvic trainers are designed to help you correctly locate and target your pelvic floor muscles, allowing you to exercise them effectively. The balls are inserted into the vagina and whilst your body holds and clenches around them, you are effectively working out and strengthening the pelvic floor muscles. The set contains 6 different weights and 2 sizes, allowing you to gradually train your pelvic floor muscles to become stronger. You can begin with the lightest trainer, and gradually progress onto larger and heavier balls.

Material
 String: 100% Silicone.
 Ball: Steel covered with ABS.
 Size: 3pcs - 89 mm x 34mm, 3pcs - 39 mm x 34 mm

18g 31g 57,6g 57,6g
 18g

33,4g 46,4g 73g 104,2g 117,2g 143,8g

EN

Instructions:

- Cleaning:** Before use, clean the Kegel ball with soapy water and rinse thoroughly.
- Usage:** Gently insert the Kegel ball against your vaginal opening until it comfortably slides in. Insertion may be easier while lying down or squatting. For removal, bear down on pelvic muscles and gently pull the ball by using the string gently until it slides free. Extra lubricant can aid in the removal process.

Beginners: While lying down, squeeze pelvic floor muscles around the balls for 5-10 seconds. Relax, take a short rest, and repeat the squeeze. Aim for 3 sets of 10 repetitions or as many as comfortable.

Intermediate/Advanced: Try holding the ball in place while standing, a gentle, prolonged hold of pelvic floor muscles for 5 to 10 minutes. If the weights slip, gently push them back into place with your fingers. With practice, they will be less prone to slipping.

- Removal:** Get comfortable in your preferred position, lubricate the vaginal opening, and gently pull the string to ease the Kegel ball out, minimizing discomfort.
- Cleaning and Storage:** Run the Kegel ball under hot water, wash with antibacterial soap, and rinse the balls thoroughly to eliminate any soap residue. Dry with a clean towel or allow them to air dry for a few hours before storage.

CAUTION!
 Muscles can fatigue quickly; 10-15 minutes use is recommended. Use high quality lubricants and make sure you are not allergic to the ingredients. Wash and clean the product before and after use.