

Materials:

- This sheet printed on 120gsm dark coloured paper
- Diffraction grating sheet
 Edmund Optics 54.509, cut long
 side in nine pieces each 34mm,
 cut short side in six pieces each
 25mm
- Plastic sheet 28 x 38 mm
- Sticky tape
- Black electrical tape

Tools:

- Scissors
- Smartphone or tablet, with built-in camera

Method:

- Make a 1mm slit from plastic sheet and two pieces of black electrical tape, 1mm apart
- Cut out spectroscope on solid line
- Cut out two rectangular holes
- Fold up spectroscope such that black surface is inside and fix with tape
- Fix diffraction grating over hole with sticky tape
- Fix slit over hole with sticky tape
- Place grating onto smartphone camera lens

Use:

- Point slit toward light source (fluorescent tube etc.)
- Adjust camera zoom
- View spectrum on smartphone
- Warning: Do **not** view the sun with your eye