
**BT103S**

**Variable-Speed Peristaltic Pump**

**Operating Manual**

# Contents

Safety Precautions .....	1
1 Description .....	6
2 Functions and Features.....	6
3 Components and Connectors.....	8
4 Display Panel and Operating Keypad.....	9
4.1 Keypad.....	9
4.2 LCD Screen Display.....	10
5 Parameter Settings .....	12
6 External Control Interface.....	15
7 Operation Instructions .....	15
7.1 Before Operating.....	15
7.2 Power Connection.....	16
7.3 Change Control Mode .....	16
7.4 Internal Control Mode.....	16
7.5 External Control Mode.....	17
7.6 Time Dispense Mode.....	19
7.7 Logic Level 1 Control Mode (footswitch).....	20
7.8 Logic Level 2 Control Mode (footswitch and direction switch) .	21
7.9 Communication Mode .....	23
7.10 Set Speed .....	24
7.11 Wi-Fi Network Setting.....	24
8 Maintenance .....	25
8.1 Warranty.....	25
8.2 Regular Maintenance .....	25
8.3 Malfunction Solutions .....	26
9 Dimensions .....	27
10 Naming Rule .....	28
11 Specifications .....	28

## Safety Precautions

### **Danger**

- Use the correct voltage indicated on the rating plate label of the pump to avoid any damage.
- Do not make any unauthorized dismantling, changes or modifications to the pump which could result in malfunctions or even potential accidents.
- Turn off the pump drive before installing or removing tubing, attempting any maintenance, cleaning or repair of the drive, connecting or disconnecting external control devices or a communication interface. Fingers or loose clothing could get caught in the drive mechanism.

### **Warning**

- Ensure no chemical reactions occur between the handled fluid with the material of the pump head and tubing before use.
- Tubing should be checked regularly to avoid breakage. Tubing breakage may result in fluid being sprayed from the pump. Use appropriate measures to protect the operator and equipment. The operator is solely liable for damages resulting from tubing breakage, particularly the leakage of toxic or valuable liquids.
- The pump is provided with a grounded plug which must be well grounded at all times.
- This device is not designed for nor intended for usage in patient-connected applications, including but not limited to medical and dental use.
- Observe all other applicable regulations concerning working safety, operational safety, environmental protection and relevant local regulations.

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# 1 Description

The BT103S variable-speed peristaltic pump offers basic functions such as start/stop, speed control, prime and reversible direction. It was upgraded with an LCD display, a time dispense mode and an anti-drip function. With the MODBUS RS 485 interface, it can be easily controlled by external devices such as PC, HMI or PLC.

**BT103S flow range 0.0001-480 mL/min, speed 0.1-100 rpm**

# 2 Functions and Features

Peristaltic pumps can handle fluid that is particularly abrasive, corrosive or viscous. There are no seals in contact with the medium pumped and no valves to clog. The inner surfaces are smooth and easy to clean. The fluid contacts only the tubing or tube material. Suction lift and priming can be up to 8m water column at sea level. It can handle the most shear sensitive of fluids like latex or firefighting foam with low shearing. It is capable of running dry and pumping fluids with high quantities of entrained air, such as black liquor soap. The high volumetric efficiency allows operation in metering or dosing applications where high accuracy is required. Tubing and tube materials are available for food and pharmaceutical use.

- LCD displays speed and working mode
- Membrane operating keypad.
- Time dispense mode: automatic dispense by setting the duration for each dose, lag time between doses and the number of cycles.
- High-precision speed control
- Optional external analog signal controls speed, logic level signal controls start-stop, and directions. Control signals are

optically isolated.

- The RS485 MODBUS interface allows easy control through communication with external devices.
- Internal double-layer isolation structure; circuit board with conformal coating makes the pump dust- and moisture-proof.
- Anti-electromagnetic interference feature; wide input voltage range for complex power environment.
- Plastic housing, easy to clean

## 3 Components and Connectors


Figure 1. Components and Connectors

## 4 Display Panel and Operating Keypad


Figure 2. Display Panel

### 4.1 Keypad

 **UP Key.** When pressing it shortly, the last digit of the value will increase by 1. Hold the key to increase the value quickly.

 **DOWN Key.** When pressing it shortly, the last digit of the value will decrease by 1. Hold the key to decrease the value quickly.

 **MENU Key.** When on the main screen, press the **MENU** key to enter the setting menu. When on the setting menu, press the **MENU** key to return to the main screen. When the drive is running, this key is disabled.

 **MODE key.** When the drive is not running, use the **MODE** key to change the working mode: Internal Control mode, External

Control mode, Time Dispense mode, Logic Level 1 control mode or Logic Level 2 Control mode.

**PRIME** key. Press the key to run the pump at the maximum allowed speed in the direction shown on the display. Press again to return to the previous state.

**DIRECTION** Key. Press to key to change the drive rotating direction, clockwise or counterclockwise.

**START/STOP** key. Control the start and stop of the motor.

When on the setting menu, press this key to enter the submenu.

## 4.2 LCD Screen Display

**Data entry:** When the drive is not running, press the number to input desired value in the pop-up window. See the picture below.


Figure 3. Display screen

### 4.2.1 A. Control Mode

It displays the current control mode. Press the **MODE** key to switch the control mode: Internal Control mode, External Control mode, Time Dispense mode, Logic Level 1 control mode or Logic Level 2 Control mode.

- Internal Control Mode: Use the keypad to operate the pump. Use an optional external logic signal to control the start and

stop.

- External Control Mode: Use an external analog signal to control rotation speed. Use an external logic signal to control direction, start and stop. The keypad is disabled.
- Time Dispense Mode: Dispense fluid automatically by setting the duration for each dose, the time between doses and the number of cycles.
- Logic Level 1 Control Mode (footswitch): Use an external logic signal to control the start and stop. Use the keypad to control direction and speed.
- Logic Level 2 Control Mode (footswitch and direction switch): Use an external logic signal to control start, stop and direction. Use the keypad to adjust speed.

#### 4.2.2 B. Running State

It shows the current drive running state.

*Figure 4. Running State*

#### 4.2.3 C. Direction State

It shows the current drive direction setting.

*Figure 5. Direction*

#### 4.2.4 D. Wi-Fi Signal Strength

Displays the current Wi-Fi signal strength (Wi-Fi function not available for US and European markets).

Press the **START/STOP** key to confirm and back to the setting menu screen. To return to the main screen, press and hold the **MENU** key or press the **MENU** key a few times until it goes back to the main screen.


Figure 9. Parameter Setting Flow Chart

**Note 1. Run Time Unit:** Time unit for the dispense duration of Time Dispense mode. The unit can be in days, hours, minutes or seconds.

**Note 2. Disp Time:** Dispense time for Time Dispense mode. It is the dispense duration for each dose. The range is 0.1-999 seconds/minutes/hours/days.

**Note 3. Pause Time Unit:** Time unit for the interval time of Time Dispense mode. The unit can be in days, hours, minutes or seconds.

**Note 4. Pause time for dispense mode:** It is the interval time between doses when the number of cycles is set to more than 1. The range is 0.1-999 seconds/minutes/hours/days.

**Note 5. Cycle:** Number of cycles under the Time Dispense mode. The range is 0-999 cycles. If it is set to 0, the dispensing process will keep running until **START/STOP** is pressed. Otherwise, the pump will only stop when it finishes dispensing the set number of cycles (*Figure 17*).

**Note 6. Reverse Angle:** To minimize the drip after a dispense, the drive can reverse direction to draw the fluid back at the end of the tubing. The reverse angle range is 0-720 degrees. When it is set to 0, the anti-drip function is disabled.

**Note 7. Max Speed:** Maximum speed for External Control mode. It is the maximum speed that the external analog signal can control.

**Note 8. Pump No:** It is the pump communication address for communication mode. When the number is changed, please restart the pump to apply the setting.

**Note 9. Beep Setting:** Set the key tone on or off.

**Note 10. Language:** System language setting, English or Chinese.


**Note 11. WIFI Reset:** Reset WIFI to link the drive to APP.

**Note 12. Contrast:** LCD backlight contrast setting.

**Note 13. Reset tubing life:** after replacing a tubing, the tubing life needs to be reset and the tubing life will be recounted.

## 6 External Control Interface

DB15	Mark	Note
1	ADC_W	Positive of external analog input
2	B	Communication interface, B pole of RS485
3	A	Communication interface, A pole of RS485
4	VCC_W	External DC power input
5		
6	CW_W	External input signal to control the direction
7	PWM	Pulse output
8	COM	Ground of external power
9	AGND	Negative of analog signal input
10	+12V	Positive of internal +12V power source
11	GND	Ground of Internal power source
12	CW	Internal direction signal output
13	RS_W	External start/stop signal input terminal
14	PWM_W	External pulse signal input
15	RS	Internal start/stop signal output

## 7 Operation Instructions

### 7.1 Before Operating

- 1) Please check the packing slip to ensure all parts are included and intact in the package. If there is a problem, please contact the manufacturer or distributor.
- 2) Read through the instruction.
- 3) A minimum of 200 mm space from the back of the pump should be maintained when it is in operation.

## 7.2 Power Connection

The voltage of the power supply should match what is indicated on the rating plate label of the pump. Please plug the power cord into the IEC Power Connector on the rear of the pump and plug the opposite end of the power cord into an electrical outlet. Flip the power switch located on the rear of the pump.

## 7.3 Change Control Mode

Turn on the power switch. The display will show a welcome message and then show the main screen. Press the **MODE** key to change the working mode (*Figure 10*).


Figure 10. Change Working Mode

## 7.4 Internal Control Mode

Use the keypad to operate the pump. Use an optional external pulse signal to control the start and stop.

- 1) Turn on the power switch. The pump will display the main screen.

- 2) Press the **MODE** key to change the mode to Internal Control mode (Int shown on the screen).
- 3) Press the **UP** or **DOWN** key to adjust the speed.
- 4) Press the **DIRECTION** key to change the rotating direction.
- 5) Press the **START/STOP** key to start or stop the pump.
- 6) Press the **PRIME** key, the pump will run at the maximum allowed speed.


Figure 11. Normal vs Full Speed

## 7.5 External Control Mode

In this mode, use an external logic level signal to control direction, start and stop; use an external analog signal to control rotation speed. The keypad is disabled.


Figure 12. External Control Mode

To control the pump by external signal:

- 1) Switch the power off. Wire the DB15 connector as shown in [Figure 13](#) or [Figure 14](#). Connect it to the DB15 port on the rear of the pump.


Figure 13. DB15 Wiring with External 12VDC Power Source


Figure 14. DB15 Wiring with Internal 12VDC Power Source

- 2) Turn on the power switch. The pump will display the main screen.
- 3) Press the **MODE** key to change the mode to External Control mode (ExtV or ExtI shown on the screen).
- 4) Close the external RS\_W switch, and turn on the external analog signal power source. The speed will change according to the intensity of the input signal. Open the RS\_W to stop the drive.
- 5) Open the CW\_W switch, the pump will run in a clockwise

direction; close the CW\_W switch, the pump will run in a counterclockwise direction.

**Note:** The external DC power source can be 5V or 12V. If it is 24V, a 1.5K resistor is needed to protect the internal circuit.



Figure 15. DB15 Wiring with External 24VDC Power Source

## 7.6 Time Dispense Mode

The pump will dispense fluid automatically after setting the duration for each dose, pause time between doses and number of cycles. When dispensing, the display will show the dispensing time or lag time that has passed (1.2s on the figure below), and the total cycles (4 on the figure) that has dispensed.

Dispense Stopped

Dispense Running

Figure 16. Time Dispense Mode

To set Time Dispense mode:

- 1) Turn on the power switch. The pump will display the main screen.
- 2) Press the **MODE** key to change the mode to Time Dispense

mode (Disp shown on the screen).

- 3) Press the **MENU** key to enter the Setting menu.
- 4) Set time unit and duration for each dose, pause time between doses and set the number of cycles.
- 5) Return to the main screen.


Figure 17. Dispense Cycle Setting

To run Time Dispense:

- 1) Press the **DIRECTION** key to set the running direction, clockwise or counterclockwise.
- 2) Press the **START/STOP** key to start dispensing.
- 3) During operating, press the **START/STOP** key to stop the drive anytime when necessary.
- 4) A footswitch can be used to start/stop the pump.

## 7.7 Logic Level 1 Control Mode (footswitch)

Use an external logic level signal to control the start/stop. Use the keypad to control direction and speed.

Figure 18. Logic Level 1 Control Mode

- 1) Switch the power off. Wire the DB15 connector as shown in [Figure 19](#) or [Figure 20](#), and connect it to the DB15 port on

the rear of the pump.


Figure 19. Logic Level 1 Control with Internal 12V Power Source


Figure 20. Logic Level 1 Control with External 12V Power Source

- 2) Turn on the power switch. The pump will display the main screen.
- 3) Press the **MODE** key to change the mode to Logic Level 1 control mode (shows lev1 on the screen).
- 4) Press the **UP** or **DOWN** key to adjust the speed.
- 5) Press the **DIRECTION** key to change the rotating direction.
- 6) When the switch is closed, the drive will start running; when the switch is open, the drive will stop.

Note: Use this mode to work with a TIME CONTROLLER.

## 7.8 Logic Level 2 Control Mode (footswitch and direction switch)

Use an external logic level signal to control start, stop and direction. Use the keypad to control speed.

Figure 21. Logic Level 2 Control Mode

- 1) Switch the power off. Wire the DB15 connector as shown in [Figure 22](#) or [Figure 23](#), and connect it to the DB15 port on the rear of the pump.


Figure 22. Logic Level 2 Control with Internal 12V Power Source



Figure 23. Logic Level 2 Control with External 12V Power Source


- 2) Turn on the power switch. The pump will display the main screen.
- 3) Press the **MODE** key to change the mode to Logic Level 2 control mode (shows lev2 on the screen).
- 4) Press the **UP** or **DOWN** key to adjust the speed.
- 5) When the switch RS\_W is closed, the drive will be running at the set speed; when the switch RS\_W is open, the drive will stop.
- 6) When the switch CW\_W is open, the drive will be running in a clockwise direction. When the switch CW\_W is closed, the drive will be running in a counterclockwise direction.

## 7.9 Communication Mode

The RS485 interface supports standard MODBUS protocol. The pump can communicate with external devices via the communication port. Please refer to the [Communication Instruction manual](#) for the parameters and supported commands.

Communication Disconnected

Communication Connected  
The pump number is set to 1.

Figure 24. Communication Mode

- 1) When the power is off, wire the DB15 connector as shown in [Figure 25](#), and connect it to the DB15 port on the rear of the pump. An external DC power source is recommended to avoid electrical interference.


Figure 25. RS485 MODBUS Wiring

- 2) Turn on the power switch. The pump will display the main screen.
- 3) Press the **MODE** key to change the mode to Internal Control mode or Time Dispense mode.
- 4) Control pump with the communication interface.
- 5) Press the **START/STOP** key to stop the drive anytime.

## 7.10 Speed setting

On the main screen, the speed resolution is 0.1 rpm when the speed is between 0.1 and 100 rpm; the speed resolution is 1 rpm when the speed is above 100 rpm. The speed can be set by pressing the **UP** or **DOWN** key. Press the **UP** or **DOWN** key shortly, the last digit of the value will change by 1. Hold the **UP** or **DOWN** key to change the speed continuously and rapidly.

Hold the **PRIME** key, then press the **UP** key to set the speed to maximum directly; hold the **PRIME** key, then press the **DOWN** key to set the speed to the minimum (0.1 rpm) directly.

## 7.11 Wi-Fi Network Setting

When the device is powered on, the display shows that Wi-Fi is not connected (Wi-Fi function is not available for US/ Europe Markets).


Figure 26. Wi-Fi Disconnected

To set up the Wi-Fi connection

- 1) Connect your smartphone/pad to the current Wi-Fi network and log in to your cloud control APP account.
- 2) In the upper right corner of the cloud control home page, click the "+Add Device" button and select "Wi-Fi Configuration" to add devices. The app will automatically search for available Wi-Fi. The added device will show in about 20 seconds after the Wi-Fi password was given and "Confirm" was pressed.


Figure 27. Wi-Fi connected

## 8 Maintenance

### 8.1 Warranty

The product comes with one-year labor and parts warranty. The limited warranty does not cover any damage that is caused by improper usage and handling.

### 8.2 Regular Maintenance

- Always check the tubing and connections to avoid leakage.
- Do not cover the fan on the rear of the pump.
- Do not use water to wash the pump. Keep the pump head dry.
- Do not use chemical solvents to clean the pump and pump head.

## 8.3 Malfunction Solutions

No.	Malfunction	Description	Solution
1	Hardware	No display	<ol style="list-style-type: none"> <li>1. Check the power cord</li> <li>2. Check the fuse. If it was blown, replace it with a 0.5A slow-blow fuse</li> <li>3. Check the internal power cord connection inside the pump.</li> </ol>
2	Hardware	Motor does not work	<ol style="list-style-type: none"> <li>1. Check the indicator of the driver board.</li> <li>2. Check the wire connection between the motor and the driver board.</li> <li>3. Check the wire connection between the driver and the main board.</li> <li>4. Check the power voltage for the pump.</li> </ol>
3	Hardware	Motor only runs in one direction	Check the connection between the drive board and the main control board.
4	Hardware	Keypad does not work	<ol style="list-style-type: none"> <li>1. Check the wire connection between the keypad and the main board.</li> <li>2. Check if the key is broken.</li> </ol>
5	Hardware	External control does not work	<ol style="list-style-type: none"> <li>1. Check the wiring of the connector.</li> <li>2. Check if the external control power voltage is provided.</li> <li>3. Check the connections of the external control board.</li> </ol>
6	Hardware	RS485 com does not	<ol style="list-style-type: none"> <li>1. Check the wiring of the connector.</li> </ol>

		work	<ol style="list-style-type: none"> <li>2. Check if the external control power voltage is provided.</li> <li>3. Check the connections of the communication board.</li> </ol>
7	Hardware	Noisy when running	Check the screws and level on the pump head to make sure they are secure.
9	Software	External control does not work	Check if the pump is on External Control Mode.
10	Software	RS485 does not work right	<ol style="list-style-type: none"> <li>1. Check if the display shows the communication is ready.</li> <li>2. Reset the address of the pump.</li> <li>3. Check whether on the bus there are two pumps using the same address</li> </ol>

If the problem cannot be solved, please contact the manufacturer or distributor.

## 9 Dimensions



Figure 28. Dimensions (mm)

## 10 Naming Rule

## 11 Specifications

Suitable pump heads	YZ15, YZ25, YT15, YT25, DG1, DG2, DT10-18, DT10-28
Speed resolution	0.1-100 rpm, 0.1 rpm resolution; 100-600 rpm, 1 rpm resolution
Speed accuracy	0.5%
Power supply	AC 90-240V, 50/60Hz
Power consumption	40W
External logic control signal	5V, 12V (standard), 24V (optional)
External analog control signal	0-5V, 0-10V, 4-20mA
Communication interface	RS485 MODBUS, Wi-Fi
Operating condition	Temperature 0~40°C Relative humidity <80%
IP grade	IP31
Display	132x32 LCD
Dimensions (LxWxH)	232x140x145mm
Weight	2.9 kg

**Golander LLC**

4405 International Blvd  
Ste B117, Norcross, GA 30093  
USA

Tel: +1 678-587-8806

[info@golanderpump.com](mailto:info@golanderpump.com)

[www.golanderpump.com](http://www.golanderpump.com)

**Golander GmbH**

Dechant-Heimbach-Str. 29  
53177 Bonn  
Germany

Tel: +49 228 50446952

[info@golander.de](mailto:info@golander.de)

[www.golander.de](http://www.golander.de)