PROSTATE AND PSA

The Prostate Specific Antigen (PSA) is a protein produced by the epithelial cells of the prostate gland. Its function is to fluidify the ejaculate and therefore increase sperm mobility. PSA also plays a role in dissolving cervical mucus, allowing sperm to enter the uterus.

Small quantities of PSA are normally found in the blood of men with a healthy prostate. Prostate pathologies (such as prostatitis, benign prostatic hyperplasia, cancer) lead to an increase in the blood levels of PSA.

PSA levels also increase naturally with age, after prostate surgery or after specific diagnostic tests (rectal examination, prostate biopsy, prostate ultrasound scan, etc.).


• For men who do not have symptoms, are aged over 45, who wish to check their prostate;

• For men with urinary tract symptoms (difficulty urinating, pain, increased frequency).

- REFERENCES**
1. Catalona WJ, Southwick PC, Slawin K.M., Partin A.W., Brauer M.K., Flanigan R.C., Patel A., Richie J.P., Walsh R.C., Scardino R.T., Lange R.H., Gostor G.H., Loveland K.G., Bray K.R.: Comparison of percent free PSA, PSA density, and age-specific PSA cutoffs for prostate cancer detection and staging. *Urology*, 2000
 2. Fritsche H.A., Babayan R.J.: Analytical performance goals for measuring prostate-specific antigen. *Clin Chem*, 1993 Jul;39(7):1525-29.
 3. Catalona WJ., "Prostate Cancer Screening" *Med Clin North Am*, 2018 Mar;102(2):199-214.

WHO IS THE PROSTATE PSA TEST FOR?

• The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

* Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).

* Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.

* The desiccant packet must not be used. Dispose of it with household waste without opening it.

PRECAUTIONS

1. Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
2. Keep the Test out of the reach of children.
3. Do not use the Test after the expiry date or if the package has been damaged.
4. Follow the Test components at a temperature between +4°C and +30°C. Do not freeze.
6. Use the test and lancing device once only.
7. The test is for external use only. DO NOT SWALLOW.
8. In vitro diagnostic device for individual use.
9. Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
10. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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Phoenix Innovative Healthcare Manufacturers Pvt. Ltd.
El-209, Shri Mahape Road,
TC Industrial Area, Mahape,
Navi Mumbai 400 710 MH | India

Hualan Tianda Medical Instruments Co., Ltd.
No. 106, East Songjiang Road, Huaiyin
Economic & Technological Development Zone
223002 Hualan City, Jiangsu - China

Shanghai International
Holding Corp.GmbH
(Europe), Eiffelstrasse 80,
20537 Hamburg, Germany









Advena Ltd.,
2nd Flr,
Tower Street,
Swatar, BKR 4013, Malta

Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
+41 91 61 605 1030
primahometest.com

800080_FU_11-PH_3.0 04/2021 REF 800080

Registered trademarks.

ENGLISH TEST PROCEDURE

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag.
3) Carefully rotate the protective cap of the sterile lancet 360° without pulling it.
4) Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended).
5) Holding the hand down, massage the finger until a large drop of blood forms.
6) Take the pipette without pressing the bulb.
7) Place the blood collected with the pipette into the well indicated on the cassette (S) by pressing the pipette bulb.
8) Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on).
9) Wait 10 minutes and read the results referring to the next section for result interpretation.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES

The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line. This result means that the value of Prostate-Specific Antigen (PSA) in the blood is higher than normal, so you must consult a physician.

NEGATIVE RESULTS

A coloured band appears only under the C (control) sign. This means that the value of Prostate-Specific Antigen (PSA) in the blood is normal. A negative result does not exclude the possibility of prostatic pathologies.

NON VALID RESULT

No bands appear or there is a line only under the T (test) sign and not under the C (control) sign. In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test with a new blood sample.

F.A.Q. - QUESTIONS AND ANSWERS

HOW DOES THE PROSTATE PSA TEST WORK?

The Prostate-Specific Antigen is a protein secreted normally by the prostate and circulating in our body in small quantities. Values above normal may be a sign of the presence of pathological conditions in the prostate, which may be of infectious or inflammatory origin (prostatitis), benign (benign prostatic hyperplasia) or malignant (cancer). The PROSTATE PSA test is a rapid qualitative test for measuring the levels of Total PSA in whole blood. The test uses a unique combination of specific antibodies that can selectively identify PSA.

The diluent is accidentally spilt, if the amount of blood is not enough or the number of drops dispensed in the well is more than 3.

HOW DO I INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT? The colour and the intensity of the bands are not important for the interpretation of the result. The bands must be full and homogeneous. The test is positive regardless of the intensity of the colour of the test band.

IS THE RESULT RELIABLE IF READ AFTER 10 MINUTES? No. The test must be read 5 minutes after the procedure is completed.

WHAT SHOULD I DO IF THE RESULT IS POSITIVE? If the result is positive, the concentration of PSA in the blood is above normal levels and you must consult a physician, who will establish what to do next. PSA IS A SPECIFIC PROSTATE MARKER BUT IT DOES NOT NECESSARILY IMPLY THE PRESENCE OF NEOPLASIA.

WHAT SHOULD I DO IF THE RESULT IS NEGATIVE? If the result is negative, the concentration of PSA is normal. A negative result does not exclude the possibility of prostatic pathologies. If the result is negative but you have urinary tract symptoms (difficulty urinating, pain, increased frequency) you should talk to your physician.

HOW ACCURATE IS THE PROSTATE PSA TEST? The Test is very accurate. Assessment reports show a concordance rate above 90% (CI 95%: 85.67-95.67) with reference methods.

CAN THE RESULT BE INCORRECT? The result is correct if the instructions are followed carefully. However, the result may not be correct if: the device comes into contact with other liquids before being used, if part of

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