
Vac Star

DENTAL VACUUM SYSTEM

PART NUMBERS: VS20, VS40, VS50, VS50H, VS80 AND VS80H



USER'S MANUAL





TABLE OF CONTENTS

SECTION		PAGE
Congratula	itions	3
Purpose o	f this Manual	3
Warranty .		4
On-Line W	arranty Registration	4
Safety Ins	ructions	4
Sizing Gui	de	5
Optional Ir	-Line Filter Kit	5
Operating	Information	6
Key Parts	Identification - Single Pumps	6
Key Parts	Identification - Dual Pumps	7
Installation	Information	8
Product S	pecifications/Dimensions	11
Site Requi	rements	11
Replacem	ent/Reorder	12
Troublesh	poting	13
Accessorie	es/Options	13
Maintenan	ce	15
LIST OF	ILLUSTRATIONS	
FIGURE	TITLE	PAGE
1	VacStar 20, 40 Parts Location	
2	VacStar 50, 50H, 80 and 80H Parts Location	
2A		
3	Plumbing Connections Detail View	/
4	Plumbing Connections Detail View	
	Bypass Valve Educator Assembly	7
5	Bypass Valve Educator Assembly	
5 6	Bypass Valve Educator Assembly	
	Bypass Valve Educator Assembly	
6	Bypass Valve Educator Assembly	
6 7	Bypass Valve Educator Assembly	
6 7 8	Bypass Valve Educator Assembly	
6 7 8 9	Bypass Valve Educator Assembly	
6 7 8 9 10	Bypass Valve Educator Assembly	

Congratulations on the purchase of your new **VacStar**[™] Dental Vacuum System hereafter referred to as **VacStar**[™] in this manual.

Your VacStar™ has been engineered to deliver maximum air flow at the ideal vacuum level without creating traumatic suction pressure that could harm patients' delicate tissue. The VacStar™ is a water ring pump that produces consistent high-volume air flow, even with multiple users on-line. The balanced, corrosion free bronze impeller minimizes noise and a patented vacuum relief valve monitors and maintains constant uniform vacuum pressure. A capacitor-start type motor, with a highly reliable electrical contactor and powerful transformer can be depended on to start every time. The VacStar™ is designed with everything accessible from the front, including the easy to replace solids collector.

If your **VacStar**[™] comes with an integral HydroMiser (VS50H or VS80H), water consumption will be reduced by up to 75%. If not, a HydroMiser can be integrated into your **VacStar**[™] at a later date. The HydroMiser separates the liquid and gas discharge from the operatories. The gases are vented out and the liquid and its particulates are directed down the drain. The clean water extracted during this separation process is directed back toward the **VacStar**[™] where it is mixed with fresh water and then directed into the pump chamber to create vacuum. This efficient reuse of water reduces the VacStar's fresh water consumption.

Thousands of dentists have depended on the **VacStar**[™] since 1987. Now that your practice has a **VacStar**[™], or a **VacStar**[™] with the water saving HydroMiser, you too can depend on constant, uniform delivery of vacuum to your operatories and proven trouble-free operation.

PURPOSE OF THIS MANUAL

This manual provides installation, operation and maintenance instructions for the support of the six available **VacStar**™ systems listed below. Review and follow the guidelines included in this User Manual to ensure that the system provides the highest level of service.

VS20 2 User Capability VS40 3 User Capability VS80 VS80

Note: VS50H and VS80H models include an integral HydroMiser that reduces fresh water consumption

The **VacStar**[™] is warranted to be free from defects in material and workmanship from the date of installation for a period of twenty-four (24) months.

Any item returned to our factory through an Air Techniques Authorized Dealer, will be repaired or replaced at our option at no charge provided that our inspection shall indicate it to have been defective. Dealer labor, shipping and handling charges are not covered by this warranty.

This warranty does not apply to damage due to shipping, misuse, careless handling or repairs by other than authorized service personnel. Warranty is void if equipment is installed or serviced by other than dealer service personnel authorized by Air Techniques. Air Techniques, Inc. is not liable for indirect or consequential damages or loss of any nature in connection with this equipment.

This warranty is in lieu of all other warranties expressed or implied. No representative or person is authorized to assume for us any liability in connection with the sale of our equipment.

ON-LINE WARRANTY REGISTRATION

Quickly and easily register your new **VacStar**[™] on-line. Just have your product model and serial numbers available. Then go to the Air Techniques website, **www.airtechniques.com**, click the **register a product link** and complete the registration form. This on-line registration ensures a record for the warranty period and helps Air Techniques keep you informed of product updates and other valuable information.

SAFETY INSTRUCTIONS

Use of the **VacStar**™ not in conformance with the instructions specified in this manual may result in permanent failure of the unit.

WARNING: To prevent fire or electrical shock, do not expose this

appliance to rain in or moisture.

All user serviceable items are described in the maintenance section.

Manufacturing date code on serial number label is in the format Month YYYY.

ATTENTION USERS:



Alerts users to important Operating and Maintenance instructions. Read carefully to avoid any problems.



Warns users that uninsulated voltage within the unit may be of sufficient magnitude to cause electric shock.



Indicates the ON and OFF position for the Equipment power switch.



Indicates that the equipment complies with the Medical Device Directive 93/42/EEC.



MEDICAL ELECTRICAL EQUIPMENT

WITH RESPECT TO ELECTRICAL SHOCK, FIRE, MECHANICAL AND OTHER SPECIFIED HAZARDS ONLY IN ACCORDANCE WITH UL-60601-1, CAN/CSA C22.2 NO.601.1 66CA



Indicates the ON and OFF position for the Equipment power switch.



Medical Device Safety Service Schiffgraben 41 30175 Hannover, Germany Choosing the correct size **VacStar**[™] for your practice depends on the number of HVE (High Volume Evacuator) and SE (Saliva Ejector) users anticipated. To assure optimum vacuum, the vacuum demands should not exceed the number of HVE and SE users shown in the chart below:

Recommended Number of Simultaneous Users

VacStar 20 HVE's + SE's	VacStar 40 HVE's + SE's	*VacStar 50 & 50H HVE's + SE's	*VacStar 80 & 80H HVE's + SE's
2 + 0	3 + 0	4 + 0	7 + 0
1 + 1	2 + 2	3 + 2	6 + 1
0 + 4	1 + 4	2 + 4	5 + 3
	0 + 6	1 + 5	4 + 4
NOTES:	3 + 6		
HVE = High Volume E SE = Saliva Ejector	2 + 8		
* These combinations If only one pump is rur	1 + 10		
in only one pump is rui	0 + 13		

OPTIONAL IN-LINE FILTER KIT

Since larger quantities of particulates may occur initially when a VacStar is replacing another vacuum pump, an optional In-Line Filter is recommended to be installed at the intake connection (see Key Parts). This In-Line Filter is designed to collect larger quantities of particulates from the discharge BEFORE it flows into the VacStar. The larger quantities of debris is mainly due to the VacStar's increased pulling power and the effectiveness of the CleanStream Evacuation System Cleaner's ability to break down proteinaceous deposits and synthetic debris that have accumulated in the existing vacuum lines.

Use the In-Line Filter Kit P/N 55078 for single vacuum pump units VS20 and VS40. The kit part number for twin pump units (VS50, VS50H, VS80 and VS80H) is 55079. Refer to the Maintenance Section for recommended maintenance requirements.

□ AT THE START OF THE DAY

Always TURN ON THE WATER before TURNING ON THE POWER.

- ☐ The VacStar may be turned on/off from a single, convenient location within the dental suite using a Remote Control Panel (See Optional Accessories).
- □ The vacuum level is factory preset at 10 In Hg (inches of mercury). This is the reading on the gauge when all HVE's (High Volume Evacuator) and SE's (Saliva Ejector) are CLOSED. Should this setting be too high for your needs, contact your dealer to readjust the setting.
- ☐ It is recommended that the system run continuously during the day. However, the VacStar can be turned off if suction is not required for a period of 15 minutes or longer.
- ☐ If one pump is being operated at a time, it is important to alternate pumps on an every other day schedule so that the pumps are used evenly.

□ AT THE END OF THE DAY

Always TURN THE POWER OFF, then TURN THE WATER OFF.

KEY PARTS IDENTIFICATION - SINGLE PUMPS




Figure 1. VacStar 20 and 40 Parts Location

NOTE: VACSTAR 20 shown. VACSTAR 40 is similar except main power connection is made via provided hospital grade NEMA 6-15P line cord.




Figure 2. VacStar 50, 50H, 80 and 80H Parts Location NOTE: VACSTAR 50H Shown, Other Models are Similar




Figure 3.
Bypass Valve Educator Assembly




Figure 2A
Plumbing Connections Detail View

□ Plumbing (water) lines

- To assure that the **VacStar**™ provides optimum vacuum, incoming water pressure must be maintained between 20 and 100 psi.
- If heavy combinations of particulates exist in the incoming water, an in-line filter should be installed. (See Accessories/Options for the Remote Control Water Valve.) This will prevent the VacStar's water inlet filter from clogging too frequently.
- Incoming water temperature should be between 40° and 75°F.
- Water connection location is shown in Fig. 1 and 2a (water inlet connection).

☐ Suction

- For VacStar 20 and 40, suction hose is connected at suction intake, found on intake solids collector assembly. See Fig. 1.
- For VacStar twin pump units, suction hose is connected at suction intake, found on intake solids collector assembly. See Fig. 2.

Drain lines

- For VacStar 20 and 40 without a HydroMiser or an Air/Water Separator, see Fig. 4.
- For VacStars without a HydroMiser or an Air/Water Separator, the effluent should be discharged into an open drain or a closed vented drain. See Fig. 5.
- For VacStars with a HydroMiser (see Fig. 6) or an Air/Water Separator (see Fig. 7), gases should be vented out according to code. The waste water (with particulates) from the operatories can be discharged via an open drain or a closed vented drain.

Note: For VacStars without a HydroMiser, the drain may be up to 36" above the unit.




Figure 4. VacStar 20, 40 without a HydroMiser or Air/Water Separator

Figure 5. VacStar 50, 80 without a HydroMiser or Air/Water Separator

Figures 6, 7 and 8 show the typical installations of the dual **VacStar**[™] models 50, 50H, 80 and 80H. Install the **VacStar**[™] by referring to the figure corresponding to the system to be installed. Single **VacStar**[™] models 20 and 40 are installed in the same manner. Make sure that all notes, warnings and requirements are followed.





Figure 6. VacStar with Built-In HydroMiser

Figure 7. VacStar with Wall Mounted Air/Water Separator



Note: For all installations, vent to outside with 2-inch schedule 40 pipe.

WARNING:

CONDENSATION OF WATER WILL OCCUR IN VENT PIPING. AVOID ACCUMULATION OF WATER IN VENT, SLOPE PIPING TOWARD SEPARATOR.

■ Wall-mounted HydroMiser

If the existing drain is higher than the HydroMiser outlet, the HydroMiser must be mounted so that its outlet is above the drain. The HydroMiser can be installed up to 36" above the base of the VacStar with the HydroMiser Wall Mount Kit (#55087).

IMPORTANT NOTE: ALL INSTALLATIONS

Ambient temperature for all VacStar installations should be 40°- 104°F (5°- 40°C).

The liquid drain from the HydroMiser or an Air/Water Separator must slope downward at least 1/4" for every 10 feet of run toward the drain.

(Avoid local low sections, avoid creating traps in the line.)




Figure 8. VacStar with Wall-Mounted HydroMiser

□ Electric

- VacStar 20, 50, 50H, 80, 80H must be wired directly from an electrical box that complies with local electrical codes to the VacStar's Electrical Connection Box. See Figure. 10.
- All VS40 VacStars are wired with a provided hospital grade NEMA 6-15P line cord and requires a hospital-grade 6-15R receptacle.
- If the voltage falls below the minimum 105V or 205V or above the maximum 125V or 253V during operation, a Buck/Boost Transformer must be installed.






Figure 9. VacStar Electrical Junction Box - Interior View VacStar 20 (Factory Set for 230 V)

ALL INSTALLATIONS MUST CONFORM TO LOCAL CODES





Dual Circuit *				
Pumps	Power Leads			
RIGHT	(L1) BLACK (L2) WHITE			
LEFT	(L1) RED (L2) BLUE			

For Single Circuit, connect Black and Red wires together (L1) and White and Blue wires together (L2).

VacStar 20 Handy Box Connection

VacStar 50, 50H, 80, 80H Handy Box Connection

Figure 10. VacStar Electrical Box Connection

Spec/Dimension	VacStar Models						
ELECTRICAL	VS20	VS40	VS50	VS50H	VS80	VS80H	
Voltage Rating	*115/230	230	230	230	230	230	
Voltage Minimum/Maximum	*105/125 205/253	205/253	205/253	205/253	205/253	205/253	
Full Load Amps	*16/8	13.4	16	16	26.8	26.8	
WATER Inlet Water	20 - 100	20 - 100	20 - 100	20 - 100	20 - 100	20 - 100	
Pressure (PSI)							
Flow Rate Per Pump (gal/min) w/HydroMiser	0.12	0.18	N/A	0.12	N/A	0.18	
Flow Rate Per Pump (gal/min) w/o HydroMiser	0.50	0.75	0.50	N/A	0.75	N/A	
Water Temperature (°F)	40 - 75	40 - 75	40 - 75	40 - 75	40 - 75	40 - 75	
VACUUM LEVEL Preset at Factory (In Hg)	10	10	10	10	10	10	
SHIPPING WEIGHT (lbs)	68	85	160	170	200	210	
DIMENSIONS Inches (HxWxD)	14 x 11x 11	17 x 11 x 11	22 x 28 x 16	25 x 28 x 16	22 x 28 16	25 x 28 x 16	

^{*} VacStar 20 may be converted from 230 Volts (Factory Set) to 115 Volts at installation site. See Figure 9.

SITE REQUIREMENTS

Environment Conditions:

Operating Conditions Indoor use at altitudes up to 2000m

Temperature 5 to 40° C (41 to 104° F).

Maximum relative humidity 80% for temperatures up to 31° C,

decreasing linearly to 50% relative humidity at 40° C.

Supply voltage fluctuation of +/- 10% of nominal voltage.

ouppi) rollings intollinus or i, rollings.

IEC 60601-1 Classification Not suitable for use in the presence of a flammable anes-

thetics mixture with air or with oxygen or nitrous oxide.

Class I Installation Category

Ordinary equipment (IPXO). Does not protect against

ingress of water.

Unit is suitable for continuous operation.

Requirement	VacStar Models					
ELECTRICAL	VS20	VS40	VS50	VS50H	VS80	VS80H
Minimum Circuit Breaker Rating	30A(115V) 20A(230V	20A	2 ea. 15A or 1@. 30A	2 ea. 15A or 1 @ 30A	2 ea. 20A or 1 @ 40A	2 ea. 20A or 1 @ 40A
Wire Size AWG (Minimum Gauge)	10(115V) 12(230V)	12	2 ea. #14 or 1 @ #10	2 ea. #14 or 1 @ #10	2 ea. 12 or 1@ #8	2 ea. 12 or 1@ #8
PLUMBING						
Minimum CFM @ 0" Hg	16	22	32	32	44	44
Air Exhaust	2" schedule 40 pipe					
Ambient Temperature	40 - 104°F (5 - 40°C)					
Overhead Plumbing Main Line Minimum/Maximum Inside Diameter (inches)	1 / 1½	11/4 / 2	11/4 / 11/2	11/4 / 11/2	1½/2	1½/2
End Fitting Maximum	3/4" FNPT	3/4" FNPT	1" FNPT	1" FNPT	1" FNPT	1" FNPT
Riser Diameter Overhead Main Line	½" ID					
Floor Plumbing Main Line Minimum/Maximum Inside Diameter (inches)	1 / 1½	11/4 / 2	11/4 / 11/2	11/4 / 11/2	1½/2	1½/2
End Fitting Maximum	3/4" FNPT	3/4" FNPT	1" FNPT	1" FNPT	1" FNPT	1" FNPT
Branch Line Diameter Minimum/Maximum Inside Diameter (inches)	3/4 / 1½	1 / 1½	1 / 1½	1 / 1½	1 / 1½	1 / 1½

NOTE: Suction piping must slope at least a ¼" for each 10 feet of run towards the pump.

Use PVC Schedule 40 or Copper Type M.

ALL INSTALLATIONS MUST CONFORM TO LOCAL CODES

REPLACEMENT/REDRDER

DESCRIPTION	MODEL	PART NUMBER
Solids Collector Replacement Kit	VacStar 20, 40	55880 (3/4 inch)
	VacStar 50, 50H, 80, 80H	55094 (1 inch)
In-Line Filter Replacement Kit	All VacStar Models	55094

PROBLEM	POSSIBLE CAUSE	POSSIBLE SOLUTIONS
1. Low suction.	a. Water filter or solids collector clogged. b. Check valves are stuck.	 a. Clean/replace filter and/or collector. b. Use a system cleaner like CleanStream; turn vacuum on and off to free check valve. If valve remains stuck, call your authorized Air Techniques dealer for repair service.
	c. Low water pressure. d. HydroMiser water recycler is clogged.	c. Raise water pressure. d. Open bypass valve to run VacStar. Call your authorized Air Techniques dealer for repair service.
	e. HydroMiser clogged. f. Solenoids not operating.	e. Call your authorized Air Techniques dealer for repair service. f. Open bypass valve to run VacStar. Call your authorized Air Techniques dealer for repair service.
	g. Restricted air exhaust.	g. Check air exhaust pipe size to make sure it conforms to specification; check for and clear possible restrictions in air exhaust system.
2 No suction.	a. Pumps off.b. Pumps not running.c. Inlet check valves stuck closed.d. Water inlet filter and/or solids collector clogged.	 a. Turn pumps on. b. Call your authorized Air Techniques dealer for repair service. c. Call your authorized Air Techniques dealer for repair service. d. Clean/replace filter.
	e. Suction hose collapsed. f. Solenoids not operating. g. Water off.	 e. Hose needs to be replaced, call your authorized Air Techniques dealer for repair service. f. Call your authorized Air Techniques dealer for repair service. g. Turn water on via water inlet valve.
3. Excessive suction.	a. Relief valve stuck closed.b. Relief valve filter clogged.c. Relief valve set too high.	 a. Call your authorized Air Techniques dealer for repair service. b. Call your authorized Air Techniques dealer for repair service. c. Lower Relief valve settting.
4 Pumps do not run	a. Main circuit breakers off. b. Electrical problem.	a. Turn main circuit breakers on. b. Call your authorized Air Techniques dealer for repair service.
5 Noisy pumps.	a. Inadequate water supply.b. HydroMiser eductor clogged.c. Drain line collapsed.d. Solenoids not operating.	 a. Call plumber to improve water supply system. b. Call your authorized Air Techniques dealer for repair service. c. Hose needs to be replaced. Call your authorized Air Techniques dealer for repair service. d. Call your authorized Air Techniques dealer for repair service.

ACCESSORIES/OPTIONS

The following lists the description, part number, the applicable model and description for accessory components and options available to maintain and expand the **VacStar™** System to meet your professional needs. Contact your Air Techniques Dealer for information.

Description	Model	Part Number	
HydroMiser Wall Mount Kit	VacStar 50H, 80H	55087	
Remote Control Panels with 24V switches	VacStar 20, 40	53250 or 53251	
	VacStar 50, 50H, 80, 80H	53113 or 53149	
Remote Control Water Valve, with filter	All VacStar Models	53020 (24V) - 3 4" pipe	53020-1 (115V) - 3 4" pipe
		53170 (24V) - 1" pipe	53171 (115V) - 1" pipe
	VacStar 20	H-2	, , , ,
HudraMigar Kit	VacStar 40	H-4	
HydroMiser Kit	VacStar 50	56041	
	VacStar 80	56042	
Air/Water Separator	VacStar 20, 40, 50, 80	55540	
In-Line Filter Kit	VacStar 20, 40	55078 - 3 4" pipe	
III-LIIIE FIILEI KIL	VacStar 50, 50H, 80, 80H	55079 - 1" pipe	
CleanStream Evacuation System Cleaner	All VacStar Models	57660 Starter Kit	57640 1 Box of 32 Packets

Like all precision products, your **VacStar**[™] requires a certain amount of care on a regularly scheduled basis. A well-organized maintenance program aids dependable equipment operation and reduces problems to a minimum. Routine checks help to detect general overall wear, and replacement of parts can often be made before a problem occurs.

Consequently, we have established minimum maintenance requirements listed below that include routine inspections and the replacement of filters. Adherence to this recommended maintenance schedule will ensure that the equipment will continue performing at its best with uninterrupted service.

Daily Maintenance - Clean vacuum lines

Flush all vacuum lines and tubing in the dental system with CleanStream Evacuation System Cleaner.

□ Routine Inspection - Monthly

- 1. Check tubing for kinks or cracks.
- 2. Check for abnormal noises and leaks.
- 3. Dual Units: Make sure both motors are running
- 4. Check exterior surfaces for dirt and debris, clean if necessary.
- 5. Make sure that no flammable, corrosive, or combustible materials are stored in the equipment room (especially in the area around the equipment).
- 6. Refer to Figure 11 and check the vacuum relief valve filter, clean if necessary.




Figure 11. Vacuum Relief Valve Filter Location

Caution: Solids collector may contain biologically hazardous material. Wear protective gloves. Dispose of

waste in approved bio-hazard container.

Important: When a VacStar is replacing another vacuum pump, clean the collector DAILY during the first week

of operation since larger quantities of particulates may initially occur.

A worn or missing gasket and/or failure to tightly screw the bowl to the solids collector body will cause

poor suction due to air leakage.

DO NOT OPERATE THE VACSTAR WITHOUT THE SCREEN INSIDE THE FILTER BOWL.

Intake Solids Collector Replacement - Monthly

Refer to Figure 12 and using the replacement kit listed below for the specific **VacStar**™ models, replace the solids bowl, screen and gasket. Do the same replacement if using an optional in-line filter.

VacStar Model	Solids Collector Kit Part No.	Optional In-Line Filter Kit Part No.
VS20 & VS40	55880	55094
VS50, VS50H, VS80 & VS80H	55094	55094

Replacement Procedure

- 1. Turn off the power and water supply.
- 2. Unscrew the solids bowl (counter clock-wise) and remove the screen and gasket. Dispose of all three items.
- 3. Assemble a new bowl, screen and gasket included in the Solids Collector Replacement Kit.
- 4. Install the new solids collector by screwing the bowl into the solids collector body.




Figure 12. Intake Solids Collector Location

☐ Check/Clean Solenoid Water Inlet Strainer - Semi-Annually

Check the Inlet Strainer for dirt and debris by performing the following steps.

- 1. Turn off the power and water supply to the equipment.
- 2. Use a 1 3/16 inch wrench to unscrew (turn counter clockwise) the cover nut.
- 3. Remove the cover nut and strainer.







Figure 13. Water Inlet Strainer Location

Assembling the Water Inlet Strainer (All Models)

- 1. Orienting the assembly with the cover nut facing down as shown, seat strainer into the cover nut.
- 2. Insert the strainer up into the strainer body and tighten the cover nut.
- 3. Make sure the strainer stays perpendicular to the strainer body.
- 4. Push up and tighten the cover nut making sure not over tighten.





NOTES		



For over 50 years, Air Techniques has been a leading innovator and manufacturer of dental products. Our priority is ensuring complete satisfaction by manufacturing reliable products and providing excellent customer and technical support. Whether the need is digital imaging, utility room equipment or merchandise, Air Techniques can provide the solution via our network of authorized professional dealers. Proudly designed, tested and manufactured in the U.S., our products are helping dental professionals take their practices to the next level.

Air Techniques' family of quality products for the dental professional include:

Digital Imaging

- Digital Radiography
- Intraoral Camera
- Caries Detection Aid
- Intraoral X-ray
- Film Processors

Utility Room

- Dry Vacuums
- Wet Vacuums
- Air Compressors
- Amalgam Separator
- Utility Accessories
- Utility Packages

☐ Merchandise

- Surface Disinfectant
- Enzymatic Cleaner
- Hand Sanitizer and Lotion
- Waterline Cleaner
- Evacuation System Cleaner
- Imaging Accessories
- Chemistry
- Processor Accessories

Corporate Headquarters

1295 Walt Whitman Road | Melville, New York 11747- 3062 | Phone: 800-247-8324 | Fax: 888-247-8481

Western Facility

291 Bonnie Lane, Suite 101 | Corona, CA 92880 - 2804 | Phone: 800-247-8324 | Fax: 951-898-7646

www.airtechniques.com





