

Medical Grade CO<sub>2</sub> gas is used for insufflation. Gas is delivered to the insufflator from CO<sub>2</sub> cylinders or from a central supply wall source located in each operating room. A high pressure hose is easily connected to the supply source and to the rear of the insufflator.



After connecting the hospital grade power cord to rear of unit and to a properly grounded power source, turn on the main power switch at the rear of the unit as shown.




Activate power to the unit by touching the front of the LCD screen. Prior to activation, the LED on-off indicator located at the lower left of the screen is "blue".




Once the LCD screen is activated, the LED on-off indicator located at the lower left of the screen turns "green". The insufflator will now complete a very brief self-diagnostic check.



After the self-diagnostic check is completed, the "Select Operating Mode" screen will appear. The user will need to choose from 5 different modes: - General Laparoscopic - Bariatric - Pediatric - Vessel Harvest - Custom



Once the appropriate operating mode is selected, the user must confirm the selected mode.



The operating mode screen appears displaying preselected Pressure and Flow rates. Factory set defaults can be used or new user preferences can easily be stored using the "Set-Up Menu." Easy toggling through Verres, High and Maximum Flow settings are made using the "Flow Preset" icon.



Using the "Up and Down" arrows allows for pressure and flow adjustments at any time.




Prior to use, the disposable, filtered insufflation tubing set needs to be connected as shown.



If "Continuous Pressure Sensing" is required, the disposable, filtered CPS™ tubing set is connected. CPS™ is essential to insure rapid recovery of any lost intra-abdominal pressure, reducing the breathing effect of the abdominal surface and for totally uninterrupted pressure monitoring. CPS™ is highly recommended for pediatric laparoscopy.\*



Warming the CO<sub>2</sub> gas can help to reduce post-operative patient pain, reduce PACU recovery time and reduce fogging of the laparoscope.\* A reusable In-Line Gas Warmer is available and is easily connected to the front panel of the insufflator as shown.



The In-Line Warmer is connected 3-4 inches from the distal end of the tubing.




To initiate gas flow for insufflation, depress the green arrow as shown. The insufflator will deliver gas to maintain the desired pressure.




To stop the flow of gas, depress the orange "stop" button. The flow of gas is now paused. To start gas flow again, simply depress the green arrow.




The user can access the set-up menu at any time by depressing the SET-UP MENU icon located in the lower left corner of the screen.



The "Set-Up Menu" allows adjustments to: - 5 Alert Volume levels - 14 Language choices - Viewing Service Contact Information - Operating Mode Default - Gas supply options



Default settings can be changed and stored while in the "Program Defaults" screen. Restoring the factory defaults can be done by depressing the "Restore Factory Default" icon as shown.



When finished using the GS2000 insufflator product, depress the POWER DOWN icon located at the upper left side of the screen. To re-start the unit, see step #1.

### Specifications

#### 50L Insufflator Console (GS2000)

**Dimensions:** 12 in. (w) x 6 in. (h) x 12 in. (d)  
(305 x 153 x 305 mm)  
**Weight:** 11.5 lbs. (4.3 kg)  
**Input Voltage Range:** 90-240 VAC, 50-60 Hz  
**Protection Class:** Class I, Type B Applied part  
**Electrical Leakage:** 100 micro amps maximum

#### CO<sub>2</sub> Warmer (GS1012)

**Gas Temperature:** 95-106°F (non-adjustable) (35-40°C)  
**Dimensions:** 1.0 in. (d) x 3.4 in. (l)  
**Weight:** 2.0 oz. (56g)

### Certifications

UL 60601-1  
CAN/CSA C22.2 NO. 601.1  
IEC (EN) 60601-1 IEC (EN) 60601-1-2

CAUTION: This guide is not a replacement or substitute for the ConMed GS2000 operator's manual. Please be sure to read and understand the GS2000 operator's manual thoroughly before using the insufflator.

### Ordering Information

| Description                                       | Catalog Number | Description                                              | Catalog Number |
|---------------------------------------------------|----------------|----------------------------------------------------------|----------------|
| 50L Insufflator                                   | GS2000         | Insufflation Tubing Set, disposable, CPC/Luer, w/filter  | GS1016         |
| Hose and Yoke for CO <sub>2</sub> , E-Tank, 4 ft. | GS1052         | Insufflation Tubing Set, disposable, Luer/Luer, w/filter | C7010**        |
| Dual Tank Shuttle Valve System                    | GS2050         | Insufflation Tubing Set, disposable, IEC/Luer, w/filter  | C7011**        |
| CO2 Warmer, in-line, reusable                     | GS1012         | CPS Tubing Set, disposable, w/filter                     | GS1013         |
| Cleaning Brushes for Warmer assembly, 10/pkg      | GS1054         |                                                          |                |
| Adapter, CPC to IEC (for C7011 tubing)            | GS1050         |                                                          |                |
| Adapter, CPC to Luer (for C7010 tubing)           | GS1051         |                                                          |                |

\*Data on file

\*\*C7010 tubing requires a GS1051 adapter. C7011 tubing requires a GS1050 adapter.