## BM 900®

Spaltlampe  
Lampe à fente  
Slit lamp

Gebrauchsanweisung  
Mode d'emploi  
Instruction manual

Zusätzliche Informationen  
Informations supplémentaires  
Additional informations

## Vorwort

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## Introduction

We would like to thank you for your decision to purchase this Haag-Streit product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.

## Purpose of use

The **BM 900**<sup>®</sup> is an AC-powered slitlamp biomicroscope intended for use in eye examination of the anterior eye segment, from the cornea epithelium to the posterior capsule. It is used to aid in the diagnosis of diseases or trauma which affect the structural properties of the anterior eye segment.

## Contents

<b>1</b>	<b>Safety</b> .....	5
<b>2</b>	<b>Overview</b>	
2.1	Headrest .....	8
2.2	Slit lamp .....	8
<b>3</b>	<b>Operating</b>	
3.1	Adjusting the eyepieces .....	11
3.2	Prepare the patient .....	11
3.3	Operating the instrument .....	11
3.4	Filters .....	13
3.5	Apertures .....	13
<b>4</b>	<b>Equipment maintenance routines</b>	
4.1	Replacement of the bulb .....	15
4.2	Check electrical contacts of bulb and lamp housing .....	15
4.3	Use of contact medium .....	15
<b>5</b>	<b>Technical specifications</b> .....	17
	<b>Additional informations</b> .....	21

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# 1 Safety

## Ambient conditions

- Transportation Temperature -40 °C to +70 °C  
Air pressure 500 hPa to 1060 hPa  
Relative humidity 10% to 95%
- Storage Temperature -10 °C to +55 °C  
Air pressure 700 hPa to 1060 hPa  
Relative humidity 10% to 95%
- Working Temperature +10 °C to +35 °C  
Air pressure 800 hPa to 1060 hPa  
Relative humidity 30% to 75%

## Installation

- To avoid condensation, allow the instrument to adjust to room temperature for several hours before unpacking.
- Always mount the slit lamp and the headrest on an electric insulated and fire resistant table top.
- The rail covers (33) prevent from a tip over of the slit lamp.
- Check: do the connection pieces of the accessories fit tightly (like screwed connections and quick fit)?

## Operation and surrounding

- Only qualified and trained personnel should operate the equipment, the training is at the owner's responsibility.
- Use only Haag-Streit accessories.
- Shut down after every use. In case the dust cover is used: risk of overheating.
- Never use the ocular to look at the sun.
- Never operate the instrument in the same room with combustible gases, volatile solvents (alcohol, benzol) or flammable anesthetic agents.
- Avoid humidity.

## Changing the light bulbs (not for LI 900®)

- Switch off the main switch - disconnect the mains connector.
- Allow the bulb to cool down.
- Do not open the slit lamp for any other reason than to change the light bulb.

## Electrical


- Only a hospital grade 3-conductor electrical power supply cable must be used. *For USA and Canada: Detachable Power Supply Cord Set, UL Listed, type SJE, SJT or SJO, 3-conductor, not smaller than 18 AWG.*
- Plug, cable and ground lead connection of the socket have to be in perfect condition.

## Cleaning

- Clean the slit lamp housing only with a slightly water dampened cloth.
- No liquids, no alcohol or corrosive agents.
- Clean only the exterior surfaces of the optical parts.
- Use chin paper and clean the forehead band with a cotton pad soaked in alcohol.

## Light toxicity

- Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. This device should be used with filters that eliminate UV radiation (< 400 nm) and, whenever possible, filters that eliminate short-wavelength blue light (< 420 nm).
- The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

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• While no acute optical radiation hazards have been identified for slit lamps, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk.

• The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

- See caution notes on page 17

#### Warranty / product liability

• The instrument should be operated in accordance with the chapter 'Safety'. Incorrect operation can damage the instrument. Thus no warranty claims can be accepted.

• Continued use of an instrument that has been damaged by incorrect operation can lead to personal injury. The manufacturer cannot accept liability in this case.

• Repairs and alterations on this equipment should only be carried out by Haag-Streit service technicians or by authorized persons.

#### Statutory requirements

• The slit lamp BM 900® has been designed and constructed to conform with the IEC / EN 60 601-1-2 and EN ISO 15004-2 standards. Manufacturing procedures, testing, commissioning, maintenance and repair are conducted under the observance of Swiss and international regulations.

• When combining different medical and/or nonmedical electrical equipment standard IEC / EN 60 601-1 applies.

• The 'CE' marking confirms compliance of the slit lamp BM 900® with the directive 93/42/EEC.

• The slit lamp BM 900® meets the electromagnetic compatibility requirements of IEC / EN 60601-1-2. The appliance is designed to keep the generation and emission of electromagnetic interference to a level that does not affect the normal operation of other items of equipment. The appliance is likewise designed to provide a certain level of resistance to outside sources of electromagnetic interference. • All statutory accident prevention regulations are to be observed.

#### Classification

Standard IEC / EN 60 601-1 Slit Lamp BM 900® equipment safety class I.  
Application part Type B.  
Operation mode: continuous operation

CE-Regulation 93/42 EEC Class I  
FDA Class II

#### Pictograms

- a Read very carefully (in the instruction manual)!
- b Warning for hot surfaces!
- c Disconnect the electrical supply plug!
- d Manufacturer
- e Serial number
- f HS-Part Number
- g Application part Type B
- h Date of manufacture

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#### Environment

• Electric and electronic equipment must be separated from house waste! This equipment has been introduced into the market after the 13th of August 2005.

• Disposal via your local collecting point or your Haag-Streit distributor.


• Thus it is ensured, that no harmful materials get into the environment and that valuable raw materials can be used again.

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
## 2 Overview

### 2.1 Headrest


- 1 Lamp cable
- 2 Headrest
- 3 Forehead band
- 4 Headrest height marker (patient's eye)
- 5 Adjustable focusing target
- 6 Chin rest
- 7 Chin rest height adjustment

### 2.2 Slit Lamp BM 900®

- 8 Lamp cover tungstene and halogen (LED LI 900® see separate instruction manual)
- 9 Lever for grey - and redfree filter
- 10 Scale for angled position of the slit image
- 11 Interchangeable illumination mirror
- 12 Breath shield
- 13 Protective cover
- 14 Illumination unit / microscope angle scale
- 15 Illumination arm locking screw
- 16 Microscope arm locking screw
- 17 Slit width controls
- 18 Registered trade mark BM 900®
- 19 Slit / diaphragm scale
- 20 Slit length, slit rotation, blue filter and fixation star control
- 21 Cap screw for accessories base
- 22 Stereo microscope with eyepieces
- 23 Lens-changing lever
- 24 Mounting screw for stereo microscope
- 25 Centering screw
- 26 5° to 20° inclination angle latch
- 27 Screw to block horizontal movement
- 28 Rail covers
- 29 Control lever
- 30 Gliding plate

\* LED Beleuchtung LI 900® siehe separate Anleitung  
Illumination LED LI 900® voir mode d'emploi séparé  
LED illumination LI 900® see separate instruction manual

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
## 3 Operating

### 3.1 Adjusting the eyepieces

Before the first examination, the eyepieces should be adjusted singly to correspond with the refraction divergence of the examiner. Insert the provided focusing rod (31) in place of the protective cover (32) with the flat projection surface facing the microscope. Line up the illumination unit and the microscope in the central position (0°).

Adjust each eyepiece separately, by turning the knurled ring with the diopter scale separately, until the projected slit is seen in focus. Adjust from (+) to (-) at low magnification.

Younger examiners are recommended to make an additional correction on both eyepieces of - 1 to - 2 diopters each, to compensate for the accommodation that occurs as a result of convergence when conducting binocular examinations. The requisite pupil distance on the microscope can then be adjusted.


### 3.2 Prepare the patient

- The height of the table must be selected to enable the patient to obtain a firm rest for the forehead and chin.
- In order to ensure that only that part of the eye to be examined is illuminated and to avoid disturbing glare, the height of the light beam should be appropriately adjusted.
- Those parts of the equipment coming into contact with patients should be cleaned with a dry cloth before each examination.
- The light intensity should be reduced to a minimum on completion of each examination. This extends the life of the light bulb, slit lamp and protects the patient from severe glare at the beginning of an examination (not for LI 900®).

### 3.3 Operating the instrument

1. Move the chin rest (34) up or down with the adjustment screw (35) until the **eyes** of the patient are **level** with the black mark on the headrest column (33).
2. Focus the eyepieces (36) to suit **your refraction** by turning the knurled rings and then set to your **interpupillary distance**.
3. Switch on the **illumination** unit by turning the power supply switch.

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4. Rotate the control lever (45) until the **light beam** is at **eye level**.

5. Holding the control lever (45), inclined towards yourself, move the instrument base until the **slit** appears to be approximately in focus on the cornea. This coarse setting is achieved with the naked eye. Fine adjustments are obtained by rotating and tilting the lever while observing the image through the microscope (40).

6. The **slit width** is adjusted to the left or right by means of the knob (44), which also adjusts the angle between the microscope and the illumination unit.

7. The **slit image** is made vertical, horizontal, or given any desired **angle** by means of control knob (39) (notches at 45°, 90° and 135°; stops at 0° and 180°, markers in 10° steps).


8. In order to ensure the possibility of conducting a **binocular fundus examination** in optical cross section without difficulty from side angles between 3° and 10°, the shorter mirror (11) should be installed, the illumination unit is rotated through 90° with the knurled control knob (39) and with the latch (43) in stages of 5°, the illumination unit and the microscope are turned to the center position (0°).

9. The **magnification** of the microscope in the Slit Lamp BM 900® is altered by changing the lens with the lever (41) or by exchanging the eyepieces.

10. Loosening the centering screw (42) will allow the slit image to be moved off center for **scleral illumination**. Tightening the screw will recenter the slit image in the visual field of the microscope.

11. For **general observations** over a wide field the grey filter is placed in position by means of lever (38), and the slit is opened wide.

12. For examination of the **posterior of the eye**, preset lenses or diagnostic contact lenses have to be used.

### 3.4 Filters (38)


- 46 Open aperture
- 47 Grey filter (10%)
- 48 Redfree filter
- 49 Empty space for accommodating an add. filter

### 3.5 Apertures (37) operating (39)

- 50 Fixation star\*
- 51 Apertures producing areas of 8, 5, 3, 2, 1 and 0.2 mm ø
- 52 Indication of wedge-shaped diaphragm for continuous slit length adjustment in mm
- 53 Blue filter

\* The target with the fixation star is particularly suitable for examining the fixation in amblyopia.

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## 4 Equipment maintenance routines

The slit lamp BM 900® in combination with the illumination set LI 900® may be used without any maintenance during the entire life cycle.

### 4.1 Replacement of the bulb (not for LI 900®)

The lamp cover can be removed by depressing both pressure buttons (54) to expose the bulb.

#### CAUTION

Before changing the bulb, turn off the switch on the power supply and disconnect the power cord!

#### CAUTION

Increased temperatures can occur on the upper surface of the lamp housing! Allow the light bulb to cool down before removing it or use protective gloves to remove. Use only original HS (55) lamp bulbs!

When installing a new bulb, care should be taken to place the notch in the centering base into the positioning catch. Replace the lamp cover and push it down until the two pressure buttons engage securely.

### 4.2 Check electrical contacts of bulb and lamp housing

With suspicion of bad electrical contacts (flickering or loss of light):

- Clean contact surface (with sandpaper, screwdriver etc.) of bulb (56), the contact springs in the lamp housing (58) and (59)
- If necessary bend slightly the contact surfaces in the lamp housing (58) and (59) in direction of the bulb
- Bend slightly both of the retaining clips (57) outwards




### 4.3 Use of contact medium on the electrical contacts of the bulb and lamp housing

At each change of light bulb, lubricate the contact points (A) and (B) with enough contact medium. The contacts of the bulb and the lid of the lamp should not show any scorch marks.

- Fit the lamp cover and remove it to check if there is enough contact medium on the points of contact (C) and (D). Reapply some medium if needed.
- The product does not require labeling. (Perfluorocarbons, Polyetheroil PTFE).

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## 5 Technical specifications

### Slit illumination

Slit image width	0 - 8 mm continuous
Slit image length	1 - 8 mm continuous
Illumination fields	Circle $\varnothing$ : 8 / 5 / 3 / 2 / 1 / 0.2 mm Test mark with fixation star
Slit image radial range	$\pm 90^\circ$
Radial movement of the slit illumination relative to the microscope axis	Horizontal $\pm 90^\circ$ , vertical 0 - 20°
Filters	Blue, redfree (green), grey (10%). The UV filter and the heat absorption filter are permanently mounted.
Light source	Tungsten / halogen incandescent light bulb 6 V, 4.5 A / 7.5 V, 38 W
Illumination intensity	0 V 0 Lux at position 0 5 V ~150 000 Lux at position 1/2 6 V ~300 000 Lux at position 1 7.5 V ~600 000 Lux at position 2

LED illumination LI 900®, see separate instruction manual


### Slit illumination tungsten

Maximum examination times according to ISO 15004-2 and ISO 10939

**CAUTION** – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 3.2 minutes with slit illumination or 2.5 minutes with both slit and peripheral illumination.

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
### Slit illumination halogen

Maximum examination times according to ISO 15004-2 and ISO 10939

**CAUTION** – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 5.5 minutes with slit illumination or 3.5 minutes with both slit and peripheral illumination.

### Spektrum Halogen Beleuchtung / Spectre éclairage halogène / Spektrum halogen illumination

X= Wellenlänge / Longueur d'onde / Wavelength [nm]  
Y= Rel. Intensität / Intensité relative / Rel. intensity

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**Stereo-microscope**

Stereo angle	13°
Magnification changer	1x and 1.6x
Magnification	
with ocular 10x (standard)	10x and 16x
with ocular 25x	25x and 40x
Range of adjusting eyepieces	+ 8 to - 8 diopters
Inter pupillary distance	54 - 94 mm

**Instrument base**

Operation	Single handed 3-dimensional operation of the guide lever
Adjustment of the instrument base	80 mm (length) 30 mm (height) 100 mm (side)

**Weight**

11.4 kg (without  
power supply, headrest  
and options)

**Classification**

Against electric shock	Class I, type B
Against water	IPX0
Sterilisation and disinfection	
• for the slit lamp	not necessary
• for the headrest	use chin paper and clean the forehead band with a cotton pad soaked in alcohol
Operation mode	continuous

Certified Quality Management System

PRODUCTS CERTIFIED FOR BOTH THE U.S  
AND CANADIAN MARKETS, TO THE APPLI-  
CABLE U.S. AND CANADIAN STANDARDS

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## Additional informations

### A Indications perating the equipment

- A.1 Lamp bulb ..... 21
- A.2 Regulating the brightness ..... 21
- A.3 Microscope and eyepiece ..... 23
- A.4 Adjustable focusing target ..... 23

### B Mounting of the equipment

- B.1 Microscope and illumination unit ..... 25
- B.2 Headrest and power supply ..... 25

### C Equipment maintenance routines

- C.1 Replacement of the mirror ..... 27
- C.2 Cleaning of the optical parts ..... 27
- C.3 Adjustment of friction  
of the slit width control ..... 27
- C.4 Cleaning of the gliding plate  
and the rails ..... 27
- C.5 Cleaning of the axle ..... 27
- C.6 Plastic dust cover ..... 27


### D Accessories ..... 29

## A Indications instrument and operating

### A.1 Lamp bulb (not for LI 900®)

The tungsten / halogen lamp bulb is centered ex-works to obtain an optimal illumination. The 'H-S' marking on the centering base-cap indicates a genuine Haag-Streit spare part.

The light beam from the pre-centered bulb passes through the condenser, the slit diaphragm and the illumination lens, before it is directed to the patient's eye by the oblique mirror.

### A.2 Regulating the brightness (not for LI 900®)

The brightness of the slit light illumination can be varied by means of

- the knob on the power supply which adjusts the intensity of the slit illumination
- the 10 % grey filter on the upper part of the illumination unit of the slit lamp

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
### A.3 Microscope and eyepiece

The magnification of the microscope can be altered by changing the lens mounted in a revolving turret piece, or by exchanging the eyepieces. The lever (60) is to be adjusted from 1x to 1.6x when changing the lens. In this instance, the examination need not be interrupted; nor is a new adjustment of the microscope needed. This change will alter the overall magnification from 10x to 16x (on the 10x eyepiece), and from 16x to 25x (on the 16x eyepiece).

The separate focussing of each eyepiece serves not only to compensate the refraction divergence of the examiner, but also to achieve the best possible focussing, in view of the often slightly unsymmetrical passage of both microscope beam paths in fundus examinations.

- 60 Lens-changing lever
- 61 Prism housing
- 62 Eyepiece mounting
- 63 Diopter scale for eyepiece 16x for setting the refraction of the examiner ( $\pm 8$  D)
- 64 Diopter scale for eyepiece 10x for setting the refraction of the examiner ( $\pm 8$  D)
- 65 Index (white spot)
- 66 Eyepiece 10x or 25x

Lens	1x	1.6x	1x	1.6x
Eyepiece	10x	10x	25x	25x
Total magnification	10x	16x	25x	40x
Object image section $\varnothing$	18 mm	11.3 mm	8 mm	5 mm


### A.4 Adjustable focusing target

If the lamp cable (67) is connected, the fixation lamp burns as soon as the power supply switch is activated, then the focusing target (72) is visible.

The fixation target can be brought into any position by moving the bearer's arm (68). This permits uninterrupted observation through the microscope and simultaneous guidance of the patient's eye to any desired position. This is especially advantageous in examining the fundus.

It's a condition, that the patient's eye is correctly positioned at the headrest's height marker (71).

Ametropia of the patient from +10 D to -15 D (69) can be corrected (70). Steady fixation of the eye is maintained as accommodation and convergence are avoided.

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
## B Installation of the equipment

### B.1 Microscope and illumination unit


Place the stereo microscope on the saddle bearing of the microscope arm. Move it to the rear up against the stop piece and secure with the screw (74) appearing from obliquely below (use the large screwdriver).

Secure the breath shield (75) by means of the knurled screw (73) on the inner side of the bearer arm.


Before attaching the illumination unit to the column of the instrument base, the set screw (77) in the hub of the bearer arm is to be sufficiently loosened so that its tip no longer protrudes into the screw hole.

Grasp the illumination unit arm (76) as shown and place in position.


Tighten the set screw (77) in the hub of the bearer arm with the small screwdriver.

Small screwdriver  
HS-Part number 1001436

Large screwdriver  
HS-Part number 1001435

### B.2 Power supply

The corresponding Haag-Streit operating instructions should be observed!

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HS-Part number 1001434

## C Additional maintenance routines

### C.1 Replacement of the mirror


Removal of the mirror is easiest when the microscope and illumination unit are well separated and the latter inclined by 10° or so.

### C.2 Cleaning of the optical parts

A dust brush is to be employed for dusting the exposed glass surfaces.


**Not for LI 900®:** From time to time it is necessary to clean the upper surface of the condenser. The condenser surface is now visible and can be cleaned with a brush for cleaning the optics. Surfaces which have become really dirty should first be dusted, then wiped carefully with a soft dry clean cloth, washed linen, chamois leather, or some such material which will not scratch the surface.

Brush for cleaning the optics  
HS-Part number 1001398

### C.3 Adjustment of friction of the slit width control


The friction of the slit width control (79) can be regulated by means of the small screw in the centre of the right hand control. The movement can be stiffened by turning the screw clockwise or loosened by turning it counter-clockwise. The screw should be at least tight enough to ensure that the slit does not close spontaneously.

### C.4 Cleaning of the gliding plate and the rails

If it becomes difficult to move the slit lamp around on the gliding plate (83), the plate should be cleaned with a slightly oily rag.

After removing the rail covers (82), the rails (80) should be cleaned with a stiff brush.

### C.5 Cleaning of the axle (81)

Clean only with dry, lint-free cloths!

### C.6 Plastic dust cover

It is recommended to protect the slit lamp with a dust cover in case it is not in use.

Plastic dust cover, small (for slit lamp)

HS-Part number 1001395

Plastic dust cover, large (for several instruments)

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**Haag-Streit AG**