



**VETERINARY ENDOSCOPY**  
**- SMALL ANIMALS -**

**4th EDITION 2/2012**

Chapter			Pages
<b>1</b>	<b>OTOSCOPY AND RHINOSCOPY</b>	VET-S-RH	1-10
<b>2</b>	<b>BRONCHOSCOPY AND GASTROSCOPY, RIGID INSTRUMENT SETS</b>	VET-S-BRO	11-18
<b>3</b>	<b>FLEXIBLE ENDOSCOPES FOR RESPIRATORY AND URINARY TRACTS</b>	VET-S-BRO-FIB	19-28
<b>4</b>	<b>VIDEO GASTROSCOPES</b>	VET-S-VE-GE, VET-S-VE-ACC	29-46
<b>5</b>	<b>LAPAROSCOPES AND ACCESSORIES</b>	VET-S-LAP-TEL, VET-S-LAP-TROC-INTRO, VET-S-LAP-TROC, VET-S-LAP-INST, VET-S-LAP-NH, VET-S-LAP-SET	47-88
<b>6</b>	<b>HANDLES AND INSTRUMENTS FOR SUCTION AND IRRIGATION</b>	VET-S-LAP	89-92
<b>7</b>	<b>ARTHROSCOPES AND ACCESSORIES</b>	VET-S-ART, VET-S-ART-INST	93-118
<b>8</b>	<b>CYSTOSCOPY AND TRANSCERVICAL INSEMINATION</b>	VET-S-CYST	119-132
<b>9</b>	<b>ENDOSCOPY IN BIRDS AND SMALL EXOTIC ANIMALS</b>	VET-S-B	133-140
<b>10</b>	<b>VITOM® 25 SYSTEM – VISUALIZATION SYSTEM FOR OPEN SURGERY WITH MINIMAL ACCESS</b>	VET-S-VITOM	141-156
<b>11</b>	<b>UNITS AND ACCESSORIES</b>	VET-S-UNITS-INTRO, VET-S-UNITS	U 1-52
<b>12</b>	<b>EXCERPT FROM TELEPRESENCE, IMAGING SYSTEMS – DOCUMENTATION, ILLUMINATION – EQUIPMENT CARTS</b>	VET-S-DOK-HD, VET-S-DOK-3D, VET-S-DOK-TRI, VET-S-DOK-TELE, VET-S-TP-AIDA, VET-S-DOK-MEDIA, VET-S-DOK-L, VET-S-DOK-V, VET-S-TP-SET	TP 1-140
<b>13</b>	<b>EXCERPT FROM HYGIENE, MAINTENANCE, STERILIZATION, STORAGE TECHNIQUES</b>	VET-S-CL	CL 1-56
<b>14</b>	<b>COMPONENTS / SPARE PARTS</b>	VET-S-SP	SP 1-44



# **THE WORLD OF ENDOSCOPY**

Not all the products listed in this document are certified according to Regulation 2017/745/EU.

For this reason, some products requiring certification under this regulation may not be available in these countries.

## **VETERINARY ENDOSCOPY - SMALL ANIMALS -**

**4th EDITION 2/2012**

## Important Notes:

Endoscopes and accessories contained in this catalog have been designed in part with the cooperation of physicians and are manufactured by the KARL STORZ group. If subcontractors are hired to manufacture individual components, these are made according to proprietary KARL STORZ plans or drawings. Furthermore, these products are subject to strict quality and control guidelines of the KARL STORZ group. Both contractual and general legal provisions prohibit subcontractors from supplying components manufactured by order of KARL STORZ to competitors.

Any assumptions that competitors' endoscopes and accessories are acquired from the same suppliers as the KARL STORZ products are not correct. Moreover, endoscopes and instruments provided by competitors are not manufactured according to the design specifications of KARL STORZ. This means it cannot be assumed that these endoscopes and accessories – even if they look identical on the outside – are constructed in the same manner and have been tested according to the same criteria.

### Standardized Design and Labeling

KARL STORZ participates both in national and international bodies involved in the development of standards for endoscopes and endoscopic accessories. Standardized design and development therefore have long been implemented consistently by KARL STORZ. The user can rest assured that all products by the KARL STORZ group have been designed and constructed not only in compliance with strict internal quality guidelines, but also with international standards. All data relevant for safe use, such as viewing direction, sizes and diameters, or notes regarding sterilization of telescopes, are applied to the instruments, have been formulated according to international standards, and therefore provide reliable information.

As we constantly seek to improve and modify our products, we reserve the right to make changes in design that vary from catalog descriptions.

### Original or Counterfeit

KARL STORZ products are name brand articles renowned around the world and represent the state of the art in important areas of healthcare. A large number of “copy cat” products are currently being offered in many markets. These products are designed intentionally to resemble KARL STORZ products and use marketing strategies that at least point out their compatibility with KARL STORZ products. These products are by no means genuine products, since genuine KARL STORZ products are sold worldwide exclusively under the name of KARL STORZ, which appears on the packaging and the product. In the absence of such labeling, the product is not from KARL STORZ.

KARL STORZ, therefore, is unable to ensure that such products are actually compatible with genuine KARL STORZ products or can be used with them without injury to the patient.

# OTOSCOPY AND RHINOSCOPY





## Ear Diagnostic Set

Diagnostic otoscopy in veterinary medicine is one of the most important fields of application in endoscopy. Disorders of the external ear are common in dogs, therefore, the easy accessibility afforded by this technique makes otoendoscopy an ideal examination method.

A rigid endoscope can also be utilized to examine the middle ear.

KARL STORZ now offers a new veterinary otoscopy set which consists of an otoscope with a 0° forward-oblique telescope and a 3-way stopcock.

The shape of the telescope and its length facilitate thorough examination even in large breed dogs.

The instrument set can also be used for treatment, i.e. flushing the ear canal, suction, removal of foreign bodies, etc.

Used in conjunction with a KARL STORZ camera system and a video recorder or printer, it provides a precise means of assessing clinical results and documenting findings.



67260 OSA

67260 OSA **Veterinary Otoscope Set** including:

**Otoscope**, for dogs and cats, **autoclavable**, diameter 5 mm, length 8.5 cm, with working channel diameter 5 Fr.

**Stopcock Attachment**, with integrated working channel and LUER-Lock adaptor for irrigation

### Optional Accessories



60651 BK

**Cleaning Brush**, flexible, for use with working channels up to 2.3 mm, length 25 cm



**NEW** 69003 HH

**Otoscope Camera Head Holder**, for combining Otoscope 67260 AA with **non-autoclavable** TELECAM camera head and Veterinary Video Camera III camera head



**NEW** 15006 D

**“ULTRA STOP” Antifog Solution**, 15 ml, atomizer bottle

**Optional accessories for ear flushing and suction** see chapter 11, UNITS, VETPUMP® II

**Camera systems and accessories for documentation** see chapter 12, EXCERPT FROM TELEPRESENCE

**Components/Spare Parts** see chapter 14

# Otoscopy

## Instruments



67261 K

**Ear Curette**, small, 5 Fr., length 30 cm



67261 KL

**Ear Curette**, large, for insertion in the working channel from the distal end, length 30 cm



67161 T

**Grasping Forceps**, flexible, double action jaws, 5 Fr., length 34 cm



67161 Z

**Biopsy Forceps**, flexible, oval, double action jaws, 5 Fr., length 34 cm




67071 XS

**Miringotomy Needle**, with outer tube including:  
**Flexible Needle**, outer diameter 0.7 mm  
**Outer Tube**, 1.5 mm



67159 L


**Polypectomy Snare**, 5 Fr., length 32 cm including:  
**Handle**  
**Tube**  
**Snare**  
**Adaptor**



67479 S

67479 S

**Outer Tube**, with LUER-Lock adaptor, outer diameter 1.5 mm, working length 18 cm, for use with Flexible Needle 67071 S, or alone for suction and irrigation



62992



**NEW** 62992

**OTEX® Brush**, 5 Fr., bristle diameter 5 mm, length 19 cm



**NEW** 62993

**OTEX® Brush**, 5 Fr., bristle diameter 7 mm, length 19 cm

**Components/Spare Parts** see chapter 14






11301 D1/D3/D4/DE/DF


- 11301 D1 **Battery Light Source LED for Endoscopes**, with fine screw thread, brightness > 100 lm / > 150 klx, burning time > 120 min, weight approx. 150 g, waterproof and fully immersible for manual cleaning and disinfection, with 2 Photo Batteries 121306 P
- 11301 D3 **Same**, with coarse thread
- 11301 D4 **Same**, with fast screw thread, **suitable for wipe disinfection**
- 121306 P **Photo Battery**, lithium, 3 V, CR 123 A
- 11301 DE **Battery Light Source LED for Endoscopes**, rechargeable, with click connection, brightness > 110 lm / > 150 klx, color temperature 5500 K, lithium-ion batteries, charging time 60 min, burning time at 100% brightness 40 min, weight approx. 150 g ready for use, **suitable for wipe disinfection**
- 11301 DF **Battery Light Source LED for Endoscopes**, rechargeable, with fast screw thread, brightness > 110 lm / 150 klx, color temperature 5500 K, lithium-ion batteries, charging time 60 min, burning time at 100% brightness 40 min, weight approx. 150 g ready for use, **suitable for wipe disinfection**
- 094129 **Battery Charger Li-Ion**, for charging the rechargeable battery box, power supply 100 – 240 VAC, 50/60 Hz
- 094127 **Mains Cord**, for Battery Charger 094129, length 150 cm

# Rhinocopy


1.9 mm Telescopes, for rhinoscopy, pharyngoscopy and laryngoscopy, especially for cats and small dogs




67030 BA





64301 AA



61029 C



61029 D



67030 BA

**Operating Telescope 30°**, 9.5 Fr., working length 14 cm, **autoclavable**, fiber optic light transmission incorporated, with working channel 3 Fr. and 2 LUER-Lock cones  
color code: red

64301 AA

**HOPKINS® Straight Forward Telescope 0°**, diameter 1.9 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated,  
color code: green

64301 BA

**HOPKINS® Forward-Oblique Telescope 30°**, diameter 1.9 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated,  
color code: red

61029 C

**Protection and Examination Sheath**, working length 8.5 cm, with Obturator 61029 CO and 2 LUER-Lock cones,  
color code: blue

61029 D

**Operating Sheath**, 9 Fr., diameter 2.5 cm, working length 7 cm, with Obturator 61029 DO and working channel 3 Fr.,  
color code: red



61071 TJ

**Grasping Forceps**, for removing foreign bodies, double-action jaws, flexible, 3 Fr., length 28 cm




61071 ZJ



**Biopsy Forceps**, double-action jaws, flexible, 3 Fr., length 28 cm

# Rhinoscropy

2.7 mm Telescopes, for universal use



Telescope and sheath with working channel for biopsy and the removal of foreign bodies



64019 AA

**HOPKINS® II Straight Forward Telescope 0°**, diameter 2.7 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

64019 BA

**HOPKINS® II Forward-Oblique Telescope 30°**, diameter 2.7 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

67065 C

**Sheath**, 14.5 Fr., working length 15 cm, with working channel 5 Fr. and two stopcocks, including blunt Obturator 67065 CO

67161 T

**Grasping Forceps**, flexible, double action jaws, 5 Fr., length 34 cm

67161 Z

**Biopsy Forceps**, flexible, oval, double action jaws, 5 Fr., length 34 cm





64018 US

**Examination and Protection Sheath**, diameter 3.5 mm, working length 17 cm, for Telescope 64019 BA

10-002


Additional instruments for the working channel of Sheath 67065 C see chapter 9, page 138

# Nasal Forceps



## Nasopharyngeal Endoscopes

### Postrhinoscope with Uvula Retractor



723014

7230 EA


723014

**Uvula Retractor**, for use with HOPKINS® II Telescopes 7230 CWA/DA/EA




7230 CWA

**HOPKINS® II Wide Angle Lateral Telescope 70°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: yellow



7230 DA

**HOPKINS® II Lateral Telescope 90°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: blue



7230 EA

**HOPKINS® II Retrograde Telescope 120°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: white



650227

**Double Spoon Forceps**, horizontal opening, 110° curved backwards, spoon diameter 4 mm, with cleaning connector, working length 12 cm




**NEW** 723029

**Biopsy and Grasping Forceps**, vertical opening, malleable, cupped jaws diameter 4 mm, working length 18 cm



# BRONCHOSCOPY AND GASTROSCOPY RIGID INSTRUMENT SETS







Rigid endoscopes have certain advantages in bronchoscopy over flexible endoscopes. Superior image size, brightness and resolution guarantee high quality photographic and video documentation. A thorough examination of the bronchial tree can be performed

using endoscopes with various viewing angles. Samples for cytology, histopathology and micropathology are easily obtained. A variety of instruments for removal of foreign bodies is available.

#### **Indications for Rigid Bronchoscopy**

- **Cytology**
- **Retrieval of foreign bodies**
- **Removal of fibrinous casts and mucous plugs**
- **Bronchoalveolar lavage**
- **Diagnosis and staging of tumors**
- **Diagnosis and staging of chronic bronchitis**
- **Documentation**

# Standard Bronchoscopy Esophago-Gastroscope

Rigid Instrument Set, diameter 5.5/7 mm, working length 54 cm

61315 M 61315 N 61338 S

12015 AA

61924 D

61924 E

61015 US

61016 L

61016 N

61016 K

15006 B

61318 B

61315 M

61315 N

61338 S

61924 D

61924 E

61016 L - N

61016 K 61016 L 61016 N

61318 B

Basic set for rigid bronchoscopy in medium-sized and large dogs. Ideal for the retrieval of large and deeply embedded foreign bodies. This set can also be used for esophagoscopy.

12015 AA **HOPKINS® Straight Forward Telescope 0°**, diameter 5.5 mm, length 53 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

61015 US **Protection and Examination Sheath**, for HOPKINS® Telescope 12015 AA, with insufflation stopcock, diameter 7 mm, working length 53 cm (not illustrated)

61016 L **Optical Biopsy Forceps**, working length 53 cm

61016 N **Optical Biopsy Forceps**, universal jaws, for the removal of foreign bodies and denatured tissue, working length 53 cm

61016 K **Optical Forceps**, for soft foreign bodies

15006 B **“ULTRA STOP” Antifog Solution**, 25 ml, pipette bottle (not illustrated)

61318 B **Bronchoscope Tube Universal**, with Fiber Optic Light Carrier 61318 L, size 8.5, length 43 cm

61315 M **Glass Window Plug**

61315 N **Rubber Telescope Guide**

61338 S **Instrument Guide**, for suction catheter


61924 D **Adaptor**, for respirator

61924 E **Sealing Plug**, for respiration connector

Bronchoscopy may also be performed without the aid of a tracheal tube in the case of difficult anatomical conditions. While using the bronchoscope tube the examination can be performed under mechanical ventilation. The use of an anti-fogging solution to prevent fogging of the telescope lens is recommended.

# Standard Bronchoscopy Set

for use in small animals less than 10 kg,  
diameter 2.9/3.4 mm, working length 36 cm



61324 AA **HOPKINS® Straight Forward Telescope 0°**, diameter 2.9 mm, length 36 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

61378 H **Optical Alligator Forceps**, for hard foreign bodies, working length 35 cm, for use with Telescope 61324 AA

61378 L **Same**, for biopsy

61378 U **Optical Forceps**, universal, biopsy and grasping, working length 35 cm, for use with Telescope 61324 AA

61378 CF **Optical Forceps**, for broncho-esophagoscopy, 2 x 2 teeth, for controlled grasping of coins and flat foreign bodies, for use with Telescope 61324 AA

61324 US **Protection and Examination Sheath**, for HOPKINS® Telescope 61324 AA, with insufflation stopcock, diameter 4 mm, working length 36 cm (not illustrated)

61339 CD **Bronchoscope Tube**, size 3.5, length 30 cm

61338 M **Glass Window Plug**

61338 N **Rubber Telescope Guide**, for telescopes with diameter 2.7 and 2.9 mm and optical forceps

10101 FA **Prismatic Light Deflector, autoclavable**, with connection to fiber optic light cable

61338 S **Instrument Guide**, for suction catheter

61924 D **Adaptor**, for respirator

61924 E **Sealing Plug**, for respiration connector


The combined use of a HOPKINS® telescope with various optical forceps provide optimal conditions for performing biopsies or for removing foreign bodies, even in the case of very small lumina. The protection and examination sheath may provide an additional oxygen supply.




61378 CF 61378 H 61378 L 61378 U

# Accessories for Bronchoscopy

Instruments for Aspiration and Irrigation




- 10432 N **Aspirator**, for collecting secretions, with cut-off hole, with Spare Glass 10432 E, can be connected to suction tubes or suction catheters
- 10432 ED **Collection Device**, plastic, for single use, with lid, package of 100, for use with HUZLY Aspirator and Bronchus Irrigator 10434 and Intermediate Ring 10434 Z
- 10432 EE **Shipping Container**, for Collection Device 10432 ED, plastic, for single use, with lid, package of 100
- 10434 **HUZLY Aspirator and Bronchus Irrigator**, with Graduated Container 10455 A and Rubber Bulb 437000. After switching manifold from "Irrigation" to "Suction", the injected irrigation fluid is immediately withdrawn again together with the bronchial secretion



- 10434 A **Suction Tube**, diameter 4 mm, length 50 cm, for use with HUZLY Aspirator and Bronchus Irrigator 10434
- 10434 B **Same**, diameter 3 mm, length 35 cm
- 10434 V **Adaptor**, for catheters
- 10434 Z **Intermediate Ring**, for use with Collection Device 10432 ED with HUZLY Aspirator and Bronchus Irrigator 10434

For use with Telescope 12015 AA



61016 E **Esophagoscope Tube**, size 11, outer diameter 12 mm, inner diameter 8.6 x 10.3 mm, with suction channel and air inflation device, including Glass Window Plug with LUER-Lock 12016 EW, with Fiber Optic Light Carrier 12016 EM

61016 C **A. F. (anti-fog) Tube**, with working channel and anti-fog device, for use with Straight Forward Telescope 12015 AA

12016 U **Forceps**, spoon-shaped, flexible, for biopsy, width of jaws 3.2 mm, working length 61 cm

40924 **Rubber Insufflation Bulb**, with detachable tube, with LUER-Lock including:  
**Rubber Bulb**  
**Rubber Tube**

40924

61016 C


61016 E

12016 U

**Components/Spare Parts** see chapter 14

# Accessories for Gastroscopy and Bronchoscopy

Foreign Body Retrieval Forceps



10370 UX



10370 UX

**Forceps**, universal, for biopsy and foreign bodies removal, double action jaws, sheath diameter 2.5 mm, width 4 mm, working length 60 cm




**NEW** 60370 FB

**Foreign Body Forceps**, double action jaws, jaws 30 mm, jaw opening 3 cm, for tracheal, oesophageal and gastric foreign body removal, positioned parallel to an esophago-gastroscope, length 70 cm



## Instruments for airway sample collection



10436


10436

**Aspiration Biopsy Needle**, for transbronchial biopsy, distal straight, diameter 1.6 mm, working length 47 cm



10438

**Same**, distally angled



12016 S

**Suction Catheter**, 10 Fr., length 60 cm, with Adaptor 10479 D

# FLEXIBLE ENDOSCOPES FOR RESPIRATORY AND URINARY TRACTS



Fine diameter flexible endoscopes are ideal for numerous respiratory and urinary tract applications in dogs, cats and small exotic animals. KARL STORZ offers several options in this category that always include unsurpassed optics in a robust design that is easy to handle and includes a durable working channel for flushing or passage of instrumentation and catheters. The integral working channels facilitates the collection of biopsy,

fluid and cytology samples. The retrieval of foreign bodies under visual control is also possible.

KARL STORZ offers a complete range of accessories, i.e. catheters for cytology as well as various biopsy and grasping forceps. Since these flexible endoscopes are fully immersible and available in a variety of sizes, the applications are virtually unlimited.

#### **Indications for Flexible Bronchoscopy**

- **Cytology**
- **Retrieval of foreign bodies**
- **Bronchoalveolar lavage**
- **Diagnosis and staging of tumors**
- **Diagnosis and staging of chronic bronchitis**
- **Endoscopic placement of tracheal stents**
- **Documentation**

#### **Indications for Urethrocystoscopy**

- **Procurement of biopsies**
- **Identification of bleeding source**
- **Diagnosis of ectopic ureter**
- **Diagnosis and staging of tumors**
- **Diagnosis and staging of chronic cystitis**
- **Removal of uroliths**
- **Endoscopic placement of urethral and ureteral stents**
- **Documentation**




# The KARL STORZ Video-Uretero-Cystoscope

When our vision meets your most challenging requests: KARL STORZ as a pioneer of Endoscopy opens a new

door in terms of diagnosis and treatment in narrow anatomies



## Special features:

- CMOS technology for outstanding quality of view
- LED technology ensuring homogenous illumination
- Ergonomic handle due to the integrated light technology
- Extreme thinness and maneuverability
- Maximal deflection of 270° up- and downwards
- 1 to 1 response of the tip by manipulating the handle due to the high torque stability
- Protection against thermal effects by LASERITE ceramic liner



# The KARL STORZ Video-Uretero-Cystoscope

8.5 Fr.

Order No.	Deflection	Direction of view	Field of view	Working length	Working channel inner diameter	Sheath size
<b>NEW</b> 11278 V		0°	90°	70 cm	3.6 Fr.	8.5 Fr.
<b>NEW</b> 11278 VU*		0°	90°	70 cm	3.6 Fr.	8.5 Fr.

\* with contrapositive deflection mechanism

The flexible Video Uretero-Renoscope 11278 V/11278 VU is fully compatible with the IMAGE 1 HUB™ HD camera control unit (CCU)

Set 11278 VK/VUK

**Accessories included:**

- 27677 X **Case**
- 13242 XL **Leakage Tester**, with bulb and manometer
- 27651 AL **Cleaning Brush**, flexible, outer diameter 2 mm, for working channel diameter 1.2 – 1.8 mm, length 150 cm
- 27014 Y **LUER-Adaptor**, with seal

**Additional accessories:**




60275 FE

**Grasping Forceps**, double-action jaws, 3 Fr., length 120 cm



60275 ZE

**Biopsy Forceps**, oval cupped jaws, 3 Fr., length 120 cm




27023 VB

**Stone Basket**, sterile, for single use, 2.5 Fr., length 120 cm



**NEW** 11770 T

**Coagulating Electrode**, unipolar, 3 Fr., length 110 cm



**NEW** 27025 P

**Guide Wire**, with ball end, 3 Fr., package of 10



39405 AS

**Plastic Container for Flexible Endoscopes**, specially suited for gas and hydrogen peroxide (Sterrad®) sterilization and storage, for use with one flexible endoscope, external dimensions (w x d x h): 550 x 260 x 90 mm

**Note: Compatible with IMAGE 1 HUB™ HD Camera Control Units only**  
see chapter 12, EXCERPT FROM TELEPRESENCE

# The KARL STORZ Video Bronchoscope

The Compact Solution

KARL STORZ is expanding its line of flexible bronchoscopes with the new video bronchoscope with distal chip. Users who do not want to miss out on the high image quality of videoscopes any longer, can now take advantage of the quality products of KARL STORZ in this field too.

Whereas in the past comparatively large and heavy endoscopy towers were used in video endoscopy, the TELE PACK VET X control unit now provides a compact and, at the same time, multifunctional documentation unit. With its Hi-Lux high-performance light source and the high-resolution 15" flat screen monitor, the system provides an image with excellent color fidelity. Via the integrated recording button on the handle of the video bronchoscope, stills and videos can be stored on USB drive or SD card.

In combination with the TELE PACK VET X control unit, the new video bronchoscope offers outstanding image sharpness, as well as excellent brightness and representation of detail.

Other features of the video bronchoscope are:

- Comfortable handling due to ergonomic design and light weight
- Removable suction unit
- Overview due to full screen display
- Working length of 61 cm
- Full automatic brightness control at maximum and minimum range

## Very good hygiene due to removable suction unit



## Sensitive and stable control via ergonomic handle





## Short insertion tip with smooth-running and non-twisting angle mechanism




## Ergonomic video connector with sealing cap ensures comfortable handling for the personnel

## Optimally matched, non-rotational insertion sheath



# The KARL STORZ Video Bronchoscope

The Compact Solution



- 11900 BP **Video Bronchoscope, color system PAL**, including mounted Suction Unit 11900 AE, for use with TELECAM SL II SCB 20 2130 11, Veterinary Video Camera III 69 2360 11 and TELE PACK VET X 69 0450 01
- Working channel: 2.3 mm  
Direction of view: 0°  
Angle of view: 120°  
Deflection: up 180°, down 100°  
Depth of view: 3 – 50 mm  
Working length: 61 cm  
Outer diameter: 6.2 mm

11900 BN **Same, NTSC**

**Following accessories are included (for color systems PAL and NTSC):**

- 27677 KV **Case**  
11001 KL **Biopsy Forceps**, flexible, spoon-shaped, round, double action jaws, diameter 1.8 mm, working length 120 cm
- 20 2130 70 **Video Connecting Cable**, for use with TELE CAM, Veterinary Video Camera and TELE PACK VET X camera control units
- 11025 XE **Pressure Compensation Cap**, for ventilation during gas sterilization  
13242 LX **Leakage Tester**, with bulb and manometer  
27651 B **Cleaning Brush**, flexible, length 100 cm  
13272 **Mouth Piece**  
13991 SS **Irrigation Tube**  
27656 B **Instrument Oil**, 50 ml, silicone-free  
13990 SF **Protective Film**

**Please note:** To connect KARL STORZ video bronchoscopes with IMAGE1 camera control unit, IMAGE1 Video Connecting Cable 22 2000 77 is required.


Further suitable instruments with a diameter  $\leq 1.8$  can be found on page 28.

**Camera control unit** see chapter 12, EXCERPT FROM TELEPRESENCE

## Fiberscopes

### Special features:




- Veterinary-specific dimensions include extended length and reduced diameter facilitating use in a wide variety of patient sizes
- Large viewing angle and moveable distal tip facilitates orientation
- Working channel for the use of operating instruments, suction, insufflation or irrigation
- Waterproof, fully immersible for cleaning and disinfection
- Simple leakage test prevents costly repairs
- Sterilizable with EtO and FO gas, STERIS® and STERRAD®
- Exceptional optical quality of both lens and fiber optic image transmission bundle
- Resistant construction and robust mechanical components ensure long life of the instrument






# Fiberscopes

Broncho-Fiberscopes	Order No.	Distal Tip	Deflection	Direction of view	Field of view
---------------------	-----------	------------	------------	-------------------	---------------




**For routine bronchoscopy and postrhinocopy of medium-sized and large dogs**

 <p>5.2 x 85</p>	60001 VL1			0°	110°
---	-----------	---	--	----	------


**For routine bronchoscopy and postrhinocopy of cats and small animal less than 10 kg**

 <p>3.7 x 54</p>	60002 VB1			0°	90°
---	-----------	---	--	----	-----

**For deep lung examination and urethrocystoscopy in male dogs**

 <p>3.0 x 100</p>	60003 VB1			0°	90°
--	-----------	---	--	----	-----


**Accessories included:**



60025 E



13242 XL



6927691



11301 CA



11301 CB



11301 CD

60025 E

**Pressure Compensation Cap,**  
for ventilation during gas sterilization

13242 XL

**Leakage Tester,** with bulb and manometer

6927691

**Adaptor for Two-Way Stopcock,**  
LUER-Lock, with O<sub>2</sub> tube connection  
(60001 VL1 and 60002 VB1 only)

11301 CA

**Leaflet Valve,** for single use, package of 20  
(60001 VL1 and 60002 VB1 only)

11301 CB

**Suction Valve,** reusable (60001 VL1  
and 60002 VB1 only)

11301 CD

**Irrigation Adaptor,** for machine cleaning, for  
fiberscopes (60001 VL1 and 60002 VB1 only)

Working Length	Total Length	Working Channel Inner Diameter	Distal End Outer Diameter	Accessories included:			
				Case	Pressure Compensation Cap	Leakage Tester	Cleaning Brush
85 cm	113 cm	2.3 mm	5.2 mm	27677 A	60025 E	13242 XL	60651 BL
54 cm	82 cm	1.5 mm	3.7 mm	27677 A	60025 E	13242 XL	60651 AL
100 cm	128 cm	1.2 mm	3 mm	27677 A	60025 E	13242 XL	60651 AL

**Instrument-specific accessories, see table:**



60651 AL/BL


60651 AL

**Cleaning Brush**, flexible, for working channels of Fiberscope 60003 VB1 and 60002 VB1, length 150 cm

60651 BL

**Cleaning Brush**, flexible, for working channel of for Broncho-Fiberscope 60001 VL1, length 120 cm

**Optional accessory:**



11020 D

11020 D

**Ventilation Attachment**, attaches to endotracheal tube, allowing continuous administration of oxygen and anesthetic gas during tracheobronchoscopy



**Additional accessories and instruments for fiberscopes** see page 28  
**Fiber optic light cable for use with fiberscopes and Cold Light Fountains**  
see chapter 12, EXCERPT FROM TELEPRESENCE

**All fiberscopes and accessories are delivered in a sturdy case for storage and transportation.**

# Accessories for Fiberscopes and Fine Diameter Video Endoscopes

Biopsy Forceps, Grasping Forceps, Stone Baskets

## Accessories for fiberscopes with a working channel of at least 1.2 mm




**NEW** 60003 MA

**Biopsy Forceps**, oval, flexible, double action jaws, diameter 1 mm, length 110 cm



**NEW** 60003 MB

**Grasping Forceps**, flexible, double action jaws, diameter 1 mm, length 110 cm



**NEW** 27023 Y

**Brush for Cytology**, unsterile, for single use, 3 Fr., package of 5, length 100 cm, (not for use with endoscopes longer than 85 cm)



27023 VU

**Stone Basket**, sterile, for single use, 3 Fr., length 115 cm



**NEW** 11770 T

**Coagulating Electrode**, unipolar, 3 Fr., length 110 cm

## Accessories for fiberscopes and video endoscopes with working channel of at least 2 mm



60001 KL

**Biopsy Forceps**, oval cupped jaws, diameter 1.7 mm, length 120 cm



60002 KS

**Grasping Forceps**, double-action jaws, diameter 1.7 mm, length 120 cm



**NEW** 60002 KX

**Grasping Forceps**, 2 x 3 teeth, flexible, diameter 1.7 mm, working length 120 cm



11008

**Suction Catheter**, 5 Fr., length 100 cm, with Adaptor 10479 F



60007 BL

**Cytology Brush**, diameter 1.8 mm, length 100 cm

including:

**Outer Tube**

**Handle**

Components/Spare Parts see chapter 14



# VIDEO GASTROSCOPES





# Video Endoscopes

KARL STORZ's new video endoscope range sets new standards for image quality. In digital processing technology, high resolution CCD chips and an enhanced optical system results in endoscopes that produce images with unique clarity. The images are distinguished

by extreme brightness, resolution, size and depth of field. By incorporating two diametrically opposed light bundles and intelligent digital processor control, evenly illuminated images are provided, even under difficult lighting conditions.

### Special features:

- One camera control unit for use with video endoscopes, fiberscopes and rigid telescopes
- High resolution, true colour video images
- Large depth of field
- Evenly illuminated, bright images – even under difficult lighting conditions
- 120 – 140° wide angle telescopes with high light transmitting capacity
- New, sensitive control unit
- Optimal location of controls
- Anti-rotational insertion sheath
- Non-twisting supply tube with standard connection



Camera Control Unit 69 2360 11 see chapter 12, EXCERPT FROM TELEPRESENCE



# Video Endoscopes

Flexible Video Gastrosopes and Bronchoscopes	Order No.		Sheath Diameter	Working Channel Diameter	Working Length
	PAL	NTSC			
* 9.7 mm x 140 cm	60914 PKS	60914 NKS	9.7 mm	2.8 mm	140 cm
<b>NEW</b> * 7.9 mm x 140 cm	60714 PKS	60714 NKS	7.9 mm	2.8 mm	140 cm
<b>NEW</b> * 5.9 mm x 110 cm	60511 PKS	60511 NKS	5.9 mm	2 mm	110 cm

\* These video gastrosopes are also available with a working length of 110 cm.

## Accessories included:

	13991 SS		13991 YS		11025 XE		13242 LX		13991 DA
	13991 SRV		13991 RV		13242 LX		13991 DA		13991 RV
	13991 YS		13991 SRV		13991 LV		13991 SAV		13991 SV
	60991 SSV		13991 SV/SSV		13991 SSV		60991 SSV		27677 KV
					13990 SF				

Accessories for flexible endoscopes see page 35-45

Camera Control Units see chapter 12, EXCERPT FROM TELEPRESENCE

To connect KARL STORZ video endoscopes with IMAGE1 camera control unit, IMAGE1 Video Connecting Cable 22 2000 77 is required.

Deflection up/down	Deflection left/right	Field of Vision	Depth of Field	Accessories including
				Cleaning Brush included
		140°	2 – 200 mm	27651 H
		120°	2 – 200 mm	27651 H
		120°	2 – 200 mm	27651 F

**Accessories included:**



20 2130 70



27656 B



13992 BS



20 1000 31

- 27656 B
- 20 2130 70
- 13992 BS
- 20 1000 31

- Instrument Oil**, 50 ml, silicone-free
- Video Connecting Cable**, for use with TELECAM SL II, Veterinary Video Camera III and TELE PACK VET X Camera Control Units
- Water Bottle**, 250 ml
- Bottle Holder**

**Instrument-specific accessories, see table:**



27651 H

- 27651 H
- 27651 F

- Cleaning brush**, bristle diameter 5 mm, length 230 cm
- Cleaning brush**, bristle diameter 2.5 mm, length 180 cm

**Accessories not included:**



610 VE

- 610 VE

- Silicone Tubing Set**, for insufflation when using an external airpump, in combination with Adaptor 487 VE or 69487 KVE



60991 LA


- 60991 LA

- Videoendoscope LUER-Adapter**, irrigation adapter to connect to the working channel of KARL STORZ videoendoscopes

# Video Endoscopes

**Specifications:**


- Field of View: 120°/140°
- Depth of View: 2 – 200 mm
- Resolution: 480,000 pixels
- 2 push-button controls – programmable



**Accessories** see pages 35-45

# Accessories for Flexible Endoscopes

Biopsy Forceps, reusable



	round cupped jaws	round cupped jaws, with spike, fenestrated	oval cupped jaws, fenestrated	oval cupped jaws, with spike, fenestrated	oval cupped jaws, serrated
Working length					

### Diameter 1.7 mm – for use with working channel of at least size 2 mm

200 cm	60001 LL	60001 LM	60001 LN	60001 LP	–
<b>NEW</b> 160 cm	11001 ML	11001 MM	11001 MN	11001 MP	–

### Diameter 1.8 mm – for use with working channel of at least size 2.2 mm

180 cm	–	–	–	–	60180 LT
--------	---	---	---	---	----------


### Diameter 2.2 mm – for use with working channel of at least size 2.8 mm

200 cm	60250 LL	60250 LM	60250 LN	60250 LP	60250 LT
--------	----------	----------	----------	----------	----------

**NEW** 110285-10\* **Biopsy Forceps**, for single use, sterile, with oval jaws, coated, working length 160 cm, package of 10, for use with working channels as of diameter 2 mm





5-06

\*  mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# Accessories for Flexible Endoscopes

Biopsy Forceps, insulated



Working Length	round cupped jaws
	

## Diameter 1.7 mm

200 cm	60751 LL	-
<b>NEW</b> 160 cm	13751 M	-


## Diameter 2.2 mm






200 cm	-	60751 LN
--------	---	----------



# Accessories for Flexible Endoscopes

Grasping Forceps, reusable



	alligator jaws	alligator jaws, round	universal, round cupped jaws, serrated	alligator jaws, serrated, with 2 x 3/3 x 4 teeth	serrated, with 1 x 2 teeth
Working Length					

**Diameter 1.8 mm – for use with working channel at least size 2 mm**


200 cm	60002 LS	60002 LT	60002 LU	60002 LX	60002 LZ
<b>NEW</b> 160 cm	–	11002 MT	11002 MU	11002 MX	11002 MZ

**Diameter 2.3 mm – for use with working channel at least size 2.8 mm**

200 cm	60252 LS	60252 LT	60252 LU	60252 LX	60252 LZ
--------	----------	----------	----------	----------	----------

# Accessories for Flexible Endoscopes

Grasping Forceps, Foreign Body Baskets



Diameter 1.8 mm – for use with working channel of at least size 2.2 mm

Working length	Grasping Forceps, flexible, with 2 spring-action jaws, 1 x 2 teeth, sharp
	
180 cm	60185 LA

Diameter 2.4 and 1.8 mm

Working length	Foreign Body Basket
	

Diameter 2.4 mm (with 6 wires), large opening: 14 x 10 cm


200 cm	60278 FK
--------	----------

Diameter 1.8 mm (with 4 wires)

200 cm	60185 FK
--------	----------

# Accessories for Flexible Endoscopes

Cytology Brushes



Diameter 2.2 mm – for use with working channel of at least size 2.8 mm

Working length	Cytology Brush, complete	Outer Tubes, only, package of 10	Brushes, only, package of 10
200 cm	60266 L	60266 SL	60266 LE



10479 C

10479 C **Adaptor, for outer tubes**



13266 H

13266 H **Handle, for brushes**

Diameter 2.2 mm – for use with working channel of at least size 2.8 mm, for single use, package of 10




Working length	Cytology Brushes, complete
180 cm	110240-10*
<b>NEW</b> 230 cm	110241-10*

5-061


\* mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# Accessories for Flexible Endoscopes

Suction catheters for broncho-alveolar lavage (BAL), reusable, autoclavable,  
Injection Needle



**Diameter 2.3 and 1.8 mm – for use with working channel of at least size 2.8 mm**

Catheter for Aspiration Biopsy	
Working length	


**Diameter 2.3 mm – for use with working channel of at least size 2.8 mm**

500 cm	60469 C
--------	---------

**Diameter 1.8 mm – for use with working channel of at least size 2 mm**

180 cm	60725 TM
--------	----------


**Diameter 2.2 mm – for use with working channel of at least size 2.8 mm**

Injection Needle	
Working length	

200 cm	60260 L
--------	---------

# Accessories for Flexible Endoscopes

Coagulating Electrodes for use with High Frequency Surgery Units



## KARL STORZ HF Connector, 4 mm, instrument side


Diameter 1.8 mm,  
for use with working channel  
of at least size 2 mm



Diameter 2.2 mm,  
for use with working channel  
of at least size 2.8 mm



	unipolar
Working length	
200 cm	60765 L

	bipolar
Working length	
200 cm	60267 LB


High Frequency Surgery Units, see chapter 11, UNITS

# Accessories for Flexible Endoscopes

Polypectomy Snares, insulated, reusable,  
for use with High Frequency Surgery Units

**Special features:**

- Easy-to-use mechanism
- Forms snare automatically
- Multifilament cutting wire
- Asymmetric opening



Diameter 2.4 mm – for use with working channels of at least size 2.8 mm  
KARL STORZ HF connector, 4 mm, instrument side

Working length	Complete Instrument	Spare Snares, reusable	Spare Tubes
<b>Polypectomy Snare, complete, with Handle 60773 H, snare size 25 mm</b>			
200 cm	60773 LD	60273 TD	60273 TW
<b>Polypectomy Snare, complete, with Handle 60774 H, snare size 40 mm</b>			
200 cm	60773 LE	60273 TE	13273 TX
<b>Polypectomy Snare, complete, with Handle 60774 H, snare size 60 mm</b>			
200 cm	60773 LF	60273 TF	13273 TX



60774 H

60774 H **Handle**, only, reusable



60773 H

60773 H **Handle**, only, reusable


High Frequency Surgery Units, see chapter 11, UNITS

# Accessories for Flexible Endoscopes

Polypectomy Snares, insulated, reusable,  
for use with High Frequency Surgery Units

**Special features:**

- Easy-to-use mechanism
- Multifilament cutting wire



**Oval Diathermy Loops**

Diameter 1.8 mm – for use with working channels of at least size 2 mm  
KARL STORZ HF Connector, 4 mm, instrument side

<b>Working length</b>	 25 mm oval	<b>Spare Tubes</b>
	<b>Polypectomy Snare, complete, with Handle 60773 H</b>	
180 cm	60725 LA	
	<b>Spare Snares, reusable</b>	
180 cm	60725 TA	60725 TM



60773 H

60773 H

**Handle, only, reusable**

5-06


**High Frequency Surgery Units**, see chapter 11, UNITS

# Accessories for Flexible Endoscopes

Polypectomy Snares, insulated, reusable,  
for use with High Frequency Surgery Units

**Special features:**

- Easy-to-use mechanism
- Multifilament cutting wire



**Oval Diathermy Loops**

Diameter 2.4 mm – for use with working channels of at least size 2.8 mm  
KARL STORZ HF Connector, 4 mm, instrument side

<b>Working length</b>	 35 mm oval	 30 mm oval	 30 mm hexagonal	<b>Spare Tubes</b>
-----------------------	-------------------	-------------------	------------------------	--------------------

**Polypectomy Snares, complete, with Handle 60773 H**

200 cm	60773 LA	60773 LB	60773 LC
--------	----------	----------	----------

**Spare Snares**

200 cm	60273 TA	60273 TB	60273 TC	60273 TT
--------	----------	----------	----------	----------



60773 H

60773 H

**Handle, only, reusable**

High Frequency Surgery Units, see chapter 11, UNITS



# Accessories

High Frequency Cords,  
for use with unipolar and bipolar instruments

## Unipolar High Frequency Cords



KARL STORZ Instrument High Frequency Surgery Units

	26002 M	<b>Unipolar High Frequency Cord</b> , with 4 mm plug for models KARL STORZ, Erbe type T, older models and Ellman, length 300 cm
	26004 M	<b>Unipolar High Frequency Cord</b> , with 4 mm plug for models Berchtold and Martin, length 300 cm
	26005 M	<b>Unipolar High Frequency Cord</b> , with 5 mm plug for AUTOCON® II 400 SCB system (111, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ AUTOCON® system (50, 200, 350) and Erbe type ICC, length 300 cm
	26006 M	<b>Unipolar High Frequency Cord</b> , with 8 mm plug for models KARL STORZ AUTOCON® II 400 SCB system (112, 116) and Valleylab, length 300 cm

## Bipolar High Frequency Cords



KARL STORZ Instrument High Frequency Surgery Units


	26176 LE	<b>Bipolar High Frequency Cord</b> , for AUTOCON® II 400 SCB system (111, 113, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ Coagulator 26021 B/C/D, 860021 B/C/D, 27810 B/C/D, 28810 B/C/D, AUTOCON® system (50, 200, 350), Erbe-Coagulator, T and ICC series, length 300 cm
	26176 L	<b>Bipolar High Frequency Cord</b> , with 2x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Ellman, length 300 cm
	26176 LA	<b>Bipolar High Frequency Cord</b> , with 2x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Valleylab, length 300 cm
	26176 LM	<b>Bipolar High Frequency Cord</b> , for Martin and Berchtold coagulators, length 300 cm
	26176 LV	<b>Bipolar High Frequency Cord</b> , for KARL STORZ AUTOCON® II 400 SCB system (112, 114, 116, 122, 125) and Valleylab coagulators, length 300 cm

**Please note:**

All high frequency cords on this page are delivered with a length of 300 cm. If a length of 500 cm is required, please add the letter **L**, i.e. 26002 **ML**, 26176 **LVL** to the order number.



# LAPAROSCOPES AND ACCESSORIES





New laparoscopic techniques provide new alternatives in veterinary medicine. KARL STORZ offers telescopes and instrument sets for these techniques which allow diagnostic interventions and biopsies as well as minimally invasive operations such as laparoscopic ovariectomies of the female dog.

User-friendly CO<sub>2</sub> insufflators ensure optimal working conditions. The operative procedure can be displayed on a monitor by means of a video camera attached to the endoscope. Compared to conventional open surgery, laparoscopic techniques allow a less traumatic surgical intervention in most cases and a more extensive view of the entire abdomen.

**Indications for Laparoscopy and Thoracoscopy:**

- **Diagnostic laparoscopy**
- **Cryptorchidectomy**
- **Prophylactic gastropexy**
- **Ovariectomy/ovariohysterectomy**
- **Jejunostomy/gastrostomy tube placement**
- **Adrenalectomy**
- **Pericardectomy**
- **Partial lobectomy of the lung**
- **Nephrectomy**
- **Insulinoma removal**
- **Partial pancreatectomy**
- **Correction of gastric dilatation with volvulus**
- **Laparoscopic assisted cystotomy**
- **Laparoscopic assisted cystopexy**



## Telescope with variable direction of view

Until now, surgeons had to choose in advance which telescope or direction of view to use in a procedure. Moreover, surgeons were restricted to the selected direction of view throughout the surgery or had to make an intraoperative telescope change.


To prevent this predicament in the future, we developed the ENDOCAMELEON®: a telescope that allows you to adjust the desired direction of view – also during surgery – between 0° and 120°.

The ENDOCAMELEON® combines the user comfort of the proven 0° HOPKINS® telescope with the advantages and potential of a telescope featuring a variable direction of view – offering you the quality you expect from KARL STORZ telescopes.

The innovative ENDOCAMELEON® technology is not difficult to use and, due to the external moving parts, does not take up extra intracorporal space. Handling remains straightforward and ergonomic. Image alignment is the same as rigid telescopes: the direction of view is selected by simply turning the adjustment knob, making the system very intuitive to use. As the ENDOCAMELEON® is equipped with a standard eyepiece, the variable direction of view benefits all standard camera systems. Thanks to the HOPKINS® rod lens system, ENDOCAMELEON® also offers image quality that enables a useful application of three-chip cameras or HD camera systems.

To have the direction of view best suited for each situation available at all times offers the surgeon a higher degree of safety. With the ENDOCAMELEON®, visual inspection of the entire surgical field is easily achieved. Instrument movement can be controlled throughout the entire procedure and hemorrhages in

previously inaccessible areas can be detected and controlled. With a simple turn of the adjusting knob, the ENDOCAMELEON® enables the user to easily select the direction of view between 0° and 120° to suit all OR requirements.



ENDOCAMELEON® with variable direction of view, lateral view



ENDOCAMELEON® with variable direction of view, isometric view




Diameter 10 mm, lengths 32 and 42 cm

Trocar size 11 mm

**Special Features:**

- Variable direction of view 0° – 120°
- HOPKINS® telescope with unique rod lens system
- Easy-to-use adjusting knob for selecting the direction of view
- Rigid sheath with a diameter of 10 mm



26003 AE

26003 AE


**ENDOCAMELEON® HOPKINS® Telescope,** diameter 10 mm, length 32 cm, autoclavable, variable direction of view from 0° – 120°, adjustment knob for selecting the desired direction of view, fiber optic light transmission incorporated, color code: gold

26003 AEE


**Same, length 42 cm**



**Fig. 1:** Choledochus duct in SILS cholecystectomy, angle of view, 70°



**Fig. 2:** Sigma anastomosis, angle of view, 90°



**Fig. 3:** Rectal anastomosis, angle of view, 30°



**Fig. 4:** Rectal anastomosis, angle of view, 80°

Source: Thomas Carus, M. D., Cuxhaven

# HOPKINS® Laparoscopes

enlarged view, diameter 10 mm

## Multiple puncture approach Trocar size 11 mm



26003 AA




26003 AA

**HOPKINS® II Straight Forward Telescope 0°**, enlarged view, diameter 10 mm, length 31 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green



26003 BA


**HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 10 mm, length 31 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red




26003 FA

**HOPKINS® II Telescope 45°**, enlarged view, diameter 10 mm, length 31 cm, **autoclavable**, fiber optic light transmission incorporated, color code: black

## Single puncture approach Trocar size 11 mm



26034 AAV



62039 AA



**NEW** 26034 AAV

**HOPKINS® Straight Forward Telescope 0°**, with parallel eyepiece, diameter 10 mm, length 23 cm, **autoclavable**, fiber optic light transmission incorporated, **with 6 mm working channel**, including adaptor with Silicone Leaflet Valve 26034 V



**NEW** 62039 AA

**HOPKINS® Straight Forward Telescope 0°**, enlarged view, with angled eyepiece, diameter 10 mm, length 25.5 cm, **autoclavable**, fiber optic light transmission incorporated and 5 mm working channel for rigid operating instruments



Trocars for use with 10 mm telescopes, see pages 56, 61-62 and 67



# HOPKINS® II Laparoscopes

Telescopes, diameter 5/3 mm

Telescopes, diameter 5 mm  
Trocar size 6 mm  
Multiple puncture approach



62046 AA

**HOPKINS® II Straight Forward Telescope 0°**, enlarged view, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green



62046 BA



**HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red



26046 FA

**HOPKINS® II Telescope 45°**, enlarged view, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: black

Telescopes, diameter 3 mm  
Trocar size 3.5 mm  
Multiple puncture approach



7220 AA

**HOPKINS® II Straight Forward Telescope 0°**, enlarged view, diameter 3 mm, length 14 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green



7220 BA


**HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 3 mm, length 14 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

5-061

Trocars for use with 3 and 5 mm telescopes see pages 57-58, 63-65 and 67-68

The TERNAMIAN **EndoTIP** System (**Endoscopic Threaded Imaging Port**) is a new KARL STORZ instrument for optically controllable access to the abdominal cavity. It replaces conventional trocars and, due to its unique design, offers a range of notable advantages for both surgeon and patient.

- The tissue is not cut through as it is when punctured with a conventional trocar, but is merely displaced. The integrity of the fasciae is preserved, and, therefore, that of the fascial closure mechanism.
- As sharp points or cutting trocars are not used, the risk of unintentional injuries to the patient is eliminated.
- Since the TERNAMIAN **EndoTIP** System is not pushed into the peritoneum by applying pressure, but is introduced by controlled rotation, injuries to organs caused by uncontrolled penetration of the abdominal wall are virtually impossible.
- In conjunction with a HOPKINS® Telescope, the surgeon has an excellent, enlarged image of the point of




penetration and is, therefore, able to carefully monitor the penetration under optical control. The various tissue layers can be identified precisely.


- The special design of the TERNAMIAN **EndoTIP** System ensures a secure hold in the abdominal wall and prevents gas loss at the point of penetration.




A small horizontal anterior rectus fascial incision is made under magnified video control. Blunt cannula tip (T).




The TERNAMIAN **EndoTIP** cannula is engaged in the anterior rectus fascial window and rotated clockwise.




Under clear visual control the cannula advances atraumatically without applying axial penetration force.




Rotation stretches radially, transposes and then lifts the anterior rectus fascia (F) and muscle (M) on to the cannula's outer thread.



The rectus muscle (M) is transposed laterally to expose the posterior rectus fascia (PF).




Further rotation will lift the posterior rectus fascia (PF) and advance the cannula to the pre-peritoneal space (P).



The thin peritoneal membrane transilluminates and appears blue-gray in color (P).




Continuous video control demonstrates different entry and exit layers clearly. This unique feature of the TERNAMIAN EndoTIP System allows the surgeon to stop advancing the cannula if bowel, vessels or adhesions are observed.




Access is under magnified visual control and atraumatic. Fascia and muscle fibers are not transected. The cannula's thread acts as a stabilizer to prevent slippage and discourage CO<sub>2</sub> leakage.




When the surgery is completed the cannula is rotated counter-clockwise. The different tissue layers disengage from the outer thread to converge in perfect apposition.




The peritoneum appears first (P), followed by the posterior rectus fascia (PF), rectus muscle (M), and then the anterior rectus fascia.




The rectus muscle fibers meet vertically to seal the peritoneal window, while the anterior rectus fascia shuts horizontally, to restore the natural tissue fiber orientation.



Since access is under visual control, hiatal or omental herniation along the cannula's tract is avoided. Lower half of anterior rectus fascia (F).



As sharp instruments or trocars are not used, muscle and fascial fibers are not transected and no tripolar fascial rent is sustained.







Tissue integrity is preserved at the access and exit site, maintaining the natural muscle shutter-mechanism.

# TERNAMIAN EndoTIP Cannulas


Size 11 and 13.5 mm



without insufflation adaptor  
for use with telescopes diameter 10 mm



<b>Size:</b>		<b>11 mm</b>
working length:		8.5 cm
color code:		green-white
	TERNAMIAN EndoTIP Cannula including: <b>Cannula</b> , with thread, without insufflation stopcock <b>Silicone Leaflet Valve</b>	<b>60103 MS</b> 30103 T7 30103 L1
<b>Accessories</b> 	<b>Telescope Stopper</b> , sterile, for single use, package of 12	30103 CS

with insufflation adaptor  
for use with telescopes diameter 10 mm




<b>Size:</b>		<b>11 mm</b>	<b>13.5 mm</b>
working length:		8.5 cm	11.5 cm
color code:		green-white	blue
	TERNAMIAN EndoTIP Cannula including: <b>Cannula</b> , with thread and rotating insufflation stopcock <b>Multifunctional Valve</b>	<b>60103 MVR</b> 30103 T5 30103 M1	<b>NEW 30108 MTR</b> 30108 T4 30108 M1
<b>Accessories</b> 	<b>Telescope Stopper</b> , sterile, for single use, package of 12	30103 CS	

Telescopes for use with TERNAMIAN EndoTIP Cannulas, size 11 mm see page 52  
Reducers see page 70



# TERNAMIAN EndoTIP Cannulas


Size 6 mm

**with and without insufflation adaptor**  
for use with instruments and telescopes  
diameter 5 mm





## with insufflation adaptor

<b>Size:</b>		<b>6 mm</b>
working length:		6.5 cm
color code:		black-white
	TERNAMIAN <b>EndoTIP Cannula</b> including: <b>Cannula</b> , with thread, with insufflation stopcock <b>Silicone Leaflet Valve</b>	<b>60120 MSS</b>  30120 T9  30120 L1
<b>Accessories</b> 	<b>Telescope Stopper</b> , sterile, for single use, package of 12	30160 CS



## without insufflation adaptor


<b>Size:</b>		<b>6 mm</b>
working length:		6.5 cm
color code:		black-white
	TERNAMIAN <b>EndoTIP Cannula</b> including: <b>Cannula</b> , with thread, without insufflation stopcock <b>Silicone Leaflet Valve</b>	<b>60120 MS</b>  30120 T7  30120 L1
<b>Accessories</b> 	<b>Telescope Stopper</b> , sterile, for single use, package of 12	30160 CS



Telescopes for use with TERNAMIAN EndoTIP Cannulas, size 6 mm see page 53

# TERNAMIAN EndoTIP Cannulas

Size 6 mm

with insufflation adaptor  
for use with instruments  
and telescopes diameter 5 mm




Size:		6 mm	
Working length:		8.5 cm	10.5 cm
Color code:		black-white	black
	<b>TERNAMIAN EndoTIP Cannula</b> including: <b>Cannula</b> , with thread and rotatable insufflation stopcock <b>Multifunctional Valve</b>	<b>60160 MVR</b>  30160 T5  30160 M1	<b>60160 MTR</b>  30160 T4  30160 M1
<b>Accessory</b> 	<b>Telescope Stopper</b> , sterile, for single use, package of 12	30160 CS	



Telescopes for use with TERNAMIAN EndoTIP Cannulas, size 6 mm see page 53

# TERNAMIAN EndoTIP Cannulas

Size 3.9 mm

with connector for insufflation  
for use with instruments, size 3/3.5 mm  
and telescopes, diameter 3.3 mm




<b>Size:</b>		<b>3.9 mm</b>
Working length:		5 cm
Color code:		red
	<b>TERNAMIAN EndoTIP Cannula</b> including: <b>Cannula</b> , with thread and insufflation stopcock <b>Automatic Valve</b>	<b>30117 MT</b>  30117 T2  30117 A1
<b>Accessory</b> 	<b>Telescope Stopper</b> , sterile, for single use, package of 12	30117 CS



Telescopes for use with TERNAMIAN EndoTIP Cannulas, size 3.9 mm see page 53

# Trocars

Size 13.5 mm

with insufflation adaptor  
for use with linear staplers from the company Ethicon (J&J)




<b>Size:</b>		<b>13.5 mm</b>
working length:		11.5 cm
color code:		blue
	<b>Trocar, with pyramidal tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Multifunctional Valve</b>	<b>30108 MP</b>  30108 H2 30108 P 30108 M1
	<b>Trocar, with blunt tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Multifunctional Valve</b>	<b>30108 MA</b>  30108 H2 30108 A 30108 M1






# Trocars

Size 11 mm

with insufflation adaptor



Size:		11 mm			
		8.5 cm green-white	10.5 cm green	8.5 cm green-white	10.5 cm green
working length:					
color code:					
		Multifunctional Valve		Automatic Valve	
	<b>Trocar, with conical tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Valve</b>	<b>62103 WY</b> 30103 H1 30103 Y 30103 M1	<b>62103 MC</b> 30103 H2 30103 C 30103 M1	<b>62103 FY</b> 30103 H1 30103 Y 30103 A1	<b>62103 AC</b> 30103 H2 30103 C 30103 A1
	<b>Trocar, with pyramidal tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Valve</b>	<b>62103 WZ</b> 30103 H1 30103 Z 30103 M1	<b>62103 MP</b> 30103 H2 30103 P 30103 M1	<b>62103 FZ</b> 30103 H1 30103 Z 30103 A1	<b>62103 AP</b> 30103 H2 30103 P 30103 A1
	<b>Trocar, with blunt tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Valve</b>	<b>62103 WX</b> 30103 H1 30103 X 30103 M1	<b>62103 MA</b> 30103 H2 30103 A 30103 M1	<b>62103 FX</b> 30103 H1 30103 X 30103 A1	<b>62103 AA</b> 30103 H2 30103 A 30103 A1



**Reducers** see page 70




# Trocars

Size 11 mm

without insufflation adaptor

color code: green-white, size 11 mm, working length 6.5 cm




<b>Size:</b>		<b>11 mm</b>
working length:		6.5 cm
color code:		green-white
	<b>Trocar, with conical tip</b> including: <b>Cannula, with thread</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>60123 EW</b>  30123 E1 30123 W 30123 L1
	<b>Trocar, with pyramidal tip</b> including: <b>Cannula, with thread</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>60123 EZ</b>  30123 E1 30123 Z 30123 L1
	<b>Trocar, with blunt tip</b> including: <b>Cannula, with thread</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>60123 EV</b>  30123 E1 30123 V 30123 L1

# Trocars

Size 6 mm

with connector for insufflation




Size:		6 mm				
		8.5 cm black-white	10.5 cm black	8.5 cm black-white	10.5 cm black	
Working length:						
Color code:						
		Multifunctional Valve		Automatic Valve		
	<b>Trocar, with conical tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Valve</b>	<b>62160 WY</b> 30160 H1 30160 Y 30160 M1	<b>62160 MC</b> 30160 H2 30160 C 30160 M1	<b>62160 FY</b> 30160 H1 30160 Y 30160 A1	<b>62160 AC</b> 30160 H2 30160 C 30160 A1	
		<b>Trocar, with pyramidal tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Valve</b>	<b>62160 WZ</b> 30160 H1 30160 Z 30160 M1	<b>62160 MP</b> 30160 H2 30160 P 30160 M1	<b>62160 FZ</b> 30160 H1 30160 Z 30160 A1	<b>62160 AP</b> 30160 H2 30160 P 30160 A1
			<b>Trocar, with blunt tip</b> including: <b>Cannula</b> <b>Trocar only</b> <b>Valve</b>	<b>62160 WX</b> 30160 H1 30160 X 30160 M1	<b>62160 MA</b> 30160 H2 30160 A 30160 M1	<b>62160 FX</b> 30160 H1 30160 X 30160 A1





Telescopes for use with trocars size 6 mm see page 53

# Trocars <sup>NEW</sup>

Size 6 mm

with insufflation adaptor  
for use with instruments size 5 mm





Size:		6 mm		
		5 cm black-green	8.5 cm black-white	10.5 cm black
working length:				
color code:				
	<b>Trocar, with conical tip</b> including: <b>Cannula, with LUER-Lock Connector</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>30160 GYK</b>  30160 G3 30160 YK 30160 L1	<b>30160 GYM</b>  30160 G2 30160 Y 30160 L1	<b>30160 GC</b>  30160 G1 30160 C 30160 L1
	<b>Trocar, with pyramidal tip</b> including: <b>Cannula, with LUER-Lock Connector</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>30160 GZK</b>  30160 G3 30160 ZK 30160 L1	<b>30160 GZM</b>  30160 G2 30160 Z 30160 L1	<b>30160 GP</b>  30160 G1 30160 P 30160 L1
	<b>Trocar, with blunt tip</b> including: <b>Cannula, with LUER-Lock Connector</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>30160 GXK</b>  30160 G3 30160 XK 30160 L1	<b>30160 GXM</b>  30160 G2 30160 X 30160 L1	<b>30160 GA</b>  30160 G1 30160 A 30160 L1
	<b>Trocar</b> including: <b>Cannula, with LUER-Lock Connector</b> <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>30160 GRK</b>  30160 G3 30160 RK 30160 L1	<b>30160 GRM</b>  30160 G2 30160 R 30160 L1	<b>30160 GB</b>  30160 G1 30160 B 30160 L1

Telescopes for use with trocars size 6 mm see page 53




# Trocars

Size 6 mm

color code: black-white, working length 6 cm



### without insufflation adaptor

Size:		6 mm
working length:		6 cm
color code:		black-white
	<b>Trocar</b> , with conical tip including: <b>Cannula</b> , with thread <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>62120 EH</b>  30120 E1 30120 H 30120 L1
	<b>Trocar</b> , with pyramidal tip including: <b>Cannula</b> , with thread <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>62120 EK</b>  30120 E1 30120 K 30120 L1
	<b>Trocar</b> , with blunt tip including: <b>Cannula</b> , with thread <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>62120 EO</b>  30120 E1 30120 O 30120 L1

### with insufflation adaptor

Size:		6 mm
working length:		6 cm
color code:		black-white
	<b>Trocar</b> , with blunt tip including: <b>Cannula</b> , with thread <b>Trocar only</b> <b>Silicone Leaflet Valve</b>	<b>30120 EOX</b>  30120 EX1 30120 O 30120 L1

5-06

# Open Laparoscopy

Size 6 mm

## For open laparoscopy

to establish pneumoperitoneum without the use of a VERESS needle or sharp trocar





30160 C3


- 30160 C3     **Ribbed Cone**, with 2 discs for fixation of sutures, for trocars size 6 mm, color code: black
- 30103 C3     **Ribbed Cone**, with 2 discs for fixation of suture, for trocars size 11 mm, color code: green

# Trocars for Thoracoscopy


Size 6 and 11 mm

**without insufflation adaptor**  
**for use with curved surgical instruments**



<b>Size:</b>		<b>6 mm</b>	<b>11 mm</b>
working length:		8.5 cm	8.5 cm
color code:		black	green
	<b>Trocar, with blunt tip</b> including: <b>Cannula</b> , flexible, with thread and silicone leaflet valve <b>Trocar only</b>	<b>30120 NOL</b>	<b>30123 NOL</b>
		30120 NL	30123 NL
		30120 NO	30123 O
<b>Accessories</b>	<b>Plastic Cannula, autoclavable,</b> for flexible trocars, package of 5	30120 X	30123 X




<b>Size:</b>		<b>6 mm</b>	
working length:		4 cm	6 cm
	<b>Trocar, with blunt tip</b> including: <b>Cannula</b> , with thread <b>Trocar only</b>	<b>40160 AB</b>	<b>40160 CD</b>
		40160 A	40160 C
		40160 B	40160 D


# Trocars

Size 3.9 and 3.5 mm

For use in small animals less than 10 kg


For use with 2.7 mm telescope in conjunction with Protection Sheath 64018 US



<b>Size:</b>		<b>3.9 mm</b>
working length:		5 cm
color code:		red
	<b>Trocar</b> , with pyramidal tip including:	<b>NEW 62117 GK</b>
	<b>Cannula</b> , with outer thread, with LUER-lock connector for insufflation	30117 G4
	<b>Trocar only</b>	30117 PK
	<b>Silicone Leaflet Valve</b>	30117 L1

Lightweight trocar model for use with instruments and HOPKINS® telescopes 3 mm and 2 mm



<b>Size:</b>		<b>3.5 mm</b>
working length:		5 cm
color code:		green
	<b>Trocar</b> , with pyramidal tip including:	<b>NEW 62114 GK</b>
	<b>Cannula</b> , with LUER-Lock cone, with thread	62114 G1
	<b>Trocar only</b>	30114 K
	<b>Silicone Leaflet Valve</b>	30114 L1


HOPKINS® Telescopes for use with Trocars, size 3.5 mm see page 53





# Trocar and VERESS Needle

For use in small animals less than 10 kg

For use with 2 mm instruments



<b>Size:</b>	<b>2.5 mm</b>
working length:	4 cm
color code:	blue
 <p><b>Trocar</b>, with pyramidal tip including: <b>Cannula</b>, with outer thread, with LUER-lock connector for insufflation <b>Trocar only</b> <b>Silicone Leaflet Valve</b></p>	<p><b>NEW 62603 GK</b></p> <p>62603 G1</p> <p>11603 K</p> <p>11603 L1</p>



62120 J

62120 J

**VERESS Pneumoperitoneum Needle**,  
with spring-loaded blunt stylet, LUER-Lock,  
**autoclavable**, diameter 2.1 mm, length 10 cm

# Accessories

## Reducers and Reduction Sleeves

### Reduction Sleeves



30140 DB

Instrument Diameter	Cannula Outer Diameter – Color Code:	
	6 mm black	11 mm green
3 mm	30140 KA	30140 DA
5 mm	–	30140 DB

### Reducers



62141 AA



62141 DB

62141 AA

**Reducer, 6/3 mm**

62141 DB

**Reducer, 11/5 mm, flip style**  
(optional for larger ovaries)

### Reducers for use with trocars and TERNAMIAN EndoTIP cannulas with silicone leaflet valve



30160 RE

**NEW**

30160 RE

**Reducer, 6/3 mm**

**NEW**

30103 RE

**Reducer, 11/5 mm**

**NEW**

30141 HB

**Reducer, 13/5 mm and 13.5/5 mm**

### Thread Sleeves for Trocars



30160 S

30160 S

**Thread Sleeve, for trocar size 6 mm,**  
color code: black

30103 S

**Thread Sleeve, for trocar size 11 mm,**  
color code: green

### Distance holders to limit the depth of puncture



30160 D

30160 D

**Distance Holder, for trocars size 6 mm,**  
color code: black

30103 D

**Distance Holder, for trocars size 11 mm,**  
color code: green

### Sealing Caps, with color code



30100 XA


30100 XA


**Sealing Cap, size 5 mm, autoclavable,** package of 5,  
for use with trocars size 6 mm and extractors/reduction  
sleeves in conjunction with instruments,  
color code: black

30100 XB

**Sealing Cap, size 10 mm, autoclavable,** package of 5,  
for use with trocars size 11 mm and extractors/reduction  
sleeves in conjunction with instruments,  
color code: green

# DISSECTING AND GRASPING FORCEPS BIOPSY FORCEPS AND SCISSORS





**CLICK<sup>line</sup> Instruments, consisting of:**

- Handle
- Outer sheath
- Working insert

**Nothing could be simpler!**

The **CLICK<sup>line</sup>** series continues the development of instruments which have proven their value for years. The simplicity of handling and ease with which the instruments can be cleaned has been improved, in particular with respect to assembly and disassembly. When the rear handle is positioned horizontally, the

handle can be separated from the outer sheath and the working insert at the press of a button. Reassembly is equally reliable and quick. The high frequency connection is mounted in a 45° angle on the upper side of the handle, thereby ergonomically guiding the high frequency cable away from the field of operation.

**Available in size 2 to 10 mm, CLICK<sup>line</sup> instruments can be completely disassembled into separate components:**

- Handle
- Outer sheath/outer sheath with working insert
- Working insert

**This unique, reusable two/three-piece design offers the surgeon the following benefits:**

- Available in size 2 mm, length 20 cm; size 3/3.5 mm, length 20 or 30 cm; size 5 and 10 mm, length 30, 36 and 43 cm
- Choice of handle styles
- Fully rotational 360° sheath facilitates easy access in all clinical situations
- No hidden spaces that can trap blood or tissue debris
- Can be dismantled at the press of a button – reducing instrument cleaning time considerably
- Completely autoclavable design
- Cost-effective – reusable instruments reduce O.R. costs per case and simplify inventory management, eliminating the need to store large quantities of disposable instruments
- Environmentally correct, i.e. if damage occurs, only the component with the defect needs to be replaced – not the entire instrument
- Convenient and ergonomic handling
- Cleaning port allows the instrument to be cleaned without disassembly

**Due to the modular CLICK<sup>line</sup> system, the user can individually assemble the desired instrument at any time.**


# User Information

For easy location of the correct catalog number for the required CLICK'line<sup>®</sup> Instrument

## How to find the required instrument:

- (a) Select the required **instrument size** and type of **jaws**
- (c) Select the required **handle style**
- (d) The catalog number for the completely assembled instrument can be found at the point where the selected parameters intersect (see table below).

### Example:



## Dissecting and Grasping Forceps

CLICK'line – rotating, dismantling, with or without connector pin for unipolar coagulation

Size 5 mm, length 36 cm

↓  
unipolar

Length	Handle		
	(b) 33121	33122	33123
36 cm			
Working Insert No.	Complete Instrument		
33310 MF	(c) 66321 MF	66322 MF	66323 MF

CLICK'line Dissecting and Grasping Forceps, atraumatic, especially designed for grasping and removing hollow organs

— 24 —

— — Length of jaws in mm (approx. measurements)

### Please note:

Among CLICK'line instruments only the **single component parts** are numbered. The number for the **complete instrument** is not etched on it. Instruments with **insulated handles with connector** for unipolar coagulation are shown against the **red background**. Instruments with **handles without connector pin** for unipolar coagulation are shown against the **blue background**. The colour **green** indicates the **working inserts**.


# Handles, Outer Sheaths

for Dissecting and Grasping Forceps,

CLICK'liné – rotating, dismantling, with connector pin for unipolar coagulation

Plastic Handles, with and without ratchets

Made of high-quality plastics, meeting even the highest of safety requirements



33121 **CLICK'liné Plastic Handle, without ratchet**



33122 **CLICK'liné Plastic Handle, with ratchet**



33123 **CLICK'liné Plastic Handle, with hemostat style ratchet**

Metal Handles, with and without ratchets



33131 **CLICK'liné Metal Handle, rotating, without ratchet**




33132 **CLICK'liné Metal Handle, rotating, with ratchet**



33133 **CLICK'liné Metal Handle, rotating, with hemostat style ratchet**

Size 5 mm

Metal Outer Sheath, insulated



33300

**CLICK'liné Metal Outer Sheath, insulated, with LUER-Lock connector for cleaning, size 5 mm, length 36 cm**

33400

**Same, length 43 cm**

Working Inserts for insulated instruments, size 5 mm see pages 75-76

# Dissecting and Grasping Forceps

CLICK'line – rotating, dismantling, with connector pin for unipolar coagulation

Size 5 mm, length 36 and 43 cm



Length	Handle		
	33121	33122	33123
36 cm			
43 cm			


Working Insert No.	Complete Instrument		
33310 MF   — 24 —	66321 MF	66322 MF	66323 MF
<b>CLICK'line Dissecting and Grasping Forceps</b> , atraumatic, specially designed for grasping and removing hollow organs			
33310 ML   — 22 —	66321 ML	66322 ML	66323 ML
<b>CLICK'line Dissecting and Grasping Forceps</b> , curved			
33310 ME   — 14 —	66321 ME	66322 ME	66323 ME
<b>CLICK'line Grasping Forceps</b> , jaws with multiple teeth, width of jaws 4.8 mm, for atraumatic and accurate grasping			
33310 MG	66321 MG	66322 MG	66323 MG
33410 MG	–	<b>NEW</b> 60422 MG	<b>NEW</b> 60423 MG
  — 14 —	<b>CLICK'line Grasping Forceps</b> , “tiger jaws”, 2 x 4 teeth, for grasping and removing solid organs, in particular during adhesiolysis		
33310 BA	66321 BA	66322 BA	66323 BA
33410 BA	<b>NEW</b> 60421 BA	<b>NEW</b> 60422 BA	<b>NEW</b> 60423 BA
  — 16 —	<b>CLICK'line Grasping Forceps</b> , atraumatic, with multiple teeth, fenestrated, single action jaws		
33310 R	<b>NEW</b> 33321 R	<b>NEW</b> 33322 R	<b>NEW</b> 33323 R
  — 16 —	<b>CLICK'line Dissecting and Grasping Forceps</b> , right angled		









10-002

# Scissors, Biopsy Forceps

CLICK'line – rotating, insulated, with connector pin for unipolar coagulation

Size 5 mm, length 36 and 43 cm



Length	Handle	
	33121	33123
36 cm		
43 cm		
Working Insert No.	Complete Instrument	
34310 DZ	66321 DZ	–
34410 DZ	<b>NEW</b> 60421 DZ	–
  — 13 —	<b>CLICK'line Biopsy Forceps</b> , with 2 teeth	
34310 MA	66321 MA	–
  — 20 —	<b>CLICK'line Scissors</b> , serrated, curved, spoon-shaped jaws	
34310 MW	66321 MW	–
  — 15 —	<b>CLICK'line Scissors</b> , serrated, curved, conical	
34310 EH	66321 EH	–
  — 10 —	<b>CLICK'line Hook Scissors</b>	
34310 DB	66321 DB	66323 DB
34410 DB	<b>NEW</b> 60421 DB	–
  — 14 —	<b>CLICK'line Biopsy Forceps</b>	
34310 DH	66321 DH	–
33410 DH	<b>NEW</b> 60421 DH	–
  — 10 —	<b>CLICK'line Biopsy Punch</b> , through-cutting	



# Mini-Laparoscopy Instruments

for animals less than 10 kg, CLICK'line® – rotating, dismantling,  
with connector pin for unipolar coagulation

Size 3 mm, length 20 cm

For use with trocars size 3.5 mm



Length	Handle		
	33121	33122	33123
20 cm			

## Double Action Jaws

Outer Tube No. with Working Insert	Complete Instrument		
30310 MDS	30321 MDS	30322 MDS	30323 MDS
   10	<b>CLICK'line® KELLY Dissecting and Grasping Forceps</b>		
30310 ULS	30321 ULS	30322 ULS	30323 ULS
   11	<b>CLICK'line® REDDICK-OLSEN Dissecting and Grasping Forceps, robust</b>		
30310 MWS	30321 MWS	–	–
   10	<b>CLICK'line® Scissors, serrated, curved, conical</b>		

## Single Action Jaws

30310 DBS	30321 DBS	–	–
   10	<b>CLICK'line® BLAKESLEY Biopsy Forceps</b>		

Instruments also available with a diameter of 2 mm.

## Size 10 and 5 mm, length 36 cm



33533 DU

**CLICK'line® DUVAL Grasping Forceps**, rotating, double action jaws, size 10 mm, length 36 cm including:

**Metal Handle**, with hemostat style ratchet

**Outer Sheath**, insulated

**Forceps Insert**



**NEW** 62046 GS

**Optical Biopsy Forceps**, outer diameter 10 mm, capable of procuring biopsy samples 1 gram or more, for use with telescope 62046 AA



62175 T



62175 T

**Palpation Probe**, with cm-markings, size 5 mm, length 36 cm



62175 P



62175 P


**Injection Cannula**, LUER-Lock, needle diameter 1.2 mm, length 36 cm

62175 PD


**Same**, diameter 0.8 mm

62175 R

**Same**, diameter 1.6 mm



26173 BN



26173 BN

**Suction and Irrigation Tube**, with lateral holes, with two-way stopcock for single hand control, size 5 mm, length 36 cm

**NEW** 30623 F

**Fan Retractor**, dismantling, distendable, size 5 mm, length 36 cm

**NEW** 30623 FP

**Same**, size 10 mm

## Size 3 mm, length 20/30 cm



26167 TS



26167 TS

**Palpation Probe**, with cm-markings, size 3 mm, length 20 cm

26167 T

**Same**, length 30 cm

**High Frequency Cords** see chapter 11, UNITS  
**Components/Spare Parts** see chapter 14


# TAKE-APART® Instruments

for Bipolar Coagulation

Size 5 mm

Multiple puncture approach

Operating instruments for use with trocars, size 6 mm, length 33 cm



66176 HA



66176 HA

**TAKE-APART® Bipolar Coagulating Forceps,**  
insulated, with connector pin for bipolar  
coagulation, size 5 mm

including:

**Handle**

**Outer Sheath**

**Forceps Insert**



66176 HB

**TAKE-APART® Bipolar Micro Forceps,**  
insulated, with connector pin for bipolar  
coagulation, size 5 mm

including:

**Handle**

**Outer Sheath**

**Forceps Insert**



66176 HL

**TAKE-APART® Bipolar Coagulating Forceps,**  
insulated, with connector pin for bipolar  
coagulation, size 5 mm

including:

**Handle**

**Outer Sheath**

**Forceps Insert**

10-00:


**Components/Spare Parts** see chapter 14

**TAKE-APART® Instruments** *NEW*  
for Bipolar Coagulation

Size 3 mm

Multiple puncture approach

Operating instruments for use with trocars, size 3.5 mm, length 20 cm



**TAKE-APART® Bipolar Coagulating Forceps**,  
width of jaws 1 mm, size 3 mm, length 20 cm  
including:  
**Bipolar Ring Handle**  
**Outer Sheath**  
**Forceps Insert**

Components/Spare Parts see chapter 14


# Instruments

for Bipolar Coagulation and Cutting

Size 5 mm

Multiple puncture approach

Operating instruments for use with trocars, size 6 mm, length 36 cm



38321 MZ



38321 MZ

**RoBi® Scissors**, CLERMONT-FERRAND model, straight, serrated, for cutting and bipolar coagulation, size 5 mm, length 36 cm

including:

**RoBi® Ring Handle**

**RoBi® Metal Outer Sheath**, insulated

**RoBi® Scissors Insert**

**Components/Spare Parts** see chapter 14

## KOH Macro Needle Holders

The ergonomic properties of the new 5 mm KOH macro needle holder allow for precise holding of the needles and the suture. The precision jaws and the new locking mechanism enable easy and precise guidance of the

needle with suture of size 0/0 to 7/0. Specially designed to meet the demanding requirements of laparoscopic micro surgery, the KOH ultra micro instrument set enables suturing structures of even 0.5 mm.


### Special Features:

- **Tungsten carbide inserts for firm and reliable needle positioning, as well as long durability**
- **Ease of operation, precise adjustable ratchet for easy and safe positioning of the needle**
- **Ergonomically designed handle for relaxed working**

### Size 5 mm

#### Multiple puncture approach

Operating instruments, **length 33 and 43 cm**, for use with trocars size 6 mm



**NEW**

26173 KPF

**KOH Macro Needle Holder**, with tungsten carbide insert, ergonomic pistol handle, with disengageable ratchet, ratchet position left, jaws straight, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1

26178 KPF

**Same**, length 43 cm



**NEW**

26173 KPL

**KOH Macro Needle Holder**, with tungsten carbide insert, ergonomic pistol handle, with disengageable ratchet, ratchet position left, jaws curved to left, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1

26178 KPL

**Same**, length 43 cm



**NEW**

26173 KPR

**KOH Macro Needle Holder**, with tungsten carbide insert, ergonomic pistol handle with disengageable ratchet, ratchet position right, jaws curved to right, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1

26178 KPR


**Same**, length 43 cm

## KOH Macro Needle Holders

### Size 5 mm

#### Multiple puncture approach

Operating instruments, **length 33 and 43 cm**,  
for use with trocars size 6 or 11 mm with reduction sleeve



26173 KC




26173 KL

**KOH Macro Needle Holder**, with tungsten carbide insert, ergonomic handle with ratchet, ratchet position on top, jaws curved to left, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1



26173 KC

**KOH Macro Needle Holder**, with tungsten carbide insert, ergonomic handle with ratchet, ratchet position underneath, jaws curved right, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1



26173 ML



26173 ML

**KOH Macro Needle Holder**, straight handle, with ratchet, jaws curved to left, with tungsten carbide inserts, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0

26178 ML

**Same**, length 43 cm



26173 MC

**KOH Macro Needle Holder**, straight handle, with ratchet, jaws curved to right, with tungsten carbide inserts, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0

26178 MC

**Same**, length 43 cm


### Size 3 mm



**NEW** 26167 FNS

**KOH Ultramicro Needle Holder**, with tungsten carbide insert, jaws curved to left, straight handle, with disengageable ratchet, length 20 cm, for use with suture material 7/0, 8/0 (Ethicon) and needle size BV 175-6


# Endo-Suture, Endo-Ligature, Ovariectomy Hook



30140 KA

30140 KA

**Reduction Sleeve**, reusable, instrument diameter 3 mm, trocar cannula outer diameter 6 mm, color code: black



26172 AE

26172 AE

**Endo-Loop Ligature**, with ROEDER knot, for bleeding stumps, with absorbable synthetic suture, for single use, sterile, USP 0, size 3 mm, length 33 cm, package of 12

**Knot Tier**  
for use with trocar 6 mm



**NEW**

26596 SK

KOECKERLING **Knot Tier**, for extracorporeal knotting, size 5 mm, length 36 cm


**Ovariectomy Hook**  
secures ovary transabdominally during dissection and removal



62175 H

62175 H

TANKERSLEY **Ovariectomy Hook**, with heavy handle, for percutaneous use, working length 8.5 cm



62175 HK

**NEW**

62175 HK

**Ovariectomy Hook**, with heavy handle, length 8.5 cm, with small hook for percutaneous use in small dogs and cats



# Diagnostic Laparoscopy

## Standard Set

62120 J	<b>VERESS Pneumoperitoneum Needle</b> , with spring-loaded blunt stylet, LUER-Lock, <b>autoclavable</b> , diameter 2.1 mm, length 10 cm
60160 MTR	<b>TERNAMIAN EndoTIP Cannula</b> , size 6 mm, color code: black including: <b>Cannula</b> , with thread and rotating insufflation stopcock, working length 10.5 cm <b>Multifunctional Valve</b>
62160 WZ	<b>Trocar</b> , size 6 mm including: <b>Trocar only</b> , with pyramidal tip <b>Cannula</b> , without valve, with insufflation stopcock, working length 8.5 cm <b>Multifunctional Valve</b>
62046 BA	<b>HOPKINS® II Forward-Oblique Telescope 30°</b> , enlarged view, diameter 5 mm, length 29 cm, <b>autoclavable</b> , fiber optic light transmission incorporated, color code: red
62175 T	<b>Palpation Probe</b> , with cm-markings, size 5 mm, length 36 cm
66321 DB	<b>CLICK'line® Dissecting and Biopsy Forceps</b> , rotating, single action jaws, with connector pin for unipolar coagulation, size 5 mm, length 36 cm including: <b>Plastic Handle</b> , without ratchet <b>Metal Outer Sheath</b> , insulated <b>Forceps Insert</b>
66323 ML	<b>CLICK'line® KELLY Dissecting and Grasping Forceps</b> , rotating, dismantling, with connector pin for unipolar coagulation, with irrigation connector for cleaning, double action jaws, long, size 5 mm, length 36 cm including: <b>Plastic Handle</b> , with hemostat style ratchet <b>Metal Outer Sheath</b> , insulated <b>Forceps Insert</b>
66321 DH	<b>CLICK'line® Biopsy Punch</b> , single action jaws, through-cutting, size 5 mm, length 36 cm including: <b>Plastic Handle</b> , without ratchet <b>Metal Outer Sheath</b> , insulated <b>Forceps Insert</b>

62160 MTR	<b>TERNAMIAN EndoTIP Cannula</b>
2x 62160 WZ	<b>Trocar</b> , size 6 mm including: <b>Trocar only</b> , with pyramidal tip <b>Cannula</b> , without valve, with insufflation stopcock, working length 8.5 cm <b>Multifunctional Valve</b>
62046 BA	<b>HOPKINS® II Forward-Oblique Telescope 30°</b> , enlarged view, diameter 5 mm, length 29 cm, <b>autoclavable</b> , fiber optic light transmission incorporated, color code: red
66323 ML	<b>CLICK'linē KELLY Dissecting and Grasping Forceps</b> , rotating, dismantling, with connector pin for unipolar coagulation, with irrigation connector for cleaning, double action jaws, long, size 5 mm, length 36 cm including: <b>Plastic Handle</b> , with hemostat style ratchet <b>Metal Outer Sheath</b> , insulated <b>Forceps Insert</b>
66323 BA	<b>CLICK'linē BABCOCK Grasping Forceps</b> , rotating, with connector pin for unipolar coagulation, atraumatic, with multiple teeth, fenestrated, single action jaws, size 5 mm, length 36 cm including: <b>Plastic Handle</b> , with hemostat style ratchet <b>Metal Outer Sheath</b> , insulated <b>Forceps Insert</b>
66321 MW	<b>CLICK'linē Scissors</b> , rotating, insulated, with connector pin for unipolar coagulation, serrated, curved, conical, double action jaws, size 5 mm, length 36 cm including: <b>Plastic Handle</b> , without ratchet <b>Metal Outer Sheath</b> , insulated <b>Scissors Insert</b>
66176 HA	<b>TAKE-APART® Bipolar Coagulating Forceps</b> , insulated, with connector pin for bipolar coagulation, size 5 mm including: <b>Handle</b> <b>Outer Sheath</b> <b>Forceps Insert</b>
26173 KL	<b>KOH Macro Needle Holder</b> , with tungsten carbide insert, ergonomic handle with ratchet, ratchet position on top, jaws curved to left, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1
26173 KC	<b>KOH Macro Needle Holder</b> , with tungsten carbide insert, ergonomic handle with ratchet, ratchet position underneath, jaws curved right, size 5 mm, length 33 cm, for use with suture material size 0/0 to 7/0 and needle sizes BV, SH or CT-1
26596 SK	<b>KOECKERLING Knot Tier</b> , for extracorporeal knotting, size 5 mm, length 36 cm
26172 AE	<b>Endo-Loop Ligature</b> , with ROEDER knot, for bleeding stumps, with absorbable synthetic suture, for single use, sterile, USP 0, size 3 mm, length 33 cm, package of 12

## Multi-Port LAPspay Instrument Set

- 62046 AA **HOPKINS® II Straight Forward Telescope 0°**, enlarged view, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green
- 60120 MS TERNAMIAN **EndoTIP Cannula**, size 6 mm, color code: black-white including:  
**Cannula**, with thread, working length 6.5 cm  
**Silicone Leaflet Valve**
- 60120 MSS **Same**, with insufflation stopcock
- 66322 MG **CLICKline® Grasping Forceps**, "tiger jaws", rotating, with connector pin for unipolar coagulation, single action jaws, 2 x 4 teeth, size 5 mm, length 36 cm including:  
**CLICKline® Plastic Handle**, with ratchet  
**Metal Outer Sheath, insulated**  
**Forceps Insert**
- 62120 J VERESS **Pneumoperitoneum Needle**, with spring-loaded blunt stylet, LUER-Lock, **autoclavable**, diameter 2.1 mm, length 10 cm
- 62175 H TANKERSLEY **Ovariectomy Hook**, with heavy handle, for percutaneous use, working length 8.5 cm
- 62103 FZ **Trocar**, with pyramidal tip, automatic valve (optional for larger ovaries), size 11 mm, working length 8.5 cm
- 62141 DB **Reducer**, 11/5 mm, flip style (optional for larger ovaries)

## Single Port LAPspay Instrument Set

- 26034 AAV **HOPKINS® Straight Forward Telescope 0°**, with parallel eyepiece, diameter 10 mm, length 23 cm, **autoclavable**, fiber optic light transmission incorporated, **with 6 mm working channel**, including adaptor with Silicone Leaflet Valve 26034 V
- 62103 FZ **Trocar**, with pyramidal tip, automatic valve (optional for larger ovaries), size 11 mm, working length 8.5 cm
- 60120 MSS **Same**, with insufflation stopcock
- 60422 MG **CLICKline® Grasping Forceps**, "tiger jaws", rotating, with connector pin for unipolar coagulation, single action jaws, 2 x 4 teeth, size 5 mm, length 43 cm including:  
**CLICKline® Plastic Handle**, with ratchet  
**Outer Sheath, insulated**  
**Forceps Insert**
- 62120 J VERESS **Pneumoperitoneum Needle**, with spring-loaded blunt stylet, LUER-Lock, **autoclavable**, diameter 2.1 mm, length 10 cm
- 62175 H TANKERSLEY **Ovariectomy Hook**, with heavy handle, for percutaneous use, working length 8.5 cm

# Mini-Laparoscopic Surgery

for animals less than 10 kg

7220 AA	<b>HOPKINS® II Straight Forward Telescope 0°</b> , enlarged view, diameter 3 mm, length 14 cm, <b>autoclavable</b> , fiber optic light transmission incorporated, color code: green
3x 62114 GK	<b>Trocar</b> , size 3.5 mm, color code: green including: <b>Cannula</b> , with LUER-Lock cone, with thread, working length 5 cm <b>Trocar only</b> , with pyramidal tip <b>Silicone Leaflet Valve</b>
26167 TS	<b>Palpation Probe</b> , with cm-markings, size 3 mm, length 20 cm
30321 DBS	<b>CLICK<sup>line</sup> BLAKESLEY Dissecting and Biopsy Forceps</b> , rotating, with connector pin for unipolar coagulation, with LUER-Lock adaptor for cleaning, single action jaws, size 3 mm, length 20 cm including: <b>Plastic Handle</b> , insulated, without ratchet <b>Outer Sheath</b> , with working insert, insulated
30323 ULS	<b>CLICK<sup>line</sup> REDDICK-OLSEN Dissecting and Grasping Forceps</b> , rotating, with connector pin for unipolar coagulation, with LUER-Lock adaptor for cleaning, double action jaws, robust, size 3 mm, length 20 cm including: <b>Plastic Handle</b> , with hemostat style ratchet <b>Outer Sheath</b> , with working insert, insulated
30321 MWS	<b>CLICK<sup>line</sup> Scissors</b> , rotating, with connector pin for unipolar coagulation, with LUER-Lock adaptor for cleaning, double action jaws, serrated, curved, conical, size 3 mm, length 20 cm including: <b>Plastic Handle</b> , without ratchet <b>Outer Sheath</b> , with scissors insert, insulated
26184 HCS	<b>TAKE-APART® Bipolar Coagulating Forceps</b> , width of jaws 1 mm, size 3 mm, length 20 cm including: <b>Bipolar Ring Handle</b> <b>Outer Sheath</b> <b>Forceps Insert</b>
62120 J	<b>VERESS Pneumoperitoneum Needle</b> , with spring-loaded blunt stylet, LUER-Lock, <b>autoclavable</b> , diameter 2.1 mm, length 10 cm

# HANDLES AND INSTRUMENTS FOR SUCTION AND IRRIGATION




# Modular Handles for Suction and Irrigation




30803

**Handle**, without irrigation and suction connector, **autoclavable**, for use with 5 mm coagulating and dissecting electrodes with suction channel



30804

**Handle with Trumpet Valve**, for irrigation or suction, **autoclavable**, for use with 5 mm suction-coagulation cannulas and 5 mm suction/irrigation tubes




30805

**Handle with Two-Way Stopcock**, for irrigation and suction, **autoclavable**, for use with 5 mm suction and irrigation tubes



69810

**Handle**, for irrigation and suction, **autoclavable**, for use with Tubing Set 69811, in conjunction with 5 mm suction and irrigation tubes and catheters any size with LUER-Lock connection



69811

**Instrument Tubing Set**, autoclavable, connects bottles to Handle 69810, for use with VETPUMP® II




9896890




**LUER-Lock Connector**, for use with Handle 69810

**Suction and irrigation system** see chapter 11, UNITS

# Suction and Irrigation Tubes

Size 5 mm

Length	Instrument
30 cm	
36 cm	
43 cm	

Distal End	Instrument
	37260 LH
	37360 LH
	37460 LH
	37260 SC
	37360 SC
	37460 SC
	37260 CP
	37360 CP
	37460 CP






Accessories for use with suction/irrigation tubes see chapter 11, UNITS

# Coagulating Suction Tubes

for unipolar coagulation, size 5 mm

For use with handle

Length	Instrument
36 cm	

Distal End	Instrument
	62370 SC <b>Coagulating Suction Tube</b>
	62370 DW <b>Coagulating Suction Tube, W-shaped</b>
	62370 DB <b>Coagulating Suction Tube, blunt</b>
	62370 DL <b>Coagulating Suction Tube, L-shaped</b>
	62370 DU <b>Coagulating Suction Tube, U-shaped</b>

**Accessories for use with suction/irrigation tubes** see chapter 11, UNITS  
**High frequency cables and high frequency surgical units** see chapter 11, UNITS



# ARTHROSCOPES AND ACCESSORIES





# Canine Arthroscopy

Most articular interventions in humans and horses are performed by means of arthroscopy. Due to recent advances in instruments and technology, arthroscopy is gaining widespread acceptance in the articular surgery of dogs. The advantages of arthroscopy over conventional arthrotomy are improved access to joint struc-

tures, magnification of lesions, arthroscopic lavage, and a more rapid recovery of the animal patient. KARL STORZ has developed the instruments described here in conjunction with leading veterinarians over the course of several years.

## Indications for Arthroscopy in the Canine Patient – Forelimb

### Shoulder

- Osteochondral lesions (osteochondritis dissecans (OCD), articular fractures of the border of the glenoid)
- Subluxation of the shoulder
- Evaluation of shoulder luxations
- Biceps tendon lesions: tendinitis, partial or complete tears
- Tenosynovitis of the biceps tendon
- Synovial biopsies
- Arthroliths
- Septic arthritis
- Lesions of the glenoid labrum
- Tendinitis, rupture of the subscapular muscle tendon

### Elbow

- Excision of the fragmented coronoid process
- OCD of the humeral condyle
- Treatment of arthritis of the elbow
- Treatment of the ununited anconeal process
- Removal of arthroliths
- Synovectomy
- Articular trauma
- Septic arthritis
- Elbow pain of unknown origin

### Carpus

- Cartilage fractures
- Carpal fractures without displacement (distal radius, radial carpal bone)
- Carpal trauma
- Synovial biopsies
- Septic arthritis

## Indications for Arthroscopy in the Canine Patient – Hind Limb

### Hip

- Labral lesions
- Arthroliths
- Articular trauma
- Preoperative evaluation of hip dysplasia before Triple Pelvic Osteotomy (TPO)
- Postoperative evaluation of TPO
- Postoperative evaluation of treatment of hip luxations and fractures
- Degenerative joint disease
- Septic arthritis
- Osteochondromatosis
- Tear of the ligament of the femoral head

### Stifle (Knee)


- Preoperative evaluation of torn cruciate ligaments
- The stable, osteoarthritic, non-inflamed stifle
- Lesions of the common digital extensor
- Lesions of the popliteal ligament
- Partial tear of the LCA
- Postoperative evaluation of cruciate ligament repair
- Differentiation between torn caudal cruciate ligament versus cranial cruciate ligament
- Cartilaginous lesions (fractures, OCD)
- Diagnosis of lesions associated with traumatic tears of the cruciate ligament
- Septic arthritis
- Synovial biopsies (polyarthritis, synovitis)
- Articular tumors

### Tarsus


- Excision of OCD fragments of the talus
- Tarsal arthritis
- Tarsocrural instability

# Slender Arthroscopes


for the arthroscopy of small joints, diameter 1.9 mm



64301 AA




64302 BN




64301 AA **HOPKINS® Straight Forward Telescope 0°**, diameter 1.9 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

64301 BA **HOPKINS® Forward-Oblique Telescope 30°**, diameter 1.9 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red


64302 BN **High-Flow Arthroscope Sheath**, diameter 2.8 mm, working length 8.5 cm




28305 BA



28306 BN



28306 S



**NEW** 28305 BA **HOPKINS® II Forward-Oblique Telescope 30°**, diameter 1.9 mm, length 6.5 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

**NEW** 28306 BN **High-Flow Arthroscope Sheath**, with traditional coupling mechanism, diameter 2.5 mm, working length 4 cm, one stopcock, **fixed**, for use with HOPKINS® II Telescope 30° 28305 BA and Obturators 28306 BS/BT color code: red

**NEW** 28306 S **High-Flow Arthroscope Sheath**, with snap-in coupling mechanism, diameter 2.5 mm, working length 4 cm, one stopcock, **fixed**, for use with HOPKINS® II Telescope 30° 28305 BA and Obturators 28306 BS/BT


# Slender Arthroscopes

for the arthroscopy of small joints, diameter 1.9 mm

For use with Arthroscopy Sheath 64302 BN



64302 BS



64302 BT

64302 BS     **Obturator, sharp**

64302 BT     **Obturator, blunt**

For use with Arthroscopy Sheath 28306 BN and 28306 S



28306 BT



28306 BS


**NEW** 28306 BT     **Obturator, blunt**

**NEW** 28306 BS     **Obturator, sharp**


# Slender Arthroscopes

for the arthroscopy of small joints, diameter 2.4 mm


Diameter 2.4 mm,  
sheath diameter 3.2 mm,  
working length 8.5 cm




64300 BA



64303 BN



28303 BS



30°

64300 BA

**HOPKINS® Wide Angle Forward-Oblique Telescope 30°**, diameter 2.4 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

64303 BN

**High Flow Arthroscope Sheath**, diameter 3.2 mm, working length 8.5 cm, color code: red

**NEW** 28303 BS

**High-Flow Arthroscope Sheath**, with snap-in coupling mechanism, diameter 3.2 mm, working length 8.5 cm, one stopcock, **fixed**, for use with HOPKINS® II Telescope 30° 28300 BA, 64300 BA and Obturators 28302 BU/BV

# Slender Arthroscopes

for the arthroscopy of small joints, diameter 2.4 mm

For use with Arthroscopy Sheath 64303 BN



64303 BU




64303 BV

64303 BU **Obturator, sharp**

64303 BV **Obturator, blunt**

For use with Arthroscopy Sheath 28303 BS



28302 BU



28302 BV


**NEW** 28302 BU **Obturator, sharp**

**NEW** 28302 BV **Obturator, blunt**


# Slender Arthroscope

with HOPKINS® II Telescope, diameter 2.7 mm


Sheath diameter 4 mm, working length 13 cm,  
with increased inflow and outflow rate



64019 BA



64128 AR



30°


64019 BA

**HOPKINS® II Forward-Oblique Telescope 30°**, diameter 2.7 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

64128 AR

**Arthroscope Sheath**, diameter 4 mm, working length 13 cm, with rotating stopcock, automatic lock, including LUER-Lock adaptor, with Stopcock 27502, color code: red


## For use with Arthroscope Sheath 64128 AR



64122 AS

64122 AS


**Obturator, sharp**



64122 AB

64122 AB

**Obturator, blunt**




64122 BS

**NEW**

64122 BS

**Obturator, sharp, short tip**



64122 BT

**NEW**

64122 BT


**Obturator, blunt, short tip**



# Arthroscope Sheath with Obturators

for HOPKINS® II Telescope, diameter 2.7 mm

Sheath diameter 4 mm, working length 13 cm,  
with increased inflow and outflow rate



28132 S


**NEW** 28132 S

**Arthroscope Sheath**, with snap-in coupling mechanism, diameter 4 mm, working length 13.5 cm, one stopcocks, **rotatable**, for use with HOPKINS® II telescopes 30°, 70° and Obturators 28132 BC/BT

**NEW** 28132 DS

**Same**, with two stopcocks


For use with Arthroscope Sheath 28132 S and 28132 DS



28132 BC

**NEW** 28132 BC

**Obturator**, blunt, for use with Arthroscope Sheath 28132 S/DS



28132 BT


**NEW** 28132 BT

**Obturator**, semisharp, for use with Arthroscope Sheaths 28132 CR/S/DS


# Slender Arthroscope

with HOPKINS® II Telescope, diameter 2.7 mm


Sheath diameter 4 mm, working length 10 cm,  
with increased inflow and outflow rate




67208 BA



64147 BN



28133 DS



30°

67208 BA

**HOPKINS® Forward-Oblique Telescope 30°**, diameter 2.7 mm, length 11 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red


64147 BN

**Arthroscope Sheath**, diameter 4 mm, working length 10 cm, for use with HOPKINS® Telescope 30°, 67208 BA, color code: red


**NEW** 28133 DS

**Arthroscope Sheath**, with snap-in coupling mechanism, diameter 4 mm, working length 7.5 cm, two stopcocks, **rotatable**, for use with HOPKINS® II telescope 30° and Obturator 28133 BC

## For use with Arthroscope Sheath 64147 BN



64147 BS



64147 BT


64147 BS

**Obturator**, sharp

64147 BT

**Obturator**, blunt

## For use with Arthroscope Sheath 28133 DS



28133 BC



**NEW** 28133 BC

**Obturator**, blunt

# Arthroscope Recommended for Canine Stifle

HOPKINS® II Telescope, diameter 4 mm

Sheath diameter 5.5 mm, working length 13 cm, with increased inflow and outflow rate



**NEW** 64230 BWA

**HOPKINS® II Wide Angle Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

**NEW** 64124 AR

**Arthroscope Sheath**, with one stopcock, rotating, for use with HOPKINS® II Telescopes 30°, diameter 5.5 mm, working length 13 cm color code: red

For use with Arthroscope Sheath 64124 AR



**NEW** 64127 BS

**Obturator**, sharp, color code: red

**NEW** 64127 BT

**Obturator**, blunt, color code: red

Short version, sheath diameter 5mm, working length 6.5 cm, with two stopcocks



**NEW** 64728 BWA

**HOPKINS® II Wide Angle Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 12 cm, **autoclavable**, fiber optic light transmission incorporated, for use with Arthroscope Sheath 28129 RK and Elbow Arthroscopy Set 28128 S, color code: red

**NEW** 64126 KR

**Arthroscope Sheath**, diameter 5.5 mm, working length 6.5 cm, with snap-in coupling mechanism, with one stopcock, **rotating**, for use with HOPKINS® II Telescope 30° and Obturator 64129 BT, diameter color code: blue




For use with High Flow Arthroscope Sheath 64126 KR

**NEW** 64129 BT

**Obturator**, semisharp, for Arthroscope Sheath 64129 RK, color code: red

# Forceps

Working length 11 cm



# ARTHROFORCE® III TAKE-APART® Instruments

dismantling

## ARTHROFORCE® III TAKE-APART®

### Punches, Scissors, Biopsy and Grasping Forceps – dismantling

ARTHROFORCE® III TAKE-APART® hand instruments feature a modular design. These instruments can be dismantled into a working attachment and a handle so that various working attachments can be exchanged.

The two-part design of the ARTHROFORCE® III TAKE-APART® instrument simplifies cleaning and sterilization. This reduces the risk of blood and tissue debris becoming trapped and prolongs instrument life. The specially designed cleaning adaptor enables cleaning during the main sterilization process.

The instruments are well balanced and precision-ground. The agile, flat yet robust cutting edge construction allows easy access in narrow areas.

ARTHROFORCE® III TAKE-APART® offers the following benefits:


- Instruments can be quickly and easily dismantled and reassembled
- Interchangeable working attachments provide perfect function and sharp cutting area
- Few handles necessary
- No hidden spaces that can trap blood and tissue debris
- Improved instrument cleaning
- Fully autoclavable
- Ergonomic design



# ARTHROFORCE® III TAKE-APART® Forceps

dismantling

Sheath diameter 3.5 mm, working length 13 cm



64168 P

**ARTHROFORCE® III TAKE-APART® Forceps**,  
straight, sheath diameter 3.5 mm, working length 13 cm  
including:  
**Handle**, without ratchet  
**Working Attachment**



64166 A

**BLAKESLEY ARTHROFORCE® III TAKE-APART® Forceps**,  
straight, sheath diameter 3.5 mm, working length 13 cm  
including:  
**Handle**, without ratchet  
**Working Attachment**



**NEW** 64168 B


**ARTHROFORCE® III TAKE-APART® Scissors**,  
hook-shaped, straight, sheath diameter 3.5 mm,  
working length 13 cm  
including:  
**Handle**, without ratchet  
**Working Attachment**

**Components/Spare Parts** see chapter 14

# ARTHROFORCE® III TAKE-APART® Grasping Forceps

dismantling

Sheath diameter 3.5 mm, working length 13 cm



64169 L

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, with hook, sheath diameter 3.5 mm, working length 13 cm including:

**Handle**, without ratchet

**Working Attachment**

64169 LS

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, with hook, sheath diameter 3.5 mm, working length 13 cm including:

**Handle**, with ratchet

**Working Attachment**



64169 G

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, for foreign bodies, without hook, sheath diameter 3.5 mm, working length 13 cm including:

**Handle**, without ratchet

**Working Attachment**

64169 GS

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, without hook, sheath diameter 3.5 mm, working length 13 cm including:

**Handle**, with ratchet

**Working Attachment**



64166 LL


**Cleaning Adaptor**, for **ARTHROFORCE® TAKE-APART®** instruments

**Components/Spare Parts** see chapter 14

# ARTHROFORCE® III TAKE-APART® Grasping Forceps

dismantling

Sheath diameter 2.3 mm, working length 10 cm



64183 L

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, with hook, sheath diameter 2.3 mm, working length 10 cm including:

**Handle**, without ratchet

**Working Attachment**




64183 LS

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, with hook, sheath diameter 2.3 mm, working length 10 cm including:

**Handle**, with ratchet

**Working Attachment**



64183 G

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, for arthroliths, without hook, sheath diameter 2.3 mm, working length 10 cm including:

**Handle**, without ratchet

**Working Attachment**



64183 GS

**ARTHROFORCE® III TAKE-APART® Grasping Forceps**, without hook, sheath diameter 2.3 mm, working length 10 cm including:

**Handle**, with ratchet

**Working Attachment**



64183 P

**ARTHROFORCE® III TAKE-APART® Forceps**, straight, sheath diameter 2.3 mm, working length 10 cm including:

**Handle**, without ratchet

**Working Insert**



64166 LL

**Cleaning Adaptor**, for **ARTHROFORCE® TAKE-APART®** instruments


Components/Spare Parts see chapter 14



# ARTHROFORCE® III TAKE-APART® Forceps

dismantling

Cutting width 1.4 mm, sheath diameter 2.8 mm, working length 13 cm



64167 A

**BLAKESLEY ARTHROFORCE® III TAKE-APART® Forceps**,  
straight, sheath diameter 2.8 mm, working length 13 cm

including:


**Handle**, without ratchet

**Working Attachment**



64166 LL

**Cleaning Adaptor**, for **ARTHROFORCE® TAKE-APART®**  
instruments



28171 HC



**NEW** 28171 HC

**SilGrasp® Grasping Forceps**, spoon-shaped jaws, jaw width 2.7 mm, straight jaws, sheath diameter 3 mm, straight, working length 13 cm



**NEW** 28171 AG

**SilGrasp® Alligator Grasping Forceps**, surgical, straight jaws, sheath diameter 3 mm, straight, working length 13 cm

**Diameter 2 mm**



64302 DH **Punch**, through-cutting,  
working length 10 cm



64302 EK **Scissors**, downbiting,  
working length 10 cm




64302 EO **Same**, upbiting



64302 L **Forceps with Double Spoon**,  
working length 10 cm



64302 U **Grasping Forceps**,  
working length 10 cm



64302 M/N



64302 M **Sickle Knife**, bayonet-shaped,  
diameter 1.5 mm, working length 7.5 cm




64302 N **Sickle Knife**, straight, diameter 1.5 mm,  
working length 7.5 cm



64302 S

64302 S **Palpation Hook**, graduated, diameter 1.5 mm,  
working length 7.5 cm, length of hook 1 mm



64302 W

64302 W **Changing Rod**, length 15 cm,  
for use with sheaths and trocars




64302 X

64302 X **Trocar**, diameter 2.5 mm, for use with  
Instruments 64302 L/DH/EK/U/M/N/S

64302 XS **Obturator**, sharp, for Trocar 64302 X

64302 XT **Same**, blunt



64145 S



64145 S **Hook and Retractor**, graduated, diameter 3.5 mm, length of hook 2 mm, working length 8 cm



28145 SN **Probe**, graduated, length of hook 2 mm, diameter 1 mm, working length 4 cm



64146 M **Knife**, curved, pointed, cuts both sides, width 3.5 mm



64146 N **Knife**, sickle-shaped



64146 F **Knife**, round, width 3.5 mm



64146 L **Knife**, hooked, width 3.5 mm



64728 CF **Micro Fracturing Chisel**, straight, working length 9 cm



64728 CG **Micro Fracturing Chisel**, angled, 30°, working length 9 cm



64728 CH **Micro Fracturing Chisel**, angled, 70°, working length 9 cm



64146 K **Curette**, spoon-shaped, round, straight, diameter 2.3 mm, for use with Trocar 64183 X



**NEW** 28147 R **Knife**, retractable, backward cutting, with exchangeable knife attachment, diameter 3.5 mm, working length 9 cm



64611



64611 **Curette**, size 0, length 21 cm



64620 **Same**, size 00




64630 **Same**, size 000



64640 **Same**, size 0000

## Operating Instruments for Arthroscopy



64906



64906

**Elevator**, sharp, curved,  
width 6 mm



64800



64800

**Elevator**, sharp, slight curved,  
width 4.5 mm



64100



64100

**Osteotome**, flat, graduated, curved,  
width 6 mm, length 18.5 cm




64104



64104

**Osteotome**, flat, cuts both sides,  
width 4 mm, length 19 cm



64527 S



64527 S

**Hand Burr**, with sleeve, diameter 2.7 mm,  
length 12 cm

including:

**Hand Burr**

**Hand Burr Sleeve**



64535 S

**Same**, diameter 3.5 mm

5-061

**Components/Spare Parts** see chapter 14

For straight operating instruments up to diameter 4.8 mm



64146 X/XS



64146 XT

64146 X      **Trocar**, diameter 5.5 mm

**For use with Trocar 64146 X**

64146 XS      **Obturator**, sharp

64146 XT      **Obturator**, blunt

For straight operating instruments up to diameter 3.5 mm



64169 X


64169 X      **Trocar**, diameter 4.5 mm

**For use with Trocar 64169 X**

64169 XS      **Obturator**, sharp

64169 XT      **Obturator**, blunt

For straight operating instruments up to diameter 2.8 mm



64183 X

64183 X      **Trocar**, diameter 3.5 mm


**For use with Trocar 64183 X**

64183 XS      **Obturator**, sharp

64183 XT      **Obturator**, blunt


## Cannulas, Changing Rods

### Operating Cannula



64146 VA/VT

64146 VA **Operating Cannula**, with stopcock, inner diameter 5.6 mm, outer diameter 6 mm, length 7.5 cm, for use with arthroscope sheaths, shavers, operating instruments




64146 VS

64146 VS **Obturator**, sharp

64146 VT **Obturator**, blunt

### Changing Rod Systems




64124 BX

64124 BX **Changing Rod**, length 23 cm, for use with arthroscope sheaths with inner diameter 4 mm

64124 BZ **Same**, for use with arthroscope sheaths with inner diameter 2.8 mm

### WISSINGER Rods, for large joints





64124 BY

64124 BY **Rod**, with pressure plate, length 35 cm, for use with arthroscope sheaths with inner diameter 4 mm

64124 BW **Same**, for use with arthroscope sheaths with inner diameter 2.8 mm

The changing rod can be used in arthroscopy to localize the anterior portal to the joint by means of the arthroscope sheath. The rod is introduced through the arthroscope sheath and advanced through the joint. A second arthroscope sheath can then be placed over the rod.

# LEIPZIG Stifle Distractor <sup>NEW</sup>



64820 LSD LEIPZIG **Stifle Distractor**  
including:  
**Distractor**  
**Fixation Screws**, package of 6  
**Lever Arm**  
**Screwdriver**

**Components/Spare Parts** see chapter 14



# Egress Cannulas




64146 R

64146 R **Egress Cannula**, LUER-Lock, diameter 2.2 mm, length 7 cm



64146 Q


64146 Q **Egress Cannula**, with Sharp Obturator 64146 QO, LUER-Lock, diameter 3.2 mm, length 7 cm



64146 T

64146 T **Egress Cannula**, with Sharp Obturator 64146 QO, LUER-Lock, with stopcock, diameter 3.2 mm, length 7 cm

28146 QB **Obturator**, blunt, for use with Egress Cannula 64146 Q



26 3100 38  
26 3101 38

26 3100 38 **Pressure Infusion Cuff**, 1.5 l,  
max. 400 mmHg, sterilizable

26 3101 38 **Same**, 3 l



20 3100 91

20 3100 91 **Manometer**, for use with Pressure Cuffs  
26 3100 38 and 26 3101 38



20 3100 90


20 3100 90 **Rubber Foot Pump**, with silicone tubing,  
length 200 cm, with pressure relief valve,  
fits Pressure Cuffs 26 3100 38 and 26 3101 38



20 3100 41

20 3100 41 **Silicone Tubing Set**, sterilizable,  
length 250 cm, for use with Pressure  
Cuffs 26 3100 38 and 26 3101 38 with  
Foot Pump 20 3100 90

# CYSTOSCOPY AND TRANSCERVICAL INSEMINATION





Uroendoscopy is increasingly performed by veterinary practitioners at all levels. Besides radiographic imaging, endoscopy is the most useful method of diagnosing urinary tract pathology. It is easy to perform and offers several therapeutic modalities as well as diagnostic, which are much less invasive than conventional cystotomy.

Rigid cystoscopes, available in various sizes, are introduced into the bitch transurethrally. In male dogs, a rigid cystoscope must be introduced into the bladder either transabdominally, or through a perineal urethrotomy. Fluid distension of the bladder and urethra are maintained by gravity flow of sterile saline or ringers lactate, from a bag that is suspended above the patient.

**Indications for Cystoscopy:**

- **Chronic cystitis**
- **Pollakiuria**
- **Hematuria**
- **Stranguria and tenesmus**
- **Stones in the bladder and urethra**

The endoscopist controls the degree of distension with the inflow and outflow valves on the cystoscope sheath.


The superior visualization facilitates procurement of biopsies, as well as the introduction of lasers, stone baskets and other instrumentation through the sheath. In cases of neoplasia, recurring urolithiasis and ectopic ureters, uroendoscopy constitutes the least traumatic therapeutic intervention available.

For removal of large stones, or to gain access to the bladder in male dogs, a rigid cystoscope can be introduced onto the bladder directly through the abdominal wall or with a laparoscope-assisted approach.

- **Neoplastic cells in the urinary sediment**
- **Abnormal radiological findings**
- **Incontinence**
- **Trauma**

# HOPKINS® II Telescopes and Cystoscope Sheaths

HOPKINS® II Telescopes for use with sheaths 9.5 – 14.5 Fr.



67065 C

64019 BA



64019 BA

**HOPKINS® II Forward-Oblique Telescope 30°**, diameter 2.7 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red



67065 C

**Sheath**, 14.5 Fr., working length 15 cm, with working channel 5 Fr. and two stopcocks, including blunt Obturator 67065 CO




67031 E

64017 B




64017 BA

**HOPKINS® II Forward-Oblique Telescope 30°**, diameter 1.9/2.1 mm, length 19 cm, **autoclavable** fiber optic light transmission incorporated, color code: red



67031 E

**Cystoscope-Urethroscope Sheath**, 10 Fr., working length 14 cm, with working channel 3 Fr., with Obturator 63031 EO and 2 LUER-Lock cones, color code: white-blue



67030 BA

67030 BA

**Operating Telescope 30°**, 9.5 Fr., working length 14 cm, **autoclavable**, fiber optic light transmission incorporated, with working channel 3 Fr. and 2 LUER-Lock cones color code: red

Components/Spare Parts see chapter 14

# Cysto-Urethroscope


8 Fr.

## Special Features

- Atraumatic distal sheath with rounded tip
- Minimal shaft diameter
- Maximum inflow due to large central instrument channel with 2 lateral irrigation ports
- The large working channel allows the use of rigid instruments and probes up to 4 Fr.
- 2 right angled lateral irrigation connections
- autoclavable

## Specifications

Instrument sheath:	8 Fr., one-step 8 – 12 Fr.
Central channel:	5 Fr. for use with instruments up to 4 Fr.
Telescope:	KARL STORZ Fiberoptic system, direction of view 6°
Length:	13 cm
Eyepiece:	angled, rigid



**NEW** 27030 KA

**Cysto-Urethroscope**, 8 Fr., 6°, 1 step, 8 – 11 Fr., length 13 cm, **autoclavable**, with angled eyepiece, fiber optic light transmission incorporated, with 2 lateral irrigation ports and 1 working channel, 5 Fr. for operating instruments 4 Fr.

## Following accessories are included:

27001 G	<b>Instrument Port with Sealing System and Quick Release Lock</b> , 1 channel
27550 N	<b>Seal</b> , for working channel, package of 10, single use recommended
27500	<b>LUER-Lock Tube Connector</b> , male, tube diameter 9 mm
27502	<b>LUER-Lock Tube Connector</b> , with stopcock, dismantling
27001 RA	<b>Cleaning Adaptor</b>
27001 E	<b>Insertion Aid</b> , for guide wires in ureterorenoscopy

## Optional accessory:

<b>NEW</b> 27001 GF	<b>Instrument Port with Sealing System and Quick Release Lock</b> , 2 channels
<b>NEW</b> 27001 GH	<b>Instrument Port with Sealing System and Quick Release Lock</b> , 2 channels, 1 straight channel, 1 lateral channel
<b>NEW</b> 39501 XK	<b>Wire Tray for Cleaning, Sterilization and Storage of pediatric ureterorenoscopes</b> and compact miniature hysteroscopes with a maximum working length of 27 cm, including cleaning adaptor for washer-disinfector, with lid, spare parts basket 39501 XS and silicone telescope holders, external dimensions (w x d x h): 460 x 150 x 80 mm, for use with instruments 270xx and 260x2 K with up to 27 cm working length

**Components/Spare Parts** see chapter 14

# Cysto-Urethroscope


9.5 Fr.

## Special Features

- Atraumatic distal sheath with rounded tip
- Minimal shaft diameter
- Maximum inflow due to large central instrument channel with 2 lateral irrigation connections
- The large working channel allows the use of rigid instruments and probes up to 5 Fr.
- 2 right angled lateral irrigation connections

## Specifications

- Instrument sheath:** 9.5 Fr., one-step 9.5 – 12 Fr.
- Central channel:** 6 Fr. for use with instruments up to 5 Fr.
- Telescope:** KARL STORZ Fiberoptic system, direction of view 6°
- Length:** 13 cm
- Eyepiece:** angled, rigid



**NEW** 27030 KB

**Cysto-Urethroscope**, 9.5 Fr., 6°, 1 step, 9.5 – 11 Fr., length 13 cm, **autoclavable**, with angled eyepiece, fiber optic light transmission incorporated, with 2 lateral irrigation ports and 1 working channel, 6 Fr. for operating instruments 5 Fr.

## Following accessories are included:

- 27001 G **Instrument Port with Sealing System and Quick Release Lock**, 1 channel
- 27550 N **Seal**, for working channel, package of 10, single use recommended
- 27500 **LUER-Lock Tube Connector**, male, tube diameter 9 mm
- 27502 **LUER-Lock Tube Connector**, with stopcock, dismantling
- 27001 RA **Cleaning Adaptor**
- 27001 E **Insertion Aid**, for guide wires in ureterorenoscopy

## Optional accessory:



- NEW** 27001 GF **Instrument Port with Sealing System and Quick Release Lock**, 2 channels
- NEW** 27001 GH **Instrument Port with Sealing System and Quick Release Lock**, 2 channels, 1 straight channel, 1 lateral channel
- NEW** 39501 XK **Wire Tray for Cleaning, Sterilization and Storage of pediatric ureterorenoscopes** and compact miniature hysteroscopes with a maximum working length of 27 cm, including cleaning adaptor for washer-disinfector, with lid, spare parts basket 39501 XS and silicone telescope holders, external dimensions (w x d x h): 460 x 150 x 80 mm, for use with instruments 270xx and 260x2 K with up to 27 cm working length

**Components/Spare Parts** see chapter 14



# HOPKINS® II Telescopes

Diameter 4 mm, for use with sheaths 17 – 25 Fr.  
and Miniature Straight Forward Telescope, diameter 1 mm



63005 AA


**HOPKINS® II Straight Forward Telescope 0°**,  
enlarged view, diameter 4 mm, length 30 cm,  
**autoclavable**, fiber optic light transmission  
incorporated,  
color code: green



63005 BA

**HOPKINS® II Forward-Oblique Telescope 30°**,  
enlarged view, diameter 4 mm, length 30 cm,  
**autoclavable**, fiber optic light transmission  
incorporated,  
color code: red

## For cystoscopy of the male cat



26120 JX


26120 JX

**Optical VERESS Pneumoperitoneum Needle**,  
with spring-loaded blunt inner cannula and  
insufflation stopcock, diameter 2.1 mm,  
length 17.5 cm, with telescope lock,  
for use with Miniature Telescope 62512


62512

**Miniature Straight Forward Telescope 0°**,  
semi-rigid, diameter 1 mm, length 20 cm,  
fiber optic light transmission incorporated,  
color code: green

# Cystoscope Sheaths, Deflecting Mechanism



63026 A



63163 V



63026 A

**Cystoscope-Urethroscope Sheath, 25 Fr.,** working length 22 cm, with Obturator 63026 AO and 2 LUER-Lock cones, color code: white

63026 B

**Same, 22 Fr.,** with Obturator 63026 BO, color code: blue

27026 C

**Same, 20 Fr.,** with Obturator 27026 CO, color code: red

27026 D


**Same, 19 Fr.,** with Obturator 27026 DO, color code: green

63026 U

**Same, 17 Fr.,** with Obturator 63026 UO, color code: yellow

63163 V


**Examination Sheath,** with stopcock, diameter 5 mm, working length 28 cm



63026 E

63026 E

**Catheter Deflecting Mechanism,** with ALBARRAN lever, with 1 lockable channel, with ratchet



63025 G

63025 G


**Telescope Bridge,** for use with Telescope 63005 BA and Sheaths 63026 with 1 lockable channel

Sheath with Obturator	Catheter Capacities Deflecting Mechanism 63026 E	Catheter Capacities Telescope Bridge 63025 G	Color Code
63026 A, 25 Fr.	10 Fr.	12 Fr.	white
63026 B, 22 Fr.	9 Fr.	10 Fr.	blue
27026 C, 20 Fr.	6 Fr.	7 Fr.	red
27026 D, 19 Fr.	5 Fr.	6 Fr.	green
63026 U, 17 Fr.	–	5 Fr.	yellow

**Components/Spare Parts** see chapter 14

# Extra Long Cystoscope-Urethroscopes

for urethro-cystoscopy and artificial insemination in bitches,  
working length 29 cm



63325 BA



63027 KL/UL



30°

63325 BA


**HOPKINS® Forward-Oblique Telescope 30°**, diameter 3.5 mm, working length 36 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

63027 KL

**Cystoscope-Urethroscope Sheath**, working length 29 cm, 22 Fr., accepts instruments and catheters up to 9 Fr., with Obturator 63027 KO and 2 LUER-Lock cones, color code: blue

63027 UL


**Cystoscope-Urethroscope Sheath**, working length 29 cm, 17 Fr., accepts instruments and catheters up to 6 Fr., with Obturator 63027 UO and 2 LUER-Lock cones, color code: yellow



63027 E

63027 E

**Catheter Deflecting Mechanism**, with ALBARRAN lever, with 1 lockable channel, with ratchet, only for use with 63027 KL and with instruments 9 Fr.



63027 NL

63027 NL

**Telescope Bridge**, with 1 lockable channel, for use with Telescope 63325 and Sheaths 63027

**Please note:**

Extra Long Cystoscope-Urethroscope Sheaths series 63027 require instruments 40 cm or longer.

**Components/Spare Parts** see chapter 14

# TCl Endoscope

for transcervical insemination and hysteroscopy in dogs,  
working length 43 cm



**NEW** 63002 LA

**Minitube TCl Endoscope**, 9.5 Fr., 6°, 1 step, conical, 9.5 – 13.5 Fr., length 43 cm, **autoclavable**, with angled eyepiece, fiber optic light transmission incorporated, 2 LUER-Lock adaptors and 1 working channel 6 Fr., for use with instruments up to 5 Fr., with Instrument Port 27001 GK, Seal 27550 N and Cleaning Adaptor 27001



**NEW** 63002 MK

**Minitube TCl Catheter**, 5 Fr., length 70 cm




**NEW** 27425 F

**Grasping Forceps**, for stone fragments, double action jaws, 5 Fr., length 60 cm  
color code: red



**NEW** 27425 P

**Grasping Forceps**, for larger stones and fragments, double action jaws, 5 Fr., length 60 cm,  
color code: red



**NEW** 27425 Z

**Biopsy Forceps**, double action jaws, 5 Fr., length 60 cm,  
color code: red

**Components/Spare Parts** see chapter 14

For use with Cystoscope Sheaths 67065 C, 63026 U, 27026 D, 63027 UL and Operating Cystoscope 27030 KB




67161 T

**Grasping Forceps**, flexible, double action jaws, 5 Fr., length 34 cm



67161 Z

**Biopsy Forceps**, flexible, oval, double action jaws, 5 Fr., length 34 cm



**NEW** 27095 P


**Grasping Forceps**, for larger stone fragments, double action jaws, 5 Fr., length 30 cm

For use with Cystoscope Sheaths 63026 U, 27026 D and 63027 UL



27034 FK

**Forceps for Foreign Body Removal**, flexible, double action jaws, 5 Fr., length 40 cm



27034 FL

**Biopsy Forceps**, flexible, double action jaws, 5 Fr., length 40 cm



67023 VK


**Stone Basket**, flexible, 5 Fr., length 60 cm including:  
**3-Ring Handle**  
**Basket**  
**Coil**

For use with Cystoscope Sheaths 63026 A – C, 63027 KL and 63027 NL



63175 A

**Grasping Forceps**, for removing foreign bodies, double-action jaws, 7 Fr., length 40 cm



63177 A

**Biopsy Forceps**, double action jaws, 7 Fr., length 40 cm, for use with Sheaths 63026 A/B/C/D



63182 A

**Injection Cannula**, 6 Fr., length 50 cm

For use with Cystoscope-Urethroscope Sheath 67031 E and Operating Cystoscopes 67030 BA and 27030 KA



61071 TJ

**Grasping Forceps**, for removing foreign bodies, double-action jaws, flexible, 3 Fr., length 28 cm



61071 ZJ



**Biopsy Forceps**, double-action jaws, flexible, 3 Fr., length 28 cm



27770 AA

**Ball Electrode**, unipolar, 3 Fr., length 53 cm

**Components/Spare Parts** see chapter 14



- 27770 AA **Ball Electrode**, unipolar, 3 Fr., length 53 cm
- 27770 A **Same**, 4 Fr.
- 27770 B **Same**, 5 Fr.
- 27770 C **Same**, 6 Fr.
- 27770 D **Same**, 7 Fr.
- 27770 E **Same**, 8 Fr.
- 27770 F **Same**, 10 Fr.




- 27772 AA **Needle Electrode**, unipolar, 3 Fr., length 53 cm
- 27772 A **Same**, 5 Fr.
- 27772 B **Same**, 7 Fr.



- 27778 A **Loop Electrode**, unipolar, 7 Fr., length 53 cm
- 27778 B **Same**, 10 Fr.


### Unipolar High Frequency Cords

KARL STORZ High Frequency  
Instrument Surgery Units

	26002 M	<b>Unipolar High Frequency Cord</b> , with 4 mm plug for models KARL STORZ, Erbe type T, older models and Ellman, length 300 cm
	26004 M	<b>Unipolar High Frequency Cord</b> , with 4 mm plug for models Berchtold and Martin, length 300 cm
	26005 M	<b>Unipolar High Frequency Cord</b> , with 5 mm plug for AUTOCON® II 400 SCB system (111, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ AUTOCON® system (50, 200, 350) and Erbe type ICC, length 300 cm
	26006 M	<b>Unipolar High Frequency Cord</b> , with 8 mm plug for models KARL STORZ AUTOCON® II 400 SCB system (112, 116) and Valleylab, length 300 cm

**Please note:**

All high frequency cables on this page are delivered with a length of 300 cm only. If a 500 cm cable is required, please add the letter **L**, e.g. 26002 ML to the order number.





- 63957      **Vaginoscope**, diameter 20 mm, working length 15 cm, with Obturator 63957 O
- 63959      **Vaginoscope**, diameter 15 mm, working length 20 cm, with Obturator 63959 O
- 63958 E    **Vaginoscope**, diameter 13 mm, working length 11 cm, with Obturator 63958 EO
- 63960      **Illumination Rod**, with connector pin for fiber optic light carrier
- 63338 T    **Vaginoscope Attachment**, with adjustable magnifier, magnification 6x
- 63956 J    **Sheath Insert**, for Telescope 63005 B, for use with Vaginoscopes 63957, 63959 and 63958 E
- 63591 T    **Biopsy Forceps**, with round cupped jaws, diameter 3.9 mm, working length 25 cm

# Instrument Set for Transurethral Viscous Fluid Injection

for the treatment of incontinence resulting from sterilization

## Special Features

- Safe positioning of the injection needle
- Precise and gentle control of the needle
- Working Element 63054 EJ ensures maximum stability of semi-rigid injection needles
- Rotatable LUER-Lock connection of the 4 Fr. needle enables easy handling by both right- and left-handed users
- Slender 20 Fr. operating sheath reduces the risk of urethral trauma to a minimum
- HOPKINS® II telescopes with 12° and 30° directions of view complete the instrument set, delivering high light transmission, outstanding detail resolution and high contrast



63005 BA

**HOPKINS® II Forward-Oblique Telescope 30°**, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red



63005 FA

**HOPKINS® II Telescope 12°**, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: black




63054 EJ

**Working Element**, with channel for Injection Needle 63054 NJ

63054 SJ

**Operating Sheath**, 20 Fr., with Obturator 27054 JO and 2 LUER-Lock cones, color code: red




63054 NJ

63054 NJ

**Injection Needle**, semi-rigid, 4 Fr., distal end 2 Fr., with rotatable LUER-Lock connection, for use with Working Element 63054 EJ (single use recommended)



# ENDOSCOPY IN BIRDS AND SMALL EXOTIC ANIMALS





**Due to the anatomical and physiological characteristics peculiar to avian and small exotic animal patients, endoscopy is an ideal method for examination.**

**Diagnosis of illnesses in the following organs and body parts:**

- Heart
- Lung
- Airsac
- Liver
- Proventriculus, ventriculus
- Intestines
- Pancreas
- Spleen
- Adrenal glands
- Kidneys
- Urethra
- Reproductive organs
- Esophagus
- Crop
- Cloaca
- Trachea
- Syrinx
- Fat bodies
- Swim bladder
- Gills
- Teeth
- Pharynx


**Endoscopic interventions:**

- Coelioscopy
- Biopsy
- Gender identification
- Sampling of debris
- Retrieval of foreign bodies
- Cloacoscopy
- Thoracoscopy
- Arthroscopy
- Minimally invasive surgery
- Laser and HF surgery
- Fine needle injection or aspiration
- Gastrointestinal endoscopy
- Cystoscopy
- Rhinoscopy
- Otoscopy
- Vaginoscopy
- Intubation aid
- Therapeutic curettage and lavage of lesions
- Salpingectomy, ovariectomy and orchidectomy
- Examination of oral cavities and trimming teeth

# HOPKINS® Telescopes

Diameter 1.9 and 2.7 mm

## Telescope, diameter 2.7 mm, long version




64019 BA



64019 BA

**HOPKINS® II Forward-Oblique Telescope 30°**, diameter 2.7 mm, length 18 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

## Telescopes, diameter 1.9 mm




64301 BA




64301 BA

**HOPKINS® Forward-Oblique Telescope 30°**, diameter 1.9 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red



64301 AA

**HOPKINS® Straight Forward Telescope 0°**, diameter 1.9 mm, length 10 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green



67030 BA




67030 BA

**Operating Telescope 30°**, 9.5 Fr., working length 14 cm, **autoclavable**, fiber optic light transmission incorporated, with working channel 3 Fr. and 2 LUER-Lock cones color code: red

Sheaths see page 137


## Examination Sheaths

For use with Telescope 64019 BA



64018 US

64018 US **Examination and Protection Sheath**, diameter 3.5 mm, working length 17 cm, for Telescope 64019 BA




67065 CV

64019 BA


67065 CV **Examination Sheath**, 14.5 Fr., working length 14 cm, with working channel 5 Fr. and Obturator 67065 CB, for use with 5 Fr. instruments

For use with Telescope 64301 BA



61029 C


61029 C **Protection and Examination Sheath**, working length 8.5 cm, with Obturator 61029 CO and 2 LUER-Lock cones, color code: blue



61029 D

61029 D **Operating Sheath**, 9 Fr., diameter 2.5 cm, working length 7 cm, with Obturator 61029 DO and working channel 3 Fr., color code: red

5-061



69117 Z



69117 Z

**Biopsy Forceps**, rigid, round cupped jaws, diameter 1.7 mm, working length 15 cm



67161 Z

**Biopsy Forceps**, flexible, oval, double action jaws, 5 Fr., length 34 cm



67161 T

**Grasping Forceps**, flexible, double action jaws, 5 Fr., length 34 cm



62501 EK

**Scissors**, flexible, single action jaws, 4 Fr., length 34 cm



67071 X

**Needle**, flexible, with outer tube, diameter 0.7 mm, length 31 cm



**NEW** 67161 M

**Measurement Instrument**, flexible, scale in mm, diameter 5 Fr., length 40 cm



67023 VK

**Stone Basket**, flexible, 5 Fr., length 60 cm

including:

**3-Ring Handle**

**Basket**

**Coil**



67159 L

**Polypectomy Snare**, 5 Fr., length 32 cm

including:

**Handle**

**Tube**

**Snare**

**Adaptor**



67772 A

**Needle Electrode**, unipolar, 5 Fr., length 53 cm



61071 TJ

**Grasping Forceps**, for removing foreign bodies, double-action jaws, flexible, 3 Fr., length 28 cm




61071 ZJ

**Biopsy Forceps**, double-action jaws, flexible, 3 Fr., length 28 cm

Components/Spare Parts see chapter 14

# Miniature Telescope and Optical VERESS Pneumoperitoneum Needle

For tracheoscopy, cystoscopy and rhinoscopy in smaller patients



- 26120 JX    **Optical VERESS Pneumoperitoneum Needle,**  
with insufflation stopcock, diameter 2.1 mm,  
length 17.5 cm, with telescope lock for use with  
Miniature Telescope 62512
- 62512    **Miniature Straight Forward Telescope 0°,**  
semi-rigid, diameter 1 mm, length 20 cm,  
fiber optic light transmission incorporated,  
color code: green






**VITOM® 25 SYSTEM**  
**VISUALIZATION SYSTEM FOR OPEN**  
**SURGERY WITH MINIMAL ACCESS**



# VITOM® 25 <sup>NEW</sup>

Visualization System for open Surgery with Minimal Access

**STORZ**  
KARL STORZ — ENDOSKOPE



3-12

# VITOM® 25 <sup>NEW</sup>

Visualization System for open Surgery with Minimal Access

VITOM® 25 from KARL STORZ represents a revolutionary and innovative way of displaying open surgery with minimal access in a high quality and ergonomic manner.

VITOM® 25 is based on the renowned HOPKINS® II rod lens system from KARL STORZ. With the help of a holding system, the VITOM® 25 is placed at a working distance of 25 – 60 cm above the surgical field. This gives the surgeon more room to work.

The small size of the VITOM® 25 reduces space requirements in the OR to a minimum. Due to its slim and compact design, the surgical field is not obstructed and even long instruments can be used with ease. The VITOM® 25 provides excellent depth of field, magnification, contrast and color reproduction, which are the requisites for FULL HD display and recording.

#### The VITOM® 25 system offers:

- Excellent FULL HD image quality
- Great depth of view
- Large working distance
- Ergonomic monitor work

The first-class enhanced imaging can be observed via a FULL HD monitor from a convenient distance by the surgeon, the assistant as well as the entire OR team.

The VITOM® 25 system has proven to be an excellent alternative to OR illumination cameras, loupes or operating microscopes in various surgical disciplines.

The VITOM® 25 is equipped with an integrated illumination, which can be supplemented with an optional illuminator.

The VITOM® 25 exoscope adds new capabilities to any existing KARL STORZ video system. A FULL HD endoscope imaging solution from KARL STORZ can be used with the VITOM® 25 exoscope.

- Compact design requiring minimal space in the OR
- Use of existing KARL STORZ FULL HD endoscopy system possible



**Exoscope and Illumination**



**Camera System and Display**



**Holding System**



**Documentation System**




KARL STORZ AIDA® compact NEO advanced

**Units and Transport Cart**


Equipment Cart 29005 DRB  
Transport Cart for POINT SETTER® 20 0200 83

**Cold Light Fontains XENON 300 SCB, Cameras and Monitors** see chapter EXCERPT FROM TELEPRESENCE  
**KARL STORZ AIDA®** see catalog KARL STORZ OR1 NEO™

**Exoscope and Illumination**  
Diameter 10 mm, length 11 cm



60 9160 20



60 9160 20

**VITOM® 25 HOPKINS® II Straight Forward Telescope 0°**, working distance 25 – 60 cm, diameter 10 mm, length 11 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

**Fiber Optic Light Cables 69495 NE, 495 NVB or 495 NVC recommended**



69495 NE



495 NVB



495 NVC

69495 NE


**Fiber Optic Light Cable**, diameter 4.8 mm, length 300 cm

495 NVB

**Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 4.8 mm, length 300 cm

495 NVC

**Fiber Optic Light Cable**, with 90° deflection to the instrument, very narrow radius of curvature, diameter 4.8 mm, length 300 cm



20 9180 20

20 9180 20

**VITOM® 25 Distance Rod**, length 25 cm

optional:



20 9170 00


20 9170 00 **VITOM® 25 Illuminator**, for straight forward telescopes 0°, 2 adjustable lenses with holding device for VITOM® 25, autoclavable

**Y-Fiber Optic Light Cable 495 UV recommended**



495 UV

495 UV **Y-Fiber Optic Light Cable**, 2x diameter 3.5 mm, length 230 cm, for simultaneous connection of two instruments



39501 A2

39501 A2 **Wire Tray for Cleaning, Sterilization and Storage** of two rigid endoscopes and one light guide cable, including holder for light post adaptors, silicone telescope holders and lid, external dimensions (w x d x h): 352 x 125 x 54 mm, for rigid endoscopes up to diameter 10 mm and working length 20 cm

**Cold Light Fontains XENON 300 SCB, Cameras and Monitors** see chapter EXCERPT FROM TELEPRESENCE

<b>Working distance:</b>	25 – 60 cm	
<b>Depth of view at working distance of: Depth of view:</b>	<b>25 cm</b> approx. 3.5 cm	<b>50 cm</b> approx. 7 cm
<b>Field of view at working distance of:</b> H3-Z camera zoom 1x H3-Z camera zoom 2x	25 cm 5 cm 3.5 cm	50 cm 10 cm 7 cm
<b>Reproduction scale at working distance of:</b>	25 cm	50 cm
<b>26" Monitor:</b> H3-Z camera zoom 1x H3-Z camera zoom 2x	approx. 8x approx. 16x	approx. 4x approx. 8x
<b>42" Monitor:</b> H3-Z camera zoom 1x H3-Z camera zoom 2x	approx. 14x approx. 28x	approx. 7x approx. 14x
<b>52" Monitor:</b> H3-Z camera zoom 1x H3-Z camera zoom 2x	approx. 17x approx. 34x	approx. 8x approx. 16x



POINT SETTER® Holding System



28172 WKS

- 28172 WKS **POINT SETTER®**, pneumatic holding arm including:  
**POINT SETTER® Holding Arm**, with OR table adaptor  
**KSLOCK Adaptor**, for clamping jaws  
**Clamping Jaw**, large  
**Clamping Jaw**, small  
**Clamping Jaw**, for fiberscopes  
**Pressure Regulator**, 7 bar  
**Drape**, for single use, sterile, package of 10
- 28272 CN **Clamping Cylinder**, folding, for flexible mounting of 10 mm telescopes on the telescope sheath, **autoclavable**. The clamping cylinder allows vertical movement and rotation of the telescope.
- 28172 WD **Drape**, for single use, sterile, package of 10

**Note:** A connecting tube is required to operate the POINT SETTER® arm. Please select an appropriate product on page 149





**Components/Spare Parts** see chapter 14



**POINT SETTER® Holding System**  
**Alternative and Spare Connecting Tubes**

	28172 WA	<b>Connecting Tube</b> , for POINT SETTER®, Dräger, length 600 cm
	28172 WB	<b>Connecting Tube</b> , for POINT SETTER®, Dräger air motor, length 600 cm
	28172 WC	<b>Connecting Tube</b> , for POINT SETTER®, compressor, length 600 cm
	28172 WN	<b>Connecting Tube</b> , for POINT SETTER®, Schrader, length 600 cm
	28172 WO	<b>Connecting Tube</b> , for POINT SETTER®, with open end, length 600 cm

**Accessories and Expendable Materials**

	20 0310 02	<b>Air Compressor</b> , for POINT SETTER®, 0 – 8 bar, power supply 230 VAC, 50/60 Hz including: <b>Fabric Reinforced Tube</b> , for CALCUSPLIT® SCB <b>Connecting Tube</b> , for POINT SETTER®, length 600 cm
	20 0310 02 C	<b>Air Compressor</b> , for POINT SETTER®, 0 – 8 bar, power supply 115 VAC, 50/60 Hz including: <b>Fabric Reinforced Tube</b> , for CALCUSPLIT® SCB <b>Connecting Tube</b> , for POINT SETTER®, length 600 cm
	2761291	<b>Pressure Regulator</b> , 7 bar, for POINT SETTER®, must be used for pressure > 8 bar on supply line
	28272 CN	<b>Clamping Cylinder</b> , folding, for flexible mounting of 10 mm telescopes on the telescope sheath, <b>autoclavable</b> . The clamping cylinder allows vertical movement and rotation of the telescope. For use with Clamping Jaw 28272 UG, 28272 UGN and POINT SETTER® universal adaptor 10 – 15 mm.
	28272 ECW	<b>Spring Balance</b> , loading capacity 25 N
	28172 WD	<b>Drape</b> , for single use, sterile, package of 10

**Components/Spare Parts** see chapter 14


Transport Cart



20 0200 83 **STOR-E® Transport Cart**, for POINT SETTER®  
holding arm  
including:  
**Subrack**, for POINT SETTER®  
**Accessory Set**, for POINT SETTER®  
**Top Cover**  
**Wire Basket**  
**Manual**

Components/Spare Parts see chapter 14


**ENDOCRANE® Holding System**



28272 EH

- 28272 EH **ENDOCRANE®**, piezoregulated holding arm, including stand including:  
**Socket**, to clamp to the OR table  
**Control Unit**  
**Cover\***, sterile, package of 25  
**Spring Balance**  
**Mains Cord**  
**Case**
- 28272 UGN **Clamping Jaw**, metal, with axial intake, for use with all square-headed KARL STORZ HOPKINS® telescopes, clamping range 16.5 up to 23 mm, with quick release coupling KSLOCK (male)
- 28272 CN **Clamping Cylinder**, folding, for flexible mounting of 10 mm telescopes on the telescope sheath, **autoclavable**. The clamping cylinder allows vertical movement and rotation of the telescope. For use with Clamping Jaw 28272 UG, 28272 UGN and POINT SETTER® universal adaptor 10 – 15 mm.

3-12

\*  mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

**Components/Spare Parts** see chapter 14




**Special Features:**

- Simple, fast and accurate positioning
- New, reinforced versions of various articulated stands
- Five joint functions enable a large number of different positions
- All five joint functions can be fixed by means of a mechanical central clamp
- Additional height adjustment by using the socket
- Additional angle adjustment by using the socket 28172 HR
- Socket for use with European and United States standard rails
- Ergonomical positioning at the operating table
- Eases the work routine of the assistant
- Instruments and telescopes are clamped securely
- Steady imaging of the operation field
- Modular composition using individual elements
- Maintenance-free solid construction
- Autoclavable




28272 RKB

Mechanical Holding System



28272 HC



28172 HR



28272 UGN











28272 CN

- 28272 HC **Articulated Stand**, L-shaped, long, reinforced version, only, especially large swivel range, with one mechanical central clamp for all five joint functions, height 48 cm, swivel range 66 cm, with quick release coupling KSLOCK (female)
- 28172 HR **Rotation Socket**, to clamp to the operating table, with one mounted Butterfly Nut 28172 HRS, for European and US standard rails, with lateral clamp for height and angle adjustment of the articulated stand
- 28272 UGN **Clamping Jaw**, metal, with axial intake, for use with all square-headed KARL STORZ HOPKINS® telescopes, clamping range 16.5 up to 23 mm, with quick release coupling KSLOCK (male)
- 28272 CN **Clamping Cylinder**, folding, for flexible mounting of 10 mm telescopes on the telescope sheath, **autoclavable**. The clamping cylinder allows vertical movement and rotation of the telescope. For use with Clamping Jaw 28272 UG, 28272 UGN and POINT SETTER® universal adaptor 10 – 15 mm.

# Mechanical Holding Systems

with KSLOCK



Socket	Articulated Stand					
28172 HK	28272 HA		28272 HB		28272 HC	
						
28172 HR						
	A 30 cm	B –	A 48 cm	B 15 cm	A 48 cm	B 15 cm
	C 20 cm	D 17 cm	C 20 cm	D 17 cm	C 27 cm	D 24 cm

Clamping Jaw	Complete instrument		
28272 UGN	28272 KGA	28272 KGB	28272 KGC
	28272 RGA	28272 RGB	28272 RGC



**Clamping Jaw**, metal, with axial intake, for use with all square-headed KARL STORZ telescopes, clamping range 16.5 up to 23 mm

28272 UKN	28272 KKA	28272 KKB	28272 KKC
	28272 RKA	28272 RKB	28272 RKC



**Clamping Jaw**, metal, with axial intake, for use with instrument, irrigation and telescope sheaths, clamping range 4.8 up to 12.5 mm

28272 UL	28272 KLA	28272 KLB	28272 KLC
	28272 RLA	28272 RLB	28272 RLC



**Clamping Jaw**, universal, clamping range 0 up to 18 mm

28272 UF	28272 KFA	28272 KFB	28272 KFC
	28272 RFA	28272 RFB	28272 RFC



**Clamping Jaw**, for use with all KARL STORZ polymer housing fiberscopes

## Clamping Cylinder



28272 CN **Clamping Cylinder**, folding, for flexible mounting of 10 mm telescopes on the telescope sheath, **autoclavable**. The clamping cylinder allows vertical movement and rotation of the telescope.

# Mobile OR-Rail <sup>NEW</sup>



28272 MR

28272 MR **Mobile OR-Rail**, for installation on tables in veterinary and laboratory settings or experimental surgery including:  
**Holding Plate**  
2x **C-Clamp**  
2x **Stop-Block**





# UNITS AND ACCESSORIES



**INSUFFLATORS**

**SUCTION AND IRRIGATION SYSTEMS**

**MOTOR SYSTEMS**

**HIGH FREQUENCY SURGERY UNITS**

**SCB**

calcusplit





- **INSUFFLATORS**
- **SUCTION AND IRRIGATION SYSTEMS**
- **MOTOR SYSTEMS**
- **HIGH FREQUENCY SURGERY UNITS**

INSUFFLATORS

SUCTION AND IRRIGATION  
SYSTEMS

MOTOR SYSTEMS

HIGH FREQUENCY  
SURGERY UNITS

The units manufactured by KARL STORZ combine long-lasting precision mechanics with state-of-the-art micro-electronic programmable controls. At KARL STORZ, the greatest emphasis is placed on user and patient safety. The quality assurance system of KARL STORZ is certified in accordance with the requirements of ISO 9001/EN 46001. It guarantees constant quality testing in the selection of materials and components. At the end of each manufacturing process, tests are carried out with automatic measuring and testing systems developed specially for this purpose. The results are recorded and logged. That gives each device a distinct “fingerprint” that can be checked at any time before and after it is delivered to the customer.

The standardized, modular design of KARL STORZ units was developed based on extensive ergonomic studies and is conceived for ease of care and cleaning and user-friendly practice, as well as to meet the demands of the special hygienic standards required in surgery. Clearly laid out adjacent function keys and displays guarantee efficient operation and make it easier to constantly monitor actual and set parameters. Acoustic and visual warning signals also assist the user. The settings can be changed manually at any time. Automatic microelectronic control systems guarantee optimum operating conditions and therefore relieve the surgeon in his work who can then fully concentrate on medical procedures.

The overall KARL STORZ product line includes the following categories of units with accessories:

- **Insufflators**
- **Suction and Irrigation Systems**
- **Motor Systems**
- **Units for High Frequency Surgery**





## ■ INSUFFLATORS

Electronic CO<sub>2</sub> ENDOFLATOR® SCB

*NEW* CO<sub>2</sub> ENDO-ARTHROFLATOR-VET

# Electronic CO<sub>2</sub> ENDOFLATOR® SCB

Recommended System Configuration

## Special Features:


- High degree of patient safety due to CF application part
- Easy to use
- Clear, adjacent displays for set value and actual value facilitate monitoring of the insufflation process
- Touch keys for precise preselection of set values
- Optical and acoustic alarm signals in the event of patient overpressure
- Fully automatic, electronically controlled gas refill (e. g. in case of gas loss when changing instruments)
- SECUVENT® Safety System: Constant monitoring of intraabdominal pressure; any overpressure is reduced immediately
- SCB model with connections to the KARL STORZ Communication Bus (KARL STORZ-SCB)



**62 4305 08-1 Electronic CO<sub>2</sub> ENDOFLATOR® SCB**,  
power supply 100 – 240 VAC, 50/60 Hz  
including:  
**Mains Cord**  
**Silicone Tubing Set**, sterilizable  
**Universal Wrench**  
**SCB Connecting Cable**, length 100 cm  
**CO<sub>2</sub>/N<sub>2</sub>O Gas Filter\***, sterile, for single use,  
package of 10

## Specifications:

Operating Mode	Operating Mode	SECUVENT® Safety System	●
Gas Flow	0-20 l/min	Power Supply	100-240 VAC, 50/60 Hz
Pressure	0-30 mmHg (3990 Pa)	Dimensions w x h x d	305 x 164 x 233 mm
Gas	CO <sub>2</sub>	Weight	6 kg
Measuring/Control System	electronic	Certified to	IEC 601-1, CE acc. to MDD
Parameter Display	- Insufflation Pressure - Intraabdominal Pressure: 0-50 mmHg (6650 Pa) - Gas Flow (l/min): 0-20 - Gas Consumption (l): 0-999		

\*  mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

**Optional Accessories for the Electronic CO<sub>2</sub> ENDOFLATOR® SCB** see page U 10  
**Components/Spare Parts** see chapter 14

# Electronic CO<sub>2</sub> ENDOFLATOR® SCB

System Components

CO<sub>2</sub> Bottle, Pin-Index connector




26 4000 91 empty  
26 4000 93 filled

CO<sub>2</sub> Bottle, German connector



26 4000 90 empty  
26 4000 92 filled

High Pressure Tube, Pin-Index connector



20 4000 22 length 55 cm  
20 4000 28 length 102 cm

High Pressure Tube, German connector



20 4000 21 length 55 cm  
20 4000 27 length 102 cm

High Pressure Tube, ISO connector



20 4002 22 length 102 cm

Low Pressure Tube, for central CO<sub>2</sub> gas supply



20 4000 26 length 300 cm

UNIT SIDE

PATIENT SIDE




SCB®

CO<sub>2</sub>/N<sub>2</sub>O Gas Filter, sterile, for single use




031122-10\*

Silicone Tube, for insufflation



20 4001 43

VERESS Needle and Trocars



62120 J

62160 MCL

62103 APL

10-002



\* mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# CO<sub>2</sub> ENDO-ARTHROFLATOR-VET <sup>NEW</sup>

Recommended System Configuration

### Special Features:

- One unit with laparoscopy mode and arthroscopy mode
- High operating safety
- Easy to use
- Continuously adjustable insufflation pressure
- Visual and acoustic alarm signals in the event of patient overpressure
- Fully automatic, electronically controlled gas refill (e.g. in case of gas loss when changing instruments)



62 4325 01 **CO<sub>2</sub> ENDO-ARTHROFLATOR-VET**,  
power supply 100 – 240 VAC, 50/60 Hz  
including:

**Mains Cord**

**Mains Cord, US-version**


**Silicone Tubing Set**, sterilizable

**Universal Wrench**

**CO<sub>2</sub>/N<sub>2</sub>O Gas Filter\***, sterile, for single use,  
package of 10

### Specifications:

Operating Mode	- ARTHRO mode - LAPARO mode	Parameter Display	- Insufflation Pressure - Intraabdominal Pressure: 0-100 mmHg (13.3kPa) - Gas Flow (l/min): 0-20 - Gas Consumption (l): 0-999
Gas Flow	- 0-20 l/min - 0-20 l/min	Power Supply	100-240 VAC, 50/60 Hz
Pressure	- 50-100 mmHg (6650 Pa-13332 Pa) - 0-15 mmHg (1999 Pa)	Dimensions, w x h x d	305 x 112 x 234 mm
Gas	CO <sub>2</sub>	Weight	5.5 kg
Measuring/Control System	electronic	Certified to	IEC 601-1, CE acc. to MDD

\*  mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

**Optional Accessories for the CO<sub>2</sub> ENDO-ARTHROFLATOR-VET** see page U 10  
**Components/Spare Parts** see chapter 14



# CO<sub>2</sub> ENDO-ARTHROFLATOR-VET <sup>NEW</sup>

System Components

CO<sub>2</sub> Bottle, Pin-Index connector




26 4000 91 empty  
26 4000 93 filled

CO<sub>2</sub> Bottle, German connector



26 4000 90 empty  
26 4000 92 filled

High Pressure Tube, Pin-Index connector



20 4000 22 length 55 cm  
20 4000 28 length 102 cm

High Pressure Tube, German connector



20 4000 21 length 55 cm  
20 4000 27 length 102 cm

High Pressure Tube, ISO connector



20 4002 22 length 102 cm

Low Pressure Tube, for central CO<sub>2</sub> gas supply



20 4000 26 length 300 cm

UNIT SIDE

PATIENT SIDE




CO<sub>2</sub>/N<sub>2</sub>O Gas Filter, sterile, for single use



031122-10\*

Silicone Tube, for insufflation



20 4001 43

VERESS Needle, Trocar and Arthroscopy Sheath



62120 J

62160 MCL

64131 CR



mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# Optional Accessories

for Insufflators


INSUFFLATORS

for use with

Electronic CO<sub>2</sub>  
ENDOFLATOR® SCB  
26 4305 20-1

CO<sub>2</sub> ENDO-ARTHRO-  
FLATOR-VET  
62 4325 20

	62120 J 62120 JK 62120 JLL 62120 JL	VERESS <b>Pneumoperitoneum Needle</b> , with spring-loaded blunt stylet, LUER-Lock, <b>autoclavable</b> , diameter 2.1 mm, length 10 cm <b>Same</b> , length 7 cm <b>Same</b> , length 15 cm <b>Same</b> , length 13 cm	•	•
	30103 HC	<b>Trocar</b> , size 11 mm	•	•
	64131 CR	<b>High-Flow Arthroscope Sheath</b> , diameter 5.5 mm, working length 13.5 cm	-	•
	20 4000 32	<b>High Pressure Inline Gas Filter</b>	•	•
	20 4001 43	<b>Silicone Tubing Set</b> , sterilizable	•	•
	26 4000 90 26 4000 92	<b>CO<sub>2</sub> Bottle</b> , empty, with German connection <b>Same</b> , filled	•	•
	26 4000 91 26 4000 93	<b>CO<sub>2</sub> Bottle</b> , empty, with Pin-Index connection <b>Same</b> , filled	•	•
	20 4000 21 20 4000 27	<b>CO<sub>2</sub> High Pressure Tube</b> , American connection/German connection, length 55 cm <b>Same</b> , length 102 cm	•	•
	20 4000 22 20 4000 28	<b>CO<sub>2</sub> High Pressure Tube</b> , American connection/ Pin-Index connection, length 55 cm <b>Same</b> , length 102 cm	•	•
	20 4002 22	<b>CO<sub>2</sub> High Pressure Tube</b> , American connection/ISO connection, length 102 cm	•	•
	* 031122-10*	<b>CO<sub>2</sub>/N<sub>2</sub>O Gas Filter</b> , sterile, for single use, package of 10	•	•

\*  mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany



## ■ SUCTION AND IRRIGATION SYSTEMS

ARTHROPUMP® SCB

VETPUMP®

VETPUMP® II

# ARTHROPUMP® SCB

Irrigation System,  
Recommended Complete System Configuration

## Special features

- Irrigation pressure and flow rate are electronically controlled
- Facing displays of operating parameter setpoints and actual performance data allow ready operator monitoring of unit operational status
- Preselectable values for irrigation pressure and flow are retained in memory while the unit is switched off and will be redisplayed when it is switched on again
- State-of-the-art technology for easy operation and high patient safety
- ARTHROPUMP® with KARL STORZ Communication Bus (KARL STORZ-SCB)



**69 3300 04-1 ARTHROPUMP® SCB**,  
power supply 100 – 240 VAC, 50/60 Hz  
including:  
**Mains Cord**  
**Tubing Set\***, for single use  
**SCB Connecting Cable**, length 100 cm

### Optional Accessories:

**20 0141 30 One-Pedal Footswitch**, digital, one-stage, for irrigation  
**26 3300 42 ARTHRO Tubing Set**, for ARTHROPUMP® SCB,  
reusable, sterilizable

### Specifications:

Irrigation	- Pressure Regulated: 0-200 mmHg (26.6 kPa) - Flow Regulated: 0-1000 ml/min	Dimensions w x h x d	305 x 110 x 260 mm
Power Supply	100-240 VAC, 50/60 Hz	Weight	7 kg
		Certified to	IEC 601-1, CE acc. MDD






mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

**Optional Accessories for ARTHROPUMP® SCB** see page U 18  
**Components/Spare Parts** see chapter 14

# ARTHROPUMP® SCB

System Components



SUCTION AND IRRIGATION SYSTEMS

# VETPUMP®

Insufflation Unit

Provides air insufflation for gastrointestinal endoscopy,  
when using a light source without integrated insufflation pump.



**69 3215 01** **VETPUMP®**, power supply 230 VAC,  
50/60 Hz  
including:  
**Mains Cord**

#### Specifications:

Flow non-regulated	max. 5 l/min	Dimensions w x h x d	305 x 110 x 260 mm
Pressure non-regulated	max. 800 mmHg (106.4 kPa)	Weight	2.5 kg
Power Supply	- <b>69 3215 20</b> 230 VAC, 50/60 Hz	Certified to	IEC 601-1, CE acc. to MDD

**Optional Accessories for the VETPUMP®** see page U 18  
**Components/Spare Parts** see chapter 14

UNIT SIDE

PATIENT SIDE



Silicone Tubing Set



610 VE

Adaptor, for use with light sources XENON NOVA® and XENON 175 SCB



69487 KVE



# VETPUMP® II

Combined, adjustable Suction and Irrigation System,  
Recommended Complete System Configuration

## Special Features:

- Powerful suction and irrigation pump for all veterinary medical applications for both rigid and flexible KARL STORZ endoscopes
- Irrigation pressure and suction vacuum can be regulated separately
- Low noise level as all pumps operate on demand



**69 3216 01 VETPUMP® II**, power supply 100/120, 230/240 VAC, 50/60 Hz including:

**Mains Cord**

**Mains Cord, US-version**

**Suction Bottle, 1.5 l, sterilizable**

**Bottle Cap**, for use with suction bottle 1.5 l, sterilizable

**Irrigation Bottle, 1 l, sterilizable**

**Bottle Cap**, for use with irrigation bottle

**Bottle Holder**, for suction bottle

**Bottle Holder**, for irrigation bottle

**Filter\***, to block fluid during suction, for single use, unsterile, package of 10

## Optional Accessories:

**69 3200 33 Silicone Tubing Set**, irrigation-suction from unit to bottles, autoclavable

**69811 Instrument Tubing Set, autoclavable**, connects bottles to Handle 69810, for use with VETPUMP® II

**69810 Handle**, for irrigation and suction, **autoclavable**, for use with Tubing Set 69811, in conjunction with 5 mm suction and irrigation tubes and catheters any size with LUER-Lock connection

**20 3000 42 Silicone Tubing Set**, for suction, sterilizable

**610 VE Silicone Tubing Set**, for insufflation when using an external airpump, in combination with Adaptor 487 VE or 69487 KVE

## Specifications:

Pressure Regulated	max. 600 mmHg	Dimensions w x h x d	305 x 110 x 260 mm
Flow Regulated	max. 3.5 l/min	Weight	3.8 kg
Suction Pressure	regulated: max. 580 mmHg	Certified to	EN 61010-1
Power Supply	100/120, 230/240 VAC, 50/60 Hz		




mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

**Optional Accessories for the VETPUMP® II** see page U 18, **Additional suction and irrigation accessories** see chapter 6, HANDLES AND INSTRUMENTS FOR SUCTION AND IRRIGATION

**Components/Spare Parts** see chapter 14





5-061











mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# Optional Accessories

for Pump/Suction or Irrigation System

SUCTION AND IRRIGATION SYSTEMS

			for use with		
			ARTHROPUMP® 69 3300 04-1	VETPUMP® 69 3215 01 69 3215 01 C 69 3215 01 U	VETPUMP® II 69 3216 01
	20 0141 30	<b>One-Pedal Footswitch</b> , digital, one-stage, for irrigation	●	-	-
	20 3000 42	<b>Silicone Tubing Set</b> , for suction, sterilizable	-	-	●
	20 3000 50 20 3000 32 20 3000 33	<b>Suction Bottle</b> , 5 l, sterilizable <b>Bottle Stand</b> , for suction bottle 5 l <b>Bottle Stand Holder</b> , for Bottle Stand 20 3000 32	-	-	●
	26173 BN	<b>Suction and Irrigation Tube</b> , with lateral holes, with two-way stopcock for single hand control, size 5 mm, length 36 cm	-	-	●
	69811	<b>Instrument Tubing Set, autoclavable</b> , connects bottles to Handle 69810, for use with VETPUMP® II	-	-	●
	69 3200 33	<b>Silicone Tubing Set</b> , irrigation-suction from unit to bottles, <b>autoclavable</b>	-	-	●
	610 VE	<b>Silicone Tubing Set</b> , for insufflation when using an external airpump, in combination with Adaptor 487 VE or 69487 KVE	-	●	●
	69 0106 20	<b>VET Switch</b> , for VETPUMP® II, consisting of shut-off device and foot pedal	-	-	●



## ■ MOTOR SYSTEMS

**NEW** UNIDRIVE® S III ARTHRO SCB

## UNIDRIVE® S III ARTHRO SCB

The aim of all arthroscopic procedures is the same: to complete the procedure quickly and with minimal trauma, thereby reducing patient time under anesthesia.

With UNIDRIVE® S III ARTHRO, KARL STORZ has developed a motor system that is customized to meet these requirements.

- Easy to use
- Automatic handpiece recognition
- Control of rotation speed similar to the multifunctional handpiece
- Up to 15,000 rpm
- Comprehensive portfolio of single-use and reusable blades
- Communication with ARTHROPUMP® Power



# UNIDRIVE® S III ARTHRO SCB <sup>NEW</sup>

Recommended System Configuration

**Special Features:**

- Arthroscopic motor system for all joints
- Stable torque throughout entire speed range
- Motors with speed ranging from 6000 rpm to max. 15,000 rpm
- Multifunctional handpiece with accessories:
  - Jacobs chucks
  - Pin drivers
  - Sagittal saws
  - Synthes-style drivers
- Various control options possible
  - Hand controls
  - Footswitch
  - Touch Screen
- Fast and easy change of blade via quick coupling
- Single-use and reusable blades available



**28 7230 01-1 UNIDRIVE® S III ARTHRO SCB**, with color display, touch screen operation, two motor outputs, integrated SCB module, power supply 100 – 120/230 – 240 VAC, 50/60 Hz including:  
**SCB Connecting Cable**, length 100 cm  
**Mains Cord**  
**Manual**

**Specifications:**

Operating Mode	- Counter Clockwise - Clockwise - Oscillating	User-Defined Memory	Programmable Physicians' Profiles
Connections	- Small Joint Handpiece	Power Supply	100-127 / 220-250 VAC, 50/60 Hz
	- Hybrid Handpiece	Dimensions w x h x d	305 x 155 x 318 mm
	- Standard Handpiece	Weight	9.8 kg
	- High-Speed Handpiece	Certified to	IEC 601-1, CE acc. MDD
	- Multifunction Handpiece		
	- Footswitch		
	- KARL STORZ-SCB		


**System requirements:**

SCB control NEO System with installed SCB control NEO Software Release **20 0900 01-44** or higher

**Components/Spare Parts** see chapter 14

# UNIDRIVE® S III ARTHRO SCB

Handpieces



28 7210 40

28 7210 40 **Small Joint Shaver Handpiece**, up to 6000 rpm, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB



28 7210 50

**NEW** 28 7210 50 **Hybrid Shaver Handpiece**, up to 7000 rpm, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB



28 7210 35

28 7210 35 **Standard Shaver Handpiece**, up to 8000 rpm, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB



28 7210 45

28 7210 45 **Standard Shaver Handpiece**, up to 8000 rpm, **without handpiece buttons**, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB



28 7210 30

28 7210 30 **High-Speed Shaver Handpiece**, up to 15000 rpm, for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL SCB



28 7210 36

28 7210 36 **Multifunction Handpiece**, for use with UNIDRIVE® S II ARTHRO SCB and POWERSHAVER SL SCB

Optional accessories for UNIDRIVE® S III ARTHRO SCB see pages U 23 and 25

# UNIDRIVE® S III ARTHRO SCB

Handpieces

**Specifications:**

Handpiece	KARL STORZ Quick Coupling	for use with Reusable Blades	Rotation (rpm)	Oscillation (rpm)	Connecting Cable	Weight
28 7210 40 <b>Small Joint Shaver Handpiece</b>	●	●	1000 to 6000	500 to 3000	3.5 m	0.27 kg
28 7210 50 <b>Hybrid Shaver Handpiece</b>	●	●	1000 to 7000	500 to 3000	3.5 m	0.43 kg
28 7210 35 <b>Standard Shaver Handpiece</b>	●	●	1000 to 8000	500 to 3000	3.5 m	0.47 kg
28 7210 45 <b>Standard Shaver Handpiece</b>	●	●	1000 to 8000	500 to 3000	3.5 m	0.44 kg
28 7210 30 <b>High-Speed Shaver Handpiece</b>	●	●	1000 to 15000	500 to 3000	3.5 m	0.58 kg
28 7210 36 <b>Multifunction Handpiece</b>	●	–	0 to 1700	–	3.5 m	0.65 kg

**Optional accessories for UNIDRIVE® S III ARTHRO SCB shaver handpieces and blades**



- 28205 SA     **Irrigation Adaptor**, for validated cleaning of shaver handpieces
- 28205 SAA     **Irrigation Adaptor**, for validated cleaning of outer blades
- 28205 SAI     **Irrigation Adaptor**, for validated cleaning of inner blades

**Optional accessories for Multifunctional Handpiece 28 7210 36:**






- 28 7210 31     **Jacobs Chuck**, up to diameter 4.5 mm
- 28 7210 37     **Jacobs Chuck**, up to diameter 6.5 mm
- 28 7210 38     **Jacobs Chuck**, keyless, up to diameter 7.4 mm, for use with Multifunction Handpiece 28 7210 36
- 28 7210 39     **Sagittal Saw Attachment**, keyless, with tilt lever, for use with Multifunction Handpiece 28 7210 36 and Saw Blades 28207 only
- 28 7210 33     **Pin Driver**, for drilling wires diameter 1.5 – 3 mm
- 28 7210 34     **Pin Driver**, for Synthes attachments
- 28 7210 71     **Connecting Cable**, for FMS pump

# UNIDRIVE® S III ARTHRO SCB

System Components



MOTOR SYSTEMS





# Optional Accessories

## for Multifunction Handpiece 28721036



**28 7210 31**    **Jacobs Chuck**, up to diameter 4.5 mm



**28 7210 37**    **Jacobs Chuck**, up to diameter 6.5 mm



**28 7210 33**    **Pin Driver**, for drilling wires diameter 1.5 – 3 mm




**28 7210 38**    **Jacobs Chuck**, keyless, up to diameter 7.4 mm, for use with Multifunction Handpiece **28 7210 36**



**28 7210 39**    **Sagittal Saw Attachment**, keyless, with tilt lever, for use with Multifunction Handpiece **28 7210 36** and Saw Blades 28207 only



**28 7210 34**    **Pin Driver**, for Synthes attachments



**28 7210 71**    **Connecting Cable**, for FMS pump

## for UNIDRIVE® S III ARTHRO SCB shaver handpieces and blades



**28205 SA**    **Irrigation Adaptor**, for validated cleaning of shaver handpieces



**28205 SAA**    **Irrigation Adaptor**, for validated cleaning of outer blades



**28205 SAI**    **Irrigation Adaptor**, for validated cleaning of inner blades

# Shaver Blades

for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL,  
for cutting and grinding

Working length 180 mm

Type	Use	Order No.	Diameter	Color Code	Application
 Aggressive Cutter	for single use	28208 BKS	4.5 mm mauve	green	soft tissue
 Full Radius Resector	for single use	28208 CCS	4.5 mm mauve	yellow	soft tissue
 Aggressive Full Radius Resector	for single use	28208 DCS	4.5 mm mauve	blue	soft tissue
 Curved Full Radius Resector	for single use	28208 EGS	4.5 mm mauve	light blue	soft tissue
 Curved Aggressive Full Radius Resector	for single use	28208 EHS	4.5 mm mauve	light green	soft tissue
 Semi Hooded Barrel Burr	for single use	28208 IDS	5.5 mm dark gray	pink	bone

# Shaver Blades

for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL,  
for cutting and grinding

Working length 120 mm

Type	Use	Order No.	Diameter	Color Code	Application
	reusable	28205 AB	3.5 mm	green	soft tissue
		28205 AC	4.2 mm		
		28205 AD	5.5 mm		
	for single use	28205 ABS	3.5 mm light gray		
		28205 ACS	4.2 mm medium gray		
		28205 ADS	5.5 mm dark gray		
	reusable	28205 AK	4.5 mm	green	soft tissue
	for single use	28205 AKS	4.5 mm mauve		
		*28208 AKS	4.5 mm dark gray		
	reusable	28205 BB	3.5 mm	red	soft tissue
		28205 BC	4.2 mm		
		28205 BD	5.5 mm		
	for single use	28205 BCS	4.2 mm medium gray		
		28205 BDS	5.5 mm dark gray		
		28205 CB	3.5 mm		
reusable	28205 CC	4.2 mm			
	28205 CD	5.5 mm			
	for single use	28205 CBS	3.5 mm light gray		
28205 CCS		4.2 mm medium gray			
28205 CDS		5.5 mm dark gray			
	reusable	28205 CK	4.5 mm	yellow	soft tissue
	for single use	28205 CKS	4.5 mm mauve		
		*28208 CKS	4.5 mm dark gray		

\*28208 blades are available with distal drill guide

# Shaver Blades

for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL,  
for cutting and grinding

Working length 120 mm

Type	Use	Order No.	Diameter	Color Code	Application
	reusable	28205 DB	3.5 mm	blue	soft tissue
		28205 DC	4.2 mm		
		28205 DD	5.5 mm		
	for single use	28205 DBS	3.5 mm light gray		
		28205 DCS	4.2 mm medium gray		
		28205 DDS	5.5 mm dark gray		
	reusable	28205 DK	4.5 mm	blue	soft tissue
	for single use	28205 DKS	4.5 mm mauve		
	reusable	28205 EG	4.2 mm	light blue	soft tissue
	for single use	28205 EGS	4.2 mm medium gray		
	reusable	28205 EH	4.2 mm	light green	soft tissue
	for single use	28205 EHS	4.2 mm medium gray		
	reusable	28205 FC	4.2 mm	orange	bone
		28205 FD	5.5 mm		
	for single use	28205 FCS	4.2 mm medium gray		
		28205 FDS	5.5 mm dark gray		
		28205 FES	6.5 mm black		

# Shaver Blades

for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL,  
for cutting and grinding

Working length 120 mm

Type	Use	Order No.	Diameter	Color Code	Application
	reusable	28205 GD	4.2 mm	purple	bone
		28205 GE	5.5 mm		
		28205 GF	6.5 mm		
	for single use	28205 GDS	4.2 mm medium gray		
		28205 GES	5.5 mm dark gray		
		28205 GFS	6.5 mm black		
	reusable	28205 HC	4.2 mm	pink	bone
		28205 HD	5.5 mm		
		28205 HE	6.5 mm		
	for single use	28205 HCS	4.2 mm medium gray		
		28205 HDS	5.5 mm dark gray		
		28205 HES	6.5 mm black		
	for single use	28205 KDS	5.5 mm medium gray	light purple	bone
	reusable	28205 MK	4.5 mm	pink	bone, soft tissue
		28205 MD	5.5 mm		
	for single use	28205 MKS	4.5 mm mauve		
		28205 MDS	5.5 mm dark gray		
	reusable	28205 NK	4.5 mm	light blue	bone, soft tissue
		28205 ND	5.5 mm		
	for single use	28205 NKS	4.5 mm mauve		
		28205 NDS	5.5 mm dark gray		

# Shaver Blades

for use with UNIDRIVE® S III ARTHRO SCB and POWERSHAVER SL,  
for cutting and grinding for small joints

Working length 70 mm

Type	Use	Order No.	Diameter	Color Code	Application
	reusable	28206 AA	2.7 mm	green	soft tissue
	for single use	28206 AAS	2.7 mm white		
		28206 ABS	3.5 mm light gray		
	reusable	28206 CA	2.7 mm	yellow	soft tissue
	for single use	28206 CAS	2.7 mm white		
		28206 CBS	3.5 mm light gray		
	reusable	28206 DA	2.7 mm	blue	soft tissue
		28206 DB	3.5 mm		
	for single use	28206 DAS	2.7 mm white		
		28206 DBS	3.5 mm light gray		
	for single use	28206 FAS	2.7 mm white	orange	bone
		28206 FBS	3.5 mm light gray		

# Accessories

Saw Blades for use with Multifunction Handpiece 28 7210 36 of UNIDRIVE® S III ARTHRO SCB and Saw Handpiece 28 7200 38, for use with Power Unit S1 and S2



28207 SLA

28207 SLA **Saw Blade**, blade width 10 mm, blade thickness 0.4 mm, cutting thickness 0.6 mm, working length 50 mm



28207 SLB

28207 SLB **Same**, blade width 14 mm



28207 SLC

28207 SLC **Saw Blade**, blade width 10 mm, blade thickness 0.4 mm, cutting thickness 0.6 mm, with stopper at 15 mm, working length 25 mm



28207 SLD

28207 SLD **Same**, blade width 6 mm



28207 SLE

28207 SLE **Same**, blade width 8 mm



28207 SLF

28207 SLF **Saw Blade**, blade width 12 mm, blade thickness 0.4 mm, cutting thickness 0.6 mm, working length 22 mm



28207 SLG

28207 SLG **Saw Blade**, blade width 6 mm, blade thickness 0.4 mm, cutting thickness 0.6 mm, working length 31 mm



28207 SLH

28207 SLH **Same**, blade width 10 mm







## ■ HIGH FREQUENCY SURGERY UNITS

*NEW* AUTOCON® II 400 SCB


*NEW* AUTOCON® II 200

*NEW* AUTOCON® II 80



**Special Features:**

- For interdisciplinary use
- Cutting-edge next-generation unit with a convenient, easy-to-disinfect touch screen
- Bi-Vascular-Safe mode for bipolar coagulation and thermofusion of large-lumen vessels
- Equipped with 2 bipolar or 2 unipolar HF outputs depending on unit
- Defibrillator-safe CF outputs for maximum patient and user safety
- Permanent safety due to continuous monitoring of the contact between the neutral electrode and the patient during unipolar use
- 2 freely programmable foot pedals can be connected simultaneously
- Automatic activation of HF energy or via handswitch or footswitch depending on mode used
- Self-test for maximum patient and user safety
- Unit versions for standard use, bipolar resection and thermofusion of large-lumen vessels according to individual user requirements
- Up to 28 pre-programmed procedures for a wide range of disciplines make the unit extremely simple to operate. 100 program memories are available for individual programming.



### **Integrated voltage stability control or arc control**

The two state-of-the-art generators both guarantee optimal surgical cutting and coagulation power that adapts continuously to the particular indication, especially when there are wide variations in the tissue structures and tissue impedances. At the press of a button, the operator can switch between arc control (TOP Cut mode) and voltage stability control (POWER Cut mode).

### **TOP Cut mode**

In this mode, the HF energy required for a cutting effect is automatically reduced to the necessary physical minimum in each individual case. The electric arc always remains at a constant level, thus ensuring a uniform surgical effect. Such precise work offers a clear safety advantage and produces a cut that conserves tissue and reduces stress for the patient.

This increased safety with the AUTOCON®II 400 SCB unit is achieved by using the latest and fastest micro-processor and sensor technology, which enables the unit to capture all the important parameters, such as variable incision speed, geometry of the active electrode, different impedance behavior of biological tissue types and fluids, and transition and contact resistances. This data is then used to adjust the HF power output and HF voltage. This means that the operator is not restricted in his work or obliged to adapt to the HF unit. On the contrary, the HF unit adapts optimally to all the user's application and operating techniques.

### **POWER Cut mode with constant HF voltage and power**

This mode ensures a uniform surgical effect and consistent cutting efficiency over a wide impedance range and many different tissue types.

### **RAM system – return electrode application monitoring**

This safety system continuously monitors the contact quality between the neutral electrode and the patient's skin, and additionally indicates it with symbols. If the contact surface area decreases, the safety system gives an early visual and acoustic alarm, thus preventing a burn under the application site of the patient electrode. To increase the contact reliability of the neutral electrode, the user can prevent the application of single-faced electrodes.

### **Gastro Cut and Papillo Cut**

These two new resection modes, which were specially developed by KARL STORZ for use in flexible endoscopy, permit a fractionated and controlled cut with no bleeding. The special HF generator technology allows controlled output of cutting and coagulation pulse current. Both the pulse sequence and the pulse speed can be set separately and specifically for each mode.

### **LF/HF leakage currents**

Stray currents and the associated risk of burns are minimized by design measures.

### **SCB- and OR1™-compatible**

The units are designed for integration in the KARL STORZ Communication Bus (SCB). Full integratability of the HF unit in the networked, and even speech-controlled, operating room of the future is already a standard feature (system requirement: RUI software release 2009001-26 or later).

### **C-Cut® mode and LAP-C-Cut mode – the intermittent coagulation-cutting mode with AUTOCON®II 400 SCB**

Designed for blood-free cutting during laparoscopy and when using irrigation liquids, the electric current is specially modulated and offers a reproducible coagulation current with high cutting efficiency. This makes time-consuming subsequent coagulation a thing of the past. Using the C-Cut® mode thus reduces the need for blood transfusions and saves surgery time. The result is an overall reduction in operating costs, and for the patient it means additional protection.

### **Bipolar generator with 370 W HF power**

This outstanding performance range allows the unit to be used in conjunction with the newly developed bipolar special and standard accessories. This high-end unit is even suitable for indications using irrigation liquids, which in the past could only be performed with special-purpose HF systems. The safety of the bipolar technology from KARL STORZ eliminates the need for applying a neutral electrode, including for interventions which used to be standard unipolar procedures.

### **Unipolar generator with 300 W HF power**

With a peak power of 300 watts, the AUTOCON®II 400 SCB is ideally equipped for interventions in all areas of use.

# AUTOCON® II 400 SCB <sup>NEW</sup>

Safety Aspects and Main Features

**STORZ**  
KARL STORZ — ENDOSKOPE

## Precise power setting and power limitation

Exact fine tuning in 1 W steps is provided for interventions requiring maximum precision with very low power.

## Up to 8 hemostatic effects

Individual selection of up to 8 hemostatic effects for unipolar and bipolar cutting, each with up to 370 W output, permits optimal control of coagulation and the surgical effect in every situation.

## Forceps auto-start function

When the forceps tips contact the tissue, bipolar coagulation is activated automatically after a freely adjustable delay of up to 9.9 seconds.

## 6.5" touch screen

The color touch screen makes the AUTOCON®II 400 SCB the world's first HF unit with this new, user-friendly operating technology. It also offers the important advantage of very easy cleaning and wipe-down disinfection.

## 100 program memories

Simple programming of the indication-related unit parameters makes the AUTOCON®II 400 SCB easy and intuitive to operate because all the programs can be stored in numeric order or text-based with user name and indication. The stored programs can be called up in the indication list at the touch of a button.

## Bipolar coagulation auto-stop function

Automatic power shutdown when the coagulation procedure has been ended.

## Self-test program

A comprehensive software safety concept ensures smooth, safe use after switching on. Detected component faults are indicated by an error code display, enabling rapid troubleshooting. The self-test also includes the connected accessories for the specific purpose of minimizing waiting times in preoperative work-up.

## Software upgrade

The service port on the back of the unit allows the HF functions of the AUTOCON®II 400 SCB to be expanded economically for future forms of HF treatment. This means that the AUTOCON®II 400 SCB is always up-to-date.



# AUTOCON® II 400 SCB, <sup>NEW</sup> AUTOCON® II 200/80

Special Features

Special Features:	AUTOCON® II 400 SCB	AUTOCON® II 200	AUTOCON® II 80
Degree of coagulation (effect) can be preselected in several steps: The degree of coagulation measures extent of coagulation depth	●	●	-
Bi-Vascular-Safe mode for bipolar coagulation and thermofusion of large-lumen vessels	●	-	-
Easy to use due to automatic mode selection thanks to recognition of instrument-cable connection	●	-	-
Spray coagulation: Coagulation with modulated HF voltage (Up > 500 V); very long arcs enable coagulation of large and bleeding areas of tissue without contact to tissue	●	-	-
Autostart function: Manual adjustment of operating time limit for bipolar coagulation	●	●	-
Voltage-regulated cutting	●	●	-
Separate papillo-cut and gastro-cut functions enable fractionated cutting with regulated HF current at different cutting speeds for flexible endoscopy	●	●	-
Autostart function for bipolar coagulation: Automatic activation of coagulation current as soon as coagulation electrode touches tissue with both branches	●	●	-
Activation of HF functions possible via footswitch or manual control switch for unipolar or bipolar	●	●	-
Bipolar resection with KARL STORZ bipolar electrotomes	●	-	-
Bipolar application with NaCl irrigation solution	●	-	-
Modular connecting sockets for unipolar and bipolar applications can be selected according to individual requirements	●	●	●
100 applications with text can be stored	●	-	-
Convenient use via 6.5" Touch Screen	●	-	-
Switchover function enables switching between two modes within a user program via a footswitch from the sterile area	●	-	-
Arc-controlled cutting, unipolar	●	-	-
Compatible with KARL STORZ Communication Bus (KARL STORZ-SCB)	●	-	-
Service port for software updates and HF functionality upgrades	●	-	-

HIGH FREQUENCY SURGERY UNITS

# AUTOCON® II 400 SCB <sup>NEW</sup>

Specifications

HF Modes	Effects	P max. at 500 Ohm	V <sub>p</sub> max. at 500 Ohm	Crest Factor	Arc Control	Voltage Control
<b>Unipolar</b>						
TOP-Cut	8	300	950	1.4	●	–
POWER-Cut	8	300	740	1.4	–	●
C-Cut®	8	200	1450	3.0	–	●
LAP-C-Cut	8	200	1450	3.0 – 3.8	–	●
Gastro-Cut	4	200	880	1.4	–	●
Papillo-Cut	4	200	880	1.4	–	●
Standard Coag	8	200 (at 50 Ohm)	190	1.4	–	●
Forced Coag	4	120	1800	6.0	–	●
Spray Coag	2	120	4300	7.4	–	●
<b>Bipolar</b>						
Bipolar-Cut	8	100	740	1.4	–	●
Saline-C-Cut	8	370	770	1.4	–	●
Saline-C-Cut ++*	8	300 (at 750 Ohm)	490	1.4	–	●
Saline-Time-C-Cut	8 time 0.1-1 sec.	370	770	1.4	–	●
Saline-Time-C-Cut ++*	8	300 (at 750 Ohm)	490	1.4	–	●
Saline Coag	8	200 (at 50 Ohm)	190	1.4	–	●
Saline Coag ++*	8	200 (at 50 Ohm)	190	1.4	–	●
Saline-Time-Coag	8 time 0.1-1 sec.	200 (at 50 Ohm)	190	1.4	–	●
Saline-Time-Coag ++*	8 time 0.1-1 sec.	200 (at 50 Ohm)	190	1.4	–	●
Bipolar Soft Coag	8	120 (at 75 Ohm)	190	1.4	–	●
Bipolar Soft with Auto-Stop	8	120 (at 75 Ohm)	190	1.4	–	●
Bi-Vascular-Safe**	8	300 (at 25 Ohm)	220	1.4	–	●

\*Only for units with additional resection module

\*\*with software package “Bi-Vascular-Safe”

**Specifications:**

Safety systems	<ul style="list-style-type: none"> <li>- Automatic self-test</li> <li>- Maldosage</li> <li>- Neutral electrode safety system (dynamic, two-part, one- and two-part NE)</li> <li>- LF/HF leakage current monitor</li> <li>- Activation time</li> <li>- Deactivable HF</li> </ul>	Power supply	<b>20 5352 2x-12x:</b> 220-240 VAC, 50/60 Hz <b>20 5352 2xU12x:</b> 100-120 VAC, 50/60 Hz
		Dimensions w x h x d	448 x 164 x 345 mm
		Weight	10 kg
		Certified to	IEC 60601-1, CE acc. to MDD

# AUTOCON® II 400 SCB <sup>NEW</sup>

High Frequency Surgery Unit,  
Recommended System Configurations



**AUTOCON® II 400 SCB,**  
power supply 220 – 240 VAC, 50/60 Hz  
including:  
**Mains Cord**  
**SCB Connecting Cable,** length 100 cm

**AUTOCON® II 400 SCB,**  
power supply 100 – 120 VAC, 50/60 Hz  
including:  
**Mains Cord**  
**SCB Connecting Cable,** length 100 cm

Application	Standard: Unipolar/Bipolar	High-End
Unit version	-122 (220 – 240 VAC) U122 (100 – 120 VAC)	-125 (220 – 240 VAC) U125 (100 – 120 VAC)
Product No.	<b>20 5352 01-122</b> <b>20 5352 01U122</b>	<b>20 5352 01-125</b> <b>20 5352 01U125</b> basic unit
	–	<b>20 5352 02-125</b> <b>20 5352 02U125</b> basic unit, incl. additional resection module
	–	<b>20 5352 03-125</b> <b>20 5352 03U125</b> basic unit, incl. Bi-Vascular-Safe mode
	–	<b>20 5352 04-125</b> <b>20 5352 04U125</b> basic unit, incl. additional resection module + Bi-Vascular-Safe mode

### Socket Position








1	Bipolar Combination 	Bipolar Combination 
2	Bipolar Combination 	Bipolar Multifunction 
3	Unipolar 3-pin and Erbe 	Unipolar 3-pin and Erbe 
4	NE 6.3 mm jack and 2-pin 	NE 6.3 mm jack and 2-pin 

**Optional Accessories for AUTOCON® II 400 SCB** see pages U 46-51  
**Components/Spare Parts** see chapter 14



# AUTOCON® II 400 SCB <sup>NEW</sup>

System Components



2-12

HIGH FREQUENCY SURGERY UNITS

# AUTOCON® II 200 <sup>NEW</sup>

High Frequency Surgery Unit,  
Recommended System Configuration




- 20 5322 01 **AUTOCON® II 200**, power supply  
200 – 240 VAC, 50/60 Hz  
including:  
**Mains Cord**
- 20 5322 01 C **AUTOCON® II 200**, power supply  
100 – 120 VAC, 50/60 Hz  
including:  
**Mains Cord, US version**
- 20 5322 01-010 **AUTOCON® II 200**, with GASTRO-Cut,  
power supply 200 – 240 VAC, 50/60 Hz  
including:  
**Mains Cord**
- 20 5322 01C010 **AUTOCON® II 200**, with GASTRO-Cut,  
power supply 100 – 120 VAC, 50/60 Hz  
including:  
**Mains Cord, US version**

**Specifications:**

HF rated power	<ul style="list-style-type: none"> <li>- <b>Cutting</b> Power-Cut (unipolar): 220 Watt/500 Ohm GASTRO-Cut (unipolar): 200 Watt/75 Ohm</li> <li>- <b>Coagulation</b> Standard-Coag (unipolar): 120 Watt/50 Ohm Forced-Coag (unipolar): 120 Watt/500 Ohm Standard-Coag (bipolar): 120 Watt/75 Ohm</li> </ul>
Max. voltage	<ul style="list-style-type: none"> <li>- <b>Cutting</b> Power-Cut (unipolar): 740Vp GASTRO-Cut (unipolar): 550Vp</li> <li>- <b>Coagulation</b> Standard-Coag (unipolar): 190Vp Forced-Coag (unipolar): 1800Vp Standard-Coag (bipolar): 190Vp</li> </ul>
Degree of coagulation	<ul style="list-style-type: none"> <li>Power-Cut: 8</li> <li>GASTRO-Cut: 4</li> <li>Standard-Coag: 8</li> <li>Forced-Coag: 4</li> <li>Bipolar Standard-Coag: 8</li> </ul>

Autostart	<b>Coagulation</b> bipolar
Autostop	<b>Coagulation</b> unipolar standard
Safety systems	<ul style="list-style-type: none"> <li>- permanent power control</li> <li>- maldosage</li> <li>- neutral electrode safety system</li> <li>- automatic Self-Test</li> </ul>
Power supply	100-120 VAC or 200-240 VAC, 50/60 Hz
Dimensions w x h x d	410 x 165 x 385 mm
Weight	10 kg
Certified to	IEC 60601-1, CE acc. to MDD

**Optional Accessories for AUTOCON® II 200** see pages U 46-51  
**Components/Spare Parts** see chapter 14



2-12

HIGH FREQUENCY  
SURGERY UNITS

**AUTOCON® II 80** <sup>NEW</sup>  
 High Frequency Electrosurgical Unit,  
 Recommended System Configurations



**20 5308 01**     **AUTOCON® II 80,**  
 power supply 100 – 240 VAC, 50/60 Hz  
 including:  
**Mains Cord**


**Specifications:**

HF rated power	<ul style="list-style-type: none"> <li>- Dry Cut: 80 Watt/500 Ohm</li> <li>- Forced Coag: 50 Watt/500 Ohm</li> <li>- Soft Coag: 80 Watt/100 Ohm</li> <li>- Auto Cut: 80 Watt/500 Ohm</li> <li>- Bipolar Soft: 80 Watt/100 Ohm</li> </ul>	Safety systems	<ul style="list-style-type: none"> <li>- permanent power control</li> <li>- maldosage</li> <li>- neutral electrode safety system</li> <li>- automatic self-test</li> </ul>
Max. voltage	<ul style="list-style-type: none"> <li>- Dry Cut: 830Vp</li> <li>- Forced Coag: 1200Vp</li> <li>- Soft Coag: 180Vp</li> <li>- Auto Cut: 500Vp</li> <li>- Bipolar Soft: 180Vp</li> </ul>	Power supply	100-240 VAC, 50/60 Hz
		Dimensions w x h x d	280 x 135 x 300 mm
		Weight	4 kg
		Certified to	IEC 60601-1, CE acc. to MDD

**Optional Accessories for AUTOCON® 80** see pages U 46-51  
**Components/Spare Parts** see chapter 14

# AUTOCON® II 80 <sup>NEW</sup>

System Components



2-12

HIGH FREQUENCY SURGERY UNITS

# Optional Accessories

for AUTOCON® II 400 SCB and AUTOCON® II 200/80

			for use with		
			AUTOCON® II 400 SCB	AUTOCON® II 200	AUTOCON® II 80
	20017831	<b>Three-Pedal Footswitch</b> , with switch function for switching between mode levels A and B, for use with AUTOCON® II 400 SCB	●	-	-
	20017830	<b>Two-Pedal Footswitch</b> , for use with AUTOCON® II 400 SCB and AUTOCON® II 200	●	●	-
	20017832	<b>One-Pedal Footswitch</b> , for activating coagulation, for use with AUTOCON® II 400 SCB and AUTOCON® II 200	●	●	-
	20017834	<b>Two-Pedal Footswitch</b> , digital, one-stage	-	-	●
	20017833	<b>One-Pedal Footswitch</b> , digital, one-stage	-	-	●
	27805	<b>Neutral Electrode</b> , of conductive silicone with 2 rubber ties for fastening, contact surface A = 500 cm <sup>2</sup> , for use with Connecting Cable 27806	●	●	●
	860021 E	<b>Neutral Electrode</b> , of conductive silicone, with 1 rubber tie for fastening, contact surface A = 187 cm <sup>2</sup> , for use with Connecting Cable 27806	-	-	●
	27806	<b>Neutral Electrode Connecting Cable</b> , for Neutral Electrodes 27805 and 860021 E, length 400 cm	-111 -115 -122 -125	●	●
	27806 UR	<b>Neutral Electrode Connecting Cable</b> , for Neutral Electrode 27805	-112 -116 -122 -125	●	●
	27806 US	<b>Neutral Electrode Connecting Cable</b> , for Neutral Electrode 27802	-112 -116 -122 -125	●	●
	27802	<b>Neutral Electrode</b> , for single use, contact surface divided into two, A = 169 cm <sup>2</sup> , package of 50, Connecting Cable 27801 required	●	●	●
	27801	<b>Connecting Cable</b> , for connecting Neutral Electrode 27802, length 500 cm	-111 -115 -122 -125	●	●
	26520043	<b>Electrode Handle</b> , with 2 buttons for activating the unipolar generator, yellow button: unipolar cutting, blue button: unipolar coagulation (Cable 26520045 required)	●	●	●
	26520045	<b>High Frequency Cable</b> , for Electrode Handle 26520043, length 400 cm	-111 -115 -122 -125	●	●
	26520046	<b>Electrode Handle</b> , without buttons, with integrated connecting cable, length 300 cm	-111 -115 -122 -125	●	●
	26176 BP <i>NEW</i>	<b>Bipolar Adaptor</b> , for AUTOCON® to 66177 CC/CCL and 26520048	●	●	●

HIGH FREQUENCY  
SURGERY UNITS

# Surgery Electrodes Set

Accessories



20 5300 08

**Surgery Electrodes Set**  
including:

**Container with Lid and Sterilizing Insert**, for 16 electrodes with diameter 4 mm



**Wire Snare**, 5 mm



**Same**, 10 mm



**Ribbon Snare**, 10 mm



**KIRSCHNER Spatula Electrode**, straight



**MAGENAU Knife Electrode**, angled



**Same**, lancet-shaped



**Ball Electrode**, 2 mm



**Same**, 4 mm



**Same**, 6 mm



**Needle Electrode**



**Flat Electrode**, 8 x 10 mm




**Same**, 10 x 15 mm

For use with Electrode Handles 26 5200 43 and 26 5200 46  
Components/Spare Parts see chapter 14





# Accessories

## Unipolar High Frequency Cords



### Unipolar High Frequency Cords

KARL STORZ High Frequency  
Instrument Surgery Units

	26002 M	<b>Unipolar High Frequency Cord</b> , with 4 mm plug for models KARL STORZ, Erbe type T, older models and Ellman, length 300 cm
	26004 M	<b>Unipolar High Frequency Cord</b> , with 4 mm plug for models Berchtold and Martin, length 300 cm
	26005 M	<b>Unipolar High Frequency Cord</b> , with 5 mm plug for AUTOCON® II 400 SCB system (111, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ AUTOCON® system (50, 200, 350) and Erbe type ICC, length 300 cm
	26006 M	<b>Unipolar High Frequency Cord</b> , with 8 mm plug for models KARL STORZ AUTOCON® II 400 SCB system (112, 116) and Valleylab, length 300 cm

**Please note:** All high frequency cords of this page are delivered with a length of 300 cm. If a length of 500 cm is requested please add letter **L** to the part number, e. g. 26002 **ML**, 26176 **LVL**.






# Accessories

## Bipolar High Frequency Cords


### Bipolar High Frequency Cords



KARL STORZ High Frequency  
Instrument Surgery Units

	<p>26176 LE <b>Bipolar High Frequency Cord</b>, for AUTOCON® II 400 SCB system (111, 113, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ Coagulator 26021 B/C/D, 860021 B/C/D, 27810 B/C/D, 28810 B/C/D, AUTOCON® system (50, 200, 350), Erbe-Coagulator, T and ICC series, length 300 cm</p>
	<p>26176 LM <b>Bipolar High Frequency Cord</b>, for Martin and Berchtold coagulators, length 300 cm</p>
	<p>26176 L <b>Bipolar High Frequency Cord</b>, with 2x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Ellman, length 300 cm</p>
	<p>26176 LA <b>Bipolar High Frequency Cord</b>, with 2x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Valleylab, length 300 cm</p>
	<p>26176 LV <b>Bipolar High Frequency Cord</b>, for KARL STORZ AUTOCON® II 400 SCB system (112, 114, 116, 122, 125) and Valleylab coagulators, length 300 cm</p>




	<p>205400 21 <b>BiVascularClamp</b>, total length 20 cm, with thermoinsulation, for open surgery, integrated MF cable length 4 m</p>
	<p>205400 22 <b>Same</b>, total length 21 cm, with corrugated jaws</p>
	<p>205400 23 <b>Same</b>, total length 27 cm</p>

# Accessories

## Bipolar Coagulating Forceps

Bipolar Coagulating Forceps, insulated  
for use with Bipolar High Frequency Cords 847000 or 847000 A/E/M/V



	841045	<b>Bipolar Coagulating Forceps</b> , insulated, angled tip, sharp, very delicate, length 10.5 cm
	841046	<b>Same</b> , blunt, tip 1 mm
	841111	<b>Bipolar Coagulating Forceps</b> , insulated, straight, blunt, length 11 cm
	841116	<b>Bipolar Coagulating Forceps</b> , insulated, straight, pointed, length 16 cm
	841216	<b>Same</b> , blunt, tip 1 mm, length 16 cm
	841219	<b>Same</b> , blunt, tip 1 mm, length 19 cm
	842016	<b>Bipolar Coagulating Forceps</b> , insulated, angled tip, very delicate, tip 0.5 mm, length 16 cm
	842116	<b>Bipolar Coagulating Forceps</b> , insulated, angled tip, tip 0.5 mm, length 16 cm
	842111	<b>Bipolar Coagulating Forceps</b> , insulated, angled tip, pointed, length 11 cm
	842216	<b>Same</b> , blunt, tip 1 mm, length 16 cm
	842219	<b>Same</b> , blunt, tip 2 mm, length 19 cm
	842319	<b>Same</b> , blunt, tip 1 mm, length 19 cm






# Accessories

High Frequency Cords  
for use with Bipolar Coagulation Instruments



## Bipolar High Frequency Cords

KARL STORZ Instrument High Frequency Electrosurgery Units

	<p>847000 E <b>Bipolar High Frequency Cord</b>, for AUTOCON® II 400 SCB systems (111, 113, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ Coagulator 26021 B/C/D, 860021 B/C/D, 27810 B/C/D, 28810 B/C/D, AUTOCON® systems (50, 200, 350), Erbe-Coagulator, T and ICC series, length 300 cm</p>
	<p>847000 <b>Bipolar High Frequency Cord</b>, with 2x 4 mm banana plug to KARL STORZ Coagulator 26020 XA/XB for KARL STORZ bipolar coagulation forceps, length 300 cm</p>
	<p>847000 A <b>Bipolar High Frequency Cord</b>, with 2x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Valleylab, length 300 cm</p>
	<p>847000 M <b>Bipolar High Frequency Cord</b>, for Martin and Berchtold coagulators, length 300 cm</p>
	<p>847000 V <b>Bipolar High Frequency Cord</b>, for KARL STORZ AUTOCON® II 400 SCB systems (112, 114, 116, 122, 125) and Valleylab coagulators, length 300 cm</p>

**Please note:**

All high frequency cords on this page are delivered with a length of 300 cm only.

# Compatibility

High Frequency Cords to AUTOCON® HF Electrosurgical Units

<b>AUTOCON® II 80</b>	20 5308 20	–	20 5308 20	–	–
<b>AUTOCON® II 200</b>	20 5225 20 20 5225 20C 20 5225 20-010 20 5225 20C010	–	20 5225 20 20 5225 20C 20 5225 20-010 20 5225 20C010	–	–
<b>AUTOCON® II 400 SCB</b>	20 5352 20-111 20 5352 20-115	20 5352 20-112 20 5352 20-116	20 5352 20-122	20 5352 20-125 20 5352 21-125	20 5352 22-125 20 5352 23-125

## Unipolar High Frequency Cords

27806	●	–	●	●	●
27801	●	–	●	●	●
27806 UR	–	●	●	●	●
27806 US	–	●	●	●	●
26 5200 45	●	–	●	●	●
26 5200 46	●	–	●	●	●
26002 M	●	●	●	●	●
26005 M	●	–	●	●	●
26006 M	–	●	–	–	–


## Bipolar High Frequency Cords

26176 LE	●	–	●	●	●
26176 LV	–	●	●	●	●
847000 E	●	–	●	●	●
847000 V	–	●	●	●	●

## Bipolar High Frequency Cords/Instruments to Multifunction Socket

<b>AUTOCON® II 400 SCB</b>	20 5352 20-115	20 5352 20-116	–	20 5352 20-125 20 5352 21-125	20 5352 22-125 20 5352 23-125
20 5400 21	–	–	–	–	●
20 5400 22	–	–	–	–	●
20 5400 23	–	–	–	–	●

**EXCERPT FROM**  
**TELEPRESENCE**  
**IMAGING SYSTEMS - DOCUMENTATION**  
**ILLUMINATION - EQUIPMENT CARTS**





# Excerpt from TELEPRESENCE

Imaging Systems – Documentation – Illumination – Equipment Carts



*NEW*

*NEW*

*NEW*

## ■ Camera Systems

- **IMAGE 1 HUB™ HD**  
Brilliant Imaging FULL HD
- **KARL STORZ 3D System**  
Excellent Depth Perception in 3D
- **TRICAM®**
- **TELE PACK VET X**
- **TELECAM**
- **Veterinary Video Camera III**
- **C-HUB™ and C-CAM™**

## ■ HD Imaging with Operating Microscopes

## ■ Data Management and Documentation

## ■ Monitors

*NEW*

- **KARL STORZ FULL HD Monitors**
- **KARL STORZ HD and HD WIDEVIEW™ Monitors**
- **15" and 19" Touch Screens**

*NEW*

- **ZeroWire**
- **VET Monitors**

## ■ Light Sources

*NEW*

- **HALOGEN**
- **LED**
- **XENON**

## ■ SMARTCART Equipment Cart Series

## ■ Components Spare Parts

CAMERA SYSTEMS

DATA MANAGEMENT  
AND DOCUMENTATION

MONITORS

LIGHT SOURCES

EQUIPMENT CARTS

COMPONENTS  
SPARE PARTS





# Table of Contents

Excerpt from TELEPRESENCE

## TABLE OF CONTENTS

Table of Contents	I-II
Index	III-V
Numerical Index	VI-VIII

## IMAGE 1 HUB™ HD – BRILLIANT IMAGING IN FULL HD

Brilliant Imaging in FULL HD	TP 7
KARL STORZ FULL HD – Highest Standards for Image Input and Image Display	TP 8
IMAGE 1 HUB™ HD – FULL HD Camera Head, FULL HD Camera Control Unit, FULL HD Monitor	TP 9
IMAGE 1 HUB™ HD – FULL HD Camera Control Unit	TP 10-TP 12
IMAGE1 HD – FULL HD Camera Control Unit	TP 13-TP 14
IMAGE1 FULL HD Camera Heads – for Endoscopy	TP 15-TP 16
IMAGE1 FULL HD Camera Heads – for Microscopy and ICG Application	TP 17
IMAGE1 Camera Heads	TP 18
IMAGE 1 HUB™ HD, IMAGE1 – Accessories	TP 19-TP 20

## KARL STORZ 3D SYSTEM

The KARL STORZ 3D System – Discover a New World in 3D	TP 22
KARL STORZ 3D System – 3D Camera Control Unit, 3D TIPCAM®1	TP 23
KARL STORZ 3D System – 3D Monitors	TP 24
Accessories – for KARL STORZ 3D System	TP 25

## TRICAM® – THREE-CHIP CAMERA SYSTEM

TRICAM® SL II – Three-Chip Camera System, Color Systems PAL/NTSC, Digital Imaging Processor	TP 28
TRICAM® SL II – Camera Control Unit SCB	TP 29
TRICAM® SL II – Camera Heads	TP 30-TP 31
TRICAM® SL II – Optional Accessories	TP 32

## TELE PACK VET X COMPACT SYSTEM AND TELECAM ONE-CHIP CAMERA SYSTEM

TELE PACK VET X – A multipurpose unit for diagnostic and surgical procedures including still image and video capture software.	TP 34-TP 35
--	-------------

## TELECAM

TELECAM – One-Chip Camera System, Color Systems PAL/NTSC, Digital Imaging Processor	TP 38
TELECAM SL II – Camera Control Unit	TP 39
Veterinary Video Camera III – One-Chip Camera System, Color Systems PAL/NTSC	TP 40
Veterinary Video Camera III – Camera Control Unit, Camera Head	TP 41
TELE PACK VET X, TELECAM and Veterinary Video Camera III – Camera Heads	TP 42
TELE PACK VET X, TELECAM and Veterinary Video Camera III – Accessories	TP 43-TP 44

## OVERVIEW AND ACCESSORIES

Camera Heads, Camera Control Unit – Matrix	TP 46-TP 47
Video Endoscopes, Camera Control Unit – Matrix	TP 48-TP 49
Video Cable and SCB Cable	TP 50
Sterile Covers – for all KARL STORZ Camera Heads	TP 51

## C-HUB™ AND C-CAM™

C-HUB™ and C-CAM™	TP 54
C-MAC® Monitor System	TP 55

# Table of Contents

Excerpt from TELEPRESENCE

## **DATA MANAGEMENT AND DOCUMENTATION – KARL STORZ AIDA®**

Data Management and Documentation – KARL STORZ AIDA® compact NEO advanced/standard	TP 58-TP 61
Data Management and Documentation – Options – Network Integration, DICOM, HL7	TP 62
Accessories – for AIDA compact NEO advanced/standard	TP 63-TP 68

## **MONITORS**

KARL STORZ FULL HD Monitors	TP 70-TP 71
KARL STORZ HD, HD WIDEVIEW™ and HD LED Backlight Monitors	TP 72-TP 75
15" and 19" Touch Screens	TP 76-TP 77
ZeroWire – Wireless FULL HD Video Transmission System	TP 78-TP 79
Color Monitors – for Video Cameras, Color Systems PAL/NTSC, without IEC 60601-1 certification	TP 80
Accessories and Optional Accessories – for HD WIDEVIEW Monitors, TFT Flat Screens and Touch Screens	TP 81

## **LIGHT SOURCES**

Light Sources	TP 84
Cold Light Fountain HALOGEN	TP 85
Cold Light Fountain HALOGEN – Overview	TP 86
Cold Light Fountain VETERINARY	TP 87
Cold Light Fountain HALOGEN 250 twin	TP 88
Cold Light Fountains – Accessories	TP 89
Cold Light Fountain LED	TP 91-TP 92
Cold Light Fountain LED NOVA 150	TP 93
Cold Light Fountain LED NOVA 100 twin	TP 94
Cold Light Fountain Power LED 175 SCB	TP 95
Cold Light Fountains XENON	TP 97
Cold Light Fountains XENON – Overview	TP 99
Cold Light Fountain XENON 100 SCB – with Integrated Insufflation Pump	TP 100
Cold Light Fountain XENON NOVA® 175	TP 101
Cold Light Fountain XENON 175 SCB	TP 102
Cold Light Fountain XENON NOVA® 300	TP 103
Cold Light Fountain XENON 300 SCB	TP 104
Cold Light Fountains XENON – Spare Lamps and Spare Lamp Modules	TP 105
Fiber Optic Light Cables – for Cold Light Fountains	TP 106-TP 107
Adaptors – for Light Sources and Light Cables	TP 108
Adaptor Series – for HALOGEN and XENON Light Sources	TP 109

## **SMARTCART EQUIPMENT CART SERIES**

SMARTCART 2 – The New Generation of KARL STORZ Equipment Carts	TP 112
Equipment Carts	TP 113-TP 117, TP 125
Equipment Carts – Accessories	TP 118-TP 124
Equipment Carts E Series	TP 126
Equipment Carts E Series – Accessories	TP 127

## **COMPONENTS, SPARE PARTS**

TP 129-TP 139

# Index

Excerpt from TELEPRESENCE

## 1

15" KARL STORZ Touch Screen	TP 65, TP 66, TP 76, TP 77, TP 135
15" Monitor	TP 73, TP 75, TP 136
19" HD Monitor	TP 73, TP 75, TP 135, TP 136
19" KARL STORZ HD Monitor	TP 71, TP 134
19" KARL STORZ TFT Flat Screen	TP 63, TP 64, TP 135, TP 136
19" KARL STORZ Touch Screen	TP 65, TP 66, TP 76, TP 77, TP 135

## 2

24" 3D Monitor	TP 24
24" HD WIDEVIEW Monitor	TP 73, TP 75, TP 136
26" HD Monitor	TP 136
26" HD Monitor (Highbright with LED Backlight)	TP 73, TP 75, TP 136
26" HD Monitor with LED Backlight	TP 73, TP 75, TP 136
26" KARL STORZ HD Monitor	TP 71, TP 134

## 3

32" 3D Monitor	TP 24
3D Camera Control Unit SCB	TP 131
3D Camera Control Unit SCB, with ICM Module	TP 23, TP 131
3D Clip-on Glasses	TP 25
3D Glasses	TP 25, TP 131
3D TIPCAM®1	TP 23

## A

Accessory Side Shelf	TP 118
Adaptor	TP 35, TP 43, TP 81, TP 100, TP 108, TP 109, TP 135
Adaptor ACC 1	TP 20, TP 32, TP 43
Adaptor Plate	TP 24
AIDA compact HD Frame Grabber Board	TP 134
Air Filter	TP 100

## B

Bar Code Scanner	TP 67
Basic Equipment Cart	TP 116
BNC/BNC Video Cable	TP 50, TP 130, TP 131, TP 132, TP 133, TP 134, TP 135, TP 136
BNC Connecting Cable	TP 134
Bottle Holder	TP 122, TP 125

## C

Cable Conduit Horizontal	TP 119
Cable Manager	TP 121
Camera Cover	TP 51
Camera Holder	TP 122
Camera Mount	TP 122
Carrier Set	TP 122

Case	TP 132
C-CAM™ Camera Head, 6-pin	TP 54
C-CAM™ Camera Head, 8-pin	TP 54
CD-R	TP 68
C-HUB™ Camera Control Unit	TP 54, TP 133
C-MOUNT Lens	TP 43, TP 132
Cold Light Fountain HALOGEN 250 twin	TP 88, TP 137
Cold Light Fountain LED NOVA 100 twin	TP 94, TP 138
Cold Light Fountain LED NOVA 150	TP 93, TP 138
Cold Light Fountain Power LED 175 SCB	TP 95, TP 138
Cold Light Fountain VETERINARY	TP 87, TP 137
Cold Light Fountain XENON 100 SCB	TP 100, TP 139
Cold Light Fountain XENON 175 SCB	TP 102, TP 139
Cold Light Fountain XENON 300 SCB	TP 104, TP 139
Cold Light Fountain XENON NOVA® 175	TP 101, TP 139
Cold Light Fountain XENON NOVA® 300	TP 103, TP 139
Color Ink Jet Printer	TP 67
Color Laser Printer	TP 67
Connecting Cable	TP 50, TP 130, TP 131, TP 132, TP 133, TP 134

## D

Double Pedal Footswitch	TP 67
Drawer Unit with Lock	TP 121
DV Connecting Cable	TP 132
DVD+R	TP 68
DVD+RW	TP 68
DVD-R	TP 68
DVI-D Connecting Cable	TP 50, TP 130, TP 131, TP 134, TP 135, TP 136

## E

Earth Leakage Monitor	TP 119
Equipment Cart	TP 113, TP 114, TP 115, TP 117, TP 125
Equipment Cart Protective Cover	TP 124
Equipment Rail	TP 127
Ethernet Network Board	TP 68
Europlug Adaptor	TP 133
Extension Cable	TP 44, TP 50

## F

Fiber Optic Light Cable	TP 25, TP 106, TP 107
Fiber Optic Network Board	TP 68
Fiber Optic Network Cable	TP 68
Fiberscope Carrier Set	TP 125
Fixation for Power Supply	TP 124

## H

Halogen Spare Lamp	TP 89
Holder	TP 122
Holder for Irrigation Liquid Bottle	TP 125

Holder for One-Pedal Footswitch	TP 124, TP 127	Media Converter	TP 68
Holder for Single-Use Container	TP 121	Mobile Stand	TP 125
Holder for Two-Pedal Footswitch	TP 123, TP 124, TP 127	Monitor Arm	TP 118
<b>I</b>			
IMAGE1 F3 Three-Chip ENT and Arthroscopy Camera Head	TP 18	Monitor for CMOS Endoscopes	TP 55, TP 133
IMAGE1 H3-FA Three-Chip FULL HD Camera Head	TP 16	Monitor Holder	TP 118
IMAGE1 H3-M COVIEW® Three-Chip FULL HD Microscope Camera Head	TP 17	Monitor Holding Arm	TP 113, TP 118
IMAGE1 H3-ZA Three-Chip FULL HD Camera Head	TP 16	Monitor Power Supply	TP 134
IMAGE1 H3-Z FI Three-Chip ICG FULL HD Camera Head	TP 17	Monitor Swivel Arm	TP 67, TP 113, TP 118
IMAGE1 H3-Z Three-Chip FULL HD Camera Head	TP 15	Multifunctional Holder	TP 122
IMAGE1 HD Camera Control Unit SCB	TP 13, TP 131	<b>O</b>	
IMAGE1 HD Camera Control Unit SCB, with ICM module	TP 13, TP 131	One-Pedal Footswitch	TP 20, TP 32, TP 35, TP 43
IMAGE1 HD Camera Head Cable	TP 19	OR1™ Network Isolator	TP 68
IMAGE 1 HUB™ HD Camera Control Unit SCB	TP 11, TP 130	Otoscope Camera Head Holder	TP 122
IMAGE 1 HUB™ HD Camera Control Unit SCB, with ICM module	TP 10, TP 130	<b>P</b>	
IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI/ICM module	TP 10, TP 130	Pedestal	TP 24, TP 66, TP 71, TP 73, TP 81, TP 135
IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI module	TP 11, TP 131	Plastic Container for Sterilization and Storage	TP 20, TP 32, TP 44
IMAGE 1 HUB™ HD Extension Cable	TP 19	Plug Adaptor Australia	TP 133
IMAGE1 SD Extension Cable	TP 19	Plug Adaptor UK Standard	TP 133
Instruction Manual	TP 137	Plug Adaptor USA/Japan	TP 133
Irrigator Rod	TP 122, TP 125	Potential Equalization Cabling	TP 120
Irrigator Rod, height adjustable	TP 122	Potential Equipotential Cable	TP 68, TP 120
Isolation Transformer	TP 119	Power Box	TP 119
<b>K</b>			
KARL STORZ AIDA® compact NEO advanced	TP 61, TP 134	Power Supply	TP 81, TP 133, TP 135, TP 136
KARL STORZ AIDA® compact NEO Software	TP 134	Power Supply C-HUB™	TP 133
KARL STORZ AIDA® compact NEO standard	TP 61, TP 134	Printer Cable	TP 68
KARL STORZ AIDA® control NEO	TP 134	PS/2 Compact Keyboard	TP 67, TP 134
KARL STORZ USB Flash Drive	TP 20, TP 25, TP 35, TP 44, TP 130, TP 131, TP 132, TP 134	Pump Holder	TP 123
Keyboard	TP 19, TP 32, TP 44, TP 130, TP 131, TP 132	<b>R</b>	
Keyboard Shelf	TP 127	Receiver	TP 137
Keyboard Tray	TP 121	Receptacle	TP 108
<b>L</b>			
LCD Color Monitor	TP 80	RGBS Connecting Cable	TP 50
Light Adaptor	TP 106	RJ45 Network Cable	TP 68
<b>M</b>			
Mains Adaptor	TP 66, TP 81	RS 232 Connecting Cable	TP 68, TP 81
Mains Cord	TP 81, TP 120, TP 130, TP 131, TP 132, TP 133, TP 134, TP 135, TP 136, TP 137, TP 138, TP 139	<b>S</b>	
Mains Cord Adaptor	TP 120	SCB Connecting Cable	TP 50, TP 130, TP 131, TP 132, TP 138, TP 139
		Screw Base	TP 108
		Shelf	TP 121, TP 127
		Shelf with Handles	TP 121
		SHORT-E HF Cart	TP 126
		Silicone Keyboard	TP 19, TP 25, TP 32, TP 35, TP 44, TP 67, TP 131
		Silicone Tubing Set	TP 139
		Single Pedal Holder	TP 124, TP 127
		Spare Lamp	TP 35, TP 43
		Special RGBS Connecting Cable	TP 50, TP 130, TP 131, TP 132

# Index

Excerpt from TELEPRESENCE

Standard Equipment Rail	TP 119
Sterile Cover	TP 66, TP 81
SVGA Connecting Cable	TP 68, TP 81
SVGA Monitor Cable	TP 50, TP 135
S-Video (Y/C) Connecting Cable	TP 50, TP 130, TP 131, TP 132, TP 133, TP 135, TP 136
S-Video (Y/C) Extension Cable	TP 50

## T

TELECAM One-Chip Camera Head	TP 42
TELECAM SL II Camera Control Unit	TP 39, TP 132
TELE PACK VET X	TP 35, TP 132
Touch Pen	TP 66, TP 81, TP 135
Transmitter	TP 137
Tray	TP 125
TRICAM® PDD Three-Chip Camera Head	TP 31
TRICAM® PDD Three-Chip Pendulum Camera Head	TP 31
TRICAM® SL II Camera Control Unit SCB	TP 29, TP 132
TRICAM® SL II Extension Cable	TP 32
TRICAM® Three-Chip Camera Head	TP 30
TROLL-E Mobile Stand	TP 25, TP 126
Tube Clip	TP 123
Two-Pedal Footswitch	TP 43

## U

USB Connecting Cable	TP 133
----------------------	--------

## V

Veterinary Video Camera C-MOUNT One-Chip Camera Head	TP 41, TP 42, TP 132
Veterinary Video Camera III	TP 40, TP 132
Veterinary Video Camera III Camera Control Unit	TP 41, TP 132, TP 133
Video Connecting Cable	TP 44
Video Endoscope Adaptor	TP 19
Video Endoscope Holder	TP 123
Video Endoscope Plug Intake	TP 123

## W

Wall Mount MOUNT-E™	TP 126
Wire Basket	TP 124, TP 127
Wire Tray for Cleaning, Sterilization and Storage	TP 25

## X

XENON Spare Lamp	TP 105
XENON Spare Lamp Module	TP 105

## Z

ZeroWire 2ft Spare DVI-D Cable	TP 79, TP 137
ZeroWire 42"/52"/55"	TP 79
ZeroWire Power Supply	TP 79
ZeroWire Radiance Y-Cable 19"/26"	TP 79, TP 137
ZeroWire Table Stand	TP 79
ZeroWire VESA Holder	TP 79, TP 137

# Numerical Index

Excerpt from TELEPRESENCE

040112-40	TP 51	20040281	TP 20, TP 25, TP 130, TP 131	20161101	TP 94, TP 138
040113-50	TP 51	20040372	TP 50, TP 135	20161120	TP 138
040114-15	TP 51	20040374	TP 50	20161201	TP 93, TP 138
040115-40	TP 51	20040402-20	TP 134	20161220	TP 138
040169-40	TP 51	20040577-V02	TP 134	20161401-1	TP 95, TP 138
041165-20	TP 66, TP 81	20040912-EN	TP 61, TP 134	20161420	TP 138
041265-20	TP 66, TP 81	20040913-EN	TP 61, TP 134	20200042	TP 43
103	TP 89	20046120-A	TP 134	20200043	TP 43, TP 132
105	TP 89	20090034	TP 68	20200073	TP 133
13991 DET	TP 123	20090036	TP 68	20200231	TP 19, TP 32, TP 44
13991 ET	TP 123	20090037	TP 68	20200231F	TP 19, TP 32, TP 44
20010170	TP 120	20090070	TP 50	20200231I	TP 19, TP 32, TP 44
20010270	TP 68, TP 120	20090170	TP 50, TP 130, TP 131, TP 132, TP 138, TP 139	20200231P	TP 19, TP 32, TP 44
20010370	TP 68, TP 120	20090171	TP 68	20200231S	TP 19, TP 32, TP 44
20010570	TP 120	20090234	TP 67	20200231U	TP 19, TP 32, TP 44, TP 130, TP 131, TP 132
20010670	TP 120	20090234U	TP 134	20200270	TP 44
20014130	TP 20, TP 32, TP 35, TP 43	20090270	TP 50	20203270	TP 50, TP 130, TP 131, TP 132
20014330	TP 43	20090331	TP 135	20204080	TP 68
20014430	TP 67	20090370	TP 50	20204081	TP 68
20020023	TP 125	20090386	TP 66, TP 81, TP 135	20204082	TP 68
20020031	TP 125	20090401	TP 65, TP 66, TP 76, TP 77, TP 135	20204083	TP 68
20020032	TP 125	20090403	TP 65, TP 66, TP 76, TP 77, TP 135	20212030	TP 42
20020046	TP 127	20090405	TP 65, TP 66, TP 76, TP 77, TP 135	20212040	TP 42
20020047	TP 127	20090406	TP 65, TP 66, TP 76, TP 77, TP 135	20212130	TP 42
20020049	TP 127	20090437	TP 135	20212140	TP 42
20020080	TP 126	20090470	TP 50	20213011U	TP 39, TP 132
20020081	TP 25, TP 126	20090481	TP 66, TP 81	20213020	TP 132
20020082	TP 126	20090483	TP 81, TP 135	20213070	TP 44
20030030	TP 100	20090485	TP 66, TP 81	20220270	TP 32
20040073	TP 68, TP 81	20090486	TP 66, TP 81, TP 135	20221030	TP 30
20040074	TP 68	20090487	TP 66, TP 81, TP 135	20221037	TP 31
20040076	TP 68	20090488	TP 81	20221039	TP 31
20040082	TP 132	20090570	TP 50	20221040	TP 30
20040086	TP 134	20090670	TP 50	20221070	TP 50, TP 130, TP 131, TP 132, TP 133, TP 134
20040087	TP 134	20090770	TP 50	20221075	TP 50
20040089	TP 50, TP 130, TP 131, TP 134, TP 135, TP 136	20113301	TP 88, TP 137	20221076	TP 50
20040174	TP 68	20113320	TP 137	20221077	TP 50
20040240DE	TP 35, TP 44	20131520	TP 139	20221130	TP 30
20040240ES	TP 35, TP 44	20132025	TP 105	20221137	TP 31
20040240FR	TP 35, TP 44	20132026	TP 105	20221139	TP 31
20040240IT	TP 35, TP 44	20132101-1	TP 102, TP 139	20221140	TP 30
20040240PT	TP 35, TP 44	20132120-1	TP 139	20223011U1	TP 29, TP 132
20040240RU	TP 35, TP 44	20132526	TP 105	20223020-1	TP 132
20040240US	TP 35, TP 44, TP 67	20132580	TP 100	20230141	TP 43
20040245DE	TP 19, TP 25, TP 32	20132620-1	TP 139	20230145	TP 43
20040245ES	TP 19, TP 25, TP 32	20133027	TP 105	20233070	TP 20, TP 32, TP 43
20040245FR	TP 19, TP 25, TP 32	20133028	TP 105	20290101	TP 54, TP 133
20040245IT	TP 19, TP 25, TP 32	20133101-1	TP 104, TP 139	20290120-PS	TP 133
20040245PT	TP 19, TP 25, TP 32	20133120-1	TP 139	20290131	TP 54
20040245US	TP 19, TP 25, TP 32, TP 131	20134001	TP 103, TP 139	20290132	TP 54
20040272	TP 68, TP 81	20134020	TP 139	22200077	TP 19
20040273	TP 67			22200170	TP 19
20040275	TP 134			22201011U1	TP 11, TP 130
20040280	TP 35, TP 44, TP 132			22201011U102	TP 11, TP 131

# Numerical Index

Excerpt from TELEPRESENCE

22201011U110	TP 10, TP 130	29005 HFE	TP 124, TP 127	487 O	TP 109
22201011U112	TP 10, TP 130	29005 HFH	TP 124, TP 127	487 P	TP 109
22201020-1	TP 130	29005 HFS	TP 123	487 U	TP 109
22201020-102	TP 131	29005 KH	TP 122	487 UO	TP 109
22201020-110	TP 130	29005 KHN	TP 122	487 W	TP 109
22201020-112	TP 130	29005 KKK	TP 120	487 WI	TP 109
22202011U1	TP 13, TP 131	29005 KKL	TP 120	487 Z	TP 109
22202011U110	TP 13, TP 131	29005 KKM	TP 120	495 AD	TP 107
22202020-1	TP 131	29005 KSK	TP 120	495 BD	TP 107
22202020-110	TP 131	29005 KW	TP 121	495 C	TP 108
22204011U114	TP 23, TP 131	29005 LAP	TP 114	495 D	TP 108
22204020	TP 131	29005 MH	TP 118	495 DC	TP 107
22220033-3	TP 18	29005 MHB	TP 118	495 DV	TP 107
22220054-3	TP 17	29005 MHS	TP 118	495 E	TP 108
22220055-3	TP 15	29005 MSK	TP 118	495 EW	TP 106
22220060-3	TP 16	29005 MSL	TP 118	495 F	TP 108
22220061-3	TP 16	29005 MZ	TP 118	495 FU	TP 108
22220070	TP 19	29005 MZD	TP 113	495 G	TP 108
22220071	TP 19	29005 NBB	TP 124	495 GO	TP 108
22220072	TP 19	29005 NBS	TP 124	495 KS	TP 35, TP 100
22220085-3	TP 17	29005 PA-AN	TP 120	495 M	TP 108
22220133-3	TP 18	29005 PA-GI	TP 120	495 N	TP 108
26604 AA	TP 23	29005 PA-GU	TP 120	495 NAC	TP 106
26614 AA	TP 23	29005 PA-HNO	TP 120	495 NAS	TP 107
28003 TE	TP 120	29005 PA-LAP	TP 120	495 NCSC	TP 106
28003 TO	TP 120	29005 PB	TP 119	495 NTAS	TP 107
29003 EW	TP 121	29005 PH	TP 123	495 NTW	TP 106
29003 FH	TP 122, TP 125	29005 SG	TP 121	495 NTXS	TP 106, TP 107
29003 FHG	TP 122	29005 SK	TP 123	495 NV	TP 106
29003 FL	TP 122	29005 SKN	TP 123	495 NVB	TP 106
29003 FT	TP 125	29005 SKO	TP 121	495 NVC	TP 106
29003 GK	TP 118	29005 SKS	TP 121	495 NVL	TP 106
29003 GS	TP 119	29005 SU	TP 124	495 NW	TP 106
29003 IW	TP 119	29005 SZ	TP 67, TP 118	495 NWL	TP 106
29003 KK	TP 119	29005 SZD	TP 113	495 NWM	TP 106
29003 LC	TP 125	29005 TAG	TP 121	495 NWMW	TP 107
29003 MH	TP 122	29005 TAK	TP 121	495 PA	TP 108
29003 PBK	TP 119	29005 TBG	TP 119	495 Q	TP 108
29003 S	TP 122, TP 125	29020 GM	TP 43	495 S	TP 108
29003 SV	TP 122	29020 GN	TP 43	495 SM	TP 108
29003 TBK	TP 119	29020 GO	TP 43	495 SON	TP 108
29005 AK	TP 124, TP 127	29020 GP	TP 43	495 ST	TP 108
29005 AVS	TP 123	39301 ACT	TP 20, TP 32, TP 44	495 T	TP 108
29005 BAS	TP 115	39301 PHTS	TP 20	495 TIP	TP 25, TP 106
29005 BG	TP 121	39301 Z3TS	TP 20	495 UD	TP 107
29005 BGG	TP 121	39501 TC	TP 25	495 W	TP 108
29005 BGI	TP 117	400 A	TP 81, TP 130, TP 131, TP 132, TP 133, TP 134, TP 135, TP 136, TP 137, TP 138, TP 139	495 X	TP 108
29005 BGK	TP 121	400 D	TP 120	536 ME	TP 50
29005 BK	TP 121	400 DM	TP 120	536 MF	TP 50
29005 BWK	TP 116	400 F	TP 132	536 MK	TP 50, TP 130, TP 131, TP 132, TP 133, TP 134, TP 135, TP 136
29005 DFH	TP 124, TP 127	400 FUS	TP 132	536 MKD	TP 134
29005 DRB	TP 113	40150031	TP 133	536 MN	TP 50
29005 DRS	TP 113	487 A	TP 109	547 RG	TP 50
29005 EFH	TP 124, TP 127	487 M	TP 109		
29005 GI	TP 117				
29005 GU	TP 115				

# Numerical Index

Excerpt from TELEPRESENCE

547 S	TP 50, TP 130, TP 131, TP 132, TP 133, TP 135, TP 136	69495 NE	TP 106	9524 NB	TP 73, TP 75, TP 136
547 SE	TP 50	69495 NL	TP 106	9524 NG	TP 136
547 SN	TP 50	69495 NT	TP 106	9526 NBL	TP 73, TP 75, TP 136
547 SR	TP 50	69495 NTA	TP 106	9526 NGL	TP 136
610 AFT	TP 139	81131021	TP 35, TP 43	9526 PS	TP 136
69003 FT	TP 122	8402 ZX	TP 55, TP 133	9526 SF	TP 24, TP 73, TP 81
69003 HH	TP 122	9104250	TP 132	9619 NB	TP 71, TP 134
69045001	TP 35, TP 132	9370 TQ	TP 67	96206173E	TP 137
69045002	TP 35, TP 132	9372 LA	TP 67	9626 NB	TP 71, TP 134
69045020	TP 132	9415 NG	TP 136	9626 NB-2	TP 71, TP 134
69111001	TP 87, TP 137	9415 NSF	TP 81	9626 P	TP 134
69111001 C	TP 87, TP 137	9419 NG	TP 136	9626 PS	TP 134
69111020	TP 137	9419 NSF	TP 81	9626 SF	TP 71
69131501	TP 101, TP 139	9419 PS	TP 81, TP 135, TP 136	9726 NB	TP 73, TP 75, TP 136
69132601-1	TP 100, TP 139	9500 RXG	TP 137	9726 NG	TP 136
69215034	TP 41, TP 42, TP 132	9500 TXG	TP 137	9800 C	TP 25
69215134	TP 41, TP 42, TP 132	9500 ZW	TP 79, TP 137	9800 G	TP 25, TP 131
69236001	TP 40, TP 132	9500 ZW-1	TP 79	9824 NB-3D	TP 24
69236011	TP 41, TP 133	9500 ZW-2	TP 79, TP 137	9832 A	TP 24
69236020	TP 132, TP 133	9500 ZW-4	TP 79	9832 NB-3D	TP 24
69236101	TP 40, TP 132	9500 ZW-5	TP 79	9832 SF	TP 24
69274 N	TP 80	9500 ZW-6	TP 79, TP 137	ET38-1717618	TP 133
69487 KVE	TP 109	9500 ZW-7	TP 79, TP 137	ET38-1717707	TP 133
69495 NA	TP 106	9500 ZW-8	TP 79, TP 137	ET38-1717715	TP 133
69495 NB	TP 106	9500 ZW-9	TP 79	ET38-1800496	TP 133
69495 NCS	TP 106	9515 NB	TP 73, TP 75, TP 136	W20475	TP 68
69495 ND	TP 106	9519 NB	TP 73, TP 75, TP 136	W29519 N	TP 63, TP 64
		9519 NG	TP 135	W29519 NB	TP 63, TP 64, TP 135
		9523 PS	TP 81		



**IMAGE 1 HUB™ HD**  
**Brilliant Imaging in FULL HD**



2-12



## Brilliant Imaging in FULL HD

IMAGE1  
HD 1

CAMERA SYSTEMS

The 16:9 three-chip CCD system in the IMAGE1 FULL HD camera head captures images that can be displayed in natural format without being converted. The implemented progressive scan technology, therefore, enables natural and stress-free viewing of the endoscopic image.

To provide veterinarians with the ideal conditions for diagnosis, great emphasis was placed on both image quality and an optimal signal-to-noise ratio when designing the IMAGE1 camera head family. A low signal-to-noise ratio is guaranteed even with high light sensitivity, resulting in bright, crystal-clear and ultra-sharp images. The unique combination of easy handling and excellent image quality ensures time efficiency and high patient safety. Consequently, IMAGE1 FULL HD camera heads represent a further milestone in the evolution of endoscopic imaging.

- With a resolution of 1920 x 1080 pixels, IMAGE1 FULL HD camera heads capture ultra-sharp imaging with unsurpassed brilliance.
- The integrated Parfocal Zoom Lens ensures optical magnification and excellent rendition of details without information loss. Digital image processing delivers excellent image quality, even under difficult viewing conditions.
- Also available:  
Autoclavable FULL HD camera heads as well as FULL HD camera heads for perfusion diagnosis of tissue and organs with indocyanine green (ICG)

IMAGE1 FULL HD camera heads enable direct control of KARL STORZ units with the KARL STORZ Communication Bus (SCB). The camera, therefore, is the cornerstone of the endoscopic operating environment.



# KARL STORZ FULL HD

Highest Standards for Image Input and Image Display



Maximum resolution and the consistent use of the 16:9 aspect ratio guarantee FULL HD.

Endoscopic and microscopic camera systems are equipped with three CCD sensors which support the 16:9 input format as well as capturing images with a resolution of 1920 x 1080 pixels.

These factors provide users with significantly better image definition and a correspondingly larger display window.

The benefits of High Definition Technology (HD) for medical applications are:

- Up to 6 times\* higher input resolution of the camera delivers more detail and depth of focus.
- Using 16:9 format during image acquisition enlarges the field of vision and supports ergonomic viewing.
- The brilliant color rendition enables optimal diagnosis.
- Lateral view is enhanced by 32% when the endoscope is withdrawn slightly, providing the same image enhancement as a standard system. Any vertical information loss is restored and the lens remains clean.

The KARL STORZ FULL HD full-screen standard has the following resolution:

1080p (progressive); 1920 x 1080 pixels, 16:9 format  
**16:9 format for input and output!**

Further formats available on the market with lower resolution or half-frames:


720p (progressive); 1280 x 720 pixels, 16:9 format  
1080i (interlaced); 1920 x 1080 pixels, 16:9 format

### The KARL STORZ standard for endoscopic imaging

In contrast to the frequently used interlaced method, the KARL STORZ 1080p HD progressive scan system delivers an extremely stable image, with no flickering or interference.

For medical users this improved image quality virtually eliminates the eye fatigue which can result from working with a monitor.

Image acquisition is a determining factor for resolution and a larger display window.




### Image rendering in other systems

Interlaced scan (half-frame) displays 25/30\* "half-frames" per second.



### KARL STORZ image rendering

Progressive scan (full frame) displays 50/60\* "full frames" per second.



\* depending on the color system (PAL/NTSC) used

# IMAGE 1 HUB™ HD

FULL HD Camera Head, FULL HD Camera Control Unit, FULL HD Monitor

IMAGE  
HD 1

CAMERA SYSTEMS

## Ready for the future

Continuous improvement in standardization is a cornerstone of KARL STORZ camera systems. All IMAGE 1 HUB™ HD products are compatible with previous products, ensuring upward and downward compatibility and securing your investment in the future.

## Image acquisition with KARL STORZ IMAGE1 HD

- The High Definition three-chip CCD in the IMAGE1 FULL HD camera head captures images in 16:9 format that can be displayed on WIDEVIEW monitors in the same format without being converted.
- Digital input guarantees the best possible image quality. This eliminates the need to convert analog input signals to digital output formats which can result in a loss of image quality.

## IMAGE 1 HUB™ HD Camera Control Unit

- The update for the IMAGE1 CCU, combined with a new front panel and an extension of the plug-in slots on the back of the unit, enable all IMAGE1 users to experience the advantages of FULL HD.
- The full range of IMAGE1 functions are also available in the FULL HD system.
- **IMAGE 1 HUB™ HD with backward compatibility**  
All existing systems can continue to be used after upgrading the IMAGE1 system to FULL HD. Equipment modules and adaptations to existing infrastructures remain completely compatible.

## FULL HD Monitor

### Image format and image rendering determine quality

- Using 16:9 aspect ratio in place of the conventional 4:3 format provides a larger display window and demonstrably improves the viewing ergonomics.
- KARL STORZ FULL HD monitor produces the highest possible resolution of 1920 x 1080 pixels, providing medical users with improved depth of focus and color contrast for enhanced image quality.



IMAGE1 FULL HD Camera Head H3-Z with 16:9 three-CCD chip



IMAGE 1 HUB™ HD Camera Control Unit – Update of IMAGE1 CCU



KARL STORZ FULL HD Monitor

# IMAGE 1 HUB™ HD

FULL HD Camera Control Unit

**STORZ**  
KARL STORZ — ENDOSKOPE

CAMERA SYSTEMS

IMAGE  
HD 1



## 222010 11U110 IMAGE 1 HUB™ HD Camera Control Unit SCB, with ICM module

for use with IMAGE1 FULL HD and IMAGE1 standard one- and three-chip camera heads, max. resolution 1920 x 1080 pixels, **with integrated ICM (Image Capture Module), KARL STORZ SCB** and digital Image Processing Module, color systems PAL/NTSC, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**BNC/BNC Cable**

**S-Video (Y/C) Connecting Cable**

**Special RGBS Connecting Cable**

2x **Connecting Cable**, for controlling peripheral units

**DVI-D Connecting Cable**

**SCB Connecting Cable**

**Keyboard**, with US English character set

2x **KARL STORZ USB Flash Drive**, 4 GB

## 222010 11U112 IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI/ICM module

for use with IMAGE1 FULL HD and IMAGE1 standard one- and three-chip camera heads, max. resolution 1920 x 1080 pixels, **with integrated ICM (Image Capture Module) and SDI (Serial Digital Interface) modules, KARL STORZ-SCB** and digital Image Processing Module, color systems PAL/NTSC, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**BNC/BNC Video Cable**

**S-Video (Y/C) Connecting Cable**

**Special RGBS Connecting Cable**

2x **Connecting Cable**, for controlling peripheral units

**DVI-D Connecting Cable**

**SCB Connecting Cable**

**Keyboard**, with US English character set

2x **KARL STORZ USB Flash Drive**, 4 GB

**Further accessories for IMAGE 1 HUB™ HD** see pages TP 19-20

**Component Parts** see COMPONENTS, SPARE PARTS

2-12

# IMAGE 1 HUB™ HD

FULL HD Camera Control Unit

IMAGE  
HD 1



## 222010 11U1 IMAGE 1 HUB™ HD Camera Control Unit SCB

for use with IMAGE1 FULL HD and IMAGE1 standard one- and three-chip camera heads, max. resolution 1920 x 1080 pixels, with integrated **KARL STORZ-SCB** and integrated digital Image Processing Module, color systems PAL/NTSC, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**BNC/BNC Video Cable**

**S-Video (Y/C) Connecting Cable**

**Special RGBS Connecting Cable**

2x **Connecting Cable**, for controlling peripheral units

**DVI-D Connecting Cable**

**SCB Connecting Cable**

**Keyboard**, with US English character set

## 222010 11U102 IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI Module

for use with IMAGE1 FULL HD and IMAGE1 standard one- and three-chip camera heads, max. resolution 1920 x 1080 pixels, **with integrated SDI (Serial Digital Interface) module, KARL STORZ-SCB** and integrated Image Processing Module, color systems PAL/NTSC, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**BNC/BNC Video Cable**

**S-Video (Y/C) Connecting Cable**

**Special RGBS Connecting Cable**

2x **Connecting Cables**, for controlling peripheral units

**DVI-D Connecting Cable**

**SCB Connecting Cable**

**Keyboard**, with US English character set

**Further accessories for IMAGE 1 HUB™ HD** see pages TP 19-20

**Component Parts** see COMPONENTS, SPARE PARTS

# IMAGE 1 HUB™ HD

FULL HD Camera Control Unit



### Specifications:


Signal-to-Noise Ratio	IMAGE 1 HUB™ HD, three-chip camera systems $\geq 60$ dB
AGC	Microprocessor-controlled
Video Output	- FULL HD signal to DVI-D socket (2x) - SDI signal to BNC socket (only IMAGE 1 HUB™ HD with SDI module) (2x) - RGBS signal to D-Sub socket - S-Video an 4-pol. Mini-DIN-Buchse (2x) - Composite signal to BNC socket
Input	Keyboard for title generator, 5-pin DIN socket

Control Output /Input	- KARL STORZ-SCB at 6-pin Mini-DIN socket (2x) - 3.5 mm stereo jack plug (ACC 1, ACC 2), - Serial port at RJ-11 - USB port (only IMAGE 1 HUB™ HD with ICM) (2x)
Dimensions w x h x d	305 x 89 x 335 mm
Weight	3.35 kg
Power Supply	100-240 VAC, 50/60 Hz
Certified to	IEC 601-1, 601-2-18, CSA 22.2 No. 601, UL 2601-1 and CE acc. to MDD, protection class 1/CF defibrillation-safe

**ICM:** 2x USB outputs (1x printer output on rear panel, 1x storage medium on the front panel)

**SDI:** 2x outputs

### Additional Modules:



**SDI: Serial Digital Interface**  
optimized for medical display on flat screens, routing with OR1™ and digital recording with AIDA-DVD-M

**ICM: ICM**  
USB connection for recording still images and video sequences on a USB storage medium or to connect to a USB printer for printing recorded still images directly.  
Compatible printers:  
KARL STORZ USB Printer 549 M

### Recommended storage media:

KARL STORZ 4 GB USB Flash Drive **20** 0402 81

Further accessories for IMAGE 1 HUB™ HD see pages TP 19-20





**22 2020 11U110 IMAGE1 HD Camera Control Unit SCB,  
with ICM module**

for use with IMAGE1 FULL HD three-chip camera heads, max. resolution 1920 x 1080 pixels, **with integrated ICM (Image Capture Module), KARL STORZ-SCB** and digital Image Processing Module, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

2x **Connecting Cable**, for controlling peripheral units

**DVI-D Connecting Cable**

**SCB Connecting Cable**

**Keyboard**, with US English character set

2x **KARL STORZ USB Flash Drive**, 4 GB

**22 2020 11U1 IMAGE1 HD Camera Control Unit SCB**

for use with IMAGE1 FULL HD three-chip camera heads, max. resolution 1920 x 1080 pixels, with integrated **KARL STORZ-SCB** and digital Image Processing Modul, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

2x **SCB Connecting Cable**

**DVI-D Connecting Cable**

**Connecting Cable**, for controlling peripheral units

**Keyboard**, with US English character set

**Further accessories for IMAGE1 HD** see pages TP 19-20

**Component Parts** see COMPONENTS, SPARE PARTS

# IMAGE1 HD <sup>NEW</sup>

FULL HD Camera Control Unit



**Specifications:**

Signal-to-noise ratio	IMAGE1 HD, three-chip camera systems $\geq 60$ dB
AGC	Microprocessor-controlled
Video output	FULL HD signal to DVI-D socket
Input	Keyboard for title generator, 5-pin DIN socket
Control output /input	- USB Port (only IMAGE1 HD with ICM) (2x) - Serial port at RJ-11 - 3.5 mm stereo jack plug (ACC 1, ACC 2) - KARL STORZ-SCB at 6-pin Mini-DIN socket (2x)

Dimensions w x h x d	305 x 89 x 335 mm
Weight	3.35 kg
Power supply	100-240 VAC, 50/60 Hz
Certified to	IEC 601-1, 601-2-18, CSA 22.2 No. 601, UL 2601-1 and CE acc. to MDD, protection class 1/CF defibrillation-safe

**ICM:** 2x USB outputs (1x printer output on rear panel, 1x storage medium on the front panel)

**Additional Modules :**



**ICM:**

**ICM**

USB connection for recording still images and video sequences on a USB storage medium or to connect to a USB printer for printing recorded still images directly.

Compatible printers:

KARL STORZ USB Printer 549 M

**Recommended storage media:**

KARL STORZ 4 GB USB Flash Drive 20 0402 81

Further accessories for IMAGE1 HD see pages TP 19-20

# IMAGE1 FULL HD Camera Heads <sup>NEW</sup>

for Endoscopy



For use with IMAGE 1 HUB™ HD Camera Control Unit SCB 22 2010 11U1xx  
and IMAGE1 HD Camera Control Unit SCB 22 2020 11U1xx



22 2200 55-3

22 2200 55-3	50 Hz 60 Hz	<b>IMAGE1 H3-Z</b> <b>Three-Chip FULL HD Camera Head</b>
--------------	----------------	---

max. resolution 1920 x 1080 pixels, progressive scan, soakable, gas- and plasma-sterilizable, with integrated Parfocal Zoom Lens, focal length  $f = 15 - 31$  mm (2x), 2 freely programmable camera head buttons, for use with color systems **PAL/NTSC**

#### Specifications:

IMAGE1 FULL HD Camera Head	H3-Z
50/60 Hz	22 2200 55-3
Image sensor	3x 1/8" CCD chip
Pixel output signal H x V	1920 x 1080
Dimensions (w x h x d)	39 x 49 x 114 mm
Weight	270 g
Optical interface	integrated Parfocal Zoom Lens, $f = 15-31$ mm (2x)
Min. sensitivity	F 1.4/1.17 Lux
Grip mechanism	standard eyepiece adaptor
Cable	non-detachable
Cable length	300 cm

Further accessories for IMAGE1 FULL HD see pages TP 19-20

# IMAGE1 FULL HD Camera Heads <sup>NEW</sup>

for Endoscopy



For use with IMAGE 1 HUB™ HD Camera Control Unit SCB 22 2010 11U1xx and IMAGE1 HD Camera Control Unit SCB 22 2020 11U1xx



222200 61-3

**222200 61-3**    **50 Hz**    **IMAGE1 H3-ZA**  
**60 Hz**    **Three-Chip FULL HD Camera Head**

**autoclavable**, max. resolution 1920 x 1080 pixels, progressive scan, soakable, gas- and plasma-sterilizable, with integrated Parfocal Zoom Lens, focal length  $f = 15 - 31$  mm (2x), 2 freely programmable camera head buttons, for use with color systems **PAL/NTSC**



222200 60-3

**222200 60-3**    **50 Hz**    **IMAGE1 H3-FA**  
**60 Hz**    **Three-Chip FULL HD Camera Head**

**autoclavable**, max. resolution 1920 x 1080 pixels, progressive scan, soakable, gas- and plasma-sterilizable, **fixed focus**, focal length  $f = 17$  mm, 2 freely programmable camera head buttons, for use with color systems **PAL/NTSC**

**Specifications:**

IMAGE1 FULL HD Camera Heads	H3-ZA	H3-FA
50/60 Hz	222200 61-3 (PAL/NTSC)	222200 60-3 (PAL/NTSC)
Image sensor	3x 1/8" CCD chip	3x 1/8" CCD chip
Pixel Output signal H x V	1920 x 1080	1920 x 1080
Dimensions (w x h x d)	39 x 49 x 100 mm	39 x 49 x 93 mm
Weight	299 g	261 g
Optical interface	Integrated Parfocal Zoom Lens, $f = 15-31$ mm	fixed focus $f = 17$ mm
Min. sensitivity	F 1.4/1.17 Lux	F 1.4/1.17 Lux
Grip mechanism	Standard eyepiece adaptor	standard eyepiece adaptor
Cable	non-detachable	non-detachable
Cable length	300 cm	300 cm

Further accessories for IMAGE1 FULL HD see pages TP 19-20

# IMAGE1 FULL HD Camera Heads

NEW

for Microscopy and ICG Application



For use with IMAGE 1 HUB™ HD Camera Control Unit SCB 22 2010 11U1xx and IMAGE1 HD Camera Control Unit SCB 22 2020 11U1xx



22 2200 54-3

**22 2200 54-3**    **50 Hz**    **IMAGE1 H3-M COVIEW® Three-Chip**  
**60 Hz**    **FULL HD Microscope Camera Head**

max. resolution 1920 x 1080 pixels, progressive scan, with C-MOUNT thread for coupling to microscopes, 2 freely programmable camera head buttons, with detachable camera head cable, length 900 cm, for use with color systems **PAL/NTSC**



22 2200 85-3

**22 2200 85-3**    **50 Hz**    **IMAGE1 H3-Z FI Three-Chip ICG**  
**60 Hz**    **FULL HD Camera Head**

for perfusion diagnosis of tissue and organs with indocyanine green (ICG) in conjunction with light source D-LIGHT P, max. resolution 1920 x 1080 pixels, progressive scan, with integrated Parfocal Zoom Lens, focal length f = 15 – 31 mm (2x), with 2 freely programmable camera head buttons, for use with color systems **PAL/NTSC**

**Specifications:**

IMAGE1 FULL HD Camera Heads	H3-M COVIEW®	H3-Z FI
50/60 Hz	22 2200 54-3 (PAL/NTSC)	22 2200 85-3 (PAL/NTSC)
Image sensor	3x 1/8" CCD chip	3x 1/8" CCD chip
Pixel Output signal H x V	1920 x 1080	1920 x 1080
Dimensions (w x h x d)	45 x 50 x 60 mm	39 x 49 x 114 mm
Weight	240 g	246 g
Optical interface	C-MOUNT connection	Integrated Parfocal Zoom Lens, f = 15-31 mm (2x)
Min. sensitivity	F 1.9/1.4 Lux	F 1.4/1.17 Lux
Grip mechanism	C-MOUNT connection	standard eyepiece adaptor
Cable	detachable	non-detachable
Cable length	900 cm	300 cm

**Components for connecting the H3-M COVIEW® camera head to operating microscopes**  
see complete TELEPRESENCE catalog

**Further accessories for IMAGE1 FULL HD** see pages TP 19-20

# IMAGE1 Camera Heads

For use with IMAGE 1 HUB™ HD Camera Control Unit SCB 22 2010 11U1xx



22 2200 33-3 / 22 2201 33-3

22 2200 33-3	<b>PAL</b>	<b>IMAGE1 F3 Three-Chip ENT and Arthroscopy Camera Head</b>
22 2201 33-3	<b>NTSC</b>	

color systems **PAL/NTSC**, soakable, gas-sterilizable, fixed focus  
f = 16.8 mm, 2 freely programmable camera head buttons

Further accessories for **IMAGE1** see pages TP 19-20

# IMAGE 1 HUB™ HD, IMAGE1

Accessories



- 20 2002 31** **Keyboard**, with German character set, PS/2 connector, including AT-adaptor and protective cover
- 20 2002 31F** **Same**, with French character set
- 20 2002 31I** **Same**, with Italian character set
- 20 2002 31P** **Same**, with Portuguese character set
- 20 2002 31S** **Same**, with Spanish character set
- 20 2002 31U** **Same**, with US English character set



- 20 0402 45DE** **Silicone Keyboard**, with German character set, DIN AT connector, IP 68
- 20 0402 45FR** **Same**, with French character set
- 20 0402 45IT** **Same**, with Italian character set
- 20 0402 45PT** **Same**, with Portuguese character set
- 20 0402 45ES** **Same**, with Spanish character set
- 20 0402 45US** **Same**, with US English character set



- 22 2200 72** **IMAGE 1 HUB™ HD Extension Cable**, length 520 cm, for use with IMAGE1 HD camera heads (except H3-M/H3-M COVIEW® FULL HD microscope camera heads)



- 22 2001 70** **IMAGE1 SD Extension Cable**, length 760 cm, for use with IMAGE1 standard camera heads



- 22 2200 70** **IMAGE1 HD Camera Head Cable**, length 300 cm, for use with H3-M/H3-M COVIEW® FULL HD Microscope Camera Heads **22 2200 54-3** and **22 2201 54-3** (spare part)



- 22 2200 71** **IMAGE1 HD Camera Head Cable**, length 900 cm, for use with H3-M/H3-M COVIEW® FULL HD Microscope Camera Heads **22 2200 54-3** and **22 2201 54-3** (spare part)



- 22 2000 77** **Video Endoscope Adaptor**, color systems **PAL/NTSC**, length 90 cm, for use with GI CCU and IMAGE1

# IMAGE 1 HUB™ HD, IMAGE1

Accessories



39301 Z3TS **Plastic Container for Sterilization and Storage** of camera heads **IMAGE1 H3-Z, H3-ZA and H3-FA, autoclavable**, suitable for use with steam, gas and hydrogen peroxide sterilization, Sterrad® compatible, external dimensions (w x d x h): 385 x 255 x 75 mm

**Please note:** The instrument displayed is not included in the plastic container.

Only camera heads marked “autoclave” can be placed in the tray for steam sterilization.



**NEW** 39301 PHTS **Plastic Container for Sterilization and Storage** of camera heads **IMAGE1 H3-P and H3-ZI, autoclavable**, suitable for use with steam, gas and hydrogen peroxide sterilization, Sterrad® compatible, external dimensions (w x d x h): 385 x 255 x 75 mm

**Please note:** The instrument displayed is not included in the plastic container.

Only camera heads marked “autoclave” can be placed in the tray for steam sterilization.



39301 ACT **Plastic Container for Sterilization and Storage**, of camera heads **IMAGE1 P3, TRICAM® and TELECAM, autoclavable**, for use with steam, gas and plasma sterilization

**Please note:** The instrument displayed is not included in the plastic container.

Only camera heads marked “autoclave” can be placed in the tray for steam sterilization.

(not usable with IMAGE1 FULL HD camera heads)



20 0141 30 **One-Pedal Footswitch**, one-stage, for freeze frame or video recording in combination with IMAGE 1 HUB™ 22 2010 20-110/-112 and 22 2000 20-112

20 2330 70 **Adaptor ACC 1**, for Footswitch 20 0141 30




**NEW** 20 0402 81 **KARL STORZ USB Flash Drive, 4 GB**



KARL STORZ 3D System *NEW*

**STORZ**  
KARL STORZ — ENDOSKOPE

CAMERA SYSTEMS



2-12

# The KARL STORZ 3D System <sup>NEW</sup>

Discover a New World in 3D

Exact depth perception inside the patient's body is essential for any endoscopic procedure. The KARL STORZ 3D system offers the surgeon excellent depth of field combined with outstanding hand-eye coordination. This first-class stereoscopic imaging system is particularly valuable for activities such as suturing and knotting that demand a high degree of spatial perception. The KARL STORZ 3D system ensures accurate positioning with minimal fatigue and, therefore, leads to significant improvements in surgical efficiency and patient safety in the OR.

The 3D TIPCAM®1 design is based on the 0° endoscope with mounted camera head. The system can be used like a normal endoscope, which simplifies work for the user. Light and data cables are placed at the proximal end of the endoscope to provide the surgeon with additional flexibility.

Its slender and lightweight design makes the 3D endoscope an excellent instrument. Combined with an ergonomic handle, the 3D TIPCAM®1 with integrated endoscope opens up a whole new dimension in stereoscopic imaging without compromising ease of use.

Two image sensors at the distal end of the 3D TIPCAM®1 are precisely aligned mechanically. Complemented by precise electronics in the endoscope and camera control unit, the system generates endoscopic 3D images that are true to life. Combined with the latest generation of 3D monitors, this system offers the surgeon a complete system that reduces fatigue to a minimum even after long hours of surgery.



KARL STORZ 3D glasses



distal sensor chips



3D TIPCAM®1

# KARL STORZ 3D System <sup>NEW</sup>

3D Camera Control Unit, 3D TIPCAM®1



## 222040 11U114 3D Camera Control Unit SCB, with ICM Module

max. resolution 1920 x 1080 pixels, with integrated ICM (Image Capture Module) and KARL STORZ-SCB, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**BNC/BNC Video Cable**

**S-Video (Y/C) Connecting Cable**

2x **Connecting Cable**, for controlling peripheral units

**DVI-D Connecting Cable**

**SCB Connecting Cable**

**Silicone Keyboard**, with US English character set

2x **KARL STORZ Flash Drive**, 4 GB

8x **KARL STORZ 3D Glasses**

For use with KARL STORZ 3D Camera Control Unit SCB 222040 11-114



26604/26614 AA

26604 AA	<b>PAL</b>	<b>3D TIPCAM®1</b>
26614 AA	<b>NTSC</b>	

with two distal CCD image sensors, direction of view 0°, diameter 10 mm, length 31 cm, for use with KARL STORZ 3D Camera Control Unit SCB 222040 11U114

Further accessories for KARL STORZ 3D System see page TP 25

Component Parts see COMPONENTS, SPARE PARTS



9832 NB-3D

9832 NB-3D **32" 3D Monitor**

wall-mounted with VESA 200 adaption, with integrated power supply, color systems **PAL/NTSC**, max. screen resolution 1920 x 1080, image format 16:9, power supply 100 – 240 VAC, 50/60 Hz



9824 NB-3D

9824 NB-3D **24" 3D Monitor**

wall-mounted with VESA 100 adaption, color systems **PAL/NTSC**, max. screen resolution 1920 x 1200, image format 16:10, power supply 100 – 240 VAC, 50/60 Hz


9832 SF **Pedestal**, for 3D Monitor 9832 NB-3D

9526 SF **Pedestal**, for 95xx and 97xx series and 3D Monitor 9824 NB-3D

9832 A **Adaptor Plate**, VESA 200 to VESA 100, for 3D Monitor 9832 NB-3D, for fixation to the VESA 100 holding systems



- 20040245DE **Silicone Keyboard**, with German character set, DIN AT connector, IP 68
- 20040245FR **Same**, with French character set
- 20040245IT **Same**, with Italian character set
- 20040245PT **Same**, with Portuguese character set
- 20040245ES **Same**, with Spanish character set
- 20040245US **Same**, with US English character set



- 39501 TC **Wire Tray for Cleaning, Sterilization and Storage** of 3D TIPCAM®1 266x4 AA, **autoclavable**, incl. compartment for light cable, external dimensions (w x d x h): 480 x 320 x 75 mm



- 20 0402 81 **KARL STORZ USB Flash Drive**, 4 GB



- 9800 G **3D Glasses**, passive, circularly polarized, package of 8



- 9800 C **3D Clip-on Glasses**, circularly polarized



- 495 TIP **Fiber Optic Light Cable**, highly heat resistant, diameter 4.8 mm, length 300 cm, optimized for 3D TIPCAM®1 and TIPCAM®



- 20 0200 81 **TROLL-E Mobile Stand**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, for mounting TFT monitors with VESA 75/100 connection, integrated cable conduit in vertical beam, 1 shelf and cable manager, load capacity for monitor: max. 15 kg, Dimensions:  
Mobile stand: 670 x 1660 x 670 mm (w x h x d),  
Shelf: 455 x 350 mm (w x d), caster diameter: 100 mm  
Trolley is delivered unassembled.




# TRICAM® – Three-Chip Camera System

**STORZ**  
KARL STORZ – ENDOSKOPE



CAMERA SYSTEMS



2-12

# TRICAM® SL II

Three-Chip Camera System, Color Systems PAL/NTSC,  
Digital Imaging Processor

**TRICAM® SL II** is the endoscopic video platform offered by KARL STORZ for all endoscopically-assisted applications and video display of the highest image quality.

Thanks to its ergonomic, compact design, this CCD sensor technology separates the optical image of the 3-CCD image sensor into the three primary colors by means of a specially coated color prism.

Each individual CCD sensor receives a specific color with the maximum line representation.

The separate processing of the primary color information red, green, and blue means that all colors can be optimally captured with this process. Furthermore, the use of all three CCD sensors offers up to three times more detail enhancement.

#### Special Features:


- **Best image quality and display due to 3-CCD sensor chip technology**
- **Enhanced signal quality through Digital Signal Processing (contrast and filter enhancement for rigid, semirigid and flexible endoscopes)**
- **Optimal image size display, regardless of the telescope used, thanks to patented Parfocal Optical Zoom Lens**
- **Simplified operation reduces camera movement to, e.g., prevent damage to or soiling of the lens**
- **Consistent illumination due to the auto-exposure system ensures seamless transitions between images, regardless of imaging area or available light**
- **Functionality at the touch of a button with easy-to-program operating elements via intuitive control of all functions in the sterile area**
- **Parameters of the OR equipment used, light sources, etc. can be displayed and controlled (KARL STORZ Communication Bus) on the monitor**
- **Digital video output (DV) for digital recording**
- **Digital zoom 1 – 2x can be adjusted to 4 levels**
- **“Freeze Frame” function can be activated via camera head buttons or an external footswitch**
- **Fluorescence imaging mode available**

**TRICAM® SL II** offers a high degree of optoelectronic accuracy and precision.

Signals are processed by means of a digital image processor. The endoscopic image is produced true to the original image. The image processor also assumes automatic image optimization during continuously changing illumination and color conditions. Manual adjustment is therefore no longer necessary.

The user has several pre-assigned digital-optimization-filters at his disposal for adjustment to rigid, semi-rigid and flexible endoscopes at the touch of a button.

Video outputs such as Digital-Video (DV) and all analog signal variations are an integral feature of **TRICAM® SL II** standard equipment.





# TRICAM® SL II

Camera Control Unit SCB



## 20 2230 11U1 TRICAM® SL II Camera Control Unit SCB

color systems **PAL/NTSC**, with **KARL STORZ-SCB** and integrated digital Image Processing Module, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**Keyboard**, with US English character set

2x **Connecting Cable**, for controlling peripheral units



**BNC/BNC Video Cable**

**S-Video (Y/C) Connecting Cable**

**Special RGBS Connecting Cable**

**SCB Connecting Cable**

**DV Connecting Cable**, 6 pin to 4 pin



*Photodynamic Diagnosis (PDD) with D-LIGHT light sources, left: white light, right: blue light*

### Specifications:

Video output	- RGBS signal to 4x BNC sockets - S-Video signal to 4 pin Mini-DIN socket (2x) - Composite signal at BNC socket - DV signal to 6 pin DV socket	Camera Control Unit (CCU)	dimensions: 305 x 88 x 254 mm (w x h x d)
Input	keyboard input for character generator and control functions, 6-pin Mini-DIN socket	Weight	2.7 kg
Control out/input	KARL STORZ-SCB to 6 pin Mini-DIN socket (2x), 3.5 mm stereo jack plug (ACC 1, ACC 2)	Power supply	100-240 VAC, 50/60 Hz
		Certified to	IEC 601-1, 601-2-18, CSA 22.2 No. 601, UL 2601, and CE acc. to MDD, protection class 1/BF

**Further accessories for TRICAM® SL II** see page TP 32

**Component Parts** see COMPONENTS, SPARE PARTS

# TRICAM® SL II

Camera Heads

For use with TRICAM® SL II Camera Control Unit SCB 20 2230 11U1



20 2210 40 / 20 2211 40

20 2210 40	<b>PAL</b>	<b>TRICAM®</b>
20 2211 40	<b>NTSC</b>	<b>Three-Chip Camera Head</b>

color systems **PAL/NTSC**, autoclavable, soakable, gas-sterilizable, with integrated Parfocal Zoom Lens, f = 14 – 28 mm (2x), 2 freely programmable camera head buttons, including Plastic Container for Sterilization and Storage 39301 ACT



20 2210 30 / 20 2211 30

20 2210 30	<b>PAL</b>	<b>TRICAM®</b>
20 2211 30	<b>NTSC</b>	<b>Three-Chip Camera Head</b>

color systems **PAL/NTSC**, soakable, gas-sterilizable, integrated Parfocal Zoom Lens, f = 25 – 50 mm (2x), 2 freely programmable camera head buttons

**Note:** All camera heads can be disinfected by soaking in suitable solution, or sterilized in gas. The use of a sterile cover is recommended (see page TP 51). In addition, camera heads **20 2210 40** and **20 2211 40** are autoclavable up to 134 °C.

**Further accessories for TRICAM® SL II** see page TP 32

# TRICAM® SL II

Camera Heads

For use with TRICAM® SL II Camera Control Unit SCB 20 2230 11U1



20 2210 39 / 20 2211 39

20 2210 39	<b>PAL</b>	<b>TRICAM® PDD Three-Chip</b>
20 2211 39	<b>NTSC</b>	<b>Pendulum Camera Head</b>

color systems **PAL/NTSC**, soakable, gas-sterilizable, **with pendulum system and fixed focus, f = 14 mm**, 2 freely programmable camera head buttons, **for photodynamic early diagnosis**



20 2210 37 / 20 2211 37

20 2210 37	<b>PAL</b>	<b>TRICAM® PDD</b>
20 2211 37	<b>NTSC</b>	<b>Three-Chip Camera Head</b>

color systems **PAL/NTSC**, soakable, gas-sterilizable, with integrated Parfocal Zoom Lens, f = 25 – 50 mm (2x), 2 freely programmable camera head buttons, **for photodynamic early diagnosis**

**Specifications:**

TRICAM® Camera Heads	Three-Chip Camera Head, autoclavable	Three-Chip Camera Head	PDD Three-Chip Pendulum Camera Head	PDD Three-Chip Camera Head
<b>PAL</b>	20 2210 40	20 2210 30	20 2210 39	20 2210 37
<b>NTSC</b>	20 2211 40	20 2211 30	20 2211 39	20 2211 37
Image sensor	3x ¼" CCD chip	3x ½" CCD chip	3x ¼" CCD chip	3x ½" CCD chip
Pixels output signal H x V	752 x 582 (PAL) 768 x 494 (NTSC)	752 x 582 (PAL) 768 x 494 (NTSC)	752 x 582 (PAL) 768 x 494 (NTSC)	752 x 582 (PAL) 768 x 494 (NTSC)
Resolution (horizontal)	min. 750 lines	min. 750 lines	min. 750 lines	min. 750 lines
Dimensions	diameter 28-33 mm, length 97 mm	37 x 47 mm, length 124 mm	diameter 29-37 mm, length 117 mm	37 x 47 mm, length 124 mm
Weight	195 g	240 g	170 g	240 g
Min. sensitivity	F 1.4/3 Lux	F 1.4/3 Lux	F 1.4/3 Lux (standard)	F 1.4/3 Lux (standard)
Lens	integrated Parfocal Zoom Lens, f = 14-28 mm	integrated Parfocal Zoom Lens, f = 25-50 mm	pendulum System, fixed Focus, f = 14 mm	integrated Parfocal Zoom Lens, f = 25-50 mm
Grip mechanism	standard eyepiece adaptor	standard eyepiece adaptor	standard eyepiece adaptor	standard eyepiece adaptor

**Further accessories for TRICAM® SL II** see page TP 32

# TRICAM® SL II

Optional Accessories



- 20 2002 31 **Keyboard**, with German character set, PS/2 connector, including AT-adaptor and protective cover
- 20 2002 31F **Same**, with French character set
- 20 2002 31I **Same**, with Italian character set
- 20 2002 31P **Same**, with Portuguese character set
- 20 2002 31S **Same**, with Spanish character set
- 20 2002 31U **Same**, with US English character set



- 20 0402 45DE **Silicone Keyboard**, with German character set, DIN AT connector, IP 68
- 20 0402 45FR **Same**, with French character set
- 20 0402 45IT **Same**, with Italian character set
- 20 0402 45PT **Same**, with Portuguese character set
- 20 0402 45ES **Same**, with Spanish character set
- 20 0402 45US **Same**, with US English character set



- 20 2202 70 **TRICAM® SL II Extension Cable**, length 600 cm



- 39301 ACT **Plastic Container for Sterilization and Storage**, of camera heads IMAGE1 P3, TRICAM® and TELECAM, **autoclavable**, for use with steam, gas and plasma sterilization

**Please note:** The instrument displayed is not included in the plastic container.  
Only camera heads marked “autoclave” can be placed in the tray for steam sterilization.  
(except IMAGE1 FULL HD camera heads)



- 20 2330 70 **Adaptor ACC 1**, for Footswitch 20 0141 30, for freeze function with TELECAM DX II, TELECAM SL II and TRICAM® SL II
- 20 0141 30 **One-Pedal Footswitch**, digital, one-stage

# TELE PACK VET X Compact System and TELECAM One-Chip Camera System



2-12

# TELE PACK VET X <sup>NEW</sup>

A multipurpose unit for diagnostic and surgical procedures including still image and video capture software.

**STORZ**  
KARL STORZ — ENDOSKOPE

TELE PACK VET X is a compact, portable and flexible system that has been developed for use in a large number of fields. It can be used in veterinary hospitals as well as for mobile practice. TELE PACK VET X is widely used for various applications such as gastrointestinal endoscopy, bronchoscopy, laparoscopy, arthroscopy, otoscopy and exotic animal endoscopy among others.

To enable swift and easy work, TELE PACK VET X combines all that is needed: monitor, camera, light source and insufflation pump. Consideration has also been given to documentation: Integrated data management enables comprehensive recording of examinations or surgical interventions. Multiple USB ports and an SD card slot are available to store the data.

## Crystal clear display

- 15" LCD display
- Image rotation
- 24-bit color intensity for natural color rendition
- DVI-D video output for connecting HD monitors

## Flexible storage possibilities

- SD card slot for high storage capacity
- USB ports for external hard drives, USB flash drives and post-script printers
- Picture gallery for records
- Playback of saved videos
- Print-ready patient report documentation

## Natural illumination

- HiLux 50 Watt high-performance light source
- Natural colour rendition close to daylight with a color temperature of 5700 K
- Up to 1000 hours lamp operating time

## Easy control combined with highest safety

- Membrane keyboard for wipe-down disinfection
- Hot keys for rapid and direct adjustment
- Arrow keys for intuitive control
- Connection socket for pedal control

## Additional information

- Sturdy, portable housing
- Ergonomically designed handle for convenient transport
- Universal power supply unit: 100 – 240 VAC, 50/60 Hz
- Dimensions (w x h x d): 450 x 350 x 150 mm
- Weight: 7 kg



# TELE PACK VET X <sup>NEW</sup>

A multipurpose unit for diagnostic and surgical procedures including still image and video capture software.



## 69 0450 01 TELE PACK VET X

endoscopic video unit for use with all KARL STORZ TELECAM one-chip camera heads, Veterinary Video Camera Head and video endoscopes, incl. 50 W Hi-Lux light source, with integrated insufflation pump and digital Image Processing Module, 15" LCD monitor, USB/SD memory module, color systems PAL/NTSC, power supply 100 – 240 VAC, 50/60 Hz

including:

**USB Flash Drive**

**Mains Cord**

**Mains Cord, US plug**

**Case**

## 69 0450 02 Same

without **Case**

### Optional Accessories:

**20 0141 30 One-Pedal Footswitch**, digital, one-stage, for irrigation

**20 0402 80 KARL STORZ USB Flash Drive**, 4 GB

**495 KS Adaptor**, for standard KARL STORZ fiber optic light cables (not 495 NCS) to KARL STORZ XENON 100 SCB

**81131021 Spare Lamp**, 50 W

**20 0402 40DE Silicone Keyboard**, with German character set, USB plug, with touchpad

**20 0402 40US Same**, with US English character set

**20 0402 40PT Same**, with Portuguese character set

**20 0402 40ES Same**, with Spanish character set

**20 0402 40FR Same**, with French character set

**20 0402 40RU Same**, with Russian character set

**20 0402 40IT Same**, with Italian character set

### Specifications:

Power supply	100-240 VAC, 50/60 Hz/12VDC	Video out- /input	DVI-D in/out (XGA)
Dimensions w x h x d	450 x 350 x 150 mm	Audio	3,5mm plug, 2 x Microphone, Line-in, Line-out
Weight	7.3 kg	Interface for printing	USB 2.0
Ergonomic carry handle	USB 2.0, SD- memory card	Flow non-regulated:	max. 3.5l/min
Documentation interface		Pressure non -regulated:	max. 760 mmHG

**Component Parts** see COMPONENTS, SPARE PARTS






TELECAM

STORZ  
KARL STORZ — ENDOSKOPE



CAMERA SYSTEMS



2-12

# TELECAM

One-Chip Camera System, Color Systems PAL/NTSC,  
Digital Imaging Processor

## Multifunctional, Cost-Effective Camera System for Veterinary Clinics

The cost-effective **TELECAM SL II** and **Veterinary Video Camera III** camera platform offers functionality, image quality and high performance, ensuring a high degree of patient safety combined with cost-efficiency.

The new KARL STORZ camera system provides such a wide range of applications that it is also suitable for routine use in the OR.



Special Features	TELECAM SL II	Veterinary Video Camera III
<b>Compatible with KARL STORZ Video Endoscopes</b>	●	●
Automatic exposure control for superior image quality and maximum control	●	●
Anti-moiré filter	●	●
Automatic white balance with memory function	●	●
Horizontal resolution of more than 450 lines	●	●
Programmable function keys for four camera functions or peripheral devices	●	●
Freeze Frame to “freeze” images for the purpose of documenting clinical findings, for diagnostic consultations with colleagues or patients and for patient documentation	●	●
Lightweight design allows ergonomic use	●	●
Composite and S-Video output	●	●
Up to 60 times higher light sensitivity due to integrated digital processing technology (IPM)	●	–
Optical Parfocal Zoom Lens, integrated in all camera heads used	●	–
Digital video output (DV)	●	–
Keyboard Input, for entering patient data	●	–
Electronic zoom 1-2x can be adjusted to 4 levels	●	–
Manual or fully automatic exposure control	●	●
Contrast enhancement and digital filter for improved image quality	●	–
Patient Information Display	●	–

# TELECAM SL II

Camera Control Unit



## 202130 11U TELECAM SL II Camera Control Unit

color systems **PAL/NTSC**, with integrated digital Image Processing Module, power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

2x **Connecting Cable**, for controlling peripheral devices

**BNC/BNC Video Cable**

**S-Video (Y/C) Connecting Cable**

**DV Connecting Cable**, 6 pin to 4 pin

**Keyboard**, with US English character set

### Specifications:

Video output	- S-Video signal to 4 pin Mini-DIN socket (2x) - Composite signal at BNC socket - DV signal to 6 pin DV socket	Control Unit (CCU)	Dimensions: 305 x 88 x 254 mm (w x h x d)
Input	keyboard input for title generator and camera functions to 5 pin DIN socket	Control unit	2.7 kg
Control output	3.5 mm stereo jack plug (ACC 1, ACC 2)	Power supply	100-240 VAC, 50/60 Hz
		Certified to	IEC 601-1, 601-2-18, CSA 22.2 No. 601, UL 2601 and CE acc. to MDD, protection class 1/BF

**Further accessories for TELECAM SL II** see pages TP 43-44

**Component Parts** see COMPONENTS, SPARE PARTS

# Veterinary Video Camera III



One-Chip Camera System, Color Systems PAL/NTSC

## Special Features:

All previously compatible TELECAM camera heads can be operated with the new Veterinary Video Camera III platform.

The new platform features the following installed functions:

- Compatible with the PAL and NTSC color systems
- Integrated freeze frame function
- Integrated electronic fiberscope filter
- Patented autoexposure system
- Menu available in 6 languages: (German, English, French, Italian, Portuguese and Spanish)
- Freely programmable camera head buttons



### 69 2360 01 PAL Veterinary Video Camera III

color system **PAL**, power supply 100 – 240 VAC, 50/60 Hz including:

**Veterinary Video Camera III Camera Control Unit**, color systems **PAL/NTSC**, with integrated digital fiberscope filters

#### Mains Cord

2x **Connecting Cable**, for controlling peripheral units, length 180 cm

**BNC/BNC Video Cable**, length 180 cm

**S-Video (Y/C) Connecting Cable**, length 180 cm

**C-MOUNT Lens**, f = 38 mm, soakable

**VET C-MOUNT Camera Head, PAL**

### 69 2361 01 NTSC Veterinary Video Camera III

color system **NTSC**, power supply 100 – 240 VAC, 50/60 Hz including:

**Veterinary Video Camera III Camera Control Unit**, color systems **PAL/NTSC**, with integrated digital fiberscope filters

#### Mains Cord

2x **Connecting Cable**, for controlling peripheral units, length 180 cm

**BNC/BNC Video Cable**, length 180 cm

**S-Video (Y/C) Connecting Cable**, length 180 cm

**C-MOUNT Lens**, f = 38 mm, soakable

**VET C-MOUNT Camera Head, NTSC**

**Optional accessories for Veterinary Video Camera III** see pages TP 43-44

**Component Parts** see COMPONENTS, SPARE PARTS

# Veterinary Video Camera III

Camera Control Unit, Camera Head

## Camera Control Unit



69 2360 20

### 69 2360 11 Veterinary Video Camera III Camera Control Unit

color systems **PAL/NTSC**, with integrated digital fiberscope filters including:

#### Veterinary Video Camera III Camera Control Unit Mains Cord

2x **Connecting Cable**, for controlling peripheral units, length 180 cm

**BNC/BNC Video Cable**, length 180 cm

**S-Video (Y/C) Connecting Cable**, length 180 cm

#### Specifications:

Video output	- Composite signal to BNC socket - S-Video signal to 4-pin Mini DIN socket (2x)	Power supply	100-240 VAC, 50/60 Hz
Control output	3.5 mm stereo jack plug (ACC 1, ACC 2)	Certified to:	IEC 601-1, 601-2-18, CSA 22.2 No. 601, UL 2601, and CE acc. to MDD, protection class 1/BF
Control Unit (CCU)	- Dimensions: 305 x 88 x 254 mm (w x h x d) - Weight: 2.7 kg		

## Veterinary Video Camera III C-MOUNT One-Chip Camera Head

### Veterinary Video Camera C-MOUNT One-Chip Camera Head



69 2150 34 / 69 2151 34

69 2150 34	<b>PAL</b>	<b>Veterinary Video Camera</b>
69 2151 34	<b>NTSC</b>	<b>C-MOUNT One-Chip Camera Head</b>

color systems **PAL/NTSC**, 2 freely programmable camera head buttons

For use with Veterinary Video Camera III Camera Control Unit 69 2360 11

#### Specifications:

Image Sensor	- 1/2" CCD chip	AGC	Microprocessor-controlled
Pixels	- 752 (H) x 582 (V) pixels per chip (PAL) - 768 (H) x 494 (V) pixels per chip (NTSC)	Min. Sensitivity	3 Lux (f = 1.4 mm)
Resolution	> 450 lines (horizontal)	Exposure Control	- 1/50 sec.-1/10000 sec. (PAL) - 1/60 sec.-1/10000 sec. (NTSC)
Signal-to-Noise Ratio	> 50 dB (PAL) > 52 dB (NTSC)	Lens	C-MOUNT: integrated special lens, C-MOUNT

**Optional Accessories for Veterinary Video Camera III** see pages TP 43-44

**Component Parts** see COMPONENTS, SPARE PARTS

# TELE PACK VET X, TELECAM and Veterinary Video Camera III

Camera Heads

For use with TELECAM SL II Camera Control Unit 20 2130 11U, Veterinary Video Camera III Control Unit 69 2360 20, TELE PACK VET X 69 0450 20



20 2120 40 / 20 2121 40

20 2120 40	<b>PAL</b>	<b>TELECAM</b>
20 2121 40	<b>NTSC</b>	<b>One-Chip Camera Head</b>

color systems **PAL/NTSC**, autoclavable, soakable, gas-sterilizable, with integrated Parfocal Zoom Lens, f = 14 – 28 mm (2x), 2 freely programmable camera head buttons, including plastic container 39301 ACT for sterilization



20 2120 30 / 20 2121 30

20 2120 30	<b>PAL</b>	<b>TELECAM</b>
20 2121 30	<b>NTSC</b>	<b>One-Chip Camera Head</b>

color systems **PAL/NTSC**, soakable, gas-sterilizable, with integrated Parfocal Zoom Lens, f = 25 – 50 mm (2x), 2 freely programmable camera head buttons



69 2150 34 / 69 2151 34

69 2150 34	<b>PAL</b>	<b>Veterinary Video Camera</b>
69 2151 34	<b>NTSC</b>	<b>C-MOUNT One-Chip Camera Head</b>

color systems **PAL/NTSC**, 2 freely programmable camera head buttons

Further accessories for TELECAM SL II and Veterinary Video Camera III see pages TP 43-44

# TELE PACK VET X, TELECAM and Veterinary Video Camera III

Accessories




- 20 2000 42 **C-MOUNT Lens**, f = 30 mm, soakable
- 20 2000 43 **C-MOUNT Lens**, f = 38 mm, soakable
- 20 2301 41 **C-MOUNT Lens**, f = 25 mm, soakable
- 20 2301 45 **C-MOUNT Lens**, f = 12 mm, soakable



- 81 1310 21 **Spare Lamp**, 50 W, for TELE PACK VET X



- 20 0141 30 **One-Pedal Footswitch**, digital, one-stage, for irrigation
- 20 2330 70 **Adaptor ACC 1**, for Footswitch 20 0141 30



- 20 0143 30 **Two-Pedal Footswitch**, one-stage



- 20 2330 70 **Adaptor ACC 1**, for Footswitch 20 0101 30

- 29020 GM **Adaptor**, for Machida fiberscopes



- 29020 GN **Adaptor**, for Olympus fiberscopes, new type



- 29020 GO **Adaptor**, for Olympus fiberscopes, old type



- 29020 GP **Adaptor**, for Pentax and Fujinon fiberscopes

# TELE PACK VET X, TELECAM and Veterinary Video Camera III

Accessories



- 20 2002 31 **Keyboard**, with German character set, PS/2 connector, including AT-adaptor and protective cover
- 20 2002 31F **Same**, with French character set
- 20 2002 31I **Same**, with Italian character set
- 20 2002 31P **Same**, with Portuguese character set
- 20 2002 31S **Same**, with Spanish character set
- 20 2002 31U **Same**, with US English character set



- 20 0402 40DE **Silicone Keyboard**, with German character set, USB plug, with touchpad
- 20 0402 40FR **Same**, with French character set
- 20 0402 40IT **Same**, with Italian character set
- 20 0402 40PT **Same**, with Portuguese character set
- 20 0402 40ES **Same**, with Spanish character set
- 20 0402 40RU **Same**, with Russian character set
- 20 0402 40US **Same**, with US English character set



- 20 0402 80 **KARL STORZ USB Flash Drive**, 4 GB



- 20 2002 70 **Extension Cable**, for TELECAM SL II and Veterinary Video Camera III, length 600 cm



- 20 2130 70 **Video Connecting Cable**, between KARL STORZ video endoscopes and TELECAM SL II, Veterinary Video Camera III or TELE PACK VET X, length 60 cm



- 39301 ACT **Plastic Container for Sterilization and Storage**, of camera heads IMAGE1 P3, TRICAM® and TELECAM, **autoclavable**, for use with steam, gas and plasma sterilization

**Please note:** The instrument displayed is not included in the plastic container.  
Only camera heads marked “autoclave” can be placed in the tray for steam sterilization.  
(except IMAGE1 FULL HD camera heads)





# Camera Heads, Camera Control Unit

Matrix

CAMERA SYSTEMS

Camera Control Units		IMAGE 1 HUB™ HD	IMAGE1 HD	
Camera Heads	Cat. No.	22 2010 11-1xx	22 2020 11-1xx	
<b>IMAGE1 FULL HD Camera Heads</b>				
H3-Z	22 2200 55-3	●	●	
H3-FA	22 2200 60-3	●	●	
H3-ZA	22 2200 61-3	●	●	
H3-M COVIEW®	22 2200 54-3	●	●	
H3-Z FI	22 220085-3	●	●	
<b>IMAGE1 SD Camera Heads</b>				
F3	22 2200 33-3 (PAL) 22 2201 33-3 (NTSC)	●	-	
<b>3D TIPCAM®1</b>				
3D TIPCAM®1	26604 AA (PAL) 26614 AA (NTSC)	-	-	
<b>TRICAM® Camera Heads</b>				
TRICAM®, autoclavable	20 2210 40 (PAL) 20 2211 40 (NTSC)	-	-	
TRICAM®	20 2210 30 (PAL) 20 2211 30 (NTSC)	-	-	
TRICAM® PDD Pendel	20 2210 39 (PAL) 20 2211 39 (NTSC)	-	-	
TRICAM® PDD	20 2210 37 (PAL) 20 2211 37 (NTSC)	-	-	
<b>TELECAM-compatible Camera Heads</b>				
VET C-MOUNT Camera Head	69 2150 34 69 2151 34	-	-	
TELECAM, autoclavable	20 2120 40 (PAL) 20 2121 40 (NTSC)	-	-	
TELECAM C-MOUNT	20 2120 34 (PAL) 20 2121 34 (NTSC)	-	-	
TELECAM Beamsplitter (f = 25 mm)	20 2120 32 (PAL) 20 2121 32 (NTSC)	-	-	
TELECAM DCI® II	20 2620 30 (PAL) 20 2621 30 (NTSC)	-	-	
<b>C-CAM™</b>				
C-CAM™	20 2901 32	-	-	
<b>Videoendoscope Adapters</b>				
Video Endoscope Adaptor (for connecting video endoscopes with TELECAM, Veterinary Video Camera III and TELE PACK VET X)	20 2130 70	-	-	
Video Endoscope Adaptor (suitable for fluorescence applications with PDD/AF video endoscopes)		-	-	
Video Endoscope Adaptor (for connecting video endoscopes with IMAGE1)	22 2000 77 (PAL)	●	-	2-12

	3D System	TRICAM® SL II	TELECAM SL II	Veterinary Video Camera III	TELE PACK VET X	C-HUB™
	22 2040 11-114	20 2230 11U1xx	20 2130 11	69 2360 11	69 0450 01 69 0450 02	20 2901 01
	-	-	-	-	-	-
	-	-	-	-	-	-
	-	-	-	-	-	-
	-	-	-	-	-	-
	-	-	-	-	-	-
	-	-	-	-	-	-
	●	-	-	-	-	-
	-	●	-	-	-	-
	-	●	-	-	-	-
	-	●	-	-	-	-
	-	●	-	-	-	-
	-	-	●	●	●	-
	-	-	●	●	●	-
	-	-	●	●	●	-
	-	-	●	●	●	-
	-	-	●	●	●	-
	-	-	-	-	-	●
	-	-	●	●	●	-
	-	-	-	-	-	-
2-12	-	-	-	-	-	-

# Video Endoscopes, Camera Control Unit

## Matrix







CAMERA SYSTEMS

Camera Control Units		IMAGE 1 HUB™ HD	GI	
Video Endoscopes	Cat. No.	22 2010 11U1xx	22 2030 11-1xx	
<b>Bronchoscopy</b>				
Video Bronchoscope	11900 BP 11900 BN	*22 2000 77 *22 2001 77	-	
<b>Gastroenterology</b>				
Feline Videoendoscope	60511 PKS/NKS	*22 2000 77 *22 2001 77	*22 2000 77 *22 2001 77	
Multipurpose Video Endoscopes	60714 PKS/NKS 60914 PKS/NKS	*22 2000 77 *22 2001 77	-	
<b>Urology</b>				
Video Uretero-Renoscope FLEX-X <sup>c</sup>	11278 V 11278 VU	●	-	

\* video endoscope adaptor

	TELECAM SL II FI	TELECAM SL II	Veterinary Video Camera III	TELE PACK VET X	C-HUB™
	20 2130 11-1	20 2130 11	20 2330 11	69 0450 20	20 2901 01
	*20 2130 70	*20 2130 70	-	*20 2130 70	-
	*20 2130 70	*20 2130 70	-	*20 2130 70	-
	*20 2130 70	*20 2130 70	-	*20 2130 70	-
	-	-	-	-	-

# Video Cable and SCB Cable

	<p>20 0400 89    <b>DVI-D Connecting Cable</b>, length 300 cm</p>
	<p>536 MK    <b>BNC/BNC Video Cable</b>, length 200 cm                      536 MF    <b>BNC/BNC Video Cable</b>, length 450 cm                      536 ME    <b>BNC/BNC Video Cable</b>, length 760 cm                      536 MN    <b>BNC/BNC Video Cable</b>, length 15 m</p>
	<p>547 S    <b>S-Video (Y/C) Connecting Cable</b>, length 180 cm                      547 SR    <b>S-Video (Y/C) Connecting Cable</b>, length 500 cm                      547 SN    <b>S-Video (Y/C) Connecting Cable</b>, length 15 m                      547 SE    <b>S-Video (Y/C) Extension Cable</b>, length 500 cm</p>
	<p>547 RG    <b>RGBS Connecting Cable</b>, length 180 cm                      20 2032 70    <b>Special RGB Connecting Cable</b>, length 180 cm</p>
	<p>20 0403 72    <b>SVGA Monitor Cable</b>, length 180 cm                      20 0403 74    <b>SVGA Monitor Cable</b>, length 300 cm</p>
	<p>20 2210 70    <b>Connecting Cable</b>, for controlling peripheral units, length 180 cm                      20 2210 75    <b>Connecting Cable</b>, for controlling video printers or video recorders, length 15 m                      20 2210 76    <b>Extension Cable</b>, for controlling video printers or video recorders, length 500 cm                      20 2210 77    <b>Extension Cable</b>, for controlling video printers or video recorders, length 15 m</p>
	<p>20 0900 70    <b>SCB Connecting Cable</b>, length 30 cm                      20 0903 70    <b>SCB Connecting Cable</b>, length 60 cm                      20 0901 70    <b>SCB Connecting Cable</b>, length 100 cm                      20 0902 70    <b>SCB Connecting Cable</b>, length 500 cm                      20 0904 70    <b>SCB Connecting Cable</b>, length 750 cm                      20 0905 70    <b>SCB Connecting Cable</b>, length 15 m                      20 0906 70    <b>SCB Connecting Cable</b>, length 23 m                      20 0907 70    <b>SCB Connecting Cable</b>, length 30 m</p>

# Sterile Covers

for all KARL STORZ Camera Heads

- 040112-40\* **Camera Cover**, folded inside out, 13 x 242 cm, sterile, package of 25, for use with all standard endoscopic video camera systems
- 040113-50\* **Camera Cover**, for single use, sterile, 13 x 242 cm, telescopic folding with paper insertion aid, package of 50, for use with all currently available endoscopic video camera systems
- 040114-15\* **Camera Cover**, for single use, sterile, 17 x 242 cm, rolled with plastic insertion aid, package of 15, for use with all currently available endoscopic video camera systems
- 040115-40\* **Camera Cover**, for single use, sterile, 13 x 242 cm, telescopic folding with paper insertion aid, package of 40, for use with all currently available endoscopic video camera systems
- 040169-40\* **Camera Cover**, for single use, sterile, 13 x 242 cm, with elastic telescope opening, package of 40, for use with all currently available endoscopic video camera systems as well as for KARL STORZ DCI® camera systems



## Note:

To ensure a safe and convenient use of KARL STORZ ENDOVISION® and IMAGE1 cameras, we recommend the use of sterile covers by:



mtp medical technical promotion gmbh  
Take-Off GewerbePark 46  
D-78579 Neuhausen ob Eck  
Germany  
Telephone: +49/74 67/9 45 04-0  
Fax: +49/74 67/9 45 04-99






# C-HUB™ and C-CAM™ *NEW*

**STORZ**  
KARL STORZ — ENDOSKOPE

CAMERA SYSTEMS



2-12

### Nothing could be easier!

How can a standard exam room with table, chair and computer be converted into an effective examination room with an endoscope and a camera? Quite simple: The KARL STORZ C-CAM™ can be directly connected to a USB interface on a Medical PC via C-HUB™.

Thanks to its CMOS sensor, the C-CAM™ camera head is able to provide a resolution of 640 x 480 pixels. The endoscope can be attached as usual with the proven KARL STORZ grip mechanism.

C-HUB™ features an USB interface that can be connected to all Medical PCs, whether desktops or laptops. It also features an S-Video connection for direct transmission of image data to an external monitor.

If the C-HUB™ is connected to a Medical PC, image data or videos can be directly archived on the hard drive and processed with the relevant programs.

All in all, the KARL STORZ C-CAM™ in conjunction with C-HUB™ offers an effective possibility to perform fast and easy endoscopic examinations in any practice.



20 2901 01

#### 20 2901 01 C-HUB™ Camera Control Unit

for use with C-CAM™ 20 2901 32, Electronic Module 8402 X or compatible CMOS video endoscopes, interfaces: USB 2.0, S-Video output (NTSC), power socket

including:

**C-HUB™ Power Supply**

**S-Video (Y/C) Connecting Cable**

**USB Connecting Cable**



20 2901 32/20 2901 31

#### 20 2901 32 C-CAM™ Camera Head, 8-pin

one-chip CMOS camera head, resolution 640 x 480, focal length  $f = 20$  mm, compatible with C-HUB™ 20 2901 01 and C-MAC® 8402 ZX

#### 20 2901 31 C-CAM™ Camera Head, 6-pin


one-chip CMOS camera head, resolution 640 x 480, focal length  $f = 20$  mm, compatible with C-MAC® 8401 ZX

**Component Parts** see COMPONENTS, SPARE PARTS

## Monitor/Electronic Module

### Special Features:

- Resistant ABS plastic housing
- Splash-proof according to IP54
- 7" TFT wide view angle display with resolution of 800 x 480 pixels
- Ready for use within seconds
- Documentation of still images (JPEG) and videos (MPEG4) on SD memory card
- VESA 75 norm for connecting and attaching racks
- Soft keys enable use within seconds
- Cinch video output for connecting external monitor
- System open for further components
- Operating time with lithium-ion batteries of about 2 hours
- World power supply 100 – 240 VAC, 50/60 Hz
- Operation with line voltage and rechargeable lithium-ion batteries
- Processing of the electronic module: Suitable and validated for the following low-temperature reprocessing methods up to bis max. 60 °C: manual/machine cleaning and disinfection, sterilization with Steris® AMSCO V-PRO 1, Sterrad® (50S, 100S, 200S, NX, 100NX) and EtO gas; High-Level Disinfection (HLD) acc. to US standards
- Additional standards: RTCA/DO-160F, EMI Test Report (German air rescue service DRF Luftrettung)




8402 ZX

8402 ZX

**Monitor for CMOS Endoscopes**, screen size 7", documentation can be stored directly on SD card, rechargeable Li-Ion batteries, power adaptor for EU, UK, USA and Australia, power supply 110 – 240 VAC, 50/60 Hz, additional standards: RTCA/DO-160F, EMI Test Report (German air rescue service DRF Luftrettung), **suitable for wipe disinfection**

**Component Parts** see COMPONENTS, SPARE PARTS





**DATA MANAGEMENT  
AND DOCUMENTATION**

# Data Management and Documentation

KARL STORZ AIDA® compact NEO advanced/standard

	AIDA compact NEO advanced	AIDA compact NEO standard
Digital storage	●	●
Recording of still images, video sequences, voice commentaries	●	●
Live mode via flat/touch screen	●	●
Sterile use of recording functions via:	touch screen	●
	voice control	○
	footswitch	●
	camera head buttons	●
Choice of 3 video inputs on user interface	●	●
File system storage	●	●
Storage on:	CD-ROM, DVD	●
	USB flash drive	●
	internal or external hard drive	●
	network path	●
Multipatient storage	●	●
Cache	●	●
Standard report	●	●
Communication package (DICOM/HL7)	●	●
Integrated image editor	●	●
Integration of image material in applications, i. e. training materials and presentations	●	●
Printer connection via USB/network	●	●
Image formats	BMP, JPG	BMP, JPG
Video formats	MPEG2	MPEG2
Video signal inputs:	analog	Y/C (S-Video), Composite
	digital	DVI-D
Still image resolution	1920 x 1080, 720 x 576	1920 x 1080, 720 x 576
Video sequence resolution	1920 x 1080 (1080p) 1280 x 720 (720p) 720 x 576 (PAL) 720 x 480 (NTSC)	1280 x 720 (720p) 720 x 576 (PAL) 720 x 480 (NTSC)
Adjustable video bit rate	●	●
Compatible with:	SCB control NEO, as of Release 41	●
	OR1™ AV system NEO	●


- Standard
- Options

## The Compact Documentation Solution


The new KARL STORZ AIDA® compact NEO advanced/standard system combines all the required functions for integrated and precise documentation in a single system. All functions are centrally controlled via touch screen, footswitch, camera head buttons or voice control directly from the sterile area.



**AIDA compact NEO advanced/standard:** Review screen



**AIDA compact NEO advanced/standard:** Automatic creation of standard reports



**AIDA compact NEO advanced/standard:** Efficient archiving

### Data Acquisition

AIDA compact NEO advanced/standard records still images, video sequences (in HD quality) and spoken comments of findings and intraoperative procedures directly from the sterile area. Recordings are activated via touch screen, footswitch or camera head buttons.

Live display of camera images on the touch screen enables immediate monitoring and selection of the recorded data.

### Flexible Review

Before final archiving, the saved data can be viewed or listened to on the review screen. Data no longer required can be simply deleted.

Individual images, video and audio sequences can be renamed and given more meaningful names. A pre-defined selection list with keywords simplifies and speeds up data entry.

A voice entry of the case report can even be recorded while viewing video and image files.

### Automated Data Archiving

Once a treatment is completed, AIDA compact NEO advanced/standard automatically stores the data on a DVD, creates a standard report and prints it as an overview, if required.

### Multipatient

Efficient data archiving is assured as the data of several patients can be saved on a DVD.

# Data Management and Documentation

KARL STORZ AIDA® compact NEO advanced/standard

## Special Features:

- Digital storage of still images with a resolution of 1920 x 1080 in both systems (advanced/standard) possible. Video sequences to 1080p with advanced and to 720p with standard variant possible (via HD-SDI)
- Audio files can be archived in both systems
- DICOM/HL7 interface package included
- Printing from review area (individual images with meta-data)
- DVD burning, Blu-ray reading capabilities
- AIDA Restore Configuration supports simple import and export of system settings
- Reference screen with new QuickView (favorites folder)
- Video settings (e.g., contrast, ...) separate for all channels
- Adaptation options for GUI
- Compressed DICOM
- Supports OR1™ Checklist V1.1
- Enhanced support of 15" monitors with OR1™ checklist installation
- Scalable watermark
- High-quality function and switching of image/video quality without additional settings
- Sterile, ergonomic operation via touch screen, camera head buttons and/or footswitch
- Data export to DVD, CD-ROM or USB flash drive, multi-session and multi-patient
- Network storage possible
- Automatic creation of standard reports
- Systems for use in the OR area certified according to EN 60601-1
- Compatible with KARL STORZ Communication Bus (SCB) and KARL STORZ OR1™ AV NEO
- KARL STORZ AIDA® compact advanced/standard as an attractive, digital alternative to video printers, video recorders and dictating machines





# Data Management and Documentation

KARL STORZ AIDA® compact NEO advanced/standard

## 20 0409 13-EN\* **KARL STORZ AIDA® compact NEO advanced**

Documentation system for digital storage of still images, video sequences and audio files, power supply 115/230 VAC, 50/60 Hz including:

**KARL STORZ AIDA® control NEO**, with integrated DVD/CD writer, Blu-ray reading possible

**AIDA compact HD Frame Grabber Card**

**PS/2 Compact Keyboard**, US English, with cover

**KARL STORZ AIDA® compact NEO Software**, with voice control and software protection

**KARL STORZ USB Flash Drive**

2x **Connecting Cables**, for controlling peripheral units

**BNC/BNC Video Cable**

**BNC Connecting Cable**

**DVI-D Connecting Cable**

**Mains Cord**

**Connecting Cable**, Mini-DIN plug to BNC socket

## 20 0409 12-EN\* **KARL STORZ AIDA® compact NEO standard**

Documentation system for digital storage of still images, video sequences and audio files, power supply 115/230 VAC, 50/60 Hz


### Specifications:

Video Systems	- PAL, NTSC	Video Formats	- MPEG2
Signal Inputs	- S-Video (Y/C) - Composite - RGB (only standard) - SD-SDI (only standard) - HD-SDI (only standard) - DVI-D (only advanced)	Audio Formats	- WAV
Image Formats	- JPG - BMP	Storage Media	- DVD+R - DVD+RW - DVD-R - DVD-RW - CD-R - CD-RW - USB flash drive

\* Available in the following languages: DE, ES, FR, IT, PT, PL, RU, DK, SE, JP, CN

**Component Parts** see COMPONENTS, SPARE PARTS

## AIDA – Link to PACS and HIS



### DICOM

DICOM is an acronym for Digital Imaging and Communications in Medicine. It describes a file format accepted worldwide for exchanging image data within clinics. In contrast to pure image formats, the DICOM format also includes patient data. It ensures that information is clearly allocated during data processing. Radiological, ultrasound images and CT scans for example contain meta-information such as demographic patient data, admission dates, etc. which can be transferred as a single unit within the clinic network. The DICOM standard enables straightforward and secure communication of medical imaging equipment and servers to archive image data.

### HL7

HL7 is the acronym for Health Level Seven. It describes an international standard for exchanging medical data. HL7 originated in the USA and is generally used in the Hospital Information System (HIS) for transferring patient-related data which do not contain image data. The objective of HL7 is to improve communication channels within clinic networks.

# Accessories

for AIDA compact NEO advanced/standard

## KARL STORZ TFT Flat Screens



W29519 N



W29519 NB

KARL STORZ TFT Flat Screens	19"
Desktop with pedestal	W29519 N
Wall mounted with VESA 100 adaption	W29519 NB
<b>Inputs:</b>	
SDI	●
HD-SDI	●
RGBS	●
S-Video	●
Composite	●
SOG	●
DVI-D	●
Fiber Optic	—
VGA	●
<b>Outputs:</b>	
SDI	●
HD-SDI	●
RGBS	●
S-Video	●
Composite	●
DVI-D	●
<b>Signal Format Display:</b>	
4:3	●
5:4	●
16:9	●
16:10	●
Picture-in-Picture (PiP, HD – SD)	●
Picture-in-Picture (PiP, HD – HD)	●
PAL/NTSC compatible	●
Picture-Perfect Technology	●

**The following accessories are included:**

**Mains Cord**, length 300 cm

**Power Supply**

**Pedestal**, for 15", 19", 24" and 26" monitors for the 95xx or 97xx series

# Accessories

for AIDA compact NEO advanced/standard

## KARL STORZ TFT Flat Screens

### W29519 N 19" KARL STORZ TFT Flat Screen

desktop with pedestal, color systems **PAL/NTSC**, max. screen resolution 1280 x 1024, power supply 100 – 240 VAC, 50/60 Hz

### W29519 NB 19" KARL STORZ TFT Flat Screen

wall-mounted with VESA 100 adaption, color systems **PAL/NTSC**, max. screen resolution 1280 x 1024, power supply 100 – 240 VAC, 50/60 Hz

including:

- Power Supply**
- Mains Cord**
- BNC/BNC Video Cable**
- S-Video (Y/C) Connecting Cable**
- DVI-D Connecting Cable**
- SVGA Monitor Cable**

#### Specifications:

<b>KARL STORZ TFT Flat Screens</b>	<b>19"</b>
<b>Desktop with pedestal</b>	W29519 N
<b>Wall mounted with VESA 100 adaption</b>	W29519 NB
Brightness	300 cd/m <sup>2</sup>
Max. viewing angle	178° vertical
Pixel distance	0.294 mm
Reaction time	10-16 ms
Contrast ratio	600:1
Mount	100 mm VESA
Weight	6.8 kg
Rated power	65 Watt
Operating conditions	0-38 °C
Storage	-20-60 °C
Rel. humidity	5-85%, non-condensing
Dimensions w x h x d	465 x 400 x 98 mm
Power supply	100-240 VAC
Certified to	EN 60601-1, UL 60601-1, MDD 93/42/EEC, protection class IPX1

**Component Parts** see COMPONENTS, SPARE PARTS

# Accessories

for AIDA compact NEO advanced/standard

## KARL STORZ Touch Screens



20 0904 03/20 0904 06



20 0904 01/20 0904 05

KARL STORZ Touch Screens	15"	19"
<b>Desktop with pedestal</b>	20 0904 03	20 0904 06
<b>Wall mounted with VESA 75/100 adaption</b>	20 0904 01	20 0904 05
<b>Inputs:</b>		
VGA	●	●
SVGA	●	●
XGA	●	●
SXGA	-	●
<b>Screen Diagonal:</b>		
38 cm	●	-
48 cm	-	●
<b>Max. Screen Resolution:</b>		
1024 x 768	●	-
1280 x 1024	-	●

### Specifications:

KARL STORZ Touch Screens	15"	19"
<b>Desktop with pedestal</b>	20 0904 03	20 0904 06
<b>Wall mounted with VESA 75/100 adaption</b>	20 0904 01	20 0904 05
Monitor type	TFT-LCD monitor with active matrix	TFT-LCD monitor with active matrix
Display	304 x 228 mm	376 x 301 mm
Touch screen	AccuTouch	AccuTouch
Colors	16 million with pixel shading	16 million with pixel shading
Brightness	250 cd/m <sup>2</sup> (typ.)	250 cd/m <sup>2</sup> (typ.)
Viewing angle	horizontal total 80° vertical total 80°	horizontal total 85° vertical total 85°
Contrast ratio	400:1 (typ.)	500:1 (typ.)
Monitor reaction time	25 ms	25 ms
Pixel distance	0.30 mm	0.29 mm
Input signal level	0.7 Vpp	0.7 Vpp
Rated power	48 VA	90 VA
Connections	- D-Sub - VGA-compatible analog-RGBS	- D-Sub - VGA-compatible analog-RGBS
Power supply	100-240 VAC, 50/60 Hz	100-240 VAC, 50/60 Hz
Certified to	EN 60601-1, CSA 22.2 Nr. 601.1, UL 2601-1	EN 60601-1, CSA 22.2 Nr. 601.1, UL 2601-1
Mechanical values	- weight: max. 6.5 kg - dimensions w x h x d: 375 x 356 x 205 mm	- weight: max. 10.5 kg - dimensions w x h x d: 432 x 467 x 251 mm

# Accessories

for AIDA compact NEO advanced/standard

## KARL STORZ Touch Screens

### 20 0904 01 15" KARL STORZ Touch Screen

wall-mounted, VGA/SVGA/XGA, resolution max. 1024 x 768 (XGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

#### Touch Pen

Adaptor, VESA 75

Signal cables: RS 232, SVGA

### 20 0904 03 15" KARL STORZ Touch Screen

desktop, VGA/SVGA/XGA, resolution max. 1024 x 768 (XGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

#### Touch Pen

#### Pedestal

Signal cables: RS 232, SVGA

### 20 0904 05 19" KARL STORZ Touch Screen

wall-mounted, VGA/SVGA/XGA/SXGA, resolution max. 1280 x 1024 (SXGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

#### Touch Pen

Signal cables: RS 232, SVGA

### 20 0904 06 19" KARL STORZ Touch Screen

desktop, VGA/SVGA/XGA/SXGA, resolution max. 1280 x 1024 (SXGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

#### Touch Pen

#### Pedestal

Signal cables: RS 232, SVGA

### Accessories:

- 20 0904 81 **Mains Adaptor**, for 15" KARL STORZ Touch Screen 20 0903 31, with 300 cm low voltage cable
- 20 0904 85 **Mains Adaptor**, for 19" KARL STORZ Touch Screen 20 0904 37, with 300 cm low voltage cable
- 20 0903 86 **Touch Pen**, for KARL STORZ touch screens
- 20 0904 86 **Pedestal**, for 15" KARL STORZ Touch Screen 20 0903 31
- 20 0904 87 **Pedestal**, for 19" KARL STORZ Touch Screen 20 0903 37 and 20 0904 37
- 041165-20\* **Sterile Cover**, for 15" KARL STORZ touch screen
- 041265-20\* **Sterile Cover**, for 19" KARL STORZ touch screen








mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

**Component Parts** see COMPONENTS, SPARE PARTS

# Accessories

for AIDA compact NEO advanced/standard

	<p>29005 SZ    <b>Monitor Swivel Arm</b>, height and side adjustable, can be mounted on the left or on the right side, swivel range 180°, reach 600 mm, for side boom as of 800 mm, loading capacity max. 14 kg, with monitor holder VESA 75/100, for Equipment Cart 29005 LAP/GU/BAS</p>
	<p>20 0402 73    <b>Bar Code Scanner</b></p>
	<p>9370 TQ    <b>Color Ink Jet Printer</b>, with horizontal paper feed, incl. mains cord, USB printer cable, length 200 cm, black and color ink jet cartridges</p>
	<p>9372 LA    <b>Color Laser Printer</b>, including mains cord, USB printer cable, length 200 cm</p>
	<p>20 0902 34    <b>PS/2 Compact Keyboard</b>, German version, with cover</p>
	<p>20 0402 40US    <b>Silicone Keyboard</b>, with US English character set, USB plug, with touchpad</p>
	<p>20 0144 30    <b>Two-Pedal Footswitch</b>, two-stage</p>


# Accessories

for AIDA compact NEO advanced/standard

DATA MANAGEMENT  
AND DOCUMENTATION

	20 0402 72	<b>SVGA Connecting Cable</b> , length 600 cm
	20 0400 73	<b>RS 232 Connecting Cable</b> , length 600 cm
	20 0400 74	<b>Printer Cable</b> , length 300 cm
	20 0401 74	<b>Printer Cable</b> , length 500 cm
	20 0103 70	<b>Potential Equalization Cable</b> , length 200 cm
	20 0102 70	<b>Potential Equalization Cable</b> , length 500 cm
	W20475	<b>OR1™ Network Isolator</b> , for galvanic decoupling of units with RJ45 network connection (100 Mbit)
	20 0900 34	<b>Ethernet Network Board</b> , 10/100 Mbit/s
	20 0900 37	<b>Fiber Optic Network Board</b> , 100 Base-FX (SC)
	20 0901 71	<b>Fiber Optic Network Cable</b> , LWL-patch cable with ST and SC plugs, configured, length 100 cm
	20 0900 36	<b>Media Converter</b> , converts a ST-Port 100 Base-FX to a RJ45 100 Base-TX, including power supply
	20 0400 76	<b>RJ45 Network Cable</b> , length 500 cm
	20 2040 80	<b>CD-R</b> , 700 MB, package of 10
	20 2040 81	<b>DVD+RW</b> , 4.7 GB, package of 5
	20 2040 82	<b>DVD+R</b> , 4.7 GB, package of 10
	20 2040 83	<b>DVD-R</b> , 4.7 GB, package of 10





MONITORS

# KARL STORZ FULL HD Monitors <sup>NEW</sup>

## Sophisticated Design in Excellent Quality

The new KARL STORZ monitor series offers the user superior quality combined with outstanding product characteristics that are optimally matched to FULL HD applications in endoscopy. The high-quality aluminium housing in the exclusive KARL STORZ design has been given a silicium dioxide nanocoating finish. The closed housing enables easy cleaning of the monitor through wipe disinfection and thus guarantees optimal hygiene.


The latest electronic technology provides a Picture-in-Picture function, 90° image rotation as well as horizontal and vertical mirror imaging. Input signals are captured directly and immediately transmitted to the display, saving time in the OR. DVI inputs and outputs

guarantee flawless endoscopic FULL HD images with a resolution of 1920 x 1080 pixels in 16:9 format. The progressive scan technology of the KARL STORZ FULL HD system provides imaging with minimal latency, even when there is rapid movement.

The innovative LED backlight technology reduces power consumption and expands the service life of the monitors. Easy and intuitive handling via a ring slider provides immediate feedback on the display. The inbuilt RS 232 interface permits remote control of the monitor, enabling integration of the monitor in the KARL STORZ OR1™ system.

## Special Features:

- 1920 x 1080 pixels, 50/60 Hz in progressive mode, 16:9 format
- High-quality, fully enclosed housing with nanocoating finish for optimal cleaning and disinfection requirements
- Excellent hygiene characteristics with minimum thermal radiation
- Intuitive use without function keys via ring slider behind closed glass front
- Optional fiber connector available
- Easy integration into KARL STORZ OR1™ system
- High level of flexibility due to diverse connection possibilities
- Low voltage protection via external 24VDC mains power supply
- VESA 100 adaptor for flexible installation options
- Fast response time for imaging with minimal latency
- Premium model with dual video inputs available
- Natural color rendition due to optimal integration into the KARL STORZ FULL HD video chain
- 26" monitor with LED backlight for low energy consumption and longer service life for the monitors
- Picture-in-Picture (PiP) and Picture-by-Picture (PbP) display
- Special settings for mirror imaging and image rotation
- Medical grade FDA, UL and CE approval





9619 NB



9626 NB/NB-2

KARL STORZ HD and FULL HD Monitors	19"	26"	26"
Wall mounted with VESA 100 adaption	9619 NB	9626 NB	9626 NB-2
<b>Inputs:</b>			
DVI-D	1x	1x	2x
Fiber Optic	optional	optional	optional
RGBS/VGA	1x	1x	2x
S-Video	1x	1x	2x
Composite/FBAS	1x	1x	2x
<b>Outputs:</b>			
DVI-D	●	●	●
S-Video	●	●	●
Composite/FBAS	●	●	●
<b>Signal Format Display:</b>			
4:3	●	●	●
5:4	●	●	●
16:9	●	●	●
Picture-in-Picture	●	●	●
PAL/NTSC compatible	●	●	●

**The following accessories are included:**

**Mains Cord**

**External 24VDC Power Supply**

**Signal cables: DVI-D, BNC**

**Optional accessories:**

9626 SF

**Pedestal, for 96xx monitor series**

**Specifications:**

KARL STORZ FULL HD Monitors	19"	26"
Desktop with pedestal	optional	optional
Wall mounted with 100 adaption	9619 NB	9626 NB/NB-2
Brightness	280 cd/m <sup>2</sup>	400 cd/m <sup>2</sup>
Max. viewing angle	178° vertical	178° vertical
Pixel distance	0.29 mm	0.30 mm
Reaction time	12 ms	12 ms
Contrast ratio	700:1	700:1
Mount	100 mm VESA	100 mm VESA
Weight	10 kg	14 kg
Rated power	120 W	120 W
Operating conditions	0-40°C	0-40°C
Storage	-20-60°C	-20-60°C
Rel. humidity	max. 80%	max. 80%
Dimensions w x h x d	469.5 x 416 x 75.5 mm	699 x 445.6 x 87.5 mm
Power supply	85-264 VAC	85-264 VAC


**Larger KARL STORZ FULL HD Monitors** see catalog KARL STORZ OR1 NEO™

**Signal cables** see page TP 50

# KARL STORZ HD, HD WIDEVIEW™ and HD LED Backlight Monitors

## Special Features:

- FULL HD display in 16:9 HD format, max. resolution of 1920 x 1200 pixels
- 1080p/50 and 1080p/60 displays possible, PAL/NTSC compatible
- Anti-reflection coated front glass
- Compact and lightweight design
- Easy-to-reach control buttons on the housing front
- Dripwater protected, dustproof housing
- Low voltage protection via external 24VDC mains power supply
- VESA 100 adaptor for flexible installation options
- Fast response time for virtually lag-free imaging
- Optional fiber optic input for OR1™ integration
- Picture-Perfect technology for true-to-life image reproduction that incorporates deinterlacing algorithms
- Picture-in-Picture display (PiP)
- Up to 5 different user profiles can be stored
- Meets the medical requirements for FDA, UL and CE approval
- 26" monitor with LED backlight for low energy consumption and longer service life for the monitors





9515 NB, 9519 NB  
9524 NB/NBO, 9526 NBL, 9726 NB

MONITORS

	TFT Flat Screens		HD WIDEVIEW Monitors		
	15"	19"	24"	26"	
<b>Wall mounted with VESA 100 adaption</b>	9515 NB	9519 NB	9524 NB	9526 NBL	9726 NB
<b>Inputs:</b>					
SDI	●	●	—	—	—
HD-SDI	●	●	●	●	●
3G-SDI	—	—	—	●	●
RGBS	●	●	●	●	●
S-Video	●	●	●	●	●
Composite/FBAS	●	●	●	●	●
SOG	—	●	●	●	●
DVI-D	●	●	●	●	●
Fiber Optic	—	—	—	—	●
VGA	●	●	●	●	●
<b>Outputs:</b>					
SDI	—	●	●	●	●
HD-SDI	—	●	●	●	●
3G-SDI	—	—	—	●	●
RGBS	—	●	●	●	●
S-Video	—	●	●	●	●
Composite/FBAS	—	●	●	●	●
DVI-D	—	●	●	●	●
<b>Signal Format Display:</b>					
4:3	●	●	●	●	●
5:4	—	●	●	●	●
16:9	—	—	●	●	●
16:10	—	—	●	—	—
Picture-in-Picture	—	●	●	●	●
PAL/NTSC compatible	●	●	●	●	●

**The following accessories are included:**

**Mains Cord**

**External 24VDC Power Supply**

**Signal cables**

**Optional accessories:**

9526 SF      **Pedestal**, for 15", 19", 24" and 26" monitors for the 95xx or 97xx series

**Larger KARL STORZ FULL HD Monitors** see catalog KARL STORZ OR1 NEO™

**Signal cables** see page TP 50

# KARL STORZ HD, HD WIDEVIEW™ and HD LED Backlight Monitors

## Specifications:

Wall mounted with VESA 100 adaption	15"	19"	24"	26"	26"
Brightness	9515 NB	9519 NB	9524 NB	9526 NBL	9726 NB
Brightness	430 cd/m <sup>2</sup>	300 cd/m <sup>2</sup>	400 cd/m <sup>2</sup>	400 cd/m <sup>2</sup>	800 cd/m <sup>2</sup>
Max. viewing angle	178° vertical	178° vertical	178° vertical	178° vertical	178° vertical
Pixel distance	0.297 mm	0.294 mm	0.270 mm	0.3 mm	0.3 mm
Reaction time	10-15 ms	10-16 ms	5-12 ms	8 ms	8 ms
Contrast ratio	500:1	600:1	1000:1	1000:1	1000:1
Mount	100 mm VESA	100 mm VESA	100 mm VESA	100 mm VESA	100 mm VESA
Weight	4.8 kg	6.8 kg	7.3 kg	7.3 kg	8.2 kg
Rated power	40 Watt	65 Watt	115 Watt	60 Watt	130 Watt
Operating conditions	0-40 °C	0-38 °C	0-40 °C	0-40 °C	0-40 °C
Storage	-20-60 °C	-20-60 °C	-20-60 °C	-20-60 °C	-20-60 °C
Rel. humidity	5-85%, non-condensing	5-85%, non-condensing	20-85%, non-condensing	20-85%, non-condensing	20-85%, non-condensing
Dimensions w x h x d	386 x 302 x 81 mm	465 x 400 x 98 mm	597 x 401 x 100 mm	73 x 418 x 88 mm	673 x 418 x 88 mm
Power supply	100-240 VAC	100-240 VAC	100-240 VAC	100-240 VAC	100-240 VAC

9515 NB **15" Monitor**

wall-mounted with VESA 100 adaption, color systems **PAL/NTSC**, max. screen resolution 1024 x 768, power supply 100 – 240 VAC, 50/60 Hz

including:

**Power Supply**, flat screens

Signal cables: S-Video (Y/C), DVI-D

9519 NB **19" HD Monitor**

wall-mounted with VESA 100 adaption, color systems **PAL/NTSC**, max. screen resolution 1280 x 1024, power supply 100 – 240 VAC, 50/60 Hz

including:

**Power Supply**

**Mains Cord**

Signal cables: DVI-D, BNC

9524 NB **24" HD WIDEVIEW Monitor**

wall-mounted with VESA 100 adaption, color systems **PAL/NTSC**, max. screen resolution 1920 x 1200, image format 16:10, power supply 100 – 240 VAC, 50/60 Hz

including:

**Power Supply**

**Mains Cord**

Signal cables: DVI-D, BNC

**LED Backlight Monitors** *NEW*

9526 NBL **26" HD Monitor with LED Backlight**

wall-mounted with VESA 100 adaption, color systems **PAL/NTSC**, max. screen resolution 1920 x 1080, image format 16:9, power supply 100 – 240 VAC, 50/60 Hz

including:

**Power Supply**

**Mains Cord**

Signal cables: DVI-D, BNC

9726 NB **26" HD Monitor  
 (Highbright with LED Backlight)**

wall-mounted with VESA 100 adaption, optical input, color systems **PAL/NTSC**, max. screen resolution 1920 x 1080, image format 16:9, power supply 100 – 240 VAC, 50/60 Hz

including:

**Power Supply**

**Mains Cord**

Signal cables: DVI-D, BNC

# 15" and 19" Touch Screens



20 0904 03/20 0904 06



20 0904 01/20 0904 05

MONITORS

KARL STORZ Touch Screens	15"	19"
<b>Desktop with pedestal</b>	20 0904 03	20 0904 06
<b>Wall mounted with VESA 75/100 adaption</b>	20 0904 01	20 0904 05
<b>Inputs:</b>		
VGA	●	●
SVGA	●	●
XGA	●	●
SXGA	-	●
<b>Screen Diagonal:</b>		
38 cm	●	-
48 cm	-	●
<b>Max. Screen Resolution:</b>		
1024 x 768	●	-
1280 x 1024	-	●

## Specifications:

KARL STORZ Touch Screens	15"	19"
<b>Desktop with pedestal</b>	20 0904 03	20 0904 06
<b>Wall mounted with VESA 75/100 adaption</b>	20 0904 01	20 0904 05
Monitor type	TFT-LCD monitor with active matrix	TFT-LCD monitor with active matrix
Display	304 x 228 mm	376 x 301 mm
Touch screen	AccuTouch	AccuTouch
Colors	16 million with pixel shading	16 million with pixel shading
Brightness	250 cd/m <sup>2</sup> (typ.)	250 cd/m <sup>2</sup> (typ.)
Viewing angle	horizontal total 80° vertical insg. 80°	horizontal total 85° vertical insg. 85°
Contrast ratio	400:1 (typ.)	500:1 (typ.)
Monitor reaction time	25 ms	25 ms
Pixel distance	0.30 mm	0.29 mm
Input signal level	0.7 Vpp	0.7 Vpp
Rated power	48 VA	90 VA
Connections	- D-Sub - VGA-compatible analog-RGBS	- D-Sub - VGA-compatible analog-RGBS
Power supply	100-240 VAC, 50/60 Hz	100-240 VAC, 50/60 Hz
Certified to	EN 60601-1, CSA 22.2 No. 601.1, UL 2601-1	EN 60601-1, CSA 22.2 No. 601.1, UL 2601-1
Mechanical values	- weight: max. 6.5 kg - dimensions w x h x d: 375 x 356 x 205 mm	- weight: max. 10.5 kg - dimensions w x h x d: 432 x 467 x 251 mm

2-12



**20 0904 01    15" KARL STORZ Touch Screen**

wall-mounted, VGA/SVGA/XGA, resolution max. 1024 x 768 (XGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

**Touch Pen**

**Adaptor**, VESA 75

Signal cables: RS 232, SVGA

**20 0904 03    15" KARL STORZ Touch Screen**

desktop, VGA/SVGA/XGA, resolution max. 1024 x 768 (XGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

**Touch Pen**

**Pedestal**

Signal cables: RS 232, SVGA

**20 0904 05    19" KARL STORZ Touch Screen**

wall-mounted, VGA/SVGA/XGA/SXGA, resolution max. 1280 x 1024 (SXGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

**Touch Pen**

Signal cables: RS 232, SVGA

**20 0904 06    19" KARL STORZ Touch Screen**


desktop, VGA/SVGA/XGA/SXGA, resolution max. 1280 x 1024 (SXGA mode), power supply 100 – 240 VAC, 50/60 Hz, with mains adaptor, CD driver and 3 touch screen covers

including:

**Touch Pen**

**Pedestal**

Signal cables: RS 232, SVGA




The new ZeroWire video transmission system from KARL STORZ enables the wireless integration of secondary monitors into the camera system.


Developed specifically for use in the OR area, ZeroWire replaces the DVI cable between the Camera Control Unit (CCU) and the monitor. The unit simplifies the installation of additional monitors in the OR. ZeroWire optimizes the clinic workflow and reduces preparation time in the operating rooms. ZeroWire offers the user the best FULL HD quality for wireless connections available on the market.

KARL STORZ ZeroWire comprises a transmitter and receiver that transmit FULL HD video signals from IMAGE1 HD CCU to a secondary monitor. Costly and sensitive light fibers or DVI cables are no longer required. ZeroWire delivers 1080p FULL HD imagery with minimal latency which is ensured by the H.264 compression method. This also provides razor-sharp image quality that can barely be distinguished from a conventional DVI connection.

MONITORS



**Image quality of a competitor**



**Very high image quality with ZeroWire**

Comparison of the image quality of a competitor with ZeroWire from KARL STORZ. ZeroWire image quality of a transmitted video can barely be distinguished from a DVI connection.

2-12



9500 ZW

- 9500 ZW     **ZeroWire Video Transmission Set,**  
transmitter and receiver  
including:  
2x **Radiance Y-Cable 19"/26"**  
2x **VESA Holder 24"/26"**  
2x **VESA Holder 19"/23"**  
2x **2ft DVI-D Cable**, for CCU-Tx and  
Monitor-Rx connection  
**Instruction Manual**

**Further Accessories:**

- 9500 ZW-1     **ZeroWire 42"/52"/55"**, accessory for video  
transmission set, with 2ft DVI-D cable,  
power supply and VESA holder 42"/52"/55"
- 9500 ZW-2     **ZeroWire 2ft Spare DVI-D Cable**,  
for CCU-Tx and Monitor-Rx connection
- 9500 ZW-4     **ZeroWire Power Supply**
- 9500 ZW-5     **ZeroWire Table Stand**
- 9500 ZW-6     **ZeroWire Radiance Y-Cable 19"/26"**,  
for looping into the monitor power supply
- 9500 ZW-7     **ZeroWire VESA Holder**, for mounting a transmitter  
or receiver to a 24" or 26" monitor
- 9500 ZW-8     **ZeroWire VESA Holder**, for mounting a transmitter  
or receiver to a 19" or 23" monitor
- 9500 ZW-9     **ZeroWire VESA Holder**, for mounting a transmitter  
or receiver to 42", 52" and 55" monitors

**Component Parts** see COMPONENTS, SPARE PARTS

# Color Monitors

for Video Cameras, Color Systems PAL/NTSC,  
without IEC 60601-1 certification



69274 N      **LCD Color Monitor**


NeovoX-17B, desktop version, color systems **PAL/NTSC**, resolution max. 1280 x 1024, power supply 100 – 240 V, screen 17", 1x S-Video in, 1x DVI-D, VGA input

**Specifications:**


Monitors	69273 N
Color Systems	PAL, NTSC
Screen Diagonal	35 cm
Inputs/Outputs	- Composite signal to BNC socket - S-Video signal to 4 pin Mini-DIN socket - VGA input
Power Input	48 Watt
Dimensions w x h x d	343 x 354 x 264 mm
Weight	6.5 kg
Power Supply	100-240 VAC, 50 Hz

# Accessories and Optional Accessories

for HD WIDEVIEW Monitors,  
TFT Flat Screens and Touch Screens

			for use with		
			HD WIDEVIEW Monitors	TFT Flat Screen	Touch Screens
<b>Accessories</b>					
	400 A	<b>Mains Cord</b> , length 300 cm	●	●	-
	9523 PS	<b>External 24VDC Power Supply</b>	●	-	-
	9415 NSF	<b>Pedestal</b> (only 9415 NN)	-	●	-
	9419 NSF	<b>Pedestal</b> (only 9419 N, 9519 N, 9523 N, 9524 N/NO)	●	●	-
	9526 SF	<b>Pedestal</b> , for 15", 19", 24" and 26" monitors for the 95xx or 97xx series	●	-	-
	20 0904 86	<b>Pedestal</b> , for 15" KARL STORZ Touch Screen 20 0903 31	-	-	●
	20 0903 86	<b>Touch Pen</b> , for KARL STORZ touch screens	-	-	●
	20 0400 73	<b>RS 232 Connecting Cable</b> , length 600 cm	-	-	●
	20 0402 72	<b>SVGA Connecting Cable</b> , length 600 cm	-	-	●
	20 0904 83	<b>Adaptor</b> , VESA 75	-	-	●
	20 0904 87	<b>Pedestal</b> , for 19" KARL STORZ Touch Screen 20 0903 37 and 20 0904 37	-	-	●
<b>Optional Accessories</b>					
	9419 PS	<b>Power Supply</b>	-	●	-
	20 0904 81	<b>Mains Adaptor</b> , for 15" KARL STORZ Touch Screen 20 0903 31, with 300 cm low voltage cable	-	-	●
	20 0904 88	<b>Same</b> , with 10 m low voltage cable	-	-	●
	20 0904 85	<b>Mains Adaptor</b> , for 19" KARL STORZ Touch Screen 20 0904 37, with 300 cm low voltage cable	-	-	●
	041165-20*	<b>Sterile Cover</b> , for 15" KARL STORZ touch screen	-	-	●
	041265-20*	<b>Sterile Cover</b> , for 19" KARL STORZ touch screen	-	-	●





LIGHT SOURCES

2-12

# Light Sources

## Cold Light

In 1960, Karl Storz discovered that it was possible to transmit light with fiber optic light cables. This discovery marked the birth of cold light endoscopy. From a light source outside the body, light is transmitted via a fiber optic light cable through an endoscope to the examination site. A patent was granted for this new type of light transmission.

Only specific and particularly powerful Halogen, LED or Xenon light sources are used in today's cold light projectors from KARL STORZ.

A widely recognized feature of the KARL STORZ illumination system is the perfectly matched light system, consisting of a lamp and a concave mirror which project light evenly and effectively onto the surface of the light guide.

A vaporized concave mirror, or the so-called cold light mirror, is designed to reflect visible light while allowing the transmission of infrared radiation.

## Light and Radiation

“Light” often refers to the portion of electromagnetic radiation visible to the human eye. The sun emits some of the energy powered by nuclear fusion at 4 and 11 million Kelvin from the solar photosphere which has a temperature of approximately 5800 Kelvin. Only those rays with wavelengths exceeding 295 nm reach the Earth's surface. Other rays are absorbed by the ozone layer in the upper atmosphere so that we are provided with relatively good protection against ultraviolet solar radiation.

There are no exact bounds to the optical spectrum as this can only be determined according to the amount of radiation reaching the retina as well as the eye sensitivity of the individual observer. The lower limit is generally considered to be between 360 nm and 400 nm and the upper limit between 760 nm and 830 nm.

Most solar radiation is located in the green part of the spectrum and it is no coincidence that the human eye is most sensitive to light at this wavelength. On the one hand, solar radiation is our most important source of radiation, light and heat and, on the other, it has influenced the development of a highly sensitive light sensor system in humans.

To measure the eye's response to light, all light processes are measured in accordance with the sensitivity curve of the human eye.

Light output or luminous flux is derived from this radiant power, which is expressed in lumens. The lumen is the basic parameter in lighting technology for light intensity, illumination, and luminance.

Although the light output of a lamp, i.e. the ratio between the radiated luminous flux and generated

electric power expressed in lumens/watts (light output = lumen/watt), is only one of many features to be considered when choosing a particular lamp for a specific application, it is nevertheless one of the most important parameters. A high light output offers the commercial benefit of low energy costs as well as an ecological advantage of resource conservation. These aspects are becoming increasingly important due to growing environmental consciousness. KARL STORZ also has a strong focus on these environmental values and continues to lead the market in this area.

To achieve high light output, it is necessary that:

- Energy in the lamp is generated by means of highly efficient processes
- Radiance should have a high amount of visible light
- Heat emission to the environment is low by means of good conduction and convection.

## Color Temperature

Light output and color temperature are closely related as the human perception of light has undergone a long evolutionary process under the influence of solar radiation. The human eye is most sensitive to light at about 555 nm, i.e. in the green area of the spectrum. Sensitivity drops off at the red and at the blue ends and disappears completely at about 380 nm and 750 nm. Sunshine has a temperature of 5800 Kelvin.

In contrast, light sources subject to the laws of physics have a completely different spectral distribution. Therefore such processes are required that transmit as much radiation as possible to the area of the eye governing sensitivity. The generation of light through thermal radiation is an important feature of illumination technology. This applies to both incandescent and gas discharge lamps so that the color spectrum, i.e. color temperature, can be assigned a single numerical measurement.

In accordance with the latest technology, LED (light-emitting diodes) are also used. In contrast to **incandescent light bulbs**, light-emitting diodes are not **heat generators**. LEDs emit narrow-spectrum light, which is virtually **monochrome**.

The following table provides a definition and a classification of color temperatures:

Color Temperature	Color of Light
< 3300 Kelvin	Warm white
3300 – 5000 Kelvin	Neutral white
6000 Kelvin	Daylight white
6400 Kelvin	Daylight white LED



# Cold Light Fountain HALOGEN



LIGHT SOURCES

2-12

# Cold Light Fountain HALOGEN

Overview

## Cold Light Fountain VETERINARY



## Cold Light Fountain HALOGEN 250 twin



## Halogen Reflector Lamps

The light bulb owes its comeback not just to halogen technology but to a larger extent to the glass reflectors associated with halogen lamps. The surface of these lamps are specially designed for optimal light focus. A further technical feature is the multi-layered coating on the surface which reduces thermal stress in the light bundle by 60%. This is why these lamps are also called cold light mirror reflector lamps. Further technical benefits include a halogen light bulb with an axial coil for optimal light directability, which is especially useful for light output in cold light sources, and a special quartz glass bulb which retains UV radiation. It is possible to raise the temperature of the coil in these lamps to 3400 Kelvin to produce a whiter light than that of standard light bulbs (2800 Kelvin).

A special feature of these lamps is that light can be directed where required, i.e., onto a particular point.

All KARL STORZ halogen cold light sources of at least 150 Watt are suitable for all types of endoscopic procedures. Integrated light intensity adjustment allows illumination to be adapted to a specific OR situation. The halogen cold light models have a color temperature of approximately 3400 Kelvin.

The lamp systems optimized by KARL STORZ provide high light output and their modern design ensures ease of use as well as unproblematic service and exchange of light sources.

## Cold Light Fountain VETERINARY

### Special Features:

- Compact design
- Optimum lighting power
- Manual light intensity control, 3 settings
- Integrated insufflation pump



**69 1110 01    Cold Light Fountain VETERINARY**

with integrated insufflation pump, power supply 230 VAC, 50/60 Hz

including:

**Mains Cord**

**69 1110 01 C    Same**

Power supply 115 VAC, 50/60 Hz

LIGHT SOURCES

### Specifications:

Lamp type	Halogen 24 V, 150 Watt	Dimensions w x h x d	305 x 84 x 265 mm
Color temperature	approx. 3400 K	Weight	6 kg
Light outlets	1	Certified to	IEC 601-1, CE acc. to MDD, protection class 1/CF
Light intensity adjustment	in 3 steps		

**Spare Lamps for Cold Light Fountains HALOGEN** see page TP 89

**Component Parts** see COMPONENTS, SPARE PARTS

## Cold Light Fountain HALOGEN 250 twin

### Special Features:

- Compact and lightweight design
- Automatic backup operation in case of lamp failure
- Manual light intensity control



### 20 1133 01 Cold Light Fountain HALOGEN 250 twin

power supply 100/120/230/240 VAC, 50/60 Hz

including:

**Mains Cord**

### Specifications:

Lamp type	HALOGEN 24 V, 250 Watt	Dimensions w x h x d	305 x 109 x 265 mm
Color temperature	approximately 3400 K	Weight	7 kg
Light outlets	1	Certified to	IEC 601-1, CE acc. to MDD, protection class 1/CF
Light intensity adjustment	adjustable in 3 steps		

**Spare Lamps for Cold Light Fountains HALOGEN** see page TP 89

**Component Parts** see COMPONENTS, SPARE PARTS

# Cold Light Fountains

Accessories

for use with

*Cold Light Fountain  
VETERINARY*

*Cold Light Fountain  
HALOGEN 250 twin*

	103	<b>Halogen Spare Lamp,</b> 150 W, 15 V	●	-
	105	<b>Halogen Spare Lamp,</b> 250 W, 24 V	-	●



# Cold Light Fountain LED



LIGHT SOURCES

2-12

## Cold Light Fountain LED

### Cold Light Fountain LED NOVA 150



**NEW**

### Cold Light Fountain LED NOVA 100 twin



**NEW**

### Cold Light Fountain Power LED 175 SCB



**NEW**

### Cold Light Fountains LED

The advantages of LED (Light Emitting Diode) light sources over conventional light sources are clear to see: LEDs combine high luminance with a very long service life. The average operational life of LED lamps is 30,000 hours: this results in considerable savings in maintenance costs as the need to change lamps is virtually eliminated.

Minimal heat generation also offers a number of advantages. As ventilation is no longer required to cool lamps, disturbing fan noises are a thing of the past. With an optimized ventilation system, it is possible to construct LED cold light fountains with a more compact housing. This frees up valuable space on the PC tower or in the OR.

With a color temperature of approx. 6400 K, LEDs provide natural daylight coloration. This enables exact color reproduction and, therefore, a more accurate diagnosis of the examined tissue.

A self-evident yet nevertheless extremely important factor is that all KARL STORZ light cables can easily be connected to all LED cold light sources.

#### Key Features and Benefits:

- High energy efficiency
- Quiet operation: No fan noise
- High environmental performance: No lead or mercury
- Economic benefits

#### The following LED light sources are available:

- LED NOVA 150
- LED NOVA 100 twin
- Power LED 175



# Cold Light Fountain LED NOVA 150 <sup>NEW</sup>

## Special Features:

- The best of both worlds: High cost-effectiveness coupled with bright and durable LED technology
- Average lamp life of approx. 30,000 hours
- Very high energy efficiency with the latest LED technology
- Extremely quiet operation combined with minimal heat generation
- High environmental performance as neither lead nor mercury are used in the LED lamp
- Continuously variable, adjustable brightness control



### 20 1612 01 Cold Light Fountain LED NOVA 150

power supply 100 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

## Specifications:

Lamp type	high-performance LED	Weight	2.3 kg
Color temperature	approximately 6400 K	Certified to	IEC 601-1, protection class 1/CF
Light outlets	1		
Dimensions w x h x d	305 x 84.5 x 238 mm		

**Component Parts** see COMPONENTS, SPARE PARTS

# Cold Light Fountain <sup>NEW</sup> LED NOVA 100 twin

## Special Features:

- No ventilator, low noise level
- Average lamp life of approx. 23,000 hours
- Suitable for all KARL STORZ standard light cables
- Compact, lightweight construction



## 20 1611 01 Cold Light Fountain LED NOVA 100 twin

with two LED lamps and two light outlets, power supply  
110 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

## Specifications:

Lamp type	LED
Color temperature	approximately 6400 K
Light outlets	2
Dimensions w x h x d	305 x 84.5 x 243 mm

Weight	3.8 kg
Certified to	IEC 601-1, protection class 1/CF

**Component Parts** see COMPONENTS, SPARE PARTS

# Cold Light Fountain <sup>NEW</sup> Power LED 175 SCB

## Special Features:

- Low noise level
- Powerful LED light source similar to Xenon technology
- Lamp life of approx. 30,000 hours
- Energy-saving due to adjustable brightness control
- High serviceability
- SCB – light source control possible via KARL STORZ Communication Bus



### 20161401-1 Cold Light Fountain Power LED 175 SCB

power supply 110 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

**SCB Connecting Cable**


## Specifications:

Lamp type	high-performance LED	Weight	4 kg
Color temperature	approximately 6400 K	Certified to	IEC 601-1, protection class 1/CF
Light adjustment	in 20 steps	miscellaneous	two SCB ports
Light outlets	1		
Dimensions w x h x d	305 x 110 x 233 mm		

**Component Parts** see COMPONENTS, SPARE PARTS



# Cold Light Fountains XENON



LIGHT SOURCES

2-12



# Cold Light Fountains XENON

Overview

## Cold Light Fountain XENON 100



## Xenon Cold Light Sources

Xenon technology by KARL STORZ marks the next major development in illumination systems. Xenon cold light sources have two distinct advantages over standard light bulbs: the Xenon lamp produces twice the light output as a modern halogen lamp while consuming 2/3 of the energy. The improved light output provides brighter and wider illumination of the examination site. The natural daylight color temperature of the Xenon light source is accommodating towards the human eye.

## Cold Light Fountain XENON NOVA® 175



In addition to Xenon Cold Light Fountains, KARL STORZ also offer several units at various price ranges to provide optimal customer satisfaction and maximum customer benefits.

## Cold Light Fountain XENON 175 SCB



The XENON 175 SCB and XENON NOVA® 175 models are operated with a 175 Watt Xenon lamp, providing excellent illumination with a natural color temperature of 6000 Kelvin. A mechanical dimmer provides a continuous manual adjustment of brightness. The lamp current remains unaltered which ensures maximum lamp life.

## Cold Light Fountain XENON NOVA® 300



Models from the XENON 300 SCB and XENON NOVA® 300 series meet all requirements for endoscopic procedures and also provide optimal image and video documentation. With a natural color temperature of 6000 Kelvin and equipped with a 300 Watt Xenon lamp, these models are particularly suitable for the aforementioned documentation.

The XENON NOVA® 300 has a mechanical dimmer to regulate brightness. With the electronic XENON 300 SCB, light is adjusted via a microprocessor-controlled optomechanical dimmer. The unaltered lamp current also ensures maximum lamp life. Light adjustment can be controlled manually via push button or automatically with the video camera's output signal. The standard XENON 175 SCB and XENON 300 SCB models can be operated with the KARL STORZ-SCB system via touch screen or the VOICE1® voice control function.

## Cold Light Fountain XENON 300 SCB



# Cold Light Fountain XENON 100 SCB

with Integrated Insufflation Pump

for use with KARL STORZ video endoscopes

## Special Features:

- Optimal light efficiency
- Electronically adjustable light intensity
- Powerful ventilation with low noise level
- No overheating of environment
- Special diaphragm focuses the light onto the field of view
- New dielectric heat protection filter improves light efficiency by around 20%
- With KARL STORZ Communication Bus (KARL STORZ-SCB)
- Problem-free lamp replacement
- Integrated insufflation pump
- High light intensity
- Xenon lamp provides daylight coloration
- Display of lamp service life
- Insufflation pump with 3 performance levels
- 100 W Xenon lamp



**NEW** 69 1326 01-1 **Cold Light Fountain XENON 100 SCB**

incl. integrated insufflation pump, power supply 100 – 240 VAC, 50/60 Hz, for use with KARL STORZ video endoscopes

**including:**

**Mains Cord**

**Accessories:**

20 1325 80 **Air Filter**, XENON 100

495 KS **Adaptor**, for standard KARL STORZ fiber optic light cables (not 495 NCS) to KARL STORZ XENON 100 SCB

20 0300 30 **Adaptor**, for KARL STORZ fiber optic light cables to KARL STORZ XENON 100 SCB, includes LUER-Lock adaptor to integrated insufflation pump

## Specifications

Xenon lamp	- capacity: 100 W - lamp life: approx. 500 h - color temperature: approx. 6000 K	Power supply	- 230 VAC +/-10%, 50/60 Hz - 115 VAC +/-10%, 50/60 Hz
Pump	- pressure levels: 3 - max. pressure: 0.51 bar - max. flow rate: 5 l/min.	Dimensions w x h x d	305 x 110 x 380 mm
Operating conditions	- operating temperature: 0 °C – +40 °C - humidity: 30 – 90%	Weight	5.4 kg
		Certified to	- IEC 60601-1, 60601-2-18, UL 2601.1, CSA 22.2, No. 601.1-M 90 - Mode of operation: continuous

The system can be operated with KARL STORZ-SCB units.

**Spare Lamps and Spare Lamp Modules for Cold Light Fountains XENON** see page TP 105

**Component Parts** see COMPONENTS, SPARE PARTS



## Cold Light Fountain XENON NOVA® 175

### Special Features:

- High light intensity due to 175 W Xenon lamp, thus ideally suitable for applications requiring high illumination levels
- Continuous manual adjustment of light intensity
- Excellent price/performance ratio



**69 1315 01 Cold Light Fountain XENON NOVA® 175**

power supply 100 – 125/220 – 240 VAC, 50/60 Hz  
including:

**Mains Cord**

LIGHT SOURCES

### Specifications:

Lamp type	XENON 15 V, 175 Watt	Dimensions w x h x d	305 x 110 x 260 mm
Color temperature	approximately 6000 K	Weight	4 kg
Light outlets	1	Certified to	IEC 601-1, CE acc. to MDD, protection class 1/CF
Light intensity adjustment	continuous manual adjustment		

**For use with video endoscopes, Adaptor 69487 KVE is required.**

**Spare Lamps and Spare Lamp Modules for Cold Light Fountains XENON** see page TP 105

**Component Parts** see COMPONENTS, SPARE PARTS

# Cold Light Fountain XENON 175 SCB

**Special Features:**

- Built-in antifog pump
- With integrated KARL STORZ Communication Bus (KARL STORZ-SCB)



**20 1321 01-1 Cold Light Fountain XENON 175 SCB**

power supply 100 – 125/220 – 240 VAC, 50/60 Hz  
including:

**Mains Cord**

**Silicone Tubing Set**, length 250 cm

**SCB Connecting Cable**, length 100 cm

**Specifications:**

Lamp type	XENON 15 V, 175 Watt
Color temperature	approximately 6000 K
Light outlets	1
Light intensity adjustment	continuously adjustable, either via a membrane keyboard or KARL STORZ Communication Bus Signal

Dimensions w x h x d	305 x 165 x 260 mm
Weight	7.7 kg
Certified to	IEC 601-1 and UL 544, protection class 1/CF

**For use with video endoscopes, Adaptor 69487 KVE is required.**

**Spare Lamps and Spare Lamp Modules for Cold Light Fountains XENON** see page TP 105

**Component Parts** see COMPONENTS, SPARE PARTS

## Cold Light Fountain XENON NOVA® 300

### Special features:

- Extremely high light intensity due to 300 Watt Xenon lamp
- Continuous manual adjustment of light intensity
- Excellent price/performance ratio



### 20 1340 01 Cold Light Fountain XENON NOVA® 300

power supply 100 – 125/220 – 240 VAC, 50/60 Hz

including:

**Mains Cord**

### Specifications:

Lamp type	XENON 15 V, 300 Watt	Dimensions w x h x d	305 x 110 x 260 mm
Color temperature	approximately 6000 K	Weight	4.2 kg
Light outlets	1	Certified to	IEC 601-1, CE acc. to MDD, protection class 1/CF
Light intensity adjustment	continuous manual adjustment		

**For use with video endoscopes, Adaptor 69487 KVE is required.**

**Spare Lamps and Spare Lamp Modules for Cold Light Fountains** see page TP 105

**Component Parts** see COMPONENTS, SPARE PARTS

# Cold Light Fountain XENON 300 SCB

**Special Features:**

- Extremely high light intensity due to 300 Watt Xenon lamp
- Built-in antifog pump
- With integrated KARL STORZ Communication Bus (KARL STORZ-SCB)



**20 1331 01-1 Cold Light Fountain XENON 300 SCB**

power supply 100 – 125/220 – 240 VAC, 50/60 Hz  
including:

**Mains Cord**

**Silicone Tubing Set**, length 250 cm

**SCB Connecting Cable**, length 100 cm

**Specifications:**

Lamp type	XENON 15 V, 300 Watt	Dimensions w x h x d	305 x 165 x 335 mm
Color temperature	approximately 6000 K	Weight	7.96 kg
Light outlets	1	Certified to	IEC 601-1 and UL 544, protection class 1/CF
Light intensity adjustment	continuously adjustable, either via a membrane keyboard or KARL STORZ Communication Bus Signal		






**For use with video endoscopes, Adaptor 69487 KVE is required.**

**Spare Lamps and Spare Lamp Modules for Cold Light Fountains** see page TP 105

**Component Parts** see COMPONENTS, SPARE PARTS

# Cold Light Fountains XENON

Spare Lamps and Spare Lamp Modules

			for use with				
			XENON 100 SCB	XENON NOVA® 175	XENON 175 SCB	XENON NOVA® 300	XENON 300 SCB
	20 1320 26	<b>XENON Spare Lamp,</b> 175 W, 15 V	-	●	●	-	-
	20 1330 28	<b>XENON Spare Lamp,</b> 300 W, 15 V	-	-	-	●	●
	20 1325 26	<b>XENON Spare Lamp,</b> 100 W, 15 V	●	-	-	-	-
	20 1320 25	<b>XENON Spare Lamp Module,</b> 175 W, 15 V	-	-	●	-	-
	20 1330 27	<b>XENON Spare Lamp Module,</b> 300 W, 15 V	-	-	-	-	●

LIGHT SOURCES

# Fiber Optic Light Cables

for Cold Light Fountains



## Fiber Optic Light Cable with Straight Connector

- 69495 NE **Fiber Optic Light Cable**, diameter 4.8 mm, length 300 cm
- 69495 NCS **Fiber Optic Light Cable**, extremely heat-resistant, diameter 4.8 mm, length 250 cm
- NEW** 495 NCSC **Fiber Optic Light Cable**, extremely heat-resistant, with safety locking device, diameter 4.8 mm, length 250 cm
- 69495 NB **Fiber Optic Light Cable**, diameter 4.8 mm, length 180 cm
- 69495 ND **Fiber Optic Light Cable**, diameter 3.5 mm, length 300 cm
- 69495 NA **Fiber Optic Light Cable**, diameter 3.5 mm, length 230 cm
- NEW** 495 NAC **Fiber Optic Light Cable**, with safety locking device, diameter 3.5 mm, length 230 cm
- 69495 NL **Fiber Optic Light Cable**, diameter 3.5 mm, length 180 cm
- 69495 NTA **Fiber Optic Light Cable**, diameter 2.5 mm, length 230 cm
- 69495 NT **Fiber Optic Light Cable**, diameter 2.5 mm, length 180 cm
- NEW** 495 TIP **Fiber Optic Light Cable**, highly heat resistant, diameter 4.8 mm, length 300 cm, optimized for 3D TIPCAM®1 and TIPCAM®



## Fiber Optic Light Cable with 90° Deflection to the Light Source

- 495 NWL **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 300 cm
- 495 NWM **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 230 cm
- 495 NW **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 180 cm
- 495 NTXS **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain, diameter 2.5 mm, length 230 cm
- 495 NTW **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 2.5 mm, length 180 cm



495 NVC



495 NV/NVL/NVB

## Fiber Optic Light Cable with 90° Deflection to the Instrument

- NEW** 495 NVC **Fiber Optic Light Cable**, with 90° deflection to the instrument, very narrow radius of curvature, diameter 4.8 mm, length 300 cm
- 495 NVB **Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 4.8 mm, length 300 cm
- 495 NVL **Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 3.5 mm, length 300 cm
- 495 NV **Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 3.5 mm, length 230 cm



## Light Adaptor, angled 90°

- 495 EW **Light Adaptor**, angled 90°, diameter 4.8 mm, free rotatable, to connect with standard telescopes


Adaptor for connecting KARL STORZ fiber optic light cables with endoscopes and light sources from other manufacturers see pages TP 108-109

# Fiber Optic Light Cables

for Cold Light Fountains



495 AD/BD



495 UD



495 DC/DV

## Fiber Optic Light Cable with Straight Connector

- 495 AD** **Fiber Optic Light Cable**, diameter 3.5 mm, length 180 cm, for double incidental light radiation in film and TV and when used as demonstration devices
- 495 BD** **Fiber Optic Light Cable**, diameter 4.8 mm, length 180 cm, for double incidental light radiation in film and TV and when used as demonstration devices
- 495 UD** **Fiber Optic Light Cable**, 2x diameter 3.5 mm, length 180 cm, for simultaneous light transmission for 2 instruments
- 495 DC** **Fiber Optic Light Cable**, diameter 3.5 mm, length 320 cm, for use with DCI® Camera Heads **20 2620 30**, **20 2621 30**, **22 2600 31-3** and **22 2601 31-3**
- 495 DV** **Fiber Optic Light Cable**, diameter 2.5 mm, length 320 cm, for use with DCI® Camera Heads **20 2620 30**, **20 2621 31**, **22 2600 31-3** and **22 2601 31-3**

The given diameter is the diameter of the built-in glass fibers on the instrument side.



## Fiber Optic Light Cable with Straight Connector

- 495 NAS** **Fiber Optic Light Cable**, diameter 3.5 mm, length 230 cm
- 495 NTAS** **Fiber Optic Light Cable**, diameter 2.5 mm, length 230 cm
- 495 NWMS** **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain, diameter 3.5 mm, length 230 cm
- 495 NTXS** **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain, diameter 2.5 mm, length 230 cm

LIGHT SOURCES

### Please note:

The high light concentration at the end of the light cable causes heat to be generated in the focal point. The end of the light cable should never be placed on the patient's drape or skin as long as the cold light fountain is turned on, since the light intensity could cause burns in the patient or set the drape on fire.

### Sterilization and Disinfection:

The fiber optic light cables can be sterilized in the autoclave at 273 ° F (134 ° C). Gas sterilization is recommended and chemical disinfection is possible.

### Recommended combination: light cable with endoscope\*

Diameter light cable	Diameter endoscope
4.8 – 5.0 mm	6.5 – 12.0 mm
3.0 – 3.5 mm	3.0 – 6.5 mm
2.0 – 2.5 mm	0.8 – 2.9 mm


\*Special endoscopes may deviate.

Adaptor for connecting KARL STORZ fiber optic light cables with endoscopes and light sources from other manufacturers see pages TP 108-109

# Adaptors

for Light Sources and Light Cables

HOPKINS® Telescopes with universal connector for light cables



Manufacturer	I = Instrument CF = Cold Light Fountain LC = Light Cable	KARL STORZ				Special Applications
		Light Cable	Cold Light Fountain	Former Instrument	Current Instrument	
ACMI	I	495 W	-	-	-	ACMI connector to ACMI connector 495 Y
	CF	495 PA	-	-	-	
	LC	-	495 C	-	-	
Downs	CF	495 PA	-	-	-	Olympus Techno Fountain for KARL STORZ Cable 81495 OT, Olympus Light Cable for KARL STORZ Techno Fountain 81450/81495 OL
Dyonics	CF	495 Q	-	-	-	
Effner	CF	495 Q	-	-	-	
Fuji	CF	495 FU	-	-	-	
KLI	LC	-	495 C	-	-	
	CF	495 PA	-	-	-	
Machida	CF	495 SM	-	-	-	
Olympus	CF	495 SON	-	-	-	
	LC	-	495 D	-	-	
Pentax	CF	495 T	-	-	-	
	LC	-	495 D	-	-	
HSW	CF	495 Q	-	-	-	
	LC	-	495 E	495 N	495 N	
Stryker	I	495 ST	-	-	-	
	CF	495 PA	-	-	-	
Richard Wolf	I	495 X	-	-	-	
	CF	495 Q	-	-	-	
	LC	-	495 E	495 M	-	
Olympus/Winter & Ibe	CF	495 S	-	-	-	<b>495 S:</b> The use of KARL STORZ fiberscopes with integrated fiber optic light cable and KARL STORZ light cables with non-graduated distal connection (diameter 7 mm) requires an additional adaptor <b>495 PWS</b> to adaptor <b>495 S</b>
	LC with attached lens	-	-	495 GO	495 GO	

LIGHT SOURCES



# Adaptor Series

for HALOGEN and XENON Light Sources



- 487 A **Adaptor**, for KARL STORZ light sources in combination with ACMI light cable
- 487 M **Adaptor**, for KARL STORZ light sources in combination with Machida light cable
- 487 O **Adaptor**, for KARL STORZ light sources in combination with Olympus light cable
- 69487 KVE **Adaptor**, for KARL STORZ light sources (except TELE PACK) in combination with KARL STORZ video endoscopes or Olympus OES flexible endoscopes, with integrated heat filter, includes connector for insufflation
- 487 P **Adaptor**, for KARL STORZ light sources in combination with Pilling light cable
- 487 U **Adaptor**, for KARL STORZ light sources in combination with KARL STORZ, Wolf, Dyonics, ACMI, Pilling, V. Müller, Stryker light cable
- 487 UO **Adaptor**, for KARL STORZ light sources in combination with KARL STORZ, Wolf, Dyonics, ACMI, Olympus, V. Müller, Stryker light cable
- 487 W **Adaptor**, for KARL STORZ light sources in combination with Wolf light cable
- 487 WI **Adaptor**, for KARL STORZ light sources in combination with Winter-IBE light cable
- 487 Z **Adaptor**, for KARL STORZ light sources in combination with Zeiss light cable

Adaptor	HALOGEN LIGHT SOURCE VETERINARY	HALOGEN 250 twin	LED NOVA 150	LED NOVA 100 twin	Power LED 175 SCB	XENON 100 SCB	XENON NOVA® 175	XENON 175 SCB	XENON NOVA® 300	XENON 300 SCB	TELE PACK VET X
487 A	●	●	-	-	-	-	●	●	●	●	-
487 M	●	●	-	-	-	-	●	●	●	●	-
487 O	●	●	-	-	-	-	●	●	●	●	-
69487 KVE	-	-	-	-	-	-	●	●	●	●	-
487 P	●	●	-	-	-	-	●	●	●	●	-
487 U	●	●	-	-	-	-	●	●	●	●	-
487 UO	●	●	-	-	-	-	●	●	●	●	-
487 W	●	●	-	-	-	-	●	●	●	●	-
487 WI	●	●	-	-	-	-	●	●	●	●	-
487 Z	●	●	-	-	-	-	●	●	●	●	-



# SMARTCART Equipment Cart Series



2-12

EQUIPMENT CARTS

# SMARTCART 2 – The New Generation of KARL STORZ Equipment Carts

## Special Features:

- Choose from a broad range of sizes and features
- The carts can be easily customized to any given requirements with the help of a wide variety of useful accessories – for example, with camera holders, irrigator rod, special storage areas for keyboard or mouse, etc.
- Convenient and functional cable management within the double rear wall with integrated electrical sub-distributors through which cables and tubes can be safely passed.
- Vertical cable managers ensure tidy cabling right up to the appliance.
- Central power supply via the powerbox, which contains all the mains and grounding terminals while offering enough space to hold an isolation transformer.
- Smooth, readily accessible, easy-to-clean surfaces



TFT monitor holders with VESA 75/100 adaption



Detailed view with open rear panel



Lateral push handles for easy handling



Optional integrable isolation transformer with grounding terminal for maximum patient safety



## Equipment Carts <sup>NEW</sup>



29005 DRB

### 29005 DRB Equipment Cart

including:

**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes (rear), 3 fixed shelves, 1 with push handles, mains switch in side boom, 1 drawer unit with lock, integrated cable conduits in both vertical beams, 1 set of non-sliding stands for units, double rear panel with integrated electrical subdistributors (with 12 sockets), holder for power supplies, potential earth connectors and cable winding on the outside, 1 camera holder, 2 equipment rails sidewise

Dimensions:

Equipment Cart: 730 x 1490 x 716 mm (w x h x d),

Shelf: 630 x 480 mm (w x d),

Caster diameter: 150 mm

#### Accessories:

29005 MZD **Monitor Holding Arm**, height and side adjustable, swiveling and tilting, central mount, swivel range 360°, reach 300 mm, loading capacity max. 15 kg, with monitor holder VESA 75/100, with integrated cable conduit, for use with Equipment Cart 29005 DRB

29005 SZD **Monitor Swivel Arm**, height and side adjustable, can be mounted on the left or on the right side, swivel range 180°, reach 600 mm, for side boom as of 800 mm, loading capacity max. 14 kg, with monitor holder VESA 75/100, for Equipment Cart 29005 DRB



29005 DRS

### 29005 DRS Equipment Cart

including:


**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes (front), 3 fixed shelves, 1 with push handles, mains switch in side boom, 2 push handles on the side, 1x drawer unit with lock, integrated cable conduits in both side booms, 1x set of non-sliding stands for units, double rear panel with integrated electrical subdistributors with 12 sockets, holder for power supplies, potential earth connectors and cable manager on the outside, 1x camera mount, 2x equipment rails on the side,

Dimensions:

Equipment Cart: 530 x 1455 x 645 mm (w x h x d),

Shelf: 430 x 480 mm (w x d),

Caster diameter: 125 mm



29005 LAP

## 29005 LAP    **Equipment Cart**

including:

**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, 3 fixed shelves, 1 with push handles, mains switch in side boom, 1x drawer unit with lock, integrated cable conduit in both booms, 3x horizontal cable conduits, 2 with 4x subdistribution board for cold device sockets, 1 with cable manager, 1x set of non-sliding stands for units, 1x monitor holding arm VESA 75/100, 1x camera mount, 8x mains cords 50 cm, 2x mains cord 200 cm, 2x equipment rails on the side, isolation transformer 230 VAC (50/60 Hz) 2000 VA with 8 plug-in slots, 8 grounding plugs and earth leakage monitor, 1x CO<sub>2</sub> bottle holder, max. diameter 155 mm


Dimensions:

Equipment cart: 730 x 1490 x 716 mm (w x h x d),

Shelf: 630 x 480 mm (w x d),

Caster diameter: 150 mm

## Equipment Carts



29005 GU

### 29005 GU Equipment Cart

including:

**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, 2 fixed shelves, 1 with push handles, main switch in side boom, 1x drawer unit with lock, integrated cable conduit in both booms, 3x horizontal cable conduits, 2 with 4x subdistribution board for cold device sockets, 1 with cable manager, 1x set of non-sliding stands for units, 1x monitor holding arm VESA 75/100, 1x camera mount, 8x mains cords 50 cm, 2x mains cord 200 cm, isolation transformer 230 VAC (50/60 Hz) 2000 VA with 8 plug-in slots, 8 grounding plugs and earth leakage monitor

Dimensions:

Equipment cart: 730 x 1490 x 716 mm (w x h x d),

Shelf: 630 x 480 mm (w x d),

Caster diameter: 150 mm

### 29005 BAS Equipment Cart

including:

**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, 2 fixed shelves, 1 with push handles, main switch in side boom, 1x drawer unit with lock, integrated cable conduit in both booms, 3x horizontal cable conduits, 1 with cable manager, 1x set of non-sliding stands for units, 1x monitor holding arm VESA 75/100, 1x camera mount

Dimensions:

Equipment cart: 730 x 1490 x 716 mm (w x h x d),

Shelf: 630 x 480 mm (w x d),

Caster diameter: 150 mm



29005 BKW

### 29005 BKW **Equipment Cart**

including:

**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, main switch in side boom, 1x drawer unit with lock, integrated cable conduit in both vertical beams, 2x horizontal cable conduits, one with cable manager, 1x set of non-sliding stands for units

Dimensions:

Equipment cart: 530 x 1090 x 645 mm (w x h x d),

Shelf: 430 x 480 mm (w x d),

Caster diameter: 125 mm



## Equipment Carts



29005 GI

### 29005 GI Equipment Cart

including:

**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, 2 fixed shelves, 1 with push handles, main switch in side boom, 1x drawer unit with lock, integrated cable conduit in both vertical beams, 2x horizontal cable conduits, 1 with cable manager, 1x set of non-sliding stands for units, 1x pull-out keyboard at the front, 1x video endoscope intake, 1x video endoscope plug intake, 1x tube clip, 1x monitor holder VESA 75/100, 4x equipment rails, 6x mains cords, 3x 200 cm and 3x 100 cm, isolation transformer 230 VAC (50/60 Hz) 1200 VA with 6 plug-in slots and 6 grounding plugs

Dimensions:

Equipment cart: 530 x 1200 x 645 mm (w x h x d),  
Shelf: 430 x 480 mm (w x d),  
Caster diameter: 125 mm

### 29005 BGI Equipment Cart

including:


**basic equipment cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes, 2 fixed shelves, 1 with push handles, main switch in side boom, 1x drawer unit with lock, integrated cable conduit in both vertical beams, 2x horizontal cable conduits, 1 with cable manager, 1x set of non-sliding stands for units

Dimensions:

Equipment cart 530 x 1200 x 645 mm (w x h x d),  
Shelf: 430 x 480 mm (w x d),  
Caster diameter: 125 mm

# Equipment Carts

Accessories



29005 SZ **Monitor Swivel Arm**, height and side adjustable, can be mounted on the left or on the right side, swivel range 180°, reach 600 mm, for side boom as of 800 mm, loading capacity max. 14 kg, with monitor holder VESA 75/100, for Equipment Cart 29005 LAP/GU/BAS



29005 MZ **Monitor Holding Arm**, height and side adjustable, swiveling and tilting, central mount, swivel range 190°, reach 300 mm, loading capacity max. 15 kg, with monitor holder VESA 75/100, for Equipment Cart 29005 LAP/GU/BAS




29005 MSK **Monitor Holding Arm**, height and side adjustable, tilting, can be mounted either on the left or on the right side, swivel range 190°, reach 300 mm, loading capacity max. 15 kg, with monitor holder VESA 75/100, for Equipment Carts 29005 xx

29005 MSL **Same**, reach 600 mm, loading capacity max. 14 kg, for equipment carts 29005 LAP/GU/BAS/DRB



29005 MH **Monitor Arm**, height adjustable, swiveling and tilting, central mount, swivel range approx. 210°, loading capacity max. 18 kg, with monitor holder VESA 75/100, for Equipment Carts 29005 HNO/GI/AN/BGI/BWK/DRS

29005 MHB **Same**, for Equipment Cart 29005 BAS/DRB and 29003 LCG



29005 MHS **Monitor Holder**, lateral, with VESA 75/100 monitor holder, including monitor module and multifunctional holder, for use with Equipment Carts 29005 xx



29003 GK **Accessory Side Shelf**, for mounting on side boom, with 2 arms, load capacity 10 kg, floor space: 310 x 310 mm (w x d)

# Equipment Carts

Accessories



29003 GS **Standard Equipment Rail**, for mounting on side with equipment carts of series 29003 and 29005



29003 KK **Cable Conduit Horizontal**, width 630 mm, for Equipment Cart 29005 LAP/GU/BAS



29005 PB **Power Box**, socket board with 12 mains sockets, 12 grounding plugs, for Equipment Cart 29005 BAS

29003 PBK **Power Box**, socket board with 6 mains sockets, 6 grounding plugs for Equipment Carts 29005 HNO/AN/BGI/BWK



29005 TBG **Isolation Transformer**, 230 VAC, 2000 VA, with 8 special mains socket expulsion fuses, 8 grounding plugs, for Equipment Carts 29005 BAS/DRB/LAP/GU








29003 TBK **Isolation Transformer**, 230 VAC, 1200 VA, 6 special mains socket expulsion fuses, 6 grounding plugs, for Equipment Cart 29005 DRS/BGI/BWK/GI/HNO/AN



29003 IW **Earth Leakage Monitor**, for use with Isolation Transformer 29003 TBK and 29005 TBG

# Equipment Carts

Accessories

	<p>400 D <b>Mains Cord</b>, mains connector 90° angulated, unit connector straight, length 400 cm</p> <p>400 DM <b>Same</b>, length 250 cm</p>
	<p>29005 KKK <b>Mains Cord</b>, length 50 cm, with cold device plug and socket, with UL approval</p> <p>29005 KKM <b>Same</b>, length 100 cm</p> <p>29005 KKL <b>Same</b>, length 200 cm</p>
	<p>28003 TE <b>Mains Cord Adaptor</b>, for devices with European plugs, for use with Isolation Transformer 29003 TBG/TBK/TBR and 29005 TBG</p>
	<p>28003 TO <b>Mains Cord Adaptor</b>, for devices with fixed mains cord with shockproof plug, for use with Isolation Transformer 29003 TBG/TBK/TBR and 29005 TBG</p>
	<p>29005 KSK <b>Mains Cord Adaptor</b>, for devices with cold device sockets to power box with shockproof socket</p>
	<p>20 0101 70 <b>Potential Equipotential Cable</b>, length 10 m</p> <p>20 0102 70 <b>Potential Equalization Cable</b>, length 500 cm</p> <p>20 0103 70 <b>Potential Equalization Cable</b>, length 200 cm</p> <p>20 0105 70 <b>Potential Equipotential Cable</b>, length 150 cm</p> <p>20 0106 70 <b>Potential Equipotential Cable</b>, length 100 cm</p>
	<p>29005 PA-HNO <b>Potential Equalization Cabling</b>, for configuring equipment carts, including Potential Equipotential Cables 2x 20 0105 70, 1x 20 0106 70 and 1x 20 0102 70 for 3 units</p> <p>29005 PA-GI <b>Potential Equalization Cabling</b>, for configuring equipment carts, including Potential Equipotential Cables 3x 20 0105 70, 1x 20 0106 70 and 1x 20 0102 70 for 4 units</p> <p>29005 PA-AN <b>Potential Equalization Cabling</b>, for configuring equipment carts, including Potential Equipotential Cables 1x 20 0103 70 and 1x 20 0102 70 for 1 unit</p> <p>29005 PA-LAP <b>Potential Equalization Cabling</b>, for configuring equipment carts, including Potential Equipotential Cables 5x 20 0103 70, 2x 20 0105 70 and 1x 20 0102 70 for 7 units</p> <p>29005 PA-GU <b>Potential Equalization Cabling</b>, for configuring equipment carts, including Potential Equipotential Cables 4x 20 0103 70, 1x 20 0105 70 and 1x 20 0102 70 for 5 units</p>

# Equipment Carts

Accessories



- 29005 BK **Shelf**, load capacity max. 50 kg, dimensions: 430 x 480 mm (w x d), for Equipment Cart HNO/AN/GI/BWK/BGI/DRS
- 29005 BG **Shelf**, load capacity max. 50 kg, dimensions: 630 x 480 mm (w x d), for Equipment Cart 29005 LAP/GU/BAS/DRB



- 29005 TAK **Keyboard Tray**, for mounting under Shelf 29005 BK, load capacity 10 kg, dimensions: 430 x 480 mm (w x d)
- 29005 TAG **Keyboard Tray**, for mounting under Shelf 29005 BG, load capacity 10 kg, dimensions: 630 x 480 mm (w x d)



- 29005 BGG **Shelf with Handles**, load capacity max. 50 kg, dimensions: 630 x 480 mm (w x d), for Equipment Cart 29005 LAP/GU/BAS/DRB
- 29005 BGK **Shelf with Handles**, load capacity max. 50 kg, dimensions: 430 x 480 mm (w x d), for Equipment Cart 29005 HNO/AN/BWK/GI/BGI/DRS



- 29005 SKS **Drawer Unit with Lock**, load capacity max. 50 kg, with cable slot, dimensions: 430 x 125 x 480 mm (w x h x d), for Equipment Cart 29005 HNO/AN/BWK/GI/BGI/DRS
- 29005 SG **Drawer Unit with Lock**, load capacity max. 50 kg, with cable slot, dimensions: 630 x 125 x 480 mm (w x h x d), for Equipment Cart 29005 LAP/GU/BAS/DRB
- 29005 SKO **Drawer Unit with Lock**, load capacity max. 50 kg, dimensions: 430 x 125 x 480 mm (w x h x d), for Equipment Cart 29005 HNO/AN/BWK/GI/BGI/DRS



- 29005 KW **Cable Manager**, for camera cables, can be installed in T-grooves of leg columns, package of 2



- 29003 EW **Holder for Single-Use Container**, for mounting to Irrigator Rod 29003 S/SV

# Equipment Carts

Accessories

	<p>29005 KH <b>Camera Holder</b>, adjustable height, can be installed in T-grooves of leg columns, suitable for all KARL STORZ endoscopy cameras</p>
	<p>29005 KHN <b>Camera Mount</b>, for mounting to Standard Equipment Rail 29003 GS, suitable for all KARL STORZ endoscopy cameras</p>
	<p>29003 FH <b>Bottle Holder</b>, for CO<sub>2</sub> bottles with bottom diameter of max. 155 mm 29003 FHG <b>Bottle Holder</b>, for CO<sub>2</sub> bottles with bottom diameter of max. 200 mm</p>
	<p>69003 FT <b>Carrier Set</b>, for flexible endoscopes for mounting at equipment rail, including holding device for 2 endoscopes, 2 acrylic tubes and 2 equipment rails</p>
	<p>29003 FL <b>Holder</b>, for fitting the Bottle Stand Holder 20 3000 33</p>
	<p>69003 HH <b>Otoscope Camera Head Holder</b>, for combining Otoscope 67260 AA with <b>non-autoclavable</b> TELECAM camera head and Veterinary Video Camera III camera head</p>
	<p>29003 S <b>Irrigator Rod</b>, with 4 bottle hooks, for mounting to equipment cart series 29005, incl. two multifunctional holders, length 130 cm 29003 SV <b>Irrigator Rod with Holder</b>, height adjustable, with 4 bottle hooks, for mounting to equipment cart series 29005, incl. two multifunctional holders, length 120 – 200 cm</p>
	<p>29003 MH <b>Multifunction Holder</b>, can be installed in T-grooves of upright columns, for integrating accessories and centering pins</p>

# Equipment Carts

## Accessories



- 29005 SK **Tube Clip**, for attaching suction and irrigation tubes to Keyboard Shelf 29005 TAK/TAG
- 29005 SKN **Tube Clip**, to attach to suction and irrigation tubes, for mounting to Equipment Rail 29003 GS



- 13991 ET **Video Endoscope Holder**, for mounting to 2 Equipment Rails 29003 GS including:  
**Video Endoscope Holder**  
**Drip Container**, with clamp for fixation to the holding rod
- 13991 DET **Video Endoscope Holder**, for mounting on 2 Equipment Rails 29003 GS including:  
**Double Video Endoscope Holder**  
**2 Drip Containers**, each equipped with clamp for fixation to the holding rod



- 29005 AVS **Video Endoscope Plug Intake**, for mounting to Irrigator Rod 29003 S/SV and Video Endoscope Holder 13991 ET and 13991 DET



- 29005 PH **Pump Holder**, for attaching suction pump UNIMAT® 30 25 3200 01 to equipment cart, for mounting to 2 Equipment Rails 29003 GS



- 29005 HFS **Holder for Two-Pedal Footswitch**, for mounting to equipment cart side boom model 29005 xx, for use with Footswitch 20 0103 30, 20 0104 30, 20 0114 30, 20 0124 30, 20 0126 30, 20 0128 31 and 20 0133 30

# Equipment Carts

Accessories



29005 HFH **Holder for Two-Pedal Footswitch**, suspended, for mounting to equipment rails, for use with Footswitch **20 0138 30** and **20 0138 31**

29005 HFE **Single Pedal Holder**, suspended, for mounting to equipment rails, for use with One-Pedal Footswitch **20 0138 32**




29005 DFH **Holder for Two-Pedal Footswitch**, for mounting to equipment rail, for use with Footswitch **20 0143 30**, **20 0144 30**, **20 0154 30**, **20 0164 30**, **20 0166 30**, **20 0168 31**, **20 0178 30** and **20 0178 31**

29005 EFH **Single Pedal Holder**, for mounting to equipment rail, for use with Footswitch **20 0141 30**, **20 0142 30** and **20 0178 32**



29005 AK **Wire Basket**, for accessories, for mounting to equipment rail, dimensions: 300 x 100 x 200 mm (w x h x d)



29005 NBB **Fixation for Power Supply**, double, for mounting 2 power supplies, for attachment to Shelf **29005 BG/BGG**

29005 NBS **Fixation for Power Supply**, double, for mounting 2 power supplies, for attachment to Shelf **29005 BK/BGK**



29005 SU **Equipment Cart Protective Cover**, with 2 zippers at the front, dimensions: 740 x 1450 x 730 mm (w x h x d), for use with Equipment Carts **29003 LCG** and **29005 LAP/GU/BAS/DRB**



## Equipment Cart



29003 LC

### 29003 LC Equipment Cart

Dimensions:  
Equipment cart: 730 x 1450 x 660 mm (w x h x d),  
Shelf: 630 x 480 mm (w x d),  
Caster diameter: 125 mm

including:

**Equipment Cart**, rides on 4 antistatic dual wheels, 2 equipped with locking brakes (front), 2 fixed shelves, 1 drawer unit with lock, integrated cable conduits in both vertical beams, mains switch in vertical beam, 1 set of non-sliding stands, 1 camera holder

**Power Box**, socket board with 12 sockets and 12 potential earth connectors

#### Accessories:

29003 FH **Bottle Holder**, for CO<sub>2</sub> bottles with bottom diameter of max. 155 mm

29003 FT **Fiberscope Carrier Set**, for mounting to equipment rail, including holding device for 2 fiberscopes, 2 acrylic tubes and 2 Equipment Rails 29003 GS

29003 S **Irrigator Rod**, with 4 bottle hooks, for mounting to equipment cart series 29005, incl. two multifunctional holders, length 130 cm



20 0200 23

### 20 0200 23 Mobile Stand

height adjustable, rides on 5 antistatic casters, 2 equipped with locking brakes, tray made of plastic-coated V2A, for use with cold light fountains, laparoflators, pump and motor systems

Dimensions:  
Tray: 330 x 330 mm (w x d),  
Working height: 900 – 1500 mm

#### Accessories:

20 0200 31 **Tray**, made of V2A, plastic-coated, for use with Mobile Stand 20 0200 23, dimensions: 260 x 170 x 55 mm (w x h x d)

20 0200 32 **Holder for Irrigation Liquid Bottle**, with 4 hooks for bottles, for use with Mobile Stand 20 0200 23



**20 0200 81 TROLL-E Mobile Stand**

rides on 4 antistatic dual wheels, 2 equipped with locking brakes, for mounting TFT monitors with VESA 75/100 connection, integrated cable conduit in vertical beam, 1 shelf, cable manager, load capacity for monitor: max. 15 kg

Dimensions:  
Mobile stand: 670 x 1660 x 670 mm (w x h x d),  
Shelf: 455 x 350 mm (w x d), caster diameter: 100 mm

**Note:** Trolley is delivered unassembled.



**20 0200 80 SHORT-E HF Cart**

rides on 4 antistatic dual wheels, 2 equipped with locking brakes, integrated cable conduit in vertical beam, 1 HF shelf for AUTOCON® II 400 SCB, with hook for storing cables, cable manager, max. load capacity 51 kg,

Dimensions:  
HF cart: 670 x 1020 x 670 mm (w x h x d),  
HF shelf: 455 x 375 mm (w x d), caster diameter: 100 mm

**Note:** Trolley is delivered unassembled.



**20 0200 82 Wall Mount MOUNT-E™**

for mounting monitors with VESA 75/100 connection, integrated cable conduit in vertical beam, max. load capacity 20 kg,

Dimensions: 140 x 750 x 260 mm (w x h x d)

# Equipment Carts E Series <sup>NEW</sup>

Accessories



**20 0200 46** **Keyboard Shelf**, max. capacity 5 kg, for Mobile Stand **20 0200 81** and Wall Mount **20 0200 82**, dimensions: 465 x 185 mm (w x d)



**20 0200 47** **Shelf**, incl. installation accessories, max. capacity 12 kg, dimensions: 455 x 350 mm (w x d), for use with subracks for Mobile Stands **20 020060** and **20 020061**



**20 0200 49** **Equipment Rail**, for mounting to central boom on HF Cart **20 0200 80** and Mobile Stand **20 0200 81**, incl. installation accessories dimensions: 350 x 25 x 10 mm (w x h x d)



**29005 HFH** **Holder for Two-Pedal Footswitch**, suspended, for mounting to equipment rails, for use with Footswitch **20 0138 30** and **20 0138 31**

**29005 HFE** **Single Pedal Holder**, suspended, for mounting to equipment rails, for use with One-Pedal Footswitch **20 0138 32**



**29005 DFH** **Holder for Two-Pedal Footswitch**, for mounting to equipment rail, for use with Footswitch **20 0143 30**, **20 0144 30**, **20 0154 30**, **20 0164 30**, **20 0166 30**, **20 0168 31**, **20 0178 30** and **20 0178 31**


**29005 EFH** **Single Pedal Holder**, for mounting to equipment rail, for use with Footswitch **20 0141 30**, **20 0142 30** and **20 0178 32**



**29005 AK** **Wire Basket**, for accessories, for mounting to equipment rail, dimensions: 300 x 100 x 200 mm (w x h x d)



# Components Spare Parts



2-12

COMPONENTS  
SPARE PARTS

# IMAGE 1 HUB™ HD Camera Control Units SCB

Components / Spare Parts

Catalog page

**22201011U110 IMAGE 1 HUB™ HD Camera Control Unit SCB, with ICM module**

TP 10

22201020-110 **IMAGE 1 HUB™ HD Camera Control Unit SCB, with ICM module**

400 A **Mains Cord**

536 MK **BNC/BNC Video Cable**

547 S **S-Video (Y/C) Connecting Cable**

20203270 **Special RGBS Connecting Cable**

20221070 **Connecting Cable, for controlling peripheral units**

20040089 **DVI-D Connecting Cable**

20090170 **SCB Connecting Cable**

20200231U **Keyboard, with US English character set**

20040281 **KARL STORZ USB Flash Drive, 4 GB**



Components / Spare Parts

Catalog page

**22201011U112 IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI/ICM module**

TP 10

22201020-112 **IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI/ICM module**

400 A **Mains Cord**

536 MK **BNC/BNC Video Cable**

547 S **S-Video (Y/C) Connecting Cable**

20203270 **Special RGBS Connecting Cable**

20221070 **Connecting Cable, for controlling peripheral units**

20040089 **DVI-D Connecting Cable**

20090170 **SCB Connecting Cable**

20200231U **Keyboard, with US English character set**

20040281 **KARL STORZ USB Flash Drive, 4 GB**



Components / Spare Parts

Catalog page

**22201011U1 IMAGE 1 HUB™ HD Camera Control Unit SCB**

TP 11

22201020-1 **IMAGE 1 HUB™ HD Camera Control Unit SCB**

400 A **Mains Cord**

536 MK **BNC/BNC Video Cable**

547 S **S-Video (Y/C) Connecting Cable**

20203270 **Special RGBS Connecting Cable**

20221070 **Connecting Cable, for controlling peripheral units**

200400 89 **DVI-D Connecting Cable**

20090170 **SCB Connecting Cable**

20200231U **Keyboard, with US English character set**



# IMAGE 1 HUB™ HD, IMAGE1 HD, 3D Camera Control Units SCB

Components / Spare Parts

Catalog page

**22201011U102 IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI module**

TP 11

22201020-102 **IMAGE 1 HUB™ HD Camera Control Unit SCB, with SDI module**

400 A **Mains Cord**

536 MK **BNC/BNC Video Cable**

547 S **S-Video (Y/C) Connecting Cable**

20203270 **Special RGBS Connecting Cable**

20221070 **Connecting Cable, for controlling peripheral units**

20040089 **DVI-D Connecting Cable**

20090170 **SCB Connecting Cable**

20200231U **Keyboard, with US English character set**



Components / Spare Parts

Catalog page

**22202011U110 IMAGE1 HD Camera Control Unit SCB, with ICM module**

TP 13

22202020-110 **IMAGE1 HD Camera Control Unit, with ICM module**

400 A **Mains Cord**

20221070 **Connecting Cable, for controlling peripheral units**

20040089 **DVI-D Connecting Cable**

20090170 **SCB Connecting Cable**

20200231U **Keyboard, with US English character set**

20040281 **KARL STORZ USB Flash Drive, 4 GB**



Components / Spare Parts

Catalog page

**22202011U1 IMAGE1 HD Camera Control Unit SCB**

TP 13

22202020-1 **IMAGE1 HD Camera Control Unit SCB**

400 A **Mains Cord**

20090170 **SCB Connecting Cable**

20040089 **DVI-D Connecting Cable**

20221070 **Connecting Cable, for controlling peripheral units**

20200231U **Keyboard, with US English character set**



Components / Spare Parts

Catalog page

**22204011U114 3D Camera Control Unit SCB, with ICM Module**

TP 23

22204020 **3D Camera Control Unit SCB, with ICM module**

400 A **Mains Cord**

536 MK **BNC/BNC Video Cable**

547 S **S-Video (Y/C) Connecting Cable**

20221070 **Connecting Cable, for controlling peripheral units**

20040089 **DVI-D Connecting Cable**

20090170 **SCB Connecting Cable**

20040245US **Silicone Keyboard, with US English character set**

20040281 **KARL STORZ USB Flash Drive, 4 GB**

9800 G **KARL STORZ 3D Glasses**



# TRICAM® SL II SCB, TELE PACK VET X, TELECAM SL II, Veterinary Video Camera III

Components / Spare Parts

Catalog page

- 20223011U1**    **TRICAM® SL II Camera Control Unit SCB**
- 20223020-1    **TRICAM® SL II Camera Control Unit SCB**
- 400 A            **Mains Cord**
- 20200231U    **Keyboard, with US English character set**
- 20221070      **Connecting Cable, for controlling peripheral units**
- 536 MK         **BNC/BNC Video Cable**
- 547 S           **S-Video (Y/C) Connecting Cable**
- 20203270      **Special RGBS Connecting Cable**
- 20090170      **SCB Connecting Cable**
- 20040082      **DV Connecting Cable, 6 pin to 4 pin**



TP 29

Components / Spare Parts

Catalog page

- 69045001**    **TELE PACK VET X**
- 69045020      **TELE PACK VET X**
- 20040280      **KARL STORZ USB Flash Drive, 4 GB**
- 400 F            **Mains Cord**
- 400 FUS         **Mains Cord, US plug**
- 9104250        **Case**
  
- 69045002**    **Same, without Case**



TP 35

Components / Spare Parts

Catalog page

- 20213011U**    **TELECAM SL II Camera Control Unit**
- 20213020      **TELECAM SL II Camera Control Unit**
- 400 A            **Mains Cord**
- 20221070      **Connecting Cable, for controlling peripheral units**
- 536 MK         **BNC/BNC Video Cable**
- 547 S           **S-Video (Y/C) Connecting Cable**
- 20040082      **DV Connecting Cable, 6 pin to 4 pin**
- 20200231U    **Keyboard, with US English character set**



TP 39

Components / Spare Parts

Catalog page

- 69236001**    **Veterinary Video Camera III, PAL**
- 69236020      **Veterinary Video Camera III Camera Control Unit, color systems PAL/NTSC**
- 400 A            **Mains Cord**
- 20221070      **Connecting Cable, for controlling peripheral units**
- 536 MK         **BNC/BNC Video Cable**
- 547 S           **S-Video (Y/C) Connecting Cable**
- 20200043      **C-MOUNT Lens, f = 38 mm, soakable**
- 69215034      **VET C-MOUNT Camera Head, PAL**
  
- 69236101**    **Same, NTSC**
- 69215134      **VET C-MOUNT Camera Head, NTSC**



TP 40



# Veterinary Video Camera III, C-HUB™, C-MAC® Video Laryngoscope

Components / Spare Parts

Catalog page

- 69236011**    **Veterinary Video Camera III Camera Control Unit**
- 69236020**    **Veterinary Video Camera III Camera Control Unit**
- 400 A**        **Mains Cord**
- 20221070**    **Connecting Cable**, for controlling peripheral units
- 536 MK**      **BNC/BNC Video Cable**
- 547 S**        **S-Video (Y/C) Connecting Cable**

TP 41



Components / Spare Parts

Catalog page

- 20290101**    **C-HUB™ Camera Control Unit**
- 20290120-PS** **C-HUB™ Power Supply**
- 547 S**        **S-Video (Y/C) Connecting Cable**
- 20200073**    **USB Connecting Cable**

TP 54



Components / Spare Parts

Catalog page

- 8402 ZX**      **Monitor for CMOS endoscopes**

TP 55



Spare Parts



**40150031**  
**Power Supply**, for charging and starting up C-MAC® 8401 ZX and 8402 ZX, with 4 adaptors for Europe, Great Britain, USA/Japan and Australia, power supply 100 – 240 VAC, 50/60 Hz



**ET38-1717707**  
**Europlug Adaptor**, to mains adaptor for C-MAC® video laryngoscope



**ET38-1717715**  
**Plug Adaptor USA/Japan**, to mains adaptor for C-MAC® video laryngoscope



**ET38-1717618**  
**Plug Adaptor UK Standard**, to mains adaptor for C-MAC® video laryngoscope



**ET38-1800496**  
**Plug Adaptor Australia**, to mains adaptor for C-MAC® video laryngoscope

# AIDA compact NEO advanced/standard, 19" and 26" HD Monitors

Components / Spare Parts

Catalog page

**20040913-EN/20040912-EN KARL STORZ AIDA® compact NEO advanced/standard**

TP 61

- 20046120-A **KARL STORZ AIDA® control NEO**, with integrated DVD/CD writer, Blu-ray reading possible
- 20040577-V02 **AIDA compact HD Frame Grabber Card**
- 20090234U **PS/2 Compact Keyboard**, US version, with cover
- 20040402-20 **KARL STORZ AIDA® compact NEO Software**, with voice control and software protection
- 20040275 **KARL STORZ USB Flash Drive**
- 20221070 **Connecting Cable**, for controlling peripheral units
- 536 MK **BNC/BNC Video Cable**
- 536 MKD **BNC Connecting Cable**
- 20040086 **DVI-D Connecting Cable**
- 400 A **Mains Cord**
- 20040087 **Connecting Cable**, Mini-DIN plug to BNC socket



Components / Spare Parts

Catalog page

- 9619 NB 19" HD Monitor**
- 9626 PS **Monitor Power Supply**
- 20040089 **DVI-D Connecting Cable**
- 536 MK **BNC/BNC Video Cable**

TP 71



Components / Spare Parts

Catalog page

- 9626 NB 26" HD Monitor**
- 9626 NG **Monitor**
- 9626 PS **Monitor Power Supply**
- 20040089 **DVI-D Connecting Cable**
- 536 MK **BNC/BNC Video Cable**

TP 71



Components / Spare Parts

Catalog page

- 9626 NB-2 26" HD Monitor**
- 9626 NG **Monitor**
- 9626 P **Monitor Power Supply**
- 20040089 **DVI-D Connecting Cable**
- 536 MK **BNC/BNC Video Cable**

TP 71



COMPONENTS  
SPARE PARTS

2-12

# 19" TFT Flat Screens, 15" and 19" Touch Screens

Components / Spare Parts Catalog page

<b>W29519 NB</b>	<b>19" KARL STORZ TFT Flat Screen</b>	TP 64
9519 NG	<b>19" TFT Flat Screen</b>	
9419 PS	<b>Power Supply</b>	
400 A	<b>Mains Cord</b>	
536 MK	<b>BNC/BNC Video Cable</b>	
547 S	<b>S-Video (Y/C) Connecting Cable</b>	
20040089	<b>DVI-D Connecting Cable</b>	
20040372	<b>SVGA Monitor Cable</b>	

Components / Spare Parts Catalog page

<b>20090401</b>	<b>15" KARL STORZ Touch Screen</b>	TP 66, 77
20090331	<b>15" Touch Screen</b> , including mains adaptor, CD driver	
20090386	<b>Touch Pen</b>	
20090483	<b>Adaptor, VESA 75</b>	



Components / Spare Parts Catalog page

<b>20090403</b>	<b>15" KARL STORZ Touch Screen</b>	TP 66, 77
20090331	<b>15" Touch Screen</b> , including mains adaptor, CD driver	
20090386	<b>Touch Pen</b>	
20090486	<b>Pedestal</b>	



Components / Spare Parts Catalog page

<b>20090405</b>	<b>19" KARL STORZ Touch Screen</b>	TP 66, 77
20090437	<b>19" Touch Screen</b> , including mains adaptor, CD driver	
20090386	<b>Touch Pen</b>	



Components / Spare Parts Catalog page

<b>20090406</b>	<b>19" KARL STORZ Touch Screen</b>	TP 66, 77
20090437	<b>19" Touch Screen</b> , including mains adaptor, CD driver	
20090386	<b>Touch Pen</b>	
20090487	<b>Pedestal</b>	



# 15" and 19" HD Monitors, 24" and 26" HD WIDEVIEW Monitors

Components / Spare Parts Catalog page

<b>9515 NB</b>	<b>15" Monitor</b>	TP 75
9415 NG	<b>15" Monitor</b>	
9526 PS	<b>Power Supply, flat screens</b>	
20040089	<b>DVI-D Connecting Cable</b>	
547 S	<b>S-Video (Y/C) Connecting Cable</b>	

Components / Spare Parts Catalog page

<b>9519 NB</b>	<b>19" HD Monitor</b>	TP 75
9419 NG	<b>19" HD Monitor</b>	
9419 PS	<b>Power Supply</b>	
400 A	<b>Mains Cord</b>	
20040089	<b>DVI-D Connecting Cable</b>	
536 MK	<b>BNC/BNC Video Cable</b>	


Components / Spare Parts Catalog page

<b>9524 NB</b>	<b>24" HD WIDEVIEW Monitor</b>	TP 75
9524 NG	<b>24" HD WIDEVIEW Monitor</b>	
9526 PS	<b>Power Supply</b>	
400 A	<b>Mains Cord</b>	
20040089	<b>DVI-D Connecting Cable</b>	
536 MK	<b>BNC/BNC Video Cable</b>	




Components / Spare Parts Catalog page

<b>9526 NBL</b>	<b>26" HD Monitor with LED Backlight</b>	TP 75
9526 NGL	<b>26" HD Monitor</b>	
9526 PS	<b>Power Supply</b>	
400 A	<b>Mains Cord</b>	
20040089	<b>DVI-D Connecting Cable</b>	
536 MK	<b>BNC/BNC Video Cable</b>	



Components / Spare Parts Catalog page

<b>9726 NB</b>	<b>26" HD Monitor (Highbright with LED Backlight)</b>	TP 75
9726 NG	<b>26" HD Monitor</b>	
9526 PS	<b>Power Supply</b>	
400 A	<b>Mains Cord</b>	
20040089	<b>DVI-D Connecting Cable</b>	
536 MK	<b>BNC/BNC Video Cable</b>	



# ZeroWire Video Transmission Set, Cold Light Fountains VETERINARY and HALOGEN

Components / Spare Parts

Catalog page

**9500 ZW** ZeroWire Video Transmission Set

TP 79

9500 TXG

**Transmitter**

9500 RXG

**Receiver**

9500 ZW-6

**Radiance Y-Cable 19"/26"**

9500 ZW-7

**VESA Holder 24"/26"**

9500 ZW-8

**VESA Holder 19"/23"**

9500 ZW-2

**DVI-D Cable**, for CCU-Tx and Monitor-Rx connection

96206173E


**Instruction Manual**



9500 TXG/RXG



9500 ZW-6



9500 ZW-7/ZW-8



9500 ZW-2

Components / Spare Parts

Catalog page

**69111001** Cold Light Fountain VETERINARY

TP 87

69111020

**Cold Light Fountain VETERINARY**, power supply 230 VAC, 50/60 Hz

400 A

**Mains Cord**

**69111001 C** Same, power supply 115 VAC, 50/60 Hz



Components / Spare Parts

Catalog page

**20113301** Cold Light Fountain HALOGEN 250 twin

TP 88

20113320

**Cold Light Fountain HALOGEN 250 twin**, power supply 100/120/230/240 VAC, 50/60 Hz

400 A

**Mains Cord**



# Cold Light Fountains LED

Components / Spare Parts

Catalog page

- 20161201 Cold Light Fountain LED NOVA 150**
- 20161220 **Cold Light Fountain LED NOVA 150,**  
power supply 100 – 240 VAC, 50/60 Hz
- 400 A **Mains Cord**



TP 93

Components / Spare Parts

Catalog page

- 20161101 Cold Light Fountain LED NOVA 100 twin**
- 20161120 **Cold Light Fountain LED NOVA 100 twin,**  
with two LED lamps and two light outlets,  
power supply 110 – 240 VAC, 50/60 Hz
- 400 A **Mains Cord**



TP 94

Components / Spare Parts

Catalog page

- 20161401-1 Cold Light Fountain Power LED 175 SCB**
- 20161420 **Cold Light Fountain Power LED 175,**  
power supply 110 – 240 VAC, 50/60 Hz
- 400 A **Mains Cord**
- 20090170 **SCB Connecting Cable,** length 100 cm



TP 95

# Cold Light Fountains XENON

Components / Spare Parts

Catalog page

- 69132601-1 Cold Light Fountain XENON 100 SCB**  
 20132620-1 **Cold Light Fountain XENON 100 SCB**,  
 incl. integrated insufflation pump, power supply  
 100 – 240 VAC, 50/60 Hz, for use with  
 KARL STORZ video endoscopes  
 400 A **Mains Cord**



TP 100

Components / Spare Parts

Catalog page

- 69131501 Cold Light Fountain XENON NOVA® 175**  
 20131520 **Cold Light Fountain XENON NOVA® 175**,  
 power supply 100 – 125 VAC/  
 220 – 240 VAC, 50/60 Hz  
 400 A **Mains Cord**



TP 101

Components / Spare Parts

Catalog page

- 20132101-1 Cold Light Fountain XENON 175 SCB**  
 20132120-1 **Cold Light Fountain XENON 175 SCB**,  
 power supply 100 – 125/220 – 240 VAC, 50/60 Hz  
 400 A **Mains Cord**  
 610 AFT **Silicone Tubing Set**, length 250 cm  
 20090170 **SCB Connecting Cable**, length 100 cm



TP 102

Components / Spare Parts

Catalog page

- 20134001 Cold Light Fountain XENON NOVA® 300**  
 20134020 **Cold Light Fountain XENON NOVA® 300**,  
 power supply 100 – 125/220 – 240 VAC, 50/60 Hz  
 400 A **Mains Cord**



TP 103

Components / Spare Parts

Catalog page


- 20133101-1 Cold Light Fountain XENON 300 SCB**  
 20133120-1 **Cold Light Fountain XENON 300 SCB**, power supply  
 100 – 125/220 – 240 VAC, 50/60 Hz  
 400 A **Mains Cord**  
 610 AFT **Silicone Tubing Set**, length 250 cm  
 20090170 **SCB Connecting Cable**, length 100 cm



TP 104







**EXCERPT FROM  
HYGIENE, MAINTENANCE  
STERILIZATION, STORAGE TECHNIQUES**



## KARL STORZ Service Program

### Full Service at First Hand

KARL STORZ has developed and optimized a special service program for the protection of patients, users and investments: ENDOPROTECT1

KARL STORZ ENDOPROTECT1 provides a comprehensive range of safe and cost-effective services to meet the needs of endoscopic equipment in everyday clinical practice.



### Repair Service

- Manufacturer quality with new product guarantee
- Technology updates
- Product group and third-party product replacement
- Transparent and calculable costs

### OR1™ Service

- Service and maintenance
- Hotline support
- On-site service
- Comprehensive support services

### Set Audits

- Standardized sets
- Stock optimization
- Reduced processing time
- Reduced repair costs


### Hygiene


- Processing and maintenance
- Innovative storage systems
- Wide range of cleaning accessories
- Hygiene consulting service

### Training

- Avoidance of errors
- Success in learning ensured
- Learning by doing
- Self-initiative

ENDOPROTECT1 Service Program:  
modular – safe – cost-effective





- **CLEANING AND DISINFECTION**
- **MAINTENANCE**
- **STERILISATION**
- **STORAGE TECHNIQUES**

### **Quality Standards during the Instrument Cycle**

Hygiene must be treated as an integral responsibility and a fundamental component of the entire health care process.

This principle begins with instrument design: in addition to technical solutions, effective cleaning, disinfection and sterilization represent the central factors in the development of all instruments. The objective is to create a reusable instrument, the materials, design and quality of which allow preparation using a suitable process that takes all applicable norms and laws into consideration, while maintaining the full functionality of the instrument.

Optimal handling throughout the instrument cycle goes further still and complements the instrument with individually coordinated accessories for optimized storage during the entire process – from the operating theater to transportation to the central sterile supply department, right down to the individual preparation steps.

Intelligent storage systems cut handling expenses. A reduced amount of repacking guarantees maximum protection for the instruments, optimizes operating theater procedures and ensures the desired results of the preparation.

KARL STORZ achieves this objective by continuously developing new storage systems and improving existing ones according to the demands of medical professionals.

# Table of Contents

Excerpt from HYGIENE, MAINTENANCE,  
STERILIZATION, STORAGE TECHNIQUES

## TABLE OF CONTENTS

Table of Contents	I
Index	II-III
Numerical Index	IV-V

## CLEANING, STERILIZATION AND MAINTENANCE, STORAGE TECHNIQUES

Accessories for Maintenance and Cleaning	CL 2-CL 6
Cleaning Accessories – for Fiber- and Video Endoscopes with Sheath Diameter $\leq$ 10 mm	CL 7
Cleaning Accessories – for Video Endoscopes with Sheath Diameter $\geq$ 10 mm	CL 8
Cases for Cleaning and Disinfection – Plastic Basins for Cleaning and Disinfection Solutions	CL 9
Quiver for Disinfection and Storage of Flexible Endoscopes	CL 10
Wire Trays for Cleaning, Sterilization and Storage of Endoscopes	CL 11-CL 13
Wire Trays for Cleaning, Sterilization and Storage of Semiflexible Endoscopes	CL 14
Wire Trays for Cleaning, Sterilization and Storage of Endoscopes and Instruments	CL 15-CL 17
Wire Trays for Cleaning, Sterilization and Storage of Instruments	CL 18, CL 21-CL 22
Accessories for Cleaning, Sterilization and Storage of Instruments	CL 19
Instrument Racks for Cleaning, Sterilization and Storage of Instruments	CL 20, CL 23-CL 24
Plastic Containers for Cleaning, Sterilization and Storage of Endoscopes	CL 25
Plastic Containers for Sterilization and Storage of Endoscopes	CL 26-CL 30
Plastic Containers for Sterilization and Storage of Camera Heads	CL 31
Plastic Containers for Sterilization and Storage of Instruments	CL 32-CL 34
Plastic Containers for Sterilization and Storage	CL 35-CL 36
Stainless Steel Tray for Sterilization and Storage of Instruments	CL 37
Aluminium Tray for Sterilization and Storage of Instruments	CL 38-CL 39
Sterilization Container with MicroStop® System for Sterilization and Sterile Storage	CL 40-CL 41
Accessories for Sterilization Container with MicroStop® System for Sterilization and Sterile Storage	CL 42-CL 44
Universal Aluminium Cases for HOPKINS® Telescopes	CL 45
Carrying Cases for Instruments and Accessories	CL 46

## COMPONENTS, SPARE PARTS

Grasping Forceps, Cleaning Pistol, Cleaning Tubing Set Gastro	CL 48
Wire Trays for Cleaning, Sterilization and Storage	CL 49
Instrument Racks and Plastic Containers for Cleaning, Sterilization and Storage	CL 50
Plastic Containers for Sterilization, Trays and Plastic Containers	CL 51
Plastic Containers for Sterilization and Storage	CL 52-CL 53
Sterilization Containers with MicroStop®	CL 54
MiniSet Containers with MicroStop®	CL 55

# Index

Excerpt from HYGIENE, MAINTENANCE,  
STERILIZATION, STORAGE TECHNIQUES

<b>A</b>			
Adaptor		CL 7	
Adaptor, air/water		CL 48	
Adaptor for aspiration		CL 48	
Adaptor for BHT Cleaning Machines		CL 7	
Adaptor for Leakage Tester		CL 8	
Adaptor for MIELE Cleaning Machines		CL 7	
Aluminium Case		CL 45	
Aluminium Tray for Sterilization and Storage		CL 38	
<b>B</b>			
Base		CL 40, CL 41, CL 55	
Brush		CL 5	
<b>C</b>			
Case Insert		CL 46	
Centering		CL 54	
Cleaning Adaptor		CL 3, CL 4	
Cleaning Agent		CL 2	
Cleaning Attachment		CL 6	
Cleaning Brush		CL 5	
Cleaning Pistol		CL 6, CL 48	
Cleaning Pistol Set		CL 6	
Cleaning Set		CL 5	
Cleaning Tubing Set Gastro		CL 8, CL 48	
Coding Framelet		CL 44	
Coding Label		CL 42, CL 44	
Color Tab		CL 42	
Cotton Applicator		CL 3	
<b>D</b>			
Disinfection System		CL 9	
Double Valve Seal		CL 48	
<b>F</b>			
Fixation Pin		CL 19, CL 21, CL 22	
<b>G</b>			
Grasping Forceps		CL 3, CL 48	
<b>H</b>			
Holder		CL 10	
<b>I</b>			
Instrument Oil		CL 2	
Instrument Rack for Cleaning, Sterilization and Storage		CL 20, CL 23, CL 50	
Instrument Rack for Storage		CL 34	
Instrument Rack Set for Cleaning, Sterilization and Storage		CL 20	
Irrigation Adaptor			CL 4
Irrigation Tube			CL 48, CL 50
<b>L</b>			
Large Diamond Grid Silicone Insert		CL 19, CL 21, CL 22	
Latch		CL 50, CL 51	
Leakage Tester		CL 7, CL 8	
Leakage Tester Adaptor		CL 48	
Lid		CL 21, CL 22, CL 40, CL 41, CL 52, CL 55	
Locking Piece, with clamp		CL 49	
Log Labels		CL 42	
Long Wire Mesh Tray		CL 18	
LUER-Lock Tube Connector		CL 4	
<b>M</b>			
Metal Tube			CL 10
MicroStop® Microbiological Barrier		CL 40, CL 41, CL 54, CL 55	
MiniSet Container		CL 43, CL 55	
<b>O</b>			
Oil Dropper			CL 2
Opening Spring			CL 3
O-Rings, for nozzle			CL 48
Oval Head Screw			CL 51
<b>P</b>			
Plastic Case			CL 46
Plastic Container for Cleaning, Sterilization and Storage			CL 25, CL 50
Plastic Container for Flexible Endoscopes			CL 30
Plastic Container for Sterilization			CL 26, CL 29, CL 30, CL 51
Plastic Container for Sterilization and Storage		CL 28, CL 31, CL 32, CL 33, CL 34, CL 35, CL 52	
Plastic Container for Sterilization and Storage of Endoscopes			CL 26, CL 27, CL 51
Plastic Container for Sterilization and Storage of Variable Instrument Sets			CL 36, CL 53
<b>Q</b>			
Quiver Head			CL 10
Quiver Set			CL 10
<b>R</b>			
Rod			CL 51
Rubber Jaws Insert			CL 48
<b>S</b>			
Seal			CL 42, CL 44
Sealing Cap			CL 10, CL 50
Silicone Grid Insert			CL 19, CL 43

# Index

Excerpt from HYGIENE, MAINTENANCE,  
STERILIZATION, STORAGE TECHNIQUES

Silicone Seal	CL 54, CL 55
Silicone Sealing Ring	CL 54
Silicone Telescope Holder	CL 49, CL 51
Silicone Tie-Downs	CL 19, CL 21, CL 22, CL 53
Snap-in Clip	CL 53
Spare Part Set, for measuring out	CL 48
Special Lubricant	CL 2
Sponge Holder	CL 3
Stainless Steel Tray for Sterilization, Storage and Transport	CL 37
STEINER Aluminium Tray	CL 39
Sterilization Container	CL 40, CL 41, CL 54
Strips of Silicone Knobs	CL 51

## T

Telescope Holder	CL 50
Telescope Lock Plug	CL 4
Tool	CL 53

## U

Universal Instrument Rack with Silicone Holders and Tie-downs	CL 20
--	-------

## W

Wadding Silver Polish "duraglit"	CL 2
Wire Tray for Cleaning, Sterilization and Storage	CL 11, CL 12, CL 13, CL 14, CL 15, CL 16, CL 17, CL 18, CL 21, CL 22, CL 24, CL 49
Wire Tray for MiniSet Container for Cleaning, Sterilization and Storage	CL 43

# Numerical Index

Excerpt from HYGIENE, MAINTENANCE,  
STERILIZATION, STORAGE TECHNIQUES

10384 A	CL 3	27658	CL 2	39301 LA	CL 28
10384 B	CL 3	27660	CL 6	39301 PHTS	CL 31
10387 A	CL 3	27660 A	CL 6	39301 Z3TS	CL 31
10387 B	CL 3	27660 B	CL 6	39305 C1S	CL 25, CL 50
11301 BM	CL 7	27660 C	CL 6	39305 C1SH	CL 50
11301 FF2	CL 7	27660 D	CL 6	39305 C2S	CL 25, CL 50
11301 GG2	CL 7	27660 E	CL 6	39305 C2SH	CL 50
11301 HH	CL 7	27660 F	CL 6	39305 L1S	CL 25, CL 50
11301 KK2	CL 7	27660 G	CL 6	39305 L1SH	CL 50
13242 LX	CL 8	27660 H	CL 6	39305 L2S	CL 25, CL 50
13242 XA	CL 7	27660 P	CL 6, CL 48	39305 L2SH	CL 50
13242 XL	CL 7	27661	CL 2	39305 LS	CL 50, CL 51
13242 XT	CL 8	27662	CL 3, CL 48	39311 P	CL 28
13991 CD	CL 48	27662 A	CL 48	39312 B	CL 34
13991 CK	CL 8, CL 48	27677 A	CL 46	39312 X	CL 34
13991 CL	CL 48	27677 B	CL 46	39312 XR	CL 34
13991 CR	CL 48	27677 C	CL 46	39314 G	CL 27
13991 CS	CL 48	27695 AE	CL 51	39351 J	CL 33
13991 SS	CL 48	27717 A	CL 33	39360 AK	CL 36, CL 53
26040 BX	CL 4	27717 D	CL 35	39360 AP	CL 53
27640 A	CL 45	28166 LL	CL 4	39360 AR	CL 53
27640 B	CL 45	28205 SA	CL 4	39360 AS	CL 19, CL 21, CL 22, CL 53
27640 C	CL 45	28205 SAA	CL 4	39360 BK	CL 36, CL 53
27645	CL 9	39100 L	CL 3	39401 AS	CL 29
27646	CL 9	39100 PS	CL 19, CL 21, CL 22	39402 AS	CL 29
27647 V	CL 9	39100 S	CL 19, CL 21	39403 AS	CL 30
27648 A	CL 5	39100 SH	CL 19, CL 22	39405 AS	CL 30
27650 A	CL 5	39100 SHS	CL 19	39501 A1	CL 11, CL 49
27650 B	CL 5	39100 SL	CL 19	39501 A2	CL 12, CL 49
27650 C	CL 5	39100 SMS	CL 43	39501 B1	CL 11, CL 49
27650 D	CL 5	39105 A	CL 3	39501 B2	CL 12, CL 49
27650 E	CL 5	39105 B	CL 3	39501 BAE	CL 13
27650 F	CL 5	39219 XF	CL 23, CL 50	39501 BEC	CL 13, CL 49
27650 G	CL 5	39219 XFB	CL 23, CL 50	39501 BR	CL 15, CL 49
27651 AK	CL 5	39219 XFS	CL 50	39501 C	CL 11, CL 49
27651 AL	CL 5	39219 XFV	CL 50	39501 CEC	CL 13, CL 49
27651 B	CL 5	39219 XU	CL 20	39501 ES	CL 16, CL 49
27651 C	CL 5	39219 XV	CL 20	39501 U	CL 17, CL 49
27651 E	CL 5	39219 XX	CL 20	39501 V	CL 17, CL 49
27651 F	CL 5	39301 ACTS	CL 31	39501 X	CL 14, CL 49
27651 G	CL 5	39301 AHS	CL 51	39501 XK	CL 14, CL 49
27651 H	CL 5	39301 AS	CL 26, CL 51	39502 L	CL 21
27651 I	CL 5	39301 BS	CL 26, CL 51	39502 LH	CL 22
27651 K1	CL 5	39301 CHS	CL 51	39502 MA	CL 43
27651 K2	CL 5	39301 CS	CL 26, CL 51	39502 MB	CL 43
27651 K3	CL 5	39301 D	CL 26, CL 51	39502 V	CL 24
27651 K4	CL 5	39301 F	CL 32, CL 52	39502 X	CL 18, CL 21
27651 K5	CL 5	39301 FB	CL 52	39502 XH	CL 18, CL 22
27651 UA	CL 5	39301 G	CL 32, CL 52	39502 Z	CL 18, CL 21
27651 UB	CL 5	39301 GB	CL 52	39502 ZH	CL 18, CL 22
27652	CL 5	39301 H	CL 32, CL 52	39502 ZL	CL 18
27655	CL 2	39301 HB	CL 52	39731 A	CL 38
27656 B	CL 2	39301 HCTS	CL 31	39731 B	CL 39
27656 S	CL 2	39301 J	CL 32, CL 52	39740 A2	CL 40, CL 54
27657	CL 2	39301 JB	CL 52		



# Numerical Index


Excerpt from HYGIENE, MAINTENANCE,  
STERILIZATION, STORAGE TECHNIQUES

39740 AB2	CL 40	39762 C	CL 44	39781 A	CL 43, CL 55
39741 A2	CL 40, CL 54	39762 C2	CL 42	39781 AB	CL 55
39741 AB2	CL 40	39762 D	CL 44	39785 AA	CL 55
39742 A2	CL 40, CL 54	39762 D2	CL 42	39785 AD	CL 55
39742 AB2	CL 40	39762 E	CL 44	39910 SC	CL 37
39743 A2	CL 40, CL 54	39762 E2	CL 42	5911202	CL 51
39743 AB2	CL 40	39762 F	CL 44	600007	CL 4
39745 AA2	CL 40	39762 F2	CL 42	60652 LL	CL 5
39745 AD	CL 54	39762 G	CL 44	620 N	CL 10
39750 A2	CL 41, CL 54	39762 G2	CL 42	621 C	CL 10
39750 AB2	CL 41	39762 H	CL 44	621 CT	CL 10
39751 A2	CL 41, CL 54	39762 H2	CL 42	621 D	CL 10
39751 AB2	CL 41	39762 I	CL 44	621 DT	CL 10
39752 A2	CL 41, CL 54	39762 I2	CL 42	621 H	CL 10
39752 AB2	CL 41	39763 A	CL 44	621 HT	CL 10
39753 A2	CL 41, CL 54	39763 A2	CL 42	621 L	CL 10
39753 AB2	CL 41	39763 B	CL 44	621 M	CL 10
39755 AA2	CL 41	39763 B2	CL 42	7920690	CL 51
39755 AD	CL 54	39763 E	CL 42	7981790	CL 49
39755 X	CL 54	39763 PA	CL 42	8021390	CL 48
39760 A	CL 55	39763 PB	CL 42	8021490	CL 48
39760 A2	CL 40, CL 41, CL 54, CL 55	39763 PKA	CL 44	8021590	CL 49
39762 A	CL 44	39763 PKB	CL 44	8021690	CL 49
39762 A2	CL 42	39764 A	CL 44	9200030	CL 46
39762 B	CL 44	39764 B	CL 44	9210560	CL 46
39762 B2	CL 42	39780 A	CL 43, CL 55	9211580	CL 46
		39780 AB	CL 55		





**CLEANING  
STERILIZATION AND MAINTENANCE  
STORAGE TECHNIQUES**



27655

**Oil Dropper**, for Instrument Oil



27656 S

**Oil Dropper**, for maintenance of instruments, silicone-free



27656 B

**Instrument Oil**, 50 ml, silicone-free



27657

**Special Lubricant**, for stopcocks



27658

**Wadding Silver Polish "duraglit"**



27661

**Cleaning Agent**, for distal or proximal lenses and fiber optic surface of endoscopes, content 5 g



10384 A **Cotton Applicator**, working length 50 cm  
 10384 B **Same**, working length 35 cm



10387 A **Sponge Holder**, with spring handle, working length 50 cm  
 10387 B **Same**, working length 35 cm



27662 **Grasping Forceps**, with serrated rubber jaws, for carefully handling the instruments








39100 L **Opening Spring**, for SilCut® and similar forceps



39105 A **Cleaning Adaptor**, for manual or machine cleaning of instruments without irrigation port from 1.5 to 5 mm diameter

39105 B **Cleaning Adaptor**, for manual or machine cleaning of instruments without irrigation port from 5 to 8 mm diameter


	<p>28166 LL     <b>Cleaning Adaptor</b>, for <b>ARTHROFORCE® TAKE-APART®</b> instruments</p>
	<p>600007     <b>LUER-Lock Tube Connector</b>, male, tube diameter 6 mm</p>
	<p>26040 BX     <b>Telescope Lock Plug</b>, with LUER-Lock for aspiration of liquid</p>
	<p>28205 SA     <b>Irrigation Adaptor</b>, for validated cleaning of shaver handpieces</p>
	<p>28205 SAA     <b>Irrigation Adaptor</b>, for validated cleaning of outer blades</p>



27652




27650 A-G



27651 AK-UB

Art. No.	Flat	Round/rigid	Round/flexible	Length	Outer diameter	Channel diameter	Notes
27652	●	-	-	-	-	-	for cleaning the jaws
27648 A	-	●	-	55 cm	16 mm	4 – 14 mm	
27650 A	-	●	-	35 cm	11 mm	3 – 9 mm	
27650 B	-	●	-	35 cm	7 mm	2.5 – 5 mm	
27650 C	-	●	-	35 cm	2.5 mm	2 – 2.4 mm	
27650 D	-	●	-	50 cm	11 mm	3 – 9 mm	
27650 E	-	●	-	50 cm	7 mm	2.5 – 5 mm	
27650 F	-	●	-	50 cm	2.5 mm	2 – 2.4 mm	
27650 G	-	●	-	50 cm	2.5 mm	-	for <b>CLICK'line</b> and RoBi®
27651 AK	-	-	●	75 cm	2 mm	1.2 – 1.8 mm	
27651 AL	-	-	●	150 cm	2 mm	1.2 – 1.8 mm	
27651 B	-	-	●	100 cm	3 mm	1.8 – 2.6 mm	
27651 C	-	-	●	230 cm	4 mm	2.4 – 3.6 mm	
27651 E	-	-	●	350 cm	4 mm	2.4 – 3.6 mm	
27651 F	-	-	●	180 cm	2.5 mm	1.4 – 2.3 mm	
27651 G	-	-	●	180 cm	3 mm	1.8 – 2.6 mm	
27651 H	-	-	●	230 cm	5 mm	2.2 – 4.6 mm	
27651 I	-	-	●	270 cm	5 mm	2.2 – 4.6 mm	
27651 UA	-	-	●	60 cm	2 mm	1.4 – 1.8 mm	
27651 UB	-	-	●	60 cm	2.2 mm	0.8 – 2 mm	
60652 LL	-	-	●	400 cm	5 mm	2.2 – 4.6 mm	veterinary medicine
27651 K1	-	-	●	40 cm	-	0.4 – 0.6 mm	for single use, 10 pcs.
27651 K2	-	-	●	40 cm	-	0.6 – 0.8 mm	for single use, 10 pcs.
27651 K3	-	-	●	40 cm	-	0.8 – 1.4 mm	for single use, 10 pcs.
27651 K4	-	-	●	40 cm	-	0.2 – 0.4 mm	for single use, 10 pcs.
27651 K5	-	-	●	40 cm	-	0.25 – 0.4 mm	for single use, 10 pcs.

7-974



27660

27660

**Cleaning Pistols Set, autoclavable**, for fast and easy rinsing or drying of instruments, attachable to standard hospital 3/4" water or compressed air connections, with fabric tube, including Cleaning Attachments 27660 A – H and wall support



27660 P

**Cleaning Pistol, only, autoclavable**, with fibrous tubing, without accessories

### Cleaning Attachments, for use with Cleaning Pistol 27660 P



27660 A

**Attachment**, for cannulas



27660 B

**Attachment**, for channels



27660 C

**Attachment**, for valves



27660 D

**Attachment**, for tubes



27660 E

**Attachment**, for Luer connectors



27660 F

**Attachment**, for surfaces and hollow spaces



27660 G

**Attachment**, for surfaces and hollow spaces



27660 H

**Attachment**, as water jet pump for suction of corners and hollow spaces



# Cleaning Accessories

for Fiber- and Video Endoscopes with Sheath Diameter ≤ 10 mm



11301 BM **Adaptor**, for leakage test, for Belimed washer-disinfectors



11301 HH **Adaptor for BHT Cleaning Machines**, for automatic leakage test of flexible KARL STORZ endoscopes



13242 XA **Adaptor**, for manual leakage test with 13242 XL or automatic leakage test in washer/disinfector



11301 FF2 **Adaptor for MIELE Cleaning Machines**, with safety valve, for automatic leakage test of flexible KARL STORZ endoscopes



11301 GG2 **Adaptor**, for cleaning and disinfecting the irrigation and working channels of flexible endoscopes, for MIELE-ETD washer-disinfectors



11301 KK2 **Adaptor**, for working channel of flexible endoscopes, for MIELE-ETD 03 washer/disinfectors  
Please note: Adaptors 11301 FF2 and 11301 GG2 have to be ordered separately!



13242 XL **Leakage Tester**, with bulb and manometer

# Cleaning Accessories

for Video Endoscopes with Sheath Diameter  $\geq 10$  mm



13991 CK

**Cleaning Tubing Set Gastro**, for the cleaning and disinfection of video endoscopes in gastroenterology, particularly for use with RDG-E automatic units with one-channel monitoring  
including:  
**Seal for double valve**  
**Leakage tester adaptor**  
**Adaptor for aspiration**  
**Adaptor air/water**  
**Irrigation Tube**, reusable



13242 XT

**Adaptor for Leakage Tester**, large, with 60 cm tube for GI video endoscopes



13242 LX

**Leakage Tester**, with bulb and manometer

**Components** see COMPONENTS, SPARE PARTS

# Cases for Cleaning and Disinfection

Plastic Basins for Cleaning and Disinfection Solutions



- 27645      **Disinfection System**, for use with fiberscopes including:  
1. **Soak Basin**  
2. **Perforated Storage Tray**  
3. **Lid**  
Internal dimensions (w x d x h): 503 x 186 x 132 mm
- 27646      **Disinfection System**, for use with fiberscopes including:  
1. **Soak Basin**  
2. **Perforated Storage Tray**  
3. **Lid**  
Internal dimensions (w x d x h): 740 x 220 x 132 mm
- 27647 V      **Disinfection System**, for use with fiberscopes including:  
1. **Soak Basin**  
2. **Perforated Storage Tray**  
3. **Lid**  
Internal dimensions (w x d x h): 694 x 381 x 141 mm

**Please note:**

- Instruments must be disassembled and cleaned prior to disinfection. To disinfect, immerse instrument completely in solution.
- Follow manufacturer's instructions regarding time period for soaking and solution density.
- To wash out instruments, remove perforated storage tray from the disinfectant solution and rinse instruments thoroughly with sterile or desalinated water in a second basin.
- Tilt basin gently to remove air bubbles during disinfection.

# Quiver for Disinfection\* and Storage of Flexible Endoscopes

For fixation to storage cabinets



621 C

**Quiver Set**, for disinfection\* and storage of fiberscopes including:

**Quiver Head**, for fiberscopes

**Metal Tube**

**Sealing Cap**

**Holder**, to be screwed to a storage cabinet

Quiver Set	Working Length max.	Quiver Head	Metal Tube	Sealing Cap	Holder
621 C	30 cm	621 L	621 CT	621 M	620 N
621 D	40 cm	621 L	621 DT	621 M	620 N
621 H	76 cm	621 L	621 HT	621 M	620 N

**Please note:** Precisely follow manufacturer's specifications regarding solution density and soaking period when mixing and applying solutions. If soaked too long, instruments may be damaged. Store endoscopes in a dry place after rinsing with demineralized or sterilized water.


**Note:** Steam sterilization of quiver set is possible at a temperature of 134 °C.

\* Please also take note of the statutory regulations valid in your country regarding the disinfection of quivers.


# Wire Trays for Cleaning, Sterilization and Storage of Endoscopes

For safe storage of telescopes during transport and machine cleaning. For sterile storage the trays may either be wrapped individually or put into a container.

The trays may be used for sterilization with steam, gas or plasma.



39501 A1




39501 B1

- 39501 A1 **Wire Tray for Cleaning, Sterilization and Storage** of one rigid endoscope, including holder for light post adaptors, silicone telescope holders and lid, external dimensions (w x d x h): 290 x 60 x 52 mm, for rigid endoscopes up to diameter 5 mm and working length 20 cm
- 39501 B1 **Wire Tray for Cleaning, Sterilization and Storage** of one rigid endoscope, including holder for light post adaptors, silicone telescope holders and lid, external dimensions (w x d x h): 430 x 65 x 52 mm, for rigid endoscopes up to diameter 10 mm and working length 34 cm
- 39501 C **Wire Tray for Cleaning, Sterilization and Storage** of one rigid telescope, with silicone telescope holder and lid, external dimensions (w x d x h): 670 x 80 x 52 mm, for telescopes for bronchoscopy and esophagoscopy


**Please note:** The instruments displayed are not included in the sterilization tray.

# Wire Trays for Cleaning, Sterilization and Storage of Endoscopes



39501 A2

39501 A2 **Wire Tray for Cleaning, Sterilization and Storage** of two rigid endoscopes and one light guide cable, including holder for light post adaptors, silicone telescope holders and lid, external dimensions (w x d x h): 352 x 125 x 54 mm, for rigid endoscopes up to diameter 10 mm and working length 20 cm




39501 B2

39501 B2 **Wire Tray for Cleaning, Sterilization and Storage** of two rigid endoscopes and one light cable, including holder for light post adaptors, silicone telescope holders and lid, external dimensions (w x d x h): 487 x 125 x 54 mm, for rigid endoscopes up to diameter 10 mm diameter and working length 34 cm


**Please note:** The instruments displayed are not included in the sterilization tray.

# Wire Trays for Cleaning, Sterilization and Storage of Endoscopes <sup>NEW</sup>



39501 BAE

39501 BAE **Wire Tray for Cleaning, Sterilization and Storage** of one rigid endoscope with angled eyepiece, including holder for light post adaptors, silicone telescope holders and lid, external dimensions (w x d x h): 460 x 80 x 52 mm




39501 BEC

39501 BEC **Wire Tray for Cleaning, Sterilization and Storage** of one ENDOCAMELEON®, length 32 cm and one light cable, including holder for light post adaptor, silicone telescope holder and lid, external dimensions (w x d x h): 480 x 125 x 54 mm

39501 CEC **Wire Tray for Cleaning, Sterilization and Storage** of one ENDOCAMELEON®, length 42 cm and one light cable, including holder for light post adaptor, silicone telescope holder and lid, external dimensions (w x d x h): 604 x 125 x 54 mm

**Please note:** The instruments displayed are not included in the sterilization tray.

# Wire Trays for Cleaning, Sterilization and Storage of Semiflexible Endoscopes




39501 XK

- 39501 X **Wire Tray for Cleaning, Sterilization and Storage of uretero-renoscopes**, including cleaning adaptor for washer-disinfector, with lid, spare parts basket 39501 XS and silicone telescope holders, external dimensions (w x d x h): 644 x 150 x 80 mm, for use with URS models with working length up to 43 of the product numbers 270xx
- 39501 XK **Wire Tray for Cleaning, Sterilization and Storage of pediatric uretero-renoscopes and compact miniature hysteroscopes** with a maximum working length of 27 cm, including cleaning adaptor for washer-disinfector, with lid, spare parts basket 39501 XS and silicone telescope holders, external dimensions (w x d x h): 460 x 150 x 80 mm, for use with instruments 270xx and 260x2 K with up to 27 cm working length

**Please note:** The instruments displayed are not included in the sterilization tray.



# Wire Trays for Cleaning, Sterilization and Storage of Endoscopes and Instruments



39501 BR


39501 BR

**Wire Tray for Cleaning, Sterilization and Storage of telescopes and accessories for bronchoscopy, including silicone telescope holders and lid, external dimensions (w x d x h): 622 x 150 x 78 mm**

3-071

**Please note:** The instruments displayed are not included in the sterilization tray.

# Wire Trays for Cleaning, Sterilization and Storage of Endoscopes and Instruments




39501 ES

39501 ES

**Wire Tray for Cleaning, Sterilization and Storage of instruments and accessories for oesophagoscopy, including silicone telescope holders and lid, external dimensions (w x d x h): 640 x 216 x 78 mm**

**Please note:** The instruments displayed are not included in the sterilization tray.

# Wire Trays for Cleaning, Sterilization and Storage of Endoscopes and Instruments



39501 U

39501 U **Wire Tray for Cleaning, Sterilization and Storage of telescopes and instruments**, with silicone telescope holders, tray for small parts and lid, external dimensions (w x d x h): 510 x 185 x 82 mm, for angled telescopes and uretero-renoscopes up to working length 34 cm

39501 V **Wire Tray for Cleaning, Sterilization and Storage of telescopes and instruments**, with silicone telescope holders, tray for small parts and lid, external dimensions (w x d x h): 608 x 185 x 82 mm, for angled telescopes and uretero-renoscopes up to working length 43 cm


# Wire Trays for Cleaning, Sterilization and Storage of Instruments <sup>NEW</sup>

## Components of modular processing systems

- Trays made from electropolished stainless steel mesh wire 39502 XH/X or with hole plate walls 39502 ZH/Z/ZL
- Wide-mesh silicone knob mats “Diamond” and/or “Large Diamond” 39502 SHK/SH/S/SL
- Fixiation Pin 39100 PS for screwing instruments into position
- Silicone Tie-down 39360 AS
- The modular processing system can be ideally combined with Instrument Racks 39219 XV/XU



39502 ZH




39502 Z

**Standard wire trays for cleaning, sterilization and storage** of instruments, mesh wire or hole plate walls, stackable and foldaway handles

Wire Tray	External Dimensions (w x d x h)	Special Features	Suitable for
39502 XH	240 x 250 x 50 mm	mesh wire	div. instruments
39502 ZH	240 x 250 x 66 mm	hole plate wall	div. instruments
39502 X	480 x 250 x 50 mm	mesh wire	div. instruments
39502 Z	480 x 250 x 66 mm	hole plate wall	div. instruments
39502 ZL	535 x 250 x 66 mm	hole plate wall	div. instruments

**Please note:** The instruments displayed are not included in the sterilization tray.

# Accessories for Cleaning, Sterilization and Storage of Instruments <sup>NEW</sup>



39100 S

- 39100 S     **Large Diamond Grid Silicone Insert**, blue, extra wide mesh, for the storage of instruments in standard wire trays, plastic and sterilization containers, external dimensions: (w x d): 470 x 240 mm
- 39100 SH    **Same**, external dimensions (w x d): 230 x 240 mm
- 39100 SL    **Same**, external dimensions (w x d): 530 x 240 mm
- 39100 SHS   **Silicone Grid Insert, "Small Diamond Grid"**, blue, wide meshed, for storage of instruments in standard wire baskets and sterilization container. External dimensions: 240 x 230 mm



39100 PS


- 39100 PS     **Fixation Pin**, including screw and washer, to screw instruments into position in wire trays, height 38 mm, package of 12, for use with Silicone Tie-Downs 39360 AS



39360 AS

- 39360 AS     **Silicone Tie-Downs**, package of 12, for use with Fixation Pins 39100 PS and 39360 AP


# Instrument Racks for Cleaning, <sup>NEW</sup> Sterilization and Storage of Instruments



39219 XX

39219 XX **Instrument Rack Set for Cleaning, Sterilization and Storage** of up to 14 instruments with diameter 2.5 to 10 mm, incl. variable bars with silicone holders, rack with Tray 39502 V for drawer and Wire Tray 39502 X, external dimensions (w x d x h): 480 x 250 x 125 mm

39219 XV **Instrument Rack for Cleaning, Sterilization and Storage** of up to 14 instruments with diameter 2.5 to 10 mm, incl. variable bars with silicone holders, rack with Tray 39502 V for drawer and Wire Tray 39502 X, for use with standard wire trays, external dimensions (w x d x h): 463 x 238 x 125 mm




39219 XU

39219 XU **Universal Instrument Rack with Silicone Holders and Tie-downs**, with Tray 39502 V as drawer, for storage of instrument sets, for use with standard wire trays, external dimensions (w x d x h): 460 x 233 x 145 mm

**Please note:** The instruments displayed are not included in the instrument tray.

# Wire Trays for Cleaning, Sterilization and Storage of Instruments <sup>NEW</sup>

Modular processing system: Possible set combinations




39502 X

### Set 1

- 39502 X **Wire Tray, stackable, mesh wire**
- 39502 L **Lid**
- 39100 S **Silicone Insert**
- 39100 PS **Fixation Pins, package of 12**
- 39360 AS **Silicone Tie-Downs, package of 12**

Can be combined with 39219 XV or 39219 XU  
(cannot be used with Lid 39502 L)



39502 Z

### Set 2

- 39502 Z **Wire Tray, stackable, with hole plates**
- 39502 L **Lid**
- 39100 S **Silicone Insert**
- 39100 PS **Fixation Pins, package of 12**
- 39360 AS **Silicone Tie-Downs, package of 12**

Can be combined with 39219 XV or 39219 XU  
(cannot be used with Lid 39502 L)

**Please note:** The instruments displayed are not included in the sterilization tray.

# Wire Trays for Cleaning, Sterilization and Storage of Instruments <sup>NEW</sup>

Modular processing system: Possible set combinations



39502 XH

## Set 1B

- |          |  |
|----------|--|
| 39502 XH | <b>Wire Tray, stackable, mesh wire</b>   |
| 39502 LH | <b>Lid</b>                               |
| 39100 SH | <b>Silicone Insert</b>                   |
| 39100 PS | <b>Fixation Pins, package of 12</b>      |
| 39360 AS | <b>Silicone Tie-Downs, package of 12</b> |



39502 ZH

## Set 2B

- |          |   |
|----------|---|
| 39502 ZH | <b>Wire Tray, stackable, with hole plates</b> |
| 39502 LH | <b>Lid</b>                                    |
| 39100 SH | <b>Silicone Insert</b>                        |
| 39100 PS | <b>Fixation Pins, package of 12</b>           |
| 39360 AS | <b>Silicone Tie-Downs, package of 12</b>      |

**Please note:** The instruments displayed are not included in the sterilization tray.




# Instrument Racks for Cleaning, Sterilization and Storage of Instruments <sup>NEW</sup>

Instrument Racks, specially designed for SinuFIT® and LaryngoFIT® instruments

**Special Features:**

- Optimized instrument storage
- Separate storage for inner and outer sheaths of 3-part dismantling instruments
- Individual connection of sheaths via LUER adaptor
- Positive opening of jaws of 2-part dismantling instruments
- Systematic flushing of handles
- Direct connection of instrument rack to washer-disinfector



39219 XFB

**39219 XF** **Instrument Rack for Cleaning, Sterilization and Storage** of up to 15 **LaryngoFIT®** or 15 **SinuFIT®** instrument attachments and 6 handles, incl. variable bars with silicone holder and direct connection of the instrument attachments via a multi-port bridge to the washer-disinfector, external dimensions (w x d x h): 480 x 250 x 160 mm

**39219 XFB** **Instrument Rack for Cleaning, Sterilization and Storage** of up to 15 **LaryngoFIT®** or 30 **SinuFIT®** working inserts and 6 handles, incl. variable bars with silicone holder and direct connection of the working inserts via a multi-port bridge to the washer-disinfector, external dimensions (w x d x h): 480 x 250 x 160 mm

**Please note:** The instruments displayed are not included in the sterilization tray.

# Instrument Racks for Cleaning, <sup>NEW</sup> Sterilization and Storage of Instruments

## Optional accessories



39502 V


39502 V

**Wire Tray for Cleaning, Sterilization and Storage**  
of instruments, with hole plate walls and foldaway  
handles, external dimensions (w x d x h):  
228 x 248 x 45 mm

# Plastic Containers for Cleaning, <sup>NEW</sup> Sterilization and Storage of Endoscopes

## Special Features:

- Combines the advantages of the proven wire mesh tray for telescopes and that of a lightweight plastic tray
- Wide-mesh tray for minimal influence on the cleaning kinetics of the washer-disinfector
- Minimal instrument contact surface for maximum cleaning performance
- Vibration-absorbing silicone holders for maximum protection of optical systems during transport and machine cleaning
- Approx. 60 % weight reduction compared to the proven steel wire trays



39305 L2S


- 39305 C1S **Plastic Container for Cleaning, Sterilization and Storage** of one rigid endoscope with a diameter up 5 mm and working length 34 cm, including holders for light post adaptors, for variable positioning of silicone telescope holders and lid, external dimensions (w x d x h): 470 x 105 x 50 mm
- 39305 C2S **Plastic Container for Cleaning, Sterilization and Storage** of two rigid endoscopes with a diameter up 5 mm and working length 34 cm, including holders for light post adaptors, for variable positioning of silicone telescope holders and lid, external dimensions (w x d x h): 470 x 105 x 50 mm
- 39305 L1S **Plastic Container for Cleaning, Sterilization and Storage** of one rigid endoscope with a diameter up 10 mm and working length 34 cm, including holders for light post adaptors, for variable positioning of silicone telescope holders and lid, external dimensions (w x d x h): 470 x 105 x 50 mm
- 39305 L2S **Plastic Container for Cleaning, Sterilization and Storage** of two rigid endoscopes with a diameter up 10 mm and working length 34 cm, including holders for light post adaptors, for variable positioning of silicone telescope holders and lid, external dimensions (w x d x h): 470 x 105 x 50 mm

**Please note:** The instruments displayed are not included in the sterilization tray.

# Plastic Containers for Sterilization and Storage of Endoscopes

Plastic container with silicone telescope holders for secure fixation of rigid endoscopes during sterilization and storage.

The plastic containers are suitable for steam, gas and hydrogen peroxide (Sterrad®) sterilization.



39301 BS

- NEW** 39301 AS **Plastic Container for Sterilization**, suitable for steam, gas, and hydrogen peroxide (Sterrad®) sterilization and storage, perforated, with lid, for use with two rigid endoscopes up to max. 20 cm working length, external dimensions (w x d x h): 321 x 90 x 45 mm
- NEW** 39301 BS **Plastic Container for Sterilization**, suitable for steam, gas, and hydrogen peroxide (Sterrad®) sterilization and storage, perforated, with lid, for use with two rigid endoscopes up to max. 34 cm working length, external dimensions (w x d x h): 446 x 90 x 45 mm
- NEW** 39301 CS **Plastic Container for Sterilization**, suitable for steam, gas, and hydrogen peroxide (Sterrad®) sterilization and storage, perforated, with lid, for use with two rigid endoscopes with a working length of over 34 cm, external dimensions (w x d x h): 520 x 90 x 45 mm
- 39301 D **Plastic Container for Sterilization and Storage of Endoscopes**, perforated, with transparent lid, with silicone telescope holder, external dimensions: (w x d x h): 670 x 116 x 80 mm

**Please note:** The instruments displayed are not included in the plastic containers.

# Plastic Containers for Sterilization and Storage of Endoscopes

Plastic container with transparent lid and silicone telescope holders for secure fixation of rigid endoscopes during sterilization and storage. The

plastic containers are suitable for steam, gas and hydrogen peroxide sterilization.



39314 G

39314 G

**Plastic Container for Sterilization and Storage of Endoscopes**, with angled eyepiece, perforated, with transparent lid, with inserts for two angled telescopes of up to 10 mm, external dimensions (w x d x h): 520 x 250 x 84 mm

**Please note:** The instruments displayed are not included in the plastic containers.

# Plastic Containers for Sterilization and Storage of Endoscopes



39311 P

39301 LA

- 39311 P **Plastic Container for Sterilization and Storage**, with silicone holders and box for small parts, for use with 1 uretero-renaloscope of up to 34 cm working length, external dimensions (w x d x h): 510 x 115 x 90 mm
- 39301 LA **Plastic Container for Sterilization and Storage**, with silicone holders and box for small parts, for use with 2 uretero-renaloscopes with angled eyepiece of up to 45 cm working length, external dimensions (w x d x h): 812 x 113 x 82 mm

**Please note:** The instruments displayed are not included in the plastic containers.

# Plastic Containers for Sterilization and Storage of Endoscopes



39401 AS



39402 AS

39401 AS **Plastic Container for Sterilization**, specially suited for gas and hydrogen peroxide (Sterrad®) sterilization and storage, perforated, with lid, for use with flexible endoscopes up to max. 39 cm working length, external dimensions (w x d x h): 690 x 170 x 92 mm

39402 AS **Plastic Container for Sterilization**, specially suited for gas and hydrogen peroxide (Sterrad®) sterilization and storage, perforated, with lid, for use with flexible endoscopes up to max. 95 cm working length, external dimensions (w x d x h): 550 x 260 x 92 mm

**Please note:**

Fiberscopes stored in these containers can only be sterilized with gas or plasma. The instruments displayed are not included in the plastic containers.

# Plastic Containers for Sterilization and Storage of Endoscopes



39403 AS



39405 AS

39403 AS **Plastic Container for Sterilization**, specially suited for gas and hydrogen peroxide (Sterrad®) sterilization and storage, perforated, with lid, for use with flexible video-cysto-, rhino- and bronchoscopes, external dimensions (w x d x h): 550 x 260 x 90 mm

39405 AS **Plastic Container for Flexible Endoscopes**, specially suited for gas and hydrogen peroxide (Sterrad®) sterilization and storage, for use with one flexible endoscope, external dimensions (w x d x h): 550 x 260 x 90 mm

**Please note:** The instruments displayed are not included in the plastic containers.



## Plastic Containers for Sterilization and Storage of Camera Heads



39301 Z3TS

- 39301 ACTS **Plastic Container for Sterilization and Storage** of camera heads **IMAGE1, TRICAM® and TELECAM**, suitable for steam, gas, and hydrogen peroxide sterilization, Sterrad® compatible, external dimensions (w x d x h): 385 x 255 x 75 mm
- 39301 HCTS **Plastic Container for Sterilization and Storage** of **HD camera heads**, suitable for steam, gas, and hydrogen peroxide sterilization, Sterrad® compatible, external dimensions (w x d x h): 385 x 255 x 75 mm
- 39301 Z3TS **Plastic Container for Sterilization and Storage** of camera heads **IMAGE1 H3-Z, H3-ZA and H3-FA, autoclavable**, suitable for use with steam, gas and hydrogen peroxide sterilization, Sterrad® compatible, external dimensions (w x d x h): 385 x 255 x 75 mm
- 39301 PHTS **Plastic Container for Sterilization and Storage** of camera heads **IMAGE1 H3-P and H3-ZI, autoclavable**, suitable for use with steam, gas and hydrogen peroxide sterilization, Sterrad® compatible, external dimensions (w x d x h): 385 x 255 x 75 mm

**Please note:** The instruments displayed are not included in the plastic containers.

# Plastic Containers for Sterilization and Storage of Instruments

Plastic container with transparent lid, or secure fixation of instruments during sterilization and storage.

The plastic containers may be used for sterilization with steam, gas or hydrogen peroxide sterilization.



39301 H

- 39301 F **Plastic Container for Sterilization and Storage,** perforated, with transparent lid, with silicone mat, external dimensions (w x d x h): 260 x 233 x 65 mm
- 39301 G **Plastic Container for Sterilization and Storage,** perforated, with transparent lid, with silicone mat, external dimensions (w x d x h): 515 x 114 x 65 mm
- 39301 H **Plastic Container for Sterilization and Storage,** perforated, with transparent lid, with silicone mat, external dimensions (w x d x h): 515 x 233 x 65 mm
- 39301 J **Plastic Container for Sterilization and Storage,** perforated, with transparent lid, with silicone mat, external dimensions (w x d x h): 532 x 252 x 141 mm

# Plastic Containers for Sterilization and Storage of Instruments



39351 J




27717 A

39351 J **Plastic Container for Sterilization and Storage**, perforated, with transparent lid, for **two-level storage**, for use with forceps and instruments, external dimensions (w x d x h): 530 x 250 x 145 mm

27717 A **Plastic Container for Sterilization and Storage**, perforated, with transparent lid and silicone mat, for two-level storage, external dimensions (w x d x h): 543 x 255 x 112 mm

# Plastic Containers for Sterilization and Storage of Instruments



39312 X **Plastic Container for Sterilization and Storage**, with removable Instrument Rack 39312 XR and insert tray for storage of up to 12 instruments with diameter 2.5 to 10 mm and up to 6 trocars, perforated, with transparent lid, external dimensions (w x d x h): 532 x 254 x 165 mm


39312 XR **Instrument Rack for Storage** of up to 12 instruments with diameter 2.5 to 10 mm, including 2 bars with silicone holders



39312 B **Plastic Container for Sterilization and Storage**, perforated, with transparent lid and silicone mat, for use with up to 15 instruments with working length up to 43 cm, external dimensions (w x d x h): 680 x 325 x 155 mm

**Please note:** The displayed instruments are not included in the instrument rack.

# Plastic Containers for Sterilization and Storage




27717 D

27717 D **Plastic Container for Sterilization and Storage,** perforated, with transparent lid, for use with 12 **CLICKline** instruments with working length up to 36 cm, external dimensions (w x d x h): 585 x 255 x 145 mm

**Please note:** The instruments displayed are not included in the plastic container.

# Plastic Containers for Sterilization and Storage



39360 AK


39360 BK

- 39360 AK **Plastic Container for Sterilization and Storage of Variable Instrument Sets**, perforated, with transparent lid, with silicone mat, **two-level storage**, (1 additional insert), external dimensions (w x d x h): 545 x 260 x 115 mm including:  
2x **Snap-in Clip**, package of 12  
2x **Silicone Tie-Downs**, package of 12  
**Tool**
- 39360 BK **Plastic Container for Sterilization and Storage of Variable Instrument Sets**, perforated, with transparent lid, with silicone mat, **single-level storage**, (1 additional insert), external dimensions (w x d x h): 525 x 240 x 70 mm including:  
**Snap-in Clip**, package of 12  
**Silicone Tie-Downs**, package of 12  
**Tool**

**Components** see COMPONENTS, SPARE PARTS

**Please note:** The instruments displayed are not included in the plastic container.

## Stainless Steel Tray for Sterilization and Storage of Instruments




39910 SC

39910 SC

**Stainless Steel Tray for Sterilization, Storage and Transport** of up to 12 SilCut® or RHINOFORCE® forceps, with lid, latches and handles, removable rack with double silicone holders for max. 12 instruments, external dimensions (w x d x h): 340 x 250 x 145 mm

**Please note:** The instruments displayed are not included in the stainless tray.

# Aluminium Tray for Sterilization and Storage of Instruments



39731 A


39731 A

**Aluminium Tray for Sterilization and Storage of STAMMBERGER-FESS instrument sets, two-level storage, with instrument rack, external dimensions (w x d x h): 525 x 248 x 202 mm**

**Please note:** The instruments displayed are not included in the aluminium tray.



# Aluminium Tray for Sterilization and Storage of Instruments



39731 B

39731 B **STEINER Aluminium Tray**, for sterilization and storage of instrument sets for laser laryngoscopy, external dimensions (w x d x h): 525 x 248 x 202 mm

1-0032

**Please note:** The instruments displayed are not included in the aluminium tray.

# Sterilization Container <sup>NEW</sup> with MicroStop® System for Sterilization and Sterile Storage

**Special Features:**

- Cost-effective, time-saving and eco-friendly
- Ergonomic, standardized and stackable
- MicroStop® microbiological barrier with a separation efficiency of 99.99997%
- Secure safety lock with seal slot



39740 A2

**Sterilization Container**, with MicroStop® microbiological barrier,  
for sterilization and sterile storage

including:

**Base**

**Lid**

**Microbiological Barrier**

Container	Base	Lid	Micro- biological Barrier	External Dimensions (w x d x h)	Inner Dimensions (w x d x h)
<b>39740 A2</b>	39740 AB2	39745 AA2	39760 A2	300 x 300 x 110 mm	262 x 267 x 88 mm
<b>39741 A2</b>	39741 AB2	39745 AA2	39760 A2	300 x 300 x 140 mm	262 x 267 x 113 mm
<b>39742 A2</b>	39742 AB2	39745 AA2	39760 A2	300 x 300 x 160 mm	262 x 267 x 138 mm
<b>39743 A2</b>	39743 AB2	39745 AA2	39760 A2	300 x 300 x 210 mm	262 x 267 x 186 mm

# Sterilization Container <sup>NEW</sup> with MicroStop® System for Sterilization and Sterile Storage

**Special Features:**

- Cost-effective, time-saving and eco-friendly
- Ergonomic, standardized and stackable
- MicroStop® microbiological barrier with a separation efficiency of 99.99997%
- Secure safety lock with seal slot



39750 A2

**Sterilization Container**, with MicroStop® microbiological barrier,  
for sterilization and sterile storage

including:

**Base**

**Lid**

**Microbiological Barrier**

Container	Base	Lid	Micro-biological Barrier	External Dimensions (w x d x h)	Inner Dimensions (w x d x h)
<b>39750 A2</b>	39750 AB2	39755 AA2	39760 A2	600 x 300 x 110 mm	548 x 267 x 88 mm
<b>39751 A2</b>	39751 AB2	39755 AA2	39760 A2	600 x 300 x 140 mm	548 x 267 x 113 mm
<b>39752 A2</b>	39752 AB2	39755 AA2	39760 A2	600 x 300 x 160 mm	548 x 267 x 138 mm
<b>39753 A2</b>	39753 AB2	39755 AA2	39760 A2	600 x 300 x 210 mm	548 x 267 x 186 mm

# Accessories for Sterilization Container <sup>NEW</sup> with MicroStop® System for Sterilization and Sterile Storage

## Coding Labels



39763 B2

- 39763 A2     **Coding Label**, inscribed  
39763 B2     **Coding Label**, blank

## Color Tabs



39762 C2


Art. No.	Designation	Color
39762 A2	Color Tab	black
39762 B2	Color Tab	white
39762 C2	Color Tab	red
39762 D2	Color Tab	blue
39762 E2	Color Tab	green
39762 F2	Color Tab	yellow
39762 G2	Color Tab	violet
39762 H2	Color Tab	orange
39762 I2	Color Tab	grey

## Optional Accessories



39763 E

- 39763 E     **Log Labels**, self-adhesive,  
package of 1000



39763 PB

- 39763 PA     **Seal**, without indicator, package of 1000  
39763 PB     **Seal**, with indicator, package of 1000

# Accessories for Sterilization Container <sup>NEW</sup> with MicroStop® System for Sterilization and Sterile Storage



39780 A

39780 A **MiniSet Container**, with MicroStop®, for sterilization and sterile storage, external dimensions (w x d x h): 310 x 189 x 65 mm, internal dimensions (w x d x h): 283 x 177 x 40 mm including:

**Base**

**Lid**


**Microbiological Barrier**, with MicroStop®

39781 A **MiniSet Container**, with MicroStop®, for sterilization and sterile storage, external dimensions (w x d x h): 310 x 189 x 80 mm, internal dimensions (w x d x h): 283 x 177 x 60 mm including:

**Base**

**Lid**

**Microbiological Barrier**, with MicroStop®



39502 MA

39502 MA **Wire Tray for MiniSet Container for Cleaning, Sterilization and Storage** of instruments, external dimensions (w x d x h): 260 x 165 x 28 mm

39502 MB **Wire Tray for MiniSet Container for Cleaning, Sterilization and Storage** of instruments, external dimensions (w x d x h): 260 x 165 x 47 mm

39100 SMS **Silicone Grid Insert, "Small Diamond Grid"**, blue, wide mesh, for the storage of instruments in MiniSet wire trays and MiniSet plastic sterilization containers, external dimensions: 250 x 155 mm

# Accessories for Sterilization Container <sup>NEW</sup> with MicroStop® System for Sterilization and Sterile Storage


## Coding Labels, for lids



39763 B

- 39763 A **Coding Label**, inscribed
- 39763 B **Coding Label**, blank


## Coding Labels, for containers



39764 B

- 39764 A **Coding Label**, inscribed
- 39764 B **Coding Label**, blank


## Coding (Logistics) Framelets, for lids



39762 A

Art. No.	Designation	Color
39762 A	Coding Framelet	black
39762 B	Coding Framelet	white
39762 C	Coding Framelet	red
39762 D	Coding Framelet	blue
39762 E	Coding Framelet	green
39762 F	Coding Framelet	yellow
39762 G	Coding Framelet	violet
39762 H	Coding Framelet	orange
39762 I	Coding Framelet	grey

## Accessories




39763 PKB

- 39763 PKA **Seal**, without indicator, package of 1000
- 39763 PKB **Seal**, with indicator, package of 1000

# Universal Aluminium Cases for HOPKINS® Telescopes

Universal cases are suitable for use in steam or gas sterilization.



27640 A/B/C

- 27640 A **Aluminium Case**, with lid, for use with HOPKINS® Telescopes 27015 A/B/C/E, external dimensions (w x d x h): 405 x 73 x 55 mm
- 27640 B **Aluminium Case**, with lid, for use with HOPKINS® Telescopes 27018 A/B/C, external dimensions (w x d x h): 265 x 73 x 55 mm
- 27640 C **Aluminium Case**, with lid, for use with HOPKINS® Telescopes 10320 A/B/D, 12015 A, 26031 B, 24880 A, external dimensions (w x d x h): 620 x 78 x 55 mm

**Please note:** Before sterilization the foam insert is to be removed and the container be laid out with drapes.

## Carrying Cases for Instruments and Accessories



27677 A/B/C

- |         |   |
|---------|---|
| 27677 A | <b>Plastic Case</b> , without inserts, internal dimensions (w x d x h): 730 x 325 x 85 mm |
| 9200030 | <b>Case Insert</b> , grooved foam   |
| 27677 B | <b>Plastic Case</b> , without inserts, internal dimensions (w x d x h): 605 x 315 x 80 mm |
| 9211580 | <b>Case Insert</b> , grooved foam   |
| 27677 C | <b>Plastic Case</b> , without inserts, internal dimensions (w x d x h): 480 x 285 x 80 mm |
| 9210560 | <b>Case Insert</b> , foam material, grooved   |



# Components Spare Parts




# Grasping Forceps, Cleaning Pistol, Cleaning Tubing Set Gastro

Components /spare parts

Catalog page

**27662 Grasping Forceps**, with serrated rubber jaws

VET-S-CL 3



## Spare parts

27662 A  
**Rubber Jaws Insert**

Components /spare parts

Catalog page

**27660 P Cleaning Pistol**, only, autoclavable

VET-S-CL 6



## Spare parts

8021390  
**O-Rings, for nozzle**, package of 3

8021490  
**Spare Part Set, for measuring**

Components /spare parts

Catalog page

**13991 CK Cleaning Tubing Set Gastro**

VET-S-CL 8


13991 CD **Double Valve Seal**  
13991 CL **Leakage Tester Adaptor**  
13991 CS **Adaptor, for aspiration**  
13991 CR **Adaptor, air/water**  
13991 SS **Irrigation Tube**, reusable



13991 CD



13991 CL



13991 CS




13991 CR



13991 SS

# Wire Trays for Cleaning, Sterilization and Storage

Components /spare parts		Catalog page
39501 A1	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 11
39501 B1	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 11
39501 C	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 11



## Spare parts

7981790  
**Locking Piece**, with clamp, screw and nut


8021690  
**Silicone Telescope Holder**, for head of telescope

8021590  
**Silicone Telescope Holder**, straight

Components /spare parts		Catalog page
39501 A2	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 12
39501 B2	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 12
39501 BEC	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 13
39501 CEC	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 13
39501 X	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 14
39501 XK	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 14
39501 BR	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 15
39501 ES	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 16
39501 U	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 17
39501 V	Wire Tray for Cleaning, Sterilization and Storage	VET-S-CL 17

## Spare parts

7981790  
**Locking Piece**, with clamp, screw and nut




# Instrument Racks and Plastic Containers for Cleaning, Sterilization and Storage

Components /spare parts

Catalog page

- 39219 XF** Instrument Rack for Cleaning, Sterilization and Storage
- 39219 XFB** Instrument Rack for Cleaning, Sterilization and Storage

- VET-S-CL 23
- VET-S-CL 23



**Spare parts**

39219 XFS  
**Irrigation Tube**, with two LUER-Lock tube connectors, length 140 mm


39219 XFV  
**Sealing Cap**, package of 15, for multi-port bridge

Components /spare parts

Catalog page

- 39305 C1S** Plastic Container for Cleaning, Sterilization and Storage of Telescopes
- 39305 C2S** Plastic Container for Cleaning, Sterilization and Storage of Telescopes
- 39305 L1S** Plastic Container for Cleaning, Sterilization and Storage of Telescopes
- 39305 L2S** Plastic Container for Cleaning, Sterilization and Storage of Telescopes

- VET-S-CL 25
- VET-S-CL 25
- VET-S-CL 25
- VET-S-CL 25



**Spare parts**

Spare parts for 39305 C1S/C2S

Spare parts for 39305 L1S/L2S

39305 C1SH  
**Telescope Holder**

39305 L1SH  
**Telescope Holder**

39305 C2SH  
**Telescope Holder**


39305 L2SH  
**Telescope Holder**

Spare parts for 39305 C1S/C2S/L1S/L2S

39305 LS  
**Latch**

# Plastic Containers for Sterilization, Trays and Plastic Containers

Components /spare parts		Catalog page
39301 AS	Plastic Container for Sterilization	VET-S-CL 26
39301 BS	Plastic Container for Sterilization	VET-S-CL 26
39301 CS	Plastic Container for Sterilization	VET-S-CL 26
39301 D	Plastic Container for Sterilization and Storage of Endoscopes	VET-S-CL 26



## Spare parts

Spare parts for 39301 AS/BS/CS

39301 AHS  
**Silicone Telescope Holder**, up to diameter 5 mm

39301 CHS  
**Silicone Telescope Holder**, up to diameter 10 mm

5911202  
**Oval Head Screw**

39305 LS  
**Latch**

Spare parts for 39301 D

7920690  
**Rod**

Components /spare parts	Catalog page
<b>Trays and Plastic Containers, diverse</b>	—

## Spare parts

27695 AE  
**Strips of Silicone Knobs**, for cutting out individually sized pieces, external dimensions (w x d): 600 x 300 mm

# Plastic Containers for Sterilization and Storage

## Components /spare parts

## Catalog page

<b>39301 F</b>	<b>Plastic Container for Sterilization and Storage</b>
<b>39301 G</b>	<b>Plastic Container for Sterilization and Storage</b>
<b>39301 H</b>	<b>Plastic Container for Sterilization and Storage</b>
<b>39301 J</b>	<b>Plastic Container for Sterilization and Storage</b>

VET-S-CL 32
VET-S-CL 32
VET-S-CL 32
VET-S-CL 32



## Spare parts

Spare parts for 39301 F

39301 FB  
**Lid**

Spare Parts for 39301 H

39301 HB  
**Lid**

Spare Parts for 39301 G

39301 GB  
**Lid**

Spare Parts for 39301 J

39301 JB  
**Lid**

# Plastic Containers for Sterilization and Storage

Components /spare parts

Catalog page

**39360 AK** Plastic Container for Sterilization and Storage of Variable Instrument Sets

VET-S-CL 36

39360 AP Snap-in Clip, package of 12

39360 AS Silicone Tie-Downs, package of 12

39360 AR Tool



Components /spare parts

Catalog page


**39360 BK** Plastic Container for Sterilization and Storage of Variable Instrument Sets

VET-S-CL 36

39360 AP Snap-in Clip, package of 12

39360 AS Silicone Tie-Downs, package of 12

39360 AR Tool



# Sterilization Containers with MicroStop®

## Components /spare parts

## Catalog page

- 39740 A2 Sterilization Container, with MicroStop®**
- 39741 A2 Sterilization Container, with MicroStop®**
- 39742 A2 Sterilization Container, with MicroStop®**
- 39743 A2 Sterilization Container, with MicroStop®**

- VET-S-CL 40
- VET-S-CL 40
- VET-S-CL 40
- VET-S-CL 40



## Spare parts

- 39760 A2  
**MicroStop® Microbiological Barrier**

- 39745 AD  
**Silicone Seal, for lid 300 x 300**

## Components /spare parts

## Catalog page

- 39750 A2 Sterilization Container, with MicroStop®**
- 39751 A2 Sterilization Container, with MicroStop®**
- 39752 A2 Sterilization Container, with MicroStop®**
- 39753 A2 Sterilization Container, with MicroStop®**

- VET-S-CL 41
- VET-S-CL 41
- VET-S-CL 41
- VET-S-CL 41



## Spare parts

- 39760 A2  
**MicroStop® Microbiological Barrier**

- 39755 X  
**Centering, for lid seal in container 600 x 300**

- 39755 AD  
**Silicone Sealing Ring, for lid 600 x 300**



# MiniSet Containers with MicroStop®

Components /spare parts

Catalog page

- 39780 A**      **MiniSet Container, with MicroStop®**
- 39780 AB      **Base**
- 39785 AA      **Lid**
- 39760 A      **Microbiological Barrier, with MicroStop®**
  
- 39781 A**      **MiniSet-Container, with MicroStop®**
- 39781 AB      **Base**
- 39785 AA      **Lid**
- 39760 A      **Microbiological Barrier, with MicroStop®**

VET-S-CL 43



**Spare parts**

39760 A2  
**MicroStop® Microbiological Barrier**

39785 AD  
**Silicone Seal**





**COMPONENTS  
SPARE PARTS**




The chapter “Components / Spare Parts” contains detailed information on KARL STORZ instruments.

For easy location and reference, an index is available which lists the order number of the spare parts as well as those of the entire instrument, set or unit.

### Hotline

Queries concerning products, exchange, maintenance and cleaning can be addressed to the KARL STORZ EP1 Hotline: 07461/708-980, from Monday to Thursday from 7-18 h and Friday from 7-17 h.

### For example:

Components / Spare parts		Catalog page
<b>67260 OSA</b>	<b>Veterinary Otoscope Set</b>	3
67260 AA	<b>Otoscope</b>	
67262 S	<b>Stopcock Attachment</b>	
Components / Spare parts		Catalog page
<b>67071 XS</b>	<b>Myringotomy Needle</b>	4
67071 S	<b>Flexible Needle</b> , outer diameter 0.7 mm	
67479 S	<b>Outer Tube</b> , 1.5 mm	
Components / Spare parts		Catalog page
<b>67159 L</b>	<b>Polypectomy Snare</b>	4, 138
67159 LA	<b>Handle</b>	
67159 LB	<b>Tube</b>	
67159 LC	<b>Snare</b>	
67159 LD	<b>Adaptor</b>	

*Spare parts assigned to instrument with catalog page reference and order numbers for individual components/spare parts*

# Table of Contents, Numerical Index

## COMPONENTS/SPARE PARTS

### INSTRUMENTS

Protection Tubes, Veterinary Otoscope Set, Myringotomy Needle, Polypectomy Snare	SP 4
Operating Telescope 30°, Protection and Examination Sheath	SP 5
Operating Sheath, Sheaths, Rubber Insufflation Bulb, Cytology Brush	SP 6
Leakage Testers	SP 7
Suction Valve, Air/Water Valve, Water Bottle	SP 8
HOPKINS® Telescope	SP 9
Multifunctional Valves	SP 10
Automatic Valves	SP 11-SP 12
Cannulas	SP 13
Cannulas, Ribbed Cones	SP 14
Reducers	SP 15
Reducer, Thread Sleeve	SP 16
Thread Sleeve, Distance Holders	SP 17
VERESS Pneumoperitoneum Needle, Outer Sheaths	SP 18
Outer Sheaths, DUVAL Grasping Forceps, Fan Retractors	SP 19
TAKE APART® Bipolar Forceps	SP 20-SP 21
TAKE APART® Scissors	SP 21
Handles	SP 22
Suction and Irrigation Tube, Arthroscope Sheaths	SP 23
ARTHROFORCE® III TAKE-APART® Scissors	SP 24
ARTHROFORCE® III TAKE-APART® Forceps	SP 24-SP 26

Knife, Trocars	SP 27
Distractors	SP 28
Cystoscope-Urethroscope Sheath, Pediatric Cystoscope-Urethroscopes	SP 29
Cystoscope-Urethroscope Sheaths, Catheter Deflecting Mechanism, Telescope Bridges	SP 30
Minitube TCI Endoscope, Stone Basket, Operating Sheath	SP 31
POINT SETTER®, Air Compressors, STOR-E® Transport Cart	SP 32
ENDOCRANE®, Mobile OR-Rail	SP 33

### UNITS

Insufflators	
Electronic CO <sub>2</sub> ENDOFLATOR® SCB	SP 34-SP 35
CO <sub>2</sub> ENDO-ARTHROFLATOR-VET	SP 34-SP 35
Suction and Irrigation Systems	
ARTHROPUMP® SCB	SP 36-SP 38
VETPUMP®	SP 38
VETPUMP® II	SP 39-SP 40
Motor Systems	
UNIDRIVE® S III ARTHRO SCB	SP 41
High Frequency Surgery Units	
AUTOCON® II 400 SCB	SP 42
AUTOCON® II 200	SP 43
AUTOCON® II 80	SP 43
Surgery Electrodes Set	SP 44

Product	Page	Product	Page	Product	Page
031122-10	SP 34	20020061	SP 32	20300185	SP 38
031124-10	SP 39, SP 40	20020083	SP 32	20300580	SP 40
031128-01	SP 36	20031002	SP 32	20330088	SP 38
040146-25	SP 33	20031002 C	SP 32	20330185	SP 37
1059400	SP 39	20031020	SP 32	20330188	SP 37
1059700	SP 41	20031020 C	SP 32	20330386	SP 36, SP 37
11603 L1	SP 12	20090170	SP 34, SP 36, SP 41, SP 42	20400021	SP 35
1222890	SP 41	20100055	SP 8	20400022	SP 35
13242 LX	SP 7	2027390	SP 34, SP 36	20400027	SP 35
13242 XL	SP 7	2027590	SP 34, SP 38	20400028	SP 35
13266 H	SP 6	2027690	SP 42, SP 43	20400030	SP 34
13991 LV	SP 8	2028090	SP 42, SP 43	20400081	SP 35
13991 SSV	SP 8	20300034	SP 39, SP 40	20400143	SP 34, SP 35
13991 SV	SP 8	20300042	SP 39	20530008	SP 44
13991 SW	SP 8	20300052	SP 39	20530031	SP 44
13992 BS	SP 8	20300180	SP 39	20530801	SP 43
20020055	SP 32	20300181	SP 39	20530820	SP 43
20020057	SP 32	20300184	SP 35	20532201	SP 43

# Numerical Index

## COMPONENTS/SPARE PARTS

Product	Page	Product	Page	Product	Page
20532201 C	SP 43	2731291	SP 37	30103 S	SP 16
20532201-010	SP 43	2731391	SP 37	30103 T5	SP 13
20532201C010	SP 43	2731591	SP 37	30108 H2	SP 13
20532220	SP 43	2731691	SP 37	30108 M1	SP 10
20532220-010	SP 43	2731791	SP 37	30108 T4	SP 13
20532220C	SP 43	2732091	SP 38	30114 L1	SP 12
20532220C010	SP 43	2732191	SP 38	30117 A1	SP 11
20780020	SP 33	27500	SP 38, SP 39	30117 G4	SP 14
2386390	SP 40	27504	SP 29	30117 L1	SP 12
26034 AAV	SP 9	27550 A	SP 6, SP 29	30117 T2	SP 13
26034 V	SP 9	27550 C	SP 6, SP 30	30120 L1	SP 12
26173 BN	SP 23	27550 N	SP 29, SP 31	30120 T9	SP 13
26173 BT	SP 22, SP 23	2761291	SP 32	30123 L1	SP 11
26176 HS	SP 20	27630022	SP 32	30141 HB	SP 15
26176 HT	SP 20	27677 BV	SP 33	30160 A1	SP 11
26176 HU	SP 20	28083 G	SP 26	30160 C3	SP 14
26176 HX	SP 20	28083 L	SP 25, SP 26	30160 D	SP 17
26184 HCS	SP 21	28147 R	SP 27	30160 G1	SP 14
26184 HM	SP 21	28147 RM	SP 27	30160 G2	SP 14
26184 HSS	SP 21	28172 HG	SP 33	30160 G3	SP 14
26184 HVS	SP 21	28172 W	SP 32	30160 H1	SP 13
26196 HR	SP 20	28172 WC	SP 32	30160 H2	SP 13
26330042	SP 36	28172 WD	SP 32	30160 L1	SP 12
26430520-1	SP 34	28172 WKS	SP 32	30160 M1	SP 10
26520031	SP 44	28272 ECW	SP 33	30160 RE	SP 15
26520032	SP 44	28272 EH	SP 33	30160 S	SP 17
26520033	SP 44	28272 EHS	SP 33	30160 T4	SP 13
26520034	SP 44	28272 MC	SP 33	30160 T5	SP 13
26520035	SP 44	28272 MP	SP 33	30310 DBS	SP 19
26520036	SP 44	28272 MR	SP 33	30310 MDS	SP 19
26520037	SP 44	28272 MW	SP 33	30310 MWS	SP 19
26520038	SP 44	28272 UA	SP 32	30310 ULS	SP 19
26520039	SP 44	28272 UF	SP 32	30623 F	SP 19
26520040	SP 44	28272 UGN	SP 32	30623 FP	SP 19
26520041	SP 44	28272 UKN	SP 32	30804	SP 22
26520042	SP 44	28306 BN	SP 23	30805	SP 22
27001 G	SP 29, SP 31	28330020-1	SP 36	33133	SP 19
27001 GF	SP 29	28723001-1	SP 41	33300	SP 18
27001 GH	SP 29	28723020-1	SP 41	33400	SP 18
27001 RA	SP 29, SP 31	29005 AK	SP 32	33500	SP 18, SP 19
27021 O	SP 30	2903390	SP 35	33510 DU	SP 19
27023 VR	SP 31	2903490	SP 35	33533 DU	SP 19
27023 VS	SP 31	29100	SP 14, SP 18, SP 19, SP 21	38121	SP 21
27023 VT	SP 31	2911590	SP 40	38300	SP 21
27026 C	SP 30	30103 A1	SP 11	38310 MZ	SP 21
27026 D	SP 30	30103 C3	SP 14	38321 MZ	SP 21
27030 KA	SP 29	30103 D	SP 17	400 A	SP 33, SP 34, SP 36, SP 38, SP 39, SP 41, SP 42, SP 43
27030 KB	SP 29	30103 H1	SP 13	400 B	SP 34, SP 39, SP 43
2719791	SP 36, SP 37	30103 H2	SP 13	40924	SP 6
2720091	SP 36	30103 L1	SP 11	40924 B	SP 6
2731091	SP 37	30103 M1	SP 10	40924 S	SP 6
2731191	SP 37	30103 RE	SP 16	495 F	SP 5, SP 9, SP 29, SP 31

# Numerical Index

## COMPONENTS/SPARE PARTS

Product	Page	Product	Page	Product	Page
495 G	SP 5, SP 9, SP 29, SP 31	64069 G	SP 25	67159 LA	SP 2, SP 4
5500735	SP 7	64069 L	SP 24, SP 25	67159 LB	SP 2, SP 4
5903010	SP 27	64083 P	SP 26	67159 LC	SP 2, SP 4
5905098	SP 10	64147 BN	SP 23	67159 LD	SP 2, SP 4
5905610	SP 22	64166 A	SP 24	67260 AA	SP 2, SP 4
5907900	SP 7	64167 A	SP 26	67260 OSA	SP 2, SP 4
5925300	SP 8	64168 B	SP 24	67262 S	SP 2, SP 4
5925600	SP 8	64168 P	SP 24	67479 S	SP 2, SP 4
5925900	SP 8	64169 A	SP 24, SP 25, SP 26	6875990	SP 20
5934400	SP 8	64169 AS	SP 25, SP 26	69300133	SP 39
5936400	SP 7	64169 G	SP 25	69300233	SP 39
60007 BL	SP 6	64169 GS	SP 25	69311601	SP 39
60007 EL	SP 6	64169 L	SP 24	69311610	SP 39, SP 40
60007 SL	SP 6	64169 LS	SP 25	69321501	SP 38
6052191	SP 23	64169 X	SP 27	69321520	SP 38
6089590	SP 11	64183 G	SP 26	69321601	SP 39
61029 C	SP 5	64183 GS	SP 26	69321620	SP 39
61029 D	SP 6	64183 L	SP 25	69330004-1	SP 36
6126990	SP 27	64183 LS	SP 26	6958691	SP 31
6127090	SP 11, SP 15, SP 27	64183 P	SP 26	6979091	SP 22
6127290	SP 12	64183 X	SP 27	6985691	SP 5, SP 6, SP 9, SP 13, SP 18, SP 23, SP 29, SP 30
6127390	SP 12	64302 BN	SP 23	7072390	SP 21
6127590	SP 9, SP 10, SP 11, SP 12, SP 15, SP 16	64302 X	SP 27	7080690	SP 10
6127690	SP 14, SP 17	64303 BN	SP 23	7228790	SP 16
62114 G1	SP 14	64527	SP 27	723750 A	SP 4
62120 J	SP 18	64527 S	SP 27	723750 B	SP 4
62141 AA	SP 15	64530	SP 27	723750 E	SP 4
62141 DB	SP 15	64535	SP 27	723750 H	SP 4
62430508-1	SP 34	64535 S	SP 27	7302990	SP 10
62432501	SP 34	64820 LSD	SP 28	7342490	SP 27
62432520	SP 34	64821 C	SP 28	7351090	SP 16, SP 17
62603 G1	SP 14	64821 D	SP 28	7351190	SP 16, SP 17
6266690	SP 22	64821 H	SP 28	7351490	SP 16, SP 17
6266691	SP 22	64821 S	SP 28	7385590	SP 17
63002 LA	SP 31	6526295	SP 28	7391190	SP 14, SP 16
63025 G	SP 30	6526795	SP 28	7415690	SP 10
63026 A	SP 30	6526895	SP 28	7416090	SP 10
63026 B	SP 30	6573990	SP 5	7430790	SP 17
63026 E	SP 30	66176 HA	SP 20	7616690	SP 10, SP 11
63026 U	SP 30	66176 HB	SP 20	7720590	SP 9, SP 12, SP 15
63027 E	SP 30	66176 HL	SP 20	7909890	SP 18
63027 KL	SP 30	6626690	SP 30	7944790	SP 18
63027 NL	SP 30	6681291	SP 19	7962190	SP 11, SP 16
63027 UL	SP 30	6681391	SP 19	8259290	SP 27
63054 SJ	SP 31	67023 VK	SP 31	8407390	SP 14
6353591	SP 23	67030 BA	SP 5	8415090	SP 13
6376990	SP 9, SP 18, SP 31	67031 E	SP 29	8458190	SP 13, SP 29, SP 30, SP 31
6377091	SP 30, SP 31	67065 C	SP 6	8464790	SP 14
64066 A	SP 24	67065 CV	SP 6	8515090	SP 5, SP 6, SP 23, SP 29
64067 A	SP 26	67071 S	SP 2, SP 4	90100043	SP 8
64068 B	SP 24	67071 XS	SP 2, SP 4	96206358DF	SP 32
64068 P	SP 24	67159 L	SP 2, SP 4		

# Protection Tubes, Veterinary Otoscope Set, Myringotomy Needle, Polypectomy Snare

Components / Spare parts Catalog page

**723750 A**    **Protection Tube**, for HOPKINS® telescopes with length of 6 – 11 cm –



Components / Spare parts Catalog page

**723750 B**    **Protection Tube**, for use with HOPKINS® telescopes with length 18 cm –



Components / Spare parts Catalog page

**723750 E**    **Protection Tube**, for HOPKINS® II telescope with length 30 cm –

Components / Spare parts Catalog page

**723750 H**    **Protection Tube**, for use with HOPKINS® telescopes up to outer diameter 4 mm and working length 44 cm –

Components / Spare parts Catalog page

**67260 OSA**    **Veterinary Otoscope Set** 3

67260 AA    **Otoscope**

67262 S    **Stopcock Attachment**



Components / Spare parts Catalog page

**67071 XS**    **Myringotomy Needle** 4

67071 S    **Flexible Needle**, outer diameter 0.7 mm

67479 S    **Outer Tube**, 1.5 mm



Components / Spare parts Catalog page

**67159 L**    **Polypectomy Snare** 4, 138

67159 LA    **Handle**

67159 LB    **Tube**

67159 LC    **Snare**

67159 LD    **Adaptor**






# Operating Telescope 30°, Protection and Examination Sheath

Components / Spare Parts

Catalog page

67030 BA Operating Telescope 30°

6, 122, 136



## Spare Parts



495 F  
**Receptacle**,  
diameter 9 mm, for Wolf  
fiber optic light cable



6985691  
**Spring Cap**



495 G  
**Screw Base**,  
for KARL STORZ fiber  
optic light cable and  
Olympus/Winter & Ibe



8515090  
**Stopcock**




6573990  
**Sealing Cap**, diameter  
0.8 mm

Components / Spare Parts

Catalog page

61029 C Protection and Examination Sheath

6, 137



## Spare Parts



6985691  
**Spring Cap**



8515090  
**Stopcock**


# Operating Sheath, Sheaths, Rubber Insufflation Bulb, Cytology Brush

Components / Spare Parts

Catalog page

**61029 D**      **Operating Sheath**

6, 137



**Spare Parts**



**27550 A**  
**Sealing Cap**, drill hole  
diameter 0.8 mm



**8515090**  
**Stopcock**



**6985691**  
**Spring Cap**

Components / Spare Parts


Catalog page

**67065 C**      **Sheath**

7, 122

**67065 CV**     **Sheath**

137



**Spare Parts**



**27550 C**  
**Sealing Cap**, drill hole  
diameter 1.2 mm



**8515090**  
**Stopcock**



**6985691**  
**Spring Cap**

Components / Spare parts

Catalog page

**40924**      **Rubber Insufflation Bulb**

17

40924 B      **Rubber Bulb**

40924 S      **Rubber Tube**



Components / Spare parts

Catalog page

**60007 BL**     **Cytology Brush**

28

60007 EL     **Brush**

60007 SL     **Outer Tube**

13266 H      **Handle**



# Leakage Testers

Components / Spare Parts

Catalog page

**13242 XL**    **Leakage Tester**

22, 26-27



## Spare Parts



5500735  
**Silicone Tube**



5907900  
**O-Ring**

Components / Spare Parts

Catalog page

**13242 LX**    **Leakage Tester**

24, 32



## Spare Parts



5500735  
**Silicone Tube**



5936400  
**O-Ring, Silicone, red**


# Suction Valve, Air/Water Valve, Water Bottle

Components / Spare Parts

Catalog page

**13991 LV Suction Valve**

32



## Spare Parts



5925600  
**O-Ring**, for suction valve,  
large, red, 2.5 x 1.5 mm



5925300  
**O-Ring**, for suction valve,  
small, red, 2 x 1.5 mm

Components / Spare Parts


Catalog page

**13991 SV Air/Water Valve**

32

**13991 SSV Air/Water Valve**

32



## Spare Parts



5925900  
**O-Ring**, red, 3.1 x 1.6 mm


Components / Spare Parts

Catalog page

**13992 BS Water Bottle**

33

90100043 **Bottle Cap**  
20100055 **Irrigation Bottle**



## Spare Parts



13991 SW  
**Irrigation Adaptor**,  
for water bottle



5934400  
**O-Ring**, water bottle seal,  
6.1 x 1.6 mm


## HOPKINS® Telescope

Components / Spare Parts

Catalog page

26034 AAV HOPKINS® Straight Forward Telescope 0°

52, 87



### Spare Parts



26034 V  
**Adaptor**, for working channel, for use with thorascopes diameter 10 mm



6376990  
**Stopcock**, for irrigation connector



495 F  
**Receptacle**, diameter 9 mm, for Wolf fiber optic light cable



6985691  
**Spring Cap**



495 G  
**Screw Base**, for KARL STORZ fiber optic light cable and Olympus/Winter & Ibe



7720590  
**Silicone Leaflet Washer**



6127590  
**Sealing Cap**, (50/4)


# Multifunctional Valves

Components / Spare Parts

Catalog page

**30103 M1 Multifunctional Valve**

56, 61



**Spare Parts**



5905098  
**Sealing Cap, small**



7302990  
**Tappet**



7616690  
**Sealing Cap, (60/10)**

Components / Spare Parts

Catalog page

**30160 M1 Multifunctional Valve**

58, 63



**Spare Parts**



5905098  
**Sealing Cap, small**



7080690  
**Tappet**




6127590  
**Sealing Cap, (50/4)**

Components / Spare Parts

Catalog page

**30108 M1 Multifunctional Valve**

56, 60



**Spare Parts**



5905098  
**Sealing Cap, small**



7415690  
**Tappet**



7416090  
**Sealing Cap, (60/13)**

# Automatic Valves

Components / Spare Parts Catalog page

**30117 A1 Automatic Valve**

59



**Spare Parts**



6089590  
**Seal, mushroom-shaped**



6127090  
**Sealing Cap, olive, (40/3)**

Components / Spare Parts Catalog page

**30103 A1 Automatic Valve**

61



**Spare Parts**



7616690  
**Sealing Cap, (60/10)**

Components / Spare Parts Catalog page

**30160 A1 Automatic Valve**

63



**Spare Parts**



6127590  
**Sealing Cap, (50/4)**

Components / Spare Parts Catalog page

**30103 L1 Automatic Valve**

56

**30123 L1 Automatic Valve**


62



**Spare Parts**



7616690  
**Sealing Cap, (60/10)**



7962190  
**Silicone Leaflet Washer, large**

# Automatic Valves

Components / Spare Parts

Catalog page

**30120 L1**    **Automatic Valve**  
**30160 L1**    **Automatic Valve**


57, 65  
64



Spare Parts



6127590  
**Sealing Cap, (50/4)**



7720590  
**Silicone Leaflet Washer**

Components / Spare Parts

Catalog page

**30117 L1**    **Automatic Valve**  
**30114 L1**    **Automatic Valve**


68  
68



Spare Parts



6127390  
**Sealing Cap, (50/2.6)**



7720590  
**Silicone Leaflet Washer**

Components / Spare Parts

Catalog page

**11603 L1**    **Automatic Valve**


69



Spare Parts



6127290  
**Sealing Cap, (50/2.2)**



7720590  
**Silicone Leaflet Washer**




# Cannulas

Components / Spare Parts

Catalog page

<b>30103 T5</b>	<b>Cannula</b>		56
<b>30108 T4</b>	<b>Cannula</b>		56
<b>30120 T9</b>	<b>Cannula</b>		57
<b>30160 T4</b>	<b>Cannula</b>		58
<b>30160 T5</b>	<b>Cannula</b>		58



Spare Parts



6985691  
**Spring Cap**




8458190  
**Stopcock**

Components / Spare Parts

Catalog page

<b>30117 T2</b>	<b>Cannula</b>		59
-----------------	----------------	--	----



Spare Parts



6985691  
**Spring Cap**




8415090  
**Stopcock**

Components / Spare Parts

Catalog page

<b>30103 H1</b>	<b>Cannula</b>		61
<b>30103 H2</b>	<b>Cannula</b>		61
<b>30108 H2</b>	<b>Cannula</b>		60
<b>30160 H1</b>	<b>Cannula</b>		63
<b>30160 H2</b>	<b>Cannula</b>		63



Spare Parts



6985691  
**Spring Cap**



8458190  
**Stopcock**

# Cannulas, Ribbed Cones

Components / Spare Parts		Catalog page
30160 G1	Cannula	64
30160 G2	Cannula	64
30160 G3	Cannula	64
30117 G4	Cannula	68
62114 G1	Cannula	68
62603 G1	Cannula	69



29100 —

## Spare Parts



29100  
**Plug**, for LUER-Lock connector for cleaning, black, autoclavable, package of 10

Components / Spare Parts		Catalog page
30160 C3	Ribbed Cone	66



8464790 —

— 6127690

## Spare Parts



8464790  
**Cone**



6127690  
**Sealing Cap, (50/5)**

Components / Spare Parts		Catalog page
30103 C3	Ribbed Cone	66



8407390 —

— 7391190

## Spare Parts



8407390  
**Cone**



7391190  
**Sealing Cap, (60/14)**


## Reducers

Components / Spare Parts

Catalog page

**62141 AA Reducer**

70



### Spare Parts



6127090  
**Sealing Cap, olive, (40/3)**

Components / Spare Parts

Catalog page

**62141 DB Reducer**

70

**30141 HB Reducer**

70



### Spare Parts



6127590  
**Sealing Cap, (50/4)**

Components / Spare Parts

Catalog page

**30160 RE Reducer**

70



### Spare Parts



6127590  
**Sealing Cap, (50/4)**



7720590  
**Silicone Leaflet Washer**



6127090  
**Sealing Cap, olive, (40/3)**

# Reducer, Thread Sleeve

Components / Spare Parts

Catalog page

**30103 RE Reducer**

70



**Spare Parts**



6127590  
**Sealing Cap, (50/4)**



7962190  
**Silicone Leaflet Washer, large**




7228790  
**Sealing Cap, (60/11)**

Components / Spare Parts

Catalog page

**30103 S Thread Sleeve**

70



**Spare Parts**



7391190  
**Sealing Cap, (60/14)**



7351090  
**Metal Inner Ring**



7351190  
**Plastic Inner Ring**



7351490  
**Knurled Screw**

# Thread Sleeve, Distance Holders

Components / Spare Parts

Catalog page

**30160 S Thread Sleeve**



70

## Spare Parts



6127690  
**Sealing Cap, (50/5)**



7385590  
**Metal Inner Ring**



7430790  
**Plastic Inner Ring**



7351490  
**Knurled Screw**

Components / Spare Parts

Catalog page

**30103 D Distance Holder**



70

## Spare Parts



7351090  
**Metal Inner Ring**



7351490  
**Knurled Screw**



7351190  
**Plastic Inner Ring**

Components / Spare Parts

Catalog page

**30160 D Distance Holder**



70

## Spare Parts



7385590  
**Metal Inner Ring**



7351490  
**Knurled Screw**



7430790  
**Plastic Inner Ring**


# VERESS Pneumoperitoneum Needle, Outer Sheaths

Components / Spare Parts

Catalog page

**62120 J** VERESS Pneumoperitoneum Needle

69, 85, 87, 88, U 10



## Spare Parts



6985691  
**Spring Cap**



6376990  
**Stopcock, for irrigation connector**

Components / Spare Parts

Catalog page

**33300** CLICK'line® Metal Outer Sheath, insulated

74

**33400** CLICK'line® Metal Outer Sheath, insulated

74



## Spare Parts



29100  
**Plug, for LUER-Lock connector for cleaning, black, autoclavable, package of 10**



7909890  
**Halfround Plate**




7944790  
**Silicone Ring**

Components / Spare Parts

Catalog page

**33500** CLICK'line® Metal Outer Sheath, insulated

78



## Spare Parts



29100  
**Plug, for LUER-Lock connector for cleaning, black, autoclavable, package of 10**



7909890  
**Halfround Plate**



7944790  
**Silicone Ring**

# Outer Sheaths, DUVAL Grasping Forceps, Fan Retractors

Components / Spare Parts Catalog page

<b>30310 MWS</b>	<b>CLICK'line</b> Outer Sheath with Scissors Insert	77
<b>30310 DBS</b>	<b>CLICK'line</b> BLAKESLEY Outer Sheath with Forceps Insert	77
<b>30310 MDS</b>	<b>CLICK'line</b> KELLY Outer Sheath with Forceps Insert	77
<b>30310 ULS</b>	<b>CLICK'line</b> REDDICK-OLSEN Outer Sheath with Forceps Insert	77




**Spare Parts**



29100  
**Plug**, for LUER-Lock connector for cleaning, black, autoclavable, package of 10


Components / Spare parts Catalog page

<b>33533 DU</b>	<b>CLICK'line</b> DUVAL Grasping Forceps	78
33133	<b>Metal Handle</b> , with hemostat style ratchet	
33500	<b>Outer Sheath</b> , insulated	
33510 DU	<b>Forceps Insert</b>	




Components / Spare Parts Catalog page

<b>30623 F</b>	<b>Fan Retractor</b>	78
----------------	----------------------	----



**Spare Parts**




6681391  
**Outer Sheath**

Components / Spare Parts Catalog page

<b>30623 FP</b>	<b>Fan Retractor</b>	78
-----------------	----------------------	----



**Spare Parts**



6681291  
**Outer Sheath**

# TAKE APART® Bipolar Forceps


Components / Spare parts

Catalog page

**66176 HA TAKE-APART® Bipolar Coagulating Forceps**

79, 86

- 26196 HR **Handle**
- 26176 HS **Outer Sheath**
- 26176 HT **Forceps Insert**




Components / Spare parts

Catalog page

**66176 HB TAKE-APART® Bipolar Micro Forceps**

79

- 26196 HR **Handle**
- 26176 HS **Outer Sheath**
- 26176 HU **Forceps Insert**




Components / Spare parts

Catalog page

**66176 HL TAKE-APART® Bipolar Coagulating Forceps**

79

- 26196 HR **Handle**
- 26176 HS **Outer Sheath**
- 26176 HX **Forceps Insert**



Components / Spare parts

Catalog page

**26196 HR Bipolar Spring Handle**

79



## Spare Parts



6875990  
**Sealing Cap,**  
diameter 2.6 mm



# TAKE APART® Bipolar Forceps, Scissors

Components / Spare parts

Catalog page

**26184 HCS TAKE-APART® Bipolar Coagulating Forceps**

80, 88

26184 HM **Bipolar Ring Handle**

26184 HSS **Outer Sheath**

26184 HVS **Forceps Insert**



Components / Spare Parts

Catalog page

**26184 HM TAKE-APART® Bipolar Ring Handle**

80, 88



Spare Parts



7072390  
**Sealing Cap,**  
diameter 1.8 mm

Components / Spare parts

Catalog page


**38321 MZ RoBi® Scissors, CLERMONT-FERRAND model**

81

38121 **RoBi® Ring Handle**

38300 **RoBi® Metal Outer Sheath, insulated**

38310 MZ **RoBi® Scissors Insert**



Components / Spare Parts

Catalog page

**38300 RoBi® Metal Outer Sheath, insulated**

81



Spare Parts



29100  
**Plug,** for LUER-Lock  
connector for cleaning,  
black, autoclavable,  
package of 10


# Handles

Components / Spare Parts

Catalog page

**30804 Handle with Trumpet Valve**

90



## Spare Parts



5905610  
**Spring**



6266691  
**Piston**




6266690  
**Knurled Cap**

Components / Spare Parts

Catalog page

**30805 Handle with Two-Way Stopcock**

90



## Spare Parts



26173 BT  
**Stopcock**



6979091  
**Spring Cap**

# Suction and Irrigation Tube, Arthroscope Sheaths

Components / Spare Parts

Catalog page

**26173 BN Suction and Irrigation Tube**

78, U 18



## Spare Parts



**26173 BT  
Stopcock**



**6353591  
Spring Cap**

Components / Spare Parts

Catalog page

**64302 BN High-Flow Arthroscope Sheath**

96

**28306 BN High-Flow Arthroscope Sheath**


96

**64303 BN High-Flow Arthroscope Sheath**

98

**64147 BN Arthroscope Sheath**

102



## Spare Parts



**6052191  
Locking Lever, red**



**8515090  
Stopcock**



**6985691  
Spring Cap**

# ARTHROFORCE® III TAKE-APART® Forceps and Scissors


Components / Spare parts

Catalog page

**64168 P**      **ARTHROFORCE® III TAKE-APART® Forceps**

106

64169 A      **Handle, without ratchet**  
64068 P      **Working Attachment**




Components / Spare parts

Catalog page

**64166 A**      **BLAKESLEY ARTHROFORCE® III TAKE-APART® Forceps**

106

64169 A      **Handle, without ratchet**  
64066 A      **Working Attachment**




Components / Spare parts

Catalog page

**64168 B**      **ARTHROFORCE® III TAKE-APART® Scissors**

106

64169 A      **Handle, without ratchet**  
64068 B      **Working Attachment**




Components / Spare parts

Catalog page

**64169 L**      **ARTHROFORCE® III TAKE-APART® Grasping Forceps**

107

64169 A      **Handle, without ratchet**  
64069 L      **Working Attachment**



# ARTHROFORCE® III TAKE-APART® Forceps


Components / Spare parts

Catalog page

**64169 LS**    **ARTHROFORCE® III TAKE-APART® Grasping Forceps**

107

64169 AS    **Handle, with ratchet**  
64069 L    **Working Attachment**




Components / Spare parts

Catalog page

**64169 G**    **ARTHROFORCE® III TAKE-APART® Grasping Forceps**

107

64169 A    **Handle, without ratchet**  
64069 G    **Working Attachment**




Components / Spare parts

Catalog page

**64169 GS**    **ARTHROFORCE® III TAKE-APART® Grasping Forceps**

107

64169 AS    **Handle, with ratchet**  
64069 G    **Working Attachment**




Components / Spare parts

Catalog page

**64183 L**    **ARTHROFORCE® III TAKE-APART® Grasping Forceps**

108

64169 A    **Handle, without ratchet**  
28083 L    **Working Attachment**



# ARTHROFORCE® III TAKE-APART® Forceps

Components / Spare parts

Catalog page

**64183 LS** ARTHROFORCE® III TAKE-APART® Grasping Forceps

108

64169 AS **Handle, with ratchet**  
28083 L **Working Attachment**




Components / Spare parts

Catalog page

**64183 G** ARTHROFORCE® III TAKE-APART® Grasping Forceps

108

64169 A **Handle, without ratchet**  
28083 G **Working Attachment**




Components / Spare parts

Catalog page

**64183 GS** ARTHROFORCE® III TAKE-APART® Grasping Forceps

108

64169 AS **Handle, with ratchet**  
28083 G **Working Attachment**



Components / Spare parts

Catalog page

**64183 P** ARTHROFORCE® III TAKE-APART® Forceps

108

64169 A **Handle, without ratchet**  
64083 P **Working Insert**




Components / Spare parts

Catalog page

**64167 A** BLAKESLEY ARTHROFORCE® III TAKE-APART® Forceps

109

64169 A **Handle, without ratchet**  
64067 A **Working Attachment**



## Knife, Trocars

Components / Spare Parts

Catalog page

**28147 R Knife**



112

Spare Parts



28147 RM  
**Knife Attachment,**  
backward cutting



7342490  
**Cap**



5903010  
**Spring**

Components / Spare parts

Catalog page

**64527 S Hand Burr**



113

64527 **Hand Burr**

64530 **Hand Burr Sleeve**

Components / Spare parts

Catalog page

**64535 S Hand Burr**



113

64535 **Hand Burr**

64530 **Hand Burr Sleeve**

Components / Spare parts

Catalog page

**64302 X Trocar**



111

Spare Parts



6126990  
**Sealing Cap,**  
(40/1.5)

Components / Spare parts

Catalog page

**64169 X Trocar**



114

Spare Parts



6127090  
**Sealing Cap, olive,** (40/3)

Components / Spare parts

Catalog page

**64183 X Trocar**



114

Spare Parts




8259290  
**Sealing Cap**

# Distractors

Components / Spare parts

Catalog page

- 64820 LSD** LEIPZIG Stifle Distractor
- 64821 C **Distractor**
- 64821 S **Fixation Screw**, package of 6
- 64821 H **Lever Arm**
- 64821 D **Screwdriver**




116

Components / Spare parts

Catalog page

- 64821 C** **Distractor**

116



## Spare Parts

6526295  
**Swivel Nut**

6526895  
**Knurled Screw**

6526795  
**Drive Mechanism**



# Cystoscope-Urethroscope Sheath, Pediatric Cystoscope-Urethroscopes

Components / Spare Parts

Catalog page

**67031 E Cystoscope-Urethroscope Sheath**

122



## Spare Parts



27550 A  
**Sealing Cap**, drill hole  
diameter 0.8 mm



8515090  
**Stopcock**



6985691  
**Spring Cap**

Components / Spare parts


Catalog page

**27030 KA Pediatric Cystoscope-Urethroscope**

123

**27030 KB Pediatric Cystoscope-Urethroscope**

124



## Spare Parts



27001 G  
**Instrument Port with  
Sealing System and  
Quick Release Lock**,  
1 channel



27550 N  
**Seal**, for working channel,  
package of 10, single use  
recommended



27001 GF  
**Instrument Port with  
Sealing System and  
Quick Release Lock**,  
2 channels



495 F  
**Receptacle**,  
diameter 9 mm, for Wolf  
fiber optic light cable



27001 GH  
**Instrument Port with  
Sealing System and  
Quick Release Lock**,  
2 channels, 1 straight  
channel, 1 lateral channel



495 G  
**Screw Base**,  
for KARL STORZ fiber  
optic light cable and  
Olympus/Winter & Ibe



27001 RA  
**Cleaning Adaptor**,  
for Working Channel  
27001 G/GF/GH



6985691  
**Spring Cap**



27504  
**Flow Control Stopcock**



8458190  
**Stopcock**

# Cystoscope-Urethroscope Sheaths, Catheter Deflecting Mechanism, Telescope Bridges

## Components / Spare Parts

Catalog page

63026 A	Cystoscope-Urethroscope Sheath	126
63026 B	Cystoscope-Urethroscope Sheath	126
27026 C	Cystoscope-Urethroscope Sheath	126
27026 D	Cystoscope-Urethroscope Sheath	126
63026 U	Cystoscope-Urethroscope Sheath	126
63027 KL	Cystoscope-Urethroscope Sheath	127
63027 UL	Cystoscope-Urethroscope Sheath	127



## Spare Parts



6985691  
**Spring Cap**



8458190  
**Stopcock**


## Components / Spare Parts

Catalog page

63026 E	Catheter Deflecting Mechanism	126
63027 E	Catheter Deflecting Mechanism	127



## Spare Parts



27021 O  
**Rod**



6377091  
**Spring Cap**



27550 C  
**Sealing Cap**, drill hole  
diameter 1.2 mm




6626690  
**Stopcock**

## Components / Spare Parts

Catalog page

63025 G	Telescope Bridge	126
63027 NL	Telescope Bridge	127



## Spare Parts



27550 C  
**Sealing Cap**, drill hole  
diameter 1.2 mm



6626690  
**Stopcock**



6377091  
**Spring Cap**

# Minitube TCI Endoscope, Stone Basket, Operating Sheath

Components / Spare parts

Catalog page

**63002 LA Minitube TCI Endoscope**

128



## Spare Parts



27001 G  
**Instrument Port with Sealing System and Quick Release Lock**, 1 channel



495 G  
**Screw Base**, for KARL STORZ fiber optic light cable and Olympus/Winter & Ibe



27001 RA  
**Cleaning Adaptor**, for Working Channel 27001 G/GF/GH



6958691  
**Spring Cap**



27550 N  
**Seal**, for working channel, package of 10, single use recommended



8458190  
**Stopcock**



495 F  
**Receptacle**, diameter 9 mm, for Wolf fiber optic light cable

Components / Spare parts

Catalog page


**67023 VK Stone Basket**, flexible, 5 Fr., length 60 cm

129, 138

27023 VT **3-Ring Handle**

27023 VR **Basket**

27023 VS **Coil**




Components / Spare parts

Catalog page

**63054 SJ Operating Sheath**

132



## Spare Parts



6376990  
**Stopcock**, for irrigation connector



6377091  
**Spring Cap**

# POINT SETTER®, Air Compressors, STOR-E® Transport Cart

## Components / Spare parts

Catalog page

<b>28172 WKS</b>	<b>POINT SETTER®</b> , pneumatic holding arm
28172 W	<b>POINT SETTER® Arm</b> , with OR table adaptor
28272 UA	<b>KSLOCK Adaptor</b> , for clamping jaws
28272 UGN	<b>Clamping Jaw</b> , large
28272 UKN	<b>Clamping Jaw</b> , small
28272 UF	<b>Clamping Jaw</b> , for fiberscopes
2761291	<b>Pressure Regulator</b> , 7 bar
28172 WD	<b>Drape</b> , for single use, sterile, package of 10



148

## Components / Spare Parts

Catalog page

<b>20031002</b>	<b>Air Compressor</b> , for POINT SETTER®
20031020	<b>Compressor</b> , 0 – 8 bar, power supply 230 VAC, 50/60 Hz
27630022	<b>Fabric Reinforced Tube</b> , for CALCUSPLIT® SCB
28172 WC	<b>Connecting Tube</b> , for POINT SETTER®, length 600 cm



20031002

149

## Components / Spare Parts

Catalog page

<b>20031002 C</b>	<b>Air Compressor</b> , for POINT SETTER®
20031020 C	<b>Compressor</b> , 0 – 8 bar, power supply 115 VAC, 50/60 Hz
27630022	<b>Fabric Reinforced Tube</b> , for CALCUSPLIT® SCB
28172 WC	<b>Connecting Tube</b> , for POINT SETTER®, length 600 cm



20031002 C

149

## Components / Spare parts

Catalog page

<b>20020083</b>	<b>STOR-E® Transport Cart</b> , for POINT SETTER® holding arm
20020061	<b>Subrack</b> , for POINT SETTER®
20020055	<b>Accessory Set</b> , for POINT SETTER®
20020057	<b>Top Cover</b>
29005 AK	<b>Wire Basket</b>
96206358DF	<b>Manual</b>



20020055

150

## ENDOCRANE®, Mobile OR-Rail

Components / Spare Parts

Catalog page

- 28272 EH**    **ENDOCRANE®**, piezoregulated holding arm
- 28272 EHS    **ENDOCRANE® Arm**, including stand
- 28172 HG    **Socket**, to clamp to the OR table
- 20780020    **Control Unit**
- 040146-25\*    **Cover**, sterile, package of 25
- 28272 ECW    **Spring Balance**
- 400 A    **Mains Cord**, length 300 cm
- 27677 BV    **Case**



151


28272 EHS

Components / Spare parts

Catalog page

- 28272 MR**    **Mobile OR-Rail**
- 28272 MP    **Holding Plate**
- 28272 MC    **C-Clamp**
- 28272 MW    **Stop-Block**

155



28272 MP

28272 MW

28272 MC

# Insufflators

Electronic CO<sub>2</sub> ENDOFLATOR® SCB, CO<sub>2</sub> ENDO-ARTHROFLATOR-VET

Components / Spare parts

Catalog page

- 62430508-1 Electronic CO<sub>2</sub> ENDOFLATOR® SCB**
- 26430520-1 **Electronic CO<sub>2</sub> ENDOFLATOR® with KARL STORZ-SCB,**  
power supply 100 – 240 VAC, 50/60 Hz
- 400 A **Mains Cord**
- 20400143 **Silicone Tubing Set,** sterilizable
- 20400030 **Universal Wrench**
- 20090170 **SCB Connecting Cable,** length 100 cm
- 031122-10\* **CO<sub>2</sub>/N<sub>2</sub>O Gas Filter,** sterile, for single use,  
package of 10



62430520-1

## Spare Parts



2027590  
**Mains Fuse,**  
T 2.0 AL (SB),  
package of 10

Components / Spare parts

Catalog page

- 62432501 CO<sub>2</sub> ENDO-ARTHROFLATOR-VET**
- 62432520 **CO<sub>2</sub> ENDO-ARTHROFLATOR-VET,**  
power supply 100 – 240 VAC, 50/60 Hz
- 400 A **Mains Cord**
- 400 B **Mains Cord,** US-version
- 20400143 **Silicone Tubing Set,** sterilizable
- 20400030 **Universal Wrench**
- 031122-10\* **CO<sub>2</sub>/N<sub>2</sub>O Gas Filter,** sterile, for single use,  
package of 10



62432520

## Spare Parts



2027390  
**Mains Fuse,**  
T 1.0 A (SB),  
package of 10



mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# Insufflators

Electronic CO<sub>2</sub> ENDOFLATOR® SCB, CO<sub>2</sub> ENDO-ARTHROFLATOR-VET

## Components / Spare Parts

Catalog page

20400022 CO<sub>2</sub> High Pressure Tube

U 10

20400028 CO<sub>2</sub> High Pressure Tube

U 10



## Spare Parts



2903390  
**Seal**, for use with  
CO<sub>2</sub> bottle, Pin-Index  
connector

## Components / Spare Parts

Catalog page

20400021 CO<sub>2</sub> High Pressure Tube

U 10

20400027 CO<sub>2</sub> High Pressure Tube

U 10



## Spare Parts



2903490  
**Seal**, for use with  
CO<sub>2</sub> bottle, German  
connection

## Components / Spare Parts

Catalog page

20400143 **Silicone Tubing Set**, for insufflation

U 10



## Spare Parts



20400081  
**Plastic Connector**, male



20300184  
**LUER-Lock Tube  
Connector**, male, tube  
diameter 10 mm

# Suction and Irrigation Systems

ARTHROPUMP® SCB

Components / Spare parts

Catalog page

- 69330004-1 ARTHROPUMP® SCB**
- 28330020-1 **ARTHROPUMP® with KARL STORZ-SCB,**  
power supply 100 – 240 VAC, 50/60 Hz
- 400 A **Mains Cord**
- 031128-01\* **Tubing Set,** for single use
- 20090170 **SCB Connecting Cable,** length 100 cm



U 12

28330020-1

## Spare Parts



2027390  
**Mains Fuse,**  
T 1.0 A (SB),  
package of 10

Components / Spare Parts

Catalog page

- 26330042 ARTHRO Tubing Set**


U 12



## Spare Parts



2719791  
**Pressure Dome ARTHRO,**  
without silicone tubing set



20330386  
**Pump Tubing ARTHRO,**  
without connectors,  
package of 25

2720091  
**Silicone Tubing Set**  
**ARTHRO,** without  
pressure dome



mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany



# Suction and Irrigation Systems

ARTHROPUMP® SCB

Components / Spare Parts

Catalog page

**2719791 Pressure Dome ARTHRO**

U 12



## Spare Parts



2731091  
**Lock Ring ARTHRO**



2731591  
**Compression Spring ARTHRO/HYST**



2731191  
**Dome**



2731691  
**Plastic Disk**



2731291  
**Dome Ring**



20330185  
**Diaphragm for Pressure Dome, package of 10**



2731391  
**Sealing Disk, for pressure relief valve**



20330188  
**O-Ring, package of 10**

Components / Spare Parts

Catalog page

**20330386 Pump Tubing ARTHRO/LAP**

U 12



## Spare Parts



2731791  
**Connector Set LAP**

# Suction and Irrigation Systems

ARTHROPUMP® SCB, VETPUMP®

Components / Spare Parts

Catalog page

## Patient Tube

U 12



## Spare Parts



20300185  
**Tube Connector**



27500  
**LUER-Lock Tube Connector, male,**  
tube diameter 9 mm

Components / Spare Parts

Catalog page

## Inflow Tube, with Y-adaptor


U 12



## Spare Parts



20330088  
**Puncture Needle,**  
for irrigation bottle



2732091  
**Tube Clip, white,**  
package of 10



2732191  
**Y-Adaptor,** package of 10

Components / Spare parts

Catalog page

**69321501 VETPUMP®**

U 14

69321520 **VETPUMP®,**  
power supply 230 VAC, 50/60 Hz

400 A **Mains Cord**



69321520

## Spare Parts



2027590  
**Mains Fuse,**  
T 2.0 AL (SB),  
package of 10

# Suction and Irrigation Systems

VETPUMP® II

Components / Spare parts

Catalog page

- 69321601 VETPUMP® II**
- 69321620 **VETPUMP® II**, power supply 100/120, 230/240 VAC, 50/60 Hz
- 400 A **Mains Cord**
- 400 B **Mains Cord**, US-version
- 20300052 **Suction Bottle**, 1.5 l, sterilizable
- 20300034 **Bottle Cap**, for use with suction bottle 1.5 l, sterilizable
- 69311601 **Irrigation Bottle**, 1 l, sterilizable
- 69311610 **Bottle Cap**, for use with Irrigation Bottle 69311601
- 69300233 **Bottle Holder**, for suction bottle
- 69300133 **Bottle Holder**, for irrigation bottle
- 031124-10\* **Filter**, to block fluid during suction, for single use, unsterile, package of 10



69321620 C

U 16

## Spare Parts



1059400  
**Mains Fuse**,  
THT/T0, 800A/5 x 20 mm

Components / Spare Parts

Catalog page

- 20300042 Silicone Tubing Set**, for suction



U 16

## Spare Parts



27500  
**LUER-Lock Tube Connector**, male,  
tube diameter 9 mm



20300180  
**Tubing Connector Set**,  
for Bottle Cap 20300034



20300181  
**LUER-Lock Tube Connector**, female,  
tube diameter 8 mm

# Suction and Irrigation Systems

VETPUMP® II

Components / Spare Parts

Catalog page

20300034 **Bottle Cap**


U 16



## Spare Parts



20300580  
**Overflow Case**, gray,  
for Bottle Cap 20300034



031124-10\*  
**Filter**, to block fluid during  
suction, for single use,  
unsterile, package of 10



2911590  
**Plunger Ball**, green,  
to overflow case,  
2 pcs required

Components / Spare Parts

Catalog page

69311610 **Bottle Cap**

U 16



## Spare Parts

2386390  
**Sealing Ring**



mtp medical technical promotion gmbh,  
Take-Off GewerbePark 46, D-78579 Neuhausen ob Eck, Germany

# Motor Systems

UNIDRIVE® S III ARTHRO SCB

Components / Spare parts

Catalog page

- 28723001-1 UNIDRIVE® S III ARTHRO SCB
- 28723020-1 UNIDRIVE® S III ARTHRO SCB
- 20090170 SCB Connecting Cable, length 100 cm
- 400 A Mains Cord
- 96156013D Manual

U 21



28723020-1

## Spare Parts



1059700  
**Mains Fuse,**  
100 – 240 VAC,  
T 2.5 AL (SB),  
package of 10



1222890  
**Fuse,**  
100 – 127 VAC, T 5.0 A,  
package of 10

# High Frequency Surgery Units

AUTOCON® II 400 SCB

Components / Spare parts

Catalog page

**2053520x-12x AUTOCON® II 400 SCB**

U 40

2053522x-12x **AUTOCON® II 400 SCB,**  
power supply 220 – 240 VAC, 50/60 Hz  
400 A **Mains Cord**  
20090170 **SCB Connecting Cable,** length 100 cm



**2053520xU12x AUTOCON® II 400 SCB**

2053522xU12x **AUTOCON® II 400 SCB,**  
power supply 100 – 120 VAC, 50/60 Hz  
400 A **Mains Cord**  
20090170 **SCB Connecting Cable,** length 100 cm

2053522x-12x  
2053522xU12x

**Spare part for use at 230 V:**



2027690  
**Mains Fuse,**  
T 4.0 AL (SB),  
package of 10

**Spare part for use at 110 V:**



2028090  
**Mains Fuse,**  
T 8.0 AL (SB),  
package of 10

# High Frequency Surgery Units

AUTOCON® II 200, AUTOCON® II 80

Components / Spare parts

Catalog page

**20532201**     **AUTOCON® II 200**  
20532220     **AUTOCON® II 200**,  
power supply 200 – 240 VAC, 50/60 Hz  
400 A     **Mains Cord**, length 300 cm



U 42

**20532201-010**     **AUTOCON® II 200**  
20532220-010     **AUTOCON® II 200**, with GASTRO-Cut,  
power supply 200 – 240 VAC, 50/60 Hz  
400 A     **Mains Cord**, length 300 cm

20532220  
20532220-010  
20532220C  
20532220C010

**20532201 C**     **AUTOCON® II 200**  
20532220C     **AUTOCON® II 200**,  
power supply 100 – 120 VAC, 50/60 Hz  
400 B     **Mains Cord**, US version, length 200 cm

**20532201C010**     **AUTOCON® II 200**  
20532220C010     **AUTOCON® II 200**, with GASTRO-Cut,  
power supply 100 – 120 VAC, 50/60 Hz  
400 B     **Mains Cord**, US version, length 200 cm

**Spare part for use at 230 V:**



2027690  
**Mains Fuse**,  
T 4.0 AL (SB),  
package of 10

**Spare part for use at 110 V:**



2028090  
**Mains Fuse**,  
T 8.0 AL (SB),  
package of 10

Components / Spare parts

Catalog page

**20530801**     **AUTOCON® II 80**  
20530820     **AUTOCON® II 80**,  
power supply 100 – 240 VAC, 50/60 Hz  
400 A     **Mains Cord**, length 300 cm



U 44

**Spare part**



2027690  
**Mains Fuse**,  
T 4.0 AL (SB),  
package of 10

20530820

# Surgery Electrodes Set

Components / Spare parts

Catalog page

**20530008 Surgery Electrodes Set**

U 47

- 20530031 **Container with Lid and Sterilizing Insert,**  
for 16 electrodes with diameter 4 mm
- 26520031 **Wire Snare, 5 mm**
- 26520032 **Wire Snare, 10 mm**
- 26520033 **Ribbon Snare, 10 mm**
- 26520034 **KIRSCHNER Spatula Electrode, straight**
- 26520035 **MAGENAU Knife Electrode, angled**
- 26520036 **MAGENAU Knife Electrode, lancet-shaped**
- 26520037 **Ball Electrode, 2 mm**
- 26520038 **Ball Electrode, 4 mm**
- 26520039 **Ball Electrode, 6 mm**
- 26520040 **Needle Electrode**
- 26520041 **Flat Electrode, 8 x 10 mm**
- 26520042 **Flat Electrode, 10 x 15 mm**

