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Please call toll-free:

1-855-822-6999 M-F 9 a.m.-5 p.m. CST

E-mail: service@healthcare-manager.com

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
Easy Healthcare Corporation

360 Shore Dr., Burr Ridge, IL USA 60527

Made in China

Fabriqué en Chine

User Manual / Mode d'emploi

Contents

Easy@Home EHE080 TENS Overview	01
Intended Use	01
Features & Benefits	01
What's Included	02
Device Components	02
Let's Get Started with EHE080 (Setup)	03
Instruction for Use	03
Position Placement for Best Use	05
Practices for Best Use	06
Battery Charging	07
Cleaning and Maintenance	07
Frequently Asked Questions	08
Troubleshooting	09
Safety Warning	10
Technical Specifications	14
Product Programs	15
Warranty Information	16
Electromagnetic Compatibility and FCC Compliance Statement	16
Explanation of Symbols	20

Easy@Home EHE080 TENS Overview

Easy@Home EHE080 TENS and Muscle Stimulator is an effective and portable TENS and EMS combo unit with a rechargeable battery. EHE080 is a powerful edition to broaden Easy@Home TENS and EMS product line.

As one of the leading personal healthcare brands in pain management, Easy@Home brings years of experience and understanding from the professional area to the personal healthcare sector. We are proud to provide you the Easy@Home TENS and EMS products.

Intended Use

TENS (Transcutaneous Electrical Nerve Stimulation): To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household work activities.

EMS/PMS (Electrical/Powered Muscle Stimulation): To be used to stimulate healthy muscles in order to improve and facilitate muscle performance.

Features & Benefits

EHE080 is not only a TENS Unit (Transcutaneous Electrical Nerve Stimulation), but also EMS (Electrical Muscle Stimulation). It is also known as PMS (Powered Muscle Stimulation). This means the EHE080 unit is a combination device that provides both pain relief and muscle stimulation. As a compact and portable device, EHE080 can be used at home or on-the-go.


It is small and powerful at the same time! In fact, you can choose from 24 modes to personalize your pain relief session.

Those suffering with chronic pain related to muscles or nerves, tension, or even many kinds of arthritis, will find this unit quite helpful. The device will stimulate nerves to block pain sensations and release natural pain-relieving endorphins.



EHE080 will also stimulate muscles for those looking to tone, firm, or enhance their muscles with its EMS modes available on the device. This portable unit is extremely easy to take with you and use just about anywhere, such as the gym for post workout recovery or at work for pain relief.

What's Included

- 1 x TENS control unit
- 10 x Electrode pads
- 2 x Lead wires
- 2 x Pads holders
- 1 x Storage Case
- 1 x USB cable
- 1 x Manual

Device Components

1. ON/OFF switch
2. USB charging port
3. CH1 channel output
4. CH2 channel output
5. Mode Type Selector (TENS/EMS/MASSAGE function switch)
6. Right navigation (mode increment)
7. Left navigation (mode decrement)
8. Timer adjuster/Pause (short press to adjust time, long press for 2 seconds to pause/restart output)
9. CH1 channel intensity increase/decrease "+""-"
10. CH2 channel intensity increase/decrease "+""-"
11. CH1/CH2 switch/lock (short press to switch channels, long press for 2 seconds to lock or unlock)
12. CH1 mode display
13. CH1 intensity and frequency display
14. CH2 mode display
15. CH2 intensity and frequency display
16. Mode Type (TENS/EMS/MASSAGE) display
17. Pause icon
18. Battery icon
19. Timer icon
20. Lock icon

Let's Get Started with EHE080 (Setup)

The following steps are used to prepare the device for proper operation, and the details for each step are listed in the category sections below.

Step 1: Make sure the TENS and Muscle Stimulator is fully charged.

Step 2: Connect/Snap two electrode pads on each lead.

Step 3: Connect the other end of the leads to the lead channel output on the button of the device.

Step 4: Wash skin that requires treatment thoroughly and then dry it before applying electrode pads.

Step 5: Apply the whole surface of the electrode pads firmly on the cleaned skin, making sure to remove the plastic backing from the electrodes first.

NOTE:

1. You must use two electrode pads per channel but you can choose one or both channels; two channels allow you to receive stimulation at two different areas at the same time.

2. Make sure to keep at least 1.9cm (3/4") gap between each pad of the pair. Never have the two pads touching each other during the massage session in operation, or you can risk hurting your skin.

3. Always turn off the power before removing or re-positioning the electrode pads.

Instruction for Use

The following steps are a guide to use the device after setup, and the details for each step are listed in the table below.

Step 1: Slide "ON/OFF" switch to turn on the power.

Step 2: Select your desired output channel by short pressing the rounded button in the middle.

Step 3: Select the Mode type you need (TENS, EMS, MASSAGE) by pressing the "M" button.

Step 4: Select one of the stimulation modes by pressing "<" ">".

Step 5: Choose the stimulation time by short pressing "T" button.

Step 6: Adjust the stimulation intensity for each channel by pressing "+ " "-".

Enjoy your massage therapy!

● "ON/OFF" button to turn on the power

Slide the ON/OFF button at the top left of the device to turn on the unit.

● Selecting output channel

The default channel is CH1. You may switch the channels by short pressing the rounded button in the middle. Each channel can be separately controlled with different massage modes and intensity settings.

● Selecting stimulation modes

Press "M" to choose the desired mode type between TENS, EMS and MASSAGE.

Press ">" or "<" to select your preferred mode. For safety reasons, the default intensity is set to 0 each time after you make a new mode selection. The mode of each channel (CH1 or CH2) can be set separately to the desirable one.

● Selecting stimulation time and pause/restart

Short press the " T/|| " button down. Repeat until you have chosen the time you prefer for your session. You can choose from 10, 20, 30, 40, 50, 60 minutes. The timer indication is on the upper right part of the display. 20 minutes is the default time. The device will shut off as a safety mechanism once the timer is complete.

Long press the " T/|| " button for 2 seconds to pause/restart.

NOTE:

1. Pause for 5 minutes, timer will not change, but device will be off automatically after 5 minutes. Slide the ON/OFF button to turn on the unit again.

2. Device will be off automatically after 5 minutes without any operation.

● Adjust stimulation intensity

Press the "+" button to increase the stimulation intensity, and press the "-" button to decrease the intensity. The intensity of the two channels can be set separately to the desirable level.

NOTE:

1. With the increase of intensity, you may experience sensations like tingling, vibration, etc. Gradually increase the intensity, and stop increasing when a comfortable level is reached.

2. Make sure you are increasing the correct channel. You can see the change of the intensity represented by the change in bars on the display.

- **Lock and unlock**

Press and hold the rounded button in the middle for 2 seconds to lock the entire setting display and avoid accidental changes. Press and hold for 2 seconds again to release the lock.

- **Enjoy the stimulation massage therapy**

After the above mode, time, and intensity are set to your preference, enjoy the stimulation provided by the device. Slide the“ON/OFF” to turn off the power to stop use, or the device will automatically shut off when the timer has completed.

Position Placement for Best Use

You can place TENS Pads on many areas of the body to relieve muscle and nerve pain and to stimulate healthy muscles. Some areas are great for the TENS Mode such as shoulders,back, waist, arms,legs, feet and hips. Muscle stimulation can be applied to many muscle areas on the body, such as abdomen, buttocks,arms and legs.

LOW BACK

Sitting for too long, having weak hip flexors, or exercising can aggravate pain in our hips and waist. To alleviate tension or pain, attach the electrode pads around your waist and lower back. Space the pads from 5.08cm (2") to 15.24cm (6") apart for the best results.

WAIST

TENS is perfect for relief from pain in the waist area. Whether weak hip flexors or a sore side muscle, the aching muscles can be alleviated. The TENS can soothe the nerves and help with hip pain from issues with tendons and muscles by providing gentle movement.

BUTTOCKS

The EMS mode can be used for strengthening and tightening muscle on the buttocks or hip areas.TENS mode can be used to alleviate tension or pain on the hips and buttocks too. Attach the one pair of pads at the lower back and another pair on the hip area. Space the pads at least 5.08cm (2") apart.

ABDOMEN

The EMS feature can help stimulate the gluteus muscles for better athletic performance. With the contraction of the muscles this device provides, the muscles in this area can be strengthened. This area is sensitive so it does not need much intensity to feel the effects.

ARMS

Use EMS mode to strengthen arm muscle like biceps by attaching one pair at the upper end of the muscle next to the shoulder and another pair on the lower end of the muscle. Gradually increase the intensity until you get a maximum muscle contraction, without being painful.

THIGHS/LEGS

Bicyclists know the incredible stress and strain that your thigh muscles can experience. Since these large muscles are so crucial to our mobility, it is important to take care of them and help them relax. With TENS mode, placing the electrode pads on the top and back side of your thighs is an uncomplicated way to relax and de-stress this area.

Use EMS mode to strengthen leg muscles like hamstrings and quadriceps by attaching one pair at the upper end of the muscle and another pair on the lower end of the muscle. Gradually increase the intensity until you get a maximum muscle contraction, without being painful.

Practices for Best Use

- **Electrode pads**

The electrode pads can be used until they lose the stickiness and/or conductivity. Their lifetime greatly depends on the body skin condition and protection. Purchase new electrode pads as needed.

NOTE:

1. Only use the electrode pads of the same brand as the controller device.
2. Do not use the same electrode pads on different people. Each user must have his/her own pads.

- **Low intensity to start**

With the increase of intensity, you may experience sensations like tingling, vibration, etc. Therefore, gradually increase the intensity, and stop increasing when a comfortable level is reached.

- **Cleanliness**

Keep the device clean by wiping with a moistened cloth from time to time. This prevents growth of germs and dirt and can extend the life of the pads.

- **Storing device**

When not in use, store the device and accessory in a cool place, out of direct sunlight.

- **Areas to avoid**

Do not use your TENS Unit on your head, neck, chest near heart, or genitals. Do not apply stimulation across the chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the heart.

For instance, DO NOT use one of a TENS pad set on one arm and the other of the same set on the other arm, because the electrical current may go through your heart from one arm to another arm.

Battery Charging

The TENS and Muscle Stimulator (TENS Unit) comes with a built-in rechargeable Lithium battery. Make sure the unit is fully charged before using. The unit can be charged with the adapter and cable through a standard wall outlet or through the USB port of your computer with the cable only. The device takes approximately 2-3hours for a full charge through the wall outlet.

Charge the unit when you see the battery indicator has only one bar left, or the screen is dim, the sensation is weak, or the display does not light up when turning the device on.

Cleaning and Maintenance

Please use wipes to clean the device first, and then use a dry cloth to wipe it again. The electrode pads that come with the device are disposable and should be replaced when they lose their adhesiveness. Contact the seller for replacements. Do not let the sticky side of the pad touch anything, including the greasy finger tips.

Frequently Asked Questions

1. How does my TENS unit relieve pain or stimulate muscles? How does it work?

EHE080 relieves pain by blocking pain signals and stimulating natural painkillers known as endorphins.

This unit can also stimulate muscles, which helps with muscle performance.

2. How many modes and intensity levels does the EHE080 have?

This unit has 24 massage modes and 20 intensity levels you can choose from, so you can get benefits at your own level of comfort.

3. Is it safe for anybody to use my TENS anywhere on the body?

The EHE080 unit is an effective and drug-free device for pain relief and muscle stimulation when used properly. But, it is not suitable for some body areas of some people. Please read the Safety Warning (on page 10) before using this unit. If you have any medical condition and are unsure about using the device safely, consult your physician first.

4. What is the difference between EMS and TENS? Does this unit provide both?

The EHE080 Unit is a TENS and EMS Combo Unit. EMS stands for Electrical Muscle Stimulation and TENS stands for Transcutaneous Electrical Nerve Stimulation.

TENS modes on the device are specifically for pain and soreness relief.

The EMS/PMS modes are for muscle stimulation and better blood circulation in the lower body for healthy muscles.

Yes, this device provides both.

5. What can I do to elongate the life of my electrode pad?

Cleanse the area of skin you will be placing the electrode pads on with soap and water, or with a damp cloth. Make sure the area is dry before applying the electrodes. When storing the electrode pads, we recommend placing them back onto the plastic backing. Electrode pads that become wet from sweat or water, and do not stick anymore need to be replaced.

You can also cleanse the pads with a damp cloth. If they are too dirty or your skin feels numb after the pads are cleaned, replace the pads.

6. What type of conditions does TENS/EMS really benefit?

TENS is useful in relieving many kinds of pain and soreness. TENS is great for lower back pain, knee pain, sciatica, fibromyalgia. These are just some examples.

EMS modes that are great for improving muscle performance or increasing blood circulation in the lower extremities from a sedentary lifestyle if your muscles need toning and firmness.

Troubleshooting

Problem	Possible Causes	Suggested solution
The intensity level can not be increased	Pads are not firmly attached on the skin. Or there is a loose connection. Or the user's skin is too dry.	Make sure both pads are firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Keep the skin moisturized.
Unit is on, but no sensation is felt in the pads.	Both pads are not firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased.	Make sure both pads are firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.
One pad feels stronger than the other.	This is normal. Different areas of your body will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.
During the massage, the skin feels a painful burning sensation or the stimulation becomes weakened.	Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry.	Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application.
Adhesive gel pads do not stick to skin even after cleaning and moistening the gel pad.	Adhesive gel pads need to be replaced.	Replace the set of pads. Pads can be purchased from local distributor.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Does not provide adequate pain relief.	Not using it long enough.	Use your massager 20-30 minutes at a time, 3-6 times per day.

Safety Warning

● Contraindications

Do not use this device on users who have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device, because this may cause electric shock, burns, electrical interference, or death.

Do not use this device if your pain syndromes are undiagnosed.

● Warnings

Do not apply stimulation on your neck because this could cause severe muscles spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.

Do not apply stimulation across your chest, because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.

For instance, do not use one of a TENS pad set on one arm and the other of the same set on the other arm, because the electrical current may go through your heart from one arm to another arm.

Do not apply stimulation on your head, chest, near your heart, or genitals areas.

Do not apply stimulation over or in proximity to cancerous lesions.

Do not apply stimulation when the patient is in the bath or shower.

If you have one of the following conditions, please consult with your physician before purchasing or using this device:

Acute disease, malignant tumor, infective disease, pregnant, heart disease, high fever, abnormal blood pressure, lack of skin sensation or an abnormal skin condition, any condition requiring the active supervision of a physician.

● Precautions

Do not use this device while driving.

Do not use this device while sleeping.

Do not use this device in high humidity areas such as a bathroom.

Keep the device away from wetness, high temperatures and directly-sunlit places.

Keep this device out of the reach of children.

Stop using this device at once if you feel pain, discomfort, dizziness or nausea, and consult your physician.

Do not attempt to move the electrode pads while the device is operating.

Do not use the device around the heart or on the head, mouth, genitals, or blemished skin areas.

Do not apply stimulation of this device in the following conditions:

- (1) Across the chest. The introduction of electrical current into the chest may cause rhythm disturbances to the heart, which could be lethal;
- (2) Over painful areas. Please consult with your physician before using this device if you have painful areas;
- (3) Over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins). Apply stimulation only to normal, intact, clean, healthy skin;
- (4) In the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms). The electronic stimulator may not operate properly when the electrical stimulation device is in use;
- (5) While operating machinery, or during any activity in which electrical stimulation can put you at risk of injury;
- (6) On children.

Be aware of the following:

- (1) Consult with your physician before using this device. The stimulation with the device may:
 - i. Cause lethal rhythm disturbances to the heart in susceptible individuals;
 - ii. Disrupt the healing process after a recent surgical procedure;
- (2) The device is not effective for pain in your head area, including for headaches;
- (3) The device is not a substitute for pain medications and other pain management therapies;
- (4) The device has no curative value;
- (5) The device is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism;
- (6) The long-term effects of electrical stimulation are unknown;
- (7) The user may experience skin irritation, burns or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel);
- (8) If the user has suspected or diagnosed epilepsy, the user should follow precautions recommended by his or her physician;
- (9) Use caution if the user has a tendency to bleed internally, such as following an injury or fracture;
- (10) Use caution if stimulation is applied over the menstruating uterus;
- (11) Use caution if stimulation is applied over areas of skin that lack normal sensation;
- (12) Stop using the device if the device does not provide pain relief;
- (13) Use this device only with the leads, electrodes and accessories that the manufacturer recommends;

- (14) Do not share the use of the electrode pads with others;
- (15) Do not use the device while it is charging;
- (16) The device contains a lithium battery. If overheating of the device occurs during charging, stop charging or operation immediately and report to the seller;
- (17) Dispose of the battery-containing device according to the local, state, or federal laws.

The long-term effects of electrical stimulation are unknown.

Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrode pads should not be placed on opposite sides of the head.

The safety of electrical stimulation during pregnancy has not been established.

Some users may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel).

Users with suspected or diagnosed heart disease should follow precautions recommended by their physicians.

Users with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.

Use caution if stimulation is applied over the menstruating or pregnant uterus.

• Adverse reactions

You may experience skin irritation and burns beneath the stimulation electrodes applied to the skin;

You may experience a headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.

You should stop using the device and should consult with your physician if you experience adverse reactions from the device.

Environmental condition for normal working,transport and storage

- Normal working ambient temperature: 5~40°C(41-104°F)
- Normal working ambient humidity: 15%~90% RH
- Atmospheric pressure: 70~106kPa
- Storage and transport ambient temperature:-25~70°C(-13-158°F)
- Storage and transport ambient humidity: 0%~90% RH

Safety and essential performance test standards:

- Medical Devices Directive 93/42/EEC
- IEC 60601-1:2005 +A1:2012 /EN 60601-1:2006 +A1:2013 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- ANSI/AAMIES60601-1:2005 +C1:2009 +A2:2010+A1:2012 Medical electrical equipment — Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2:2014/EN 60601-1-2:2015 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
- IEC 60601-2-10:2012+A1:2016/EN 60601-2-10:2005+A1:2016 Medical electrical equipment - Part 2-10Particular requirements for the safety and essential performance of nerve and muscle stimulators
- IEC60601-1-11:2015/EN 60601-1-1:2015 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medicalel ectrical systems used in the home healthcare environment.
- ISO 15223-1:2021/ENISO 15223-1:2021 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements EN 1041 Information supplied by the manufacturer with medical devices
- IEC60601-1-6:2010+A1:2013/EN 60601-1-6:2010+A1:2015 Medical electrical equipment -Part1-6: General requirements for basic safety and essential performance- Collateral standard: Usability
- IEC62304:2006+A1:2015/EN 62304:2006+A:2015 Medical device software - Software life-cycle processes
- IEC 62366:2007 +A1:2014/EN 62366:2008 +A1:2015 Medical devices - Application of usability engineering to medical devices
- ISO 10993-1:2009/ENISO 10993-1:2009 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

Technical Specifications

Mode/type	EHE080	Weight	65g
Power supply	Powered by internal 3.7V li-ion battery	Default working time	20 minutes
Wave form and wave shape	Biphasic rectangular wave pulse	Degree of protection against electric shock	Type BF applied part
Pulse duration	100/150/200/240µs	Type of protection against electric shock	Internally powered equipment
Pulse frequency	1-132Hz	Grade of waterproof	IP 22
Output Voltage	RMS Output: 33mA	Product life	3 Years
Treatment time	10, 20, 30,40,50, 60 minutes	Lifetime for electrode	Storage for 2 years (no use)
Modes	24 auto modes	Software version	A0
Typical operation time of Battery	Approximately 18 hours (Note: This time is based on laboratory data with both channels operating at intensity level of 5 under mode 1.)	The time required for my equipment to warm from the minimum storage temperature between uses until it is ready for intended use	30 minutes
Behavior of my equipment while the rechargeable internal electrical power source is charging:	The battery indicator on the display will indicate the battery charging status.	The time required for my equipment to cool from the maximum storage temperature between uses until it is ready for ntended use	15 minutes
Typical service life of battery	300 times of recharging	Adapter for charging	Input:100-240V 50/60Hz/0.25A Max output: DC5V 550mA
NOTE: Not intended to be sterilized.			
Not for use in an OXYGEN-RICH ENVIRONMENT			

Class II, not applied part, not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide continuous operation

Product Programs

TENS			
Mode	Frequency(Hz)	Pulse Width(μs)	Time(min)
1	46	100	10,20,30,40,50,60
2	83	100	10,20,30,40,50,60
3	3	100	10,20,30,40,50,60
4	10-38	100	10,20,30,40,50,60
5	1-12	100	10,20,30,40,50,60
6	65	100	10,20,30,40,50,60
7	54	100	10,20,30,40,50,60
8	45	100	10,20,30,40,50,60
EMS			
Mode	Frequency(Hz)	Pulse Width(μs)	Time(min)
9	92	100	10,20,30,40,50,60
10	88	100	10,20,30,40,50,60
11	54	100	10,20,30,40,50,60
12	48	100	10,20,30,40,50,60
13	45	100	10,20,30,40,50,60
14	48	100	10,20,30,40,50,60
15	48	100	10,20,30,40,50,60
16	48	100	10,20,30,40,50,60
MESSAGE (TENS and EMS combination mode)			
Mode	Frequency(Hz)	Pulse Width(μs)	Time(min)
17	12-132	100	10,20,30,40,50,60
18	36	100	10,20,30,40,50,60
	6	100	10,20,30,40,50,60
19	8	100	10,20,30,40,50,60
	4	100	10,20,30,40,50,60
20	2	100	10,20,30,40,50,60
	58	100	10,20,30,40,50,60
21	65	100	10,20,30,40,50,60
	2	100	10,20,30,40,50,60
22	3-25	100	10,20,30,40,50,60
23	65	100	10,20,30,40,50,60
	3	100	10,20,30,40,50,60
24	2-22	100	10,20,30,40,50,60
	54	100	10,20,30,40,50,60

Warranty Information

This device carries a limited warranty of one year from the date of delivery. The warranty applies to the device only, and the accessories are not covered by this warranty. During the warranty period, defective items will be repaired or replaced at no charge. Any evidence of misuse, abuse, alternations, or externally caused damage may have this warranty invalid. For more information, please contact the manufacturer.

Electromagnetic Compatibility and FCC Compliance Statement

(1) This product needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile radio frequency (RF) communications equipment.

(2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.

(3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!

(4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacture's declaration – electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used such an environment.		
EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT-GUIDANCE
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacture's declaration -electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below, The customer or the user of device should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETICENVIRONMENT -GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines	± 0.5kV, ± 1 kV line(s) to lines	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315° 0% UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315° 0% UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	Main power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacture's declaration -electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified on following page. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETICENVIRONMENT -GUIDANCE
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur bands between 0.15 MHz and 80 MHz	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur bands between 0.15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.167 \sqrt{P}$ $d = 1.167 \sqrt{P}$ 80MHz to 800 MHz $d = 2.333 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol 
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m	
NOTE 1: At 80 MHz and 800MHZ, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations.Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 6 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer of the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter(W)	Separation distance according to frequency of transmitter(m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d=1.167 \sqrt{P}$	$d=1.167 \sqrt{P}$	$d=2.333 \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333
For transmitters rated at a maximum output power not listed above, the recommended separation distance (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

The subject device has been tested and found to comply with the limits for a Class B digital device pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

The product generates, uses, and can radiate radio frequency energy and, if not installed and used accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that the interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on or off, the user is encouraged to try to correct the interference by one or more of the following measures:

- a) Reorient or relocate the receiving antenna;
- b) Increase the separation between the product and the receiver;
- c) Consult the dealer or an experienced radio/TV technician for help.
- d) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Changes or modifications to this product not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Explanation of Symbols

Information essential for proper use shall be indicated by using the corresponding symbols. The following symbols may be seen on the device and labeling.

Symbol	Title	Symbol	Title
	Batch code		Manufacturing date
	Ingress protection rating		Type BF applied part
	Fragile, handle with care		Keep dry
	Refer to instruction manual		Recyclable
	Atmospheric pressure limitation of 70kPa~106kPa		Temperature limit of -13°F~158°F (-25°C~70°C)
	Humidity limitation of 0% ~ 90%		Radio frequency transmitters

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