

# visual **Baby Nasal Aspirator**



User Manual Model No.:ENA101

Please read the user manual carefully before using.

### **Contents**

### Overview

The Easy@Home visual baby nasal aspirator, with a camera, helps parents to safely and naturally clean baby's nasal cavity for unobstructed breathing. The visual nasal aspirator has 3-level suction intensity to clean dried mucus of different volume and excess phlegm. Its soft suction prevents from disturbing infants sleep. The camera of nasal aspirator allows parents to clearly observe baby's nasal cavity. The Easy@Home visual baby nasal aspirator is an effective childcare kit.

## **Product Components**


### Contents:

- 1 \* Visual Baby Nasal Aspirator
- 2 \* Silicone Suction Nozzles
- 1 \* Tweezer
- 1 \* Charging Cable (adaptor not included)

1 \* Brush

- 2 \* Replaceable Sponges
- 1 \* Manual





## Direction of Use

### Step1: Choose a suitable nozzle.

Please choose a suitable size nozzle for your child. If child's nasal cavity is thin or for deeper suction, then choose the thin nozzle. Choose the standard nozzle for regular suction or a shorter nose cavity.

### Step 2: Download the App.

The aspirator needs the help of Easy Nasal Care to connect with a smartphone. Search "Easy Nasal Care" in Apple Store or Google Play Store to download and install the app, or simply scan the QR code provided below to download it.





QR Code to download app

Easy Nasal Care

### Step 3: Connect the smartphone with the baby nasal aspirator.

Press the power button and hold for a while (about 3 seconds) until the indicator light shines which means the device is turning on. Enter the app and connect your smartphone with the baby nasal aspirator's WiFi according to the tips on the app. The camera fill light on the aspirator will shine if the baby nasal aspirator connects with the smartphone successfully.

If you can't connect the WiFi successfully, please check whether several

phones are open the WiFi setting at the same time. The aspirator will primarily default to connecting the smartphone previously connected with it before. It is best to keep other phones WiFi function off if you want to connect the aspirator with a phone never connected before.

### Step4: Observe nasal cavity through camera.

You can easily observe children's nasal cavities by camera. The app provides video and photo modes. Also, there is an enlargement mode for observing small, dried mucus clearly. Photos and videos will be stored in album.









Step 5: Clean the nasal secretions.

There are three-level intensities of suction: weak, medium, and strong, Click the power button to choose a suitable suction intensity and then aim into the nasal cavity to suck the secretions. The suction nozzle should be as close as possible to the secretions, and it will suck out dried mucus more easily.





DO NOT use it to clean baby's other physical body parts, such as ear cavity, eves or oral cavity.

Please clean the nozzle and storage bin to avoid bacterial infection after using every time.

### Charging Reminder

Low battery: Three indicator lights shining at the same time means the device in low battery and it will turn off automatically after 90 seconds.

Charging: Three indicator lights shining circularly means the device currently charging.

Full charged: Three indicator lights steadily shining brightly at the same time means the device has a full battery.

### Automatic Shut-down

The aspirator will shut-off automatically if no operation for 10 minutes after turning on.


01 02

## **Cleaning and Maintenance**

1.Remove the nozzle and storage bin. Gently remove the nozzle with your fingers and then unscrew the storage bin by turning counterclockwise.



2. Cleaning methods shown as the following picture.



## **Storage Condition**

Operation Temperature: 50°F~104°F (10°C~40°C)

Operation Humidity: Under 85% R.H.

Storage Temperature: 14°F~122°F  $(-10^{\circ}\text{C}\sim50^{\circ}\text{C})$ 

Storage Humidity: Under 85% (Non-condensing)

## Warranty

This device carries a limited warranty of one year from the date of delivery. The warranty applies to the device only, and the accessories are not covered by this warranty.

During the warranty period, defective items will be repaired or replaced at no charge. Any evidence of misuse, abuse, alternations, or externally caused

damage may have this warranty invalid.

For more information, please contact the manufacturer.

## **Service Life**

The service life of the device is 3 years.

## Warning

- 1. Children should use this product under the assistance of adults.
- 2. Child with nose/eyes/ears allergies, rhinitis, heart disorders, sensitive skin and acute diseases should consult a doctor before using.
- 3. This product is only used for cleaning nasal secretions, please do not use it for other purposes
- 4. If feeling unwell while using this product, please stop using immediately, and check whether there is improper operation, or the product isn't suitable for your baby.
- 5.Remove the battery if discarding the product. Remove the battery: ①Cut off the power supply. ②Open the outer casing. ③Remove the battery. Waste batteries should be disposed of to professional institutions.
- Avoid children playing with device as its accessories may cause choking hazard.
- 7. Do not microwave the product or any part of it.
- 8. Do not probe the aspirator too far into the nostril.
- 9. Do not use while charging.
- 10. Make sure to turn off the power before taking apart the suction nozzle and storage bin.

## **Troubleshooting**

Problem description	Problem analysis	Solution	
Fail to turn on	The battery is too low, or the circuit is faulty.	Try to turn on the device on again after charging it.	
Suction intensity is too weak	Improper suction nozzle installation.	Check suction nozzle and storage bin installation.	
Camera faulty	WiFi connection failed.	Try to connect WiFi again, Make sure another nearby phone is not connected.	

### **Product Parameters**

Product Name	Visual Baby Nasal Aspirator	Vacuum Intensity	≥ 55 kPa
Model Number	ENA101	Flow	≥ 2.2 LPM
Input Voltage	5V 1000mA	Battery Capacity	800 mAh
Rated Power	3.5 W	Working Time	About 60 mins
Working Current	650 mA	Weight	150g
Camera Pixel	2M	Product Size	147 x 48 x 69.4 mm
Executive Standard	GB 4706.1-2005; GB4706.15-2008		

### Q&A

Q1: What should I do if the nasal secretion is dry and stubborn?

A1: Don't compress baby's nose because the dry secretions may hurt his or her nasal cavity. It is recommended to immerse Q-tips for children in clean water until it gets wet and wipe softly in nasal cavity to soften the secretions. Otherwise, parents can use the nasal aspirator after the baby bath as the water vapor from the bath will soften secretions.

Q2: What should I pay attention to when using a nasal aspirator?
A2:

- 1) Avoid probing the suction nozzle too deep into the nostril when using to avoid accidentally hurting the nasal cavity.
- 2) The suction nozzle should be as close as possible to the secretions, and it will suck out dried mucus more easily.
- If feeling unwell during the use of this product, please stop using immediately, and check whether there is improper operation, or the product isn't suitable for your baby.

Q3: What should I do if the baby still has nasal congestion after suction?

A3: If parents confirm the nostrils are free of secretions but baby still has nasal congestion, parents can apply a warm towel on the baby's nose. This will relieve the symptoms of nasal congestion temporarily. Also try to adjust baby's position: lie on left side when the right nose is congested, lie on right side when the left nose is congested. This will have as temporary alleviation. If any concerns, consult a doctor immediately.

### **FCC Information**

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## **Explanation of symbols**

	Date of Manufacture	\$\frac{1}{2}	Recyclable
LOT	Batch Code	(( <u>(</u> ))	Radio frequency transmitters
[]i	Consult instructions for use		

## Sealthcare-Manager.com

Questions or Comments?

Please call toll-free: 1-855-822-6999 M-F 9 a.m.-5 p.m. CST

E-mail: service@healthcare-manager.com
Easy Healthcare Corporation
360 Shore Dr., Burr Ridge, IL USA 60527

Made in China

05 06 07 08