

Baby Nasal Aspirator



User Manual
Model No.:ENA102

Please read the user manual carefully before using.

Contents	
Contents.....	02
Overview	02
Product Components	02
Direction of Use	03
Cleaning and Maintenance	04
Storage Condition	05
Warranty	05
Service Life	05
Warning	06
Product Parameters	06
Troubleshooting	07
Q&A	07
FCC Information	08
Explanation of symbols.....	08




Overview

The Easy@Home baby nasal aspirator helps parents to safely and naturally clean baby's nasal cavity for unobstructed breathing. The nasal aspirator has 3-level suction intensity to clean dried mucus of different volume and excess phlegm. Its soft suction prevents from disturbing infants sleep. The Easy@Home baby nasal aspirator is an effective childcare kit.

Product Components

Contents:

- 1 * Baby Nasal Aspirator
- 2 * Silicone Suction Nozzles
- 1 * Tweezer 1 * Brush
- 1 * Charging Cable 1 * Manual
- 2 * Replaceable Sponges




Direction of Use

Step1: Choose a suitable nozzle.

Please choose a suitable size nozzle for your child. If child's nasal cavity is thin or for deeper suction, then choose the thin nozzle. Choose the standard nozzle for regular suction or a shorter nose cavity.

Step 2: Clean the nasal secretions.

- (1) There are three-level intensities of suction, weak, medium and strong. Click the power button to choose a suitable suction intensity.
- (2) Aim into the nasal cavity to suck the secretions and press steadily the suction button when sucking the nasal secretions.
- (3) If you want to pause suction just stop pressing the suction button. The suction nozzle should be as close as possible to the secretions, and it will suck out dried mucus more easily.



Note:

DO NOT use it to clean baby's other physical body parts, such as ear cavity, eyes or oral cavity.

Please clean the nozzle and storage bin to avoid bacterial infection after using every time.

Night Light Mode

The nasal aspirator has a night light mode which shines with a colorful light to attract the baby's attention and comfort the baby. Turn on the device and quickly press the power button 4 times to get this mode.

Charging Reminder

Low battery: Three indicator lights shining at the same time means the device in low battery and it will turn off automatically after 90 seconds.

Charging: Three indicator lights shining circularly means the device currently charging.

Full charged: Three indicator lights steadily shining brightly at the same time means the device has a full battery.

Automatic Shut-down


The aspirator will shut-off automatically if no operation for 10 minutes after turning on.

Cleaning and Maintenance

1. Remove the nozzle and storage bin. Gently remove the nozzle with your fingers and then unscrew the storage bin by turning counterclockwise.



2. Cleaning methods shown as the following picture.



Storage Condition

Operation Temperature: 50°F~104°F (-10°C~40°C)
 Operation Humidity: Under 85% R.H.
 Storage Temperature: 14°F~122°F (-10°C~50°C)
 Storage Humidity: Under 85% (Non-condensing)

Warranty

This device carries a limited warranty of one year from the date of delivery. The warranty applies to the device only, and the accessories are not covered by this warranty. During the warranty period, defective items will be repaired or replaced at no charge. Any evidence of misuse, abuse, alterations, or externally caused damage may have this warranty invalid. For more information, please contact the manufacturer.

Service Life

The service life of the device is 3 years.

Warning

- Children should use this product under the assistance of adults.
- Child with nose/eyes/ears allergies, rhinitis, heart disorders, sensitive skin and acute diseases should consult a doctor before using.
- This product is only used for cleaning nasal secretions, please do not use it for other purposes.
- If feeling unwell while using this product, please stop using immediately, and check whether there is improper operation, or the product isn't suitable for your baby.
- Remove the battery if discarding the product. Remove the battery:
 - Cut off the power supply.
 - Open the outer casing.
 - Remove the battery. Waste batteries should be disposed of to professional institutions.
- Avoid children playing with device as its accessories may cause choking hazard.
- Do not microwave the product or any part of it.
- Do not probe the aspirator too far into the nostril.
- Do not use while charging.
- Make sure to turn off the power before taking apart the suction nozzle and storage bin.

Product Parameters

Product Name	Baby Nasal Aspirator	Vacuum Intensity	≥ 55 kPa
Model Number	ENA102	Flow	≥ 2.2 LPM
Input Voltage	5V ===1000mA	Battery Capacity	800 mAh
Rated Power	3.5 W	Working Time	About 60 mins
Working Current	650 mA	Weight	150g
Executive Standard	GB 4706.1-2005; GB4706.15-2008	Product Size	147 x 48 x 69.4 mm

Troubleshooting

Problem description	Problem analysis	Solution
Fail to turn on	The battery is too low, or the circuit is faulty.	Try to turn on the device on again after charging it.
Suction intensity is too weak	Improper suction nozzle installation.	Check suction nozzle and storage bin installation.

Q&A

Q1: What should I do if the nasal secretion is dry and stubborn?
 A1: Don't compress baby's nose because the dry secretions may hurt his or her nasal cavity. It is recommended to immerse Q-tips for children in clean water until it gets wet and wipe softly in nasal cavity to soften the secretions. Otherwise, parents can use the nasal aspirator after the baby bath as the water vapor from the bath will soften secretions.

Q2: What should I pay attention to when using a nasal aspirator?
 A2:

- Avoid probing the suction nozzle too deep into the nostril when using to avoid accidentally hurting the nasal cavity.
- The suction nozzle should be as close as possible to the secretions, and it will suck out dried mucus more easily.
- If feeling unwell during the use of this product, please stop using immediately, and check whether there is improper operation, or the product isn't suitable for your baby.

Q3: What should I do if the baby still has nasal congestion after suction?
 A3: If parents confirm the nostrils are free of secretions but baby still has nasal congestion, parents can apply a warm towel on the baby's nose. This will relieve the symptoms of nasal congestion temporarily. Also try to adjust baby's position: lie on left side when the right nose is congested, lie on right side when the left nose is congested. This will have as temporary alleviation. If any concerns, consult a doctor immediately.

FCC Information

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received including interference that may cause undesired operation.

Explanation of symbols

	Manufacturing date		Recyclable
	Batch Code		Radio frequency transmitters
	Consult instructions for use		

Healthcare-Manager.com

Questions or Comments?

Please call toll-free:

1-855-822-6999 M-F 9 a.m.-5 p.m. CST

E-mail: service@healthcare-manager.com

Easy Healthcare Corporation

360 Shore Dr., Burr Ridge, IL USA 60527

Made in China