TENS and Muscle Stimulator

User Manual
Model: AS1080-16AB

Manual (EN)	1-28
Bedienungsanleitung (DE)	29-56
Manuale Utente(IT)	57-84
Manuel d'utilisation(FR)	86-112
Manual del Usuario(ES)	113-140

Contents

Easy@Home AS1080-16AB TENS Overview	2
Device Components	3
Let's Get Started with AS1080-16AB (Setup)	4
Features & Benefits	6
Functions (Intended for use).....	7
Instructions for Use	7
Mode Reference Chart	9
Position Placement for Best Use	9
Practices for Best Use	12
Frequently Asked Questions	13
Cleaning and Maintenance	14
Troubleshooting	15
Warranty Information AS1080-16AB	15
Safety Warning.....	24
Technical Specifications	25

Easy@Home AS1080-16AB TENS Overview

Easy@Home AS1080-16AB Electronic Pulse Stimulator is an effective and portable TENS and EMS combo unit with a rechargeable battery. AS1080-16AB is a powerful edition to the broad Easy@Home TENS and EMS product line.

TENS (Transcutaneous Electrical Nerve Stimulation) Unit is an effective, drug-free, safe and easy to use device to relieve the pain and soreness from the comfort of your home or on-the-go.


EMS (Electrical Muscle Stimulation) Unit sends electronic pulses to muscles to help with muscle toning and firmness, endurance improvement, and to lessen recovery time between workouts.

AS1080 - 16AB device can also be used for improving local blood circulation in healthy muscles of lower extremities for those who travel, are sedentary for long periods of time, or need to enhance their blood circulation due to a health condition.

As one of the leading personal healthcare brands in pain management, Easy@Home brings years of experience and understanding from the professional area to the personal healthcare sector. We are proud to present you the Easy@Home TENS and EMS product line.

Device Components

- | | |
|-------------------------------------|--------------------------|
| 1.On/off switch | 9.Intensity decrease "-" |
| 2.USB port | 10.Lead channel output |
| 3.Battery indicator | 11.Lead wires |
| 4.LCD screen | 12.Electrode pads |
| 5.Mode selector | 13.USB cable |
| 6.Intensity increase "+" | 14.Pads holder |
| 7.Timer adjuster | |
| 8.A-B channel switch & lock/release | |

Operation of the device is as simple as using the following main buttons:

- The "+" and "-" buttons. Increase or decrease the intensity level for the selected channel.
- The "T" button. Adjust the timer: 10, 20, 30, 40, 50, 60 minutes options.
- The "M" button. Select one of the sixteen auto procedure modes. Pressing this button will move you through each

massage type until you find the one you prefer.

- The round button in the middle. Switch channels between A and B with a quick press of the button; lock/release with a long press of the button for 3 seconds.

Let's Get Started with AS1080-16AB(Setup)

• Setup

Unpack the box, take out the product and accessories, connect/snap the electrodes pads to the leads, plug in the other end of the lead wire to the controller and apply the electrode pads to the area of skin where you want to receive the stimulation. Turn on the power, set the mode and intensity level to your comfort. It's that easy!

The following steps are used to prepare the device for proper operation, and the details for each step are listed in the category sections below.

Step 1: Make sure the Electronic Pulse Stimulator is fully charged

Step 2: Snap 2 electrode pads on each lead

Step 3: Connect the other end of the leads to the lead channel output on the bottom of the device.

Step 4: Apply the electrode pads on the body area that requires treatment, making sure to remove the plastic backing from the electrodes first.

• Battery charging

The Electronic Pulse Stimulator (TENS Unit) comes with a built-in rechargeable Lithium battery. Make sure the unit is fully charged before using. The unit can be charged with the adapter and cable through a standard wall outlet or through the USB port of your computer with the cable only. The device takes approximately 1.5 hours for a full charge through the wall outlet.

Charge the unit when you see the battery indicator has

only one bar left, or the screen is dim, the sensation is weak, or the display does not light up when turning the device on.

- **Electrode pads installment**

Snap the enclosed electrode pads onto the leads. Then plug the other ends of the leads into the lead channel output on the bottom of the controller.

Make sure to connect two pads to each lead that you are using (channel A and channel B). You must use two electrode pads per channel but you can choose one or both channels; two channels allow you to receive stimulation at two different areas at the same time.

- **Electrode pads placement on body areas**

After you have snapped the pads into the leads and plugged them into the controller, peel the transparent protective plastic sheet off the electrode pads. Place the pads on the body area. Make sure the area you choose is clean and dry. Make sure to keep at least a 3/4" gap between each pad of the pair. Never have the two pads touching each other during the massage session in operation, or you can risk hurting your skin.

Place the pads back on either side of the pads holder for proper storage.

1. Only use the electrode pads of the same brand as the controller device.

2. Do not use the same electrode pads on different people. Each user must have his/her own pads.

3. Always turn off the power before removing or re-positioning the electrode pads.

4. Wash skin thoroughly and then dry it before applying electrode pads.

5. Apply the whole surface of the electrodes firmly to the skin. Do not use electrodes that do not stick properly to the skin or only partially stick to the skin.

NOTE: Place electrode pads on skin area firmly before turning on device to avoid discomfort.

Features & Benefits

AS1080-16AB is not only a TENS Unit (Transcutaneous Electrical Nerve Stimulation), but also EMS (Electrical Muscle Stimulation). It's also known as PMS (Powered Muscle Stimulation). This means the AS1080-16AB unit is a combination device that provides both pain relief and muscle stimulation. As a compact and portable device, AS1080-16AB can be used at home or on-the-go.

It is small and powerful at the same time! In fact, you can choose from sixteen modes to personalize your pain relief session.

Those suffering with chronic pain related to muscles or nerves, tension, or even many kinds of arthritis, will find this unit quite helpful. The device will stimulate nerves to block pain sensations and release natural pain-relieving endorphins.

AS1080-16AB will also stimulate muscles for those looking to tone, firm, or enhance their muscles with its EMS modes available on the device. This portable unit is extremely easy to take with you and use just about anywhere, such as the gym for post workout recovery or work for pain relief.

Finally, those with any blood circulation issues (such as those who travel a lot or are sedentary often) will find the AS1080-16AB can increase the blood flow of the lower body parts such as the leg or foot. This provides much needed relief and better mobility by revitalizing body tissues and promoting blood flow recovery.

Functions (Intended for use)

Intended uses include the following:

- Temporary relief of pain associated with sore
- Aching muscles in the shoulder, waist, back, neck, upper extremities(arm)
- Lower extremities(leg) due to strain from exercise
- Lower extremities(leg) due to normal household work activities

Instruction for use

The following steps are a guide to use the device after setup, and the details for each step are listed in the table below.

- Step 1: Slide "ON/OFF" switch to turn on the power
- Step 2: Select one of the stimulation modes
- Step 3: Choose the stimulation time
- Step 4: Adjust the stimulation intensity for each channel Enjoy your massage therapy!

• "ON/OFF" button to turn on the power

Slide the On/Off button to turn on the unit.
This button is found on the top of the device.

• Selecting output channel

There are two output channels so that you can receive the massage at two different areas of the body at the same time. Each channel can be separately controlled with different massage modes and intensity settings. Switch the A-B channels by pressing the middle round button. The default channel is A.

• Selecting stimulation modes

Press the "M" button to switch between 16 massage types (indicated by the flashing image of the hand/hands on display) to the most desirable one. For safety reasons, the default intensity is set to the lowest each time after you make a new mode selection. The mode of each channel (A or B) can be set separately to the desirable one.

Selecting stimulation time

Press the "T" button down. Repeat until you have chosen the time you prefer for your session.

You can choose from 10, 20, 30, 40, 50, 60 minutes. The timer indication is on the upper left part of the display. 20 minutes is the default time. The device will shut off as a safety mechanism once the timer is complete.

• Adjust stimulation intensity

Press and release the "+" button to increase the stimulation intensity, and press and release the "-" button to decrease the intensity. The intensity of the two channels can be set separately to the desirable level. A or B will flash on the display to indicate which corresponding channel you are currently adjusting.

NOTE: 1. With the increase of intensity, you may experience sensations like tingling, vibration, etc. Gradually increase the intensity, and stop increasing when a comfortable level is reached.

NOTE: 2. Make sure you are increasing the correct channel. You can see the change of the intensity represented by the change in bars on the display. Press the Center Round Button to switch between channel A and B.

• Lock and unlock

Press and hold the middle round button for three seconds to lock the entire setting display and avoid accidental changes. Press and hold for three seconds to release the lock.

• Enjoy the stimulation massage therapy

After the above mode, time, and intensity are set to your preference, enjoy the stimulation provided by the device. Slide the "ON/OFF" to turn off the power to stop use, or the device will automatically shut off when the timer has completed.

Mode Reference Chart

MODE	TENS RESPONSE	PULSE FREQUENCY	TYPES
1	Kneading	On for 3.5s and off for 1.4s	EMS
2	Kneading Combination	On for 17.5s and off for 7s On for 26.4s	TENS & EMS
3	Tapping	On for 0.1s and off for 0.77s	TENS
4	Tapping Combination	On for 3.6s and off for 28s; On for 23.4s On for 3.68s and off for 1.88s	TENS & EMS
5	Cupping	On for 34.6s and off for 3.45s	TENS
6	Cupping Combination	On for 7s and off for 2.8s; On for 26.6s	TENS & EMS
7	Shiatsu-Deep	On for 7.3s and off for 3.1s	EMS
8	Shiatsu-Deep Combination	On for 10.5s and off for 4.2s; On for 2.7s and off for 1.35s	EMS
9	Acupuncture	On for 0.76s and off for 0.4s; On for 0.1s and off for 0.3s	TENS
10	Acupuncture Combination	On for 2.3s and off for 18.4s; On for 10.5s and off for 4.2s; On for 26.5s	EMS
11	Scraping	On for 17.5s and off for 0.7s	TENS
12	Scraping Combination	On for 0.828s and off for 3.6s; On for 4.92s and off for 2.82s	EMS
13	Palm Kneading	On for 0.138s and off for 0.6s; On for 1s and off for 0.8s	TENS
14	Palm Kneading Combination	On for 3.6s and off for 8.28s; On for 23.4s; On for 14s and off for 5.6s	TENS & EMS
15	Massage	On for 1.3s and off for 0.9s	EMS
16	Massage Combination	On for 4s and off for 2.5s; On for 17.5s and off for 7s	EMS

Position Placement for Best Use

You can place TENS Pads on many areas of the body to relieve muscle and nerve pain and to stimulate healthy muscles. Some areas are great for the TENS Mode, such as shoulders, back,

waist, arms, legs, feet and hips. Muscle stimulation can be applied to many muscle areas on the body, such as abdomen, buttocks, arms and legs.

LOW BACK

Sitting for too long, having weak hip flexors, or exercising can aggravate pain in our hips and waist. To alleviate tension or pain, attach the electrode pads around your waist and lower back. Space the pads from 5cm to 15cm apart for the best results.

WAIST

TENS is perfect for relief from pain in the waist area. Whether weak hip flexors or a sore side muscle, the aching muscles can be alleviated.

The TENS can soothe the nerves and help with hip pain from issues with tendons and muscles by providing gentle movement.

BUTTOCKS

The EMS mode can be used for strengthening and tightening muscle on the buttocks or hip areas. TENS mode can be used to alleviate tension or pain on the hip and buttock too. Attach the one pair of pads at the lower back and another pair on the hip area. Space the pads from at least 5cm apart.

ABDOMEN

The EMS feature can help stimulate the gluteus muscles for better athletic performance. With the contraction of the muscles this device provides, the muscles in this area can be strengthened. This area is sensitive so it does not need much intensity to feel the effects.

ARMS

Use EMS mode to strengthen arm muscle like biceps by attaching one pair at the upper end of the muscle next to the shoulder and another pair on the lower end of the muscle. Gradually increase the intensity until you get a maximum muscle contraction, without being painful.

THIGHS/LEGS

Bicyclists know the incredible stress and strain that your thigh muscles can experience. Since these large muscles are so crucial to our mobility, it is important to take care of them and help them relax. With TENS mode, placing the electrode pads on the top and back side of your thighs is an uncomplicated way to relax and de-stress this area.

Use EMS mode to strengthen leg muscles like hamstrings and quadriceps by attaching one pair at the upper end of the muscle and another pair on the lower end of the muscle. Gradually increase the intensity until you get a maximum muscle contraction, without being painful.

Practices for Best Use

- **Electrode pads**

The electrode pads can be used until they lose the stickiness and/or conductivity. Their lifetime greatly depends on the body skin condition and protection.

Purchase new electrode pads as needed.

- **Low intensity to start**

With the increase of intensity, you may experience sensations like tingling, vibration, etc. Therefore, gradually increase the intensity, and stop increasing when a comfortable level is reached.

- **Cleanliness**

Keep the device clean by wiping with a moistened cloth from time to time. This prevents growth of germs and dirt and can extend the life of the pads.

- **Storing device**

When not in use, store the device and accessory in a cool place, out of direct sunlight.

- **Desired results**

Most find using the TENS once a day from 10-30 minutes a day provides the best benefits. You can choose your own intensity level and sensation mode based on your own comfort level.

- **Areas to avoid**

Do not use your TENS Unit on your head, neck, chest near heart, or genitals.

Do not apply stimulation across the chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the heart. For instance, DO NOT apply one set of TENS pads on one arm and the other set of TENS Pads on the other arm because the electrical current can go through your heart from one arm to another arm.

Frequently Asked Questions

How does my TENS unit relieve pain or stimulate muscles?

How does it work?

AS1080-16AB relieves pain by blocking pain signals and stimulating natural painkillers known as endorphins. This unit can also stimulate muscles, which helps with muscle performance and blood circulation in the lower extremities.

How many modes and intensity levels does the AS1080-16AB have?

This unit has 8 massage types and 16 auto procedure modes. The device has 20 intensity levels you can choose from, so you can get benefits at your own level of comfort.

Is it safe for anybody to use my TENS anywhere on the body?

The AS1080-16AB unit is an effective and drug-free device for pain relief and muscle stimulation when used properly. But, it is not suitable for some body areas or some people. Please read the Safety Warning (on page 15) before using this unit. If you have any medical condition and are unsure about using the device safely, consult your physician first.

What is the difference between EMS and TENS? Does this unit provide both?

The AS1080-16AB Unit is a TENS and EMS Combo Unit. EMS stands for Electrical Muscle Stimulation and TENS stands for Transcutaneous Electrical Nerve Stimulation. TENS modes on the device are specifically for pain and soreness.

The EMS modes are for muscle stimulation and better blood circulation in the lower body for healthy muscles. Yes, this device provides both.

What can I do to elongate the life of my electrode pad?

Cleanse the area of skin you will be placing the electrode pads on with soap and water, or with a damp cloth. Make sure the area is dry before applying the electrodes. When storing the electrode pads, we recommend placing them back onto the plastic holders. These come as accessories in the package. Electrode pads that become wet from sweat or water, and do not stick anymore need to be replaced.

You can also cleanse the pads with a damp cloth. If they are too dirty or your skin feels numb after the pads are cleaned, replace the pads.

What type of conditions does TENS/EMS really benefit?

TENS is useful in relieving many kinds of pain and soreness. TENS is great for lower back pain, knee pain, sciatica, fibromyalgia. These are just some examples.

EMS modes are great for improving muscle performance or increasing blood circulation in the lower extremities from a sedentary lifestyle if your muscles lack tone and firmness.

Cleaning and maintenance

Please use wipes to clean the device first, and then use the dry cloth to wipe it again. The electrode pads coming with the device are disposable and should be replaced when they lose their adhesiveness. Contact the seller for replacements.

Troubleshooting

Problem	Possible Cause	Solution
One pad feels stronger than the other.	This is normal. Different areas of your body will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.

During the massage, the skin feels a painful burning sensation or the stimulation becomes weakened	Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry.	Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application.
Unit is on, but no sensation is felt in the pads.	Both pads are not firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased.	Make sure both pads are firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.
Adhesive gel pads do not stick to skin even after cleaning and moistening the gel pad.	Adhesive gel pads need to be replaced.	Replace the set of pads. Pads can be purchased from local distributor.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Does not provide adequate pain relief	Not using it long enough.	Use your massager 20-30 minutes at a time, 3-6 times per day.

Warranty Information AS1080-16AB

We are so confident that you will love this product that we offer a 1 Year PRODUCT Quality Guarantee.

Safety Warning

- **Contraindications**

Do not use this device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device, because this may cause electric shock, burns, electrical interference, or death.

Do not use this device if your pain syndromes are undiagnosed.

- **Contraindications**

Do not apply stimulation on your neck because this could cause severe muscles spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.

Do not apply stimulation across your chest, because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal. For instance, DO NOT apply one set of TENS pads on one arm and the other set of TENS Pads on the other arm because the electrical current may go through your heart from one arm to another arm.

Do not apply stimulation on your head, chest, near your heart, or genitals areas.

Do not apply stimulation over -- or in proximity to -- cancerous lesions.

Do not apply stimulation when the patient is in the bath or shower.

If you have one of the following conditions, please consult with your physician before purchasing or using this device:

Acute disease, malignant tumor, infective disease, pregnant, heart disease, high fever, abnormal blood pressure, lack of skin sensation or an abnormal skin condition, any condition requiring the active supervision of a physician.

- **Precautions**

Do not use this device while driving.

Do not use this device while sleeping.

Do not use this device in high humidity areas such as a bathroom.

Keep the device away from wetness, high temperatures and directly-sunlit places.

Keep this device out of the reach of children.

Stop using this device at once if you feel pain, discomfort, dizziness or nausea, and consult your physician.

Do not attempt to move the electrode pads while the device is operating.

Do not use the device around the heart or on the head, mouth, genitals, or blemished skin areas.

Do not apply stimulation of this device in the following conditions:

(1) Across the chest. The introduction of electrical current into the chest may cause rhythm disturbances to the heart, which could be lethal;

(2) Over painful areas. Please consult with your physician before using this device if you have painful areas;

(3) Over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombo-phlebitis, varicose veins). Apply stimulation only to normal, intact, clean, healthy skin;

(4) In the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms). The electronic stimulator may

not operate properly when the electrical stimulation device is in use;

(5) While operating machinery, or during any activity in which electrical stimulation can put you at risk of injury;

(6) On children

Be aware of the following:

(1) Consult with your physician before using this device. The simulation with the device may:

I. Cause lethal rhythm disturbances to the heart in susceptible individuals;

ii. Disrupt the healing process after a recent surgical procedure;

(2) The device is not effective for pain in your head area, including for headaches;

(3) The device is not a substitute for pain medications and other pain management therapies;

(4) The device has no curative value;

(5) The device is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism;

(6) The long-term effects of electrical stimulation are unknown;

(7) The user may experience skin irritation, burns or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel);

(8) If the user has suspected or diagnosed epilepsy, the user should follow precautions recommended by his or her physician;

(9) Use caution if the user has a tendency to bleed internally, such as following an injury or fracture;

(10) use caution if stimulation is applied over the menstruating uterus;

(11) Use caution if stimulation is applied over areas of skin that lack normal sensation;

(12) Stop using the device if the device does not provide pain relief;

(13) Use this device only with the leads, electrodes, and accessories that the manufacturer recommends;

(14) Do not share the use of the electrode pads with others;

(15) Do not use the device while it's charging;

(16) The device contains a lithium battery. If overheating of the device occurred during charging, stop charging or operation immediately and report to the seller;

(17) Dispose of the battery-containing device according to the local, state, or federal laws.

The long-term effects of electrical stimulation are unknown.

Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.

The safety of electrical stimulation during pregnancy has not been established.

Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel).

Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.

Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.

Use caution if stimulation is applied over the menstruating or pregnant uterus.

- **Adverse reactions**

You may experience skin irritation and burns beneath the stimulation electrodes applied to the skin;

You may experience a headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.

You should stop using the device and should consult with your physician if you experience adverse reactions from the device.

- **Environmental condition for normal working, transport and storage**

- Normal working ambient temperature: 5~40°C
- Normal working ambient humidity: 15%~90% RH
- Store and transport ambient temperature: -25 ~ 70°C
- Store and transport ambient humidity: 0%~90% RH
- Atmospheric pressure: 70~106kPa

Electromagnetic compatibility

The AS1080-16AB has been tested and found to comply with the electromagnetic compatibility (EMC) limits for medical devices to IEC 60601-1-2: 2014. These limits are designed to provide reasonable protection against harmful interference in

a typical medical installation

Caution:

Do not use this device simultaneously with devices having high EMI levels.

Guidance and manufacturer's declaration – electromagnetic emission

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used such an environment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT- GUIDANCE
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable (internal battery powered)	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable (internal battery powered)	

Guidance and manufacturer's declaration -electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE

Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV Contact ±8 kV air	±6 kV Contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	Not applicable (internal battery powered)	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not applicable (internal battery powered)	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Not applicable (internal battery powered)	Main power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. main voltage prior to application of the test level.

Guidance and manufacturer's declaration -electromagnetic immunity

The device is intended for use in the electromagnetic environment specified on following page. The customer or the user of the device should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 Mhz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $=1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 Ghz	3 V/m	

level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE1: At 80 MHz and 800MHz, the higher frequency range applies.
NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Technical Specifications

Mode/type	AS1080-16AB	Weight	60g
Power supply	Powered by internal 3.7V li-ion battery	Automatic shutoff	20 minutes
Wave form and wave shape	Biphasic rectangular wave pulse	Degree of protection against electric shock	Type BF applied part
Pulse duration	80us (Microse-conds)	Type of protection against electric shock	Internally powered equipment
Pulse frequency	0-100Hz (Hz=vibration per second)	Warranty	1 year
Output Voltage	Max.130Vpp±20% (at 1000ohm load)	Grade of waterproof	IP 22
Treatment time	10, 20, 30, 40, 50, 60 minutes	Lifetime for electrode	Storage for 2 years(no use), Times of reuse: 30 times

Modes	16 auto modes	Software version	A0
Typical operation time of Battery	If used at the highest level, the battery can be used for about 4 hours after fully charged.	The time required for my equipment to warm from the minimum storage temperature between uses until it is ready for intended use	30 minutes
Behavior of my equipment while the rechargeable internal electrical power source is charging:	The battery indicator on the display will indicate the battery charging status.	The time required for my equipment to warm from the minimum storage emperature between uses until it is ready for intended use	15 minutes
Typical service life of battery	300 times of recharging	Adapter for charging	Please use output DCSV and output current 0.3-2.0A adapter for charging

NOTE: Not intended to be sterilized.

Not for use in an OXYGEN-RICH ENVIRONMENT

Product programs

PROGRAM NAME	TIME MIN.	FREQUENCY (Hz)	PULSE WIDTH (s)
Mode 1	10,20,30 40,50,60	66.6	80

Mode 2	10,20,30 40,50,60	66	80
Mode 3	10,20,30 40,50,60	2	80
Mode 4	10,20,30 40,50,60	97	80
Mode 5	10,20,30 40,50,60	55	80
Mode 6	10,20,30 40,50,60	97.5	80
Mode 7	10,20,30 40,50,60	47	80
Mode 8	10,20,30 40,50,60	66	80
Mode 9	10,20,30 40,50,60	90	80
Mode 10	10,20,30 40,50,60	97	80
Mode 11	10,20,30 40,50,60	98	80
Mode 12	10,20,30 40,50,60	90	80
Mode 13	10,20,30 40,50,60	5	80
Mode 14	10,20,30 40,50,60	66	80
Mode 15	10,20,30 40,50,60	63	80
Mode 16	10,20,30 40,50,60	90	80

Symbol	Meaning
	Caution
	BF Isolated Applied Part
	Date of manufacture
	Disposal in accordance with Directive 2002/96/EC
	MANUFACTURER
EC REP	AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	Read the user manual before use
	CE marking
IP22	Protected against dripping water
LOT	Batch code
	Recyclable
	Temperature: -25°C~70°C
	Humidity: 0%~90%

Disposal of Materials

A symbol of Waste Electrical and Electronic Equipment means you should find out about follow local regulations about disposing this kind of product.

- Do not dispose this product as you would other household waste. Dispose this device in accordance with the corresponding local regulations.

- Electrical and electronic devices contain hazardous substances that can have harmful effects on the environment and /or human health and should be recycled properly.

Guarantee: 1 years

Healthcare-Manager.com

Questions or comments?

Please contact us via email:

service@healthcare-manager.com

Distributed by: Easy Healthcare Corporation

360 Shore Dr., Burr Ridge, IL USA 60527

Scan the code to get the manual in more languages

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Distribuido por: Easy Healthcare Corporation