

WHERE SCIENCE
MEETS LOVE®



Welcome to ASSISI


LOOP®

Targeting what matters most,
Healthier, happier pets

 Incubate
Animal Health
Looking after innovation

 ASSISI
ANIMAL
HEALTH

ASSISI



Congratulations on your purchase of the Assisi Loop™, and easy to use, non-pharmaceutical device for soothing joints and other body structures.

Your Assisi Loop is simple an easy to use. It contains no additional parts and the battery cannot be charged or replaced.



How to position the Loop

The Assisi Loop has an Effective Working Zone (EWZ) of 25cm for the 20cm diameter loop and 13cm for the 10cm diameter loop, each side of the applicator. It can be used on either side.

Place the centre of the Loop over the area you wish to concentrate the session on. If you are unsure please consult your veterinarian. The EWZ does not extend beyond the diameter of the loop, so ensure you are using the appropriate sized loop. For very large areas, repeat sessions may required, as you change the Loop's position.





How to activate the Loop.

Press the on/off button once. Both green lights will flash on the front of the signal generator indicating the Loop is working.



How to use the Loop.

Once the Loop has been correctly positioned and turned on, keep the Loop in the correct position for 15 minutes.

- The Loop will automatically turn off after 15 minutes.
- The Assisi Loop works directly through wraps, dressings and casts.



How often to use the Loop.

Follow the instructions of your veterinarian regarding how many times to use the Loop each day on your pet. Many veterinarians will suggest you start with 2-4 15 minute sessions per day but this could vary between individuals.

- Always use the Loop in accordance with instructions from your veterinarian.
- The Loop can be beneficial in assisting recovery following trauma or surgery. Follow the instructions of your veterinarian for frequency of use. Expect to use the Loop 3-4 times a day for the first week and then taper down over the following week until its use is no longer required.
- The Loop can help soothe joints and other body structures and aid with stiffness. Follow instructions of your veterinarian for frequency of use. Expect to use the Loop 3-4 times a day for up to 10 days or more, depending on the individual, and then reduce sessions to maintain comfort. Always refer to your veterinarian if you have questions regarding frequency of use.

The Assisi Loop helps maintain health and support recovery. It is designed to be used in a protocol as per recommendation from your veterinarian, following a thorough assessment of your pet.



How soon will I see results?

Depending on your pet and their need for using the Loop you may see a noticeable positive change in your pet after the first few sessions. Other indications for use, such as stiffness, may take up to two weeks before a change is noted. It is important to continue sessions with the Loop as recommended by your veterinarian.



Caution.

- Avoid using the Loop near large metal objects (metal crates, stainless steel tables, etc.) as doing so may reduce the beneficial effects of its use.
- Do not immerse in water - the Loop is not waterproof.
- Equipment NOT SUITABLE FOR USE IN THE PRESENCE OF FLAMMABLE ANAESTHETIC MIXTURE WITH AIR or WITH OXYGEN or NITROUS OXIDE or OXYGEN ENRICHED ATMOSPHERES.
- The Assisi Loop SHOULD NOT BE USED on patients who have ANY implanted system that may contain an electrical lead, such as a pacemaker.



Cleaning Instructions.

- Loop components are not waterproof. Do not immerse.
- Turn off to clean with a soft cloth dampened with 70% isopropyl alcohol.
- Wait 15 minutes after cleaning before operating the unit.



Battery Care.

- The Assisi Loop battery is non rechargeable and cannot be replaced.
- For optimal battery life wait at least two hours between uses.
- The Assisi loop offers 150 15-minute uses. Quick flashing lights indicate the unit is reaching the end of its usable life.

Please contact your veterinarian if you have any questions

**For more information,
0330 1743194
info@incubateanimalhealth.com
incubateanimalhealth.com**



ASSISI

LOOP[®]

Contact:
0330 1743194
info@incubateanimalhealth.com
incubateanimalhealth.com

