



## TRICHROME STAIN

### WHEATLEY AND PARASITOLOGY KIT COMPONENTS:

#### PARASITOLOGY STARTER KIT (VXF-019) REFILL (VXR-010)

SCHAUDINN'S FIXATIVE: CARBOL XYLENE:		ETHANOL 70%:	ETHANOL 95%:
VXM-016 (16 oz.)	VXZ-016 (16 oz.)	VXB-032 (32 oz.)	VXC-032 (32 oz.)
VXM-032 (32 oz.)	VXZ-032 (32 oz.)	VXB-128 (128 oz.)	VXC-128 (128 oz.)
IODINE ALCOHOL: ACID ETHANOL 90%:		TRICHROME STAIN:	XYLENE:
VXI-016 (16 oz.)	VXA-016 (16 oz.)	VXT-008 (8 oz.)	VXX-016 (16 oz.)
VXI-032 (32 oz.)	VXA-032 (32 oz.)	VXT-016 (16 oz.)	VXX-032 (32 oz.)
VXI-128 (128 oz.)	VXA-128 (128 oz.)	VXT-032 (32 oz.)	

**PVA MODIFIED FIXATIVE:** VPM-016 (16 oz.)/ VPM-032 (32 oz.)/ VPM-128 (128 oz.)

Order online at [www.volusol.com](http://www.volusol.com)

**INTENDED USE:** For use in the rapid staining and differentiation of intestinal parasites. According to the Gomori procedure it can be used to counterstain histologic sections.

**SUMMARY AND PRINCIPLES:** Claude Masson introduced the first trichrome stain in 1929. It offered a simple, quick alternative to the tissue stains previously used. In 1949, George Gomori developed a single solution trichrome stain as a counterstain for hematoxylin stained tissue sections and cytological smears. Peter Wheatley's observation of *E. Histolytica* stained with chromotrope 2R in 1951 yielded the Wheatley methodology. The modification of traditional fixation and dehydration procedures, combined with trichrome stain, resulted in a procedure which stains amoeba and flagellates rapidly and easily. Volu-Sol's Trichrome Stain is recommended for use according to the Wheatley procedure for parasites, as well as the Gomori procedure for tissues. Volu-Sol's PVA and Schaudinn's Modified Fixatives are copper based and do not contain mercuric chloride. This modification eliminates the iodine-alcohol step typically used in the trichrome stain kit.

**SPECIMEN PREPARATION: WHEATLY STAINING PROCEDURE FOR FECAL PARASITES:** Fresh fecal specimens should be collected in a clean, dry container and prepared as soon as possible. If a specimen cannot be prepared immediately, it is suggested that it be fixed in PVA modified fixative. Make a thin smear of fresh fecal material on a clean, glass slide. The specimen may be emulsified in a drop of saline if necessary. If a PVA-fixed specimen is used, allow the smear to air-dry overnight and proceed to the staining procedure (heating the slide may cause distortion of the organisms). While smear is still wet, place the slide into a coplin jar containing Schaudinn's Modified Solution for one hour. **\*NOTE:** The fixing time may be shortened to 5 minutes if Schaudinn's Modified Solution is heated to 50°C. Schaudinn's Modified fixed slides must not dry before staining.\* **GOMORI TRICHROME STAINING PROCEDURE FOR TISSUE AND**

**SMEARS:** Fix smears or spread preparation in alcohol (alcohol-ether) or imbed tissue block in paraffin and section at 3-5 microns. **\*NOTE:** Change all solutions periodically to prevent "watering down" and/or carryover. Weakened trichrome stain may be strengthened by exposing to air overnight, or replenishing with fresh stain.\*

**PRINCIPLES OF THE PROCEDURE: WHEATLEY PROCEDURE:** PVA fixed fecal smears are passed through alcohol solutions and stained to produce a permanent slide for observation and identification of parasites. After mounting, slides may be stored indefinitely.

**GOMORI PROCEDURE:** Paraffin imbedded specimens are sectioned; deparaffinized and hydrated in distilled water; stained with hematoxylin and counterstained with trichrome stain. After mounting, store slides indefinitely. **\*NOTE: SCHAUDINN'S MODIFIED SOLUTION:** On the day of use, add 1/10 volume of glacial acetic acid to 9/10 volume of Schaudinn's modified solution. Stable for only one day. All slide immersion requires occasional gentle agitation.\*

**PROCEDURE: WHEATLEY'S PROCEDURE FOR INTESTINAL PARASITES:** Set up coplin jars for the following solutions: **1) ETHANOL, 70%:** Immerse slides for 3-5 minutes. **2) ETHANOL, 95%:** Immerse slides in second solution for 3-5 minutes. **3) TRICHROME STAIN:** Immerse slides for 6-8 minutes. **4) ACID ALCOHOL DECOLORIZER:** Transfer to the next solution when stain runs from smear. Two quick dips. Carryover of stain and/or decolorizing solution can cause cloudiness on slides and/or lack of contrast. **5) ETHANOL, 95%:** Dip slides twice. **6) ETHANOL, 95%:** Immerse slides for 10 minutes. **7) CARBOL XYLENE:** Immerse slides for 3 minutes. **8) XYLENE:** Immerse slides for 3 minutes. Mount smear with coverslip, using a mounting medium of choice. Slides may be screened for parasites as soon as dry. **GOMORI PROCEDURE FOR TISSUE:** Bring smears or sections to water. Stain slides with hematoxylin modified for 5 minutes. Remove excess stain by washing slides in gentle running water. Stain slides with trichrome stain for 5-10 minutes. Wash slides with several quick dips in acetic acid, 0.2%. Dehydrate and clear, using alcohols and xylene. Mount and coverslip.

**EXPECTED RESULTS: FECAL SMEARS:** Background material and artifacts will stain green. Bacterial and red blood cells usually stain red, while yeast, vegetable fibers, etc. stain green. Organisms and cysts stain blue-green with a purplish tinge. Karyosomes of the nuclei, chromatic material will appear red-to-purple. Helminth eggs and larvae also stain red-to-purple. **TISSUE SECTIONS OR SMEARS:** Connective tissue stains green, muscle and cytoplasm red, and nuclei blue.

**STORAGE AND EXPIRATION:** Store reagents at room temperature (70-77.9 °F/ 20-25.5 °C). Maximum intended shelf life is printed on the label. If necessary to interrupt staining procedure, slides may be stored in ethanol 70% without adverse effects.

#### REFERENCES:

1. Lillie R.D., 1965 Histopathology Techniques and Practical Histochemistry. McGraw Hill 3<sup>rd</sup> Ed.
2. Gomori, G., 1950 A Rapid one step Trichrome method: Am J. of Clin. Path. Vol. 20: 661-664