NOTA INFORMATIVA cod. 390017 - mod. ASTM Blu

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INFORMATIVE NOTE cod. 390017 - mod. ASTM Blu

Please read this briefing note through carefully prior to use and prior to any maintenance work. The purpose of the information contained herein is to assist and direct the user in choosing and using the PPE. Neither the manufacturer nor the distributor can be held liable for an incorrect use of the PPE. This briefing note must be kept for the duration of the PPE. Glove model: **cod. 390017 - mod. ASTM Blu.**

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The CE 0465 mark means that these gloves are personal protective equipment conforming to the essential health and safety requirements set forth in the (EU) Regulation 2016/425 and have been certified by the notified body: **A.N.C.I. Servizi Srl - Sezione C.I.M.A.C. - C.so G. Brodolini, 19 - 27029 Vigevano (PV), Italia, n° 0465.** In addition, the manufacturer is subject to the control of CE quality guarantee of the finished product by the notified body **A.N.C.I. Servizi Srl - Sezione C.I.M.A.C. - C.so G. Brodolini, 19 - 27029 Vigevano (PV), Italia, n° 0465.**

SIZES

OUR SIZE REFERENCES	SIZES CONTEMPLATED EN 420:2003+AI:2009
S	7
M	8
L	9
XL	10

The dexterity level of the fingers is equivalent to 5 out of 5.

The gloves meet the safety characteristics only if the size is right and worn properly.

RANGE OF USE

The glove is designed to protect the user against contact with chemical products and micro-organisms (intended as fungus or bacteria, not viruses).

PICTOGRAMS AND PROTECTION LEVELS

PROTECTION AGAINST CHEMICALS

PROTECTION AGAINST MICRO-ORGANISMS

Resistance to air and water penetration using the method referred to in EN 374-2:2016 standards

CHEMICAL RESISTANCE TABLE according to EN ISO 374-1: 2016:

Chemical	Permeation index [0 - 6]	Time of Passage (min)	Degradation index
G Diethylamine	2	30	74.7%
K Sodium hydroxide 40%	3	60	29.2%
L Sulphuric acid 96%	2	30	99.8%

Time of passage.

The time a chemical takes to pass through the glove material.

The values may differ when the chemical substance is used in a compound.

Do not use the gloves for handling the above-mentioned chemical products for longer than the times stated above.

COMPLIANT WITH: REGULATION (EU) N. 10/2011 - REGULATION (EC) N. 1935/2004 - DIRECTIVES 85/572/EEC AND SUBSEQUENT AMENDMENTS.			
Ethanol 20%	Acetic acid 3%	Vegetable oil	Polyphenyl oxide [MPP0]
Yes	No	No	No
2 h / 40°C	-	-	-

Since the performance levels are based on lab test results, they may differ from the actual conditions at the workplace. Consequently it is advisable to conduct a preliminary test of the gloves to see if they are suitable for their intended use.

The protective characteristics refer to new PPE, in good condition and that has never been subject to cleaning treatments.

USE

Prior to use, make sure the gloves are in good condition by which we mean no cuts, holes or cracks, etc. If the gloves are not in perfect condition, replace them immediately. The gloves must be used solely for the risks contemplated in this briefing note. Movements, tears, frictions and degradation caused by contact with the chemical substance may significantly reduce the actual time of use. In the case of corrosive agents, degradation may be the most important factor to consider when choosing the gloves.

STORING

The gloves are packed in cardboard boxes.

Keep the gloves in their original packaging in a clean, dry place, away from sources of heat and direct sunlight. The user must always inspect the gloves to ensure they are sound prior to using them.

The glove must be used, provided the packaging is perfectly intact and kept in the preservation conditions stated on the packaging itself, **within the indicated expiry date and no later.** Since these are disposable gloves to be used just the once, they must be replaced at the end of their use – they cannot be used a second time.

NOTES

Used gloves can be contaminated with chemicals or biological substances: dispose of them in compliance with the local regulations in force on this subject [dump, incinerator]. People who are sensitive to rubber or latex should not use these gloves.

If an allergic reaction occurs, stop using the gloves and seek medical advice. In the case of manufacturing defects this PPE will be replaced.

*** In the event of discrepancies between the different translations, the Italian version will exclusively apply.**

For more information, contact:

Registered Community Trade Mark, no 016928426
c/o EUIPO - Alicante - Spain

The EU conformity declaration is available at following address:
www.nerispa.com

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