

USER GUIDE: Malem Deluxe (MO24)



PERSONALISED RECORDING

MESSAGE + SOUND

WAKES DEEP SLEEPERS

Set-Up

- Connect the sensor to the base of the alarm
- Clip the sensor onto your child's underwear, in the place most likely to receive the first drops of urine
- Attach the alarm to your child's pyjama top, near the collar bone where they will hear it sound & feel it vibrate



A 5 second message can be recorded

Microphone




Use a pointed object to press down recording button

Recording A Message

- Use a pointed object to press down the record button, on the back of the alarm
- Hold the microphone hole close to your mouth and speak into it to record your message

Changing The Alert Settings

- The alarm can be set to alert by:
MESSAGE, SOUND, LIGHT & VIBRATE (V&S)
LIGHT & VIBRATE (V)
- When the sound (V&S) setting is selected, 1 of 8 loud alarm sounds plays when the alarm is activated, followed by the customised recorded message. This is repeated until the alarm is reset.



The combination of an alarm sound and a familiar voice is very effective for waking up children, even deep sleepers!

The alarm sound cycles through 8 different sounds so that your child doesn't get used to the same noise

the **bedwetting** doctor

www.thebedwettingdoctor.com