Anodyne[®]
The Anodyne[®] Shirt

TABLE OF CONTENTS

English 	2-3
Dansk 	4-5
Deutsch   	6-7
Français 	8-9
Nederlands  	10-11
Norsk 	12-13
Svenska 	14-15
Suomi 	16-17
Español 	18-19
Italiano 	20-21

INSTRUCTION MANUAL

Do I have the right size?

At first glance, the Anodyne® Shirt may seem (far) too small, but this is a deliberate design decision. The T-shirt has been created with a "seamless" structure that creates a posture-correcting effect through its elastic material and unique design. It should fit very tightly to achieve the desired effect. If you are in doubt about whether you have received the correct size, you can verify it by looking inside the T-shirt at the bottom seam. Don't worry about stretching or damaging it during use; it's designed to withstand it.

How should the product feel?

The primary effect of the Anodyne® Shirt comes from the material's dynamic "tension" and design placed on the shoulders and upper back. Most people immediately feel their shoulders pulled slightly back, the shoulder blades coming together, and they stand up straighter.

If you feel no effect, try wearing the product for a longer time. If this doesn't work, it could be because the product is too big for you. If you experience discomfort, on the other hand, it may be because the product is too small for you, or you've used it for too long initially. Our size guide is accurate and fits most, but everyone is built differently. Read more about our recommendations for using the Anodyne® Shirt below.

Instructions of use

- The Anodyne® Shirt improves your posture by stimulating and activating the muscles, which can naturally cause slight soreness. Your body needs time to adjust to a new posture, so don't use your Anodyne® Shirt for too long initially.
- We recommend the following usage times:
 - Day 1-9: 0.5-2 hours every other day.
 - Day 10-19: 2-6 hours a day.
 - Day 20+: 6-10 hours a day.
- Remember that everyone is built differently. Listen to your body and adjust the above usage times accordingly.
 - If you experience soreness or discomfort, just remove the product for the day and reduce the usage time the following day.
 - If you don't feel any effect, on the other hand, increase the usage time or the frequency with which you wear your Anodyne® Shirt.
- How much our customers use their Anodyne® Shirt in the long run is individual and depends on their specific needs. Some use the product 1-2 days a week, while others use it 5-7 days a week.

- If you have the right size based on your body type, you will soon have a good sense of whether the product is the right choice for your needs.

Washing instructions

- For best durability, hand wash the product. However, it can also be machine washed on a gentle cycle at a maximum of 30 degrees. It's recommended to use a laundry bag.
- Do not use fabric softener.
- Your Anodyne® Shirt should be air-dried. It should not be tumble dried as this will damage the product.
- The product has been pre-washed and therefore won't shrink. The shirt can stretch slightly with use and over time, but not to the extent that it affects its function.
- If you still experience unusual odor, contact us for advice. There are special detergents called sports detergent, developed for washing this type of clothing.

Do you have other questions?

- Check out our FAQ (Frequently Asked Questions) where you can learn more about Anodyne® products, size guides, and usage instructions.
- You're also very welcome to contact our customer service for personal guidance. Both our FAQ and contact details can be found on our website under the menu item "Customer Service".
- If you experience unintended side effects, it's important to contact our MDR-hotline at +45 5389 0017.

Thank you for shopping at Anodyne®/ActivePosture®.

We hope you enjoy your new Anodyne® Shirt.

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