# ***TriScroll™ 800 Inverter Dry Scroll Vacuum Pump***

*Installation and Operation Manual*

Manual No. 699904435

Revision C.00

February 2021

# Table of Contents

|                                                                        |           |
|------------------------------------------------------------------------|-----------|
| <b>Warranty</b> .....                                                  | <b>v</b>  |
| <b>Instructions for Use</b> .....                                      | <b>1</b>  |
| General Information .....                                              | 1         |
| Installation .....                                                     | 3         |
| Use .....                                                              | 3         |
| Maintenance .....                                                      | 3         |
| Environmental Conditions, Compliance and<br>Utility Requirements ..... | 4         |
| <b>Instructions d'utilisation</b> .....                                | <b>5</b>  |
| Avertissements CEM .....                                               | 6         |
| Installation .....                                                     | 7         |
| Utilisation .....                                                      | 7         |
| Entretien .....                                                        | 8         |
| <b>Instrucciones de Uso</b> .....                                      | <b>9</b>  |
| Información general .....                                              | 9         |
| Instalación .....                                                      | 11        |
| Uso .....                                                              | 11        |
| Mantenimiento .....                                                    | 11        |
| <b>Istruzioni per l'uso</b> .....                                      | <b>12</b> |
| Avvertenze CEM .....                                                   | 13        |
| Installazione .....                                                    | 14        |
| Utilizzo .....                                                         | 14        |
| Manutenzione .....                                                     | 15        |
| <b>Anwendungshinweise</b> .....                                        | <b>16</b> |
| Hinweise zur elektromagnetischen Verträglichkeit .....                 | 17        |
| Installation .....                                                     | 18        |
| Betrieb .....                                                          | 18        |
| Wartung .....                                                          | 19        |
| <b>使用上の注意</b> .....                                                    | <b>20</b> |
| 一般情報 .....                                                             | 20        |
| 設置 .....                                                               | 22        |
| 使用 .....                                                               | 22        |
| メンテナンス .....                                                           | 22        |
| 使用说明 .....                                                             | 23        |
| <b>包装</b> .....                                                        | <b>23</b> |
| CEM 警告 .....                                                           | 24        |
| <b>安装</b> .....                                                        | <b>25</b> |
| 用途 .....                                                               | 25        |
| 维护 .....                                                               | 26        |
| Mechanical Connections .....                                           | 28        |
| Operation .....                                                        | 29        |
| TriScroll 800 Inverter Pump Components View .....                      | 30        |
| <b>Technical Information</b> .....                                     | <b>31</b> |
| Pump Electrical Controller .....                                       | 33        |
| Purge Kit .....                                                        | 36        |

|                               |    |
|-------------------------------|----|
| Exhaust Silencer Kit .....    | 36 |
| Vibration Isolation Kit ..... | 36 |
| HEPA Inlet filter .....       | 36 |
| Troubleshooting .....         | 37 |
| Maintenance .....             | 38 |

## Request for Return Health and Safety Certification

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# Warranty

Products manufactured by Seller are warranted against defects in materials and workmanship for twelve (12) months from date of shipment thereof to Customer, and Seller's liability under valid warranty claims is limited, at the option of Seller, to repair, to replace, or refund of an equitable portion of the purchase price of the Product. Items expendable in normal use are not covered by this warranty. All warranty replacement or repair of parts shall be limited to equipment malfunctions which, in the sole opinion of Seller, are due or traceable to defects in original materials or workmanship. All obligations of Seller under this warranty shall cease in the event of abuse, accident, alteration, misuse, or neglect of the equipment. In-warranty repaired or replaced parts are warranted only for the remaining unexpired portion of the original warranty period applicable to the repaired or replaced parts. After expiration of the applicable warranty period, Customer shall be charged at the then current prices for parts, labor, and transportation.

Reasonable care must be used to avoid hazards. Seller expressly disclaims responsibility for loss or damage caused by use of its Products other than in accordance with proper operating procedures.

Except as stated herein, Seller makes no warranty, express or implied (either in fact or by operation of law), statutory or otherwise; and, except as stated herein, Seller shall have no liability under any warranty, express or implied (either in fact or by operation of law), statutory or otherwise. Statements made by any person, including representatives of Seller, which are inconsistent or in conflict with the terms of this warranty shall not be binding upon Seller unless reduced to writing and approved by an officer of Seller.

## Warranty Replacement and Adjustment

All claims under warranty must be made promptly after occurrence of circumstances giving rise thereto, and must be received within the applicable warranty period by Seller or its authorized representative. Such claims should include the Product serial number, the date of shipment, and a full description of the circumstances giving rise to the claim. Before any Products are returned for repair and/or adjustment, written authorization from Seller or its authorized representative for the return and instructions as to how and where these Products should be returned must be obtained. Any Product returned to Seller for examination shall be prepaid via the means of transportation indicated as acceptable by Seller. Seller reserves the right to reject any warranty claim not promptly reported and any warranty claim on any item that has been altered or has been returned by non-acceptable means of transportation. When any Product is returned for examination and inspection, or for any other reason, Customer shall be responsible for all damage resulting from improper packing or handling, and for loss in transit, notwithstanding any defect or non-conformity in the Product. In all cases, Seller has the sole responsibility for determining the cause and nature of failure, and Seller's determination with regard thereto shall be final.

If it is found that Seller's Product has been returned without cause and is still serviceable, Customer will be notified and the Product returned at its expense; in addition, a charge for testing and examination may be made on Products so returned.

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# Instructions for Use

## General Information

This equipment is designed for use by trained personnel. The user should read this instruction manual and any other additional information supplied by Agilent before operating the equipment. Agilent will not be held responsible for any events that occur due to non-compliance with these instructions, improper use by untrained persons, non-authorized interference with the equipment, or any action contrary to that provided for by specific national standards.

The TS800 is a dry scroll vacuum pump. This pump is suitable for pumping air or inert gases. The pump is not intended to pump corrosive, explosive, or particulate-forming gases.

The following paragraphs contain all the information necessary to guarantee the safety of the operator when using the equipment. Detailed information is supplied in “Technical Specifications” on page 33.

This manual uses the following standard safety protocol:

**ADVERTENCIA** *The warning messages are for attracting the attention of the operator to a particular procedure or practice which, if not followed correctly, could lead to serious injury.*

**PRECAUCIÓN** *The caution messages are displayed before procedures, which if not followed, could cause damage to the equipment.*

**NOTA** *The notes contain important information taken from the text.*

## Storage

When transporting and storing the pump, the following environmental requirements should not be exceeded:

Temperature:  $-20^{\circ}\text{C}$  to  $+60^{\circ}\text{C}$  ( $-4^{\circ}\text{F}$  to  $140^{\circ}\text{F}$ )  
Relative humidity: 0 to 95% (non-condensing)

## Preparation for Installation

The pump is supplied in a special protective packing. If this shows signs of damage, which may have occurred during transport, contact your local sales office.

Total weight of the packing, TS800 pump included, is approximately 42 kg (93 lbs).

**ADVERTENCIA** *The TS800 pump itself weighs 32 kg (71 lbs). To avoid injury, use lifting aids and proper lifting techniques when moving the pump.*

**ADVERTENCIA** *When unpacking the pump, be sure not to drop it and avoid any kind of sudden impact or shock vibration to it.*

**NOTA** *Normal exposure to the environment cannot damage the pump. Nevertheless, it is advisable to keep the pump inlet closed until the pump is installed in the system.*

## EMC Warnings

### EN 55011 Class A Warning

This is a Class A product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

## FCC

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesirable operation.

### NOTE

*The equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generated, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is also likely to cause harmful radio communications interference in which case the user will be required to correct the interference at his own expense.*

## Installation

The TS800 is *NOT* intrinsically safe. Do not install or use the pump in an environment exposed to atmospheric agents (rain, snow, ice), dust, aggressive gases, or in explosive environments or those with a high fire risk.

During operation, the following environmental conditions must be respected:

Temperature: +5 °C to +40 °C (41 °F to 104 °F)  
Relative humidity: 0 to 95% (non-condensing)

Connect the pump to the power supply using an IEC-320 style power cord of at least 10 A capacity (see Table 1 on page 27).

### CAUTION

*When installing the vacuum pump it is recommended that the pump be positioned in a way that the power switch and power cord connection remain easily accessible to the user in order to be available as an option to disconnect power to the pump quickly if needed.*

## Use

In order to reach maximum vacuum, the pump must be left running for at least an hour with the inlet sealed.

There are no special instructions for starting the pump; it need only be switched on using the On/Off switch.

The optional inlet valve will open 30 seconds after the pump is started.

### WARNING

*The pump is designed for operation with neutral or noncorrosive fluids. It is absolutely forbidden to use it with potentially explosive or inflammable substances.*

There are no special instructions for stopping the pump; it need only be disconnected from the electric power source by the On/Off switch.

The optional inlet valve will close immediately after the pump is stopped.

## Maintenance

Personnel responsible for pump operation and maintenance must be well-trained and aware of the accident prevention rules.

### WARNING

- *Death may result from contact with high voltages. Always take extreme care and observe the accident prevention regulations in force.*
- *When machine is powered up, be careful of moving parts and high voltages.*
- *If you have to perform maintenance on the pump after a considerable time in operation, allow it to cool as the temperature of the outer surface may be in excess of 60 °C.*
- *Always disconnect your power supply to the pump before beginning maintenance work.*

### NOTE

*Before returning the pump to the factory for repair, the “Request for Return” sheet attached to this instruction manual must be completed and sent to the local sales office. A copy of the sheet must be inserted in the pump package before shipping.*

If a pump is to be discarded, it must be disposed of in accordance with specific national standards.

### Disposal

Meaning of the WEEE logo found in labels.

The following symbol is applied in accordance with the EC WEEE (Waste Electrical and Electronic Equipment) Directive. This symbol (valid only in countries of the European Community) indicates that the product it applies to must NOT be disposed of together with ordinary domestic or industrial waste but must be sent to a differentiated waste collection system.

The end user is therefore invited to contact the supplier of the device, whether the Parent Company or a retailer, to initiate the collection and disposal process after checking the contractual terms and conditions of sale.

For more information refer to: <http://www.agilent.com/environment/product/index.shtml>.

# Environmental Conditions, Compliance and Utility Requirements

## Environmental Conditions

This equipment meets the following IEC classifications:

- Pollution Degree 2 (See Note.)
- Installation Category II (See Note.)
- Equipment Class 1

### NOTE

*Pollution level describes the degree to which a solid, liquid, or gas which deteriorates dielectric strength is adhering. 2 applies to a normal indoor atmosphere.*

*Installation category implies the regulation for impulse withstand voltage. It is also called the Over voltage category. II applies to electrical equipment.*

*This equipment requires the following space for ventilation, maintenances access, and easy access to the Main Power Breaker Switch. There must be a clear space of at least 60 cm (24") on all sides of the equipment. The bench in your laboratory must be able to support the entire system and other laboratory equipment.*

## Electromagnetic Compatibility

### EN55011/CISPR11

**Group 1 ISM equipment:** Group 1 contains all ISM equipment in which there is intentionally generated and/or used conductively coupled radio-frequency energy which is necessary for the internal functioning of the equipment itself.

**Class A equipment** is equipment suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

This device complies with the requirements of CISPR11, Group 1, Class A as radiation professional equipment. Therefore, there may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.

Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try one or more of the following measures:

- Make sure that appropriate cables are used to connect the device to peripheral equipment.
- Changes or modifications not expressly approved by Agilent Technologies could void the user's authority to operate the equipment.

## South Korean Class A EMC declaration

This equipment is Class A suitable for professional use and is for use in electromagnetic environments outside of the home.

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## ICES/NMB-001

This ISM device complies with Canadian ICES-001.

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<http://www.agilent.com/environment/product/index.shtml>

## Power Cord

Several power cord options are available from your Agilent dealer. Descriptions of the available power cords and their ordering numbers are given in Table 1.

### NOTE

*The pump must be connected to the power supply using a high voltage IEC-320 type power cord of at least 10 A capacity.*

**Table 1 Power Cord Selection**

| Country       | Power Cord Specification  | Order     |
|---------------|---------------------------|-----------|
| Europe        | 10 A / 220-230 VAC, 2.5 m | 656494220 |
| Denmark       | 10 A / 220-230 VAC, 2.5 m | 656494225 |
| Switzerland   | 10 A / 230 VAC, 2.5 m     | 656494235 |
| UK/Ireland    | 13 A / 230 VAC, 2.5 m     | 656494250 |
| India         | 10 A / 220-250 VAC, 2.5 m | 656494245 |
| Israel        | 10 A / 230 VAC, 2.5 m     | 656494230 |
| North America | 10 A / 230 VAC, 2.5 m     | 656494255 |

## Grounding Instructions

This product should be grounded. In the event of an electrical short circuit, grounding reduces the risk of electric shock by providing an escape wire for the electric current. This pump is equipped with a power cord that has a grounding wire with an appropriate grounding plug. The plug must be inserted into an outlet that is properly installed and grounded in accordance with all local codes and ordinances.

### WARNING

*Check with a qualified electrician or serviceman if the grounding instructions are not completely understood, or if you are in doubt as to whether the product is properly grounded.*

*Do not modify the plug provided; if it does not fit the outlet, have the proper outlet installed by a qualified electrician.*

*Connect the product only to an outlet that has the same configuration as the plug.*

*Do not use an adapter with this product.*

## Extension Cords

If you must use an extension cord with this product:

- Agilent recommends using only extension cords with a minimum of 16-gauge wire and a maximum length of 25' (7.6 m).
- Use only a 3-wire extension cord that accepts the plug.
- Ensure the extension is in good condition.
- Ensure the extension cord is rated high enough to carry the current the products draws. An undersized cord causes a drop in line voltage, resulting in loss of power and overheating.

## Circuit Breakers

### WARNING

- *Protect against short circuits by installing a circuit breaker of the proper capacity.*
- *The circuit breaker must be suitably located and easily reached.*
- *The circuit breaker must be marked at the disconnecting device for the equipment with the on/off symbol I/O and meet the relevant requirements of IEC60947-1 and IEC60947-3.*

Table 2 lists the circuit breaker capacity required for the pump. Locate the switch or circuit breaker near the pump.

**Table 2 Circuit Breaker Requirements**

|                  |                         |
|------------------|-------------------------|
| PTS08001INV      | 200-240 VAC<br>50/60 Hz |
| Breaker capacity | 10 Amperes              |

## Mechanical Connections

### Pump Location

Locate the pump on a firm, level surface.

Clearance mounting holes provided in the frame can be used to provide permanent attachment. Attach the provided rubber mounts to the threaded holes at the frame base to minimize vibration.

Do not position the equipment so that it is difficult to operate the disconnecting device.

### Vacuum Pump Isolation Valve

Scroll pumps return to atmospheric pressure quickly when shut off. Thus, an integrated fast acting, automatic, normally closed isolation valve is standard equipment to prevent pump debris from being transported back into the vacuum chamber when the pump is turned off:

- The opening of this valve must occur  $\geq 350$  ms after pump startup.
- Valve closing must occur  $\leq 250$  ms after pump shut off.  
The VPI Valve status is controlled by motor operation. The VPI Valve:
  - Opens when the pump is running and
  - Closes when it is stopped due to switch off or overload protection.

This TriScroll 800 Inverter Scroll pump includes the necessary communication cable, needed to operate the pump.

### Pump Inlet

Use NW40, or larger, clean vacuum hardware with as short a length as practical between the pump inlet and the vacuum chamber.

Use a bellows to provide both vibration isolation and strain relief between the pump and vacuum chamber.

### Pump Exhaust

A female 3/8 National Pipe Thread exhaust fitting is located underneath the scroll module. This fitting swivels 360°. Additionally, a NW25 adapter is provided with a 3/8 National Pipe Thread male thread.

### Gas Ballast

The pump incorporates an automatic gas ballast to prevent water and other condensates from accumulating within the pump. The standard configuration is a sintered filter installed in the 1/4" National Pipe Thread gas ballast port (item 8 on Figure 1 on page 30). This allows enough atmospheric air to enter the pump in order to purge condensates while not affecting pump ultimate pressure or pumping speed.

For applications where the ingress of atmospheric air is undesirable, dry nitrogen at a flow rate of  $\approx 5$  lpm can be bled into the gas ballast port. See item 8 on Figure 1 on page 30. Blocking of the gas ballast port is not recommended.

### Bearing Purge

A 1/4" National Pipe Thread bearing purge port (item 5 on Figure 1 on page 30) protects the main crankshaft bearings in applications where large amounts of water are being pumped. In the standard configuration, this port is sealed.

To enable the bearing purge, dry nitrogen at a flow rate of  $\approx 5$  lpm can be bled into the bearing purge port. This gas supply should be maintained at 2 psig or less and must be kept below 5 psig.

See "Purge Kit" on page 36.

## Operation

In order to reach ultimate vacuum, the pump must be left running for about one hour with the inlet sealed.

Unlike conventional oil-sealed pumps, Agilent's dry scroll pumps do not have fluid to cleanse them of accumulated dust and debris. Run the pump periodically at atmosphere for a minute or two to flush out the pump. Flush the pump regularly and adjust this schedule according to your specific conditions.

### Startup Procedure

1. Be sure that the vacuum system isolation valve is closed.
2. Turn on power to the pump.

#### NOTE

*Pump operation is interlock dependent. The provided mating connector must be plugged into P1 to start the pump.*

#### NOTE

*The integrated VPI solenoid must be connected to I/O P2 with the provided communications cable prior starting the pump.*

#### NOTE

*When the controller is switched on, an electronic self-test occurs, during which an orange LED illuminates for one second and the turns off for two seconds. During this time, the pump will not start.*



*This self-test time is not required if the pump is already powered, as when remotely or serially operated.*

3. Open the isolation valve.

### Shutdown Procedure

1. Close the vacuum system isolation valve. This prevents debris in pump from being transported into the vacuum system.
2. Turn off power to the pump.

# TriScroll 800 Inverter Pump Components View


- ① VPI Plastic Cover Cowling
- ② Cowling Screws; M5 (3) (two shown)
- ③ Inlet (NW40)
- ④ Inlet Screen
- ⑤ Bearing Purge Port (1/4" NPT)
- ⑥ Gas Ballast Port (1/4" NPT)
- ⑦ Mounting Holes; 11 mm diameter thru (8)
- ⑧ Rubber Feet (4 underneath)
- ⑨ Integrated VPI Valve
- ⑩ VPI Solenoid
- ⑪ VPI Communication Cable
- ⑫ Power Connector (IEC 320)
- ⑬ On/Off Switch
- ⑭ Serial Port J1
- ⑮ Interlock Port P1 I/O-1
- ⑯ Interlock Port P2 I/O-2
- ⑰ TPS Connector
- NW25 Exhaust Adapter (Not Shown)

Figure 1 TriScroll 800 Inverter Vacuum Pump


# Technical Information

**Table 3 Specifications**

|                                            |                                                                        |
|--------------------------------------------|------------------------------------------------------------------------|
| Model                                      | TriScroll™ 800 Inverter Dry Scroll Vacuum Pump (PTS08001INV)           |
| Interface dimensions                       | See Figure 2                                                           |
| Peak pumping speed                         | 800 l/m, 48m <sup>3</sup> /hr (28.3 cfm)<br>See Figure 3               |
| Media                                      | Clean air. No toxic, corrosive, explosive or particulate forming gases |
| Ultimate pressure (Torr)                   | 10.0 x 10 <sup>-3</sup> Torr (13.3 x 10 <sup>-3</sup> mbar)            |
| Maximum inlet pressure                     | 1.0 atmosphere (0 psig)                                                |
| Maximum continuous inlet pressure          | 100 Torr                                                               |
| Maximum outlet pressure                    | 1.1 atmosphere (1.5 psig)                                              |
| Inlet connection                           | NW40                                                                   |
| Exhaust connection                         | Female 3/8" National Pipe Thread (NW25 adapter provided)               |
| Gas ballast                                | Female 1/4" National Pipe Thread (40 Micron sintered filter provided)  |
| Ambient operating temperature              | 5 °C to 40 °C (41 °F to 104 °F)                                        |
| Storage temperature                        | -20 °C to 60 °C (-4 °F to 140 °F)                                      |
| Motor rating                               | 1.0 HP (0.75 kW)                                                       |
| Operating voltage                          | 50-60 Hz:200-240 VAC Single-phase ±10%, 1200 VA                        |
| Maximum line current                       | 200 VAC; 6 A<br>240 VAC; 5A                                            |
| Motor Thermal Protection                   | Automatic                                                              |
| Operating speed                            | 35-65 Hz Factory setting: 62 Hz, 1800 RPM                              |
| Cooling system                             | Air-cooled                                                             |
| Weight                                     | Pump only: 32 kg (71 lbs)<br>Shipping weight: 42 kg (92 lbs)           |
| Noise level (per ISO 11201)                | Variable with frequency, 58-71 dB(A)                                   |
| Vibration level at inlet (per ISO 10816-1) | Variable with frequency                                                |
| Main fuse                                  | 12.5 A, 250 V                                                          |
| Installation and storage                   | Installation category 2, indoor usage and storage                      |
| Pollution degree                           | 2                                                                      |
| Altitude                                   | 2000 m                                                                 |

**Figure 2 Interface Drawing with Dimensions**

**Figure 3 Pumping Speed Curves**

# Pump Electrical Controller

The pump is powered by an inverter that converts single phase line power into 3-phase power that drives an induction motor at a user selectable frequency.

Several operational variables are modifiable through the J2 serial port when using Agilent T-Plus software. A number of additional system variables are also viewable.

Pump operation is interlocked through the P1 connector. Pins 4 and 5, and Pins 8 and 9, must be connected together to enable operation. A prewired connector is supplied with the pump and should be plugged into P1. The maximum power output of the inverter is 950 Watts. At pump inlet pressures above approximately 20 Torr, the pump will be power limited and the output frequency will be reduced. This mode is known as *autotuning*.

At pump start, the drive frequency is ramped up over several seconds. This provides for smooth starting and eliminates large in-rush currents.

## Technical Specifications

- Input voltage: 200 V-240 V / 50-60 Hz
- Max frequency: 65 Hz, factory set @ 62Hz
- Normal Operation maximum power: 850 W
- Protection level: IP 20

**Table 4 P1 I/O-1 Interlock**

| PIN # | SIGNAL NAME                                     | IN/OUT |
|-------|-------------------------------------------------|--------|
| 1     | Interlock contact 1                             | Out    |
| 2     | 0-10Vdc, with reference to pin 5 remote control | In     |

**Table 4 P1 I/O-1 Interlock**

| PIN # | SIGNAL NAME                                 | IN/OUT |
|-------|---------------------------------------------|--------|
| 3     | No Connection                               |        |
| 4     | Start/Stop (-): Connect to Pin 5            | In     |
| 5     | Groundout                                   | Out    |
| 6     | Interlock contact 2                         | Out    |
| 7     | No Connection                               |        |
| 8     | Start/Stop (+): Connect to Pin 9            | In     |
| 9     | +24V $\pm$ 5% 50mA, with reference to pin 5 | Out    |


Remote speed setting allows you to control rotational speed using a 0-10 V analog input:

- $V_{in} > 9$  V: Driving frequency will be the full speed (factory set @62 Hz)
- $V_{in} < 1$  V: Driving frequency will be the minimum speed (35 Hz)

### NOTE

*Pump operation is interlock dependent. The provided mating connector must be plugged into P1 to start the pump.*

**Figure 4 Driving Frequency vs. Analog Input**

## Serial Interface Specification RS485 and RS232

**Table 5 J1 - Serial Port**

| PIN # | SIGNAL NAME                            |
|-------|----------------------------------------|
| 1     | +5V ±5%, 50mA, with reference to pin 5 |
| 2     | RS232 TX                               |
| 3     | RS232 RX                               |
| 4     | No Connection                          |
| 5     | GND                                    |
| 6     | RS485 A                                |
| 7     | No Connection                          |
| 8     | RS485 B                                |
| 9     | No Connection                          |

Physical level: RS232 or RS485

Maximum baud: 9600, 8 data bit, no parity, 1 stop bit.

**Table 6 P2 I/O-2 Port on DB9 Male Connector**

| PIN # | SIGNAL NAME                                        |
|-------|----------------------------------------------------|
| 1     | No Connection                                      |
| 2     | No Connection                                      |
| 3     | No Connection                                      |
| 4     | No Connection                                      |
| 5     | GND                                                |
| 6     | No Connection                                      |
| 7     | No Connection                                      |
| 8     | No Connection                                      |
| 9     | VPI output +24V ±5% 400mA, with reference to pin 5 |

**Table 7 T-Plus Software Variables**

| WIN | TYPE      | DESCRIPTION                                                                                                                                                                                        |
|-----|-----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 000 | Logical   | Start/Stop<br>(1=START; 0=STOP)                                                                                                                                                                    |
| 001 | Logical   | Low Speed<br>[0=OFF / 1=ON]                                                                                                                                                                        |
| 008 | Logical   | Remote/Serial Configuration<br>(1=REMOTE; 0=SERIAL)                                                                                                                                                |
| 102 | Numerical | Set Point Value<br>(Hz)                                                                                                                                                                            |
| 105 | Numerical | Set Point Hysteresis (%)<br>[0-100]                                                                                                                                                                |
| 108 | Numerical | Baud Rate<br>(0-4)<br>[600-1200-2400-4200-9600]                                                                                                                                                    |
| 117 | Numerical | Low Speed Adjust<br>(Hz)                                                                                                                                                                           |
| 120 | Numerical | Rotational Frequency Setting<br>[Hz] (High speed adjust)                                                                                                                                           |
| 200 | Numerical | Bus Current<br>[mA]                                                                                                                                                                                |
| 201 | Numerical | 3 Phase Voltage<br>[Vrms]                                                                                                                                                                          |
| 202 | Numerical | Power<br>[W]                                                                                                                                                                                       |
| 203 | Numerical | Driving Frequency<br>[Hz]: (current driving frequency)                                                                                                                                             |
| 205 | Numerical | Status:<br><br><input type="checkbox"/> 0=stop<br><input type="checkbox"/> 2=ramp<br><input type="checkbox"/> 3=autotuning<br><input type="checkbox"/> 5=normal<br><input type="checkbox"/> 6=fail |

**Table 7 T-Plus Software Variables (Continued)**

| WIN | TYPE         | DESCRIPTION                                                                                                                                                                                                                                                                                                                                                                                   |
|-----|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 206 | Numerical    | Error Code:<br><input type="checkbox"/> Bit 7: motor block<br><input type="checkbox"/> Bit 6: short circuit<br><input type="checkbox"/> Bit 5: undervoltage<br><input type="checkbox"/> Bit 4: spare<br><input type="checkbox"/> Bit 3: power fail<br><input type="checkbox"/> Bit 2: controller overt<br><input type="checkbox"/> Bit 1: pump overt<br><input type="checkbox"/> Bit 0: spare |
| 211 | Numerical    | Controller Temperature [°C]                                                                                                                                                                                                                                                                                                                                                                   |
| 216 | Numerical    | Environment Temperature [°C]                                                                                                                                                                                                                                                                                                                                                                  |
| 300 | Numerical    | Cycle Time [min]                                                                                                                                                                                                                                                                                                                                                                              |
| 301 | Numerical    | Cycle Number                                                                                                                                                                                                                                                                                                                                                                                  |
| 302 | Numerical    | Pump Life [h]                                                                                                                                                                                                                                                                                                                                                                                 |
| 319 | Alphanumeric | Controller Model                                                                                                                                                                                                                                                                                                                                                                              |
| 323 | Alphanumeric | Controller Serial Number                                                                                                                                                                                                                                                                                                                                                                      |
| 325 | Alphanumeric | Electrical Modification Level                                                                                                                                                                                                                                                                                                                                                                 |
| 500 | Logical      | Monitor Mode (write only)                                                                                                                                                                                                                                                                                                                                                                     |
| 503 | Numerical    | RS485 Serial Address Setting [0-31]                                                                                                                                                                                                                                                                                                                                                           |
| 504 | Logical      | Serial Type Select (0=RS323; 1=RS485)                                                                                                                                                                                                                                                                                                                                                         |

**Table 8 Status LEDs**

| LED STATUS     | CONTROLLER STATUS  |
|----------------|--------------------|
| Off            | Stop               |
| Green flashing | Ramp or Autotuning |
| Green          | Normal Operation   |
| Red            | Fail               |
| Orange + Off   | Reset + selftest   |

**NOTE**


When the controller is switched on, an electronic self-test occurs, during which an orange LED illuminates for one second and the turns off for two seconds. During this time, the pump will not start.

This self-test time is not required if the pump is already powered, as when remotely or serially operated.


**Using T-Plus Software To Change The Target Frequency**

1. Install T-Plus software on your PC by running setup.exe.
2. Using a 9-pin serial cable, connect your PC to the J2 serial port on the inverter (Figure 5).
3. Start the pump using the On/Off switch. The provided Interlock Connector must be installed in P1.
4. Start the T-Plus software and Figure 5 appears.

**Figure 5 T+ Pumps List**

5. Click **Open Control Panel** and Figure 6 appears.
6. Open the *Speed Settings* directory, and click **Target Frequency**.

**Figure 6 T-Unknown Pump**

7. Enter the desired operating frequency into *Insert New Value*. This must be an integer between 35 and 65.
8. Click **Send**. The pump immediately tries to change to the new target frequency. If the pump requires more than 750 Watts to attain the target frequency, the pump is power limited and will run at a lower frequency (auto-tuning)

mode). This condition continues until the power required drops below 750 Watts.

The target frequency is stored in the inverter even when power to the pump is removed.

## Purge Kit

A purge kit (Agilent part number PTSPURGEKIT) is available to properly purge either the bearing purge or the gas ballast. This kit contains a flowmeter and all necessary valving and tubing.

## Exhaust Silencer Kit

In certain applications it is desired to lower the noise level emitted by the pump. For these cases, Agilent offers Exhaust Silencer Kits (P/N EXSLRTRISCROLL, Manual P/N: 699904380), which breaks up noise pulsation in the discharge line.

The Exhaust Silencer Kits contain two basic components, as well as installation hardware. A cylindrical resonator chamber has been selected to attenuate the frequencies associated with the exhaust gas pulsation noise. The resonator is then plumbed to an exhaust muffler/filter. This muffler provides the majority of the actual noise reduction. It also contains an integral 5-micron paper filter element to capture any tip seal debris. This paper filter is replaceable; the recommended service interval is whenever the pump tip seal or module maintenance is performed.

Replacement polyester filter element rated 99% efficient at the 5-micron particle size, P/N: REPLSLRFILTER,

## Vibration Isolation Kit

A set of vibration isolation mounts (P/N PTSVIBISOKIT, Manual P/N: 699904385) are available for Agilent scroll pumps. These are used in place of the standard feet installed on the pump for applications that require minimal vibration transmission from the pump to the installation. They can be used alone, or with inlet flexible bellows to provide pump isolation.

The vibration isolation mounts in the Vibration Isolation Kit have been designed to ensure the lowest vibration level transmitted at the base of the pump.

Vibration level in every direction (axis), transmitted to the surface the pump sits on, is reduced by a significant amount. Typical vibration reduction levels are 75% for the TriScroll 800 pump.

The TriScroll Vibration Isolation Kit consists of (4) vibration isolation mounts. They have a M6 metric threaded screw to attach to the TriScroll frame; the screw is provided.

These mounts increase the height of the TriScroll pumps by 0.75" compared to the standard mounting feet.

## HEPA Inlet filter

Agilent now offers inlet filters for the scroll pumps, which protect the pump from ingested particles, and keep particles from migrating out of the pump. These HEPA filters provide a rated 99.97% efficiency in trapping particles > 0.3 microns.

NW 40 Inlet Trap with HEPA filter insert: P/N#: SCRINTRPNW40.

Replacement HEPA filter element P/N#: RPLHEPAFILTER1

# Troubleshooting

Table 9 contains a list of possible problems, their probable causes, and corrective actions.

**Table 9 Troubleshooting Chart**

| <b>Problem</b>             | <b>Probable Cause</b>        | <b>Corrective Action</b>                                                                                                                          |
|----------------------------|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Pump won't start           | P1 interlock missing         | Plug in connector (provided).                                                                                                                     |
|                            | TPS cable unplugged          | Plug in TPS cable                                                                                                                                 |
|                            | Wrong input voltage          | Plug inverter into 200-240 VAC 50/60 Hz.                                                                                                          |
|                            | Circuit breaker open         | Close breaker.<br>Identify cause of overload.                                                                                                     |
|                            | Motor thermal protector open | Allow motor to cool.<br>Identify cause of overload.                                                                                               |
|                            | Main fuse open               | Inspect and replace fuse (Type H, 10 A, 250 V (Littlefuse 0216010.HxP or equivalent)).                                                            |
|                            | Wiring loose or cut          | Repair or replace.                                                                                                                                |
|                            | Excessive voltage drop       | Check size and length of power supply cable.                                                                                                      |
|                            | Defective motor              | Inspect. Contact Agilent.                                                                                                                         |
| Poor ultimate pressure     | System leak                  | Locate and repair leak.                                                                                                                           |
|                            | Water in pump                | Flush pump with air or dry nitrogen.                                                                                                              |
|                            | Gas ballast plugged          | Replace breather vent. Contact Agilent.                                                                                                           |
|                            | Solvent in pump              | Flush pump with air or dry nitrogen.<br>Install trap or filter.                                                                                   |
|                            | Seals worn out               | Replace tip seals. (Table 10 and Table 11 on page 38 list maintenance kits and service options.)                                                  |
|                            | Poor conductance to pump     | Replumb with shorter and/or larger diameter tubing.                                                                                               |
| Pump makes hammering noise | Pump overheated              | Check ambient temperature.<br>Check ventilation to pump.                                                                                          |
|                            | Debris in pump               | Check inlet screen.<br>Flush pump.<br>Disassemble pump and inspect. (Table 10 and Table 11 on page 38 list maintenance kits and service options.) |

# Maintenance

## General Information

Agilent TriScroll 800 Inverter pump is designed to provide years of trouble-free service if maintenance procedures and intervals are observed. Bearing grease replenishment and tip seal replacement is recommended when pump base pressure has risen to an unacceptably high level for your application. Bearings, rotary seals and O-rings should also be replaced if the pump exhibits humming or grinding noises from the bearings. Main bearing life may be shortened if your application requires the pumping of high quantities of water vapor. Use of the bearing purge kit (PTSPURGEKIT) mentioned earlier will keep this water from impacting bearing life.

Maintenance should be performed in accordance with procedures, tooling and materials specified in the manuals listed below.

## Related TriScroll Manuals

Other manuals related to tip seal replacement, pump module replacement, and major maintenance of the TriScroll 800 Inverter pump are listed in Table 10.

**Table 10 Other Related Manuals**

| Title                       | Applicable TriScroll Model      | Part Number |
|-----------------------------|---------------------------------|-------------|
| Tip Seal Replacement Manual | All TriScroll 800 Series models | 699904310   |

## Maintenance and Tooling Kits

Material and tooling required to perform maintenance on TriScroll pumps is provided in kit form. A description of each kit and ordering information is provided in Table 11.

**Table 11 Maintenance and Tooling Kits**

| Description              | Contents                                                                      | Applicable TriScroll Model      | Part Number |
|--------------------------|-------------------------------------------------------------------------------|---------------------------------|-------------|
| Tip Seal Tool Kit        | All tools required to change the tip seals on the TriScroll 800 Series pumps. | All TriScroll 800 Series models | PTSTSTKIT   |
| Replacement Tip Seal Set | Replacement tip seals and static O-rings for TriScroll 800 Series pumps.      | All TriScroll 800 Series models | PTSS0800TS  |

### NOTE

*After tip seal replacement, the TriScroll 800 Inverter pump may require up to 100 hours run time to achieve full rotational speed.*

## Factory Service Options

Table 12 lists the lists the service options that Agilent offers for the TriScroll 800 Inverter pump.

**Table 12 Factory Service Options**

| Factory Service Options                                  | Part Number |
|----------------------------------------------------------|-------------|
| Advance Exchange TriScroll 800 Inverter Pump Module Only | EXPTS0800SC |

## Accessories

The accessories listed in Table 13 are available for use with the TriScroll 800 Inverter pump. Contact your local Agilent office to place an order. A list of offices is included on the rear cover of this manual.

**Table 13 Accessories**

|                         |               |
|-------------------------|---------------|
| Purge Kit               | PTSPURGEKIT   |
| Exhaust Extension       | S4807001      |
| Exhaust Filter Kit      | PTS800EXFIL   |
| Isolation Valve Kit     | VPI40INV24DC  |
| Exhaust Silencer Kit    | EXSLRTRISROLL |
| Vibration Isolation Kit | PTSVIBISOKIT  |
| HEPA Inlet Filters      | SCRINTRPNW40  |

## Contacting Agilent

In the United States, you can contact Agilent Customer Service at 1-800-882-7426. See the back cover of this manual for a listing of our sales and service offices.

Visit our web site at: <http://www.chem.agilent.com/en-US/Products/Instruments/vacuum/pages/default.aspx>.

*This page intentionally left blank.*

**Vacuum Products Division  
Instructions for returning products**

Dear Customer:

Please follow these instructions whenever one of our products needs to be returned.

1) Complete the attached Request for Return form and send it to Agilent Technologies (see below), taking particular care to identify all products that have pumped or been exposed to any toxic or hazardous materials.

2) After evaluating the information, Agilent Technologies will provide you with a Return Authorization (RA) number via email or fax, as requested.

**Note:** Depending on the type of return, a Purchase Order may be required at the time the Request for Return is submitted. We will quote any necessary services (evaluation, repair, special cleaning, eg).

**3) Important steps for the shipment of returning product:**

- Remove all accessories from the core product (e.g. inlet screens, vent valves).
- Prior to shipment, drain any oils or other liquids, purge or flush all gasses, and wipe off any excess residue.
- If ordering an Advance Exchange product, **please use the packaging from the Advance Exchange to return the defective product.**
- Seal the product in a plastic bag, and package product carefully to avoid damage in transit. You are responsible for loss or damage in transit.
- Agilent Technologies is not responsible for returning customer provided packaging or containers.
- **Clearly label package with RA number.** Using the shipping label provided will ensure the proper address and RA number are on the package. Packages shipped to Agilent without a RA clearly written on the outside cannot be accepted and will be returned.

4) Return only products for which the RA was issued.

5) **Product being returned under a RA must be received within 15 business days.**

6) **Ship to the location specified on the printable label, which will be sent, along with the RA number, as soon as we have received all of the required information.** Customer is responsible for freight charges on returning product.

7) Return shipments must comply with all applicable **Shipping Regulations** (IATA, DOT, etc.) and carrier requirements.

RETURN THE COMPLETED **REQUEST FOR RETURN** FORM TO YOUR NEAREST LOCATION:

**EUROPE:**

Fax: 00 39 011 9979 330  
Fax Free: 00 800 345 345 00  
Toll Free: 00 800 234 234 00  
[vpt-customer@agilent.com](mailto:vpt-customer@agilent.com)

**NORTH AMERICA:**

Fax: 1 781 860 9252  
Toll Free: 800 882 7426, Option 3  
[vpl-ra@agilent.com](mailto:vpl-ra@agilent.com)

**PACIFIC RIM:**

please visit our website for individual office information  
<http://www.agilent.com>

Please read important policy information on Page 3 that applies to all returns.

**1) CUSTOMER INFORMATION**

|                                      |               |                                                                                                      |  |
|--------------------------------------|---------------|------------------------------------------------------------------------------------------------------|--|
| <b>Company Name:</b>                 |               | <b>Contact Name:</b>                                                                                 |  |
| <b>Tel:</b>                          | <b>Email:</b> | <b>Fax:</b>                                                                                          |  |
| <b>Customer Ship To:</b>             |               | <b>Customer Bill To:</b>                                                                             |  |
|                                      |               |                                                                                                      |  |
|                                      |               |                                                                                                      |  |
| Europe only: <b>VAT reg. Number:</b> |               | USA/Canada only: <input type="checkbox"/> <b>Taxable</b> <input type="checkbox"/> <b>Non-taxable</b> |  |

**2) PRODUCT IDENTIFICATION**

| Product Description | Agilent P/N | Agilent S/N | Original Purchasing Reference |
|---------------------|-------------|-------------|-------------------------------|
|                     |             |             |                               |
|                     |             |             |                               |
|                     |             |             |                               |

**3) TYPE OF RETURN** (Choose one from each row and supply Purchase Order if requesting a billable service)

- 3A.  Non-Billable  Billable **➔** New PO # (hard copy must be submitted with this form):
- 3B.  Exchange  Repair  Upgrade  Consignment/Demo  Calibration  Evaluation  Return for Credit

**4) HEALTH and SAFETY CERTIFICATION**

**AGILENT TECHNOLOGIES CANNOT ACCEPT ANY PRODUCTS CONTAMINATED WITH BIOLOGICAL OR EXPLOSIVE HAZARDS, RADIOACTIVE MATERIAL, OR MERCURY AT ITS FACILITY.**

**Call Agilent Technologies to discuss alternatives if this requirement presents a problem.**

**The equipment listed above (check one):**

- HAS NOT** pumped or been exposed to any toxic or hazardous materials. OR
- HAS** pumped or been exposed to the following toxic or hazardous materials. If this box is checked, the following information must also be filled out. Check boxes for all materials to which product(s) pumped or was exposed:

- Toxic  Corrosive  Reactive  Flammable  Explosive  Biological  Radioactive

**List all toxic/hazardous materials. Include product name, chemical name, and chemical symbol or formula:**

**NOTE:** If a product is received at Agilent which is contaminated with a toxic or hazardous material that was not disclosed, **the customer will be held responsible** for all costs incurred to ensure the safe handling of the product, and **is liable** for any harm or injury to Agilent employees as well as to any third party occurring as a result of exposure to toxic or hazardous materials present in the product.

**Print Name:** \_\_\_\_\_ **Authorized Signature:** ..... **Date:** \_\_\_\_\_

**5) FAILURE INFORMATION:**

|                                                                                |
|--------------------------------------------------------------------------------|
| Failure Mode (REQUIRED FIELD. See next page for suggestions of failure terms): |
| Detailed Description of Malfunction: (Please provide the error message)        |
| Application (system and model):                                                |

**I understand and agree to the terms of Section 6, Page 3/3.**

**Print Name:** \_\_\_\_\_ **Authorized Signature:** ..... **Date:** \_\_\_\_\_

**Vacuum Products Division  
Request for Return Form  
(Health and Safety Certification)**

Please use these Failure Mode to describe the concern about the product on Page 2.

**TURBO PUMPS and TURBO CONTROLLERS**

| APPARENT DEFECT/MALFUNCTION                                                                                                                                                                   | POSITION                                                      | PARAMETERS                                                                                                                                                                                                   |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| - Does not start<br>- Does not spin freely<br>- Does not reach full speed<br>- Mechanical Contact<br>- Cooling defective<br>- Noise<br>- Vibrations<br>-Leak<br>-Overtemperature<br>-Clogging | - Vertical<br>-Horizontal<br>-Upside-down<br>-Other:<br>..... | Power:                      Rotational Speed:<br>Current:                    Inlet Pressure:<br>Temp 1:                     Foreline Pressure:<br>Temp 2:                     Purge flow:<br>OPERATING TIME: |

**ION PUMPS/CONTROLLERS**

|                         |                        |
|-------------------------|------------------------|
| - Bad feedthrough       | - Poor vacuum          |
| - Vacuum leak           | - High voltage problem |
| - Error code on display | - Other                |

**VALVES/COMPONENTS**

|                        |                  |
|------------------------|------------------|
| - Main seal leak       | - Bellows leak   |
| - Solenoid failure     | - Damaged flange |
| - Damaged sealing area | -Other           |

**LEAK DETECTORS**

|                          |                          |
|--------------------------|--------------------------|
| - Cannot calibrate       | -No zero/high background |
| - Vacuum system unstable | - Cannot reach test mode |
| - Failed to start        | - Other                  |

**INSTRUMENTS**

|                          |                     |
|--------------------------|---------------------|
| - Gauge tube not working | - Display problem   |
| - Communication failure  | - Degas not working |
| - Error code on display  | - Other             |

**SCROLL AND ROTARY VANE PUMPS**

|                        |                         |
|------------------------|-------------------------|
| - Pump doesn't start   | - Noisy pump (describe) |
| - Doesn't reach vacuum | - Over temperature      |
| - Pump seized          | - Other                 |

**DIFFUSION PUMPS**

|                        |                       |
|------------------------|-----------------------|
| - Heater failure       | - Electrical problem  |
| - Doesn't reach vacuum | - Cooling coil damage |
| - Vacuum leak          | - Other               |

Section 6) **ADDITIONAL TERMS**

**Please read the terms and conditions below as they apply to all returns and are in addition to the Agilent Technologies Vacuum Product Division – Products and Services Terms of Sale.**

- Customer is responsible for the freight charges for the returning product. Return shipments must comply with all applicable **Shipping Regulations** (IATA, DOT, etc.) and carrier requirements.
- Customers receiving an Advance Exchange product agree to return the defective, rebuildable part to Agilent Technologies **within 15 business days**. Failure to do so, or returning a non-rebuildable part (crashed), will result in an invoice for the non-returned/non-rebuildable part.
- Returns for credit toward the purchase of new or refurbished Products are subject to prior Agilent approval and may incur a restocking fee. Please reference the original purchase order number.
- Units returned for evaluation will be evaluated, and a quote for repair will be issued. If you choose to have the unit repaired, the cost of the evaluation will be deducted from the final repair pricing. A Purchase Order for the final repair price should be issued within 3 weeks of quotation date. Units without a Purchase Order for repair will be returned to the customer, and the evaluation fee will be invoiced.
- A Special Cleaning fee will apply to all exposed products per Section 4 of this document.
- If requesting a calibration service, units must be functionally capable of being calibrated.

*This page intentionally left blank.*

# Service & Support

## North America

Agilent Technologies  
Vacuum Products Division  
121 Hartwell Avenue  
Lexington, MA 02421 USA  
Tel.: +1 781 861 7200  
Toll-Free: +1 800 882 7426  
Fax: +1 781 860 5437

## Benelux

Agilent Technologies Netherlands B.V.  
Vacuum Products Division  
Herculesweg 8  
4338 PL Middelburg  
The Netherlands  
Tel: +31 118 671570  
Fax: +31 118 671569  
Toll free: 00 800 234 234 00

## Canada

Central coordination through  
Agilent Technologies  
Vacuum Products Division  
121 Hartwell Avenue  
Lexington, MA 02421 USA  
Tel.: +1 781 861 7200  
Toll-Free: +1 800 882 7426  
Fax: +1 781 860 5437

## China

Agilent Technologies (China) Co. Ltd  
Vacuum Products Division  
No.3, Wang Jing Bei Lu, Chao Yang  
District  
Beijing, 100102 China  
Tel.: +86 (0)10 64397888  
Fax: +86 (0)10 64391318  
Toll-Free: 800 820 3278

## France

Agilent Technologies France  
Vacuum Products Division  
3 avenue du Canada, Parc Technopolis  
Z.A. Courtaboeuf - CS 90263  
91978 Les Ulis cedex  
France  
Tel.: +33 (0)1 64 53 61 15  
Fax: +33 (0)1 64 53 50 01  
Toll free: 00 800 234 234 00  
e-Mail: [vpf.sales@agilent.com](mailto:vpf.sales@agilent.com)

## Germany & Austria

Agilent Technologies  
Sales & Services GmbH & Co. KG  
Vacuum Products Division  
Lyoner Str. 20  
60 528 Frankfurt am Main  
Germany  
Tel.: +49 69 6773 43 2230  
Fax: +49 69 6773 43 2250

## India

Agilent Technologies India Pvt. Ltd.  
Vacuum Products Division  
G01. Prime corporate Park,  
230/231, Sahar Road,  
Opp. Blue Dart Centre,  
Andheri (East), Mumbai – 400 099. India  
Tel: +91 22 30648287/8200  
Fax: +91 22 30648250  
Toll Free: 1800 113037

## Italy

Agilent Technologies Italia S.p.A.  
Vacuum Products Division  
Via F.lli Varian 54  
10040 Leini, (Torino) ITALY  
Tel.: +39 011 997 9111  
Fax: +39 011 997 9350  
Toll-Free: 00 800 234 234 00

## Japan

Agilent Technologies Japan, Ltd.  
Vacuum Products Division  
8th Floor, Sumitomo Shibaura Building  
4-16-36 Shibaura Minato-ku  
Tokyo 108-0023 JAPAN  
Tel.: +81 3 5232 1253  
Toll-Free: 0120 655 040  
Fax: +81 3 5232 1710

## Korea

Agilent Technologies  
Vacuum Products Division  
Shinsa 2nd Bldg. 2F  
966-5 Daechi-dong  
Kangnam-gu, Seoul KOREA 135-280  
Tel.: +82 (0)2 3452 2455  
Toll-Free: 080 222 2452  
Fax: +82 (0)2 3452 3947

## Singapore

Agilent Technologies Singapore Pte. Ltd  
Vacuum Products Division  
Agilent Technologies Building,  
No.1 Yishun Avenue 7  
Singapore 768923  
Tel: +65 6215 8045  
Fax : +65 6754 0574

## Southeast Asia

Agilent Technologies Sales Sdn Bhd  
Vacuum Products Division  
Unit 201, Level 2 uptown 2,  
2 Jalan SS21/37, Damansara Uptown  
47400 Petaling Jaya,  
Selangor, Malaysia  
Tel : +603 7712 6106  
Fax: +603 6733 8121

## Taiwan

Agilent Technologies Taiwan Limited  
Vacuum Products Division (3F)  
20 Kao-Shuang Rd.,  
Pin-Chen City, 32450  
Taoyuan Hsien, Taiwan, R.O.C.  
Tel. +886 34959281  
Toll Free: 0800 051 342

## UK & Ireland

Agilent Technologies  
Vacuum Products Division  
6 Mead Road  
Oxford Industrial Park  
Yarnton, Oxford OX5 1QU UK  
Tel.: +44 (0) 1865 291570  
Fax: +44 (0) 1865 291571  
Toll free: 00 800 234 234 00

## Other Countries

Agilent Technologies Italia S.p.A.  
Vacuum Products Division  
Via F.lli Varian 54  
10040 Leini, (Torino) ITALY  
Tel.: +39 011 997 9111  
Fax: +39 011 997 9350  
Toll-Free: 00 800 234 234 00

## Customer Service and Support

North America  
Toll Free: 800 882 7426, Option 3  
[vpf-ra@agilent.com](mailto:vpf-ra@agilent.com)

Europe:  
Toll Free: 00 800 234 234 00  
[vpt-customer@agilent.com](mailto:vpt-customer@agilent.com)

Pacific Rim:  
Please visit our website for individual  
office information <http://www.agilent.com>

## Worldwide Web Site, Catalog and Order On-line:

[www.agilent.com](http://www.agilent.com)  
Representatives in most countries

# Agilent Technologies

This information is subject to change without notice.

© Agilent Technologies, Inc., 2021