Tesmed+

www.tesmed.com

ITALIANO

ENGLISH

FRANÇAIS

DEUTSCH

ESPAÑOL

Max 830

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Instruction manual

Contents

Symbols used	33
Warnings	34
Precautions.....	34
Possible side effects	34
Operating instructions	35
Description of the device	35
Electrode connection.....	35
Switching on.....	36
Program Selection: total 115 programs	36
Customization of programs from 1 to 20.....	36
• Using programs from 1 to 20 with parameter customization	36
• Parameter adjustment	36
• Restore factory settings for a customisable programme from 01 to 20 .	37
Channel Selection: total 4.....	37
Intensity regulation: 99 levels	37
Immediate intensity reset-program interruption	37
Error detection for disconnected electrodes or cables, only above level 4 intensity	38
Key lock function	38
Pause a programme.....	38
Backlight: adjusting brightness	38
Beep: deactivating the beep	39
Switching off the device.....	39
Disconnecting the electrodes	39
Accidental activation block for travel.....	39
Battery Charging	40
Replacing the battery	41
Electro stimulation	42
High technology.....	43
Contraindications	44
Maintenance.....	45

Disposal	45
How to choose the right treatment to suit individual needs	46
Electrode polarization.....	46
Content of the package	46
Technical characteristics of the TESMED MAX 830	47
Personalised programs from 01 to 20	48
Tesmed programs from 21 to 99	49
FITNESS Treatments.....	49
BEAUTY Treatments.....	50
SPORT Treatments	53
MASSAGES TENS** Treatments	56
TESMED VITALITY 2023 Treatments	57
Sequential Impulse Stimulation Waims system	138
Positioning of electrodes.....	139-148
Normative references.....	149
Contacts, assistance, guarantee, useful information	150

symbols used

Important note, if not respected, there may be a risk for the user and/or the equipment may be damaged. We recommend that you read the information in this manual carefully and keep it in a safe place.

Low frequency parts applied.

Internally fed equipment (DC).

Parts that must be disposed of as special waste.

Read the operator's manual

Manufacturer

Double electrical isolation

Protection against dust and liquid penetration

** These are massages without medical indications. The massage has a soothing and relaxing effect on the treated area.

warnings

- You must take into account the technical characteristics of each program and read the relevant paragraphs in the User Manual.
- The Tesmed Max 830 must not be used for purposes other than those indicated.
- You must begin each treatment with the current emission value set at zero (min channel power setting) and then gradually increase and/or reduce the power to suit your tollerability.
- Max power does not necessarily mean better performance or results.
- We recommend using original electrodes only to guarantee top performance, as these have been developed specifically for the Tesmed Max 830. Non original electrodes could burn the skin.
- The Tesmed Max 830 must be kept out of the reach of children.
- Electrocardiogram monitoring equipment may be disturbed by the current used for the electro stimulator.
- We do not recommend using the Tesmed Max 830 on children under the age of 14 or the elderly over 75.
- The electrodes must not be placed less than 10 cm from the eyeballs.
- Do not use the equipment if there are special health problems without first obtaining medical advice from a doctor.
- Should the equipment malfunction, do not try to open or repair the equipment, but send the product in its original packing to **Feldi S.r.l.**, describing the problems found on the inside of the package (see warranty conditions).
- The bipolar cables not connected to the body for the treatment must be disconnected from the device
- Device not suitable for simultaneous use with high frequency surgical devices
- The application of electrodes close to the chest may increase the risk of heart fever 
- The electrodes should never be placed on the carotid breast, along the head, directly on the eyes, covering their mouths and / or on a spastic muscle.
- Device not suitable for use in the presence of flammable anaesthetic mixed with air or oxygen or nitrous oxide.

Precautions

- Do not tamper with the equipment.
- Keep the electro stimulator away from sources of heat, powder, sprayed water, humid atmospheres and avoid direct exposure to the sun.
- Avoid violent knocks that could damage the equipment.
Operating temperature: 0 °C 45 °C
Storage temperature for a period of 1 month: -20 °C 45 °C
Storage temperature for a period of 6 months: -10 °C 35 °C
Humidity: between 20 and 90% RH (without condensation)
Atmospheric pressure: between 75 and 106 kPa


Possible side effects

- Isolated cases of skin irritation may occur in people with highly sensitive skin. Should there be an allergic reaction to the electrode gel, stop the treatment and consult a specialist doctor. It is normal for the skin to be slightly red after the treatment in the areas where the electrodes were placed: this will disappear after roughly 30 minutes. We strongly advise against treating any areas of the body with serious skin problems.

OPERATING INSTRUCTIONS


Description of the device

- 1 Output channels
- 2 Tesmed logo
- 3 Keypad lock and travel function symbol
- 4 Programme number selected
- 5 Impulse frequency
- 6 Channel number
- 7 Intensity of individual channel
- 8 Function button
- 9 Decrease button
- 10 Increase button
- 11 ON/OFF button, channel selector, rapid zero intensity and other functions
- 12 Intensity symbol
- 13 Impulse amplitude
- 14 Programme time
- 15 TIME symbol
- 16 Battery symbol
- 17 Battery charging socket

Electrode connection

The electrical impulses generated by the Tesmed Max 830 are transmitted to the muscles via self-adhesive electrodes. Before switching on the Tesmed Max 830 or, in any case, before activating intensity, always apply the electrodes to the area to be treated and use the cables to connect them to the device. Insert the plugs on the two-pole cables into each pair of electrodes and connect the other end of each cable to the device by inserting the connector into its dedicated socket on the upper part of the device (**fig. 1**), in accordance with the colour coding. If this procedure is not followed, on exceeding intensity level 5, the device will report an error. Careful study and the TMS "Tesmed Multiplexer System" used in the Tesmed Max have led to the creation of a device that is capable of operating correctly without having to apply electrodes with a defined polarity, thus facilitating application, connection and use.

The plugs on the cables can be connected to either electrode as there is no difference in polarity

Switching on

Briefly press the button to switch on. The device will automatically pass to programme selection. Each time the device is switched on, it will be set to the programme used in the previous session.

Program Selection: total 115 programs

Customizable programs from 1 to 20

Pre-set programs from 21 to 99 and from A1 to B7

(N.B.: Programs from A1 to B7 are located after program 99)

When the program number is flashing, press buttons to select the desired program. To speed up the increment and decrement of the selection, hold down either button . Press button to confirm and proceed to intensity.

Attention: After 15 seconds of inactivity, the program number will become fixed, and you will automatically switch to the intensity function. In this case, if you want to go back to the program selection, press button .

Customization of programs from 1 to 20

Using programs without modifying the parameters:

1. Select the desired program with buttons
2. Press button to directly move to intensity.

Using programs from 1 to 20 with parameter customization:

Preamble:

The parameter to be modified will flash; if no operation is performed for 15 seconds, the device will automatically switch to intensity, ready to start treatment. If you wish to resume parameter adjustment mode, press button again.

Remember: Press button to switch between adjustments.

Parameter adjustment:

After selecting the program, proceed with the following operations in sequence:

1. Press button to confirm the program and proceed to adjust the program duration,
2. The default duration value will flash, set the desired duration from 1 to 99 minutes using buttons , and confirm with button ,
3. The default frequency value will flash, adjustable from 1 to 125, adjust it using buttons and confirm with button .

4. The default pulse amplitude value will flash, adjustable from 40 to 350, adjust it using buttons  , and confirm with button .
 5. The device will switch to intensity.
- Now you can start the treatment by increasing the intensity with button .

Frequency indicates how many cycles of pulses occur in one second and is measured in hertz. Pulse amplitude is the duration of the pulse and is measured in microseconds.

Restore factory settings for a customisable programme from 01 to 20

To restore the factory settings for a programme, select it and then press and hold the buttons  and  for two seconds. The programme will return to the factory settings.

Channel Selection: total 4

The selection of channels to use can only be done at intensity 00 and only when the intensity symbol  appears on the display.

By default, all 4 channels are selected. If, based on the chosen treatment or preferences, you want to use fewer than 4 channels, press button  successively to select the quantity of channels to use.

Intensity regulation: 99 levels

The channels are only independent in terms of intensity setting.

The term intensity refers to the output value of the current applied to the treated area and therefore the level of stimulation with which the muscle contracts.

When the lower part of the display shows the intensity symbol , it is possible to use the buttons   to adjust the intensity level simultaneously on all the previously selected channels.

When the intensity begins to increase, the timer icon  begins to flash, indicating that the device is supplying energy. To increase or decrease the intensity more rapidly, press and hold one of the two buttons  .

To adjust the intensity of an individual channel, press the button  to select the channel for which intensity is to be modified, and then adjust the intensity with the buttons  .

In order to obtain the desired level of muscular contraction, the intensity level can vary in accordance with multiple factors, including the position of the electrodes and the presence of adipose tissue. Therefore, the **same current intensity level may result in different sensations** for different people, on different days, or on the left and the right of the body.

Immediate intensity reset-program interruption

If for any reason it is necessary to interrupt a programme and rapidly zero the intensity, **press and hold for one second the button** , the programme will be interrupted, the intensity will drop to zero and a **beep** sound will be played. If the button is pressed again to increase intensity, the programme will recommence from the beginning.

Error detection for disconnected electrodes or cables, only above level 4 intensity

Tesmed Max 830 will detect an error in the event of a disconnected electrode. If a disconnected electrode is detected, the device will play one or two beep sounds, intensity will be reduced to zero and the number of the channel reporting the error will flash for a few seconds. If intensity is increased again in this condition, no further beep will be played but the number of the channel reporting the error will flash for a few seconds.

The error message may mean that:

- There are no electrodes connected to the channel reporting the error.
- The electrodes are not perfectly attached to the skin; reposition the electrodes, and if they are worn, try new electrodes.
- The cable to which the electrodes are connected is defective.

Check the integrity of the cables and that the electrodes are properly attached to the body.

An error can be detected if two cables are inverted, for example:

if the channel detecting the error is **CH 1**, then the relative cable is the blue cable positioned in line with the blue triangle , disconnect this from the device and connect it to the channel **CH 4**, yellow triangle , without changing the position of the 2 electrodes it is connected to. Now increase intensity and if the blue cable or the electrodes connected to it are defective, **CH 4** will detect an error and flash.

Key lock function

By simultaneously pressing the buttons  and , the display will show the key symbol . Follow the same procedure to unlock the keys. When activated, this function disables all the buttons in order to avoid involuntary modification of parameters during a treatment..

Pause a programme

To pause a programme, for example to receive a telephone call or to do something else, press and hold the button , for two seconds, the display will flash, and the device will play a beep sound. The programme timer will stop.


Press and hold the button  for 2 seconds to recommence the treatment from the point it was paused. Two beeps will be played, and the display will stop flashing. The device will switch off after 30 minutes in pause.

Backlight: adjusting brightness

To reduce the level of the display backlight, press and hold the buttons  and  simultaneously. To increase the level of the display backlight, press and hold the buttons  and  simultaneously. The backlight will remain active for 20 seconds after the most recent pressing of the buttons, after which time it will switch off. It will come on again whenever a button is pressed.

Beep: deactivating the beep

To deactivate the beep, press and hold the buttons   simultaneously for 2 seconds, the display will show the message “BEEP OFF” (**fig. 3**). Repeat the operation to reactivate the beep; the display will show the message “BEEP ON” (**fig. 4**) and the device will play a beep sound.


Switching off the device

To switch off, press and hold the button for 3 seconds  whatever the state of the device. Tesmed Max 830 will switch off automatically if left unused for more than 5 minutes in order to save the battery.

Disconnecting the electrodes

Once treatment is complete, disconnect the plugs from the cables as shown in **figure 2**, taking care not to damage the cables.

Accidental activation block for travel

With the device off, press and hold both the buttons  and  for 2.5 seconds. The display will show the Tesmed logo, the key symbol and the remaining battery   . This procedure is recommended when the device is kept in a suitcase or bag to avoid accidental activation.

To unlock, briefly press the buttons  and  again.

When this function is active, a small amount of energy is consumed which, after a few days, may completely drain the battery. It is therefore recommended to check the remaining charge every now and then via the icon  which is always visible.

Battery Charging

Take the supplied cable and insert the end with the plug into the device (as shown in **Fig. 5**) and the other end with the USB port into a mobile phone charger or a PC (as indicated in **Fig. 6**).

When charging, the display backlight will come on for 5 seconds and then switch off. The display will show the Tesmed logo and the rightmost bar on the battery symbol will flash . This will remain steady when charging is complete.

While the device is charging, pressing the button  will switch on the backlight for 4 seconds. The full battery charge indicator (**fig. 7**) will remain visible until the device is disconnected from the battery charger.

If the battery is charging, the rightmost bar will flash .

- * If the battery is fully charged, the icon will remain visible with all the bars fixed, and the display will show the message "FULL BATT" (**fig. 7**).
- * If there is a problem with the battery (e.g., disconnected or impossible to charge), the entire battery icon will flash.
- * If the battery is disconnected, the message "NO BATT" will be displayed (**fig. 8**).
- * When the device reports a battery error, a single beep will be played.
- * When the battery is fully charged, three rapid beeps will be played.

The lifespan of the battery depends on use and storage method. In the event of a long period of inactivity, charge the battery at least once a year. The battery is to be considered as a replaceable part subject to wear and can be purchased as a spare part.

⚠ Never charge the battery while cables and electrodes are connected to a person.




Replacing the battery

Tesmed Max 830 is powered with a rechargeable lithium polymer battery located to the rear of the device.

Tesmed Max is one of the few electrostimulation devices with a lithium battery that can be substituted directly by the client. The useful lifespan of the battery depends on a number of factors, including the types of programme used, the number of charge cycles carried out, periods of inactivity, temperature, etc. If the device is unused for a long period without charging, the battery may need to be replaced. To access the battery bay, simply remove the hook for the removable belt and push the cover in the direction indicated by the arrow (**fig. 9**).

The battery is connected to the device via a connector; to change the battery, remove it from the battery bay and disconnect the connector from the device (**fig. 10**).

Pay special attention to connection when inserting the new battery; on one side, the connector has a flat surface that must be facing upwards when inserted (**fig. 11**).

electro stimulation

Electro stimulation is a technique that provokes muscular contractions very similar to voluntary contraction thanks to electrical impulses acting on the motor points of the muscles. As everyone knows, muscles contract when they receive an impulse coming from the nervous system and then return to their normal length once they are no longer stimulated thanks to their elasticity. A motor nerve fibre splits into several ramifications before entering the muscle and ends at the motor end plate. The nerve fibre and the group of muscular fibres innervated by this make up the motor unit. A muscular fibre responds to the electrical stimulus only when the intensity of the stimulus reaches a threshold value, which is specific to each fibre. The most obvious use of the electro stimulation of innervated muscles is in kinesiotherapy of atrophy, functional reanimation after a long period of immobility, in vasodilation, pain reducing and in stimulating the treated tissues to combat cellulite (PASSIVE EXERCISE). When treating a patient with electro stimulation, the following are crucial: the correct positioning of the electrodes, the stimulation program and the wave type.

Electrodes

The positioning of the electrodes is crucial. It must provide accurate stimulus parameters for the zone being treated since the conductivity of body tissues varies according to their structure, aqueous and ionic make-up. The electrodes must be in perfect contact with a thoroughly clean skin (also depilated, if necessary). This ensures better penetrability and permeability of the skin when receiving the stimulus to suit the application to be carried out. **For better electrode conductivity we recommend they be replaced when they lose adherence. They can also be cleaned with running water or wiped with a damp cloth; leave to dry before using them again for a few more applications.** NB: With electro stimulation, the current density diminishes as the tissue consistency increases.

We recommend using Tesmed electrodes exclusively, as electrodes of different sizes and composition could cause burns.

Tiredness

The high percentage of simultaneous contraction/relaxing of the muscles due to electrical stimulus often results in a drop in the relaxation of the neurotransmitters, thus resulting in muscle fatigue. The scientific explanation is as follows: the stimulus tends to provoke the recruitment of large diameter motor neurones (these have a low threshold), thus exciting the more consistent and faster muscle fibres (FIBRES 2 or white fibres). These tire more easily, but they use less effort than FIBRES 1 or red fibres. **In other words, the higher the frequency of stimulus, the more tired the muscle becomes.** For this reason, it is important to include relaxation programs and get medical advice. Having said this however, Tesmed's sports programs already contain a relaxation phase. **During treatments using the electro stimulator, every increase in intensity must be done gradually. If there is EXCESSIVE contraction, you must reduce the intensity.** In fact, stimulation must not provoke muscular tetany (i.e. long continuous contractions). **There are no standard rules for the number of contractions and the duration and frequency of the sessions. These must be set to suit each person and the condition of the tissues to be treated.** When treating areas affected by CELLULITE, you should normally start the treatment by gradually raising the intensity until the patient starts to feel a slight tingling, but without creating excessive contraction of the muscles or an unpleasant sensation of localised heat.

To sum up: the vasomotor and metabolic effects induced by the passage of an electrical current considerably increase the consistency and quality of tissues. The best level of stimulus is when the patient perceives the current passing through the tissue; this should provoke a light initial contraction of the muscles (a slight trembling), but never too intense.

There are two important advantages here:

- 1) you avoid the risk of muscular over-contraction and the accumulation of lactic acid;
- 2) the frequency of muscular contractions may be much higher than with traditional gym, to get best results faster.

Advantages of using the electro stimulator

In beauty: improvement of capillary micro circulation with more oxygen supplied to tissues, thus producing better fat mobilisation; toning up and firming up of tissues; draining of excess liquids; an increase in metabolism and cellular activities; a reduction of localised fatty tissues and cellulite.

In sports: isolation of the muscular group being treated; saving of physical and psychological energy; recruitment of a larger number of muscular fibres than is the case with voluntary contraction; training of under-used muscular fibres during normal work-out sessions; shorter recovery periods; protection of the body's system of tendons and muscles; improved blood circulation. Muscular damage due to electro stimulation is most unlikely, as long as the intensity is increased gradually. Any interference in the menstrual cycle (e.g. early, late, heavy or light periods) are extremely subjective and variable. We recommend, therefore, to avoid treatments in the abdominal area during the actual period and the days immediately before and after this. To date, no side effects have been noted concerning the use of electro stimulation during lactation. However, we recommend not treating the upper body area during this time.

We strongly advise against treating any parts of the body with serious skin problems.

Results

All results in the aesthetics field are naturally subjective. We can, in any case, sustain that in the case of TONE UP, good results can be seen after just 15 days with 5-7 regular and constant weekly sessions; cellulite treatment will require 30-40 days. Better and quicker results may be obtained if treatments are added to suitable physical activity and a correct lifestyle.

high technology

Tesmed Max 830 is designed and built with cutting-edge electronic components, a high storage microcontroller with great memory, processing and storage capacity. This means that the Tesmed Max 830 is an extremely professional product with superior performance, yet compact size.

The Tesmed Max 830 programs include a sequential pulsed treatment. Thanks to the new exclusive WAIMS-system, the Tesmed Max 830 electro stimulator offers a valid alternative to manual massage (lymphatic drainage). This stimulation is an exclusive feature of the Tesmed Max 830 that, thanks to this draining, helps blood circulation remarkably, since the extracellular liquids and built-up fats are moved and can thus be removed by the body. This type of massage leaves your legs and arms feeling lighter and gives the skin an even look, at the same time considerably reducing the "orange peel" effect.

All areas treated with the WAIMSSYSTEM gain better tone and are firmer.

contraindications

ATTENTION:

Electro stimulation is not recommended in areas with any of the following:

- inflammation, open wounds, strained or torn muscles, dermatitis, bacterial infections, eritema or skin irritation;
- muscular hernias;
- fractured bones;
- prosthesis or osteosynthesis;
- varicose veins or vascular disease in general;
- neoplasia;
- mucous membranes.
- in the thoracic area to people with a cardiac arrhythmia
- patients with metal joint prosthesis and/or osteosynthesis and/or intratissue/endotissue metal;

This device is prohibited in the following cases:

- pregnant women;
- growing children (risk of juvenile heart attack);
- patients wearing pace-makers (risk of interference with the functioning of the pace-maker);
- implanted defibrillators;
- alteration of the peripheral nerves, muscular disease and serious damage to the joints;
- serious heart problems;
- patients suffering from epilepsy;
- patients with bone tumours;
- patients with serious high blood pressure;
- patients suffering from serious systematic, debilitating or dermatological disease;
- patients undergoing treatment with an HF device, as this could cause burns and damage to the stimulator itself;
- near any short wave or microwave therapy equipment, as this could cause instability in the stimulator output;
- unconscious patients or patients unable to communicate;
- people suffering from psychophysical problems that prohibit their exposure to external factors;
- in the presence of a mix of flammable anaesthetic with air, oxygen or nitrogen protoxide.

Never place the electrodes on the carotid artery and/or on a spastic muscle.

Use of the TESMED MAX 830 for therapy purposes should be supervised and recommended by a doctor. Interrupt application immediately in the event of discomfort and obtain medical advice.

maintenance

The Tesmed Max 830 is a highly sophisticated product and must be handled with care. The maintenance tips below will help you ensure its long life.

- Do not allow it to get wet. Rain, humidity and liquids or condensation generally contain corrosive mineral substances that can damage the electronic circuits.
- Keep the Tesmed Max 830 out of the reach of children.
- Do not leave the Tesmed Max 830 in particularly hot rooms. Excessive temperatures can affect the life of the electronic circuits, damage the batteries and deform the plastic parts.
- Do not leave the Tesmed Max 830 in particularly cold rooms, as condensation could form when it heats up during use and so damage the electrical circuits.
- Do not open the Tesmed Max 830. Any tampering with the equipment by inexperienced people could damage it.
- Do not allow it to fall, be knocked or shaken, as the electronic circuits inside could be irreparably damaged.
- Do not use the battery for purposes other than those intended.
- Never use the battery charger or batteries if damaged.
- Do not short circuit the battery.
- Use only the battery charger provided. Contact Customer Service in the event of a malfunction.
- Never leave the battery in hot or cold places, such as a closed car in the height of Summer or during the Winter, as this will reduce its capacity and life.
- Never throw used batteries in the fire!
- Follow local regulations to dispose of used batteries, e.g. recycle them. Never throw them away with household waste.
- To clean the Tesmed Max 830 (with the equipment switched off), use a soft cloth and alcohol, but no solvents, as these could ruin the plastic parts.

Use only original spare-parts and products for maintenance.

Follow these maintenance times to keep the Tesmed Max 830 in top working condition:

- A)** Check the state of the electrode cables and check for signs of oxidation (monthly if used regularly; fortnightly if used intensely).
- B)** Replace the electricity conducting electrode as soon as these show signs of wear. The manufacturer's warranty is invalidated if the above maintenance is not done. The manufacturer guarantees a 5 year working life of the equipment provided it is used and maintained correctly. The equipment, batteries and its accessories must be disposed of properly (treat as special waste) to protect the environment.

disposal

The equipment, batteries and its accessories must be disposed of properly (treat as special waste) to protect the environment.

how to choose the right treatment to suit individual needs

- You can choose the most suitable treatment from those in the table of treatments below, consulting the list of photos for details of how to position the electrodes on the muscles in question.
- The photo number is shown next to each treatment, in the "Photo" column. If more than one photo is indicated for the same treatment (e.g. **1, 37, 40-44**), select photo **1, 37, 40, 41, 42, 43** or **44** as required.
- We recommend that you do not treat a group of muscles with more than one session a day.
- We recommend that you do not exceed 2 hours of electro stimulation per day, even if involving different groups of muscles.
- **The treatment times shown in the tables is not binding: actual treatment times may vary considerably from one person to another.**
- If in doubt about which treatment to use, we advise you consult a doctor or a trainer.

NB: The programs in the **Tesmed Max 830** are modulates to provide various stages* of massage (1 to max 5). Some programs have one or more Waims-system stages. During treatment (and especially with the Waims-system) we always recommend that you adapt the intensity of each channel, as your sensitivity will vary various from one group of muscles to another.

* The number of stages in a program corresponds to the number of wave types used in it.

electrode polarization

Most commonly available electro stimulators adopt a common electrode standard: they either contain continuous current (DC) or generate modified waves (different positive/negative values, with more of one or the other). Hence the need to specify the polarity. After careful study and with the use of TMS "Tesmed Multiplexer used System" technology in the Tesmed Max 830 means that this equipment can work correctly without having to give a specific polarity to the the electrodes, thus making them easier for you to apply, connect and use.

Content of the package:

- 1 Stimulator
- 1 USB cable
- 4 Bipolar cables
- 20 Electrodes
- 1 User manual

TECHNICAL CHARACTERISTICS OF THE TESMED MAX 830

Power supply	Lithium-polymer (LiPo) battery rechargeable; 3.7 V 1000 mAh
Battery charger	5 V \pm 10 %, 500 mA
Time for full charge	3.5 h
Battery charger polarity, \varnothing external, \varnothing internal	 3 mm, 1.1 mm
Pulse repetition frequency	1 \div 125 Hz
Maximum energy delivered per pulse	4.9 mJ/pulse
Maximum peak power	6.1 W
Maximum continuous power	735 mW
Average output current at max. frequency	10 mA
Average output current at min. frequency	80 μ A
Max. peak output current	120 mA
Pulse shape	Rectangular, biphasic (positive component and negative component) no component in direct current
Duration of the pulses (for each phase)	40 \div 350 μ s
Preset programmes	95
Customisable programmes	20
Number of channels	4
Dimensions (mm)	128 x 58 x 18
Weight (g)	110
Body	Polycarbonate
Applied part type	BF
IP	40

All the parameters shown in the table refer to a single channel with a typical charge of 500 Ω

PERSONALISED PROGRAMS from 01 to 20

- You can change the time, frequency and impulse amplitude of the programs in the table below at your discretion.

Prog.	Wave	Parameter	Treatment Indication
1	Waims-System	1 → 2 → 3 → 4 (pause 1) 1 → 2 → 3 → 4	Lymph draining, Relaxing, Active recovery, Tiredness
2	Waims-System	1 → 2 → 3 → 4 (pause 2) 1 → 2 → 3 → 4	
3	Waims-System	1 → 2 → 3 → 4 (pause 3) 1 → 2 → 3 → 4	
4	Waims-System	1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 1) 1 → 2 → 3 → 4 → 3 → 2 → 1	
5	Waims-System	1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 2) 1 → 2 → 3 → 4 → 3 → 2 → 1	
6	Tapping 1	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	Decontracting massage, Reactivating, Cellulite
7	Tapping 2	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
8	Tapping 3	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
9	Tapping 4	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
10	Press 1	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	Stamina, Explosive force, Development force
11	Press 2	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
12	Press 3	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
13	Kneading 1	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	Tone-up, Modelling, Cellulite, Firm-up
14	Kneading 2	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
15	Kneading 3	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
16	Kneading 4	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
17	Modulation 1	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	Draining, Stretching
18	Modulation 2	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	
19	Burst 1	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	Aerobics, Winding-down
20	Burst 2	Tpulse=40÷350µs Freq.=1÷125Hz Time=1'÷99'	

TESMED PROGRAMS from 21 to 99

- Programs 21 - 99 are used for treatment with fixed time and stages (not alterable).

FITNESS			Prog.	Photo	Time Prog.	n. stages	Days
1	Increase arm and pectoral mass		23	1,37,40÷44	23'	3	30
2	Increase leg mass		59	51÷61,73÷76,78,79	21'	3	30
3	Boost pectoral mass		59	37	21'	3	30
4	Boost triceps mass		58	41	21'	3	30
5	Boost biceps and forearm mass		58	44	21'	3	30
6	Boost abdomen mass man		76	45÷48	29'	5	30
7	Boost quadriceps mass		58	54	21'	3	30
8	Boost upper back mass		58	68	21'	3	30
9	Stretching triceps		65	64, 66	10'	1	30
10	Stretching shin muscles		65	13,31,74,75	10'	1	30
11	Stretching quadriceps		65	54	10'	1	30
12	Stretching wrist muscles		65	31,34÷36,76÷79	10'	1	30
13	Lymphatic draining legs	WAIMS	1	10,14,58,76	30'	1	15
14	Lymphatic draining arms	WAIMS	2	3	30'	1	15
15	Aerobic		62	1÷79	21'	3	15
16	Anaerobic		68	1÷79	18'	3	15
17	Warm-up		21	1÷79	10'	2	*
18	Taut muscles		60	1÷79	21'	3	30
19	Inflamed muscles		61	1÷79	21'	3	30

* Follow your trainer or doctor's advice

/ No indications provided

BEAUTY			Prog.	Photo	Time Prog.	n. stages	Days
20	Lymphatic draining legs woman	WAIMS	66	10,11,33	23'	5	15
21	Lymphatic draining arms woman	WAIMS	74	3	40'	5	15
22	Lymphatic draining legs man	WAIMS	87	58,76	40'	5	15
23	Lymphatic draining arms man	WAIMS	74	64	40'	5	15
24	Soft cellulite arms		92	20	47'	2	30
25	Hard cellulite arms		89	20	45'	2	30
26	Cellulite abdomen		95	4÷7	47'	2	30
27	Cellulite hips	WAIMS	93	8,9	54'	3	30
28	Cellulite inside thighs	WAIMS	97	12÷14	49'	4	30
29	Soft cellulite outside thighs	WAIMS	91	31	52'	3	30
30	Hard cellulite outside thighs	WAIMS	98	31	50'	3	30
31	Thick cellulite side thighs	WAIMS	90	31,33	54'	3	30
32	Hard cellulite glutei + thighs	WAIMS	99	29÷31	50'	3	30
33	Light cellulite glutei	WAIMS	94	28	52'	3	30
34	Thick cellulite glutei	WAIMS	96	28	55'	3	30
35	Hard cellulite glutei	WAIMS	88	28	50'	3	30
36	Muscle definition		57	1÷79	21'	5	30
37	Localised slimming arms		26	2,20	24'	4	60
38	Localised slimming abdomen		76	4÷7	29'	5	60
39	Localised slimming hips		77	8, 9	19'	5	60
40	Draining inside thighs woman		38	12÷14,32	20'	2	15
41	Localised slimming outside thighs		37	31	26'	4	30
42	Localised slimming glutei		37	27,28	26'	4	30
43	Localised slimming inside thighs		37	12÷14,32	26'	4	30
44	Draining arms	WAIMS	22	20	30'	4	15
45	Draining arms and pectorals	⚠	76	1,37,44	29'	5	15
46	Draining abdomen area		76	4÷7	29'	5	15
47	Draining thighs		38	13	20'	2	15
48	Draining thighs + glutei		82	26	20'	1	15
49	Draining legs		38	10,11,33,58,76	20'	2	15
50	Draining heavy legs woman		82	11	20'	1	15
51	Draining glutei + thighs		82	29÷31	20'	1	15

* Follow your trainer or doctor's advice
/ No indications provided

BEAUTY		Prog.	Photo	Time Prog.	n. stages	Days
52	Draining inside thighs man	38	54,55,56	20'	2	15
53	Draining back of legs	38	10,11,14,33	20'	2	15
54	Modelling arms	29	2,20,40÷44, 64, 66, 67	28'	5	30
55	Modelling breasts	78	1	22'	4	30
56	Modelling abdomen area	76	4÷7	29'	5	30
57	Modeling and dorsal wellness	48	15,16,19, 21÷24,62,68,70	28'	5	30
58	Modelling thighs + hips	WAIMS 71	8	28'	3	30
59	Modelling legs	73	10÷14, 32÷36	15'	3	30
60	Modelling glutei	77	28	19'	5	30
61	Firm-up	70	1÷79	21'	3	30
62	Firm-up arms + pectorals	 78	1,2,20,37,38, 40÷42,44	22'	4	30
63	Firm-up breasts	 78	1	22'	4	30
64	Firm-up slack abdomen woman	76	4÷7	29'	5	30
65	Firm-up abdomen area	76	4÷7	29'	5	30
66	Firm-up back muscles	79	22,23,62,68,70	24'	4	30
67	Firm-up thighs	77	13,32	19'	5	30
68	Firm-up inside thighs	77	12÷14,32	19'	5	30
69	Firm-up outside thighs	77	31	19'	5	30
70	Firm-up legs	79	10÷14, 32÷36	24'	4	30
71	Firm-up hips + glutei	77	26	19'	5	30
72	Firm-up glutei	70	28	21'	3	30
73	Anti-stretchmarks breasts	72	1	15'	5	/
74	Anti-stretchmarks abdomen	89	4÷7	45'	2	/
75	Anti-stretchmarks arms	72	2,20	15'	5	/
76	Anti-stretchmarks thighs	84	31	15'	5	/
77	Anti-stretchmarks glutei	84	28	15'	5	/
78	Anti-stretchmarks hips	84	8,9	15'	5	/
79	Tone-up	69	1÷79	21'	3	30
80	Tone-up arms + pectorals	 69	1,2,20,37,38, 40÷42,44	21'	3	30

* Follow your trainer or doctor's advice
/ No indications provided

BEAUTY		Prog.	Photo	Time Prog.	n. stages	Days
81	Tone-up pectorals woman	80	1	21'	3	30
82	Tone-up pectorals man	69	37,38	21'	3	30
83	Tone-up arms man	36	40,41,43,44,66	20'	2	30
84	Tone-up forearms man	36	43	20'	2	30
85	Combined tone-up upper body woman WAIMS	47	21÷24	26'	5	30
86	Tone-up abdomen woman	75	4÷7	20'	4	30
87	Tone-up abdomen post-natal	67	4÷7	21'	3	30
88	Combined tone-up abdomen + pectorals man WAIMS	56	38	26'	5	30
89	Tone-up upper body muscles	69	70	21'	3	30
90	Tone-up abdomen man WAIMS	56	45÷48	26'	5	30
91	Tone-up abdomen man (gutting)	75	46	20'	4	30
92	Tone-up flexors man	69	74	21'	3	30
93	Tone-up back muscles	64	62,68	21'	3	30
94	Combined tone-up abdomen + flanks man	75	46	20'	4	30
95	Combined tone-up upper back + shoulders man	69	68	21'	3	30
96	Combined tone-up shoulders + upper arms man	69	64	21'	3	30
97	Tone-up flanks WAIMS	86	9,50	26'	4	30
98	Tone-up hips post-natal	80	8,9	21'	3	30
99	Tone-up hips man (love handles)	76	49	29'	5	30
100	Tone-up adductors man	69	59	21'	3	30
101	Tone-up inside thighs (anterior)	80	12÷14, 53÷55,59	21'	3	30
102	Tone-up thighs + hips	80	8	21'	3	30
103	Tone-up outside thighs	80	31	21'	3	30
104	Tone-up inside thighs	80	12÷14,32	21'	3	30
105	Tone-up glutei	80	28	21'	3	30
106	Tone-up quadriceps man	69	54	21'	3	30
107	Tone-up calves	80	77	21'	3	30
108	Tone-up legs man	69	51-61,73÷79	21'	3	30
109	Tone-up legs woman	80	10÷14,30, 31,33÷36	21'	3	30

* Follow your trainer or doctor's advice
/ No indications provided

SPORT		Prog.	Photo	Time Prog.	n. stages	Days
110	Relaxing arms + pectorals 	81	1,37,38,40,43,44	20'	1	*
111	Relaxing abdominals	81	45÷48	20'	1	*
112	Relaxing back muscles	81	19,23,24,64,68	20'	1	*
113	Relaxing legs	81	51÷60,73÷79	20'	1	*
114	Wind-down legs	82	51÷60,73÷79	20'	1	*
115	Sport wind-down legs	82	51÷60,73÷79	20'	1	*
116	Strength	32	1÷79	23'	3	30
117	Explosive strength	43	1÷79	23'	3	30
118	Explosive strength arms + pectorals 	43	37,40,43	23'	3	30
119	Explosive strength legs	43	51÷60,73÷79	23'	3	30
120	Max strength arms + pectorals 	51	37,40,43	23'	3	30
121	Max strength abdominals	42	45÷48	23'	4	30
122	Max strength back muscles	51	62,68	23'	3	30
123	Max strength legs	42	51÷60,73÷79	23'	4	30
124	Stamina	25	1÷79	26'	4	30
125	Stamina arms + pectorals 	44	37,40,43	23'	3	30
126	Stamina abdominals	53	45÷48	23'	3	30
127	Stamina back muscles	44	62,68	23'	3	30
128	Stamina legs	25	51÷60,73÷79	26'	4	30
129	Active recovery	34	1÷79	20'	2	*
130	Recovery after races/training	27	1÷79	20'	2	*
131	Aerobic stamina	62	1÷79	21'	3	30
132	Warm-up	21	1÷79	10'	2	*
133	Athletic tone-up legs	69	51÷60,73÷79	21'	3	30
134	Athletic tone-up arms	69	40-44,66,67	21'	3	30
135	Football training warm-up	49	1÷79	10'	2	*
136	Football training max strength	51	1÷79	23'	3	*
137	Football training explosive strength	52	1÷79	23'	3	*
138	Football training stamina	53	1÷79	23'	3	*
139	Football training active recovery	54	1÷79	15'	1	*
140	Football training recovery after training	55	1÷79	20'	1	*

* Follow your trainer or doctor's advice
/ No indications provided

SPORT		Prog.	Photo	Time Prog.	n. stages	Days
141	Football training lymphatic draining WAIMS	87	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
142	Ski training warm-up	40	1÷79	10'	2	*
143	Ski training max strength	42	1÷79	23'	4	*
144	Ski training explosive strength	43	1÷79	23'	3	*
145	Ski training stamina	44	1÷79	23'	3	*
146	Ski training active recovery	45	1÷79	15'	5	*
147	Ski training recovery after training	46	1÷79	20'	1	*
148	Ski training lymphatic draining WAIMS	87	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
149	Cycling training warm-up	21	1÷79	10'	2	*
150	Cycling training max strength	23	1÷79	23'	3	*
151	Cycling training explosive strength	24	1÷79	23'	3	*
152	Cycling training stamina	25	1÷79	26'	4	*
153	Cycling training active recovery	63	1÷79	15'	3	*
154	Cycling training recovery after training	27	1÷79	20'	2	*
155	Cycling training lymphatic draining WAIMS	87	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
156	Swimming training warm-up	40	1÷79	10'	2	*
157	Swimming training max strength	42	1÷79	23'	4	*
158	Swimming training explosive strength	43	1÷79	23'	3	*
159	Swimming training stamina	25	1÷79	26'	4	*
160	Swimming training active recovery	83	1÷79	15'	3	*
161	Swimming training recovery after training	46	1÷79	20'	1	*
162	Swimming training lymphatic draining WAIMS	74	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*

* Follow your trainer or doctor's advice
/ No indications provided

SPORT		Prog.	Photo	Time Prog.	n. stages	Days
163	Running training warm-up	21	1÷79	10'	2	*
164	Running training max strength	23	1÷79	23'	3	*
165	Running training explosive strength	24	1÷79	23'	3	*
166	Running training stamina	33	1÷79	23'	3	*
167	Running training active recovery	34	1÷79	20'	2	*
168	Running training recovery after training	27	1÷79	20'	2	*
169	Running training lymphatic draining WAIMS	87	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
170	Tennis training warm-up	49	1÷79	10'	2	*
171	Tennis training max strength	42	1÷79	23'	4	*
172	Tennis training explosive strength	52	1÷79	23'	3	*
173	Tennis training stamina	33	1÷79	23'	3	*
174	Tennis training active recovery	83	1÷79	15'	3	*
175	Tennis training recovery after training	27	1÷79	20'	2	*
176	Tennis training lymphatic draining WAIMS	74	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
177	Volleyball training warm-up	21	1÷79	10'	2	*
178	Volleyball training max strength	51	1÷79	23'	3	*
179	Volleyball training explosive strength	43	1÷79	23'	3	*
180	Volleyball training stamina	25	1÷79	26'	4	*
181	Volleyball training active recovery	83	1÷79	15'	3	*
182	Volleyball training recovery after training	46	1÷79	20'	1	*
183	Volleyball training lymphatic draining WAIMS	74	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
184	Basketball training warm-up	21	1÷79	10'	2	*
185	Basketball training max strength	51	1÷79	23'	3	*
186	Basketball training explosive strength	52	1÷79	23'	3	*
187	Basketball training stamina	25	1÷79	26'	4	*

* Follow your trainer or doctor's advice
/ No indications provided

SPORT		Prog.	Photo	Time Prog.	n. stages	Days
188	Basketball training active recovery	45	1÷79	15'	5	*
189	Basketball training recovery after training	27	1÷79	20'	2	*
190	Basketball training lymphatic draining WAIMS	87	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*

* Follow your trainer or doctor's advice
/ No indications provided

MESSAGES TENS** (soft intensity)		Prog.	Photo	Time Prog.	n. stages	Days
191	Relaxing massage neck + shoulderse	28	15	15'	1	15
		48		28'	5	30
		39		25'	5	30
192	Massage upper + lower back woman WAIMS	48	22	28'	5	30
		85		40'	5	15
		39		25'	5	30
193	Massage upper + lower back man WAIMS	48	70	28'	5	30
		85		40'	5	15
		39		25'	5	30
194	Massage legs woman	31	11	22'	4	15
195	Relaxing massage calves (cramp)	81	77	20'	1	15
196	Anti-stress massage	41	1÷79	20'	1	15
197	Wind-down massage	50	11	10'	2	15
198	Energetic massage	31	1÷79	22'	4	15
199	Relaxing massage	30	1÷79	15'	3	15
200	Toning massage	35	1÷79	21'	3	15
201	Massage TENS** all body	48	1÷113	28'	5	30

* Follow your trainer or doctor's advice
/ No indications provided

** These are massages without medical indications. The massage has a soothing and relaxing effect on the treated area.

TESMED VITALITY 2023

T.S.R.M. (Tension Soothing and Relaxation Massage)

TESMED VITALITY 2023 is a wide range of treatments aimed at promoting relaxation, well-being, and vitality in the treated area. Thanks to the Tesmed T.S.R.M. (Tension Soothing and Relaxation Massage) technology, these treatments focus on the muscles, regions with joints, and skin tissue to provide a rejuvenating experience.

An intriguing aspect of TESMED VITALITY treatments is their versatility. Each treatment is carefully designed to help improve the well-being of a specific area. In some cases, the massage can focus on relaxation, releasing tensions and accumulated stress. At other times, the energetic stimulation provides a new vitality and a feeling of freshness.

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Prog.	Onde	Paramètres			Indication Traitements
1	Waims-System	1 → 2 → 3 → 4 (pause 1) 1 → 2 → 3 → 4 →			Lymphodrainage** Relaxant, Récupération active, Défatigant
2	Waims-System	1 → 2 → 3 → 4 (pause 2) 1 → 2 → 3 → 4 →			
3	Waims-System	1 → 2 → 3 → 4 (pause 3) 1 → 2 → 3 → 4 →			
4	Waims-System	1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 1) 1 → 2 → 3 → 4 → 3 → 2 → 1			
5	Waims-System	1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 2) 1 → 2 → 3 → 4 → 3 → 2 → 1			
6	Tapping 1	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	Décontractant Capillaire, Réactivant, Cellulite
7	Tapping 2	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
8	Tapping 3	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
9	Tapping 4	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
10	Press 1	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	Force résistante, Force explosive, Développement de la force
11	Press 2	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
12	Press 3	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
13	Kneading 1	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	Tonification, Modelage, Cellulite, Raffermissment
14	Kneading 2	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
15	Kneading 3	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
16	Kneading 4	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
17	Modulation 1	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	Drainant, Stretching
18	Modulation 2	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	
19	Burst 1	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	Aérobique
20	Burst 2	Tpulse=40÷350µs	Freq.=1÷125Hz	Time=1'÷99'	

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SPORT		Prog.	Foto	Dauer Prog.	Anzahl Phasen	Tage
163	Training Laufen Warm Up	21	1÷79	10'	2	*
164	Training Laufen maximale Kraft	23	1÷79	23'	3	*
165	Training Laufen explosive Kraft	24	1÷79	23'	3	*
166	Training Laufen Widerstandskraft	33	1÷79	23'	3	*
167	Training Laufen Aktive Erholung	34	1÷79	20'	2	*
168	Training Laufen Erholung nach dem Training	27	1÷79	20'	2	*
169	Training Laufen Lymphknotendrainage WAIMS	87	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
170	Training Tennis Warm Up	49	1÷79	10'	2	*
171	Training Tennis maximale Kraft	42	1÷79	23'	4	*
172	Training Tennis explosive Kraft	52	1÷79	23'	3	*
173	Training Tennis Widerstandskraft	33	1÷79	23'	3	*
174	Training Tennis Aktive Erholung	83	1÷79	15'	3	*
175	Training Tennis Erholung nach dem Training	27	1÷79	20'	2	*
176	Training Tennis Lymphknotendrainage WAIMS	74	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
177	Training Handball Warm Up	21	1÷79	10'	2	*
178	Training Handball maximale Kraft	51	1÷79	23'	3	*
179	Training Handball explosive Kraft	43	1÷79	23'	3	*
180	Training Handball Widerstandskraft	25	1÷79	26'	4	*
181	Training Handball Aktive Erholung	83	1÷79	15'	3	*
182	Training Handball Erholung nach dem Training	46	1÷79	20'	1	*
183	Training Handball Lymphknotendrainage WAIMS	74	3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76	40'	5	*
184	Training Basketball Warm Up	21	1÷79	10'	2	*
185	Training Basketball maximale Kraft	51	1÷79	23'	3	*
186	Training Basketball explosive Kraft	52	1÷79	23'	3	*
187	Training Basketball Widerstandskraft	25	1÷79	26'	4	*

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Complies with the following directives:

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Complies with applicable technical standards:

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Contacts, assistance, guarantee, useful information

info@tesmed.com

(please specify the date of purchase)

+ 39 035 662660

every working day from 8.30 am to 12.30 am and from 1.00 p.m. to 5.00 pm

WhatsApp Tesmed Fast Check Video

+ 39 345 6979093

If appropriate, it is possible to send a short video concerning the device, to demonstrate any operating problems - the message must be accompanied by a brief written explanation, with a screen of the order number and the date of purchase. Service active on working days and hours.

After your purchase, to activate the warranty, write to: **warranty@tesmed.com** indicating in the subject: the order number - the date of purchase - Max 830

A copy of the manual in electronic format can be obtained on the **tesmed.com** website at the INFORMATION AND MANUALS page by selecting Manuals

Contacts, assistance, garantie, informations utiles

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Feldi S.r.l. - Italy