Front Support - Rectus Abdominus

How can HEALI Bump help?

- Helps to add support and relieve pain by lifting your bump
- Assists to relieve pressure in lower belly
- Reduces irritation of itchy skin/stretched out skin
- Preventing diastasis recti (common condition in pregnancy causing the splitting of central abdominal line)

What do you need?

- Prepped skin that has been cleaned with water or alcohol to get rid of oil or residue
- ✤ 3 Strips of HEALI tape



How to Tape

Tear the center of the paper backing and slowly peel back the paper

A B C A B C A B C A B C MORKING ZONE DURING ZONE DURING ZONE DURING ZONE DURING ZONE

Apply one strip horizontally across the lower belly underneath your belly button





How to Tape

Using the two additional strips: place one end of the tape under the belly and to the right side of the centerline. Repeat on the left side of the centerline.



**Rub in the tape to activate the ingredients ** Only remove tape when fully dry