KONUS[®]
Optical & Sport Systems

10X-300X LCD DIGITAL MICROSCOPE

DIGISCIENCE

USER'S MANUAL

#5024 10X-300X DIGITAL LCD MICROSCOPE DIGISCIENCE

FIGURE1:

1. Micro-SD card slot - 2. Video-out - 3. HDMI-out - 4. Focus tube - 5. Polarizer wheel - 6. Charging indicator
 7. LED light for slice - 8. Height adjustment - 9. Capture Button - 10. 5-inch IPS LCD display - 11. Dimmer for LEDs in the lens tube - 12. Digital zoom - 13. Dimmer for LEDs in the platform - 14. USB/DC-in - 15. Power on/off
 16. Setup - 17. Up Button - 18. Confirm - 19. Down Button - 20. Mode Button (Photo/Video/Playback)
 Press on/off button (15) to turn on or off the microscope.

GET FOCUS & MAGNIFICATIONS

You can get different magnifications by rotating the height adjustment wheel (8) and focus tube (4). The lower the lens, the higher the magnification will be.

For observing at low magnification rotate the height adjustment wheel to move the lens upward to the top, and then rotate the focus tube (4) toward left to get a fine focus. This magnification rate is suitable for getting an overall observation on objects like coins, stamps, jewelry, electronics etc.

For observing at high magnification rotate the height adjustment wheel to move the lens downward to stop at the lowest height, and then rotate the focus tube (4) toward right to get a fine focus.

This magnification enables you to see the details of an object.

For observing specimen slides: at the High magnification position, turn on the LED (7) at the bottom platform, and put a slice on it, you can see the details of a specimen slice.

CAPTURE PHOTO & VIDEO & PLAYBACK MODE

Press Capture button (9) to take a photo or record video. Press Mode button (20) to switch among Photo, Video and Playback.

For deleting files, at play back mode, press Setup button (16) to enter sub-menu, then confirm delete. (FIGURE 2)

For protecting files, at play back mode, press Setup button (16) to enter sub-menu, then confirm protect pictures from delete. (FIGURE 3)

USE THE POLARIZER (5)

Polarizer, by removing light reflection, is especially useful for observing objects with shiny surface like metal, diamond, jewelry etc. By rotating the polarizer wheel (5) by 90-degree, you can activate or deactivate polarizing effect.

CONNECT TO COMPUTER

Connect the microscope to computer by USB cable. Choose Mass Storage and then you can view stored photos & videos.

Or you can quit the micro-sd card and connect the microscope to the computer directly with USB cable. Choose the option "PC-Camera" and start the "PortableCapturePlus" software. In this way you can see directly on the computer display what the microscope is observing. Using the microscope in this way you cannot use the microscope's buttons and you have to use the software options.

OUTPUT TO EXTERNAL MONITOR

Insert the 1/8-inch plug on the end of the supplied video cable into the microscopes Video port (2) Or by HDMI. Insert the RCA plug on the video cable into one of your television's yellow video input connectors. Refer to your television's manual for more information on locating these connectors. Or by HDMI. 3. Turn on the television. If necessary, switch your television to the appropriate input channel. Items that you magnify appear on the television screen.

CHARGE BATTERY

Plug the small USB end of the power adaptor into the USB/DC-in port located on the bottom side of the unit. Then plug the other end of the power adaptor into a standard AC outlet. The battery indicator LED at the side of the unit will be in red to show that the battery is charging. When the battery is finished charging, the battery indicator LED will be in off. Unplug the unit.

SCREEN MENU INTRODUCTION

You can press button 16 to enter camera setup. At photo mode, press button 16 again, you can choose to select a resolution for photo taken among 12M/10M/8M/5M/3M. (FIGURE 4).

At video mode, press button 16, you can choose to select a resolution for video recording between 1080P/720P.

- **CAPTURE MODE:** At photo mode, press button 16, you can choose Capture Mode to setup interval capture. For

example, setup the microscope to take a photo every 5 seconds. (FIGURE 5)

- **LOOP RECORDING MODE:** At video mode, press button 16 to choose to setup the length of each video clip; Off means not to turn on loop recording. (FIGURE 6)

- **DATE/TIME MODE:** Press button 16 and scroll down with button 19 to choose to setup system date and time.

At photo mode, press button 16 to choose "date stamp" to setup whether put date and time mark on capture photos and videos. (FIGURE 7)

- **TV-OUT MODE:** Press button 16, you can choose to select between NTSC and PAL for the TV output of the microscope, make it compatible with that of your TV. (FIGURE 8)

- **AUTO POWER OFF:** Press button 16 and you can scroll down to choose to select a time for the microscope to power off automatically without key operation. (FIGURE 9)

- **LANGUAGE:** Press button 16 and you can choose to setup language.

- **FORMAT SD-CARD:** Press button 16 and you can choose (selecting OK) to format the inserted Micro-SD card. (FIGURE 10)

- **DEFAULT SETTING:** Press BUTTON 16, you can choose to reset the microscope configuration to default setting. (FIGURE 11)

PORTABLE CAPTURE PLUS SOFTWARE

- **INSTALL THE SOFTWARE ON PC (WIN/MAC):** Insert the driver CD into CD-ROM Drive and then open the CD. Double click to install PortableCapture Plus and the installation wizard will guide you through the whole process.

The software PortableCapture Plus is available on the included CD and www.cameradownload.net.

- **START MICROSCOPE:** Connect your Microscope to your PC USB port, start the software by clicking the icon generated on the desktop after installation and you will see the following screen (FIGURE 12) pop up.

- **CAPTURE PHOTO:** 

- **CAPTURE VIDEO:** 

- **TIMER (PHOTO / VIDEO):** 

- **CALIBRATION AND MEASUREMENT:** 

Please make calibration before measurement. Display Crossing from 'Options>Crossing'. Place the enclosed calibration chart on the specimen holder of the microscope. Focus microscope on the attached calibration ruler  and then adjust the magnification till the picture is clear. Take a picture for the selected area when the dial is parallel to the crossing as below. (FIGURE 13)

Double-click on the picture to enter into measurement interface as below.

Click Calibration icon,  and then move mouse to the picture.

Now a dialog will popup as below and you need to enter the exact value of the measurement sample into 'Actual dimension'.

Then the software will automatically calibrate the magnification rate. Click 'OK' and the adjusted magnification rate will show in the upper right corner. (FIGURE 14)

* Remark: The measured line should be paralleled.

ow you can use 'Any angle Line' under Measurement icon  to measure the value again. If the value you measure by software is nearly the same as the exact value of the dial means the calibration finish.

Double click to open the picture, Click  > Reset magnification > Reset picture magnification' as below. And then start calibration again

If you need to change the magnification rate and adjust focus to capture picture, please calibrate the magnification rate again.

VISIT OUR WEBSITE: <http://www.konus.com/en/Catalogue/Scienza/Microscopes/Slides/> in order to know all the prepared slides series with science and biology arguments.

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