PELVIC FLOOR EXERCISE TRAINER GUIDE

Before you start

- Begin with the lowest weighted trainer.
- Be patient, it can take a 6-12 weeks of daily training before your vaginal and pelvic muscles getting stronger.
- Make sure you have privacy and time when you use the weighted trainer.
- Stop if you experience pain.

Inserting the weighted trainer

- Always wash your weighted trainer and hands before and after use with warm water and a mild soap and rinse thoroughly. Dry with a clean towel.
- Apply waterbased lubricant to the tip of the weighted trainer for easy insertion.
- Lie on your back with your knees bent and slightly apart.
- Insert a weighted trainer. Make sure that the retrieval cord is outside your vagina.

Exercise A.

Tighten your pelvic floor muscles gentle to lift and squeeze the trainer. Hold for 5 seconds and relax for 10 seconds. Repeat 10 times and do 3 sets a day. Build the movement up to 10 seconds.

Exercise B.

Tighten your pelvic floor muscles gentle to lift and squeeze the trainer. Hold for 1 second and relax for 1 second. Repeat 5 times and then rest for 10 seconds. Repeat this exercise 5 times. Do 3 sets per day.

Exercise C.

Gently tighten your pelvic floor muscles and feel if the trainer is seated properly. Start your day with normal movements such as walking, bending over, etc. Do 2 sets a day of 15 minutes.

Important

This is not a medical device. We recommend that women consult with a health professional specializing in pelvic floor disorders for the most effective physical therapeutic benefit.

This information is not intended to be a substitute for your medical treatment. The information given for this product should not be used for the diagnosis of any medical condition. A licensed physician should be consulted for diagnosis and treatment of any and

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