



Fingertip Oximeter
OLED



Fingertip Pulse Oximeter


Features:

1. The product adopts double color OLED display and can show the six different display modes.
2. Low-power consumption, continuously work for more than six hours with 2 pieces AAA batteries
3. A Battery life indicator
4. In the absence of signals, the product would power-off automatically after 8 seconds.
5. Small, light in weight, and convenient to carry.

Keys function description: In standby mode, start the key instrument into the working state, Touch this button under the working state can change the display mode.

⚠Note⚠: The machine profile picture is only for reference use, specificity in kind prevails.

OLED Display Modes:



Operation Instructions

1. Install a pair of AAA batteries into the battery compartment.
2. Plug a finger into a rubber hole of the Oximeter.
3. Press button on the front panel to start.
4. Do not shake your finger or move your body when the Oximeter is working.
5. Press the button if you want to change display direction.
6. Read the relevant value from the display screen.
7. The device will automatically power off in 8 seconds if there is no finger signal.

When plugging your finger into the Oximeter, your nail surface must be upward.

Declaration: Please use the medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the Oximeter is a medical rubber which has no toxin, causes no harm, and brings no side effect such as allergy to the skin).

Hang Rope Installation

1. Put the rope thin end through the hole.
2. Put the rope coarser end through its already wearing thin end part and tighten.

Product Accessories

One user manual

One lanyard

QC card

Parameter Setup:

Press start button(>3s), into parameter setup, As menu 1:

1. In menu 1, When the “*” signal is shown on the “Alm Setup”, Press the button (>3s) and enter into the menu 2. Press the button (<1s) can select item, then Press button (>3s) to set the on/off for alarm, beep, demo and screen brightness adjustment (optional "1", "2", "3" and "4"). When the “*” signal is shown on the “Restore”, Press the button (>3s) and all the settings are back to the factory settings.

⚠Attention⚠

- **Using devices outdoors or under strong light; please adjust the screen brightness to a higher level for proper observation.**

- **It's better for users to choose a lower brightness to conserve battery power.**

2. In menu 2, When the “*” signal is shown on the “Sounds Setup”, Press the button (>3s) and enter into the menu 1, Press the button (<1s) to select an item, and then Press the button (>3s) to setup data. Choose “+” or “-” to increase or decrease values.


Settings	
Alm Setup	*
Alm	off
Beep	off
Demo	off
Restore	ok
Brightness	3
Exit	

Settings	
Sounds Setup	*
ScP2 Alm Hi	100
ScP2 Alm Lo	100
FR Alm Hi	100
FR Alm Lo	100
+/-	+
Exit	

Detailed Descriptions of Product Functions:

1. Display Type: OLED display
 2. SpO₂: Measurement range: 70% ~ 99%
Accuracy: $\pm 2\%$ on the stage of 80% ~ 99%;
 $\pm 3\%$ on the stage of 70% ~ 79%;
Below 70% no requirement
Resolution: $\pm 1\%$
 3. PR: Measurement range: 30BPM ~ 240BPM
Accuracy: ± 1 BPM or $\pm 1\%$ (the larger one)
 4. Power: two AAA 1.5V alkaline batteries
 5. Power consumption: below 30mA
 6. Automatic power-off: The product power off by itself when no finger is in the product for ≥ 8 seconds.
 7. Operation Environment:
Temperature: 5 °C ~ 40 °C
Humidity: 15% ~ 80%
Atmospheric pressure: 70kPa ~ 106kPa
- Transport, storage environment:
Temperature: -10 °C ~ 55 °C
Humidity: $\leq 95\%$
Atmospheric pressure: 50kPa ~ 106kPa, non-corrosive gas and well-ventilated environment.

Diagram of Operation Principle



1. Infrared-ray receiving tube
2. Infrared-ray transmitting tube



1300 736 330
www.medescan.com.au

Visit our YouTube channel for our latest easy to follow instructional videos. Simply search for MedeScan on YouTube and subscribe to our channel for the latest products and instruction videos.

You can also scan the QR code for instant access.

Designed in Australia
Made in China