



WARM STEAM VAPORISER



Customer Service:
1300 736 330

WSV0124

IMPORTANT SAFEGUARDS:

When using electrical appliances, basic safety precautions should always be followed, including the following:

1. Read all instructions prior to use.

2. Always place unit on a flat, level surface. **Do not** place the vaporiser directly on wooden furniture or other surfaces that could be damaged by water. A waterproof pan, tray, mat or dropcloth is recommended. Exposure to excess water or moisture can cause damage to some surfaces.
3. **Do not** place the appliance near heating devices. Keep cord away from heavy traffic areas and do not place cord under rugs or allow cord to hang over a countertop or table where it can be pulled on or tripped over.
4. **Never** place a vaporiser where it is accessible to children.
5. **Caution:** To reduce risk of burns, supervision is recommended when a vaporiser is used near children, immobile patients, or pets.
6. Position the unit so mist stream is aimed away from children, walls and furniture.
7. **Use only** 220V AC outlet. Be sure unit is plugged in firmly and completely. As with any electrical device, plug and unplug unit with dry hands only.
8. To avoid the risk of fire, shock or personal injury, do not use an extension cord.
WARNING: Never operate unit with a damaged cord or plug.
9. Fill water tank to the bottom of the "MAX FILL" line. **Do not** overfill. **Do not** block moisture outlet. This will damage the vaporiser.
10. **Do not** add medication of any type into the moisture outlet or water tank. **Do not** cover or place anything over any opening. **NEVER** drop or insert any object into any opening when plugged in.
11. Turn the vaporiser off by pushing switch to "Off" when not in use. The vaporiser should always be unplugged and emptied when not in operation or while being cleaned. Unplug vaporiser and allow to cool before moving. **Do not** move or tilt vaporiser while it is in operation. Plug and unplug unit with dry hands. **Never** pull by cord.
12. **Do not** touch unit or steam vapor during use. Steam is hot and burns can occur.
13. Use this appliance only for its intended household use as described in this manual. Any other use not recommended by the manufacturer may cause fire, electric shock or injury to persons and property damage.
14. **Do not** remove the Warm Mist Control Unit (upper housing) during operation or within 30 minutes after it is unplugged. Serious injury may result. **DANGER:** To guard against electric shock, never immerse the Warm Mist Control Unit (upper housing) in water or other liquids.
15. Regular cleaning of the unit is needed. Please refer to the cleaning section in this manual.
16. **Do not** operate outdoors. This vaporiser is intended for indoor use.
17. **Do not** pour any liquid or oils into the mist output (see Fig.1). Pour it/them only in water tank.
18. **Do not** attempt to carry the vaporiser by the head. The vaporiser should be unplugged and carried by the tank only.

ADDITIONAL SAFEGUARDS

Consumer Safety Information

Note: This is an electrical appliance requiring attention when in use.

Note: If moisture forms on the walls or windows of the room, turn off the vaporiser. The room already has plenty of humidity and additional moisture may cause damage.

PARTS



1 Red Light Therapy On/Off

2 Power On/Off

3 Warm Steam Outlet

4 Water Tank

5 Power Cable (1.8 m)

6 Warm Steam Flow Adjustment Tab

7 Warm Mist Control Unit

8 Locked / Unlocked Indicator

HOW TO OPERATE YOUR WARM STEAM VAPORISER

WARNING:

To protect against electrical shock, always unplug your vaporiser when not in use, before cleaning, rinsing, filling or moving. When in use, the Warm Mist Control Unit must be positioned on the water tank for proper operation.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities.

Keep out of reach of Children

DIRECTIONS FOR USE:

1. Always unplug vaporiser when inserting or removing the Warm Mist Control Unit. Remove the Warm Mist Control Unit. The unit will come in locked position. turn the unit to the unlocked symbol and lift the Warm Mist Control Unit out of the water tank (See Fig.1) .

2. Always rinse and empty the water tank before filling and after use.

✔ Fill water tank with cool tap water to the "MAX FILL" line (see Fig.2)

✔ Insert Warm Mist Control Unit to the water tank and turn it to the lock symbol

✔ Place the vaporiser on a secure flat surface away from children.

✔ **Optional:** Add your favourite essential oil inside water tank (see Fig.3).

3. Ensure dry hands and plug power cord into a 220VAC wall outlet only.

4. Slide the ON / OFF switch to the "On" position showing the RED indicator light (see Fig.4).

Steam will be produced after 10 minutes. Please note: mineral content in water can vary greatly by region, please refer to and follow the troubleshooting if only a little steam is produced.

5. Push the button to "Off" position if you wish to shut the unit down (see Fig.4).

Red Light Therapy

Push the Red Light Therapy ON / OFF button and the Red light will switch on to help give you and your family a restful night's sleep.



FIG 1



FIG 2



FIG 3







FIG 4

CLEANING AND MAINTENANCE INSTRUCTIONS

IMPORTANT: Clean your vaporiser regularly (after each use) to help eliminate odors and prevent bacteria growth. Depending on the mineral content of the water and frequency of use, the PTC and water inlet hole may become clogged from mineral build-up that occurs during normal use.

1. Unplug vaporiser.
2. Remove Warm Mist Control Unit from water tank by turning it counter-clockwise and lifting it straight out.
3. To remove mineral deposits from the water inlet hole, use a wooden toothpick to unclog the hole. **Do not** use any metal object or material that could enlarge the water inlet hole (see trouble shooting section).
4. Reposition the Warm Mist Control Unit making sure the cover is securely locked before operation. **Do not** attempt to operate the unit disassembled.
5. **Do not** attempt to remove or loosen the metal PTC, as these parts are attached under factory controlled conditions and are not replaceable. Clean and dry unit thoroughly prior to storage.

TROUBLESHOOTING

<p>Spitting water</p>	<p>The water level is over the "MAX FILL" line, make sure the water level is under the "MAX FILL" line.</p>	
<p>No steam</p>	<p>Make sure the device is plugged into the power supply, and make sure the switch is in the "On" position.</p>	
<p>Too little steam</p>	<ul style="list-style-type: none"> • If using mineral or distilled water, switch to tap water. • Add one to two pinches of salt to the water and mix thoroughly. • If this does not fix the problem, please see below 	
<p>Mineral deposits build-up</p>	<p>Clean incrustation from inlet holes.</p>	

If additional assistance is needed, please contact Medescan Customer Service Team

NOTE:

During the operation of unit, the indicator light will remain on (RED) indicating that the vaporiser is plugged in and turned on.

The vaporiser will automatically shut itself off when water reaches a minimum level. The vaporiser will stop producing steam and the night light will remain on until switched off. Once the unit is turned OFF - unplug and allow the vaporiser to cool before refilling the water tank and repeat steps above to operate the unit again.

For safety and best performance, the unit will be powered off automatically when the Warm Steam Control Unit is (see Fig.1) lifted even if plugged in.

Discoloration/Black Particles Note

After use, the water residue in the bottom of the water tank may contain black particles (**below Fig.3**), and in some areas, the water and the bottom of water tank may discolor. These conditions are normal and do not effect the purity of the steam or medicine. The water tank is easily cleaned with soap and water. Rinse and dry the water tank after cleaning. For more details of **CLEANING AND MAINTENANCE**, please refer to instruction manual.



FIVE YEARS LIMITED WARRANTY

This product carries a warranty that it will be free from defects in material and workmanship for a period of five years from the date of purchase. This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Keep your original sales receipt.

IMPORTANT: This warranty does not cover damages resulting from accident, misuse or abuse, lack of reasonable care, the affixing of any attachments not provided with the product, loss of parts, or subjecting the appliance to any but the specified voltage. Read directions carefully.

For all customer service issues, please contact Medescan Customer Service on 1300 736 330 or email: info@medescan.com.au

SAVE THIS FOR YOUR RECORDS

Medescan logo is a registered trademark
Made in China





STOP!

**In the unlikely event that
you have a problem
or question about your unit,
please contact our customer
service team.**

**WE ARE ALWAYS HAPPY
TO HELP!**

1300 736 330