



NO ITCH

Insect Bite Healer

MedeScan

QUICK START GUIDE



How to ditch the itch

Step 1

Press the On / Off button to power the No Itch ON

Make sure there is no insect sting in the skin.
Your MedeScan No Itch is now powered on.

- Temperature setting 37°C
 - Sensitive skin setting ON.
- 3 seconds treatment.

Solid Blue Indicator light.

Step 2

Place the heating plate directly on the insect bite for rapid relief. Heat treatment begins on contact and once you hear the beep, the treatment is complete.

Changing the temperature setting

When you treat an insect bite and you still feel the itch, it could mean that you need to increase the temperature. Go up as much as you can comfortably accept. Higher temperatures are also suitable for thicker skin areas like the elbow and heel.

Press and hold the temperature button. The temperature will change from 37°C, to 42°C to 47°C. The indicator light will change:

Blue	= 37°C
Orange	= 42°C
Red	= 47°C (Caution. HOT!)

Changing Treatment Mode

If the temperature is just right but you need a little longer for better effect, you can change from sensitive skin mode to regular skin mode or you can repeat the treatment.

After switching the No Itch on, short press the power button and the Mode will change from Sensitive skin to Regular skin mode (3 seconds to 6 seconds treatment). You can change the temperature normally as instructed above.





Keep Out of Reach of Children!

HEATING DEVICE !



How to change the batteries

Pinch both sides of the battery compartment and pull down firmly.

Replace the batteries, 2 x AA size batteries.

Replace the battery cover .

Enjoy chemical free fast relief.

General Warnings / Safety Notes

- Do not use the device if you suffer from diabetes. Diabetics may have reduced sensitivity to pain and may burn themselves.
- Do not use the device on children under 3 years of age or people with disabilities.
- Do not use on people with reduced sensitivity to heat. This includes people with skin conditions or scar tissue in the application area, taking pain relief medication, alcohol or any drugs.

- Do not use the No Itch on injured or irritated skin (e.g. inflamed skin, whether painful or not, reddened skin, rashes, (e.g. allergies), burns, frostbite, bruises, swellings, both open and healing wounds, and post-operative scars where the healing process can be affected).
- Do not use the No Itch if you are suffering from a fever.
- Do not use the No Itch on heat sensitive skin.
- Do not use the No Itch if you have any impaired senses that may reduce the feeling of pain.

- Do not use the No Itch with any lotions, creams or gels in the same area.
- Do not use the No Itch if you suffer from irritated skin due to long-term heat application.
- Any person using it on themselves must be over 12 years old. The No Itch should only be used on children under 12 years old if supervised or assisted by an adult.
- Before using the No Itch, consult your doctor if you have any skin conditions, serious conditions or malignant growths or unexplained pain in any part of the body.
- Seek medical advice immediately if the symptoms persist despite using the No Itch and stop using it.
- Seek medical advice immediately if the symptoms worsen whilst or after using the No Itch and stop using it.
- Seek medical advice immediately if it is a tick bite.
- Seek medical advice immediately if the first signs of increased itching, skin redness, swelling, dizziness, nausea, shortness of breath, formation of wheals or a drop or increase in blood pressure occur.
- Do not use near the eyes, on the eyelids, mouth or on mucous membranes.
- Please note to sanitise the device between uses, especially when using on different people.



1300 736 330

www.medescan.com.au

Visit our YouTube channel for our latest easy to follow instructional videos. Simply search for Medescan on YouTube and subscribe to our channel for the latest products and instruction videos.

You can also scan the QR code for instant access.