Guida didattica / Didactic guide

Guía didáctica

Cod. 5222

OPTIKA S.R.L.

VIA RIGLA, 30 – 24010 PONTERANICA (BERGAMO) – ITALY
Tel. +39 035 571392 - Fax +39 035 571435

www.optikascience.com

info@optikascience.com

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
Fig. 1


Fig. 2

CATHODIC RAYS TUBE FOR MAGNETIC DEVIATION

According to experimental conditions (electric field and tube rarefaction), cathodic rays have a speed that varies from 40.000 Km/s to 100.000 Km/s.


They have the following features:

- they flow in perpendicular direction to the cathode surface;
- they have a calorific action;
- they have a mechanic action;
- they excite the fluorescence and phosphorescence of corpses hit by rays;
- they are constituted by electrons;
- they are deflected by an electric field;
- they are deflected by a magnetic field
- they are deflected by a magnetic field;


In order to demonstrate the property “g”, the tube code 5222 can be used according to the disposition indicated in picture no. 1.

OPERATIVE INSTRUCTIONS

- Link electrodes A-K to the spark gap of Rumhkorff spool (code 5207).
- Action the spool. Since at the tube bottom, inclined towards the cathode, is disposed a display that becomes fluorescent when it is hit by cathodic rays, it will be clearly visible the electrons run.
- Dispose one U magnet (code 5286) in the position indicated in picture no. 1. It will be clearly visible that the beam varies in the direction indicated in picture no. 2. A similar deviation effect it is obtained by moving closer a magnetic pole of a linear magnet. When it gets closer to the other pole, deviation is inverted.

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