

Back 2Sleep


NOTICE D'UTILISATION – INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG – INSTRUCCIONES DE USO
ISTRUZIONI PER L'USO – MANUAL DE INSTRUÇÕES
INSTRUKCJA DLA UŻYTKOWNIKA
GEBRUIKSAANWIJZING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ






→ Starter Kit

→ S, M, L OR XL

→ 1 → 2 → 3

Notice d'utilisation.....	2
Instructions for use.....	6
Gebrauchsanweisung.....	10
Instrucciones de uso.....	14
Istruzioni per l'uso.....	18
Manual de instruções.....	22
Instrukcja dla użytkownika.....	26
Gebruiksaanwijzing.....	30
Οδηγίες χρήσης.....	34

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Site web : www.back2sleep.eu

Instructions for use: Back2Sleep™

Back2Sleep™ and Back2Sleep™ Starter kit Nasal airway stent. Back2Sleep is a medical device intended for patients suffering from sleep apnea or snoring.

Back2Sleep allows the patient to breathe correctly during the night, which favors a great quality of sleep.

STARTER KIT: VERY FIRST TIME USE - FIND THE LENGTH THAT SUITS YOU BEST

The STARTER KIT contains 4 lengths: S, M, L and XL and 3 tubes of gel.

Exactly as you would with shoes, you must find the size that suits you best: the one that both stops snoring and is comfortable.

Note: use only one tube at a time. Each tube can be inserted indifferently into the left or the right nostril. In case of deviation from the nasal septum or if you have difficulty inserting it on one side, it is advisable to use the opposite nostril.

A/ Check the contents of the packaging according to the diagram (symbolized by an EYE👁️), inspect the content carefully. Do not use if packaging is damaged and refer to instructions for use.

B/ Before going to bed, start with the smaller tube (size S) and follow instructions 1 to 5 below.

C/ By opening the mouth wide in front of the mirror, check the position of the tube. Ideally, the end of the tube should extend 10 to 15mm beyond the veil of the palate (a) or even lower in case of obstruction in the tonsils or tongue base. If you cannot visualize the position of the stent in the mirror, start with the shortest stent and increase the length until you find the most comfortable and effective against snoring.

D/ If the tube is too short, it can cause great discomfort. In this case, remove it and restart the operation with the upper length. In any case, discomfort will be felt the first few nights. It disappears significantly after several uses.

E/ Go to sleep with the tube. If effective, the tube can be used for a maximum of 15 nights before being replaced by a new one. If not effective, you can use the upper length the following night (b).

F/ Record the size that suits you best (S, M, L or XL). You can now order the corresponding reference directly (c).

Unlike the Starter Kit which allows you to define the most suitable length for your morphology, Back2Sleep contains 2 stents of identical length (S, M, L or XL) and 6 tubes of gel. Each stent can be used for a maximum of 15 days.

How to use the nasal airway stent Back2Sleep™?

Before opening the pack, make sure you have chosen the right size. Back2Sleep is not a sterile medical device.

1. Take Back2Sleep™ or Back2Sleep™ Starter kit out of the bag and carefully place the tube in boiling water for 2 minutes. The tube needs to be totally immersed (1).
2. Using a strainer, carefully drain the water. **To avoid any risk of burn, wait for the tube to cool down** before touching it (2).
3. **Wash your hands** thoroughly **before touching and using the boiled Back2Sleep™ / Back2Sleep™ Starter kit** (3).
4. **Apply a thin film of lubricant over five centimeters** at the tip of the tube to be inserted. Make sure that only the circumference is lubricated by carefully avoiding the inside of the tube (4).
5. In front of a mirror, insert the lubricated part of the tube in the nose while breathing normally (5). The Back2Sleep™ metallic clip part must be horizontal, in order to clip to the central part of the nose (6).
6. After the insertion, make sure the final position of Back2Sleep™ is correct: open your mouth and look inside through a mirror. The tip of the tube should be visible and reach the soft palate, next to the uvula or longer if comfortable (a).
7. Go to sleep with Back2Sleep™ inserted (b).
8. In the morning, carefully remove Back2Sleep™ out of the nose by gently pulling on it (7).
9. Immediately after retraction from the nose, scrub the tube with lukewarm flowing drinking water for at least 30 seconds by hand until all potentially adhering mucus has been removed.
10. Carefully place the tube in boiling water for 2 minutes. The tube needs to be totally immersed (1).
11. Using a strainer, carefully drain the water. **To avoid any risk of burn, wait for the tube to cool down** before touching it (2).
12. Once the tube has cooled down, store it like define in the storage conditions below (8).
13. For the re-use of the Back2Sleep™ in the evening, proceed with steps 3 to 7 again.

⚠ Replace with a new product every 15 days regardless of whether it is used or not ⚠

⚠ Check before use, and if you see any signs of discoloration, odor, or visual cracks or damage, stop immediately and replace by a new product ⚠

For more information about Instructions for Use, visit the Back2Sleep™ website:
www.back2sleep.eu

Warnings - Precautions

- **A medical examination prior to the use of back2sleep is recommended; If in doubt, consult your doctor.**
- Read carefully these instructions before using the device and always use the device accordingly.
- Do not use the device continuously for more than 10 hours.
- Back2Sleep™ can cause discomfort the first 3-4 days of use, which will gradually disappear.
- The selection of the more suitable size for your morphology may be complex. It is important to listen to your feelings:
 - In case of a too short tube: that could cause a significant and disturbing contact feeling and a discomfort. In that case, the device won't be efficient.
 - If the extremity of the tube is at the uvula level, the uvula may stick on the tube causing discomfort at the level of the soft palate.
 - If the extremity of the tube reaches the base of the tongue, it could be pushed back.
- In all these cases, select a higher length of tube.
- Any event or events arising from non-observance of the instructions given in this leaflet, are the sole responsibility of the user.
- Make sure you are using the appropriate type and size of Back2Sleep™ that fits your morphology. Consult a physician if needed.
- Back2Sleep™ may be less effective for symptoms of narrowing or obstruction of the airways due to the tongue falling if it does not go beyond the soft palate (near the uvula).
- Make sure the nose clipper fits your morphology.
- There is a risk of nausea when the device is used incorrectly.
- If you experience respiratory disease, gastrointestinal disease, nasopharyngeal disease or tonsil hypertrophy, consult your physician before using Back2Sleep™.
- **In case of trauma or disorders of the nasal mucosa, if you feel a strong discomfort, pain, or frequent nausea when using the device, if you experience breathing difficulties, immediately discontinue use and consult your physician.**
- In case of accidental swallowing, consult a physician immediately.
- Thoroughly inspect the device packaging before opening. Do not use the device if any abnormalities are found.
- After opening the package, inspect the device carefully. Do not use the device if any abnormalities are found. Never use a device with a nose clipper widened or defective. If so, use another one.
- Do not insert in your nose without applying lubricant on the device before.
- Use only the appropriate lubricant provided with the device.
- Do not bend or deform the nasal clip. Do not detach the nasal clip from the silicone tube. Do not stretch the silicone tube.
- The user has to only use the cleaning protocol described in this Instructions for Use.
- The device should never be disassembled. If it is, never use a disassembled device.
- The user cannot use a product without the nose clipper being held in place on the peg.
- Do not try to pinch or deform the product even if the insertion is difficult.

- Do not try to force the insertion.
- The user must be careful when using boiling water during the cleaning steps described in this Instructions for Use.
- Do not put two devices in at the same time.
- Do not use Back2Sleep™ if the nasal clip does not remain firmly attached on your nose.

In case of a major incident, please contact Back2Sleep™: www.back2sleep.eu

Possible Adverse Events

Use of the device may cause the following adverse events:

- Difficulty during insertion and removal
- Discomfort during insertion and once the device is in place
- Nasal or pharyngeal foreign body sensation
- Nasal tissue trauma and light bleeding
- Transient dry cough
- Migration or involuntary expulsion of the device
- Allergies to the materials of the device

Contraindications

1. The following patients are not appropriate candidates for use of the device: Patients with disorders or illness of the nasopharynx airway, unconscious patients (excluding when sleeping), patients with agitated sleep or who are unable to use the device properly, with neuromuscular, cardiac and psychotic diseases, children, patients who have difficulties consulting a physician.
2. People under guardianship must use the device under supervision of the guardian.

Storage Conditions

Store Back2Sleep nasal airway stents in clean and dry conditions away from heat (between 4°C and 40°C) and sunlight in its original packaging. Keep out of children's reach. See specific storage conditions of the gel on its own packaging and follow them.

Expiration Date: Follow expiration date on the label.

Legal manufacturer of Back2Sleep nasal airway stents: Back2Sleep™ - 188 rue de Rivoli – 75001 – Paris – France

Pack producer of Back2Sleep nasal airway stents with gel tubes: Back2Sleep™ - 188 rue de Rivoli – 75001 – Paris – France

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A/ Controleer de inhoud van de verpakking volgens het schema (gekenmerkt door een OOG (©)), inspecteer de inhoud zorgvuldig. Gebruik niet als de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing.