INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO

myinnovo.com

RESTORE YOUR PELVIC FLOOR, YOUR CONTROL AND CONFIDENCE

Thank you for purchasing INNOVO®.

INNOVO® is an effective solution to restore your pelvic floor health and improve the quality of life of those living with a weak pelvic floor. Clinically proven, INNOVO® uses a pelvic floor fitness technology called 'Innovotherapy' which strengthens and restores the pelvic floor.

INNOVO® is a non-invasive device to safely, gently and accurately train the pelvic floor muscles, directly targeting the primary cause of stress incontinence. INNOVO® delivers 180 powerful contractions of the pelvic floor during each 30 minute session, giving you the confidence that the pelvic floor muscles are being accurately targeted with each repetition.

Innovotherapy's patented criss-cross technology sends multi-directional signals which go through the pelvic floor to generate strong contractions of the pelvic floor. Innovotherapy is easy and convenient to use, giving a comfortable treatment experience in the privacy of your own home.

Should you have any questions about your INNOVO® device, please visit www.myinnovo.com and we'll be happy to help.

INNOVO® will allow you to embrace life's freedom and get back to being you again.

TABLE OF CONTENTS

Before You Start Using INNOVO®	4
Indications for use, Contraindications, Warnings, Precautions and Adverse Reactions	4
INNOVO® Treatment Plan	7
INNOVO® Pack Contents	8
INNOVO® Controls	9
INNOVO® Sockets	9
Controller Display	10
Step 1: Charger Assembly	11
Step 2: Charging and Getting Started with INNOVO®	11
Step 3: Putting On the Shorts	12
Step 4: Get Into The Correct Position For You	13
Step 5: Using INNOVO®	14
5.1 Turning on the Controller and Starting a Session	14
5.2 Intensity Levels and What To Expect During a Treatment Session	14
5.3 How to Pause a Treatment Session	15
5.4 Finishing a Treatment Session	15
5.5 Removing the Shorts	15
Additional Product Functions	16
Information Button	17
In-Program Information	19
Problem Solving Guide	19
Electrode Detection Function	20
Caring for your INNOVO®	21
Replacing the Battery	22
Technical Specifications	23
Physician Information	26
Warranty	29

BEFORE YOU START USING INNOVO®

Important:

To ensure that you use this product safely and effectively, please follow the safety information provided in this instruction for use. INNOVO® should only be set-up and operated according to the instructions laid out in this instructions for use, the Quick Start Guide and on the website - www.myinnovo.com

INDICATIONS FOR USE, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS AND ADVERSE REACTIONS

INDICATION FOR USE:

The INNOVO® is a transcutaneous electrical stimulator indicated for treatment of stress urinary incontinence in adult females. The INNOVO® is indicated for prescription use only.

CONTRA-INDICATIONS:

- Do not use if you have an electronic implant (e.g. cardiac pacemaker or defibrillator) or suffer from any other heart problem.
- Do not use if you have any conductive intra-uterine devices (IUD/Coil) or metal implants in the abdominal or pelvic area, including hip and lumbar spine.
- Do not use if you are pregnant. The safety of powered muscle stimulators for use during pregnancy has not been established.
- Do not use in close proximity (e.g. 1 meter) to shortwave or microwave therapy equipment, as this may produce instability in the controller's output.
- Do not use if wearing the device necessitates placement over areas at which drugs/medicines are administered by injection (short term or long term) e.g. hormone treatment.

WARNINGS:

- For use by adults only.
- The long-term effects of chronic electronic stimulation are unknown.
- Simultaneous connection of a patient to high-frequency surgical equipment may result in burns at the site of the electrodes and possible damage to the controller.
- The product must not be used with any other product which delivers electrical current to the body (e.g. interferential or another muscle stimulator).
- Stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrodes placed on or across the chest and the upper back or crossing over the heart.
- Stimulation applied near the thorax may increase the risk of cardiac fibrillation (cardiac arrhythmias), due to the introduction of electrical current into the heart.
- Do not use if you have cancer in the region of application of the device or in the pelvis, epilepsy or are under medical supervision for cognitive dysfunction.
- Apply stimulation only to normal, intact, clean, healthy skin.
- Do not apply over areas of skin which lack normal sensation.
- Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g. phlebitis, thrombophlebitis, varicose veins), or areas that have undergone recent surgery.
- Avoid any recent scars, broken or inflamed skin, areas of infection or susceptibility to acne, thrombosis or other vascular problems (e.g. varicose veins), or any parts of the body where feeling is limited.
- The product should not be used while sleeping.
- Do not apply stimulation while driving, operating machinery, cycling or during any activity involving exercise in which electrical stimulation can put you at risk of injury.
- Do not apply stimulation when in or near water, including a bathtub, shower, hot tub or swimming pool.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g. cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation controller is in use.
- Wireless communication equipment such as wireless home network devices, mobile phones, cordless telephones, base stations and walkie-talkies of between 1W to 10W can affect this equipment and should be

kept a distance of least 2.3 to 7.3 meters away from the controller.

- Use this controller only with the Shorts cables and accessories recommended by the manufacturer.
- For hygiene reasons the Shorts are for single person use. To prevent cross-infection and Shorts should not be shared with other people.
- Always ensure that the charger/power supply is positioned so that it is easy to unplug the power supply from the socket.
- To isolate this equipment from the mains, the charger/power supply plug should be removed from the socket.
- Do not attach any cables or connectors to INNOVO®, other than those supplied by Atlantic Therapeutics Ltd. Never connect any computer cables to the controller.
- Avoid areas of injury or restricted movement (e.g. fractures or sprains).
- Keep the product away from children due to the risk of strangulation, and the risk of small parts which may be broken off and cause a choking hazard.
- For all users, including those with impaired sensory capabilities, the controller should not be handled at any time during the charge cycle. If the charge cycle is terminated early, please do not handle the controller for at least 1 hour, as surfaces may be hot.
- Once the battery is fully charged, let the unit cool down for a couple of minutes as the battery heats up during charging.
- Do not touch the electrodes, connectors or cable pins with your fingers while the stimulation session is in progress
- The Shorts or treatment cable must not be connected to other objects.
- This product can deliver current densities in excess of 2mA/cm² when used at a high intensity. The size, shape, and type of electrodes may affect the safety and effectiveness of electrical stimulation.
- The Shorts or treatment cable must not be connected to other objects.
- You should only run one full treatment session per day to reduce the potential for muscle fatigue. An effective treatment should not cause undue discomfort.

PRECAUTIONS:

- Although compliant with applicable EMC requirements, this device may still interfere with more sensitive equipment, please move away or switch off.
- It is advised that you empty your bladder before starting a treatment session.
- Always check the product cables before use and do not use if they are damaged in any way.
- Caution should be used over the menstruating uterus.
- Keep this device out of the reach of pets, pests and children.

Please wait before using your INNOVO® until:

- At least six months after the birth of your baby or until you stop lactating.
- The heavy days of your menstruation period have finished.

Please consult with your prescribing practitioner before using INNOVO® if:

- You have any serious illness or injury not mentioned in this guide.
- You have an underlying neurological condition e.g. Multiple Sclerosis or Parkinson's disease
- You have a tendency to hemorrhage, following acute trauma or fracture.
- You take insulin for diabetes.
- You suffer from muscle or joint problems, including pelvic pain or a paravaginal defect.
- You suffer with suspected or diagnosed heart disease.
- You underwent a recent surgical procedure and where muscle contraction may disrupt the healing process.

To reposition the Shorts or Electrodes during a session:

- To reposition the Shorts during a session: Always pause the device before adjusting the position of the Shorts

ADVERSE REACTIONS:

- A small number of isolated skin irritations have been reported by people using muscle stimulation devices, including skin allergies, burns under the electrodes, a prolonged reddening of the skin and acne.
- **NOTE** that some reddening of the skin can appear under the electrode during and for a short time after a session.
- If allergic reactions last for several hours, you should stop using the product as this product may not be suitable for you.
- On very rare occasions, first-time users of transcutaneous stimulation have reported feeling light-headed or faint. In this event, we recommend that you switch off the controller and sit down until you are feeling better. You may then resume the treatment, however it is recommended that you use position B or C to begin with.
- Use of muscle stimulation devices at a high intensity can be associated with minor and transient discomfort. Ensure you use your INNOVO® at a level which is comfortable to you. An effective treatment should not cause undue discomfort.
- You should stop using the product and should consult your prescribing practitioner if you experience adverse reactions or if in any doubt about using the device.

IMPORTANT:

- The safety information provided in these Instructions For Use must be followed.
- The product should only be set-up and operated according to these Instructions For Use.
- Atlantic Therapeutics will not accept responsibility if the guidelines and instructions supplied with this product are not followed.

VALIDITY

These instructions For Use are intended for the correct application of the product. The information and technical data contained in this document relates to the INNOVO® muscle stimulator provided with Instructions For Use which is proprietary to Atlantic Therapeutics Ltd. and may be used and disseminated only for the purpose of and to the exact specifically authorized writing by the company. Each INNOVO® controller is assigned a unique serial number, which is located on the back of the controller.

RESTRICTIONS

The sale and/or operation of this equipment is subject to law in various countries. Compliance with the legislation rests with the importer, dealer or user of the equipment as appropriate.

Contact your local careline (see page cover) if:

- Your product is not working correctly. Do not use it in the meantime.

INNOVO® Treatment Plan

INNOVO® is used to treat Stress Incontinence. Below is the suggested weekly treatment plan. It is recommended that you complete just one treatment session per day and that you follow a minimum 12-week treatment plan to achieve optimum results.

Stress Incontinence

Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Prog. 1	Prog. 1	Prog. 1	Prog. 1	Prog. 1	Off	Off

After completing your treatment plan


After completing your 12 week Innovotherapy plan, you should continue to use INNOVO® to help maintain a healthy pelvic floor. How frequently you should use INNOVO® will depend on your specific condition and improvement rate. As a general guide, we would recommend you continue to use INNOVO® 1-2 times a week. This will ensure that the results achieved during the 12 week program will be sustained. One of the major advantages of using INNOVO® is that you will have an improved awareness of your pelvic floor, so you can also maintain a healthy pelvic floor by doing daily pelvic floor exercises.

Visit www.myinnovo.com to understand how to get the best from your INNOVO® product.

Important: You should only run one full treatment session per day to reduce the potential for muscle fatigue. A treatment session is defined as 30 minutes of treatment, even if this is broken down into shorter 10- or 15- minute sessions.

INNOVO® PACK CONTENTS

When you open your INNOVO® you will find the following:

1. INNOVO® Controller:

This generates comfortable stimulation which is sent via the electrodes within the INNOVO® Shorts to activate your muscles. This controller has rechargeable batteries and takes approximately 5 hours to fully charge from flat.

2. Charger:

Attach your country-specific adapter, connect to your controller and plug into an outlet socket to charge the battery.

3. INNOVO® Shorts:

The INNOVO® Shorts are worn like a pair of cycling shorts. There are two connection points to link your INNOVO® cable to the shorts. The INNOVO® Shorts contain eight integrated electrodes.

4. Spray:

Applying this safe spray ensures a reliable and comfortable connection between your skin and the eight integrated electrodes on the INNOVO® Shorts.

5. INNOVO® Cable:

This connects your INNOVO® controller to the INNOVO® Shorts.

6. Neck Strap:

The neck strap can be clipped onto your INNOVO® controller to allow it to be worn around the neck.

7. Instructions For Use:

Read these easy to follow instructions to understand how to set up and use INNOVO®.

8. Quick Start Guide:

This is an additional guide to assist you when setting up and using INNOVO®.

INNOVO® CONTROLS

1. On/Off Button (⏻)

Press and hold this button for 2 seconds to turn the controller on or off. During your treatment session, pressing this briefly will pause or resume the treatment.

Note: For your safety, you cannot switch the controller on while the battery is charging.

2. Increase Intensity Button (+)

Press this button to increase the stimulation intensity during your treatment session.

3. Decrease Intensity Button (-)

Press this button to decrease the stimulation intensity during your treatment session.

4. Program Button (P)

The Program Button is used to select the programme you wish to use:
P1 - Stress Program

5. Information Button (i)


Press this button to see various items of information and treatment details from the sessions you are running or have completed (see page 17 for more information).

6. Scroll Up Button (△)

This button is used to scroll up through the information relating to individual treatment sessions.

7. Scroll Down Button (▽)

This button is used to scroll down through the information relating to individual treatment sessions.

INNOVO® SOCKETS

1. Battery Charger Socket:

Connect the charger to this socket. Your INNOVO® controller is fully charged when all 5 bars of the battery icon on the display are visible.

2. Cables Socket:

Connect the cable to this socket.

3. Socket Isolators:

These sliders ensure that INNOVO® cannot be connected to you and to an electrical power source at the same time. To charge the controller, these must be slid inwards to cover the Cable Socket (Fig. 1). To run a treatment Program, they must be slid outwards to cover the Battery Charger Socket and PC Connector Socket (Fig. 2).

4. PC Connector Socket:

The software option is not available. Do not insert anything into this socket.


Fig. 1



Fig. 2

Notes:

- NEVER connect cables or wires from other equipment such as computers, games, cameras etc. to the controller.
- NEVER connect the cables provided with INNOVO® to any other equipment.

CONTROLLER DISPLAY

Full screen with all icons displayed

Battery Power indicator:
Shows the remaining battery power.

Warning indicator:
Indicates that there is a problem with your product - with the cables, shorts, integrated electrodes or the controller itself.

Keylock symbol:
Appears when Keylock has been activated to indicate that the chosen program or intensity level cannot be accidentally changed.

Mute symbol:
Indicates that the sound is off.

Treatment sessions:
Number of treatment sessions completed.

Treatment timer:
Displayed when viewing the total treatment time, or the duration of a treatment session.

Maximum intensity:
This is displayed when viewing the highest intensity reached to date, or in a particular treatment session.

Average intensity indicator:
This is displayed when viewing the average of all treatment sessions run to date.

Pause symbol:
Appears whenever a treatment session has been paused.

Program indicator:
Shows which program is currently selected. Type 208 has one program installed.

Individual treatment session indicator:
Displays each session (1-90).

Intensity indicator:
Appears during each contraction phase.

Intensity/time indicator:
Displays intensity level from each channel and treatment time remaining.

mA

Controller indicator:
Shows what type of information is being displayed by the intensity/time indicator.

hr min min sec

Electrode contact indicator:
Indicates electrodes are in poor contact with the skin. See page 20 for more details.

STEP 1: CHARGER ASSEMBLY

1. Align the adapter with the charger, with adapter slots facing towards the charger pins.
2. Slide the adapter into place on the charger until it clicks into place
3. Ensure the charger is firmly assembled before plugging the adapter into an electrical outlet.
4. To release the adapter from the charger, press the "PUSH" button on the charger and slide the adapter off.

Fig. 3


Fig. 4


Fig. 5



Fig. 6

STEP 2: CHARGING AND GETTING STARTED WITH INNOVO®

1. Before you start a treatment session, please ensure that the battery is sufficiently charged. If the battery charge is low, connect the charger to the controller and plug it into an electrical socket. You will need to move the sliders inwards to enable you to do this (Fig. 7). You should charge the battery for approximately 5 hours from flat. It is fully charged when all 5 bars of the battery symbol () on the display screen are visible. A full charge should be sufficient to complete approximately 5 sessions.
2. Make sure your skin is clean and free from oils, lotions, moisturisers etc as these may prevent the Shorts from functioning correctly.
3. We recommend that you wear high-thigh underwear which leaves your buttocks exposed, as the electrodes need to be in full contact with your skin to deliver the INNOVO® treatment programs successfully (Fig. 8). Alternatively, you can also use INNOVO® without underwear. It is advised that you empty your bladder before starting a treatment session.
4. Follow these instructions and Quick Start Guide carefully to ensure you set up and apply INNOVO® correctly. We strongly recommend that you review the information on myinnovo.com.

Fig. 7



Fig. 8


STEP 3: PUTTING ON THE SHORTS

Follow the steps in the Quick Start Guide.


1. Unpack all contents from box. Turn the shorts inside out. Before applying the spray, we recommend that you hang the shorts up.

Fig. 9

2. Remove the lid from the bottle of spray. Holding the bottle approx. 15-20 cms from the shorts, spray each electrode 2 - 3 times to ensure an even coverage. Spray each of the 8 integrated electrodes. Turn the shorts the right way around, so that the electrodes are on the inside, before putting them on.


Fig. 10

Important: Put on the shorts as soon as you have sprayed the electrodes.


3. Pull on the shorts until the waistband is sitting comfortably above your hips. The red and blue connectors should be to the front, with the wash label to the rear, as per Fig 11.
4. Ensure the whole surface of each electrode is in full contact with your skin.

Fig. 11

5. When the shorts are correctly positioned, the electrodes will be in place as shown in Fig 12.

Fig. 12

STEP 4: GET INTO THE CORRECT POSITION FOR YOU


Before you turn the controller on, get into the position that is most comfortable for you. You should remain in position until your treatment session is finished.

Position A:

To optimise the results from a treatment session, it is highly recommended that you remain standing with your feet flat on the floor, shoulder width apart. You should tilt slightly forward and stand near a surface to lean on if required (Fig. 13).

If standing, you may move a little during the treatment to help you remain comfortable, but you should avoid moving too much as this may interfere with the effectiveness of the treatment. However, if you wish to move to a different location, you should pause the session, move to the new location and then resume the session.


Fig. 13

Position B:

Where standing for 30 minutes is difficult, you can carry out the treatment session in the crook lying position, which is lying on your back with your legs bent at the knee, shoulder-width apart and with your feet flat (Fig. 14). To help support your back while crook lying, you can relax your quad muscles by just letting the weight of your thighs sink gently outwards into your hip sockets. This will help you to maintain a small anterior tilt of your pelvis, which has the effect of slightly arching your lower back.


Fig. 14

Position C:

The crook lying position can also be carried out with your back slightly raised (Fig. 15). Ensure the legs are bent at the knee and that your feet are flat.

Fig. 15

STEP 5: USING INNOVO®

5.1 TURNING ON THE CONTROLLER AND STARTING A SESSION

1. Turn on your INNOVO® controller by pressing the On/Off button (⏻).
2. P1 will be displayed to indicate that the program has been loaded correctly (Fig 16).
3. Start the treatment session by pressing the Increase Intensity button (+). You will see the intensity indicator rise on screen (Fig. 17). The Intensity buttons increase and decrease the intensity in single increments and the range goes from 0-120. If you keep the Intensity button pressed it will increase/decrease continuously.
4. Once you begin to increase the intensity level the treatment timer starts to count down from 30 minutes (Fig. 18).


Fig. 16

Fig. 17


Fig. 18

5.2 INTENSITY LEVELS AND WHAT TO EXPECT DURING A TREATMENT SESSION

1. You should increase the level of intensity gradually during your first treatment session. This gives you time to become accustomed to the sensation. At first you will begin to feel the stimulation in the buttock area and as you continue to increase the intensity, this stimulation will move forward between the legs, giving you a 'lifting' sensation throughout the pelvic floor (Fig. 19). You may achieve this pelvic floor lift at an intensity level of approximately 30-40.
2. With time (usually within 1 week), you should find that you are better able to tolerate higher intensities and are well on your way to attaining the recommended intensity level. Once you can feel a strong stimulation throughout the full pelvic floor you have reached a sufficient level. However, you should never exceed your comfort level. Effective stimulation should be strong but never uncomfortable. You should try to minimize movement during the therapy session

Fig. 19


Important:

You will not feel stimulation directly under the electrodes on the front or sides of your legs.

3. If you wish, you can now clip the neck strap onto the controller and hang it around your neck (Fig. 20). For more information refer to page 16.

Fig. 20

5.3 HOW TO PAUSE A TREATMENT SESSION

It is recommended to pause the stimulation if you are changing body position during a treatment session. If you wish to pause a treatment session, simply press the On/Off button (⏻) briefly. The stimulation will stop and the display will indicate that the controller is paused (Fig. 21). To resume the session, briefly press the On/Off button again. The display will return to the timer countdown and the stimulation will begin again. However, if the session is paused for 2 minutes or more, the intensity will be halved when the session is resumed. This safety feature ensures that stimulation does not feel uncomfortable when you resume treatment.

You may switch off your controller at any time during a treatment session by pressing the On/Off button for two seconds. The controller will power off and your session will terminate.

5.4 FINISHING A TREATMENT SESSION

1. The treatment session is completed when the counter reaches zero. The display will flash "00:00" (Fig. 22), stimulation is reduced to zero and the controller will beep to indicate the end of the session.
2. Once the session is complete, switch the controller off by pressing the On/Off button (⏻) for two seconds.
3. Disconnect the cables from the controller and the shorts (Fig. 23).
4. If you are using the neck strap, you can remove the controller from around your neck.

Note: If you do not switch the controller off, it will do so automatically after 10 seconds. However, if another button is pressed within this 10-second period, the controller will remain on for 2 minutes, at which point it will switch off automatically.

5.5 REMOVING THE SHORTS

1. Remove the shorts by pulling the waistband down. The shorts will be inside out when they are removed.

2. Fold the shorts in half, ensuring that electrodes remain flat when packed away with your INNOVO®.

Fig. 21

Fig. 22

Fig. 23



Fig. 24

Fig. 25

ADDITIONAL PRODUCT FUNCTIONS

Using the Neck Strap

1. When using the neck strap, attach the clip to the connector at the top of the controller (Fig. 26). Ensure the controller is attached with the screen facing away from your body.
2. The neck strap is equipped with a safety-clip which pops open if there is a sudden, forceful pull on the neck strap. If this happens, align the two halves of the clip (Fig. 27) and simply push them back together until the oval button clicks into place. If it does not “click” together, you may have one side of the clip back-to-front. Turn it around and try again.

Fig. 26



Fig. 27

Using the Keylock Function during your treatment session

When you have set the intensity to a level you find comfortable, press and hold the Information Button () for two seconds to activate the Keylock Function (Fig. 28). This “locks” the Increase and Decrease Intensity buttons, so the intensity cannot be changed accidentally. The other buttons will still function as normal. If you wish to increase or decrease the intensity, you can deactivate the Keylock Function by again pressing the Information button for two seconds. This function only remains active during the session in which it is activated. Once the treatment session is complete, the Keylock Function deactivates automatically.

Note: Keylock cannot be activated while the controller is paused.

Fig. 28




Fig. 29

Mute Function

If you want to switch off your controller’s sound effects, press and hold both the Up () and Down () Scroll buttons for 1 second (Fig. 29). This will activate the Mute Function. Deactivate the Mute Function by again pressing both the Up and Down Scroll buttons for 1 second.

Note: The Mute Function is overruled if there is an electrode contact problem (see page 20) or if there is a controller error message (see page 10).

Reset Data


The INNOVO® controller can store the data from 90 treatment sessions to enable you to review your treatment progress. Once this is reviewed, the controller can be reset and the data cleared by pressing and holding the Decrease Intensity button () and the Information button () for 4 seconds (Fig. 30). After 2 seconds, the Total Treatment Time will begin to flash on screen and when this resets to 00:00, your data has been cleared. This function will not work while a treatment session is in progress. If you don’t reset the data after 90 sessions, you can still continue to use your product: in this instance the data from the 91st and subsequent sessions will be written over the data for session #1, #2 etc.

Fig. 30

INFORMATION BUTTON

The Information button () is used to display information about your current and previous treatment sessions. Each time you press the button, a new element of your treatment progress is displayed on screen. This information remains on screen for 3 seconds and will then revert back to the countdown timer unless you press the Information button again. The information is displayed in the following sequence:

Fig. 31

1 Press: Displays your current intensity level. Fig. 31.

Fig. 32

2 Presses: Displays the maximum intensity reached during any session to date (up to 90 sessions). Fig. 32.

Fig. 33

3 Presses: Displays the average intensity of all sessions completed to date (up to 90 sessions). Fig. 33.

Fig. 34

4 Presses: Displays the total number of sessions completed to date, up to 2000 sessions. This includes partial sessions, where the full 30 minutes were not completed. Fig. 34.

Fig. 35

5 Presses: Displays the total amount of treatment time completed to date, up to 199 hours and 59 minutes. Fig. 35.

Fig. 36

6 Presses: Displays the treatment code for this particular controller. You should check that your screen displays U208. If another code is displayed, you should contact Atlantic Therapeutics for advice. Fig. 36.

Pressing the Information button a 7th time will bring you back to the treatment timer screen.

IN-PROGRAM INFORMATION

Your INNOVO® controller stores three pieces of data from every session you run. This allows you to monitor your treatment progress if you wish:

- a) The maximum intensity for each session.
- b) The average intensity for each session.
- c) The duration of each session.

This information can be reviewed in 2 different ways - either by reviewing all 3 parameters in a single session or by reviewing one specific parameter over a number of sessions.

Note: This information cannot be accessed while a program is running.

Review by Session

1. Access the list of treatment sessions by pressing either the Up or Down Scroll buttons (Fig. 37).
 - By pressing the Scroll Up button (\triangle) you will access the list of sessions from the oldest to the most recent, starting with the first treatment session you completed, then the second, then the third etc.
 - By pressing the Scroll Down button (∇) you will access your most recently completed treatment session first, then the session completed before that, then before that again, etc.
2. When you have highlighted the session you wish to review, the maximum intensity level for that session appears on screen (Fig. 38). This information remains on screen for 3 seconds and will then revert back to the start screen, unless you press the Information button again.
3. Press the Information button once to display the average intensity of the selected session (Fig. 39).
4. Press the Information button twice to display the duration of the selected session (Fig. 38).
5. Press the Information button three times to display the program used for the selected session (Fig. 40).

Review by Parameter

1. Access the list of treatment sessions by pressing either the Up or Down Scroll button.
 - By pressing the Scroll Up button you will access the list of sessions from the oldest to the most recent, starting with the first treatment session you completed, then the second, then the third etc.
 - By pressing the Scroll Down button you will access your most recently completed treatment session first, then the session completed before that, then before that again, etc.
2. When you have highlighted the session you wish to review, press the Information button until you find the parameter you wish to review. For example, Fig. 42 shows the maximum intensity in session #1.
3. Next, press the Up Scroll button to display the same parameter in the next session (i.e. the maximum intensity in sessions #2, #3 etc.- Fig. 43 & Fig. 44). Or press the Down Scroll button to display the same parameter in the previous session. This information remains on screen for 3 seconds and will then revert back to the start screen, unless you press the Scroll Buttons or the Information button again.

Note: If you wish to review the average intensity for each session instead of the maximum intensity; when you have highlighted the session you wish to review, press the Information button twice to see the average intensity for session #1. Then use the Scroll Buttons to see the average intensity in session #2, #3 etc. Press the Information button three times to see the duration of each session you wish to review.

Fig.37

Fig. 38

Fig. 39

Fig.40

Fig. 41

Fig. 42

Fig. 43


Fig. 44

PROBLEM SOLVING GUIDE

If an issue arises with your INNOVO®, check the table below to find out how to fix it.

Description of the Problem	Possible Cause	How to Fix the Problem
The display does not come on and there is no signal	The battery is out of charge	Recharge the battery
Battery symbol () flashing	The battery is low	Recharge the battery
Screen goes blank during a session	The battery is out of charge	Recharge the battery
	There is a fault with the controller	Switch off the controller and then switch it back on again. If the error persists, visit myinnovo.com for more information
	The screen may be damaged	If the screen is damaged, do not use the device. Visit myinnovo.com for contact information
The stimulation feels uncomfortable	You may not be reaching a suitable intensity level	Ensure that you are reaching a suitable intensity level to activate the pelvic floor. See pages 14 for more details on intensity levels and what to expect
	The electrodes may not be in full contact with the skin	Ensure all electrodes are flat against the skin and that there are no folds or creases in the shorts
	You may not have applied sufficient spray	Remove the shorts and apply 2 to 3 sprays per electrode as on page 12. Put the Shorts on again and Reconnect
The warning symbol () is flashing, the controller is beeping, and E'XXX' appears on screen **	There is a fault with the controller	Switch off the controller and then switch it back on again. If the error persists, visit myinnovo.com for more information
P1 and the warning symbol are flashing	There is a poor connection between the INNOVO® controller and the INNOVO® cable	Remove the INNOVO® cable from the controller and reconnect, pushing firmly to ensure a good connection. If the error persists, visit myinnovo.com for more information
The warning symbol () is flashing, along with numbers from 1-8	The Electrode Detection Function has been enabled	Refer to ELECTRODE DETECTION FUNCTION section of manual. Page 20

* Remember that your INNOVO® Shorts will become worn over time and will need to be replaced. If you can see visible signs of wear on your shorts, visit myinnovo.com or an authorised reseller to reorder.

** XXX =001, 002, 003 etc.

ELECTRODE DETECTION FUNCTION

During a treatment session, the controller continuously monitors its connection with the electrodes and can detect if they are in proper contact with both the shorts and your skin. If the controller detects poor contact, it will alert you to the problem and you can take the necessary steps to correct the issue.

If an electrode issue is detected:

- The stimulation will stop until the problem is fixed.
- The warning symbol (⚠) will flash on screen until the problem is fixed.
- The controller will beep continuously until you switch it off and fix the problem.
- The electrode contact symbol (📶) will flash on screen along with numbers from 1-8 (Fig. 45). These numbers relate to the specific electrode(s) causing the issue on the shorts. See Fig. 46 to see where each numbered electrode is positioned on your body (there may be more than one).

Possible causes:

- Poor contact between the electrode(s) and the skin.
- Poor electrode condition (i.e. with damaged surface following too many sessions).
- The shorts are in the incorrect position.

How to fix an electrode detection fault:

1. Reconnect the cable to the shorts, ensuring they are fully clicked into place and begin the treatment again.
2. If issue persists, switch off the controller by pressing the On/Off button (⏻) for 2 seconds and disconnect the cable from the shorts.
3. Remove the shorts and reapply the spray as seen on page 12.
4. Reapply the shorts so that the waistband is sitting comfortably above your hips. Ensure that the whole surface of each electrode is pressed firmly against the skin.

Fig. 45



Fig. 46

Note: If wearing underwear, ensure that the full surface of all electrodes are in contact with the skin. If this is not possible change or remove underwear.

Note: If the screen still indicates the same issue, you may have a faulty cable or shorts and should visit myinnovo.com for further information.

CARING FOR YOUR INNOVO®

Cleaning Instructions for your INNOVO® Controller

Never let your INNOVO® controller get wet. However, you can clean your INNOVO® controller using a soft cloth, lightly dampened in soapy water. Do not use your INNOVO® controller if it gets wet.

Do not use detergents, alcohol, spray aerosols or strong solvents on your controller.

Access to the interior of the controller is not required for maintenance purposes. If your controller is damaged, or there is any indication of loose or dislodged internal components, you should not use it but should return it to Atlantic Therapeutics or your local distributor for replacement or repair. Repairs, service and modifications may not be carried out by anyone other than qualified service personnel authorised by Atlantic Therapeutics. With proper care, you should expect your INNOVO® controller to maintain its safety and performance characteristics for 2 years.

If you drop your INNOVO® controller and there is any visible damage or any indication of loose/dislodged internal components as a result, please return the controller to Atlantic Therapeutics.

Washing Instructions for the Shorts

Important: Detach the controller, cables and connectors from the shorts before cleaning. Always follow the instructions on the label when washing the shorts. Avoid immersing the connector ports in water.

You can wash the shorts in lukewarm water using a small amount of handwash detergent. Do not machine wash. Do not use softeners. Do not allow the connector ports on the shorts to become wet. Remove immediately after washing.

Do not use any bleach when washing your shorts.

Do not dry clean your shorts.

Do not tumble dry your shorts. The shorts should be hung up to dry. Do not dry over anything hot (e.g. a radiator or any other direct heat source). Do not wring the shorts out. Ensure the shorts are completely dry before next use. Do not fold while drying.

The shorts should not be ironed.

Note: If the shorts become visibly damaged, contact Atlantic Therapeutics to arrange for a replacement.

Garment Materials

Legs: 91% Polyamide, 9% Elastane. Waistband: 80% Polyamide, 20% Elastane.

REPLACING THE BATTERY

The battery can be charged and discharged hundreds of times, but it will eventually wear out. To ensure optimum battery performance you should use your controller regularly, e.g. perform 1 session a day, 5 days per week. New batteries are available from the website myinnovo.com

To install a new battery, follow the steps below:

1. Turn off the controller and unplug all connected cables.
2. There is a small rubber plug on the back of the controller. Remove this and unscrew the back of the controller (Fig. 47).
3. Carefully unplug the battery connector from its socket and remove the battery. (Fig. 48).
4. Insert the new battery pack with the connector cable towards the bottom of the controller (Fig. 49).
5. Plug the battery connector into its socket, ensuring the 2 silver contacts are facing away from the battery pack (Fig. 50).
6. Replace the back of the controller by screwing it back into place. Re-insert the small rubber plug.

Battery Power

Your INNOVO® controller is powered by a rechargeable battery. The full performance of a new battery is achieved only after two or three complete charge cycles (fully charged and discharged).

The battery symbol is displayed in the upper-right corner of the screen. When the battery is fully charged, the battery symbol will have 5 bars. If the battery depletes completely, the controller will beep and the battery symbol will begin flashing. Use the battery only for its intended purpose. Never use any charger or battery which is damaged. If the controller has not been used for a significant period of time and the battery is completely discharged, it may take a few minutes before the charging indicator appears on the display. Charging may take longer than normal in this instance.

Unplug the charger from the electrical socket and the controller when it is not in use. Do not leave a fully charged battery connected to a charger, as overcharging may shorten its lifetime. If left unused, a fully charged battery will lose its charge over time. Leaving the battery in hot or cold places, such as in a closed car in summer or winter conditions, will reduce the capacity and lifetime of the battery. Always try to keep the battery between 15°C and 25°C (59°F and 77°F). A controller with a hot or cold battery may not work temporarily, even when the battery is fully charged. Battery performance is particularly limited in temperatures well below freezing. Do not dispose of batteries in a fire as they may explode. Batteries may also explode if damaged. Dispose of batteries according to local regulations. Please recycle when possible. Do not dispose as household waste.

This product is intended for continuous operation, i.e. operation in normal use, for an unlimited period of time, without the specified limits of temperature being exceeded.

Disposal of Batteries, Spray Bottles and Shorts

Used battery packs, spray bottles and shorts must never be disposed of in a fire, but in accordance with your country's national or state laws governing the disposal/recycling of such items.

Fig. 47

Fig. 48



Fig. 49

Fig. 50

Silver Contacts

TECHNICAL SPECIFICATIONS

TREATMENT PROGRAM INFORMATION

Table 1: Treatment Parameters – Stress urinary incontinence

Prog	Frequency Hz	Pulse Width μS	Contraction (Sec.)	Relaxation (Sec.)	Treatment Time (Min.)
1	50	620	5	5	30

General Specifications

Product Type: 208

Classification: Internally powered equipment, Class II charger, and Type BF applied parts.

This controller is not suitable for use in an oxygen-rich environment.

Neuromuscular Electrical Stimulation (NMES)

NMES may be defined as the application of electrical stimulation of the peripheral nervous system to contract a muscle either through the direct activation of the motor neurons in the mixed peripheral nerve, or indirectly through reflex recruitment.

Waveform: Symmetrical bi-phasic square current waveform.

PCB Fuses:

- F1: 1 Amp, fast-acting, max. Voltage Rating (V): 50, Interrupting Rating (AC/DC): 50A@50VDC/50A@VAC.
- F2: 1 Amp, fast-acting, max. Voltage Rating (V): 50, Interrupting Rating (AC/DC): 50A@50VDC/50A@VAC.


Electrical Specifications

The output current reading displayed on the screen indicates the peak current for a 500Ω resistive load only. The actual current value delivered may be lower due to varying impedance levels within users. See table on next page for measured values at different resistive values.

Program 1 (P1)

Rated Outputs – Voltage/Currents for Stress urinary incontinence

Parameters	500 Ω	1K Ω	1.5K Ω
Output RMSA	30.3mA	24.4mA	16.6mA
Output RMSV	15.1V	24.4V	25.0V
Output Frequency	50HZ	50HZ	50HZ
DC Component	0V	0V	0V
Pulse Width t1	620 μS	620 μS	620 μS
Max Output Current (peak)	120mA	100mA	70mA

INNOVO® produces a symmetric bi-phasic constant current waveform. The frequency = $1/T$, the phase durations t_1 and t_2 are the same, the overall pulse duration is the sum of $t_1+t_2+t_3$. i = current intensity.

Environmental Specifications:

Operating Range:	Temperature:	41°F to 104°F (5°C to 40°C)
	Humidity:	15 to 93% RH (non-condensing)
	Atmospheric Pressure:	70 to 106kPa=
Transport & Storage Range:	Temperature:	-13°F to 158°F (-25°C to +70°C)
	Humidity:	10 to 93% RH (non-condensing)
	Atmospheric Pressure:	50-106kPa

Time required to reach operating temperature is 1 hour. If required, always seek assistance to set up, use or maintain your INNOVO® from Atlantic Therapeutics. You should report any unexpected operation or event to Atlantic Therapeutics.

Accessories

Under no circumstances should anything other than Atlantic Therapeutics accessories be used with your INNOVO® controller (Type 208). Any others may not be compatible with your controller and could degrade the minimum safety levels. You can purchase all accessories from Atlantic Therapeutics:

INNOVO® Shorts	1005-08XX
Neck Strap	1220-0827
Cable	1200-0807
Battery Pack	3924-0312
Conductive Spray (50ML)	1500-08XX
Conductive Spray (250ML)	1500-08XX

Important:

Under no circumstances should anything other than the correct type of battery pack (rechargeable batteries 1 x 7.2V (6 x 1.2 AAA) volt (NiMH) DC battery pack) be used with your controller. These can be purchased from the Atlantic Therapeutics careline.

Power Supply (2504-0808-02)

UE Switching Power Adaptor UES12LCP-120100SPA

Input: 100-240V ~ /50-60Hz/500mA

Output: 12V DC / 1.0A

Important: Do not use any other chargers or power supplies. The power supply is included as part of this medical electrical equipment. To isolate from the mains, the charger/power supply plug should be removed from the socket.

Expected Service Life: 2 years.

At the end of the product lifecycle, do not throw this product or the batteries into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Some product materials can be re-used if you bring them to a recycling point. By re-using some parts or raw materials from used products you make an important contribution to the protection of the environment. Please contact your local authorities if you need more information about collection points in your area.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to humans.

Description of the symbols

There are a number of technical markings on your controller and Power Supply Unit. These can be explained as follows:

INNOVO® Symbols

R_X Only **CAUTION: IN THE UNITED STATES OF AMERICA FEDERAL LAW RESTRICTS THE DEVICE TO SALE OR USE BY, OR ON THE ORDER OF A PHYSICIAN OR OTHER PRACTITIONER LICENSED BY THE LAWS OF THE STATE IN WHICH HE/SHE PRACTICES.**

The controller and shorts are manufactured for Atlantic Therapeutics Ltd., Parkmore Business Park West, Galway, H91 NHT7, Ireland.

The controller requires 1 x 7.2V (6 x 1.2 AAA) (NiMH) DC battery pack; DC is indicated by the symbol:

Power (**P**): Maximum power output measured in Watts (W) into a 500Ω load.

Frequency (**F**): Number of pulses output by the controller per second, measured in Hertz (Hz).

IP22: Ingress Protection Rating

VXX: Software version.

This symbol means type BF applied parts (Shorts)

This symbol means “Attention, read the accompanying documents”

This symbol means “Keep dry”

SGS US Listed Mark.

The product packaging (box, sleeve and insert) can be recycled

This symbol means “Use by”

SN stands for ‘serial number’. The serial number specific to this product appears on a label on the back of the controller. The letter preceding the number indicates the year of manufacture, where “X” denotes 2018, “Y” denotes 2019 etc..

Power Supply Symbols

Complies with the requirements of Low Voltage Directive 2014/35/EU

Class II Equipment

Indoor use only

International Efficiency Mark rating of 6

China Compulsory Certificate mark

UL Recognized Component Mark (Canada and the United States)

China RoHS Mark (indicates an Environment Friendly Use Period of 10 years)

Indicates conformity to WEE directive 2012/19/EU

TÜV-Mark (indicates conformity to IEC 60601-1:2005 +A1 and IEC 60601-1-6:2010 +A1).

Shorts Symbols

This symbol is to indicate conformity to the requirements of the Medical Device Directive (93/42/EEC). #### is the number of the notified body.

The Shorts’ batch number is represented on its packaging and on the shorts (rear of the socket), by the number corresponding with the **LOT** symbol.

Spray Symbols

Manufacturer’s catalogue number

Keep away from sunlight

Non-sterile

Do not use if package is damaged

PHYSICIAN INFORMATION

CLINICAL STUDIES

Pilot Study

Fourteen (14) patients with stress urinary incontinence (SUI) completed a 12-week training program using a neuromuscular electrical nerve stimulation (NMES) device which was an earlier version of the INNOVO® device. The results showed a reduction in the mean urine leakage on the 1-hour and 24-hour pad weight tests from baseline to week 12, and improvements in quality of life and the number of incontinence episodes per day.

Randomized Clinical Studies

Two randomized clinical trials were performed. The first was a sham-controlled trial conducted in Germany of fifty women with mixed incontinence that was prematurely terminated due to the commencement of the second trial. The second was an active-controlled, US, randomized, non-inferiority trial that included 180 adult women with stress urinary incontinence.

Trial 1: German Sham-Controlled Trial

Fifty (50) female German subjects with predominantly stress UI were randomized into two groups; many subjects had mixed incontinence and some subjects were enrolled who had severe urge incontinence as part of an overall severe mixed incontinence. All subjects had previously failed a 6-week volitional pelvic floor muscle training program or an equivalent lifestyle and exercise program. The mean ages were 51-53 years across treatment groups.

Subjects in the active group (n=24) used an active stimulation device, while subjects in the control group (n=26) received treatment with a modified device. The modified device was programmed to produce a lower dose of stimulation with a sensory response. The active device used in this study was an earlier version of the Innovo device.

Treatment was comprised of a 30-minute session of stimulation, in a standing position, 5 days per week, for 12 weeks.

At week 12, the mean reduction from baseline in urine leakage using the one-hour pad weight test was 22.55 grams in the active group versus 17.65 grams in the control group. A greater improvement from baseline was also observed for the I-QOL in the active group. The one-hour provocative pad weight test showed improvement in the percentage of subjects who experienced at least a 50% reduction in pad weight for the treatment over control at 12 weeks (75% treatment versus 47% control). The 24-hour pad weight measures for treatment and control at 12 weeks showed no significant difference for the same outcome measure (44% treatment versus 41% control). However, there was a greater improvement in the 24-hour pad weight results for the active group versus the control group at week 12 where the mean reduction in urine leakage was as 26.54 gms for the active device and 6.6 gms for the control.

Results of this randomized, controlled, double-blind, clinical study showed greater magnitude of change in the 1-hour pad test for subjects using the active device compared to the control device after 12 weeks of treatment. The device had few side effects.

References

1. Brown CA. Literature review: Does neuromuscular electrical stimulation increase pelvic floor muscle strength in women with urinary incontinence with an ineffective pelvic floor contraction? *Journal of the Association of Chartered Physiotherapists in Women's Health*, Spring 2014, 114, 56–62.
2. Mørkved S, Bø K, Fjørtoft T. Effect of adding biofeedback to pelvic floor muscle training to treat urodynamic stress incontinence. *Obstet Gynecol*. 2002;4:730–739.
3. Alves PG, Nunes FR, Guirro EC. Comparison between two different neuromuscular electrical stimulation protocols for the treatment of female stress urinary incontinence: a randomized controlled trial. *Rev Bras Fisioter*. 2011;5:393–398.

Trial 2: Pivotal USA Trial

Aims of study: The aim of this study was to compare the efficacy and safety of the InnoVO® external electrical muscle stimulation device with an FDA-cleared intravaginal stimulation device (comparator device) for the treatment of SUI in women.

Study design: This was a prospective, randomized, single-blind, multicenter, non-inferiority study performed at 12 sites in the USA. The study sample size was 180 patients.

Women with stress predominant UI whose condition had not improved using pelvic floor muscle training were randomized to undergo treatment with either the INNOVO® or comparator device for 12 weeks. To be enrolled in the study, women with mixed UI could not have had severe urge incontinence.

Eighty-nine (89) women were randomized into the INNOVO® group and 91 women to the comparator device group.

Baseline incontinence characteristics were similar between the groups. The mean age for subjects included in the study was 45.9 years and 47.8 years for the INNOVO® and comparator groups, respectively. The majority of subjects were premenopausal (61.1%). The mean BMI was 27.52 kg/m².

Treatment was administered by the subjects at home using the device in accordance with the relevant instructions for use, which specified that the INNOVO® device was used for 30 minutes once daily for 5 days/week, and the comparator device was used consistent with its labeling.

Study endpoints: The primary endpoint was the proportion of subjects who achieved “significant improvement” (>50% reduction in pad weight from baseline) in the provocative pad weight test at 12 weeks.

Key secondary efficacy endpoints included:

- The mean change from baseline to week 12 for urine leakage in the provocative pad weight test,
- 24-hour pad weight test,
- Number of incontinence episodes/day,
- Incontinence Quality of Life questionnaire (I-QOL) score,
- Number of pads used/day,
- The proportion of patients achieving dryness (<1 g on the provocative pad weight test) at week 12.

Safety was also assessed.

Results:

For the primary effectiveness endpoint of the proportion of patients that had significant (at least 50%) improvement in the 1-hour provocative pad weight test (Table 1), the InnoVO treatment arm failed to demonstrate non-inferiority against the active comparator device based on the pre-specified 5% non-inferiority margin for the primary outcome. Overall, the primary outcome of 56.3 percent of the patients experiencing a clinically significant (at least 50% reduction in one-hour provocative pad weight) is clinically meaningful. The InnoVO performed numerically better than the active comparator for two of the seven secondary endpoints, the Global Impression of Improvement and the 24-hour pad weight test. The INNOVO® demonstrated improvements from baseline in mean urine leakage in the provocative pad weight test, 24-hour pad weight test, number of incontinence episodes, pads used per day, and I-QOL score at week 12 (Table 1).

Primary Endpoint	INNOVO® N=89	Comparator Device N=91
Proportion of Patients who attained significant improvement (50% reduction) in Provocative Pad Test between Baseline and Week 12		
ITT Population/Multiple Imputation	56.3%	63.0%
ITT Population/Last Observation Carried Forward (LOCF)	58.4%	61.5%
Key Secondary Endpoints Outcomes for INNOVO®: Change from Baseline at Week 12 (LOCF)		
Provocative pad weight test (mean) (grams)	-8.48	-9.66
24-hour pad weight test (mean) (grams)	-13.07	-9.89
Number of incontinence episodes/day (mean)	-1.24	-1.43
Incontinence Quality of Life Questionnaire (I-QOL) (mean)	13.41	15.42
Number of pads used/day (mean)	-0.30	-0.44
Other Secondary Endpoints: Change from Baseline at Week 12 (Observed Case)		
Patients who achieved Dryness (<1g leakage on the provocative pad weight test) n(%)	17 (19.1%)	29 (31.9%)
Proportion of patients reporting improvement in SUI Symptoms (Global Impression of Improvement)	70.7%	63.0%

Table 1: Overview of Pivotal Study Results

At week 12, 87.2% of the INNOVO® group were in the dry or mild categories of stress incontinence severity, compared to 54.5% at baseline (all of whom had been in the mild category at baseline). The adverse events (AEs) experienced by subjects were predominantly mild or moderate. No serious device-related adverse events occurred. The most common device-related adverse event with the INNOVO® was device discomfort (9.0%) followed by medical device pain (4.5%) and skin irritation (3.4%). Most AEs resolved by modifying the stimulation intensity. No device related urinary tract infections or vaginal infections in the INNOVO® group were reported.

Conclusion: INNOVO® external electrical stimulation device demonstrated clinically meaningful improvements in a range of subjective and objective measures of stress urinary incontinence. The INNOVO® device was well tolerated.

References:

[1] Sand PK, Richardson DA, Staskin DR, Swift SE, Appell RA, Whitmore KE, Ostergard DR. Pelvic floor electrical stimulation in the treatment of genuine stress incontinence: a multicenter, placebo-controlled trial. Am J Obstet Gynecol. 1995 Jul;173(1):72-9.

WARRANTY

Dear Customer,

Thank you for purchasing your Innovo® product ("Product"). To be covered by this Limited Warranty, the Product must have been purchased directly from an authorised retailer in the United States. This Limited Warranty extends only to the original purchaser and cannot be transferred to anyone. Please read this carefully before using the Product.

Our Limited Warranty to you

As the legal manufacturer, Atlantic Therapeutics Limited (Atlantic Therapeutics Ltd.) warrants to you, the original purchaser, that under normal use and service, the Product will be free from defects in material and workmanship for a period of two (2) years from the date of purchase ("Warranty Period"). If you discover a defect covered by this Limited Warranty during the Warranty Period, and return the Product to us, we will repair or replace the Product, at our option, using new or remanufactured components. There is no warranty or condition of any kind as to any defects discovered after the Warranty Period expires.

Certain Product failures not covered

This Limited Warranty covers defects in manufacture that arise from correct use of the Product per the product manual and instructions. It does not cover damage caused by abuse, misuse, accidents, modification or repair by any person other than Atlantic Therapeutics Ltd, moisture, extreme heat or cold, corrosive environments, high voltage surges from external sources or an unsuitable power supply. This Limited Warranty does not cover normal wear and tear on consumable items (e.g. gel pads, Shorts, etc.). This Limited Warranty does not apply to any battery provided with the Product, or if the Product has an altered or defaced serial number.

How to obtain service under this Limited Warranty

To obtain service under this Limited Warranty, please contact Atlantic Therapeutics Ltd's Customer CareLine at 1-833-891-6369 to obtain return instructions. You must submit proof of purchase in the form of a bona fide, dated bill of sale, receipt, or invoice (or copy thereof) evidencing that your request for service is made within the Warranty Period. Except as may otherwise be required by legislation in your State, Atlantic Therapeutics Ltd will not be responsible for the cost of returning the Product or for any damage to or loss of the Product while it is in transit to us. We recommend that you use a carrier that will track your shipment and ensure the Product's safe delivery. Atlantic Therapeutics Ltd cannot accept any returns sent C.O.D.

THIS LIMITED WARRANTY IS THE ONLY WARRANTY FOR THE PRODUCT, AND THERE ARE NO OTHER EXPRESS WARRANTIES, ORAL OR WRITTEN, PROVIDED BY Atlantic Therapeutics Ltd. ANY IMPLIED WARRANTIES APPLICABLE TO THE PRODUCT, INCLUDING ANY WARRANTIES OF MERCHANTABILITY OF FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE WARRANTY PERIOD SET FORTH ABOVE.

Limitation of Remedies

YOUR EXCLUSIVE REMEDY, IN LIEU OF ALL INDIRECT, INCIDENTAL, CONSEQUENTIAL, EXEMPLARY, SPECIAL OR OTHER DAMAGES OF ANY KIND, INCLUDING FOR NEGLIGENCE, IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE PRODUCT, IF DEFECTIVE, UNDER THE TERMS AND CONDITIONS OF THIS LIMITED WARRANTY.

Atlantic Therapeutics Ltd will bear no other expenses.

Atlantic Therapeutics Ltd reserves the right to refuse to repair or replace the Product if it is returned in an unsanitary condition. Atlantic Therapeutics Ltd reserves the right to substitute a product of equal or better quality if an identical product is not available at the time of replacement. This Limited Warranty is non-transferable.

Applicability of state law

Some states do not allow limitations on the duration of implied warranties, or the exclusion or limitation of incidental or consequential damages, and so the related limitations and exclusions stated above may not apply to you.

Governing law

This Limited Warranty is valid in the United States only and is not applicable in other countries. If you acquired the Product in the United States, the laws of the State of New York will apply to this Limited Warranty.

THIS LIMITED WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.

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Copyright© 2019
All rights reserved. INNOVO® is a trade
Atlantic Therapeutics Ltd no se hará cargo de otros gastos.
mark of Atlantic Therapeutics Ltd.