


ORDER FROM YOUR DENTAL SUPPLY COMPANY

Waste.

WASTED TIME IN PROCEDURE
WASTED MATERIAL
WASTED TIME CLEANING UP



2.2 CC

Not.

ROTATING TIP
REACH WITH EASE

LOCKING HEAD,
TWIST FOR
EXTENDED STORAGE

ENOUGH FOR
THREE CROWNS

NO MESSY CLEANUPS
SAVE ON
EVERY IMPRESSION

EASY TO USE
SUPERB ERGONOMICS



Stop trying to maneuver and control bulky mixing guns.

Stop having to assemble or change intraoral tips, mixing tips, cartridges, and dispensing guns for every small impression.

Stop wasting expensive material left over inside large mixing tips.

Stop cross-contaminating entire impression material cartridges.

Introducing the single use, disposable Mojo™ Syringe

The material wasted inside a standard mixing tip (up to 2.2 cc) is enough to do three additional crowns with the Mojo Syringe.

The small, yet sensibly angled tip of our syringe gives you easy access even to the distals of second molars. It allows the material to be applied right at the sulcus circumferentially, thus minimizing air bubbles and retakes.

ORDER FROM YOUR DENTAL SUPPLY COMPANY

THE MOJO™ SYRINGE

Single use, disposable

The material wasted inside a standard mixing tip (2.2 cc) is enough to do three additional crowns with our syringe.

The small yet sensibly angled tip of our syringe gives you easy access even to the distals of second molars and allows the material to be applied right at the sulcus circumferentially, thus minimizing air bubbles and retakes. Avoid messy cleanups and save money.

Simple 3 position, fill, store, activate.

Less waste

18g standard size intraoral tip

Superb ergonomics for small hands

For use with 1:1 cartridges only

Tip rotates with no-slip feature

Avoid messy cleanups

Save with every impression

Pre-loadable and disposable



Mojo 100 pk
P/N Mojo101

Technique:

NOTE: DO NOT REMOVE OR TWIST TIP PRIOR TO FILLING

1



The "Load" position of the mixing tip has the tab in a notch as indicated by the arrow. The syringe comes assembled in this position.

2



The supply cartridge, whether 25, 50 or 75ml, is firmly butted into the barrels so that there is a tight fit at the base of the syringe.

3



The VPS or core build-up materials are now injected to a level just short of the top of the barrels. A slow positive pressure allows materials at differing viscosities to self-equalize in the barrels.

4



Allowing just enough empty space at the top to accommodate the plungers, the twin plungers are firmly seated into the rear of the syringe. Be careful not to push material into the mixing chamber.

5



This photo shows the VPS material just at the mixing tip without entering it. The tab has been rotated 90° to the position above the "Λ" (Chevron). Material can be stored in this position, the plunger will not extrude the material if it is inadvertently pushed.

6



To "arm" or dispense the material, rotate the top until you hear an audible "click." The tip will now be rotated 90° from the "store" position to the "Λ" (Delta) position.

7



The material now flows freely down the mixing elements assuring a well incorporated catalyst into base of the light body or monophas VPS. Core build-up also works well in this syringe.

8



The unique design allows even the smallest hands to easily push the plunger. Its ergonomic, compact design assures comfort and accuracy.

