NADAL® Cardiac Combo Test (test cassette)

REF 282003

DE Gebrauchsanweisung	2	PL Sposób użycia	20
EN Instruction for use	6	PT Instruções de Utilização	23
FR Instructions d'utilisation	9	FI Käyttöohje	26
ES Instrucciones de uso	13	Symbols	31
IT Istruzioni per l'uso	16	Our Teams	32

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*95% Konfidenzintervall

14. Referenzen

- Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. Ann Clin Lab Sci, 26:301-12, 1996.
- Kagen LJ. Myoglobin methods and diagnostic uses. CRC Crit. Rev. Clin. Lab. Sci., 2:273, 1978.
- Chapelle JP. et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. Eur. Heart Journal, 3:122, 1982.

4. Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunoassay*, 17:24-9, 1994.
5. Lee TH, Goldman L. Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
6. Kallner A, Sylven C, Brodin. U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
7. Adams, et al. Biochemical markers of myocardial injury, *Immunoassay Circulation* 88: 750-763, 1993.
8. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J.Biol.Chem.* 266:966, 1991.

Rev. 4, 2020-06-12 OM

1. Intended Use

The NADAL® Cardiac Combo Test is a rapid visual immunoassay for the qualitative presumptive detection of human Myoglobin, CK-MB and cardiac Troponin I in whole blood, serum or plasma. The test is intended for use as an aid in the diagnosis of myocardial infarction (MI).

2. Introduction and Clinical Significance

Myoglobin (MYO), Creatine Kinase MB (CK-MB) and cardiac Troponin I (cTnI/TRO) are proteins released into the bloodstream after cardiac injury. Myoglobin is a hemeprotein normally found in skeletal and cardiac muscle with a molecular weight of 17.8 kDa. It constitutes about 2 percent of total muscle protein and is responsible for transporting oxygen within muscle cells. If muscle cells are damaged, Myoglobin is released into blood rapidly due to its relatively small size. The level of Myoglobin increases measurably above baseline within 2-4 hours post-infarct, peaking at 9-12 hours and returning to baseline within 24-36 hours. CK-MB is an enzyme with a molecular weight of 87.0 kDa which is also present in the cardiac muscle. Creatine Kinase is a dimeric molecule formed from two subunits designated as "M" and "B", which combine to form three different isoenzymes CK-MM, CK-BB and CK-MB. CK-MB is the isoenzyme of Creatine Kinase mostly involved in the metabolism of cardiac muscle tissue. The release of CK-MB into blood following an MI can be detected within 3-8 hours after the onset of symptoms. It peaks within 9 to 30 hours and returns to baseline levels within 48 to 72 hours. Cardiac Troponin I is a protein with a molecular weight of 22.5 kDa found in cardiac muscle. Troponin I is part of three subunit complex comprised of Troponin T and Troponin C. Along with tropomyosin, this structural complex forms the main component that regulates the calcium sensitive ATPase activity of actomyosin in striated skeletal and cardiac muscle. After cardiac injury occurs, Troponin I is released into blood 4-6 hours after the onset of pain. The release pattern of Troponin I is similar to that of CK-MB, however while CK-MB levels return to normal after 72 hours, Troponin I remains elevated for 6-10 days, thus providing a longer window of detection for cardiac injury.

3. Test Principle

The NADAL® Cardiac Combo Test enables the detection of Myoglobin, CK-MB and Troponin I through visual interpretation of colour development on the test membrane. Antibodies to Myoglobin, to CK-MB and to cTnI are immobilised on the respective test line regions of the membrane. During testing, the specimen reacts with antibodies to Myoglobin, to CK-MB and to cTnI conjugated to coloured particles and precoated on the sample pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with components on the membrane. If there is a sufficient amount of certain markers in the specimen, a coloured line will develop in the corresponding test line region of the membrane. The presence of the coloured test line indicates a positive result for the corresponding marker, while its absence indicates a negative result. The appearance of a coloured line in the control line region serves as a procedural control indicating that the proper volume of specimen has been added and that membrane wicking has occurred.

4. Reagents and Materials Supplied

- 5 NADAL® Cardiac Combo test cassettes
- 1 buffer
- 5 disposable pipettes
- 1 package insert

5. Additional Materials Required

- Specimen collection container
- Centrifuge (only for serum and plasma specimens)
- Timer

6. Storage & Stability

- The test kit should be stored at 2-30°C until the expiry date printed on the sealed pouch.
- The test cassette should remain in the sealed pouch or in a closed canister until use.
- **Do not freeze.**
- Care should be taken to protect the components of the test kit from contamination. Do not use the test if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Do not use the test after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged. Do not reuse tests and disposable pipettes.
- Test kits contain products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing usual safety precautions (e.g. do not ingest or inhale).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire package insert carefully prior to testing.
- Do not eat, drink or smoke in the area where the specimens and test kits are handled.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Do not interchange or mix reagents from different lots.
- Humidity and temperature can adversely affect results.
- The test components (e.g. antibodies/chemicals) present no danger if the test is carried out according to the instructions.
- Used test materials should be discarded in accordance with local regulations.

8. Specimen Collection and Preparation

- The NADAL® Cardiac Combo Test is intended for use with human whole blood, serum or plasma specimens only.
- Only clear, non-hemolysed specimens are recommended for use with this test. Serum or plasma should be separated as soon as possible to avoid hemolysis.

- Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Containers containing anticoagulants such as EDTA, citrate, or heparin should be used for whole blood storage.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and well mixed prior to testing. Avoid repeated freezing and thawing of specimens.
- If specimens are to be shipped, pack them in compliance with all applicable regulations for transportation of etiological agents.
- Icteric, lipemic, hemolysed, heat treated and contaminated specimens may cause erroneous results.
- There is a slight possibility that some whole blood specimens with very high viscosity, or those which have been stored for more than 2 days, may not run properly on the test cassette. Repeat the test with a serum or plasma specimen collected from the same patient using a new test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

As the test begins to work, colour will migrate along the membrane.

3. Wait for the coloured line(s) to appear. The result should be read at 10 minutes. Do not interpret the result after 20 minutes.


10 min

10. Result Interpretation

Positive result

One line appears in the control line region (C) and one, two or three line(s) appear(s) in the test line regions (T): (MYO) for Myoglobin, (CKMB) for CK-MB and (cTnI) for Troponin I.

Negative result

Only one coloured line appears in the control line region (C). No coloured lines appear in the test line regions (MYO), (CKMB) and (cTnI).

Invalid result

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

NOTE:

- The intensity of colour in the test line regions (MYO), (CKMB) and (cTnI) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of colour in the test line regions should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of analytes in the specimen.
- Insufficient specimen volume, incorrect test procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

- The internal procedural control is included in the test cassette. A coloured line appearing in the control line region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

9. Test Procedure

I. Preparation

- **Bring test cassettes, specimens, buffer and/or controls to room temperature (15-30°C) before use.**
- Remove the test cassette from the sealed pouch and use it as soon as possible. For the best results, the assay should be performed within one hour.

II. Test Procedure

1. Place the test cassette on a clean and level surface.

a) Serum or plasma specimens


Hold the disposable pipette vertically and transfer 2 drops of serum or plasma (approximately 50 µL) to the specimen well (S) of the test cassette. Start the timer.

OR

b) Whole blood specimens

Transfer 3 drops of whole blood specimen (approximately 75 µL) to the specimen well (S) of the test cassette with the provided disposable pipette, then add 1 drop of buffer and start the timer.

OR

c) Fingerstick whole blood specimens

Allow 3 hanging drops of fingerstick whole blood specimen to fall into the centre of the specimen well (S) of the test cassette and then add 1 drop of buffer and start the timer.

- External controls are not supplied with the test kit. In accordance with a *Good Laboratory Practice (GLP)* it is recommended that positive and negative controls be tested to confirm the test procedure and to verify proper test performance.

12. Limitations

- The NADAL® Cardiac Combo Test is for professional *in-vitro* diagnostic use only and should be used for detection of Myoglobin, CK-MB and Troponin I. No meaning should be inferred from the color intensity or width of any apparent lines.
- The NADAL® Cardiac Combo Test only indicates the presence of Myoglobin, CK-MB and Troponin I in the specimen and should not be used as the sole criterion for the diagnosis of myocardial infarction.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. The test cannot detect less than 50 ng/mL Myoglobin, 5 ng/mL CK-MB or 0.5 ng/mL Troponin I in specimens. Thus, a negative result does not at anytime rule out the existence of these cardiac markers in blood because they may be below the minimum detection level of the test.
- Like with all diagnostic tests, a confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
- Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor (RF) may affect expected results. Even if the test results are positive, further clinical evaluation should be considered with other clinical information available to the physician.

13. Performance Characteristics

The NADAL® Cardiac Combo Test (Myoglobin) vs. EIA

		NADAL® Cardiac Combo Test (Myoglobin)		
		+	-	Total
EIA	+	60	0	60
	-	9	374	383
Total		69	374	443

Relative Sensitivity: 100% (94.0% - 100.0%)*

Relative Specificity: 97.7% (95.6% - 98.9%)*

Overall Agreement: 98.0% (99.7% - 99.9%)*

*95% Confidence Interval

The NADAL® Cardiac Combo Test (CK-MB) vs. EIA

		NADAL® Cardiac Combo Test (CK-MB)		
		+	-	Total
EIA	+	54	0	54
	-	1	422	423
Total		55	422	477

Relative Sensitivity: 100% (93.4% - 100.0%)*

Relative Specificity: 99.8% (98.7% - 99.9%)*

Overall Agreement: 99.8% (98.8% - 99.9%)*

*95% Confidence Interval

The NADAL® Cardiac Combo Test (Troponin I) vs. EIA

		NADAL® Cardiac Combo Test (Troponin I)		
		+	-	Total
EIA	+	225	3	228
	-	8	505	513
Total		233	508	741

Relative Sensitivity: 98.7% (96.2% - 99.7%)*

Relative Specificity: 98.4% (97.0% - 99.3%)*

Overall Agreement: 98.5% (97.4% - 99.3%)*

*95% Confidence Interval

14. References

- Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
- Kagen LJ. Myoglobin methods and diagnostic uses. *CRC Crit. Rev. Clin. Lab. Sci.*, 2:273, 1978.
- Chapelle JP, et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
- Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunoassay*, 17:24-9, 1994.
- Lee TH, Goldman L. Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
- Kallner A, Sylven C, Brodin U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
- Adams, et al. Biochemical markers of myocardial injury, *Immunoassay Circulation* 88: 750-763, 1993.
- Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem*. 266:966, 1991.

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*95% d'intervalle de confiance

14. Bibliographie

- Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
- Kagen LJ. Myoglobin methods and diagnostic uses. *CRC Crit. Rev. Clin. Lab. Sci.*, 2:273, 1978.
- Chapelle JP. et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
- Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunology*, 17:24-9, 1994.
- Lee TH, Goldman L. Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.

6. Kallner A, Sylven C, Brodin. U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. Scand J Clin Lab Invest, 49:633-9, 1989.
7. Adams, et al. Biochemical markers of myocardial injury, Immunoassay Circulation 88: 750-763, 1993.
8. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. J.Biol.Chem. 266:966, 1991.

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Tasa de coincidencia: 98,5% (97,4% - 99,3%)*

*95% Intervalo de confianza

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1. Wong SS: Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
2. Kagen U: Myoglobin methods and diagnostic uses. *CRS Crit Rev Clin Lab. Sci*, 2:273, 1978.
3. Chappelle JP et al.: Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
4. Apple FS, Preese LM: Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunoassay*, 17:24-9, 1994.
5. Lee TH, Goldman L: Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
6. Kallner A, Sylven C, Brodin U, et al.: Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
7. Adams, et al.: Biochemical markers of myocardial injury, *Immunoassay Circulation* 88: 750-763, 1993.
8. Mehegan JP, Tobacman LS: Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem*. 266:966, 1991.

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14. Bibliografia

- Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
- Kagen LJ. Myoglobin methods and diagnostic uses. *CRC Crit. Rev. Clin. Lab. Sci.*, 2:273, 1978.
- Chapelle JP, et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
- Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunoassay*, 17:24-9, 1994.

5. Lee TH, Goldman L. Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
6. Kallner A, Sylven C, Brodin U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
7. Adams, et al. Biochemical markers of myocardial injury, immunoassay *Circulation* 88: 750-763, 1993.
8. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J.Biol.Chem.* 266:966, 1991.

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14. Bibliografia

1. Wong SS: Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
2. Kagen LJ: Myoglobin methods and diagnostic uses. *CRC Crit Rev Clin Lab. Sci*, 2:273, 1978.
3. Chappelle JP et al.: Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
4. Apple FS, Preese LM: Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunoassay*, 17:24-9, 1994.
5. Lee TH, Goldman L: Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
6. Kallner A, Sylven C, Brodin U, et al.: Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
7. Adams, et al.: Biochemical markers of myocardial injury, *Immunoassay Circulation* 88: 750-763, 1993.
8. Mehegan JP, Tobacman LS: Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem*. 266:966, 1991.

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*95% Intervalo de confiança

14. Referências

- Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
- Kagen LJ. Myoglobin methods and diagnostic uses. *CRC Crit. Rev. Clin. Lab. Sci.*, 2:273, 1978.
- Chapelle JP. et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
- Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunossay*, 17:24-9, 1994.
- Lee TH, Goldman L. Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
- Kallner A, Sylven C, Brodin U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
- Adams, et al. Biochemical markers of myocardial injury, *Immunoassay Circulation* 88: 750-763, 1993.
- Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem*. 266:966, 1991.

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Suhteellinen spesifisyys: 98,4% (97,0% - 99,3%)*

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
14. Lähteet

1. Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci*, 26:301-12, 1996.
2. Kagen L. Myoglobin methods and diagnostic uses. *CRC Crit. Rev. Clin. Lab. Sci.*, 2:273, 1978.
3. Chapelle JP. et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. *Eur. Heart Journal*, 3:122, 1982.
4. Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. *J Clin Immunoassay*, 17:24-9, 1994.
5. Lee TH, Goldman L. Serum enzyme assays in the diagnosis of acute myocardial infarction. *Ann Intern Med*, 105:221-233, 1986.
6. Kallner A, Sylven C, Brodin. U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. *Scand J Clin Lab Invest*, 49:633-9, 1989.
7. Adams, et al. Biochemical markers of myocardial injury, *Immunoassay Circulation* 88: 750-763, 1993.
8. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem*. 266:966, 1991.

Rev. 4, 2020-06-12 OL

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Portugués	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro -diagnostikkaan tarkoitettu lääkinällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik	In-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur-begränsning	Temperatuurlimiet	Temperatur-begrænsning	Temperatur begrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merkning
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tilverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel –UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16

Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel. Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607

Free Tel: +45 80 88 87 53

Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40

Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH