


Quick Start Guide

1 Insert batteries.



2 Download Meraw Health App and set up Bluetooth connection.

1. The monitor needs to be at OFF status before connection;
2. Open Meraw Health App, make sure you have authorized Bluetooth and location permission for the App (**important!**);
3. Long press Bluetooth button on the monitor until the "Bluetooth" symbol blinks in the screen(don't press "start/stop" button);
4. Follow the instruction in the App, "Add device" - "Meraw Aspen" - "Start connection";
5. The pairing will complete automatically.



App download



Video instruction

3 Prepare to measure.



1

Put on blood pressure monitor and secure the cuff.



2

Press start/stop, then adjust the vertical position.




3

When wrist is in ideal position, the symbol ✓ flashes.

4 Sit correctly and take measurement.

- Keep palm relaxed during measurement
- Sit in a chair with your legs uncrossed and feet flat on the floor.
- Sit with your back and arm supported.



◆ Contact our customer support team ◆



support@merawlabs.com



www.merawlabs.com