

Cold Gel Tennis Elbow Brace

SKU: ARB199

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Loop strap through plastic ring and slide up the forearm.

Position gel pad on top of forearm.

Step 2. Tighten strap by wrapping snugly. Do not overtighten.

Step 3. For cold therapy, place gel pad in freezer for one hour then reinsert into brace.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

