

Foam Arm Sling Shoulder Immobilizer

SKU: ARS263

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Unfasten all straps and place arm sling in front of body. Gently pull the sling over your arm and elbow. It should fit snugly around the elbow. Your hand should be at the very end of the sling.

Step 2. Take the left strap and put over the left shoulder, across to the right side of the body and Velcro onto arm sling. Take right strap and put over right shoulder, across to the left side of the body and Velcro onto arm sling, forming a criss cross in the back. Tighten the straps to keep your hand and forearm elevated above the level of your elbow.

Step 3. Optionally, add the swathe by wrapping around upper arm and chest and securing with Velcro for additional immobilization.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

