

Low Profile Lumbar Back Support Brace

SKU: BKB266

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Fully open and extend the support. The panel with the pull tab should be on the left.

Step 2. Center the two ends of the pulley system on the lower back. Wrap the sides around the abdomen with the right side overlapping the left.

Step 3. Pull the handle to the desired position. When the desired amount of tension is reached, secure the pull tab on the Velcro strap. *Tip: Pull in stomach muscles before tightening.*



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

