

90° Comfort Padded Night Splint

SKU: FTB-CL63-105

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code




Instructions For Use

Step 1. Apply the splint with the knee flexed to maximize the angle of dorsiflexion.

Step 2. Attach instep strap first to hold the heel securely in the splint. Adjust the length of the straps and close. Then extend knee to check the tightness of the instep strap.

Step 3. Close dorsiflexion assist straps and adjust them to your tolerance level.

Step 4. Place toe wedge under toe for additional stretching in necessary.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Remove splint and take out grey liner
- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse thoroughly and Air Dry.
- Do NOT dry clean, iron, or bleach. 