

ROM Knee Brace with Hinges PDAC

SKU: KNB129

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. To set the hinge to the desired flexion/ extension settings, remove the hinge from the brace and open. Drop the metal stops into the desired setting.

Step 2. Unfasten all Velcro straps and apply seated with the knee bent. Position the patella pad on the knee cap and align the hinges with the center of the knee joint.

Step 3. Fasten the calf strap first (this strap prevents slipping) and then top strap. Fasten the popliteal strap. If slipping occurs, check that the calf strap is above the largest part of the calf muscle and tightened securely.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Remove hinges before cleaning
- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse thoroughly and Air Dry.
- Do NOT dry clean, iron, or bleach.

