

Internal Deluxe Plus Size ROM Knee Brace

SKU: KNB163

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Remove the hinges from the neoprene brace. Open the cover on the ROM hinges. Set the desired degree of flexion or extension by placing the pins in the appropriate degree position. Make sure settings are the same for both hinges and reinsert into the brace.

Step 2. Unfasten Velcro straps and tabs and open the brace. Place the brace around the knee so that the hinges line up with the center of the knee cap on both sides of the leg and secure the bottom Velcro tab.

Step 3. Wrap the top around the leg and secure the Velcro tab making sure the hinges still align with the center of the knee and the brace fits snugly.

Step 4. Feed the bottom strap through the bottom D Ring and tighten to secure. Feed the top strap through the top D Ring and secure. Readjust as needed until proper fit is achieved.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Remove hinges before cleaning.
- Hand Wash with a 50% Alcohol/Water Mixture
- Air Dry only.
- Do NOT dry clean, iron, or bleach.  

